


Presentation of the system



The leader in immediate restoration

We are

we are
One!

bredent group : **One** manufacturer in the field of implantology and prosthetics

bredent medical : **One** of the leaders in immediate restorations

HELBO : **One** bacterial infection control

SKY fast & fixed : **One** session is all it takes

BioHPP SKY elegance : **One** abutment for temporary and final restoration

Dentists, dental technicians, bredent group : **One** team

All patients : **Everyone** is satisfied

System overview

- Seven implant lines – one instrument kit
- The right implant for every bone situation
- SKY implant surfaces
- SKY implant design for maximum primary stability
- SKY implants overview
- SKY implant positioning in relation to the bones
- SKY surgical protocol
- Guided implantology
- SKY surgical protocol – guided

Prosthetics

- Implant and abutment platforms
- Implant connections
- Prosthetics overview
- Classic implantology
- CAD/CAM-manufactured restorations
- Immediate and late restoration
- Prosthesis fixation

miniSKY

The implant for narrow alveolar ridges

copaSKY

The implant for broad alveolar ridges with low height

whiteSKY

The esthetic zirconium implant

Accessories and instruments for the SKY implant lines

- Drill
- Surgical tools
- Prosthetic instruments
- Screws
- Torques

HELBO® therapy

visio.lign veneer system


- Chipping repair kit
- Full range bonding kit for all prosthetic materials

Order information

Fax order forms

Seven implant lines – one instrument kit





SKY Surgical Kit
REF SKYXOT21



A study by the University of Belgrade showed that, when using the SKY drill, only a small amount of heat was generated in the bone.


Source: Marković et al: Heat generation during implant placement in low-density bone: effect of surgical technique, insertion torque and implant macro design. Clin Oral Implants Res. 2013 Jul;24(7):798-805. DOI: 10.1111/j.1600-0501.2012.02460.x. Epub 2012 Apr 2.

System overview




The right implant for every bone situation

Narrow implants
 $\text{Ø} < 3 \text{ mm}$

Implants with reduced diameter
 $\text{Ø} 3 - 3.5 \text{ mm}$





Initial situation
Following tooth extraction and unloaded (edentulous) jaws, significant bone loss occurs even after a short period of time.



*mini*¹
SKY 2.8



*mini*²
SKY 2.8



*mini*²
SKY 3.2



narrow
SKY

Angled



blue
SKY

Short
Length < 8 mm




classic
SKY (tissue level)

copa
SKY


System overview


SKY osseo connect surface (ocs)[®]



The structure of the connective tissue consists of horizontal fibres that attach themselves to the natural tooth and therefore prevents plaque attachment. The horizontal micro-grooves on the SKY implants also support attachment of the soft tissue so that a type of soft tissue cuff is created to protect the implant.

blue SKY Iso-crestal implant position







Horizontal micro-grooves support attachment of the connective tissue

1.5 mm


0.5 mm



Etched transition structure provides the bones and soft tissue possibilities for adaptation



Blasted etched surface for the attachment of osteoblasts for rapid osseointegration



classic **SKY** Semi-transgingival implant position




Design of a soft-tissue collar

SKY implant design for maximum primary stability

Immediate restoration places particular requirements on an implant system. The design of the SKY Implants and the corresponding surgical protocol ensure high primary stability


in all bone qualities and therefore form a reliable basis for immediate restoration.



Scientifically proven

A study by Prof. Marković proves that the blueSKY implant has a very high level of primary stability in comparison to implants from competitors.

It also shows that rapid osseointegration takes place with modern implant surfaces and that no loss of stability occurs after a few weeks. Secure immediate restoration is therefore possible.




Changes to implant stability during the 12-week monitoring period.

- △-- Drilling + self-cutting blueSKY implant
- Drilling + non-self-cutting competitor implant

Source: Marković et al: Evaluation of primary stability of self-tapping and non-self-tapping dental implants. A 12-week clinical study, *Clinical Implant Dentistry and Related Research* 2013

Important:

The surgical protocol for the SKY system enables the primary stability of the SKY implants to be increased easily using condensation: See page: 18



System overview

SKY implants overview

narrowSKY 3.5 N



NP - Narrow platform

3.5 N		REF	3.5 N	
	Length 8 mm	nSKY3508		Including a cover screw
	Length 10 mm	nSKY3510		
	Length 12 mm	nSKY3512		
	Length 14 mm	nSKY3514		
	Length 16 mm	nSKY3516		

blueSKY



RP - Regular platform


4.0		REF	4.0 4.5	
	Length 8 mm	bSKY4008		Including a cover screw
	Length 10 mm	bSKY4010		
	Length 12 mm	bSKY4012		
	Length 14 mm	bSKY4014		
	Length 16 mm	bSKY4016		

SKY classic



RP - Regular platform

4.0		REF	4.0 4.5	
	Length 8 mm	kSKY4008		Including a cover screw
	Length 10 mm	kSKY4010		
	Length 12 mm	kSKY4012		
	Length 14 mm	kSKY4014		
	Length 16 mm	kSKY4016		




4.5	REF
Length 8 mm	bSKY4508
Length 10 mm	bSKY4510
Length 12 mm	bSKY4512
Length 14 mm	bSKY4514

5.5	REF
Length 8 mm	bSKY5508
Length 10 mm	bSKY5510
Length 12 mm	bSKY5512




4.5	REF
Length 8 mm	kSKY4508
Length 10 mm	kSKY4510
Length 12 mm	kSKY4512
Length 14 mm	kSKY4514

SKY implant positioning in relation to the bones



Supracrestal

Isocrestal





Isocrestal

Isocrestal


Indication




Indication




Supracrestal and isocrestal implant position



Short implants



implant position



Augmented implant site

The SKY classic implant reduces the need for bone reduction in a narrow or irregular ridge.

The 8 mm implant can be used as a short implant (6.0 mm) by way of supracrestal positioning.

The SKY classic implant is perfectly suited to flapless implantation, since the long, machined neck easily allows for a semi-transgingival implant position.

The coronally structured blueSKY implant is perfectly suited to being positioned flush with the bone.

The macro-grooves ensure a high level of bone preservation.

The blueSKY implant therefore works very well with augmentations.

Surgical protocol tailored to bone quality

Bone quality from hard to soft

Hard bone



Atraumatic tapping

Prevention of bone overloading during surgery

Medium-hard to soft bone

Bone compression


Achieving primary stability




System overview

SKY Surgical protocol

- Optimised, bone quality-oriented set of instruments and surgical protocol for unsurpassed primary stability
- Drills with detachable drill stops
- Reduction of the number of drills for increased control and dependability during surgery



Detachable drill stops




The detachable drill stops are arranged so that they can be easily taken up with the drill and fastened with one hand thanks to the cavities in the OP-Tray insert. The ruler on the compartment for small parts helps to select the correct drill stop easily and quickly. The drilling depth is approx. 0.7 mm greater than the implant length.

classic
SKY






white
SKY
IMPLANT SYSTEM




With the 3.5 diameter blueSKY and SKY classic implants, the crestal drill is only sunk up to the laser marking.

In the following implants:
 • narrowSKY
 • blueSKY 4.0 to 5.5
 • SKY classic 4.0 to 4.5
 the crestal drill is completely inserted.



Crestal drill
300 rpm





System overview

SKY Surgical protocol

If it is determined during the pilot drilling or drilling with the twist drill that the bone is very soft, the primary stability can be improved by amending the clinical protocol. In these

cases, we recommend using the final drill anticlockwise as a condensation instrument:



The pilot drill and twist drill are used as described in the SKY surgical protocol.

The final drill is used anticlockwise slowly with cooling. This way, the available bone is compressed and no bone particles are lost.

The crestal drill is used in accordance with the surgical protocol.

Guided implantology

The data for the SKY implant systems are recorded in the following planning programs: (As at: 31.12.2016).
The scope of the recorded components is dependent upon the relevant options in the planning program.

3Diagnosys
3Shape Implant Studio™
Invivo
Carestream 3D Imaging
In2Guide™
coDiagnostiX™
Simplant®
i-Dixel
iRYS

3DIEMME®
3Shape
Anatmage
Carestream Dental
Cybermed / KAVO
Dental Wings
Materialise Dental
Morita
myRay

If you do not find your program in the list, please contact us, so that we can make our data available to the manufacturer of your program for updating.

NemoScan
NNT
Romexis®
Implant-3D
SICAT/Galileos Implant
SMOP
Accu Guide
Zfx Navigator


NEMOTEC
NewTom
PLANMECA
Schütz
SIRONA
Swissmeda
UniGuide Dental
Zfx GmbH

coDiagnostiX™

bredent is a distribution partner of Dental Wings for the coDiagnostiX 3D implant planning software.

coDiagnostiX™ –
Client Version
REF SplanX120

coDiagnostiX™ –
Producer Version
REF SplanX100




SKY Surgical protocol



Removing and screwing in the implant and the cover screw without changing instruments.


- blueSKY, narrowSKY and SKY classic are packed as double sterile
- The implant support is colour-coded and marked with the length details
- The implant can be removed immediately after opening the sterile films using the SKY TK long or short insertion instrument and screwed in
- The cover screw is also removed using the same instrument




1 Open the new double sterile packaging. Remove the implant with the insertion instrument for the ratchet or the contra-angle.



2 The conical Torx® allows to hold the implant safely.




3 The implant is inserted and can be screwed in immediately.



4 The cover screw can be removed using the same instrument.



5 The cover screw is held securely by the cone.



6 Screw the screw directly into the implant. The smooth cone ensures that the screw only needs to be slightly tightened and cannot become jammed. Recommended max. torque: 10 Ncm.

Guided implantology

The data for the SKY implant systems are recorded in the following planning programs: (As at: 31.12.2016).
The scope of the recorded components is dependent upon the relevant options in the planning program.

3Diagnosys
3Shape Implant Studio
Invivo
Carestream 3D Imaging
In2Guide
coDiagnostiX
SimPlant
i-Dixel
iRYS

3DIEMME
3Shape
Anatmage
Carestream Dental
Cybermed / KAVO
Dental Wings
Materialise Dental
Morita
myRay

If you do not find your program in the list, please contact us. so that we can make our data available to the manufacturer of your program for updating.

NemoScan
NNT
Romexis
Implant-3D
SICAT/Galileos Implant
SMOP
Accu Guide
Zfx Navigator


Nemotec
NewTom
Planmeca
Schütz
SIRONA
Swissmeda
UniGuide Dental
Zfx GmbH

coDiagnostiX

bredent is a distribution partner of Dental Wings for the coDiagnostiX 3D implant planning software.

coDiagnostiX –
Client Version
REF SplanX120

coDiagnostiX –
Producer Version
REF SplanX100







System overview of guided implantology

SKY surgical protocol – guided




Optimised, bone quality-oriented set of instruments and surgical protocol for unsurpassed primary stability





Crestal drill




3.5 N



3,5



4.0



4.5




SKYplanX
Crestal drill



REF SXBHI433



Final drill
Ø 4.3 mm



REF SXBHI409



Final drill
Ø 4.06 mm

Drilling depth - Crestal drill

Drill length	Implant length				
	L8	L10	L12	L14	L16
L20	K4	K3			
L24		L1	K4	K3	
L28				L1	K4


Insertion instruments

O-ring position

H (high)


M (medium)

L (low)




Implant length

8 mm 10 mm 12 mm 14 mm 16 mm



System overview of guided implantology

SKY surgical protocol – guided



Guiding sleeves



REF	SXBHM453	SXBHM555
Description	SKYplanX guiding sleeve	SKYplanX guiding sleeve
	4.53	5.55
piece	1	1
Colour	silver	silver
Internal Ø / mm	4.53	5.55
External Ø / mm	5.5	6.5
Length / mm	10	10
Can be cut to size: to 6 mm	✓	✓
Material	Titanium	Titanium
Implant Ø	3.5/4.0	4.5
narrowSKY	✓	✓
blueSKY / SKY classic	✓	✓

Internal sleeves for twist drill



REF	SXBHI239
Description	Internal sleeve for pilot drill
piece	1
Colour	green
Internal Ø / mm	2.39
External Ø / mm	4.51
Length / mm	10
Can be cut to size: to 6 mm	✓
Material	Titanium
Guiding sleeve	SXBHM453
Implant Ø	3.5 / 4.0 / 4.5
narrowSKY	✓
blueSKY / SKY classic	✓

Sleeves for drills for medium/soft bone



REF	SXBHI309	SXBHI359	SXBHI409
Description	Internal sleeve Final drill	Internal sleeve Final drill	Internal sleeve Final drill
	3.06	3.56	4.06
piece	1	1	1
Colour	yellow	pink	grey
Internal Ø / mm	3.09	3.59	4.09
External Ø / mm	4.51	4.51	4.52
Length / mm	10	10	10
Can be cut to size: to 6 mm	✓	✓	✓
Material	Titanium	Titanium	Titanium
Guiding sleeve	SXBHM453	SXBHM453	SXBHM555
Implant Ø	3.5	4.0	4.5
narrowSKY	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓

Sleeves for drills for hard bone



REF	SXBHI333	SXBHI383	SXBHI433
Description	Internal sleeve Final drill	Internal sleeve Final drill	Internal sleeve Final drill
	3.30	3.80	4.30
piece	1	1	1
Colour	gold	blue	brown
Internal Ø / mm	3.33	3.83	4.33
External Ø / mm	4.51	4.51	4.52
Length / mm	10	10	10
Can be cut to size: to 6 mm	✓	✓	✓
Material	Titanium	Titanium	Titanium
Guiding sleeve	SXBHM453	SXBHM453	SXBHM555
Implant Ø	3.5	4.0	4.5
narrowSKY	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓



SKYplanX selection box

REF SXBHSET1

- 6 x SXBHM453 Guiding sleeve
- 3 x SXBHM555 Guiding sleeve, large
- 3 x SXBHI239 Sleeve for pilot drill
- 3 x SXBHI309 Internal sleeve for final drill 3.06
- 3 x SXBHI333 Internal sleeve for final drill 3.30
- 3 x SXBHI359 Internal sleeve for final drill 3.56
- 3 x SXBHI383 Internal sleeve for final drill 3.80
- 3 x SXBHI409 Internal sleeve for final drill 4.06
- 3 x SXBHI433 Internal sleeve for final drill 4.30



SKYplanX
Short tool for drill sleeves
1 piece
REF SplanX45

SKY surgical protocol – guided



SKYplanX OP tray
REF SplanX91

SKYplanX pilot drill



Ø	Length mm	REF
2,35 mm	20	DP20D235
2,35 mm	24	DP24D235
2,35 mm	28	DP28D235

SKYplanX final drill for hard bone

3.5 N



Ø	Length mm	REF
3.30 mm	3.5 3.5 N 20	D120D330
3.80 mm	4.0 20	D120D380
4.30 mm	4.5 20	D120D430
3.30 mm	3.5 24	D124D330
3.80 mm	4.0 24	D124D380
4.30 mm	4.5 24	D124D430
3.30 mm	3.5 28	D128D330
3.80 mm	4.0 28	D128D380
4.30 mm	4.5 28	D128D430

SKYplanX final drill for medium and soft bone



Ø	Length mm	REF
3.06 mm	3.5 3.5 N 20	D320D306
3.56 mm	4.0 20	D320D356
4.06 mm	4.5 20	D320D406
3.06 mm	3.5 24	D324D306
3.56 mm	4.0 24	D324D356
4.06 mm	4.5 24	D324D406
3.06 mm	3.5 28	D328D306
3.56 mm	4.0 28	D328D356
4.06 mm	4.5 28	D328D406

Crestal drill



Ø	Length mm	REF
3.6 mm	3.5 N 20	DCN35L20
4.1 mm	3.5 20	DCR35L20
4.6 mm	4.0 20	DCR40L20
5.2 mm	4.5 20	DCR45L20
3.6 mm	3.5 N 28	DCN35L28
4.1 mm	3.5 28	DCR35L28
4.6 mm	4.0 28	DCR40L28
5.2 mm	4.5 28	DCR45L28



SKY mounting set for guided implantology
REF SKYSMSET

The insertion instrument is screwed in for a secure and precise connection with the implant.
The flexible O ring yields a little when it reaches the drill sleeve shoulder and therefore counteracts damage to the implant bed in the event that the implant is screwed in further.

Contents:

- 1 x SKY moulder, extra short
- 1 x SKYplanX insertion instrument for guiding sleeve 4.53
- 1 x SKYplanX insertion instrument for guiding sleeve 5.55

SKY prosthetics

- Implant and abutment platforms
- Implant connections
- Prosthetics overview
- Classic implantology
 - Implant analogue
 - SKY esthetic line
 - SKY esthetic gingiva former
 - SKY impression abutments
 - SKY esthetic line temporary abutments
 - SKY esthetic line titanium abutments
 - SKY standard line titanium abutments
 - SKY castable abutment
- CAD/CAM-manufactured restorations
 - SKY uni.fit scan abutments
 - SKY uni.fit titanium base
 - SKY prefab titanium
 - SKY uni.fit titanium base for CEREC
- Immediate and late restoration
 - SKY elegance overview
 - SKY elegance abutments
 - SKY elegance prefabs
 - SKY elegance titanium base
 - SKY fast & fixed – overview
 - SKY fast & fixed – components
 - SKY fast & fixed – accessories
- Prosthesis fixation
 - SKY TiSi-snap
 - SKY Locator®

Number of treatment steps – implant treatments

	Classic implantology	Immediate restoration		Implant insertion
		Temporary + permanent	Permanent	
Implant planning	X	X	X	X
Implant surgery	X	X	X	X
Implant control	X	X	X	X
Implant opening	X			X
Modelling	X	X		
Esthetic try-in	(X)			
Permanent restoration	X	X		


Surgical procedure
 Implant is restored

Implant restoration is now an established treatment for tooth loss, however the conventional bridge on ground teeth is still the standard restoration.

The reason for this is:

- The higher cost of an implant restoration
- The longer treatment duration with several visits to the practice in the case of classic implantology
- Anxiety towards surgical procedures

By choosing the correct implant, it is possible to make optimal use of the local bone and therefore avoid surgical procedures. The treatment time can be shortened by using immediate restoration, so that the total treatment can be offered at a cheaper price, without having to reduce the treatment fee.



The deciding criteria as to whether immediate restoration is possible is the primary stability of the implants. If the primary stability is over 30 Ncm, immediate restoration should be carried out, as this stimulates the bone and promotes healing.



If the primary stability is under 30 Ncm, we recommend modelling the implant position during the operation and inserting the permanent restoration into the opening. This way, the treatment time is significantly shortened ¹⁾.

¹⁾ Literature:
Deutsche Zahnärztliche Zeitschrift [The German Dental Journal] 2014; 69 (6)
F. Beuer et al.: The Munich Implant Concept (MIC): a combination of intraoral scanning device and digital fabrication.


Implant and abutment platforms

bredent medical SKY Implant system


The SKY System is designed so that the number of prefabricated components is kept to a minimum. For all indications which cannot be restored using these components, there are various custom solutions available.



blueSKY and SKY classic
have a platform switch




Implant connections




All SKY Implants have a Torx® connection.

If screws and screw joints are involved. The Torx® is the gold standard in mechanical engineering and the automotive industry and in implant dentistry as well.



Torx®: has six large force transfer surfaces


- significantly higher torques for the same force applied
- easier insertion of the implant
- no damage to the internal geometry at a high torque either



Connection with self-retaining cone

- No definitive vertical abutment height
- Height difference between laboratory and clinic
- Passive fit of bridge constructions very difficult


- Conical connection
- Flat connection



Flat connections (SKY System)

- **Defined** abutment height
- Passive-fit of bridge and bar restorations easier to achieve

Owing to the normal position and the R variants, there are 12 positions for the orientation of the angled abutments. Therefore, the abutment can be oriented to the best position in the laboratory after the operation.



SKY prosthetics

Prosthetics overview

For narrowSKY approved prosthetic components
rose gold anodised

Prosthetic components not suitable or
not approved for narrowSKY

Immediate restoration

BioHPP SKY elegance
Immediate restoration




BioHPP SKY elegance
titanium base



BioHPP SKY elegance
prefab set



SKY uni.cone



SKY fast & fixed



Individual solutions for CAD/CAM – conventional

SKY uni.fit CAD/CAM solutions / Custom abutments



SKY prefab titanium



SKY castable abutment




Implant analog




SKY temp




SKY impression abutments



SKY esthetic gingiva former



SKY esthetic abutments




SKY titanium abutments



SKY prosthesis fixation

SKY Locator



SKY Locator anglet



SKY Locator




TiSi.snap



retention.sil




Implant analog



For pre-fabricated and custom solutions.



Just a single implant analog for all prosthetic restorations on implant level, independent of the implant platform. The implant analog is made of titanium to enable both the laboratory and the clinic/surgeon to use the same material.



REF	SKY-IA40
Description	SKY implant analogue
piece	1
Height/mm	14
Shoulder Ø / mm	4
Material	Ti*
Laboratory screw	incl.
SKY prosthetic key	✓
Torque/Ncm	25
Platform	regular
narrowSKY	✓
blueSKY / SKY classic	✓

Ti*=
Grade 4 KV titanium


SKY esthetic line



SKY esthetic line:


- Three diameters of the abutment shoulder:
 - S: 4.5 mm
 - M: 5.5 mm
 - L: 7.0 mm
- Narrow platform
 - Suitable for narrowSKY
 - Platform switch with blueSKY and SKY classic
- A concave and convex abutment shape in the gingival region ensures optimal attachment of the soft tissue
- Customisable axis compensation of up to 20°
- Particularly suitable for custom transverse screw fixation

SKY esthetic line M




All components in the SKY esthetic line are matched to one another:
Corresponding concave and convex sulcus shape.
Also available for abutment shoulders with S and L diameters.


SKY esthetic gingiva former




The SKY esthetic gingiva former gives the emergence profile the optimum shape for the subsequent use of the corresponding SKY esthetic abutments.



REF	SKYESG02	SKYESG03	SKYESG04	SKYESG06
Description	SKY aesthetic gingiva former S	SKY aesthetic gingiva former S	SKY aesthetic gingiva former S	SKY aesthetic gingiva former S
Shoulder height / mm	2	3	4	6
piece	1	1	1	1
Angulation	0°	0°	0°	0°
Shoulder Ø / mm	4.3	4.7	4.7	4.7
Material	Ti*	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓	✓
Torque/Ncm	10	10	10	10
Platform	narrow	narrow	narrow	narrow
narrowSKY	✓	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓	✓





REF	SKYEMG02	SKYEMG03	SKYEMG04	SKYEMG06
Description	SKY aesthetic gingiva former M	SKY aesthetic gingiva former M	SKY aesthetic gingiva former M	SKY aesthetic gingiva former M
Shoulder height / mm	2	3	4	6
piece	1	1	1	1
Angulation	0°	0°	0°	0°
Shoulder Ø / mm	5.2	5.7	5.7	5.7
Material	Ti*	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓	✓
Torque/Ncm	10	10	10	10
Platform	regular	regular	regular	regular
narrowSKY	-	-	-	-
blueSKY / SKY classic	✓	✓	✓	✓



REF	SKYELG02	SKYELG03	SKYELG04
Description	SKY aesthetic gingiva former L	SKY aesthetic gingiva former L	SKY aesthetic gingiva former L
Shoulder height / mm	2	3	4
piece	1	1	1
Angulation	0°	0°	0°
Shoulder Ø / mm	6.5	6.8	7.2
Material	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓
Torque/Ncm	10	10	10
Platform	regular	regular	regular
narrowSKY	-	-	-
blueSKY / SKY classic	✓	✓	✓


Ti* = Grade 4 KV titanium

SKY impression abutments



Open modelling

REF	SKYnPAS1	SKYnPAM1	SKYnPAL1
Description	SKY impression abutment S open tray	SKY impression abutment M open tray	SKY impression abutment L open tray
piece	1	1	1
Angulation	0°	0°	0°
Shoulder Ø / mm	4.5	5.5	7.0
Abutment height / mm	11	11	11
Height incl. screw / mm	24	24	24
Length of Torx®	1.2	1.2	1.2
Material	Ti*	Ti*	Ti*
Trapped screw	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓
Torque/Ncm	10	10	10
Platform	narrow	narrow	narrow
narrowSKY	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓



Closed modelling

REF	SKYnPAS2	SKYnPAM2	SKYnPAL2
Description	SKY impression abutment S closed tray	SKY impression abutment M closed tray	SKY impression abutment L closed tray
piece	1	1	1
Angulation	0°	0°	0°
Shoulder Ø / mm	4.5	5.5	7.0
Abutment height / mm	9	9	9
Length of Torx®	3.5	3.5	3.5
Material	Ti*	Ti*	Ti*
Trapped screw	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓
Torque/Ncm	10	10	10
Platform	narrow	narrow	narrow
narrowSKY	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓

Ti*=
Grade 4 KV titanium


SKY esthetic line temporary abutments



SKY temp is a purely plastic abutment and is therefore suitable for temporary restoration for a maximum of 6 months. There is a risk of loosening of the screws if it is used for longer periods.

SKY temp as a custom gingiva former

The shortened and customised SKY temp can also be quickly adapted into a custom gingiva former, either chairside or in the laboratory.




REF	SKYTEMPS	SKYTEMPM	SKYTEMPL
Description	SKY temp abutment S	SKY temp abutment M	SKY temp abutment L
piece	1	1	1
Angulation	0°	0°	0°
Shoulder Ø / mm	4.1	5.5	7.0
Shoulder height / mm	2.0	3.0	3.75
Material	POM	POM	POM
2.2 screw	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓
Torque/Ncm	18	18	18
Platform	regular	regular	regular
narrowSKY	-	-	-
blueSKY / SKY classic	✓	✓	✓

SKY esthetic line titanium abutments





The concave and convex shape of the SKY esthetic abutments allows the dental technician to customise them to a large extent, and gives the gingiva a lot of space for attachment.



REF	SKYnES00	SKY-EM00	SKY-EL00
Description	SKY esthetic abutment S 0°	SKY esthetic abutment M 0°	SKY esthetic abutment L 0°
piece	1	1	1
Angulation	0°	0°	0°
Torx® alignment	Standard	Standard	Standard
Shoulder Ø / mm	4.5	5.5	7.0
Shoulder height / mm	3.0	3.0	3.0
Material	Ti*	Ti*	Ti*
2.2 screw	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓
Torque/Ncm	25	25	25
Platform	narrow	narrow	narrow
narrowSKY	✓	-	-
blueSKY / SKY classic	✓	✓	✓

Ti*=
Grade 4 KV titanium



REF	SKYnES15	SKYnES16	SKY-EM15	SKY-EM16
Description	SKY esthetic abutment S 15°	SKY esthetic abutment S 15° R	SKY esthetic abutment M 15°	SKY esthetic abutment M 15° R
piece	1	1	1	1
Angulation	15°	15°	15°	15°
Torx® alignment	Standard	R = turned by 30°	Standard	R = turned by 30°
Shoulder Ø / mm	4.5	4.5	5.5	5.5
Shoulder height / mm	3.0	3.0	3.0	3.0
Material	Ti*	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓	✓
Torque/Ncm	25	25	25	25
Platform	narrow	narrow	narrow	narrow
narrowSKY	✓	✓	-	-
blueSKY / SKY classic	✓	✓	✓	✓

SKY standard line titanium abutments




The SKY standard line that has proven itself for years for cost-effective prosthetics, now new with narrow platform.

- You can now also restore using narrowSKY
- blueSKY and SKY classic now have a platform switch



REF	SKYnPO00
Description	SKY titanium abutment NP 0°
piece	1
Angulation	0°
Torx® alignment	Standard
Shoulder Ø / mm	4.5
Shoulder height / mm	1.4
Material	Ti*
2.2 screw	incl.
SKY prosthetic key	✓
Torque/Ncm	25
Platform	narrow
narrowSKY	✓
blueSKY / SKY classic	✓

Ti*=
Grade 4 KV titanium




REF	SKYnPO15	SKYnPO16	SKYnPO25	SKYnPO26
Description	SKY titanium abutment NP 15°	SKY titanium abutment NP 15°	SKY titanium abutment NP 25°	SKY titanium abutment NP 25° R
piece	1	1	1	1
Angulation	15°	15°	25°	25°
Torx® alignment	Standard	R = turned by 30°	Standard	R = turned by 30°
Shoulder Ø / mm	4.5	4.5	4.5	4.5
Shoulder height / mm	0.5	0.5	0.5	0.5
Material	Ti*	Ti*	Ti*	Ti*
2.2 screw	incl.	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓	✓
Torque/Ncm	25	25	25	25
Platform	narrow	narrow	narrow	narrow
narrowSKY	✓	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓	✓

SKY abutment cast-on



With cast-on abutments, difficult individual clinical situations can also be easily handled from a prosthetic point of view.

This is held securely by screwing the sleeve onto the metal base, so as to facilitate quick and reliable modelling of the custom abutment. The burn-out plastic sleeve is provided already mounted.




REF	SKYnPV00
Description	SKY abutment cast-on
piece	1
Angulation	0°
Shoulder Ø / mm	4.5
Shoulder height / mm	3.5
Sleeve height / mm	9.5
Sleeve material	PMMA
Base height	2.7
Base material	Au 60%, Pd 20%, Pt 19%, Ir 1%
Melting range	1400 - 1490°C
CTE	11.9 - 12.2 10 ⁻⁶ K ⁻¹
Weight	0.63 g
2.2 screw	incl.
SKY prosthetic key	✓
Torque/Ncm	25
Platform	narrow
narrowSKY	✓
blueSKY / SKY classic	✓

SKY prosthetics

- Implant and abutment platforms
- Implant connections
- Prosthetics overview
- Classic implantology
 - Implant analogue
 - SKY esthetic line
 - SKY esthetic gingiva former
 - SKY impression abutments
 - SKY esthetic line temporary abutments
 - SKY esthetic line titanium abutments
 - SKY standard line titanium abutments
 - SKY castable abutment
- CAD/CAM-manufactured restorations
 - SKY uni.fit scan abutments
 - SKY uni.fit titanium base
 - SKY prefab titanium
 - SKY uni.fit titanium base for CEREC
- Immediate and late restoration
 - SKY elegance abutments
 - SKY elegance prefabs
 - SKY elegance titanium base
 - SKY fast & fixed / uni.cone immediate restoration of the complete jaw
 - SKY fast & fixed / uni.cone order information
 - SKY fast & fixed – accessories
- Prosthesis fixation
 - SKY TiSi-snap
 - SKY Locator®

SKY uni.fit Scan abutments





Scan abutments

The position and orientation of the implant is transferred to the virtual model using a three-point reference system.

All SKY uni.fit components are recorded in the CAD programs:

- exoCAD
- 3shape
- Dental Wings

Download: www.caelo-dental.net





REF	SKYUSCAE	SKYUSCAI
Description	SKY uni.fit scan abutment Extraoral	SKY uni.fit scan abutment Intraoral
piece	1	1
Angulation	0°	0°
Shoulder Ø / mm	3.6	3.6
Height/mm	13.5	7.5
Material	PEEK	PEEK
2.2 screw	incl.	incl.
SKY prosthetic key	✓	✓
Torque/Ncm	10	10
Platform	narrow	narrow
narrowSKY	✓	✓
blueSKY / SKY classic	✓	✓
CAD library for: (www.caelo-dental.net)	exoCAD 3shape Dental Wings	exoCAD 3shape Dental Wings

SKY uni.fit CAD abutment



The geometry of the SKY uni.fit CAD abutments is such that the custom ceramic structure can be milled to create a perfect fit. The space for excess adhesive facilitates bonding.



REF	SKYUFCAD	UFCADMOD
Description	SKY uni.fit CAD abutment	Modelling cap
piece	1	10
Angulation	0°	0°
Shoulder Ø / mm	3.6/4.5	3.8
Height/mm	3.8	10
Material	Ti*	PMMA
2.2 screw	incl.	excl.
SKY prosthetic key	✓	✓
Torque/Ncm	25	10
Platform	narrow	narrow
narrowSKY	✓	✓
blueSKY / SKY classic	✓	✓
CAD library for: (www.caelo-dental.net)	exoCAD 3shap Dental Wings	-

Ti*=
Grade 4 KV titanium

SKY prefab titanium




Information for processing:

- SKY prefab titanium is recorded in the following CAD libraries:
 - 3shape
 - Dental Wings
 - exoCAD
 Download at: www.caelo-dental.com
- The SKY prefab titanium is held in the Medentika holder
- If the SKY prefab cannot be controlled by the CAM, the Medentika workflow is then used.
 - Scanning the implant position with the Medentika scan abutment
 - Clamping the SKY prefabs titanium in the R position





REF	SKYPFTST
Description	SKY prefab titanium
piece	1
Angulation	0°
Ø / mm	11.5
Height/mm	20
Material	Ti*
2.2 screw	incl. in the original packaging
SKY prosthetic key	✓
Torque/Ncm	25
Platform	narrow
narrowSKY	✓
blueSKY / SKY classic	✓
Connector geometry	Medentika holder

Ti*=
Grade 4 KV titanium



BioHPP prefab
see page 50


SKY uni.fit titanium base for CEREC®



The SKY uni.fit titanium base for CEREC® combines the SKY Implant abutment connection with the abutment geometry of the CEREC® TiBase S. It is therefore possible to produce custom abutments for restoration with SKY Implants using the CEREC® system.

The SKY uni.fit titanium base for CEREC® combines the SKY Implant abutment connection with the abutment geometry of the CEREC® TiBase S. It is therefore possible to produce custom abutments for restoration with SKY Implants using the CEREC® system.

- The implant position is scanned using the original Sirona® scan bodies:
 - Scan bodies for Bluecam® S, 36 Stück REF 6431329
 - Scan bodies for Omnicam® S, 36 Stück REF 6431303
- For construction in the CEREC® software, a suitable implant is selected from the library, e.g. Camlog® 3.8.
- All CEREC® blocks with the S geometry can be used for this process.
- The subsequent bonding is carried out according to the manufacturer's instructions.



REF	SKYUFCTB	SKYUFCSB
Description	SKY uni.fit titanium base for CEREC®	SKY uni.fit scan base for CEREC®
piece	1	1
Angulation	0°	0°
Ø / mm	3.0	3.0
Total height / mm	4.7	10.2
Gingiva height / mm	0.5	5.5
Material	Ti*	Ti*
2.2 screw	incl.	incl.
SKY prosthetic key	✓	✓
Torque/Ncm	10	10
Platform	narrow	narrow
narrowSKY	✓	✓
blueSKY / SKY classic	✓	✓

Ti*=
Grade 4 KV titanium

SKY prosthetics

- Implant and abutment platforms
- Implant connections
- Prosthetics overview
- Classic implantology
 - Implant analogue
 - SKY esthetic line
 - SKY esthetic gingiva former
 - SKY impression abutments
 - SKY esthetic line temporary abutments
 - SKY esthetic line titanium abutments
 - SKY standard line titanium abutments
 - SKY castable abutment
- CAD/CAM-manufactured restorations
 - SKY uni.fit scan abutments
 - SKY uni.fit titanium base
 - SKY prefab titanium
 - SKY uni.fit titanium base for CEREC
- Immediate and late restoration
 - SKY elegance overview
 - SKY elegance abutments
 - SKY elegance prefabs
 - SKY elegance titanium base
 - SKY fast & fixed – overview
 - SKY fast & fixed – components
 - SKY fast & fixed – accessories
- Prosthesis fixation
 - SKY TiSi-snap
 - SKY Locator®

SKY elegance overview

The BioHPP SKY elegance abutments are a new type of abutment that combine the properties of a temporary and a permanent abutment.

- The elasticity of the BioHPP body protects the implant during healing
- The seat of the screw made in titanium ensures long-term, stable connection of the abutment to the implant

They are therefore ideal for one-time treatment without changing the abutment in immediate restoration. These abutments can of course also be used in late restoration. The implant is protected against overloading in the long-term by the elasticity.


Pre-fabricated abutments



Pre-fabricated abutments


The pre-fabricated BioHPP SKY elegance abutments are universal abutments in the shape of the SKY esthetic line and can be used both with CAD/CAM-manufactured crowns and bridges, e.g. CEREC, as well as conventionally-manufactured crowns and bridges.

Individual abutments



Laboratory-manufactured individual abutments

Any laboratory is able to manufacture an individual abutment on the BioHPP SKY elegance titanium base using the for2press device.



CAD/CAM-manufactured individual abutments

With suitable CAD/CAM systems, an individual abutment can be manufactured using the BioHPP SKY elegance prefab. The titanium base is injected already gap-free in the pre-fab.


SKY elegance abutments



For intraoral processing, we recommend Komet:



REF H 379 Q 314 023










REF H 375 RQ 314 016

The BioHPP SKY elegance abutments are hybrid abutments in which the abutment body made of BioHPP is connected to the titanium base without a gap. These abutments are best used for One-Time Therapy for immediate restoration, since they combine the properties of a temporary and a definitive abutment, i.e. it is not necessary to change the abutment. As a result, the gingiva is not subjected to multiple traumas. In addition, the time and costs are reduced.

BioHPP can be ground in the mouth as easily as dentine using carbide milling tools.

For extraoral processing, we recommend the bredent cutter set SKY elegance REF 580ELEM8



					
REF	SKYEES00	SKYEES15	SKYEEM00	SKYEEM15	SKYEEL00
Description	BioHPP SKY elegance abutment S 0°	BioHPP SKY elegance abutment S 15°	BioHPP SKY elegance abutment S 15°	BioHPP SKY elegance abutment M 15°	BioHPP SKY elegance abutment L 0°
piece	1	1	1	1	1
Angulation	0°	15°	0°	15°	0°
Shoulder Ø / mm	4.5	4.5	5.5	5.5	7.0
Shoulder height / mm	3.4	3.4	3.4	3.3	3.4
Body material	BioHPP	BioHPP	BioHPP	BioHPP	BioHPP
Base material	Ti*	Ti*	Ti*	Ti*	Ti*
2.2 screw	incl.	incl.	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓	✓	✓
Torque/Ncm	25	25	25	25	25
Platform	narrow	narrow	narrow	narrow	narrow
narrowSKY	✓	✓	-	-	-
blueSKY / SKY classic	✓	✓	✓	✓	✓

Ti*=
Grade 4 KV titanium

BioHPP SKY elegance prefab




With the BioHPP SKY elegance prefab, the abutment body made of BioHPP is pressed onto the BioHPP SKY elegance titanium base without a gap and forms a perfect mechanical connection. The required tooth shape for the custom abutment is designed in CAD software and the corresponding data set is passed on to the machine manufacturing stage.

For CAD/CAM production



REF	SKYEPFST
Description	BioHPP SKY elegance prefab
piece	1
Angulation	0°
Ø / mm	12.0
Height / mm	18.5
Connection	Sirona® Standard
Base material	Aluminum
Material construction	BioHPP
Color construction	dentín-shade
X-ray opaque	✓
Screw2,2	incl. Original packaging
SKY prosthetic key	✓
Torque / Ncm	25
Platform	narrow
narrowSKY	✓
blueSKY / SKY classic	✓

BioHPP SKY elegance titanium base





A custom abutment in a natural tooth shape is modelled on the BioHPP SKY elegance titanium base (sand-blasted by the technician). It is then embedded and remoulded with BioHPP in the *for2press* device.

This custom BioHPP abutment can then be directly veneered with the *visio.lign* veneer system to form a crown abutment or can be restored using a crown or a bridge.

For individual abutments or crown abutments in the region of the side teeth, the reinforced BioHPP SKY elegance titanium base ML is used.

Konventionelle Fertigung



REF	SKYETB00	SKYETBML
Description	BioHPP SKY elegance titanium base	BioHPP SKY elegance titanium base ML
piece	1	1
Indication	Front tooth area	Range of teeth
Angulation	0°	0°
Ø / mm	3.1	3.5
Hight / mm	5.4	5.4
Material	Ti*	Ti*
Screw 2.2	incl.	incl.
SKY prosthetic key	✓	✓
Torque / Ncm	25	25
Platform	narrow	narrow
narrowSKY	✓	
blueSKY / SKY classic	✓	✓

Ti*=
Grade 4 KV titanium


SKY fast & fixed – overview



SKY fast & fixed


Immediate restoration for potentially edentulous jaws

- Reduced number of implants
- No extensive surgical procedures such as augmentations
- Standardised work steps make the work easier
- Reduction and prevention of errors and complications
- In most cases, immediate fixed temporary bridges after only one procedure
- At an affordable price





SKY fast & fixed / SKY uni.cone

- One-time treatment – no change of abutment required
- Two shoulder diameters:
 - SKY fast & fixed: 5.65 mm
 - SKY uni.cone: 4.5 mm
- Two types of screwing:
 - Occlusally screwed
 - Transversally screwed



SKY fast & fixed / SKY uni.cone modelling

- Modelling at abutment level
 - Abutment does not need to be removed



SKY fast & fixed / SKY uni.cone transversal screwing

Self-centering transversal screw retention


Transversal screw retention is a bolting principle. The thread for the bolt screw is located in the bridge framework. The bolt screw (A) and the cylindrical surfaces (B and C) form a compact unit. The prosthetic coping is fixed by triple point attachment with the bolt screw (A) and the short cylindrical surfaces (B and C) to avoid tilting. At the same time, the short cylindrical surfaces ensure that the coping is self-centered (passive fit) when it is placed on. Thanks to the slightly inclined position of the bolt screw, the prosthetic coping is pressed on the abutment platform and gap formation is avoided when it is tightened.

Simple application

Using the pre-fabricated transversal screwing is very simple:

- because the screw always remains in the secondary structure, therefore threading is not necessary.
- because only a few rotations are required to fix or loosen the restoration.
- because the precision groove running all the way around enables a high degree of freedom when designing the access to the screw.

Combinable SKY fast & fixed and SKY uni.cone



SKY fast & fixed abutments with narrow Platform for

- blueSKY
- SKY classic

 Gingiva height in mm


SKY uni.cone abutments with narrow Platform for

- narrowSKY
- blueSKY
- SKY classic


SKY fast & fixed - components



REF	SKYFn354	SKYFn355	SKYFn173	SKYFn175
Description	SKY fast & fixed abutment 35°	SKY fast & fixed abutment 35°	SKY fast & fixed abutment 17.5°	SKY fast & fixed abutment 17.5°
Shoulder height / mm	4.0	5.0	3.3	5.0
piece	1 set	1 set	1 set	1 set
Gingiva height / mm	1.0	1.8	2.0	3.6
Angulation	35°	35°	17.5°	17.5°
Abutment height / mm	3.6	3.6	3.6	3.6
Shoulder Ø / mm	5.65	5.65	5.65	5.65
Insertion aid	incl.	incl.	incl.	incl.
Material	Ti*	Ti*	Ti*	Ti*
2.2 screw	incl.	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓	✓
Torque/Ncm	25	25	25	25
Platform	narrow	narrow	narrow	narrow
narrowSKY	✓	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓	✓




REF	SKYFn001	SKYFn002	SKYFn004
Description	SKY fast & fixed abutment 0° NP	SKY fast & fixed abutment 0° NP	SKY fast & fixed abutment 0° NP
Shoulder height / mm	1.0	2.0	4.0
piece	1	1	1
Angulation	0°	0°	0°
Abutment height / mm	3.6	3.6	3.6
Shoulder Ø / mm	5.65	5.65	5.65
Material	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓
Torque/Ncm	25	25	25
Platform	regular	regular	regular
narrowSKY	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓




REF	SKYUC001	SKYUC002	SKYUC003
Description	SKY uni.cone abutment 0°	SKY uni.cone abutment 0°	SKY uni.cone abutment 0°
Shoulder height / mm	1.0	2.0	3.0
piece	1	1	1
Angulation	0°	0°	0°
Abutment height / mm	3.3	3.3	3.3
Shoulder Ø / mm	4.5	4.5	4.5
Material	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓
Torque/Ncm	25	25	25
Platform	narrow	narrow	narrow
narrowSKY	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓

Ti*=
Grade 4 KV titanium




REF	SKYFFPKT	SKYUCPKT	SKYFFPKC	SKYUCPKC
Description	SKY fast & fixed prosthetic coping	SKY uni.cone prosthetic coping titanium	SKY fast & fixed prosthetic coping CAD/CAM	SKY uni.cone prosthetic coping CAD/CAM
piece	1	1	1	1
Construction height / mm	10.0	10.0	10.0	10.0
Angulation	0°	0°	0°	0°
Shoulder Ø / mm	5.65	4.5	5.65	4.5
Silicone tubing	incl.	incl.	-	-
Locking pin	incl.	incl.	-	-
Material	Ti*	Ti*	Ti*	Ti*
Screw	incl.	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓	✓
Torque/Ncm	18	18	18	18



REF	SKYFFPKK	SKYUCPKK	SKYFFPKH	SKYFTPKS	SKYUCPKS
Description	SKY fast & fixed prosthetic coping plastic	SKY uni.cone prosthetic coping plastic	SKY fast & fixed prosthetic coping HSL cast-on	SKY fast & fixed prosthetic coping transversal	SKY uni.cone prosthetic coping transversal
piece	1	1	1	1	1
Construction height / mm	12.4	12.4	10.0	5.0	5.0
Angulation	0°	0°	0°	0°	0°
Shoulder Ø / mm	5.65	4.5	4.5	6.1	5.1
Modelling cap	-	-	incl.	incl.	incl.
Material	Ti*	Ti*	Pt 90%, Ir 10%	Ti*	Ti*
Melting range	-	-	1770 - 1800°C	-	-
CTE	-	-	11.9 - 12.2 10 ⁻⁶ K ⁻¹	-	-
Weight	-	v	0.59 g	-	-
Sleeve material	-	-	PMMA	-	-
Screw	incl.	incl.	incl.	-	-
SKY prosthetic key	✓	✓	✓	-	-
Inbus 0.9	-	-	-	✓	✓
Torque/Ncm	18	18	18	18	18





Ti* = Grade 4 KV titanium








REF	SKYFFSPK	SKYFFLPK	SKYUFTS9
Description	SKY fast & fixed / SKY uni.cone M 1.4 screw	SKY fast & fixed / SKY uni.cone lab screw M 1.4	SKY fast & fixed / SKY uni.cone transversal screw
piece	6	10	6
Screw length / mm	6.0	6.0	1.0
Thread	M 1.4	M 1.4	M 2.0
Head Ø / mm	2.2	2.2	1.0
Material	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	-
Inbus 0.9	-	--	✓


SKY fast & fixed – components







				
REF	SKYFFSNP	SKYUCSNP	SKYFFREG	SKYUCREG
Description	SKY fast & fixed snap coping	SKY uni.cone snap coping	SKY fast & fixed bite registration	SKY uni.cone bite registration
piece	10	10	10	10
Construction height / mm	8	8	5	5
Angulation	0°	0°	0°	0°
Shoulder Ø / mm	5.65	4.5	5.65	4.5
Material	POM	POM	POM	POM



				
REF	SKYFFAGK	SKYUCAGK	SKYFFAOL	SKYUCAOL
Description	SKY fast & fixed modelling closed tray	SKY uni.cone modelling closed tray	SKY fast & fixed modelling open tray	SKY uni.cone modelling open tray
piece	1	1	1	1
Construction height / mm	10	10	15	15
Angulation	0°	0°	0°	0°
Shoulder Ø / mm	5.65	4.5	5.65	4.5
Material	Ti*	Ti*	Ti*	Ti*
Screw	incl.	incl.	incl.	incl.
SKY prosthetic key	✓	✓	✓	✓




				
REF	SKYFSCIE	SKYUSCIE	SKYFFTLA	SKYUCTLA
Description	SKY fast & fixed scan coping intraoral / extraoral	SKY uni.cone scan coping intraoral (extraoral)	SKY fast & fixed laboratory analogue	SKY uni.cone laboratory analogue
piece	1	1	1	1
Construction height / mm	12.4	12.4	12	12
Angulation	0°	0°	0°	0°
Shoulder Ø / mm	5.65	4.5	5.65	4.5
Material	PEEK	PEEK	Ti*	Ti*
Screw	incl.	incl.	incl.	incl.
SKY prosthetic key	✓	✓	-	-




REF	SKYFFS35
Description	SKY fast & fixed angulation aid set 35°
piece	1

Ti* = Grade 4 KV titanium

SKY fast & fixed – accessories





SKY fast & fixed bridge kit




SKY fast & fixed bridge kit maxilla REF 580FFBOK

- top.lign professional crown and bridge material
- novo.lign veneers
 - I47, A3
 - L3, A3



SKY fast & fixed bridge kit mandible REF 580FFBUK

- top.lign professional crown and bridge material
- novo.lign veneers
 - D38 A3
 - L3 A3



flow.sil – microgap sealing

Sealing of the gaps between the abutment and prosthetic restoration

- Has an antimicrobial effect
- Prevents or reduces bacterial colonisation
- Prevents or reduces odour formation

flow.sil Set contains:
5 ml double cartridge
4 mixing cannulas
4 tips
REF 54001270

SKY prosthetics

- Implant and abutment platforms
- Implant connections
- Prosthetics overview
- Classic implantology
 - Implant analogue
 - SKY aesthetic line
 - SKY aesthetic gingiva former
 - SKY impression abutments
 - SKY aesthetic line temporary abutments
 - SKY aesthetic line titanium abutments
 - SKY standard line titanium abutments
 - SKY castable abutment
- CAD/CAM-manufactured restorations
 - SKY uni.fit scan abutments
 - SKY uni.fit titanium base
 - SKY prefab titanium
 - SKY uni.fit titanium base for CEREC
- Immediate and late restoration
 - SKY elegance overview
 - SKY elegance abutments
 - SKY elegance prefabs
 - SKY elegance titanium base
 - SKY fast & fixed – overview
 - SKY fast & fixed – components
 - SKY fast & fixed – accessories
- Prosthesis fixation
 - SKY TiSi.snap
 - SKY Locator®

SKY TiSi.snap



Contents:
retention.sil in 3 hardnesses
in the double-mix cartridge
REF 580 RT SET




retention.sil 200
Shore hardness 25 SH
Pull-off forces 200 g / 2 Newton
REF 580RTS25




retention.sil 400
Shore hardness 50 SH
Pull-off forces 400 g / 4 Newton
REF 580RTS50



retention.sil 600
Shore hardness 65 SH
Pull-off forces 600 g / 6 Newton
REF 580RTS65




Multisil-Primer
5 ml
REF 520 0100 4



Special Silicon Trimmer
Ø 4.1 mm
REF SKY-DR41

The special silicone cutter is suitable for angled and lab handpiece.

Cost-effective prostheses can be fixed using TiSi.snap and retention.sil, as the existing prosthesis only requires grinding. Hold is guaranteed by retention.sil, which is filled into the cavity. The restoration is very resilient and offers a high level of wearing comfort.



REF	TISIOY31	TISIOY51	TISIOY53	TISIA Y17	TISIA Y35
Description	SKY TiSi.snap abutment 3/1	SKY TiSi.snap abutment 5/1	SKY TiSi.snap abutment 5/3	SKY TiSi.snap abutment 17,5°	SKY TiSi.snap abutment 35°
pieces	1	1	1	1	1
Angulation	0°	0°	0°	0°	0°
Shoulder Ø / mm	4.93	4.93	4.93	5.65	5.65
Gingiva height / mm	1.34	1.34	3.34	3.39	3.99
Material	Ti*	Ti*	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓	✓	✓
Torque / Ncm	25	25	25	25	25
Platform	narrow	narrow	narrow	regular	regular
narrowSKY	✓	✓	✓	-	-
blueSKY / SKY classic	✓	✓	✓	✓	✓
retention.sil	✓	✓	✓	✓	✓
Retentions-Ø / mm	3.87	3.87	3.87	3.87	3.87
Locator®	0° - 10°	0° - 10°	0° - 10°	0° - 10°	0° - 10°
Retention elements	10° - 20°	10° - 20°	10° - 20°	10° - 20°	10° - 20°

Ti*=
Grade 4 KV titanium

SKY Locator®



Crown and bridge material
Qu-resin rosa
REF 540 0116 1


Due to its low structural height, the SKY Locator® offers excellent possibilities for fixing prostheses in many cases.

The long-lasting stability of the bone around the implant is supported by the built-in Platform switch.

The 5 gingival heights of 1, 2, 3, 4, and 6 mm cover all of the relevant clinical situations.



REF	LOCLAB10	LOCLAB20
Description	SKY Locator® Processing set 0°-10°	Locator® Processing set 10°-20°
pieces	2 Sets	2 Sets
Material	Titan / Teflon / Nylon	Titan / Teflon / Nylon



REF	LOCZAB01	LOCZAB02	LOCZAB03	LOCZAB04	LOCZAB06
Description	SKY Locator® abutment for SKY 1 mm	SKY Locator® abutment for SKY 2 mm	SKY Locator® abutment for SKY 3 mm	SKY Locator® abutment for SKY 4 mm	SKY Locator® abutment for SKY 6 mm
pieces	1	1	1	1	1
Angulation	0°	0°	0°	0°	0°
Gingiva height	1	2	3	4	6
Material	Ti*	Ti*	Ti*	Ti*	Ti*
Coating	TiNi	TiNi	TiNi	TiNi	TiNi
Locator-Instruments	✓	✓	✓	✓	✓
Torque / Ncm	25	25	25	25	25
Platform	narrow	narrow	narrow	narrow	narrow
narrowSKY	✓	✓	✓	✓	✓
blueSKY / SKY classic	✓	✓	✓	✓	✓
retention.sil	✓	✓	✓	✓	✓
Retention Ø / mm	3.87	3.87	3.87	3.87	3.87
Retention elements 0°-10°	✓	✓	✓	✓	✓
Retention elements 10°-20°	✓	✓	✓	✓	✓

Ti*=
Grade 4 KV titanium



In atrophied jaws, implants can frequently only be inserted obliquely, which creates increased wear on the abutments and the retention elements. The direction of insertion is corrected by the angled Locator®.

The following angulations are available:

- 17.5°
- 35°

It is therefore now also possible to support prostheses from a posterior direction using implants inserted at an angle.



REF	LOCAB172	LOCAB352	LOCLAB10	LOCLAB20
Description	SKY Locator® abutment 17,5°	SKY Locator® abutment 35°	SKY Locator® Processing set 0°-10°	SKY Locator® Processing set 10°-20°
pieces	1	1	2 Sets	2 Sets
Angulation	17.5°	35°	-	-
Sholder Ø / mm	5.65	5.65	-	-
Gingiva height / mm	3.8	3.8	-	-
Coping Material	Titan Grade V	Titan Grade V	-	-
Coping coating	TiNi	TiNi	-	-
Basi material	Titan Grade IV	Titan Grade IV	-	-
Material	-	-	Titanium / Teflon / Nylon	Titanium / Teflon / Nylon
Screw	incl.	incl.	-	-
SKY prosthetic key	✓	✓	-	-
Torque / Ncm	25	25	-	-
Platform	regular	regular	-	-
narrowSKY	-	-	-	-
blueSKY / SKY classic	✓	✓	-	-
Retention Ø / mm	3.87	3.87	-	-
Retention elements 0° - 10°	✓	✓	-	-
Retention elements 10° - 20°	✓	✓	-	-

SKY Locator® Prosthesis fixation



SKY Locator®
Instrument
1 pieces
REF LOCZINST



SKY Locator®
Angular insertion
instrument, 1 pieces
REF LOCZWED6




Impression coping
4 pieces
REF LOCZAK40



Laboratory analog
4 pieces
REF LOCZLA40



Blocking out ring
20 pieces
REF LOCblock



Processing insert
black
4 pieces
REF LOCZVA11



Angular measuring
abutment
4 pieces
REF LOCZWIMP



Angular measuring
gauge
1 pieces
REF LOCZWIML



Retention insert
blue, 6.7 N, 680 g
4 pieces
REF LOCR1006



Retention insert
red, 6.7 N, 680 g
4 pieces
REF LOCR2006



Retention insert
pink, 13.4 N, 1365 g
4 pieces
REF LOCR1013



Retention insert
orange, 9.1 N, 907 g
4 pieces
REF LOCR2009



Retention insert
transparent, 22.3 N,
2270 g
4 pieces
REF LOCR1022



Retention insert
green, 17.8 N, 1815 g
4 pieces
REF LOCR2018



Prosthesis fixation
and restoration of narrow single tooth gaps

miniSKY for narrow alveolar ridges



Prosthesis fixation


With miniSKY, bredent medical offers an implant that is reduced in diameter and perfectly suited to prosthesis fixation due to its technical properties.



Restoration of narrow single tooth gaps

mini²SKY is perfectly suited to restoration of narrow single tooth gaps. In the case of low bone availability and despite residual dentition, patients can be treated with a highly-aesthetic solution for the edentulous space.


The miniSKY implant system
Implant and surface design




osseo connect surface (ocs)[®]

The mini¹SKY and mini²SKY implants have the tried-and-tested osseo-connect surface (ocs)[®] of the blueSKY implants, which ensures optimal osseointegration.



The uniform coarse surface provides ideal prerequisites for the accumulation of the osteoblasts. This is supported by the excellent hydrophilic properties of the implants.



Competing implant abutment connection



mini²SKY implant abutment connection



- High implant abutment connection
- Minimal movement due to 5° cone
- Exceptional stability

Definitive positioning of the abutment

mini¹SKY

- Temporary implant
- Fixation of templates



Length 6 mm

REF m1SKYL06

Length 10 mm

REF m1SKYL10


Length 12 mm

REF m1SKYL12

Length 14 mm

REF m1SKYL14

mini²SKY 2.8



Length 10 mm

REF m2SKYL10

Length 12 mm

REF m2SKYL12

Length 14 mm

REF m2SKYL14

mini²SKY 3.2



Length 8 mm

REF m2SK3208

Length 10 mm

REF m2SK3210

Length 12 mm


REF m2SK3212

Length 14 mm


REF m2SK3214

Prosthesis fixation

SKY O-ring housing or SKY precious metal matrix




mini²SKY retention.loc



mini²SKY TiSi.snap with retention.sil – the economical alternative

CAD/CAM manufactured bars – directly screwed



Restoration of narrow single-tooth gaps

Standard abutments




Individual abutments – analogue and digital



Component overview

Prosthetic fixation for mini¹SKY



Prosthetic fixation for mini²SKY



Prosthetic fixation with CAD/CAM manufactured bar on mini²SKY



Restoration of narrow single-tooth gaps with mini²SKY



Impression abutments, laboratory analogues, gingiva former and screws for mini²SKY



For broad alveolar ridges with low height – copaSKY, short and good



For years, the SKY implant system with narrow and angled implants has been a pioneer of minimally-invasive surgical procedures, with the aim of making optimal use of the existing bone and avoiding augmentation as far as possible.

With the new copaSKY, we now offer an impressive alternative for the restoration of broad, flat jaws with short implants.



copaSKY also has the DNA of all SKY implant systems: The only elements that have been changed are those that were necessary for the compact construction

- Mainly the new conical, parallel-walled very stable implant-abutment connection
- The implant has a back taper so that attachment of bone chips is possible
- The implant position is iso-crestal or slightly sub-crestal
- Single-start screw thread for a better feeling

The tried-and-tested properties of the SKY system have been retained:

- One connector geometry for all diameters
- Torx as rotation protection
- Manageable compact prosthetics
- Conical cylindrical implant shape
- Tried-and-tested *osseo-connect* surface (ocs)[®]
- Self-tapping deepened compression thread
- High level of primary stability
- Same surgical protocol

At the beginning, the copaSKY implant line is limited to what is important, and is expanded little by little by the tried-and-tested prosthetic elements of the SKY implant system:

- copaSKY elegance
- copaSKY uni.fit CAD
- copaSKY uni.cone abutment




copaSKY impression abutment open tray
The proven design of the SKY abutments for open trays was adopted.
REF copaPAM1




copaSKY mounter
Only one new implant mounter is required for use.

copaSKY mounter short
REF copaCTK5


copaSKY mounter long
REF copaCTK6




copaSKY Gingiva former open tray
The proven design of the SKY aesthetic abutments was adopted.
REF copaGF04




copaSKY implant analogue
Only one implant analogue for all diameters.
REF copaIA50



copaSKY exso abutment
This multi-functional abutment is suitable for closed modelling as well as a permanent titanium abutment.
REF copaEXSO



copaSKY titanium base for CEREC
The titanium base for the CEREC workflow for chairside manufacture of the superstructure.
REF copaSCTB



copaSKY elegance titanium base ML
The titanium base for the laboratory manufacture of individual abutments made from BioHPP using the for2press system
REF copaETBM

<p>Ø 4 mm</p> 	<p>Length 5.2 mm REF copa4005</p>	<p>Ø 5 mm</p> 	<p>Length 5.2 mm REF copa5005</p>	<p>Ø 6 mm</p> 	<p>Length 5.2 mm REF copa6005</p>	 <p>Including a cover screw</p>
---	---------------------------------------	---	---------------------------------------	---	---------------------------------------	--

The aesthetic zirconium implant



11 years experience – scientifically proven

High success rates

- Amberger et. al 2015 – 94.6% after 4 years
- Grassi et al. 2015 – 96.4% after 5 years

Highly stable implant

- Kohal et al. 2011

Excellent soft tissue attachment

- Stadlinger et al. 2011

Stable bone level

- Amberger et. al 2015
- Borgonovo et al. 2013
- Grassi et al. 2015

Recommendations for using the whiteSKY zirconium implants

Very thorough implant planning required.

Individual tooth:

Correct implant position is crucial for the aesthetic result

Bridges:

- Parallel positioning of the implants for the restoration required
- +/- 10° divergence balance possible
- Use drilling templates

Make use of the surgical options to achieve optimal primary stability

Temporarily restore the implant immediately

- Elastic prosthetic materials
- Interlocking with neighbouring teeth if the primary stability < 30 Ncm

Modelling like a natural tooth stump

Permanent restoration

- Using modern high-performance polymers such as BioHPP, which are veneered with composite materials such as visio.lign
- Using ceramic

Literature:

Amberger et al.
Immediate provisional restoration of single-piece zirconia implants: 4 years follow up
Med Oral Patol Oral Cir Bucal

Borgonovo et al.
Behavior of endosseous one-piece yttrium stabilized zirconia dental implants placed in posterior areas
Minerva Stomatol 2013 62, 247-57

Calvo-Guirado et al.
Histological, radiological and histomorphometric evaluation of immediate vs. non-immediate loading of a zirconia implant with surface treatment in a dog model
Clinical Oral Implants Research, Volume 25, Issue 7, July 2014, Pages: 826-830, Jose


Grassi et al.
Immediate Occlusal Loading of One-Piece Zirconia Implants: 5 years Radiographic and Clinical Evaluation
JOMI, vol. 30, N° 3, 2015

Kohal et al.
The effects of cyclic loading and preparation on the fracture strength of zirconium-dioxide implants: an in vitro investigation.
Clin. Oral Implants Res, 2011 Aug 22 (8); S. 808-814

Payer et al.
Immediate provisional restoration of single-piece zirconia implants: a prospective case series – results after 24 months of clinical function
COIR, Volume 24, Issue 5, May 2013, Pages: 569-575

Stadlinger et al.
Comparison of zirconia and titanium implants after a short healing period. A pilot study in minipigs.
Clinical & Maxillofacial Surgery, Vol. 3, 1, 2010

whiteSKY sizes




Ø 3.5 mm

Length 10 mm	REF SKY3510C
Length 12 mm	REF bSKY3512C
Length 14 mm	REF SKY3514C
Length 16 mm	REF SKY3516C



Ø 4.0 mm

Length 8 mm	REF SKY4008C
Length 10 mm	REF SKY4010C
Length 12 mm	REF SKY4012C
Length 14 mm	REF SKY4014C
Length 16 mm	REF SKY4016C



Ø 4.5 mm

Length 8 mm	REF SKY4508C
Length 10 mm	REF SKY4510C
Length 12 mm	REF SKY4512C
Length 14 mm	REF SKY4514C



OP-Kit OT21
REF SKYXOT21






Set for grinding zirconium whiteSKY
REF 580E006C




REF	SKYCPK40	SKYCPK45
Description	whiteSKY Prosthetic cap Ø 3.5 mm and 4.0 mm	whiteSKY Prosthetic cap Ø 4.5 mm
Pieces	1	1
Height / mm	8.0	8.0
Ø / mm	5.3	5.3
Material	PEEK	PEEK

Accessories and instruments for the SKY implant lines


Drill


Drill stops			L6	L8	L10	L12	L14	L16
Twistdrill		REF	SKYXST06	SKYXST08	SKYXST10	SKYXST12	SKYXST14	SKYXST16
Drill		REF	-	SKYS0840	SKYS1040	SKYS1240	SKYS1440	SKYS1640
Drill		REF	-	SKYS0845	SKYS1045	SKYS1245	SKYS1445	-


Bone bur rpm/min rpm 800-1.000

Ø 4.1  REF SKY-DR41


Twistdrill rpm 800-1.000


Ø 1.3  REF SKYDT13L

Ø 2.25  short REF SKYDT23K


Ø 2.25  long REF SKYDT23L


Pilot drill rpm 800-1.000


Ø 3.1  short REF SKY-DP06


Ø 3.1  long REF SKY-DP08

Drill for hard bone rpm 300


Ø 3.2  REF SKYD1235 3.5 3.5 N


Ø 3.8  REF SKYD1240 4.0


Ø 4.3  REF SKYD1245 4.5


Ø 5.0  REF SKYD1255 5.5

Drill for medium hard and soft bone rpm 300


Ø 3.06  REF SKYD3435 3.5 3.5 N


Ø 3.56  REF SKYD3440 4.0


Ø 4.06  REF SKYD3445 4.5


Ø 4.76  REF SKYD3455 5.5


Crestal drill rpm 300

Ø 3.5  REF SKYCD35n 3.5 N narrowSKY


Ø 4.0  REF SKYXCD35 3.5 blueSKY, SKY classic

Ø 4.0  REF SKYXCD40 4.0 blueSKY SKY classic

Ø 4.5  REF SKYXCD45 4.5 blueSKY SKY classic


Ø 5.2  REF SKYXCD55 5.5 blueSKY

SKY Drill extension

 REF SKY-DV12

Dimensions in mm

Surgical tools



SKY TK mounter for ratchet short
REF SKY-STK5

SKY TK mounter for ratchet long
REF SKY-STK6

SKY TK mounter for contra-angle short
REF SKY-WTK5

SKY TK mounter for contra-angle long
REF SKY-WTK6

miniSKY insertion instrument for contra-angle short
REF mSKYXWM6

miniSKY insertion instrument for contra-angle long
REF mSKYXWM7


SKY Connector for contra-angle
REF SKYTWCN

Parallel indicator with conical and cylindrical side,
thicker central area with hole for protection against
accidental dropping
REF SKY-PI22

SKY fast & fixed angulation aid set 35°
REF SKYFFS35

whiteSKY mounter for ratchet
REF SKYC-SM6

whiteSKY mounter for ratchet
REF SKYC-WM6



blueSKY / narrowSKY X-ray-templates

Scale 1:1 REF bSKYMS01

Scale 1.12:1 REF bSKYMS12

Scale 1.26:1 REF bSKYMS26


SKY classic X-ray-templates

Scale 1:1 REF kSKYMS01

Scale 1.12:1 REF kSKYMS12

Scale 1.26:1 REF kSKYMS26

Prosthetic tools



SKY prosthetic key short
REF SKY-SD16

SKY prosthetic key long
REF SKY-SD25

SKY Connector for contra-angle
REF SKYTWCON

SKY prosthetic key for contra-angle short
REF SKY-SD22

SKY prosthetic key for contra-angle long
REF SKY-SD28

Screwdriver Alle 0.9 for transversal screw-retention
REF 310W0106

SKY Locator® mounter for contra-angle
REF LOCZWED6

SKY Locator® Instrument
REF LOCZINST



SKY Torque Wrench pro
Precise display of the torques from 10 to 45 Ncm
REF SKYTWPRO

SKY Laboratory handle incl. SD-22

- Work end for insertion of contra-angle handpieces
 - Work end for ball head screw (corresponds to SKY-SD21)
- REF SKY-SD80

SKY Universal forceps

- titanium nitrite-coated grip surface
 - Holding of implants and abutments
 - Oral securing of the prosthetic key
- REF SKY-SD60

SKY Key holder

- Oral securing of the prosthetic key
- REF SKY-SD65


Practice



With the SKY prosthetic case for practice, the necessary wrench is always to hand in your practice.

The prosthetic case is equipped with the Torque-Wrench pro and the long and short SKY prosthetic keys as standard.

Laboratories



With the SKY prosthetic case for laboratory, the necessary wrench is always to hand in your laboratory.

The prosthetic case is equipped with the laboratory handle and the prosthetic key for the contra-angle handpiece as standard.

SKY prosthetic assortment practice REF SKYPET10




- Content:**
- SKY torque ratchet
 - SKY prosthetic key short
 - SKY prosthetic key long

The following products can be ordered separately as desired.



- SKY prosthetic key for ball head attachments
REF SKY-SD21
- SKY Connector for contra-angle
REF SKYTWCN
- Screwdriver Alle 0.9 for transversal screw-retention
SKY 310W0106
- SKY Key holder
REF SKY-SD65

SKY prosthetic assortment laboratories REF SKYPET20



- Content:**
- SKY laboratory handle
 - SKY prosthetic key for contra-angle

The following products can be ordered separately as desired.




- SKY Locator® mounter for contra-angles
REF LOCZWED6
- Screwdriver, Allen screw 0.9 for transversal screw retention
SKY 310W0106

Accessories and instruments for the SKY implant lines


Torque Wrench pro



SKY Torque Wrench pro set
incl. SKY connector
Precise display of the torques from 10 to 45 Ncm
REF SKYTWSET



- Gingiva former and impression abutment (10 Ncm)
- SKY fast & fixed / uni.cone copings (18 Ncm)
- All SKY abutments (25 Ncm)
- Range for primary stability for immediate restoration 30 – 45 Ncm (40 Ncm for improved orientation)



- Easy to clean:
- The head separates easily from the handle using finger pressure
 - Easy to reassemble after cleaning
 - Done



- SKY connector
- For contra-angle handpiece instruments
 - Snaps firmly into the ratchet by pushing with your thumb
 - Easy to remove by pushing with your thumb

CPS Cordless Prosthodontic Screwdriver

How much time do you spend loosening and fixing implant screws?

You can save app. 50 % of the time using the CPS.

Are you sure that your screws are always the correct tightness?

With CPS there will be no screw loosening due to wrong torque.

Can you reach your palatal screws easily?

With CPS you have perfect access to all sites because it is cordless.

Perfect Handling

- cordless
- ergonomic design

Highest precision

- Precision of torque values +/- 5%
- Torque range: 8 Ncm – 40 Ncm

Best hygiene

- Non-slip handpiece and removable switch
- autoclavable
- thermo-disinfectable

Universal

With screwdrivers for the most important implant systems



CPS Cordless Prosthodontic Screwdriver
+ Full Range Driver-Kit
REF 580 CPS 40


Full Range Driver-Kit

Hexagon (Inch)	0.03	0.05			
Hexagon (mm)	0.9	1.0	1.2	1.8	2.5
Torx (mm)	SKY	5.5	6.0		
Slot (mm)	1.6	2.0			

Accessories and instruments for the SKY implant lines

Screws

In two posters presented at the SKY Meeting 2012, Dr. Wentaschek from the University of Mainz illustrated how Pretension can be significantly reduced by repeatedly tightening the screw. Pretension is the decisive crucial factor for the implant-abutment connection. It is therefore recommended that the appropriate laboratory screws are used in the laboratory, in order to prevent screw loosening. These results can be transferred to the recall. In this case, we recommend replacement of the screws when the abutment is removed for cleaning.



Source: S. Wentaschek et al: Reibungskoeffizient und Vorspannung beim Implantat-Abutment-Schraubenverbindungen [Coefficient of friction and pre-tension in implant abutment screw connections]; Scientific Book SKY Meeting 2012; S. 54-55; bredent medical GmbH & Co. KG, ISBN 978-3-00-038740-1

Repeatedly tightening the screws reduces pretension by up to 40%.



REF	SKY-PS22	SKYLPS22	SKYFFSPK	SKYFFLPK	SKYUFTS9
Description	SKY screw 2.2	SKY lab screw 2.2	SKY fast & fixed / uni.cone screw M1.4 für Proesthetic coping	SKY fast & fixed / uni.cone screw M1.4 grey	SKY fast & fixed / uni.cone Transversal screw Hexagon 0.9
pieces	6	10	6	10	6
Screw length / mm	9.2	9.2	5.3	5.3	2.5
Thread	M 1.8	M 1.8	M 1.4	M1.4	M 2x0.25
Head Ø / mm	2.2	2.2	2.2	2.2	2.2
Material	Ti*	Ti*	Ti*	Ti*	Ti*
SKY prosthetic key	✓	✓	✓	✓	✓



REF	m2SKYS22	m2SKYL22
Description	mini ² SKY screw	mini ² SKY lab screw
pieces	6	10
Screw length / mm	4.7	4.7
Thread	M 1.6	M 1.6
Head Ø / mm	2.2	2.2
Material	Titan Grade V	Titan Grade V
SKY prosthetic key	✓	✓



REF	COPASM16	COPALM16
Description	copaSKY screw	copaSKY lab screw
pieces	1	1
Screw length / mm	6.5	6.5
Thread	M 1.6	M 1.6
Head Ø / mm	2.2	2.2
Material	Ti*	Ti*
SKY prosthetic key	✓	✓

Ti*=
Grade 4 KV titanium



Torques for SKY prosthetics

SKY cover screw SKY gingiva former	25 rpm	10 Ncm	
SKY impression abutments	25 rpm	10 Ncm	
SKY temp SKY uni.cone SKY fast & fixed Prosthetic screw	25 rpm	18 Ncm	
SKY abutments	25 rpm	25 Ncm	 <p style="text-align: center;">Image on a very small scale</p>

Torques for miniSKY prosthetics

miniSKY Gingiva former impression abutments Scan abutment	25 rpm	10 Ncm	
miniSKY abutments	25 rpm	20 Ncm	







Torques for copaSKY prosthetics


copaSKY Gingiva former impression abutments	25 rpm	10 Ncm	
copaSKY abutments	25 rpm	25 Ncm	

HELBO® Treatment – modern treatment


With HELBO® Treatment you can quickly relieve your patients from these inflammations/infections or prevent disturbances to wound healing. The light also accelerates the healing process and has a proven analgesic effect. The treatment integrates exceptionally into existing treatments offered by brent medical. HELBO® therefore ensures treatment success in the context of immediate restoration of whole jaws with SKY fast & fixed or even of individual tooth gaps with SKY elegance. HELBO® can also be combined perfectly with bone replacement materials, membranes or even collagen fleeces in the case of augmentative measures and socket preservation. The HELBO® treatment can be assigned to adequately trained assistants to reduce the workload for the dentist.

Singlet oxygen destroys pathogenic bacteria. The proposed treatment is based upon the marking of the bacterial wall using light-sensitive dye molecules, which diffuse into the biofilm from the Photosensitizer. The dye molecules are then activated using laser light and transfer their energy into local oxygen. This creates highly-aggressive singlet oxygen, which destroys more than 99 percent of the bacteria in the biofilm. An efficacy rate as high as this would be unthinkable using conventional treatment methods.


-  Planktonic bacteria
-  Bacteria actively growing in the biofilm
-  Persisting bacteria in the biofilm
-  Photosensitizer (light-sensitive dye)
-  Biofilm matrix
-  Quorum Sensing signal molecules




How it works




1 Accumulation of light-sensitive photosensitizers on the bacterial membrane



2 Exposure and stimulation of the photosensitizers with the HELBO®TheraLite Laser




3 Reaction with oxygen, image of aggressive singlet oxygen



4 Damage to bacterial membrane: destruction of the micro-organisms

Areas of application




Periodontitis / Periimplantitis
(closed/surgical procedure)

HELBO® therapy is used in initial and maintenance treatment. Depending on the severity of the disease and accessibility, the procedure is carried out by means of a closed or surgical procedure. Studies confirm that the inflammation subsides, healing is promoted and the probe depth is reduced with this treatment.



Image:
Dr. Tilman Eberhard, Schwäbisch Gmünd




Disinfection of the
alveoli

Socket preservation

Teeth are most often removed due to bacterial inflammations. In terms of aesthetics and additional measures such as immediate or even delayed implantation, it is particularly important to disinfect the extraction alveolus accordingly. Collapsing of the vestibular bone lamella can be avoided when using together with suitable materials (e.g. TIXXU®PROTECT collagen fleece).



Image:
Dr. Jörg Neugebauer, Landsberg am Lech




Disinfection of the
alveoli

Immediate implantation

During immediate implantation, HELBO® therapy, as a disinfecting measure following tooth extraction, can visibly reduce the occurrence of complications for single tooth restorations (e.g. with SKY elegance) and for full arch restorations (e.g. with SKY fast & fixed).



Image:
Dr. Jörg Neugebauer, Landsberg am Lech



Disinfection

Augmentation


Chronic inflammations frequently lead to bone loss. Before augmentative measures are carried out, the hard and soft tissue infected with bacteria is decontaminated using HELBO®. The bone defect is then filled with bone replacement material (such as TIXXU® GRAFT) and covered with a membrane (e.g. TIXXU® CONTROL). By using this procedure, the risk of inflammation is reduced and controlled regeneration is enabled.



Image: Dr. Sigurd Hafner, Munich



Image: Dr. Sigurd Hafner, Munich



Bone necrosis

Intraoperative staining of the biofilm with sterile *Photosensitizer* enables contamination on bones and soft tissue to be visualised and the area to be disinfected using laser radiation. This results in improved wound healing and ensures the success of the surgical measures.



Image:
Dr. Sigurd Hafner, Munich

visio.lign®

The aesthetic and functional system



<p>Coating</p>   <p>crea.lign®</p>	<p>Milling</p>   <p>visio.CAM</p>
<p>Veneering</p>   <p>novo.lign®</p>	<p>Setting up</p>   <p>neo.lign®</p>
<p>Painting</p>   <p>visio.paint</p>	<p>Bonding</p>   <p>bond.lign</p>

Global market leader* in physiological veneers

2013 visio.CAM
L-Design (novo.lign® & neo.lign®)

2014 crea.lign® Paste

2016 bond.lign
crea.lign® Paste GUM

2017

* 6800 users
1,2 Mio. restorations

Chipping repair kit



Following repair with crea.lign® composite

Instances of chipping are restored easily and quickly

- Primer works together with all prosthetic materials
- crea.lign® incisal masses and dentine (in A2 and A3 for covering in 80 % of all cases)
- 100 treatments with a single chipping repair kit
- No hydrofluoric acid necessary

Chipping Repair Kit for Dentists REF CLIGNSETA

contains crea.lign® composite in colours A2/A3 including 4 primers and consumables.



Intraoral application

- "mouth approved": approved for intraoral application
- Chipping repair for all veneer materials
- "What you see is what you get": result is immediately visible and adjustable - no kilning required
- Rapid repair in situ: direct and simple

Permanent repairs with no loss of quality

- Material is approved for permanent repairs
- No loss of aesthetics: no colour differences and no "repair transitions" can be detected
- No loss of mechanical properties
- The result is at "supra-ceramic level"

Images of intraoral, permanent repairs
© 2016 by MDT Andreas Lüdtkke, Bayreuth, Germany

Full Range Bonding Kit





The neat solution for all prosthetic materials

Simple and safe cleaning is ensured by the application of FGP insulation varnish on the parts to be bonded. Areas of inflammation due to adhesive residues are a thing of the past.

Primer for all prosthetic materials

- Titanium, zirconium, NPM – MKZ primer
- Ceramic, lithium disilicate – K primer
- PMMA, composite, BioHPP – visio.link




DTK adhesive


- Opaque (colour A2 / A3) for bonding titanium bases with ceramic abutments or crown abutments – no dark shadows impair the result
- Transparent for bonding ceramic abutments and crowns – no impairment of the translucence.

Implants


narrowSKY

	3.5 N	L 8 mm	L 10 mm	L 12 mm	L 14 mm	L 16 mm
		nSKY3508	nSKY3510	nSKY3512	nSKY3514	nSKY3516


blueSKY

	4.0	L 8 mm	L 10 mm	L 12 mm	L 14 mm	L 16 mm
		BSKY4008	BSKY4010	BSKY4012	BSKY4014	BSKY4016
	4.5	BSKY4508	BSKY4510	BSKY4512	BSKY4514	—
	5.5	BSKY5508	BSKY5510	BSKY5512	—	—


SKY classic

	4.0	L 8 mm	L 10 mm	L 12 mm	L 14 mm	L 16 mm
		KSKY4008	KSKY4010	KSKY4012	KSKY4014	KSKY4016
	4.5	KSKY4508	KSKY4510	KSKY4512	KSKY4514	—

miniSKY

	2.8	L 6 mm	L 8 mm	L 10 mm	L 12 mm	L 14 mm
		—	—	m2SKYL10	m2SKYL12	m2SKYL14
	3.2	—	m2SK3208	m2SK3210	m2SK3212	m2SK3214
	2.8	m1SKYL06	—	m1SKYL10	m1SKYL12	m1SKYL14

whiteSKY

	3.5	L 8 mm	L 10 mm	L 12 mm	L 14 mm	L 16 mm
		—	SKY3510C	SKY3512C	SKY3514C	SKY3516C
	4.0	SKY4008C	SKY4010C	SKY4012C	SKY4014C	SKY4016C
	4.5	SKY4508C	SKY4510C	SKY4512C	SKY4514C	—

whiteSKY Prosthetic coping

3.5		4.5	
4.0	SKYCPK40		SKYCPK45

Stamp

Date Signature

Illustrations not true to scale. Subject to errors and changes. All prices in € incl. VAT, plus shipping costs, if applicable. Offers cannot be used in combination with existing offers. This order is based on our terms and conditions which can be accessed at <http://www.bredent.com/en/brand/content/terms/>

Surgical and prosthetic instruments

Bone bur / Drill extension / Parallel indicator



Pilot drill



Twistdrill



Drill for hard bone



Drill for medium hard and soft bone




Crestal drill




Depth stop

	L 6 mm	L 8 mm	L 10 mm	L 12 mm	L 14 mm	L 16 mm
	SKYXST06	SKYXST08	SKYXST10	SKYXST12	SKYXST14	SKYXST16
	—	SKYS0840	SKYS1040	SKYS1240	SKYS1440	SKYS1640
	—	SKYS0845	SKYS1045	SKYS1245	SKYS1445	—

Implant mounting tool SKY



miniSKY



whiteSKY



SKY prosthetic key




SKY Torque Wrench pro



Patient passport / X-ray-templates

	Scale	1:1	1.12:1	1.26:1
blueSKY/narrowSKY		bSKYMS01	bSKYMS12	bSKYMS26
SKY classic		kSKYMS01	kSKYMS12	kSKYMS26
miniSKY		mSKYMS01	mSKYMS12	mSKYMS26
whiteSKY		SKYMS01C	SKYMS12C	SKYMS26C



Stamp

Date Signature

Illustrations not true to scale. Subject to errors and changes. All prices in € incl. VAT, plus shipping costs, if applicable. Offers cannot be used in combination with existing offers. This order is based on our terms and conditions which can be accessed at <http://www.bredent.com/en/bredent/content/terms/>

Prosthetics – Immediate and late restoration

For narrowSKY approved prosthetic components rose gold anodised

Prosthetic components not suitable or not approved for narrowSKY

BioHPP SKY elegance Immediate restoration




BioHPP SKY elegance titanium base




BioHPP SKY elegance prefab set



SKY uni.cone



SKY fast & fixed



Angulation aid set 35°



SKY locking pin



Individual solutions for CAD/CAM – conventional

SKY uni.fit CAD/CAM solutions / individual abutments



SKY prefab titanium set



SKY castable abutment



Stamp

Date Signature

Illustrations not true to scale. Subject to errors and changes. All prices in € incl. VAT, plus shipping costs, if applicable. Offers cannot be used in combination with existing offers. This order is based on our terms and conditions which can be accessed at <http://www.bredent.com/en/bredent/content/terms/>

Classic implantology



Implant analog



SKY temp




SKY impression abutments



SKY esthetic gingiva former



SKY esthetic abutments



SKY Titanium abutments



SKY prosthesis fixation


SKY Locator



SKY Locator anglet



SKY Locator



TiSi.snap



TiSi.snap anglet



Stamp

Date Signature


Illustrations not true to scale. Subject to errors and changes. All prices in € incl. VAT, plus shipping costs, if applicable. Offers cannot be used in combination with existing offers. This order is based on our terms and conditions which can be accessed at <http://www.bredent.com/en/bredent/content/terms/>

Always keep your patients well informed

Request patient information for our treatments now free of charge. Use it to always keep your patients up to date.



Patient brochures
Smile again
REF 000 540 GB



Patient brochures
HELBO®
REF 000 484 GB

Treatments explained simply

Watch our treatments in a video on YouTube. Tell your patients on your website or in the waiting room. Request the videos from marketing-medical@bredent.com.




Late implantation and immediate restoration



Delayed implantation and immediate restoration



Immediate implantation and immediate restoration



SKY fast & fixed - 3D animation



The conventional bridge - and then...

This QR code takes you straight to the videos



Link to the videos
www.youtube.com/bredentgroup

Demonstration model

Request SKY fast & fixed demonstration model and use it to explain the treatment to your patients.



Cross-section model
SKY fast & fixed
with nerve course, gingiva mask
and temporary bridge
REF 590 ZMFF 2



Posters advertising our treatments

Request the print data for your posters in the A1 format free of charge now and recommend the therapies of bredent medical using a light-hearted approach in your practice.

Send your request to:
marketing-medical@bredent.com



...gapless...
 REF OPO 201G B



...gapless... Version 2
 REF OPO 202G B



...in the thick of it...
 REF OPO 205G B




...in the thick of it... Version 2
 REF OPO 204G B




tempting...?
 REF OPO 203G B



tempting...? Version 2
 REF OPO 200G B



...bacterial infection
 REF OPO 210G B



...bacterial infection Version 2
 REF OPO 220G B

The Compendium

A variety of practical cases, both scientifically and clinically documented, complete with illustrations. Discover new approaches to restorative treatments and take away ideas that you can use in your own practice or laboratory. Interesting and exciting contents in all four Scientific & Clinical Cases, complete with university prefaces.




available in
German REF 992 975 0D
English REF 992 975 GB



available in
German REF 992 977 0D
English REF 992 977 GB



available in
German REF 992 976 0D
English REF 992 976 GB



available in
German REF 992 978 0D
English REF 992 978 GB

- Immediate restorations – SKY® fast & fixed therapy
- Physiological Prosthetic
- Immediate single-tooth restoration
- Regeneration and Biofilm management


Get one copy for free*


Send an email to marketing-medical@bredent.com, mention the copy of your choice, your contact details and the benefit code "MA459".

Literature presentation



Immediate restoration with reduced number of implants


 REF 992 9710 D


 REF 992 971 GB

Authors: Dr. Georg Bayer
Dr. Frank Kistler
Dr. Steffen Kistler
Stephan Adler
Dr. Jörg Neugebauer



The guide to modern implant prosthetics


 REF 992 9700 D


 REF 992 970 GB

Authors: Dr. Manfred Lang



The implantation simulator


 REF 992 9690 D


 REF 992 969 GB

Authors: Dr. Manfred Lang



Sinuslift


 REF 992 9680 D


 REF 992 968 GB

Authors: Dr. Manfred Lang



Scientific Book 2012

 REF 992 9740 D

 REF 992 974 GB




www.bredent.com/literatur
On this page you will find a bibliography.

Other offers that may be of interest



Tissue related management
REF 009912GB




Implant prosthetic
REF 009913GB



Immediate restoration
REF 000200GB



BioHPP - The reference
REF 000 547 GB



BioHPP elegance hybrid abutments
REF 000 534 GB



HELBO®-Treatment
REF 000 429 GB

Mistake and subject to change reserved 000250GB-20170306



40 YEARS DENTAL INNOVATIONS
1 9 7 4
2 0 1 4

bredent
medical