# Products for the dental-technical laboratory



English

## Dental innovations for a successful laboratory

## Dear Client,



Continually striving for perfection is part of human nature. We, at the bredent group, are committed to abiding by our values to achieve this and take responsibility for our conduct, both within the company and with external parties.

The intensive interaction with reference laboratories and practices has created an exceptionally large potential for product ideas, which are developed over a total operating area of 12,000 m<sup>2</sup> in Senden/Iller in Germany.

As a client, you understand the benefits of bredent products and thereby contribute to the success of the bredent group. As a family-run company, a relationship based on partnership is incredibly important to us. We would like to thank you personally for this and say that we remain available to you as a direct contact partner.

Do share your ideas on product optimisation and suggestions for an even closer cooperation by e-mailing us at peter.brehm@bredent.com or nils. brehm@bredent.com, or send us a fax on +49 (0) 7309/872-155. We look forward to hearing from you.

We hope you will be impressed when discovering the new - relaunched - bredent dental technology catalogue.

Best regards,

More than 39 years of dental innovations

Peter Brehm

Nils Brehm

## bredent group mission statement

## The bredent group

We, at the bredent group, are an internationally-active, family-run company that develops and produces optimally coordinated stand-alone products, system solutions (e.g. in implant prosthetics) and "Made in Germany" treatment concepts. These enable dentists and dental technicians to produce and maintain high-quality, cost-effective aesthetic restorations that ensure periodontal hygiene.

We strive to be amongst the best. That is why our employees are ready to deliver excellence for our clients and their patients, with the necessary flexibility and openness.

Our expertise and innovative drive make us the paradigm for the dental market.

We base our corporate culture, our collaboration with our clients and business partners and beneficial activities on our values that we set in 1995 as our foundation for relationships and long-term objectives.

### capable

We are convinced of the benefits of our work and willingly endeavour to become more efficient.

## cooperative

We are open and fair in our collaborative partnerships. This is how we build trust.

### innovative

Our expertise, flexibility and openness to the world allow us to create beneficial solutions at an early stage.

For our clients, this means that:

Our clients have a business partner that can offer benefits thanks to our performance and solutions. All of our activities demonstrate this aspiration. We are attentive to our customers and their requirements, which in turn allows us to keep up to date about any issues and wishes that affect patients. This information forms the basis for the development of new, inspiring and cost-effective solutions – which warrants the benefits of our work for clients and patients. For us, complying with the relevant regulations, standards and legislation is a matter of course.

For our organisation and our employees, this means that:

We see amendments and new requirements as opportunities and as something that challenges us to provide a rapid and flexible response. A capacity for change is therefore a significant prerequisite for our success. Consequently, we are continuously making our organisation quicker and more cost-effective in order to meet the increasing demands of our clients, their patients and of the economic environment.

For our employees, this means that:

We demonstrate engagement and a sense of responsibility in the relevant task, as individuals and as a team. Everyone takes responsibility for their involvement. Everyone is committed to achieving the company's objectives and maintaining the company's values. In return, our employees enjoy a safe place to work, personal development opportunities and scope for professional development and improvement. By using and expanding our collective knowledge, we are laying the foundations for future success.

As far as the way in which we interact with society is concerned, this means that:

We believe that businesses have a responsibility to make a lasting positive social and environmental contribution to society.


Our income serves job security and is a basis for further pioneering investments in line with our goals and values for the benefit of our clients.

The Management *Peter Brehm, Brigitte Brehm* 



## Perfectly-fitting prostheses – a fundamental prerequisite for satisfaction

## The bredent symbiosis







## Communication as a central component on a client relationship

## What is the most convenient way to reach us?

To suit your requirements, we have over 350 employees available, working in research, development, production, administration and sales. Each and every one of them helps contribute to making your day-to-day work more productive. You are the incentive for us to continue to evolve and learn - for your and our benefit.

A team of nearly 100 trained medical device advisers are on hand across the country and internationally as part of our field staff. Having someone in your area ensures excellent communication and fast access for you.

Speak to your personal adviser and benefit from their expertise on our products, system solutions and treatment concepts for a successful laboratory.



You can reach us personally:

Monday to Friday 7.00 a.m. to 4.15 p.m. Telephone: +49) 0 73 09 / 8 72-4 40

and 24 hours a day:

Fax: (+49) 0 73 09 / 8 72-4 44

E-mail: info@bredent.com Internet: www.bredent.com

Orders received by 4 p.m. will be dispatched on the same day.



Not all products in the current catalogue have been approved or are available in all markets. If you have any questions, please do not hesitate to get in touch with your to contact partner:



# Application areas

		Page
Systems	1 Preparatory work	7 - 46
	2 Prosthetic planning / 3D implant planning	47 - 60
	3 Temporary restoration	61 - 74
	4 CAD/CAM	75 – 82
	5 Framework manufacture	83 - 130
	6 Attachments, Locks and Screw connections	131 - 232
	7 Precision casting by bredent	233 - 268
	8 Processing of ceramic and zirconium	269 - 286
	9 visio.lign veneering system / Teeth / Denture resins	287 – 348
	10 Plastic moulding technology / thermoplastic resins	349 - 362
	11 Snoring therapy	363 - 368
	12 Epithetics	369 - 374
Product groups	A Devices	375 – 390
	B Instruments	391 - 410
	C Rotating tools	411 - 470
	D Polishing brushes / Polishing pastes	471 – 488
	Index	489



**Preparatory work** 

Preparatory work represents the beginning of the fabricatio process in the laboratory. The working model is the basis of the prosthetic restoration. High precision needs to be ensured, which is accomplished through the use of high-quality, matched products. This way errors in the process chain are reduced and the daily work processes are facilitated.





he impression	
breciform D impression trays	8
brecision impression materials	8
Abdruck-Cut	
Dentaclean impression desinfectant	9
Shipping bags	9
Disinfection bath 3L	
Silicone and wax surface tension	10
Surface tension reducing agent	10
he individual tray	
System steps	11
Transblock	
Isoplast ip	
Tray material UV	
Polylux 2	
Tungsten carbide burs for processing acrylics	14
he sawcut model	
System steps	
ecovac vacuum mixing system	
Exakto-Rock S	
Master-Pin Radix-S	
Master-Pin Radix-K	
Master-Sep	
Master pin drill unit mpb 1	
Master-Pin System	22
Master-Split model systemFluid-Rock	
Arti-Rock	
Thermo-syringe	
Plaster knife	
Plaster removing agent	20
Speed plaster removing agent	29
Giflex-TR	
Tungsten carbide burs for processing of plaster	
Litebloc UV	
Undercut wax	
Die varnishes	32
diephos dentine	33
Spacer varnish gold, silver, silver-blue, blue	
Gloss and hardening agent for plasters	34
he master model	
System steps	35
Functional margin wax	
Retention pins	
Pi-Ku-Plast / Pi-Ku-Plast HP 36	37
ecovac vacuum mixing system	38
Exakto-Rock S	38
Exakto-Form	
Master-Split model system	40
he implant model	
System steps	
Multisil-Mask soft	
Multisil-Mask hard	
haptosil D	
KoEx Measuring Device	45

## breciform D impression tray

### The proper impression technique to obtain exact models is ensured by the breciform D impression tray.



Special tongue support reduces the patient's gag reflex and hence minimizes sources of errors during impression taking. Essential anatomical requirements were taken into account in the development of the impression tray. With the help of kneading silicone and click-type stops as spacers, breciform D can be easily and quickly adapted and an additional patient session can be avoided.

## **breciform D impression tray** - for single use **Starter set**

10 upper and lower trays each sizes S, M, L and XL 10 breciform D triangular stops 10 breciform D bar-shaped stops REF 580 UOTS S

## brecision impression materials

The silicone-based impression materials that reproduce details very accurately offer ideal preconditions for perfect impression taking thanks to various flow characteristics.



### brecision putty soft

With Shore A hardness of 70, brecision Putty soft is a non-sticky and kneadable base impression material based on addition-curing vinyl polysiloxanes. brecision Putty soft can be easily mixed and perfectly cut and adjusted.



250 ml base (grey), 250 ml catalyst (white), 2 measuring spoons Set 4 pieces REF 580 0002 4



#### brecision implant heavy

The shore A hardness of 70, accurate reproduction of details, reduced flowability and medium hydrophilicity are preconditions for precise impressions of the situation and fabrication of the dental restoration.

### brecision implant heavy

Impression material blue, 1 x 380 ml 5 x dynamic mixers, 1 x bayonet ring blue RFF 580 BH38 0



### brecision implant light

The very high hydrophilicity and flowability, very accurate reproduction of details and a shore A hardness of 55 ensure precise impressions of delicate implants or teeth.

### brecision implant light

Impression material orange, 2 x 50 ml cartridges 10 mixing cannulas, 10 Intra-oral tips REF 580 BL05 0  $\,$ 



### security-bite blue

The thixotropic property of the security-bite blue bite registration material based on A-silicone can be adapted to the arch without dripping. The high Shore A hardness of 90 enables dimensionally stable bite registration and reproduction of all details.

### security-bite blue

2 x 50 ml cartridges (blue) 12 contouring tips (wide) 12 mixing cannulas (pink) REF 580 0002 0

### Abdruck-Cut

Undercuts can be easily and specifically removed using the scalpel-sharp loop blade.



Abdruck-Cut 1 Stück REF 360 0114 0



The scalpel-sharp loop blade allows cutting even in areas difficult to access.

### Accessories



Loop knife 1 piece REF 360 0115 0



## Dentaclean impression and denture disinfectant

Disinfecting with Dentaclean impression and denture disinfectant avoids the transmission of viruses, bacteria and fungi – from the patient to the laboratory and hence increases the protection against infections for you.



The concentrate is mixed to obtain 10 liters of ready-to use solution which is highly effective and has a surprisingly mild odor.

### Dentaclean impression and denture disinfectant 1000 ml concentrate to obtain 10 liters ready-to-use solution incl. 25 shipping bags

REF 520 0100 6

Tested and approved by the Institute for clinical hygiene and infection control, Giessen.



Pathogens can be transmitted to the laboratory with impressions.



After the use of Dentaclean impression disinfectant, active viruses, bacteria and fungi can no longer be detected.

## Shipping bags

The shipping bags have already been labeled "disinfected".



Additionally, a separate bag holds the refte to protect them against moisture.

Shipping bags 200 pieces REF 520 0100 2

### Disinfection bath 3L

The convenient filter basin enables disinfection of up 6 impressions thanks to the flat and wide design of the disinfection bath 3L. As a result, the efficiency is increased and waiting times are reduced.



**Disinfection bath 3L**W 35 x D 26 x H 14 cm
1 piece **REF 230 0015 0** 

- The bredent disinfection bath 3L has a capacity of 3 liters
- Due to the convenient filter basin careful cleaning of impression trays and instruments is simplified
- Direct skin contact with the disinfectant is avoided thanks to the integrated dripping device
- Individualizing of the filter is possible due to moveable instrument rests



Dripping device prevents direct skin contact with the solution. This guarantees safe every day usage.



The stable bath made of polypropylene is heat resistant up to 135° C and therefore suitable for autoclave and thermodisinfecting.



## Silicone and wax surface tension reducing agent

Enhances the flow properties of plaster for silicone impressions.



Spraying on the silicone and wax tension reducing agent improves the flow properties of plaster for silicone impressions. Before pouring the arch, the impression must be dry.

Silicone and wax surface tension reducing agent 750 ml REF 540 0070 5



The spraying head of the spray bottle simplifies uniform wetting of the surface with silicone and wax surface tension reducing agent.



After the application of the agent onto the surface (left), the flow characteristics of the plaster have been clearly improved.



Silicone and wax surface tension reducing agent produces a homogeneous plaster surface. This will ensure precise dental

## Surface tension reducing agent

The surface tension reducing agent for impressions. Cleans, disinfects and improves the flow characteristics of model materials.



The use for silicone, alginate and hydrocolloid impressions reduces stockkeeping.





Refill package 750 ml REF 520 ES75 0



The fine spray head of the plastic spray bottle simplifies uniform spraying of the liquid.

Accessories for silicone, wax and surface tension reducing agents

Spray bottle, plastic sp 1 piece, 125 ml REF 540 0075 0





Spray on a thin coat of surface tension reducing agent. Allow to react for 1 to 2 minutes for alginate and hydrocolloid impressions. The blow the impression dry and cast. Condensation cured silicone impressions: the impression is rinsed with water after the reaction time and blown dry subsequently. The plaster flows without any formation of bubbles and surface segregation.

## The individual tray



The initial situation The model fabricated using a class III stone.

An accurate impression tray made of stable and deformation-free resin is required for high-precision anatomic impressions, which is ensured by the lightcuring tray material from bredent.



**Blocking out** Simple and fast blocking out of remaining teeth or undercuts with transparent Transblock.



Insulating
The model is insulated with Isoplast ip to avoid adhesion of the tray materials. As a result, the hardened tray can be easily removed.



Tray material UV

The high stability ensures fast working in a way to



### Polymerizing

The fully mirrored inner chamber of Polylux 2 offers perfect illumination and hence reliable curing of all areas of the tray material UV.



Wide range of tungsten carbide and diamond tools for fast processing of the tray material.

### Materials used

Transblock	page 12
Isoplast ip	page 12
Tray material UV	page 13
Polylux 2	page 14
Processing tools	page 14



## The individual tray

### **Transblock**

The transparent block-out material for fast and systematic working.



The stability of Transblock results in uniform layer thicknesses and can be adjusted individually by scraping.



Any desired size or shape of Transblock can be produced with the help of an instrument or scissors.

Due to its stability a uni-

during the adaptation. If

can be adapted individu-

required, the thickness

ally by scraping.

form thickness is retained



Transblock 250 g REF 540 0114 9



The high flexibility simplifies placing onto the model.



The transparency of Transblock allows to check the thickness of the area that has been blocked out. This way precisely prepared models for individual trays are obtained.

## Isoplast ip


Isoplast ip is alginate based and insulates plaster against resin whilst creating a highly lustrous resin surface.







The brush pen allows to apply the material in an economic and precise way.





Isoplast ip seals the surface and the plaster exhibits a fine luster. This way the quality of the insulating layer can be checked.

Isoplast ip 750 ml REF 540 0101 9



Brush pen pk 125 REF 390 0033 0



Brush pen pk 20 REF 540 0072 0



## Tray material UV

High-stability resin to obtain individual trays for accurate impressions.



The flexibility of the material allows easy placement onto the model without tearing. The required shape can be cut with an instrument. The pink color provides the perfect basis for the set-up. The high stability allow quick placement of the tray handle without any dimensional change until polymerization is completed.



The high flexibility of the material simplifies the placement onto the model. The material will not be damaged.



Perfect adaptation to any situation guarantees uniform wall thicknesses.



Due to the high stability the position of the handle which has been determined will not be changed during the polymerization process.

The tray material UV can

be precisely cut with any

instrument. Accordingly,

the amount of work is

reduced.



The tray material UV has hardened after only 10 minutes in the Polylux 2 unit.



The high stability of tray material UV avoids deformation during impression taking. Precise models will be obtained.



Tray material UV 50 pieces UJ REF 540 0011 0



Tray material UV band 2,5 mm x 90 mm 1350 g REF 540 0016 6

Tray material UV block 1000 g REF 540 0011 3



Tray material UV 50 pieces ⊔ REF 540 0011 1

Accessories:



Polylux 2 polymerization unit with material container REF 140 0099 0

Assortment 25 Tray material UV UJ 25 Tray material UV LJ REF 540 0011 2

## The individual tray

## Polylux 2

Powerful and universally suitable light-curing unit for dental consumables with a wavelength range of 350 to 500 Nm.



Polylux 2 with drawer, 230 V Polylux 2 with drawer, 115 V

REF 140 0099 0 REF 140 0099 1

Powerful and universally suitable light-curing unit for dental consumables with a wavelength range of 350 to 500 Nm. Two different and energy-saving special lamps guarantee excellent.

- convenient and simple operation thanks to clearly arranged buttons
- fully mirrored polymerization area for selective illumination with compact drawer
- device can be accessed from three sides to enable polymerization of large objects.
- reduced energy consumption, extended service life

### Technical data

Number of light sources Wavelength range Triple timer function Mains voltage

Dimensions of unit Dimensions interior - drawer Weight

2 fluorescent lamps, 9 watts each

350 - 500 Nm

180 sec/360 sec/continuous operation

230V, 50 Hz

approx. 250 x 120 x 90 mm. approx. 140 x 110 x 55mm.

approx. 1.5 kg.

## Tungsten carbide burs for processing acrylics

The proper selection of tools reduces the amount of work.



Diatit bur 1 piece **REF D468 GG 23** 

Diatit bur

1 piece **REF D468 GG 16** 



The triple cutting tungsten carbide burs are perfectly suitable to cut off excess tray material UV. In case of shellac the shape of he Tungsten carbide bur avoids loading of the cutting edges.



Tungsten carbide 1 piece **REF H194 SH 70** 




Aggressive cutting of the super-coarse crosscutting edge allows rough shaping in a very short time.



Tungsten carbide 1 piece **REF H274 GH 60** 



The medium-coarse cross-cutting edge smoothens the surface and allows finishing in a single working step.





Diacryl grinders 1 piece REF 340 0102 0



The margin cutter produces a uniform tray margin and creates sufficient space for lip and cheek fraenums.





#### The arch

Highly accurate fabrication of the arch with Exakto-Rock S and ecovac. Marginal stability and accurate reproduction of details reduce the amount of work.

The basis for fixed restorations is created by selecting the suitable materials

Depending on further processing – digital or conventional – Exakto-Rock S super-hard stone is used. The design of the high-precision pins facilitate the removal of the model.



### The Master-Pin System

The specially shaped inner geometry of the plastic sleeves and the metal pins ensure perfect fit to facilitate the handling and to provide optimized fit of the working dies and guarantee stress-free working.



#### The model system

The Split-Cast model system offers added comfort and facilitates the assembly of models. The low-viscosity Fluid-Rock base stone provides safe retention of pins and can be poured without the formation of bubbles.



### Segmentation

The perforated design enables to achieve controlled sawcuts. The high-performance diamond discs available in various diameters allow fast and precise preparation of dies.



### Surface processing

The special cut types enables fast and trouble-free processing of plaster. Cracks in the plaster are avoided.



### The die varnish

bredent offers varnishes in various colors and thicknesses for surface hardening or as spaces for the cement gap. Tooth-colored varnishes for the veneering technique offer excellent contrast and hence facilitate processing.

### Materials used

ecovac vacuum mixing system	page 16
Exakto-Rock S	page 17
Master-Pin Radix-S	page 18
Master-Pin Radix-K	page 19
Master-Sep	page 19
Master pin drill unit mpb 1	page 20
Master-Pin System	page 22
Master-Split model system	page 25
Fluid-Rock	page 28
Arti-Rock	page 28
Thermo-syringe	page 29
Plaster knife	page 29
Plaster removing agent	
Speed plaster removing agent	page 29
Giflex-TR	page 30
Tungsten carbide burs for processing	
of plaster	page 31
Liteboc UV	page 31
Undercut wax	page 31
Die varnishes	page 32
diephos-dentine	page 33
Spacer varnish gold, silver, silver-blue, blue	page 33
Gloss and hardening agent for plasters	page 34

bredent

## ecovac vacuum mixing system



### ecovac

## Precision-fit restorations obtained through optimal use of material properties.

The user-friendly and compact design simplifies work and reduces sources of errors. A powerful and maintenance-free vacuum pump, adjustable in two different levels (15 mbars, 200 mbars), ensures bubble-free mixing of materials and results in a perfect casting surface. Stirring time and speed can be adjusted continuously to allow correct processing of different materials.

ecovac (230 V)

REF 140 0093 0

(Wall mounting, without mixing cup and base)

- 1 mains cable
- 1 spare filter
- 1 drilling template for wall mounting
- 4 screws and plugs for wall mounting

Accessories

Base ecovac, 1 piece REF 210 0045 0



## ecovac mixing spiral

The mixing spiral takes up the components to be mixed from all areas of the mixing cup and stirrs them horizontally and vertically. No unmixed materials will remain on the bottom of the mixing cup, which may cause different expansion of the material later on.

All features and components listed provide increased reliability, lead to improved fit when preparing dental restorations and avoid time-consuming reworking.

Mixing spiral,	50 ccm	REF 140 0R94 5
Mixing spiral,	250 ccm	REF 140 0R94 0
Mixing spiral,	750 ccm	REF 140 0R94 2
Mixing spiral.	1000 ccm	REF 140 OR94 3



### ecovac mixing cups

The smooth inner surface of the stainless steel mixing cup prevents any material or liquid residues from adhering to or depositing in scratches or undercuts. The conical shape ensures that material which has been taken up will flow back to the center of the mixing cup. Accordingly, the mixing ratio is retained exactly and better results can be achieved with minimal effort.

Mixing cup,	50 ccm	REF 140 0B94 5
Mixing cup,	250 ccm	REF 140 0B94 0
Mixing cup,	750 ccm	REF 140 0B94 2
Mixing cup,	1000 ccm	REF 140 0B94 3



Mixing cup, D
(for the use in the Degussa mixing unit),

REF 140 0B94 4

### Exakto-Rock S

Exakto-Rock S is a synthetic super-hard class IV stone with distinctive thixotropy and improved flow properties and is free from formaldehyde.



Reduced expansion is completed after 2 hours and is only 0.08 %. As a result, highly accurate impressions are enabled and precision-fit restorations are ensured. Moreover, thanks to optimized reflection of light, Exakto-Rock S is suitable for scanning and available in crown and ivory.

- fomaldehyde-free stone ensures safe processing and can be used for the fabrication of holistic dentures without any problem
- synthetic components ensure consistent quality and enable the fabrication of precision-fit models
- improved flow characteristics facilitate pouring of several impressions
- optimized reflection of light thanks to special dye pigments reduce reworking/adjustments in the CAD system



### Color brown:

1 x 2 kg REF 570 0SB5 2 5 x 2 kg REF 570 0SB5 1 10 x 2 kg REF 570 0SB5 0



### Color ivory:

1 x 2 kg REF 570 0SE5 2 5 x 2 kg REF 570 0SE5 1 10 x 2 kg REF 570 0SE5 0

### Expansion of various other stones



### Technical data Exakto-Rock S

Color brown, ivory Mixing ratio 100 g / 20 ml dist. water Soaking time 20 sec Manual mixing time 20 sec Vacuum mixing time 40-60 sec Processing time 5-6 min Setting time (Vicat time) approx. 10 min Removal of model after 40 min above 60 MPa Compressive strength after 1 hr Compressive strenght after 24 hrs 85 MPa Hardness after 1 hr (Brinell) 200 MPa Hardness after 24 hrs (Brinell) 280 MPa Linear expansion after 2 hrs < 0.08 % (no further expansion)

## 



The excellent processing time span allows bubble-free pouring of numerous impressions with just a single mix.



Exakto-Rock S offers high stability on the spatula and thixotropic consistency on the vibrator. Simple and clean processing is ensured.



Absolutely accurate reproduction of dimensions of the oral situation thanks to the minimal expansion value (< 0.08 %) so that precision-fit dentures are obtained.



The arches can be cut and trimmed without the formation of chips.



Preparation limits of the dies are not damaged when grinding the dies. No breaking of edges when removing the model. Consequently, precision-fit restorations are obtained.



Perfect surface detection in the scanner facilitates the design and offers the possbility to obtain high-quality, precision-fit restorations. Final expansion is completed after 2 hours and processing can be continued quickly.

### Accessories:



KoEx Measuring
Device
1 piece including
2 contraction inserts
REF 110 0148 0

More information and reference numbers on page 46.



## Master-Pin Radix-S

The sturdy, high-tech plastic provides the required stability and the root shape offers protection against twisting.



impression very easily and safely.




Master-Pin Radix-S can be easily and safely aligned and fixed in the impression.

thin plug-type pin for reduced displacement of impression



Pouring out and preparing the base of the impression are done in the usual way.

Optimized harpoon-shaped tip for safe hold in all impression materials (silicones, alginates etc.).



modified die-shaped

highly stable, extremely smooth high-tech plastic

the root-like shape of Radix-S forms the perfect counterpart in the base plaster an alveolar socket



### Master-Pin Radix-K

The favourably-priced dowel pin solution for the production of models.



Due to the special root shape only one pin can be used per die. In addition to the retention element, glueing surfaces are integrated to ensure safe hold in the die. The smooth surface of the high-tech plastic allows easy integration and removal of the die.



Perfect glueing is ensured when a small amount of adhesive is also applied to the area of the support.



The root shape that is obtained in the model base ensures exact guidance and positioning. The dies are protected against tilting and twisting.




It is also possible to place interdental Master-Pins Radix-K which are not glued in.



Master-Pin Radix-K 1000 pieces REF 360 0123 2







## Master-Sep

Assortment
250 Master-Pin Radix-K
1 Tungsten carbide
drill
REF 360 0123 4

Special plaster against plaster separating liquid with unsurpassed separating effect for sawcut models.



Master-Sep Special separating liquid for sawcut models, 200 ml REF 520 0029 0



Master-Sep penetrates into the plaster and seals the surface. Simultaneously, Master-Sep serves as a lubricant between pin and sleeves.

Arch and base can be separated more easily. A soft gliding layer is achieved by wetting the pins.



## Master pin drill unit mpb 1

The powerful, high quality and maintenance-free motor features high true running accuracy.

Accordingly, the precision of the drilled hole and the accuracy of the models are increased. Working is simplified thanks to the easy-to-operate lifting mechanism


Master pin drill unit mpb 1 REF 140 0092 0

(without 15° base) 1 spare fuse

1 flat wrench

 Master-Pin Diatit tungsten carbide bur 1 plaster collecting tray

1 plug axle 1 power cord





### Accessories:



Adapter base 15° inclination / precious wood REF 210 0044 0

made of precious wood








Tungsten carbide drill Special drill for Master-Pin Radix-K Ø 2,0 mm 3 mm shaft REF 360 0123 3



# Master-Pin Diatit tungsten carbide step drill special/yellow REF 360 0119 3





## Master-Pin Diatit tungsten carbide step drill special/red REF 360 0119 4

If the drilled hole is too large to receive the Master-Pin, the special/red Master-Pin Diatit tungsten carbide drill can be used to prepare a smaller drillhole. The diameter of this drill is 0.01 mm smaller than the one of the standard/green Master-Pin Diatit tungsten carbide drill.

## Master pin drill unit mpb 1



The diameter of the luminous spot can be adjusted individually to ensure anti-dazzling, precise focusing.




Firmly mounted model table with shape and width adapted to the arch.



Guidelines on the model table allow specific alignment of the model for exact planning of pin holes



Integrated grooves collect plaster particles and provide the precondition for arches which rest parallelly.



Exact model table mounted at an angle of 90° to the drill subsequently ensures simple removal of the arch from the model base.



The firm hold of the arch allows precise drilling of pin holes. The drill is directed to the arch without any vibration.



Simple screw mechanism for fast and precise adjustment of the drilling depth.



Drills are exchanged without the need to open the unit.





Any resulting plaster particles are automatically collected by the projecting collecting tray.



Unit, motor and collet remain clean, the collecting tray can be removed.

## Perfect model fabrication

The combination of Master pin drill unit, Master-Pin sytem, Master-Split model system and Exakto-Rock S provide the basis for precision-fit restorations.

The low expansion of Exakto-Rock S of only 0.08 % and the perfectly-fitting Master-Pins facilitate the daily work and add new brilliance to the dental restoration. This way you can improve the reputation of your laboratory.



## Master-Pin System

The pin system for perfect sawcut models.



The small drilling depth of just 4.5 mm in the arch avoids undesired perforation of the arch. Flattening of the soft plastic sleeves is ideal for pins with small distance to each other. The types of plastic of the sleeves and the design of the inner wall ensure soft and controlled removal of the dies. Ideal for bridge models.

### Your advantages at a glance



The smallest drilling depth of all pins of only 4.5 mm.
Advantage: no perforation of the arch during drilling and enhanced stability.



The Master-Pin Diatit tungsten carbide drill is adjusted so that the boundary line for drilling of the pin is exactly on the same level as the basis of the arch.



Master-Pin and Master-Pin sleeve can be easily assembled due to the taper and the rounding of the end of the pin.



The optimized glueing tip: the adhesive is spread more uniformly in the drillhole and at the glueing shaft. Advantage: safe hold of the Masterlin in the die



sleeves allows

to obtain low

sawcut models.



The sleeve rises above the Master-Pin. All Master-Pins can be clearly recognized on the underside of the model.






The funnel-shaped design of the Master-Pin sleeve simplifies assembling of die segments and model base.



The unilateral flattening of the Master-Pin sleeves serves to protect against twisting and...



solution in case of drillholes with small distance to each other.





Due to the special surface design of the inner wall of the sleeve, soft friction between Master-Pin and Master-Pin sleeve is achieved whilst ensuring maximum precision and stability





1 Preparatory work

## Master-Pin System

The Master-Pin system simplifies daily fabrication of models since the system components have been matched with each other.



Processing is simple and no new techniques need to be learned. The main advantages of the Master-Pin system are the small drilling depth and the small diameter of the drillhole. Soft integration and removal of the Master-Pin is ensured by the design of the inner wall of the Master-Pin sleeve. This is a particular advantage for bridge restorations. Easy assembling is achieved thanks to the tapering design of the Master-Pin.



Master-Pins 1000 pieces REF 360 P122 5



Master-Pin sleeves 1000 pieces REF 360 H122 5



Master-Sep Special insulating liquid for sawcut models 200 ml REF 520 0029 0

### Assortment

402 pieces

200 Master Pins

200 Master-Pin sleeves

1 Master-Pin Diatit tungsten carbide step drill

standard/green
1 Working box

REF 360 0122 6



Master-Pin Diatit tungsten carbide step drill standard/green 3 mm shaft, 1.5/2. 1 piece REF 360 0119 2



Master-Pin Diatit tungsten carbide step drill special/yellow 3 mm shaft, 1.5/2, 1 piece

REF 360 0119 3

If glueing in of the Master-Pin is too difficult, the special/yellow Master-Pin Diatit tungsten carbide drill can be used to prepare a larger drillhole. The diameter of this drill is 0.1 larger than the one of the standard/green Master-Pin Diatit tungsten carbide drill.

### Assortment

2001 pieces 1000 Master-Pins 1000 Master-Pin sleeves 1 Working box

REF 360 0122 5





Master-Pin Diatit tungsten carbide step drill special/red 3 mm shaft, 1.5/2, 1 piece

REF 360 0119 4

If the drilled hole is too large to receive the Master-Pin, the special/red Master-Pin Diatit tungsten carbide drill can be used to prepare a smaller drillhole. The diameter of this drill is 0.01 mm smaller than the one of the standard/green Master-Pin Diatit tungsten carbide drill.



Weigh resp. measure plaster and water to obtain constant results.




A thermoforming foil is placed on the impression. Uniform thickness of the arch is obtained.



The arch is trimmed to achieve uniform low height.



## Master-Pin System



The correct height of the trimmed arch is essential.



The trimmed surface can be optimized with wet grinding paper.



The inner surface of the dry arch is ground with the plaster bur H263 SH 60 slightly conically (6°) to the base.



The drillholes are positioned with the Master-Pin Diatit tungsten carbide step drill.



Drillholes are prepared- 2 for each die - beginning from the buccal direction:

1. drillhole = center of fissure

2. drillhole = approx. 3 mm away toward the palatal or lingual direction.



The correct alignment of drillholes in the arch.



The upper course of the palatal resp. lingual 6° ground edge is marked with a red pen.



The Master-Pins are precisely glued in the drillholes using cyanoacrylic adhesive.




Arch with Master-Pins glued in.



The arch as well as the Master-Pins are separated with Master-Sep.



The thicker end of the Master-Pin sleeves is put on the Master-Pins.



Even in case of Master-Pins that have only very little distance to each other, the Master-Pin sleeve can be easily used due to the lateral flattening.



The Master-Pin sleeves rise from the Master-Pins by approx. 0.5 mm so that uniform, constant height of the arch is always ensured



The Master-Split system is used to prepare base for the arch.



Place the prepared arch into the Master-Split model former and align it.



Base plaster is filled up to 1 mm below the deepest point of the red marking line (fig. 10).



After the plaster has hardened, remove the model by pressing it out of the Master-Split model former.



The removed model will receive a Split-Cast separation: the Master-Split during the preparation of the arch without any additional



Prior to trimming the model, the Master-Split base former is removed.



The sawcut model is trimmed to the smallest size possible.



The trimmed and dried working model.



The arch is removed from the model base towards the pins - parallely and without tilting.



The base of the arch and the model base must be thoroughly cleaned after trimming to ensure high precision and perfect aesthetics.



The green Master-Pin sleeves are all on the same level and can be clearly recognized on the underside of the model.



The die segments are separated using a Giflex-TR diamond disc.



Perfect fit of the working dies on the model base



It is also possible to place interdental Master-Pins that are not glued in.



Aesthetically appealing and functional models simplify daily work.



A fine and precise dental restoration is created on a fine model.



## Master-Split model system

A universal model system for economical model fabrication for all dental-technical indications.



Matched with the Master-Pin system, Master x-tray and Master pls 44. Simple and precise fabrication of the base with integrated Split-Cast which requires little space due to its shape. Three different model formers for crown and bridgework, implant prosthetics, CoCr restorations, full dentures and repairs.

### Your advantages at a glance

• helps to save time Production of the model (secondary base) in a single working step.

• plaster can be saved The respective impression size determines which of the three Master-Split model former sizes is used.

The plaster consumption is reduced to the minimum.

• high precision Since the model is produced directly on the secondary base (Master-Split base former), a perfectly

smooth, precisely fitting model underside is achieved.

extended reusability
 All individual components of the Master-Split model system are reusable and durable.

• excellent cost/benefit ratio Since time and plaster are saved, the favourably priced Master-Split model system pays for itself already

after it has been used a few times.

• optimized handling Each model will automatically obtain a Split-Cast separation. Due to the model-articulator separation,

working is performed on a small, easy-to-use and functional model.

• small height Even in cases of limited space (face-bow model assembly, etc.) the small height of the Master-Split


base former allows the use of the Master-Split model system.

• increased safety Due to the additional octagon platform the model is safely and precisely fixed on the Master-Split


base former even in the case of lateral movement in the articulator.

• perfect aesthetics Models produced with the Master-Split model system excel by their aesthetic appearance.

### Underside



## Upper side



## Master-Split model system

A universal model system for economic model fabrication for all dentaltechnical indications.



Each size of the Master-Split model system consists of two elements. Thanks to the three different Master-Split model formers, the correct size is always available for any size of arches or impressions. Saving of plaster is possible thanks to the range of different sizes. When mounting in the articulator, sufficient space is always ensured due to the small height of the Split-Cast. The surface of the material allows easy cleaning.



Master-Split model former small 2 pieces REF 360 0118 K



Master-Split model former medium 2 pieces REF 360 0118 M

## Assortment

- 1 Model former
- 3 Basis former
- 3 Metal magnetic plates

REF 360 0124 K



Master-Split model former large 2 pieces REF 360 0118 G



Master-Split base former 10 pieces REF 360 0118 0

### Assortment medium

- 1 Model former
- 3 Basis former
- 3 Metal magnetic

REF 360 0124 M



Metal magnetic plates 50 pieces REF 360 0118 1

### Assortment large

- 1 Model former
- 3 Basis former
- 3 Metal magnetic plates

REF 360 0124 G

### Application examples



bridges



Implants combined prosthetic work



dentures and CoCr work




Situation models, repairs

Tip







To ensure exact contact of the model on the Master-Split base former, the completed model is smoothed and cleaned with sanding paper 2 to 3 times. Wax or dirt that will deposit on the four model skids later on will not affect the precision.

### Care and cleaning



Plaster and wax residues can be easily recognized on the signal-green plate so that precise working is simplified.





The Master-Split base formers and the Master-Split model formers only need to be cleaned under running water since their surfaces are extremely smooth and self-separating. They are intended to be used for plaster and matched with this material. Additional separating is not



## Master-Split Modellsystem

### **Processing**



Regardless which arch or impression size is used



the Master-Split model formers fit in every case.



The green Master-Split base former – the matching counterpart to the underside of the model.



The Master-Split model former is selected according to the size of the impression or the arch.



The Master-Split base former is first inserted at the rear edge.



After fitting, the plate is positioned and pressed in using both hands.



Only then the plate is pressed in again on the table.



The plate is only then inserted properly, when there is a 0.1 mm high step at the edge.



The Master-Split metal magnetic plate is centered on the Master-Split base former.



The arch is aligned according to the markings of the Master-Split model former.



In the case of sawcut models the model base is generally prepared with a liquid base plaster, for example Fluid-Rock.



After the base plaster has hardened, the model is pressed out.



During the preparation of the base the removed die model obtains a Split-Cast separation - the Master-Split - without additional work.



Due to the special shape of the sleeve, an indentation is obtained at the model base which simplifies the removal of the plate.



The Master-Split base former is removed before the model is trimmed.



The model is trimmed with the plaster trimmer to obtain a perfect size.



The trimmed and dry working model.



If the arch is adequately prepared, trimming is no longer required after preparing the base.

### Split-Cast check



The position of the model can be easily checked despite the fixed magnet.



### The solution for a familiar problem



The plaster Split-Cast must be trimmed to be integrated into the articulator.



An individual plaster control base is much thicker than



Plaster control bases often cause problems when they are integrated into articulators



The Master-Split base former is the better choice



the specially shaped Master-Split base former



whereas the Master-Split base former always provides sufficient space.

### Fluid-Rock

### Fluid-Rock is a smoothly flowing class IV super-hard stone to prepare bases for models.



The light-blue color can be easily combined with all colors for the arch. The extended processing time allows to pour several bases at the same time. The thin consistency results in perfect flow characteristics and allows to obtain bubble-free models.

### Color - blue:

1 x 2 kg REF 570 0FB5 2 5 x 2 kg REF 570 0FB5 1 10 x 2 kg REF 570 0FB5 0

### Technical Data - Fluid-Rock

### Mixing ratio

Processing time

Setting time (Vicat time)

Comp. strength aft. 1 hr Comp. strength aft. 24 hrs Setting expansion distilled water approx. 6 min at 18° to 20° C approx. 11 min at 18° bis 20° C 48 N/mm<sup>2</sup> < 0,06 % (no further expansion after 2 hours)

### Processing in the ecovac unit:

Vacuum level 1, mixing speed: 390 rpm



Mix Fluid-Rock base stone in the ratio of 100 g powder and 25 ml distilled water to achieve a highly fluid consistency.



Fluid-Rock base stone is directly poured into the model former without using a vibrator. Perfect flow characteristics allow to obtain models without any bubbles.



Low expansion ensures constant quality when producing models. Perfectly matched with Thixo-Rock super-hard stones.

### Arti-Rock

### Low-expansion articulating stone for precision-fit restorations.



Low expansion of only 0.02 % ensures accurate position of the model when aligning according to the anatomic situation. Accurate restorations and reduced grinding time are obtained. Perfect stability and special adhesive capacity simplify mounting in the articulator and ensure safe retention of the models.

### Color white:

1 x 4 kg REF 570 0AR0 4 1 x 18 kg REF 570 0AR1 8

### Technical Data - Arti-Rock

Color Mixing ratio

Processing time span
Setting time (Vicat time)
Compressive strength
according to DIN
Expansion

white 100 g / 40 ml dist. water approx. 3 min. 5 min.

7.2 MPa

0.01 % after 20 min. 0.02 % after 48 hrs.



The smooth consistency allows trouble-free and precise mounting of models in the articulator.



When using keys, accurate reproduction of details is achieved thanks to smooth processing of the stone. The final hardness allows easy processing.



The short setting time and low expansion are perfect preconditions for accurate rebasing.



## Thermo-syringe



Fixing and glueing, that can be dissolved quickly without any residues, for any type of model situation.

The adhesive resin wax can be moulded by heating and easily placed on the models.

Thermo-syringe 1 Piece REF 110 0121 1



After heating, the adhesive resin wax is directly applied onto the glueing point using the thermosyringe. Firm bonding is ensured.



The adhesive resin wax can be applied onto any type of material. Afterwards it can be removed from the objects without leaving any residues.

### Accessories:



Adhesive resin wax Pack cont. 250 g Bucket cont. 1000 g

The extra long and nar-

row blade is perfectly

suited for cutting off

excess plaster in the

lingual region.

REF 510 0070 1 REF 510 0070 0

## Plaster knife



Multi-purpose knife with ergonomically designed plastic handle for optimum transfer of force, facilitates your daily work.

- Long blade made of stainless hardened steel.
- Dimensionally stable, easy-to-clean hard plastic handle. Ergonomic shape for right and left hand use.
- Multi-purpose element for easy removal of impression tray. Features impact surface with opposing chisel.



The special cones on the multi-purpose element simplify removal of the impression tray from the model.

Plaster knife
1 piece
REF 310 0011 4



perfectly trimmed with the permanently sharp and stable blade.

Plaster edges can be



When opening flasks, the lateral chisel ensures improved transmission of force thanks to the high leverage effect of the knife handle.



A separate impact surface has been added opposite the chisel to protect the back and the blade of the knife.

## Dentaclean plaster removing agent / speed



Ready-to-use solution to remove plaster residues from all surfaces.

The Dentaclean plaster removing agent is available in two types: normal and Speed. The ready-to-use solution removes plaster residues from all surfaces. If no time is to be wasted, Dentaclean Speed should be used.

Plaster removing agent 1000 ml REF 520 0011 9 2500 ml REF 520 0099 3 Speed plaster removing agent 1000 ml REF 520 0101 0 2500 ml REF 520 0099 4



Hard plaster particles are carefully reduced from the mixing bowl (cup) without any damage.



Gentle and fast removal of plaster protects the resin surface and the color.

### Giflex-TR


### Perforated design for controlled sawcuts.





Giflex-TR is a disc that features diamond-coating on both sides and is particularly suitable for cutting plaster and resin dies. Calculated chip spaces in the area of the diamond coating ensure quick removal of the grinding dust and increase the cutting performance of the disc. Giflex-TR allows quick, smooth and reliable cutting even of very hard plaster and resin. Troublesome chattering and jamming of the disc is avoided.

Larger holes in the diamond-free section reduce the friction heat. The disc will not overheat even if deep cuts are carried out. The high running transparency allows a better view onto the saw cut.




Giflex-TR diamond discs are coated on both sides and ready mounted.

Shaft diameter:	Standard 2.35 mm	Standard 2.35 mm	Standard 2.35 mm	Standard 2.35 mm
REF	340 0002 5	340 0012 0	340 0002 0	340 0011 0
ISO-No	806 104 377514 250	806 104 377514 300	806 104 377514 370	806 104 377514 450
Diameter (D):	25 mm	30 mm	37 mm	45 mm
Length (I):	0.3 mm	0.3 mm	0.3 mm	0.3 mm
Recommended speed:	20,000 rpm	15,000 - 20,000 rpm	15,000 - 18,000 rpm	10,000 - 15,000 rpm

## Giflex-TR Master x-tray

### Special diamond disc for processing resins.



Diamond grinding disc Giflex-TR Master x-tray REF 340 00M2 5

The Giflex-TR Master x-tray features a coarse diamond grit to achieve a cooling effect already in the diamond-coated section when separating resins.

## Tungsten carbide burs for processing of plaster

## Quick shaping and smooth surfaces for all types of plaster.

The relief protects the sharp blade against breakage of edges. This way the service life of the relief tools is three times longer than the one of comparable burs. Additionally, the processed surface is smoother and a luster is added.



Tungsten carbide 1 piece REF H263 SH 60



The super-coarse crosscutting edge allows removal of large quantities of any type of plaster.



Tungsten carbide 1 piece REF H263 GH 60



The coarse cutting edge allows finer cuts and avoids splintering of the plaster.

### For exact determination of the preparation margin for all die materials.



Rapidy microbur 1 piece REF H001 NH 31



The cross-cutting edge produces smooth and precise ditches.



Preparation bur 1 piece REF H263 GH 30



The cylindrical round shape allows to prepare an oblique ditch so that the preparation margin can be recognized more easily.

## Litebloc UV

Light-curing resin for blocking-out cavities and building up dies.





Litebloc UV 3 g REF 520 0098 0

The screwable tube allows application of the desired quantity.

The high adhesive capac-

ity of the undercut wax

offers reliable hold in the

cavity.




The fine dimensional accuracy allows perfect filling of the cavity.



After a short setting time, Litebloc UV can be coated with any die varnish.

### Undercut wax

Precise blocking out of all cavities on the die.



The undercut wax has a high melting point and is therefore perfectly suited for blocking out cavities. No bond with the dipping wax is formed.



**Undercut wax** 25 g **REF 510 0048 0** 



Low shrinkage and optimum scraping capacity simplify blocking out.



The high melting temperature also allows the use of the wax below immersion wax copings.



## Light-curing die varnish

## For smoothening and hardening the plaster



Depending on the plaster and modelling wax, different colors are available. The desired layer thickness can be achieved by applying the varnish several times and can be checked with the help of the color intensity.

### Light-curing die varnish

🕨 red, 20 ml REF 540 0100 3 yellow, 20 ml REF 540 0100 4 REF 540 0100 5 green, 20 ml blue, 20 ml REF 540 0100 0 O transparent, 20 ml REF 540 0100 6



The disposable brush allows precise application. The layer thickness can be varied by applying the material several times.



The light-curing die varnishes produce a particulary hard surface which protects the die against damage when fitting on the crowns.

Five different colours to ensure contrast to any type of modelling wax.



The varnishes are translucent. If they are applied several times, the colour becomes more intense so that the layer thickness can be controlled.



To produce a cement gap, the varnish must be cured immediately after applying. For hardening of preparation margins: Allow die to soak into the plaster, then polymerize. The varnish hardens the surface without layering.



Brush holder, bent 12 pieces REF 330 0114 1



Brush holder, straight 12 pieces REF 330 0114 9



Disposable brush 100 pieces REF 330 0114 2



Mixing block 10 pieces REF 330 0114 4

## Die varnish, light-curing, opaque


### Swift application thanks to good masking capacity.



The opaque die varnishes simplify uniform coloring of the varnish coat. The brush is already integrated in the lid. When cleaning the die with steam, the varnish coats remain intact.



Die varnish light-curing opaque is available in three different colors. The fine masking capacity allows to obtain a uniform color of the varnish layer.



red, 20 ml green, 20 ml

blue, 20 ml

REF 540 0010 4

REF 540 0010 3

REF 540 0010 1



During the application the die varnish diffuses into the plaster surface. Depth polymerization leads to abrasion-resistant bonding to the die. Light-curing die varnish resists high mechanical stress. Even steam-blasting units do not affect the strong bonding.



## diephos dentine

Light-curing varnish for abrasion-resistant surface on dies.



The use of diephos dentine allows simple assessment of the esthetic design and color of all-ceramic restorations.

- quick application
- · abrasion-resistant surface by light-curing
- assessment of esthetic appearance on the working model
- resistant to steam

diephos dentine die varnish, tooth-colored 10 ml REF 540 0010 0



Two coats of diephos dentine are applied. The first coat is applied very thinly and light-cured for at least 90 seconds. The second coat is applied in a way to ensure complete coverage and light-cured again for at least 90 seconds.



The thickness of the coat can be varied by applying several layers. This way a cement gap can be defined for all-ceramic restoration. The thicknes of a single layer is approx. 12 micro-meters. Light-curing is required after the application of each layer.



Separating liquid can be easily applied and waxups can be prepared on diephos dentine.



The modelling material will not falsify the tooth color. An opaque-like effect is achieved on the subsurface.

## Spacer varnish gold, silver, silver-blue, blue

Air-drying varnishes with metal components for scratch-resistant surfaces.



Specific layer thicknesses of 5  $\mu m$  or more can be achieved with these spacer varnishes. Each additional coat increases the layer thickness accordingly. The metal components of the gold, silver and silverblue spacer varnishes produce highly abrasion-resistant surfaces and hence protect the die. The blue spacer varnish can be used to determine and eliminate premature and undesired contacts when fitting a framework on the model.



Spacer varnish gold 20 ml REF 550 0000 5

Spacer varnish silver 20 ml REF 540 0071 7

Spacer varnish silver-blue micro, 20 ml REF 550 0000 6 Thinner for spacer varnish silver and gold 20 ml REF 540 0070 1

Thinner for spacer varnish blue and silver-blue 20 ml REF 540 0069 0 2


so that the amount of work is reduced.

The spacer varnishes can be applied easily and dried quickly

The spacer varnishes gold and silver produce a layer thickness of approx. 10 µ; the spacer varnish silver-blue produces a layer thickness of approx. 5 µ.



The spacer varnishes contain metal components. They produce a particularly abrasion-resistant surface which protects the die against damage.



spacer varnish blue 20 ml REF 550 0000 7

thinner for spacer varnish blue and silver-blue 20 ml REF 540 0069 0



The area of the cement gap of 8 – 10 µm in the inside of the crown can be easily recognized thanks to the clear color contrast.



Since the blue spacer varnish is well suited to detect undesired contact points, it can also be used as an alternative to occlusion spray.



The blue spacer varnish can be applied selectively to avoid overlaps which may result from non-uniform application of spray.



Accordingly, early contact points can be quickly eliminated.



## The sawcut model

## Gloss and hardening agent for plasters

Scratch-resistant surfaces for all plaster types without layering.

The gloss and hardening agent for plasters renders the model or die resistant to scratches. Simultaneously, a lustrous surface is achieved with a layer



thickness of only 2  $\mu$ m.

Gloss and hardening agent for plasters 20 ml REF 550 0000 1 100 ml REF 550 0000 2



Without the hardening agent models can be dam-aged when the restoration is placed on the model.



The specially adjusted consistency leads to the diffusion into the plaster surface. The high edge stability and scratch resistance avoid damage of any kind.



The gloss and hardening agent for plaster has hardened after only 2 minutes.



Gloss and hardening agent for plaster diffuses into the plaster so that it can also be used on the preparation margin.





### The functional margin

The use of functional margin wax enables exact preparation of the functional margins of the model.

Only high-quality products are used for the fabrication of the master model. The precision of fit to be achieved is ensured by the low-expansion and thixotropic Exakto-Rock S stone. The ecovac vacuum mixing unit with its abrasion-resistant mixing cups is recommended for proper and homogeneous mixing of the materials. The matched range of products facilitates the fabrication of the model and enables convenient working.



### The resin dies

The stable resin dies are quickly and easily produced using Pi-Ku-Plast modelling resin and retention pins to ensure reliable retention.



### The proper mixture

The special mixing spiral and abrasion-resistant stainless steel mixing cups guarantee a homogeneous mixing process.

One of the most important components for your success!



### The model materials

With a final expansion of only 0.08% after 2 hours, Exakto-Rock S offers utmost precision. Both stones impress with rapid hardening and good processing properties. 5 different colors for all requirements.



### The Split-Cast system

The model formers with integrated Split-Cast facilitate the fabrication and help to avoid compromising the accuracy of fit when assembling models.



### The control

The KoEx measuring device enables the control of expansion and contraction of various materials to verify the precision so that materials can be processed in a very reliable manner.

### Materials used:


Functional margin wax	page 36
Retention pins	page 36
Pi-Ku-Plast / Pi-Ku-Plast HP 36	page 37
ecovac vacuum mixing system	page 38
Exakto-Rock S	page 38
Exakto-Form	page 39
Master-Split model system	page 40



# The master model

# Functional margin wax

To produce perfect functional margins.



To produce perfect functional margins. The slightly sticky, flexible functional margin wax allows simple and safe positioning to each impression material. Final fixation is achieved by waxing up. Accordingly, uniform design of functional



Functional margin wax REF 430 0150 0



Uniform and ideal functional margins in the model guarantee perfect fit of the denture.

# Retention pins

margins is possible.

The retention pins feature retentions to guarantee safe hold when fabricating resin dies.



Also perfectly suited for milled models.



The particularly stable retention pins can be used for all resin dies.

tray represents the best prerequisite for precise models with a perfect

functional margin.



Retention pins 100 pieces REF 360 0000 1 500 pieces REF 360 0000 2



# Pi-Ku-Plast / Pi-Ku-Plast HP 36

### For precision-fit and stable resin dies in next to no time.



Brush resin in 5 different colors. Both resins only differ in their contraction values. HP 36 features a contraction of 0.036 %. Since the resin hardens quickly, it is perfectly suitable for the fabrication of resin dies or resin copings in the double crown technique.



Apply a thin coat of vaseline to the inside of the crowns



The excellent modelling characteristics allow precise filling of the crowns within a very short time.



Pi-Ku-Plast HP 36 features a very short setting time. Therefore the retention pins can be placed directly into the resin.



This way, Pi-Ku-Plast HP 36 allows to produce accurate and particularly stable resin dies within a very short time.



Resin dies are the perfect basis for precision-fit dentures.



The high stability of Pi-Ku-Plast HP 36 allows to obtain a stable basis for all types of milling work.

### Assortments, large Pi-Ku-Plast

3 vessels brush size each A+B 1 brush holder

100 ml cleaner 100 ml monomer 85 g polymer

blue REF 540 0017 3 REF 540 0017 4 yellow REF 540 0017 5 orange • red REF 540 0017 6 transparent REF 540 0017 7

### Assortments Pi-Ku-Plast HP 36

3 vessels 1 brush size each A+B 1 brush holder 100 ml cleaner 100 ml monomer 85 g polymer

blue yellow orange • red transparent REF 540 0021 9 REF 540 0021 7 REF 540 0021 8 REF 540 0022 0 REF 540 0021 6

### Refill packages

100 ml cleaner 85 g polymer 100 ml monomer blue yellow orange

red transparent REF 540 0017 8 REF 540 0017 9 REF 540 0018 0 REF 540 0018 1

REF 540 0019 6

REF 540 0016 9

REF 540 0016 7

REF 540 0016 8

Refill packages 100 ml cleaner

85 g polymer 100 ml monomer

blue yellow orange red transparent REF 540 0021 5 REF 540 0021 3 REF 540 0021 1 REF 540 0021 2 REF 540 0021 4

REF 540 0021 0

REF 540 0022 4

### Refill packages

REF 540 0017 2 Vessel cleaner, 8 ml REF 540 0017 1 Vessel monomer, 8 ml Vessel polymer, 8 ml REF 540 0017 0 brush size A and brush holder, pack. cont. 3 pcs RFF 330 0114 6 brush size B and brush holder, pack. cont. 3 pcs REF 330 0114 7

### Refill packages

Vessel cleaner, 8 ml REF 540 0020 9 REF 540 0020 7 Vessel monomer, 8 ml REF 540 0020 8 Vessel polymer, 8 ml REF 330 0114 6 brush size A and brush holder, pack. cont. 3 pcs brush size B and brush holder, pack. cont. 3 pcs REF 330 0114 7

# Assortments, small

Pi-Ku-Plast

20 ml cleaner 2 mixing trays, silicone, red 20 ml monomer red

1 brush size B + brush holder

12 g polymer



Pi-Ku-Plast separating varnish 10 ml REF 540 0018 2

# ecovac vacuum mixing system



### ecovac

### Precision-fit restorations obtained through optimal use of material properties.

The user-friendly and compact design simplifies work and reduces sources of errors. A powerful and maintenance-free vacuum pump, adjustable in two different levels (15 mbars, 200 mbars), ensures bubble-free mixing of materials and results in a perfect casting surface. Stirring time and speed can be adjusted continuously to allow correct processing of different materials.

ecovac (230 V)

REF 140 0093 0

(Wall mounting, without mixing cup and base) 1 mains cable

- 1 spare filter
- 1 drilling template for wall mounting
- 4 screws and plugs for wall mounting

Accessories:

Base ecovac, 1 piece

REF 210 0045 0







### ecovac mixing spiral

The mixing spiral takes up the components to be mixed from all areas of the mixing cup and stirrs them horizontally and vertically. No unmixed materials will remain on the bottom of the mixing cup, which may cause different expansion of the material later on.

All features and components listed provide increased reliability, lead to improved fit when preparing dental restorations and avoid time-consuming

More information and reference numbers on page 16.

### Exakto-Rock S


Exakto-Rock S is a synthetic super-hard class IV stone with distinctive thixotropy and improved flow properties and is free from formaldehyde.



Reduced expansion is completed after 2 hours and is only 0.08 %. As a result, highly accurate impressions are enabled and precision-fit restorations are ensured. Moreover, thanks to optimized reflection of light, Exakto-Rock S is suitable for scanning and available in crown and ivory.

- fomaldehyde-free stone ensures safe processing and can be used for the fabrication of holistic dentures without any problem
- synthetic components ensure consistent quality and enable the fabrication of precision-fit models
- improved flow characteristics facilitate pouring of several impressions
- optimized reflection of light thanks to special dye pigments reduce reworking/adjustments in the CAD system

### Expansion of various other stones





The excellent processing time span allows bubblefree pouring of numerous impressions with just a single mix.



Exakto-Rock S offers high stability on the spatula and thixotropic consistency on the vibrator. Simple and clean processing is ensured.



Absolutely accurate reproduction of dimensions of the oral situation thanks to the minimal expansion value (< 0.08 %) so that precision-fit dentures are obtained.

More information and reference number on page 17.



### Exakto-Form

Model resin for accurate reproduction and maximum edge stability in five different colors.



Processing does not require to change familiar working processes.



Component A
light-ivory
1 x 50 g
REF 520 0017 6



Component A olive green 1 x 50 g REF 520 0017 4

Component A

REF 520 0017 8

yellow

1 x 50 g



Component A grey 1 x 50 g REF 520 0017 5



Component A signal blue 1 x 50 g REF 520 0017 7



Component B 1 x 50 g REF 520 0017 3

Assortments cont. 240 g each

6 x 20 g A yellow

6 x 20 g B REF 520 2028 4

6 x 20 g A light-ivory

6 x 20 g B REF 520 2028 2

6 x 20 g A olive green

6 x 20 g B REF 520 2028 0

Assortments cont. 600 g each

6 x 50 g A yellow

6 x 50 g B REF 520 0028 4

6 x 50 g A light-ivory

6 x 50 g B **REF 520 0028 2** 

6 x 50 g A olive green

6 x 50 g B REF 520 0028 0

6 x 50 g A grey

6 x 50 g B REF 520 0028 3

6 x 50 g A signal-blue

6 x 50 g B REF 520 0028 1

Accessoires:



Stirring sticks 250 mm long, 100 pcs REF 390 0031 0

Mixing cups 120 ml, 100 pcs REF 390 0030 0



Measuring syringes 20 ml, 50 pcs REF 390 0036 0



Exakto-Form Insulating liquid 125 ml REF 520 0021 0



Prior to mixing, each component must be stirred so that a homogeneous mixture is obtained. Mix the sediment completely.



Add component A to component B; empty tin completely.



Mix Exakto-Form approx. 30 sec. until a uniform colour is obtained.



Two tins of Exakto-Form are sufficient to produce 2-3 complete dental arches.



The material can be removed after just 30 minutes. Final hardness is achieved after 90 minutes. Then the material can be trimmed.



If a base for the model is to be produced with Exakto-Form, the model must be previously insulated with Exakto-Form insulating liquid.



Due to its high edge stability Exakto-Form is perfectly suitable for precision-fit bridge and crown work.



Any technique can be used for sawing Exakto-Form models. Familiar working processes do not need to be changed.



If smaller quantities are used, fill component A and B into a separate syringe.



Fill equal portions of Exakto-Form into a silicone cup (approx. 2 ml each for one die) and mix to obtain a homogeneous consistency. Please note: material in the syringes must be processed within 5 days.



Pour Exakto-Form into the impression. The excellent flow properties avoid the formation of bubbles even in impressions with thin edges.




The hardened resin can be drilled and trimmed. The stability avoids dimensional changes and guarantees precise models.



# Master-Split model system

The universal model system for economic fabrication of all areas relevant in dental techniques.



Matched with the Master-Pin system. Simple and precise fabrication of the base with integrated Split-Cast which requires little space due to its shape. Three different model formers for crown and bridgework, implant prosthetics, CoCr restorations, full dentures and repairs.

### Your advantages at a glance

• helps to save time Production of the model (secondary base) in a single working step.

• plaster can be saved The respective impression size determines which of the three Master-Split model former sizes is used.

The plaster consumption is reduced to the minimum.

• high precision Since the model is produced directly on the secondary base (Master-Split base former), a perfectly

smooth, precisely fitting model underside is achieved.

• extended reusability All individual components of the Master-Split model system are reusable and durable.

• excellent cost/benefit ratio Since time and plaster are saved, the favourably priced Master-Split model system pays for itself already

after it has been used a few times.

optimized handling
 Each model will automatically obtain a Split-Cast separation. Due to the model-articulator separation,

working is performed on a small, easy-to-use and functional model.

• small height Even in cases of limited space (face-bow model assembly, etc.) the small height of the Master-Split

base former allows the use of the Master-Split model system.

• increased safety Due to the additional octagon platform the model is safely and precisely fixed on the Master-Split

base former even in the case of lateral movement in the articulator.

• perfect aesthetics Models produced with the Master-Split model system excel by their aesthetic appearance.

# Underside Sloping for enhanced mounting in the articulator Sloping for enhanced mounting in the articulator

# Easy removal of the model from the plate Only one plate size for all model sizes Raised octagon platform for more safety during plastic surface for smooth models Raised octagon platform for more safety during lateral movement in the articulator

### Application examples



Crowns and bridges



Implants and combined prosthetic work



dentures and CoCr work



Situation models, repairs

More information and reference numbers from page 25.



# The implant model



### The soft gingival mask

After fixation of the laboratory analogs, the soft gingival mask is simply placed into the impression using the dispensing device. This way the soft tissue (gingiva) is exactly reproduced.

To avoid any change of the implant position caused by expansion while the plaster hardens, bredent offers a stone (plaster) with an expansion of only 0.08% after 2 hours. As a result, the oral situation can be accurately reproduced.



### Model fabrication

The thixotropic super-hard stone is used to fabricate the model in a way to save material.



### The model system

The Master-Split model systenm facilitates the fabrication of high-quality models and reduces the mounting height of the models thanks to the integrated Split-Cast.



### The working model

The use of high-quality materials creates the basis for high precision of fit of dental restorations. A harmonized system – from the impression to low-expansion stone and the casting technique!

### Materials used

Multisil-Mask soft	page 42
Multisil-Mask hard	page 43
haptosil D	page 44
KoEx Measuring Device	nage 45



# The implant model

# Multisil-Mask soft



### Accurate reproduction of gingival tissue.

Quick and economical processing with the cartridge system and the especially adjusted silicone allow trouble-free direct application into the impression or the matrix. The natural color of the gingival mask supports perfect shade determination of the veneer. Overdimensioning of margins is recognized immediately.



The gingival situation on the unsawed sawcut model ...



... is reproduced using haptosil D kneading silicone and then the arch is sawed.



The sawcuts are coated with wax.



esthetic



Mixing cannulas Size 1 / blue REF 320 0045 0



Openings (inlet and outlet) are drilled into the matrix using the locating matrix drill and Multi-Sep is applied.



informative



Multisil-Sep 10 ml insulating liquid REF 520 0100 3




The dispensing device with cartridge and can-nula is held to the open-ing. Whilst applying the material from the dosing device, the matrix is fixed on the model ...







... to obtain the correct position of the gingival mask.



Accessories:



Dispensing device 1 piece REF 320 0044 0

### Multisil-Mask hard



# Special resin for hard gingival masks featuring stable consistency and ideal processing characteristics.

The hardness allows torsion-free and accurate placement on the model. The Vario-Stud-Snap vks-oc system is used for fixation.

Divergent implants are aligned using the implant compensating cones developed by bredent.



Multisil-Mask hard 50 ml cartridge 1 piece REF 540 0113 3

Mixing cannula blue 12 pieces REF 320 0045 0



Ø 3.5 mm, 12 pieces **REF 430 0703 5** Ø 4.0 mm, 12 pieces

**Implant** 

Ø 4.0 mm, 12 pieces REF 430 0704 0 Ø 4.5 mm, 12 pieces

compensating cones

REF 430 0704 5 Ø 5.0 mm, 12 pieces REF 430 0705 0

eces each Ø 5.5 mm, 12 pieces 0; 5.5; **REF 430 0705 5** 





Matrix vks-oc 1.7 mm 8 pieces REF 430 0659 0

# Assortment

2 x 50 ml Multisil-Mask hard in cartridges 24 pieces mixing cannulas

- 1 Assortment implant compensating cones 8 pieces patrixes vks-oc 1.7 mm
- 8 pieces matrixes vks-oc 1.7 mm

REF 540 0113 4



Implant compensating cones 20 pieces, 4 pieces each 3.5; 4.0; 4.5; 5.0; 5.5;

REF 430 0739 2



Metal transfer patrixes 1.7 mm 8 pieces REF 430 0662 0

Dispensing device 1 piece

REF 320 0044 0

### Accessories:



Light-curing die varnish transparent 20 ml REF 540 0100 6



The marginal fit of the individual abutment to the implant can always be checked.



Multisil-Mask hard permits reliable adapting of individual attachments and framework designs.



Accurate placement of pontics can be easily achieved with Multisil-Mask hard.

### Processing



Initial situation of the implant restoration with laboratory abutments.



Place the implant compensating cones on the laboratory abutments in a way to ensure that the wide side is in the angulated area.



Fill Multisil-Mask hard around the laboratory abutments at the height of the compensating cones.



Use tweezers to insert the matrixes vks-oc into the soft resin immediately after injecting the resin.



Trim the gingval mask from the basal side to obtain a straight margin.



Apply vaseline to separate the gingival mask against plaster.



Snap in of the metal transfer patrixes in the matrixes.



Fill the impression with Thixo-Rock and then ...



... box the impression with the Master-Split model system.



Use an instrument to lift the gingigal mask carefully off after boiling out the compensating cones.



The gingival mask is safely retained by the vks-oc matrixes and can always be repositioned in an accurate manner.



The completed gingival mask. Apply transparent die varnish to protect the gingival mask against scratches and to improve the aesthetic appearance.



# The implant model

# haptosil D






Addition-curing kneading silicone with a Shore A hardness of 90 for the fabrication of stable and exact keys and models.

Models for repairs and extensions can also be fabricated within a short time. Consequently, an enormous amount of time can be saved compared to conventional model fabrication.



### haptosil D Component A and B 1300 g each REF 540 0118 0

haptosil D Component A and B 7500 g each REF 540 0119 0



Precise reproduction of details with haptosil D reduces reworking time since highly accurate keys can be produced.



Equal quantities of haptosil D are removed using the enclosed portioning spoon.



Processing time is 90 - 120 seconds after mixing both components. Both components are kneaded to obtain a homogeneous mixture.



Make sure to achieve a uniform color when mixing the components; only after a uniform color is obtained, haptosil D has been properly mixed and hardens completely and evenly.



Exert uniform pressure to press haptosil D onto the desired spot/area. The softness allows accurate impressions also in areas difficult to access.



The Shore A hardness of 90 ensures a stable and safe key, which will not tear when it is removed.

### Technical data haptosil D

Addition-curing kneading silicone

Hardness: 90 Shore A
Tear strength: 4.86 N/mm²
Deformation under pressure: 1.24%
Elongation: 15%
Processing time: 90-120 seco

Processing time: 90-120 seconds
Hardening time: approx. 5-6 minutes

# KoEx Measuring Device

Why are there differences in fit between model and oral situation? For the first time, the KoEx measuring device allows to measure contraction and expansion of materials.



Facilitating contraction and expanding measurements for the first time

Why do discrepancies in fit exist between the cast and the intraoral situation?

KoEx Measuring Device 1 piece including 2 contraction inserts REF 110 0148 0

### Impression Materials, Contraction

Studies have indicated that impression materials differ greatly in their contraction (shrinkage) behavior, reproducing the oral situation inaccurately. The brecision impression material provides stable values after two hours, permitting




further processing to be performed rapidly.

### Dental stone, Expansion

Check the expansion values for your dental stone materials and compare them to those of bredent's Exakto-Rock S. Exakto-Rock S expands by a maximum of 0.06 % after two hours; after 48 hours, the expansion is still less than 0.08%.



# Comparison of contraction values for A silicones and polyether materials 10.140 0.120 0.0000 0.0000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0




### Silicone duplicating materials, Contraction

Contraction measurements of different silicone duplicating materials have shown substantial differences between these. Exaktosil N15 was stable at 30 minutes, at 0.02 %. The values for other duplicating silicones changed after



24 hours, adversely affecting the fit of the restoration.



### Investment compound, Expansion

Investment compounds that can be controlled exactly and individually are a prerequisite for non-precious alloy precision one-piece attachment casting as well as for K+B

plastic injection molding using thermopress 400.



# THE GUIDE FOR MODERN IMPLANT **PROSTHETICS**


A comprehensive guide on the rapid development of implant prosthetics with numerous patient cases - also by guest authors - and an exciting overview regarding the topic of starting out in dental implantology under the slogan "KISS" - keep it simple and successful.

**REF 9929700D (German)** REF 992970GB (English)

# Prosthetic planning / 3D implant planning

Prosthetic planning / 3D implant planning is the first stp in the fabrication process of an esthetic and functionally appealing dental restoration. Thorough and careful planning is required to ensure perfect fit of the patient's restoration.

The use of the materials listed on the following pages enables you to meet any demands and requirements on planning of the dental restoration in a perfect manner.



Bite registration with mushbite	
security-bite blue	48
Bite registration with bite blocks	
Bite blocks	49
Diagnostic wax up/Set-up	
Aesthetic wax veneers	50
Life-Color-Wax	50
beauty setup	50
Set-up wax	
Modelling wax pink Standard	
Conventional implant planning	
X-resin OPG	52
OPG-ruler	52
SKY X-ray template	52
Navigated implant planning	
System steps X-ray template	53
mini <sup>1</sup> SKY	
X-resin	55
The radiopaque silicone varnish	
X-resin flow	57
Fabrication of drilling templates	
op .	

# Bite registration with mushbite

# security-bite blue



### The solution for perfect bite registration

In order to satisfy the requirements of daily dental routines, bredent involved dentists and dental technicians for the development of the dimensionally stable bite registration material security-bite blue.

Precise bite registration is the key to accurate articulation and occlusion and avoids time-consuming grinding for the user.

- High dimensional stability
- Optimized final hardness (hardness 92 Shore A / 40 Shore D)
- Does not flow into the interdental spaces
- Can be easily processed by cutting or milling
- Neutral taste and odor

### Technical data security-bite blue

processing time 30 sek.
Residence time in the mouth 60 sek.

### Assortment

26 pieces security-bite blue

2 x 50 ml cartridges 12 mixing cannulas, pink

12 contouring tips, wide

REF 580 0002 0



security-bite blue is placed on the arch.



The patient closes the jaws in the occlusion position until the material has hardened.

### Accessories:



Dosing device 1 piece REF 320 0044 0 Mixing cannulas size 1, blue 12 pieces REF 320 0045 0

contouring tip wide, 40 pieces REF 580 0002 3

# Bite blocks




Prefabricated wax bite blocks in the shape of rods and jaws and in 4 different degrees of hardness.



medium, red, UJ/LJ 74 pieces REF 430 0022 0 medium, red, UJ 74 pieces REF 430 0020 0 medium, red, LJ 74 pieces REF 430 0021 0



soft, pink, UJ/LJ 74 pieces REF 430 0027 0 soft, pink, UJ 74 pieces REF 430 0025 0 soft, pink, ∐ 74 pieces REF 430 0026 0



The basal profile of the bite blocks allows time-saving adaptation to the base.



REF 430 0017 0 hard, yellow, UJ 74 pieces REF 430 0015 0 hard, yellow, ∐ 74 pieces REF 430 0016 0



hard, yellow, UJ/LJ 74 pieces

super-hard, white, UJ/LJ 74 pieces REF 430 0012 0 super-hard, white, UJ 74 pieces REF 430 0010 0 super-hard, white, ⊔ 74 pieces REF 430 0011 0



Bite blocks bw rods medium, red 104 pieces 14 x 8 x140 mm REF 430 0023 0



104 pieces 14 x 8 x140 mm REF 430 0028 0

soft, pink



hard, yellow 104 pieces 14 x 8 x140 mm REF 430 0018 0



super-hard, white 104 pieces 14 x 8 x140 mm REF 430 0013 0

The height and the width of prefabricated bite blocks are suitable for the use on partial dentures.



Prepare situation model in the usual way.



The basal profile of the bite block simplifies adapting on the base



The consistency of the bite blocks allows simple reduction of the height and width using the wax



The prefabricated wax bite block can be easily integrated.



Since additional application of wax to the buccal and lingual area is no longer required, a considerable amount of time and material can be saved.



The high stability and functional processing of the bite blocks ensures precise bite-taking.



# Diagnostic wax up/Set-up

# Aesthetic wax veneers





### Aesthetic wax veneers

Prefabricated wax veneers enable simple working and help to save time. The tooth shade of the wax veneers enable the patient to get an impression of the final restoration.

Assortment 240 pieces Aesthetic wax veneers (24 different shapes with 10 individual units each) REF C13 2401 0

# Life-Color-Wax



Tooth-colored wax in two consistencies. Particularly low-shrinkage wax especially developed for the waxing-up technique according to M. A. Polz; hence perfectly suited for visualization.

### Life-Color-Wax

25 g

dentine color, medium REF 510 0080 0 dentine color, hard REF 510 0081 0

100 g

dentine color, medium REF 510 0079 0 dentine color, hard REF 510 0078 0



Precise application and superior scraping properties are the distinctive features of this wax.



# beauty setup



Dentine-colored wax developed for the visio.lign veneering system for fixation of the veneers for the esthetic try-in. Two different colors that can be mixed cover the classic A–D range of shades.

beauty setup bright, 25g REF 430 0030 0

dark, 25g REF 430 0031 0



The veneer is supported with beauty setup wax and set up.



Modelling the set-up using modelling wax pink Standard.



# Set-up wax



For setting up and changing the position of acrylic teeth without heating.



Set-up wax asw 4 pink 220 g REF 430 0157 4



Set-up wax asw 5 pink 220 g REF 430 0152 0



Set-up wax asw 3 pink 220 g REF 430 0151 0



Three different sizes of the pink set-up wax allow the individual use.



Thanks to its consistency the set-up wax can be perfectly processed without being heated.





The set-up wax allows quick adapting on the base plate.



No additional wax is required for flushing of the set-up wax.



Due to the adhesive capacity of the set-up wax, acrylic teeth are fixed prior to waxing on.



Even after waxing on, acrylic teeth can be brought into any individual position.

# Modelling wax pink Standard



Modelling waxes in sheets are used for a large number of applications in denture work.

Modelling wax pink Standard sheets.

Two thicknesses and three different qualities provide the technician with individual processing options.

Sheet thickness 1.25 mm quantity 1000 g 75 x 150 x 1.25 mm soft, pink medium, pink hard, pink

REF 430 0164 3 REF 430 0164 2 REF 430 0164 1

Sheet thickness 1.50 mm quantity 1000 g 75 x 150 x 1.5 mm soft, pink medium, pink hard, pink

REF 430 0164 6 REF 430 0164 5 REF 430 0164 4



Due to the particular stability of the pink modelling wax sheets, sufficient stability for the base plates is provided.



By rolling up the wax sheets and waxing them to the base plate, acrylic teeth can be set up immediately.



Bite blocks can be easily produced by rolling up and kneading this modelling wax.



The structure of the pink modelling wax sheets allows easy blocking-out for individual trays.



Even during extended try-in, the original stability of this modelling wax is maintained.



# Conventional implant planning

### X-resin OPG



### X-resin OPG

A special radiopaque resin matched with the requirements of the OPG. In these cases a highly radiopaque material is required to indicate clearly recognizable structures, which facilitates and accelerates the analysis of the X-rays. The cartridge resin can be applied in a precise and economic way to save material and reduce the working time.



### Assortment

- 14 pieces
- 1 cartridge 50 ml
- 1 X-connector 10 ml
- 12 Mixing cannulas size 1, blue

REF 540 0115 8

### Accessories:



Dosing device 1 piece REF 320 0044 0 Mixing cannulas size 1, blue 12 pieces REF 320 0045 0

### OPG-ruler



### **OPG-ruler**

transfers the original scale of 1.25:1 from the OPG to the implant.

The ruler is made of sterilizable material and can be easily cleaned. OPG radiographs can be measured in the scale of 1.25:1.

Converting is no longer required. Thanks to the OPG-ruler, X-ray templates based on implants can be omitted.

# OPG-ruler 1 piece REF 580 0127 6

# SKY X-ray template




### X-ray template

The suitable X-ray template is put on the X-ray and it can be immediately recognized which implant can be placed. The selection of the suitable implant is facilitated, accelerated and supported.



X-ray template X-ray template SKY whiteSKY Scale 1:1 Scale 1:1 **REF SKY-MS01** REF SKYMS01C Scale 1,12:1 Scale 1,12:1 REF SKY-MS12 **REF SKYMS12C** Scale 1,26:1 Scale 1,26:1 **REF SKY-MS26 REF SKYMS26C** 



The X-ray templates for SKY, blueSKY and whiteSKY are placed on the X-ray to see at a glance which implant diameter and implant length are required.



The suitable implant could be placed thanks to the X-ray template. Photos: Dr. Alius / Dr. Gresskowski, Nürnberg



# Navigated implant planning



The initial situation
Planning model fabricated using class III stone.

### Fabricating an x-ray template

The x-ray template is fabricated to visualize the diagnostic set-up/esthetic set-up of the teeth in the x-ray picture (CT/DVT). Subsequent prosthetic/navigated implant planning is based on the scan result.



The esthetic set-up

was made with the neo.lign denture teeth from the visio.lign veneering system.



Preparing a matrix

Use Putty soft matrix silicone.



Filling up with X-resin:

Inject X-resin into the matrix, place on the planning material and allow to harden.



Repositioning the X-resin arch

The X-resin arch is removed from the matrix, finished and fixed on the model.



### Coating with crystal-clear resin

The X-resin arch is coated with powder resin to ensure stability of the x-ray template. Then it is finished and polished.

Materials used

mini¹SKYl page 54 X-resin page 55



# Navigated implant planning

### mini<sup>1</sup>SKY



### mini¹SKY Implantat

The mini<sup>1</sup>SKY implants provide the basis to transfer and carry out implant placement based on 3D planning in a minimally invasive and reliable manner. Increased wearing comfort is accomplished until the final restoration is inserted. The mini<sup>1</sup>SKY Implantate ensure reliable and precise fixation of the scan and drilling templates. As a result, the implants are accurately placed in the planned positions.



Initial situation - jaw before treatment.



The required depth in the bone was drilled using the Twistdrill from the SKY OP-Tray .



The the mini<sup>1</sup>SKY implants were screwed in.

Length 6 mm, REF m1SKYL06

Length 10 mm, REF m1SKYL10



The SKY-OR55 O-ring is placed on the mini1 SKY implant.



Then the impression is

### Accessories:



miniSKY OP-Tray Clear arrangement of the instruments and drills facilitates the surgical procedure. **REF MSKYOT41** 



brecision implant heavy 380 ml impression material, blue 5 dynamic mixers 1 bayonet ring, blue REF 580 BH38 0



brecision implant light 2 x 50 ml cartridges, impression material, orange 10 Mixing cannulas 10 Intra-oral tips REF 580 BL05 0



breciform D impression tray - for single use Starter set 10 upper and lower trays each Sizes S, M, L and XL 10 breciform D stops, triangular shape 10 breciform D stops, bar shape REF 580 UOTS S



0-ring 6 pieces **REF SKY-OR55** 



### X-resin

### Clear vision - for all concerned



### Laboratory

Using x-ray templates and reference pins, the relevant prosthetic areas can be visualized in the X-ray picture or in the planning program (software).



### Dentis

The dentist is enabled to integrate additional prosthetic information in his treatment plan.



### **Patient**

Implant planning can be outlined to the patient using easily understandable illustrations/pictures and graphics. This way visualization of the necessary treatment steps for patients and their families is achieved.

Patients who understand the therapy will be more willing to spend the money that is necessary.

### X-resin

# Radiopaque paste material for diagnoses, preoperative planning and prosthetically oriented planning

To ensure long-term success of the implants, exact placement of implants is essential in order to protect them against overloading. Consequently, prosthetic planning should be the first step in the implant planning process. The x-ray picture needs to illustrate prosthetic planning for later implant placement.

The X-resin product family includes radiopaque resins which can be used with every standard planning system to achieve prosthetically oriented implant planning in a simple, safe and fast manner.

- prosthetic planning is carried out on the basis of a situation model
- prosthetic planning is transferred to the planning model with radiopaque teeth
- radiopaque teeth are made from X-resin
- the arch made from X-resin is reinforced using transparent resin in order to be able to provide sufficient stability for the drill sleeves.





X-resin CT, DVT 50 ml cartridge REF 540 0116 8

### Assortment

- 14 pieces
- 1 cartridge 50 ml
- 1 X-connector
- 12 Mixing cannulas size 1, blue REF 540 0115 9

X-resin CT DVT has been especially matched with the requirement of CT and DVTs. In these cases a material with reduced radiopacity is required to avoid excess radiation and to recognize clear structures.





X-connector Bonding agent 10 ml REF 540 0116 0

### Fabrication of x-ray templates with X-resin CT, DVT



The planning model with the mini1 SKY implants is made from a class III stone.



The esthetic set-up is prepared using the neo. lign denture teeth from the visio.lign system.



The set-up is fixed with Putty soft matrix silicone.



The mini1 SKY planning matrixes are placed on the mini1 SKY implants and undercuts are blocked out with wax



The Putty soft matrix is filled with X-resin and then placed on the model



The teeth filled with X-resin are fixed on the model.



The X-resin arch is coated with crystal-clear resin powder.



The finished scan template.



In this case the 3D resin was injected around the X-resin arch and a combination of a scan and drilling templated was fabricated.



In der 3D picture, the scannned teeth are located in the predetermined positions and prosthetically oriented implant planning is carried out.



# Navigated implant planning

### X-resin

### Accessories:



neo.lign denture teeth anterior and posterior are the consequent supplements within the visio.lign veneering system as far as shade, design and layering are concerned.



brecision Putty soft matrix silicone Hardness 65 Shore A 250 ml base (grey) 250 ml catalyst (white) 2 measuring spoons REF 580 0002 4



mini<sup>1</sup>SKY planning matrix 3 pieces REF m1SKYPLM



mini<sup>1</sup>SKY laboratory analog titanium REF m1SKYXIA



Dentaplast KFO

1000 ml Liquid

100 g Powder REF 540 0018 3 100 ml Liquid 500 g Powder 500 ml Liquid 1000 g Powder

REF 540 0018 4 REF 540 0018 5 REF 540 0018 6 REF 540 0018 7

REF 540 0018 8



Dosing device 1 piece REF 320 0044 0

Mixing cannulas size 1, blue 12 pieces REF 320 0045 0



Haptosil D silicone for keys Hardness 90 Shore A Components A and B each 1300 g REF 540 0118 0



### X-resin flow





X-resin flow is perfectly visible after the scan (DT/

### Radiopaque silicone varnish for diagnoses, preoperative planning and prosthetically oriented planning

X-resin flow allows to fabricate an x-ray template from the existing denture in just a few minutes.

This radiopaque silicone varnish is applied to the teeth of a full denture and thinly spread using a disposable brush.

Allow the material to harden for a short moment – that's it.

The contours of the teeth can be clearly recognized in the x-ray picture (DVT or CT) and thus prosthetically oriented alignment of the implants in the implant planning software is enabled. This facilitates planning considerably. If X-resin flow is applied to the basal area, the thickness of the mucosa can be clearly determined in the x-ray picture. The silicone varnish can be easily removed after the taking the x-ray.





X-resin flow
X-resin flow is
supplied in the
double-mix cartridge
to ensure a perfect
mixing result and fast
application.
50 ml

REF 580 0115 9

### Changing the existing denture into a scan template.



X-resin flow is applied to the dry and clean denture.



X-resin flow silicone varnish is evenly spread using a brush.



X-resin flow is applied to the basal area, the denture is inserted and pressure is exerted to achieve uniform distribution of the silicone varnish. This way the mucosal thickness can be determined.



The contours of the radiopaque silicone varnish X-resin flow can be clearly recognized.





Dosing device 1 piece REF 320 0044 0

Mixing cannulas, size 1, blue 12 pieces REF 320 0045 0



# Fabrication of drilling templates

### 3D-resin

### Perfect View - For all concerned



### Laborator

Reliabe and simple processing of 3D-resin enables the laboratory to prepare precise drilling templates in an uncomplicated manner.



### Dentist

Highly accurate guidance of the drill is ensured by the use of 3D-resin drilling templates. Slipping on the ridge is avoided and hence the implant is placed precisely in accordance with the criteria.



### **Patient**

Overloading is avoided thanks to careful placement of the implants. Consequently, reliable and comfortable wearing of the implant-supported restoration is ensured for the patient.

# Special resin for the fabrication of drilling templates in implant dentistry.

Drilling templates are indispensable tools for accurate placement of implants in implant dentistry. One of the key prerequisites during surgery is the sterilization of the surgical instruments which naturally also includes the drilling template. 3D-resin was developed especially for this purpose can be sterilized in the autoclave.

- Self-curing resin in the cartridge
- Heat-resistant and dimensionally stable up to 138° C
- Easy to use and to process
- Short setting time

### 3D-resin

is a two-component, self-curing special resin supplied in a cartridge which retains dimensional stability if it is heated up to 138° for a short period. 3D-resin is used to fabricate drilling templates in implant dentistry.



3Dresin

3D-resin 50 ml cartridge REF 540 0116 9

### Assortment

- 14 pieces
- 1 cartridge 50 ml
- 1 X-connector 10 ml
- 12 Mixing cannulas size 1, blue

REF 540 0116 4



**M** connector

X-connector Bonding agent 10 ml REF 540 0116 0



The planning model is insulated with Isoplast ip.



The 3D-resin drilling template material is applied to the model.



A piece of cling film is placed on and the material is evenly spread with the fingers.



The drilling template is finished in the usual way.



The model is placed on the transfer table, the coordinates are set and the holes are drilled.



The transfer pin is used to insert the drill sleeves into the drill holes; then SERACOLL UV is used to glue the sleeves.



# 3D-resin

### Accessories:



Isoplast ip 750 ml REF 540 0101 9



Dosing device 1 piece REF 320 0044 0

Mixing cannulas size 1, blue 12 pieces REF 320 0045 0



SKYplanX drilling sleeves 5 pieces each

Ø 2.39 REF SXBHI239

Ø 3.09 **REF SXBHI309** 

Ø 3.33 REF SXBHI333

Ø 3.59 **REF SXBHI359** Ø 3.83 **REF SXBHI383** 

Ø 4.09 REF SXBHI409

Ø 4.33 REF SXBHI433



SERACOLL UV light-curing wax adhesive 2 x 3 ml 2 dosing dishes REF 540 0115 1



SKYplanX Guiding sleeve 5 pieces each Ø 4.53 REF SXBHM453 Ø 5.55 REF SXBHM555



SKYplanX Sleeve drill 1 piece each for sleeve Ø 4.53 REF SplanX47 for sleeve Ø 5.55 REF SplanX55



SKYplanX Sleeve seating tool 1 piece each Ø 2.39 REF SplanX52 Ø 4.53 REF SplanX46 Ø 5.55 REF SplanX51 SKY Analogs REF SplanX53

# **NEW OPPORTUNITIES FOR SUCCESS**

# FURTHER TRAINING AT AND WITH THE BREDENT GROUP ACADEMY



We strive to be amongst the best. In this spirit, we have reviewed our course programme and the results are now being introduced under the bredent group academy banner.

The bredent group academy course programme is aligned to the strategic direction of the bredent group. On the one hand, we have the part that is in line with the market, oriented towards tailored solutions to meet the needs of the market. Such solutions have been developed according to the specific segment. By contrast, there is the market-oriented segment. This is geared towards a changing market, such as economic cycles and trends.

However, it is far more important to recognise users' wishes from the outset. This means being intuitive to what users want and need, and knowing what the future holds - in terms of materials, systems and treatment concepts. Only this way can we as a manufacturer ensure that our company progresses and understands the associated obligations to inform, support and help you, as our client, on the path to commercial success.

Further training is a key component of this duty, since the cornerstone for new technology and a smooth transfer of knowledge are being laid. The instructors

are the life of this course – all of them experts in their field – equipped with the vision and the necessary drive to develop innovations and put them into practice. This programme unites these visions and appeals to those who never cease to develop themselves further.

The individual development and optimisation of your practice, your laboratory and your marketing has a decisive effect on realising success in your practice or laboratory, your income and therefore your abilities, ideas, wishes and dreams.

A new range of events in the areas of patient communication, hygiene, conflict management, practice marketing, data protection and social media, in addition to the exciting beginners' and improvers' courses in implantology and dental technology, can help you become even more successful. Not only is further training at the forefront of your professional qualifications, but so too is positioning your business as a successful brand on the market.

Join us and strive to be amongst the best – for the benefit of you and your patients!

Find out more about our course programme at www.bredent.com.



The temporary restoration is the first restoration the patient comes into contact with. Therefore the material used for it should have excellent quality to be accepted by the patient.

These temporary materials form the basis of the future permanent restoration.

Suitable materials for the fabrication of interim dentures for hypersensitive patients can be found in chapter 10.





Photo1 Dr. Frank und Dr. Steffen Kistler, Dr. Georg Bayer, Landsberg am Lech



The interim denture	
System steps	63
Qu-base UV	64
Thermoplastic clasp	60
The temporary	
System steps	
top.lign breformance	68
The milled temporary	
System steps	7
breCAM.resin milling blanks	72
breCAM.cutter	
Exakto-Rock S	7
Fluid-Rock	7

# **IMPLANT PROSTHETICS**

# ONE-STOP INTERDISCIPLINARY COMPETENCE – WHAT THE BREDENT GROUP STANDS FOR



In 1974, the approach of the bredent group was to ensure that they manufactured more than just a few products for lasting success for customers.

With comprehensive products, system solutions and treatment concepts for dental technology and dentistry, which are interconnected in a professionally competent manner, the bredent group is amongst the few companies that are held in high esteem equally by both groups.

In close cooperation with reference laboratories and practices, well-matched, patient-centred products, system solutions and treatment concepts are developed, which can be implemented efficiently as well as rapidly in the laboratory and the practice.

The positive effect of ground-breaking development!

# The interim denture

Thanks to the use of Qu-base UV, laborious and

and fast working.

time-consuming fabrication of interim dentures is now a thing of the past. Qu-base UV allows efficient



The initial situation
A class III stone was
used to fabricate the
model. The undercuts
were blocked out with
Biotec blocking out wax,
insulated with Isoplast Ip
and then the clasps were
attached.



Qu-base UV plate material is adapted. The high stability of Qu-base UV plate material enables convenient working and helps to save time.



Setting up the teeth neo.lign teeth from the visio.lign system are used for setting up the teeth. Polymerization is carried out in the bre.Lux Power unit in a way to save time.



Adding up gingiva Qu-resin pink in the double-mix cartridge enables controlled and efficient application.



Processing and polishing A wide range of tungsten carbide tools and polishing brushes ensure fast processing and polishing.


Materials used

Qu-base UVI Thermoplastic clasp page 64 page 66



### The interim denture

### Qu-base UV






is the new, light-curing plate material for fast fabrication of temporary dentures and drilling tem-plates.

- 50 min. time advantage compared to conventional fabrication of temporary dentures
- prefabricated plates in upper and lower jaw shapes eases the make of temporary dentures and drilling templates
- no wax set-up, no silicone key and no pressure pot required anymore, saves material and money
- excess material is used for teeth setup hence allowing for economical work



Qu-base UV 7 x UJ / 7 x LJ REF 540 0117 1



Qu-base UV 14 x UJ REF 540 0117 2



Qu-base UV 14 x LJ REF 540 0117 3

### Indication



Clasp retained temporary



Implant retained temporary dentures



Addition of temporary



Repairing of a temporary dentures



Drilling-template manufacturing

# Fabrication of an interim denture with Qu-base UV



The prefabricated Qubase UV upper plate is pressed on the model (insulated with Isoplast ip) and cut out. The thermoplastic clasps are attached.



Excess of Qu-base UV plate is placed on the ridge to set up the teeth. A layer thickness of 3 mm must not be exceeded.



Prior to the set-up, the neo.lign denture teeth are sandblasted with 110 μm aluminium oxide, a thin coat of visio.link or Qu-connector is applied and activated in the bre. lux Power unit (light curing unit) for 90 sec.



Once the teeth have been set up, the model is placed into the bre.lux Power unit (light curing unit) and polymerized. To keep polymerization shrinkage on a minimum level, the following procedure is recommended:

- 1. 90 sec. with 20% power
- 2. 90 sec. with 70% power 3. 90 sec. with 100% power
- 4. Remove the denture form the model and polymerize the bottom side for 90 sec. with 100 % power.



Qu-resin pink is used to apply the missing gingival proportions. The dispersive layer of Qu-base UV must not be removed before applying Qu-resin or Qu-connector or visio.link must be applied thinly and polymerized in the light curing unit for 90 sec before Qu-resin



The use of tooth-colored thermoplastic clasps leads to a better esthetic appearance.



### Qu-base UV

### Fabrication of a drilling template with Qu-base UV



Model with mini<sup>1</sup>SKY lab implants for retaining the drilling template.



Attach metal matrix housing, block out and insulate with Isoplast ip.



Press Qu-base UV against the metal matrixes and cure for 90 sec.



Adapt UV plate, cut out and replace into SKY5X transfer table.



Perforate the plate at the point where the drill sleeve is to be placed. Coat with Qu-base UV and carry out light curing (local) to fix the drill sleeve. This step must be repeated for additional drill sleeves.



Once all drill sleeves have been attached, the model is placed into the bre.lux Power unit (light curing unit) and polymerized. To keep polymerization shrinkage on a minimum level, the following procedure is recommended:

- 1. 90 sec. with 20 % power
- 2. 90 sec. with 70 % power
- 3. 90 sec. with 100 % power
- 4. Remove the denture form the model and polymerize the bottom side for 90 sec. with 100% power.

### Asseccories:



Biotec blocking out wax Special wax for blocking out undercuts. No discoloration on the plaster model after boiling out or cleaning with steam. 28 g, pink REF 510 0061 5



Isoplast ip alginate based and insulates plaster against resin whilst creating a highly lustrous resin surface. 750 ml REF 540 0101 9



neo.lign anterior and posterior are the logical addition to the visio. lign veneering system (see chapter 9).



### bre.Lux Power Unit

- LED technology ensures extended service life
- 370 500 nm covers the entire wavelength range and the hand lamp is also suitable for the same wavelength range
- fixation/prepolymerization/intermediate and final polymerization with a single unit directly at the working place
- high capacity for 2 models, perfectly and evenly illuminated

REF 140 0097 0



### Ou-resin

is a quick-setting, self-curing denture repair resin in two colors (pink / dentin). It is suitable for direct and indirect use.

## Qu-resin pink

Assortment 14-pieces 1 Cartridge 50 ml

- 1 Qu-connector 10 ml 12 self mixing tips
- Size 1, blue

REF 540 0116 1

### Qu-resin dentin

Assortment 14-pieces 1 Cartridge 50 ml

- 1 Qu-connector 10 ml
- 12 self mixing tips

Size 1, blue REF 540 0116 2

### Qu-resin pink 50 ml Cartridge REF 540 0116 5

Qu-resin dentin 50 ml Cartridge REF 540 0116 6

Tools for perfect surface processing can be found in chapters C and D.



# The interim denture

# Thermoplastic clasp

### 1. New fabrication of a clasp denture



**Invisible clasp in different tooth colors.**Aesthetically appealing clasp for a happy smile.



Thermoplastic clasp, left (2./4. quadrant) Package cont. 8 pieces REF 430 00LA 2 REF 430 00LA 3 REF 430 00LB 2 REF 430 00LB 3



Thermoplastic clasp, right (1./3. quadrant) Package cont. 8 pieces REF 430 00RA 2 REF 430 00RA 3 REF 430 00RB 2 REF 430 00RB 3

Available shades: A2, A3, B2, B3



In the visible area metal clasps are not cast together with the CoCr structure. Attach metal clasp only in areas that do not affect the esthetic appearance.



Select the suitable color of the prefabricated clasp profile and carry out coarse grinding using a bur.

Fix the clasp profile in Translock, heat



and press the clasp against the tooth mould. Grind the adapted clasp to obtain the final shape.



Each clasp design can be produced individually.

### 2. Repair or extension of a clasp



Time-saving fabrication of an esthetic clasp.



Invisible, thanks to a matching tooth color – this way the patient's self esteem is increased.



Thermo-Pen
1 piece
REF 110 0147 0

### Accessories:




Transblock 500 g **REF 540 0035 0** 



# The temporary



The initial situation Fabrication of the model with a class III stone. The gingival mask is prepared with Multisil-Mask soft.





Setting up the teeth Before setting up the novo.lign veneers from the visio.lign system, the model is insulated with Isoplast ip. The titanium prosthetic caps are shortened correspondingly.



A haptosil D matrix is used for fixation of the set-up.



Covering the prosthetic caps
Silicone tubes are attached as spacers to all prosthetic caps except of one in order to lute them in situ without any stress later on. Then a suitable shade of top.lign

breformance is used to fill the matrix.



Processing and polishing A wide range of tungsten carbide tools and polishing brushes ensure fast finishing and polishing.



Bonding in the patient's mouth

Qu-resin dentin is used for stress-free bonding in the patient's mouth.

Materials used

top.lign breformance

page 68

# top.lign breformance



top.lign breformance is a tooth-colored, hot- and cold-curing (top.lign hot / top.lign cold) crown and bridge resin for temporary restorations. top. lign cold features good flow properties and is distinguished by its elasticity, facture strength and shade stability.






top.lign breformance liquid heat 50 ml REF bre Ihq0 1



top.lign breformance REF bre lcq0 1 100 ml REF bre Icq0 2



1 x 50 ml top.lign breformance liquid heat 1 x 50 ml top.lign breformance liquid cold separating liquid 25 g Polymer each REF bre set0 1



Shades according to VITA classical / Vita is a registered trademark of Vita Zahnfabrik, Bad Säckingen

### Clinical use in implantoloy with SKY fast & fixed and top.lign cold



Model with screwed-on SKY fast & fixed prosthetic copings made of titanium.



novo.lign veneers were used for the set-up. The set-up is retained using the haptosil D matrix.



The novo.lign veneers were sandblasted with 110 um aluminiun oxide and placed in the matrix.



Before attaching the matrix, silicone tubes are placed on all prosthetic caps except of one.



The bridge was lined and the non-coated prosthetic cap has been fixed in place (polymerized).



The bridge was finished and space was created for direct bonding of the remaining prosthetic copings.



The remaining prosthetic caps are screwed on the abutments and bonded stress-free in the patient's mouth using Qu-resin dentin.



li is essential to polish the entire bridge structure to high gloss to avoid accumlation of plaque.



The completed temporary bridge in the patient's

Photos: Dr. Frank and Dr. Steffen Kistler, Dr. Georg Bayer and DT Stephan Adler, Lands-



# top.lign breformance

top.lign breformance hot

Easy to use heat-curing resin for temporary crown and bridge technology. Monomer reduced production is possible with breformance LiquidHeatCuring.




Deburr any sharp edges. Use Isoplast ip (REF 540 0101 9) to insulate the plaster while it is still warm. Make sure that no puddles form behind the dies. Isoplast ip facilitates removal after the polymerization process.



Embed the model so that the cover of the mould can be easily separated.



Heat the mould with boiling water to facilitate opening. To prevent the dies from fracturing, open the two halves of the mould carefully. After opening, boil out the remaining wax.



The expanded resin is now inserted into the mould taking care to avoid bubble formation. More resin is inserted than necessary so that the resin is compressed during pressing.



For the trial pressing a foil is placed between the two mould halves. After the trial pressing, individual treatment can begin.



The enamel area is cut back using a sharp instrument. The pre-reacted enamel resin is now added and a new trial pressing is made. After checking the trial pressing, the mould is closed and heated.

approx. 6 minutes.



The completed bridge with layered incisal edge. If the shade needs to be modified, individualization can be easily carried out with top.lign cold.

### Accessories



SKY fast & fixed Laboratory analog, titanium REF SKYFFOLA



SKY fast & fixed Prosthetic cap, titanium REF SKYFFPKT In bridge restorations with SKY fast & fixed, the prosthetic cap is used as a basis for the temporary restoration to enable fast and simple fabrication. The prosthetic cap – transversal screw retention – can also be used for the temporary restoration. A silicone tube is enclosed to each SKY fast & fixed titanium prosthetic cap.



Set 2 x 50 ml Multisil-Mask soft 24 mixing cannulas 10 ml Multisil-Sep REF 540 0104 1



Isoplast ip 750 ml REF 540 0101 9



Set-up waxl asw 3 pink, 220 g REF 430 0151 0 Set-up waxl asw 4 REF 430 0157 4 Set-up waxl asw 5 REF 430 0152 0



Modelling wax pink standard Sheet thickness 1.50 mm, 1000 g 75 x 150 x 1.5 mm soft, pink REF 430 0164 6 medium, pink REF 430 0164 5 hard, pink

REF 430 0164 4



neo.lign anterior and posterior are the logical addition to the visio. lign veneering system. (see chapter 9)



haptosil D
Addition-curing
kneading silicone with
a Shore A hardness of
90 for the fabrication
of stable and exact
keys and models. For
reliable fixation of the
set-up and the teeth
in the key.
Component A and B,
1300 g each
REF 540 0118 0

# The temporary

### Accessories:



### Qu-resin

is a quick-setting, self-curing denture repair resin in two colors (pink / dentin). It is suitable for direct and indirect use.

### Qu-resin pink 1 cartridge 50 ml 1 Qu-connector 10 ml 12 self mixing tips Size 1, blue REF 540 0116 1

### Qu-resin dentin Assortment 14 pieces Assortment 14 pieces 1 cartridge 50 ml 1 Qu-connector 10 ml 12 self mixing tips

Size 1, blue

REF 540 0116 2

Qu-resin pink 50 ml cartridge REF 540 0116 5

Qu-resin dentin 50 ml cartridge REF 540 0116 6

Tools for perfect surface processing can be found in chapters C and D.

# The milled temporary



The initial situation Sawcut model fabricated with Exakto-Rock S

breCAM.resin are PMMA milling blanks for the fabrication of high-quality long-term temporaries in the basic colors A, B and C as well as transparent for splints. The good mechanical values also allow the use of the milled temporaries for fully anatomical restorations.



**Der Scan**The model is scanned.



The design
The bridge is designed on the computer - in this case, fully anatomical.



The milling process
The milling system mills
the bridge from the
breCAM.resin milling
blank.



The result
The fully anatomical,
milled bridge on the
model.

#### Materials used

breCAM.resin milling blanks	page 72
breCAM.cutter	page 73
Exakto-Rock S	page 74
Fluid-Rock	page 74



# The milled temporary

# breCAM.resin milling blanks

**breCAM-Blanks for the CNC milling technique** High-quality basic materials provide ideal preconditions for high-quality restorations.

- Blanks with a diameter of 98 mm, suitable for all standard milling machines
- The integrated step ensures reliable and simple fixation of the blanks in the holder
- Different heights offer more flexibility, hence also suitable for higher implant restorations
- Different thicknesses of the blanks guarantee efficient and faster milling
- High-quality resin blanks for temporary restorations
- Wax blanks with particularly high edge stability allow high milling speeds



		16 mm	20 mm	24 mm
<b>breCAM.resin A</b> l 1 blank	REF	540 0201 0	540 0201 1	540 0201 2
<b>breCAM.resin B</b> 1 blank	REF	540 0201 3	540 0201 4	540 0201 5
breCAM.resin C 1 blank	REF	540 0201 6	540 0201 7	540 0201 8
breCAM.resin transparent 1 blank	REF	540 0201 9	540 0202 0	540 0202 1



PMMA blanks for high-quality long-term temporaries in the basic colors A, B and C and as transparent blanks for splints. Thanks to a special manufacturing process, the blanks feature superior material properties, which result in increased stability of the restoration.

Technical data - breCAM.resin:

Modulus of elasticity: 2760 MPa Flexural strength: 114 MPa Elongation at break: 7 %



The integrated step of the blanks enables safe fixation in the holder. The different thicknesses of the blanks provide high flexibility in the use and hence the blanks can also be used for higher implant restorations.



#### breCAM.resin

The use of breCAM.resin as substructure material – veneered with visio.lign veneers – enables economic fabrication of long-term temporaries. Economic fabrication of anatomical temporary restorations is also possible. The high surface density of breCAM.resin guarantees extremely low plaque affinity.

### breCAM.cutter



#### The breCAM.cutter

was developed especially for dry processing of PMMA and PEEK in CNC milling systems. Thanks to its patented geometry, it can also be used for milling thermoplastic materials, which have a tendency to rapid smearing and clogging of the tool, without cooling with water.

specially desigend for dry-milling purpose, usable for wet milling as well. Cutter stays clear, even by milling materials with the tendancy clogging the bur



- Patented blade geometry, avoiding overheating
- Single blade cutter
- Usable for popular machines
- Enhanced cutter for dry-milling ( thermoplastic material, PEEK,wax)

REF	Manufacturer	Туре	Shank diameter	 Diameter of working	Total length	— Working section length
breCAMX47	Wieland, Imes Icore, Coritec, 340i, 450i, 440i	Radius cutter	3.0	1.0	38.2	15.0
breCAMX48	Wieland, Imes Icore, Coritec, 340i, 450i, 440i	Radius cutter	3.0	2.0	38.2	15.0
breCAMX49	3M Espe, LAVA Form System 400, 500, Charly Robot	Radius cutter	3.0	1.0	38.0	15.0
breCAMX50	3M Espe, LAVA Form System 400, 500, Charly Robot	Radius cutter	3.0	2.0	38.0	15.0
breCAMX53	Roland DWX 40, DWX 50, Calidia, TDS, DMG U-Serie, Yenadent, Orgien	Radius cutter	4.0	1.0	50.0	15.0
breCAMX54	Roland DWX 40, DWX 50, Calidia, TDS, DMG U-Serie, Yenadent, Orgien	Radius cutter	4.0	2.0	50.0	15.0
breCAMX67	VHF, FinoCAM, Jeneric Pentron, Schütz, Trendgold	Radius cutter	3.0	1.0	34.0	15.0
breCAMX69	VHF, FinoCAM, Jeneric Pentron, Schütz, Trendgold	Radius cutter	3.0	2.0	34.0	15.0

# The milled temporary

### Exakto-Rock S

Exakto-Rock S is a synthetic super-hard class IV stone with distinctive thixotropy and improved flow properties and is free from formaldehyde.



Reduced expansion is completed after 2 hours and is only 0.08 %. As a result, highly accurate impressions are enabled and precision-fit restorations are ensured. Moreover, thanks to optimized reflection of light, Exakto-Rock S is suitable for scanning and available in brown and ivory.

- fomaldehyde-free stone ensures safe processing and can be used for the fabrication of holistic dentures without any problem
- synthetic components ensure consistent quality and enable the fabrication of precision-fit models
- improved flow characteristics facilitate pouring of several impressions
- optimized reflection of light thanks to special dye pigments reduce reworking/adjustments in the CAD system



#### Color - brown:


1 x 2 kg REF 570 0SB5 2 5 x 2 kg REF 570 0SB5 1 10 x 2 kg REF 570 0SB5 0



#### Color - ivory:

1 x 2 kg REF 570 0SE5 2 5 x 2 kg REF 570 0SE5 1 10 x 2 kg REF 570 0SE5 0

#### Expansion of various other stones



### Fluid-Rock

Fluid-Rock is a smoothly flowing class IV super-hard stone to prepare bases of arches.



The light-blue color can be easily combined with all colors for the arch. The extended processing time allows to pour several bases at the same time. The thin consistency results in perfect flow characteristics and allows to obtain bubble-free models.

#### Color - blue:

1 x 2 kg REF 570 0FB5 2 5 x 2 kg REF 570 0FB5 1 10 x 2 kg REF 570 0FB5 0

#### Technical Data - Fluid-Rock

Color
Mixing ratio
Processing time

Setting time (Vicat time)

Comp. strength aft. 1 hr Comp. strength aft. 24 hrs Setting expansion 100 g / 25 ml distilled water approx. 6 min at 18° to 20° C approx. 11 min at 18° bis 20° C 48 N/mm<sup>2</sup> 55 N/mm<sup>2</sup>

blue

< 0.06 % (no further expansion after 2 hours)

Processing in the ecovac unit:

Vacuum level 1, mixing speed: 390 rpm

Tools for perfect surface processing can be found in chapters C and D.



In the area of digital manufacturing, it is important to achieve optimal and reproducible results using system components of the "digital workflow" that are attuned to each other.

This includes conventional CAD software, scanners and CNC machines and the corresponding grades of materials and the tools and milling strategies optimised for them.

It also encompasses conceptual design, the assembling of independent, open hardware and software components in accordance with an individual production design and the possibility to expand these further in the future.

All that means that the choice and decisions made for CAD/CAM systems today are subject to questions of system compatibility, open interfaces and a free choice of materials from multiple manufacturers for processing in the system.

In future, the communication accompanying an order, the involvement of physicians and a free exchange of suitable documentation will be increasingly important for planning and manufacturing. The "digital workflow" will consequently be supplemented by digital production and communication management.



CAELO system group CAD/CAM	
e.order	76
e.cad	77
e.scan	77
CAD/CAM blanks	
breCAM.wax milling blanks	78
breCAM.BioHPP milling blanks	79
breCAM.resin milling blanks	80
CAD/CAM mills	
breCAM.cutter	81
CAD/CAM plasters	
Exakto-Rock S	82
Examed flock S	02

### CAELO system group CAD/CAM

### www.caelo-dental.net

The digital system group CAD/CAM. Open system range of products including hardware and software components. Modular construction; independent, compatible with other companies' software and hardware for expansion.



CAELO is a digital product range that is constantly being developed and improved. All solutions from the CAELO system group generally work with open hardware and software products from other manufacturers on the basis of processing and output of open STL files.

To let you know about the latest developments in the CAELO family, the corresponding product descriptions are available on our website:

### www.caelo-dental.net



Is a dental management and communication application, installed locally or controllable via the web using a web browser for the most important functions.

e.order is the optimal and sensible addition to the "digital workflow" to coordinate communication and data exchange between all persons and parties active in the planning or manufacturing process. **REF CAE TOOO 0** 

#### e.order user groups

#### Physicians

and other specialists involved in further processing in the areas of dental applications, including research and development.

- For exchange/communication in the area of diagnostics and treatment
- To discuss current patient cases and seek other expert opinions where necessary
- To define and implement shared actions, developments, treatments and manufacturing processes (digital order slip)
- To monitor and guide these through the implementation process
- Interested parties can expand the diagnostic circle at any time and invite particular specialists/participants
- Documents and data prepared in the cloud can be used for evaluation, processing, change and addition
- Documents saved in the cloud (in accordance with the German Medical Products Act) (flexible office) allow access independent of location to all possible digital data formats and contents

of the system (communication tools)

#### Dental technicians

stay connected to their clients during the whole of the order's manufacturing process and can access centrally prepared and saved data and communication tools together or individually with the location-independent and system-independent software e.order.

# Milling centre/service providers and dental service receivers

are offered various functions by e.order to support the digital manufacturing process:

- Network with your clients and partners
- Communicate using independent, free software
- Digital data (STL, DICOM, PDF, WORD, TXT, etc.) pertaining to an order or case can be made available to the user groups.
- The manufacturing process can be monitored and guided from any place at any time, independent
- Brokerage function, choice of required service providers by manufacturing process (milling, polishing, laser sintering, stereolithography)
- Installation and allocation of an "order thread" for the purpose of monitoring the manufacturing process and documentation exchange of associated communications
- Using additional functions, e.order supports the digital workflow in the area of communication/ monitoring, through threads, push e-mail, SMS, Skype
- ".dcm-.stl Converter", converts a .dcm file into an open .stl format
- "Digital order slip" for quick digital placing of orders by your dentist, with accompanying digital documents as an attachment in the cloud or threads
- Integrated DICOM Viewer and STL Viewer
- Cloud-storing, central data storage in accordance with MPG, also available as web-based service



# CAELO system group CAD/CAM



"Just right" is the motto of our modular CAD system.



Optimally adjusted for our system, e.scan from bredent is our open STL scanner.

Determine exactly what you need. No more, no less!

Is your weak point in crowns, bridges, telescopic and 2-component restorations, or do you require additional functionality for implantology, attachments, bridges or the visio.lign veneer system?

Our business model is visionary: a modular CAD, with selectable or deselectable functionality and structural modules, a system that adapts itself to you and your needs and grows with them.

Naturally, e.cad has already been linked to e.order. As a logical consequence of its user-friendly operation, e.order controls the scanner (e.scan) and conveys the model scan directly to e.cad for further processing. **REF CAE COOO 0** 

An unbeatable price-performance relationship:

Dental scanner for the 3D scanning of physical plaster models. Our high-speed scanner offers accuracy and precision, both of which are required for full bridges and implant bridges. Full integration with e.scan software and e.order for the optimal "scan – construction – work" process.

REF CAE SOOO O



### breCAM.wax milling blanks

breCAM blanks for the CNC milling technique
High-quality basic materials provide perfect prerequisites for top-quality dental

- Blanks with a diameter of 98 mm, suitable for all standard milling machines
- The integrated step ensures reliable and simple fixation of the blanks in the holder
- Wax blanks with particularly high edge stability allow high milling speeds



20 mm

breCAM.waxl 2 blanks REF **540 0201 1** 

The breCAM.wax blanks are made from a special milling wax (micro-crystal-line hydrocarbon wax with paraffin waxes and polyethylene) which is easy to process. Because of this, the blank is suited to use in open CAM systems, for the digital manufacture of fully-anatomical or reduced frameworks, crowns and bridges for metal casting or ceramic pressing technology and the use of high implant work.

#### Technical specifications of breCAM.wax

Size: round 98.4 x 20 mm with cervical step

with cervical step 10 x 2 mm

Material: microcrystalline hydrocarbon wax with

Colour: green

hard paraffins and polyethylene sections Dropping point: 120°C, material burns without trace







#### breCAM.wax

The significant advantage of the cutting manufacture of wax structures for the casting technique lies in the fact that the wax blank has been melted and cooled in a controlled manner and therefore has homogeneous and very limited tension behaviour, entirely independent of the various framework and connector strengths or spans of the structure. The conventional wax technique of fusing, dipping and combining various waxes quickly leads to uncontrollable results.


Overheating destroys important catalysts of the wax and different wax temperatures and applied strengths can cause strong tensions in the structure of the wax, which may lead to poorly seated and poorly fitting casts.



# breCAM.BioHPP milling blanks

**breCAM blanks for the CNC milling technique** High-quality raw materials offer the ideal basis for high-quality dentures.

- Blanks with a diameter of 98 mm, suitable for all current milling machines
- The incorporated level offers secure and simple fixing of blanks in the retainer




		16 mm	20 mm	24 mm
breCAM.BioHPP	REF	540 0203 0	540 0203 1	540 0203 2
1 blank				

BioHPP is a high-performance polymer (PEEK) material with ceramic fillers. For more than 30 years, PEEK has been used in human medicine as an implant material.

breCAM.BioHPP has great potential as a framework material in dental prosthetics, is significantly cheaper than gold and is easier and better to process in dental laboratories than non-precious metals, titanium or ceramic.

The extraordinary elasticity of the material also comes into play here, as it is close to that of bone, and is therefore capable of compensating for the torsion of the bone, particularly in larger implant work.

For processing breCAM.BioHPP we expressly recommend our breCAM.cutter mill, which has been specially developed for dry and wet processing of thermoplastic materials.



#### Technical specifications of breCAM.BioHPP

E-modulus: 4000 MPa Flexural strength: 150 MPa

(no material failure)



16 mm 20 mm 24 mm

4 CAD/CAM



#### Processing breCAM.BioHPP

Thanks to its excellent stability, its optimal polishable properties and its low plaque affinity, BioHPP is particularly well suited for producing high-quality prosthetic restorations, whereby the aesthetically-pleasing white colour "White Shade" of the material is eminently suited to use in a tooth-coloured veneer, e.g. by means of the visio. lign veneer system from bredent.



### Processing reference:

Veneer technique visio.lign Mill, breCAM.cutter Chapter 9 page 79

For the perfect high gloss Plastic polishing set REF 350 0081 0

Chapter D

# breCAM.resin milling blanks

breCAM blanks for the CNC milling technique High-quality basic materials provide perfect prerequisites for top-quality dental

- Blanks with a diameter of 98 mm, suitable for all standard milling machines
- The integrated step ensures reliable and simple fixation of the blanks in the holder
- Different heights offer more flexibility, hence also suitable for higher implant restorations
- Different thicknesses of the blanks guarantee economical and faster milling
- High-quality resin blanks for temporary restorations




	16 mm	20 mm	24 mm
<b>breCAM.resin A</b> l 1 Blank	REF <b>540 0201 0</b>	540 0201 1	540 0201 2
<b>breCAM.resin B</b> 1 Blank	REF <b>540 0201 3</b>	540 0201 4	540 0201 5
<b>breCAM.resin C</b> 1 Blank	REF <b>540 0201 6</b>	540 0201 7	540 0201 8
breCAM.resin transparent 1 Blank	REF <b>540 0201 9</b>	540 0202 0	540 0202 1

The breCAM.resin blank is produced using polymethylmethacrylate in the colours A, B and C (corresponding to Vita-colour A2, Vita-colour A3/B2, Vita-colour C2/D3) and transparent, in the industrial thermoplastic manufacturing process.

The advantage compared to the standard chemoplastic PMMA blanks on the market is a denser networking of hydrocarbon chains which leads to a high fracture resistance and flexural strength.

Thermoplastic blanks therefore exhibit a higher plaque resistance, and are generally more durable and resistant to discolouration as a result of cigarettes, tea and coffee, for example.



#### Technical data - breCAM.resin:

Modulus of elasticity: 2760 MPa Flexural strength: 114 MPa Elongation at break: 7 %



16 mm 20 mm 24 mm



#### breCAM.resin

The use of breCAM.resin as substructure material – veneered with visio.lign veneers – enables economic fabrication of long-term temporaries. Economic fabrication of anatomical temporary restorations is also possible. The high surface density of breCAM. resin guarantees extremely low plaque affinity.

#### Accessories



Generation M chapter C

breCAM.resin is free of dibenzoyl peroxide and tertiary amines. Due to its material class and processing, breCAM.resin has an extremely limited residual monomer content is therefore highly biocompatible.

The transparent material is approved for medical use and can be used for splints, try-ins, tabletops, reduction caps or framework material for burning out for the casting technique.



### breCAM.cutter



#### The breCAM.cutter

was developed especially for dry processing of PMMA and PEEK in CNC milling systems. Thanks to its patented geometry, it can also be used for milling thermoplastic materials, which have a tendency to rapid smearing and clogging of the tool, without cooling with water.

specially desigend for dry-milling purpose, usable for wet milling as well. Cutter stays clear, even by milling materials with the tendancy clogging the bur



- Patented blade geometry, avoiding overheating
- Single blade cutter
- Usable for popular machines
- Enhanced cutter for dry-milling ( thermoplastic material, PEEK,wax)

REF	Manufacturer	Туре	Shank diameter	Diameter of working se	Total length	Working section length
breCAMX47	Wieland, Imes Icore, Coritec, 340i, 450i, 440i	Radius cutter	3.0	1.0	38.2	15.0
breCAMX48	Wieland, Imes Icore, Coritec, 340i, 450i, 440i	Radius cutter	3.0	2.0	38.2	15.0
breCAMX49	3M Espe, LAVA Form System 400, 500, Charly Robot	Radius cutter	3.0	1.0	38.0	15.0
breCAMX50	3M Espe, LAVA Form System 400, 500, Charly Robot	Radius cutter	3.0	2.0	38.0	15.0
breCAMX53	Roland DWX 40, DWX 50, Calidia, TDS, DMG U-Serie, Yenadent, Orgien	Radius cutter	4.0	1.0	50.0	15.0
breCAMX54	Roland DWX 40, DWX 50, Calidia, TDS, DMG U-Serie, Yenadent, Orgien	Radius cutter	4.0	2.0	50.0	15.0
breCAMX67	VHF, FinoCAM, Jeneric Pentron, Schütz, Trendgold	Radius cutter	3.0	1.0	34.0	15.0
breCAMX69	VHF, FinoCAM, Jeneric Pentron, Schütz, Trendgold	Radius cutter	3.0	2.0	34.0	15.0

with chip channel

# **CAD/CAM plasters**

### Exakto-Rock S

Exakto-Rock S is a synthetic super-hard class IV stone with distinctive thixotropy and improved flow properties and is free from formaldehyde.



Reduced expansion is completed after 2 hours and is only 0.08 %. As a result, highly accurate impressions are enabled and precision-fit restorations are ensured. Moreover, thanks to optimized reflection of light, Exakto-Rock S is suitable for scanning and available in crown and ivory.

- fomaldehyde-free stone ensures safe processing and can be used for the fabrication of holistic dentures without any problem
- synthetic components ensure consistent quality and enable the fabrication of precision-fit models
- improved flow characteristics facilitate pouring of several impressions
- optimized reflection of light thanks to special dye pigments reduce reworking/adjustments in the CAD system



#### Color brown:


1 x 2 kg REF 570 0SB5 2 5 x 2 kg REF 570 0SB5 1 10 x 2 kg REF 570 0SB5 0



#### Color ivory:

1 x 2 kg REF 570 0SE5 2 5 x 2 kg REF 570 0SE5 1 10 x 2 kg REF 570 0SE5 0

#### Expansion of various other stones





www.caelo-dental.net

### Framework manufacture

Customised framework manufacture requires the highest level of precision. This is facilitated by the use of high-quality modelling waxes, special waxes and successful casting technology in everyday application. The investment material for exact-fit precision casting reduces the need for repetitions and therefore saves time.



bı	redent casting technique	
	Manual	84
_		
C	asting channel restoration	
	System steps	
	Isobre wax insulating liquid	86
	Elaflex	86
	Visio-Dip	86
	Dipping wax	87
	Cervical disc	87
	Waxpool duo	88
	Wax knife bwm 3	89
	Quick Change	90
	Wax for outer copings	.91
	Biotec-Wax for outer copings	.91
	K2 exact	.91
	Standard Modelling Wax	.91
	Gecko	92
	Splendido	92
	KBI-Wax	92
	Life-Color-Wax	92
	Biotec modellling wax	93
	Cervical wax	
	Biotec-Cervical wax	93

Biotec milling wax	93
SERACOLL UV	
compoForm UV	
Biotec metal-ceramic blocks without collar	
In-between pontics bwg	
In-between hollow pontics bwhg	97
Betweenblocs bwbl	97
Between hollow blocks bwhbl	98
Hollow pontic blocks hpbl	98
Metal-ceramic blocks mkbl	99
Metal-ceramic blocks with shallow collar	
fg-mkbl	99
Aesthetic and ergonomic metal-ceramic	
blocks äe-mkbl	100
Aesthetic wax veneers ä-wv	100
Hollow pontics hpg	10
Aesthetic-Gnathoflex	
Gnathoflex Premium	
Optiguss	
Crystals and beads	
vesting and casting	
System steps	11
Sprues for vacuum pressure casting	
Sprues for centrifugal casting	111
Rinsing heads	
Rinsing heads for voluminous castings	
Double rinsing heads	
Double rinsing heads for voluminous castings	
Protek wax patterns cut to size	111
Biotec reels of wax pattern	
Reels of wax pattern	
Wax pattern sticks	111
Quadro wax profile	
Quadrosticks	117
Casting pears	117
amework modelling	110
System steps	
Wax-Lite surface tension reducing agent	120
Silicone and wax surface tension reducing	400
agent	
Mould release agent	
Metal muffle rings	
Fleece liners	
Wash primer for fleece	12
Silicone muffle rings	122
Vakuum Anrührsystem ecovac	
Brevest C+B Speed	
Brevest Rapid 1	
Transfuser	
Investment marker	126
Casting ring marker	
Golden booklet	126
Brealloy C+B 270	
Brealloy MK	127
Brealloy Lot	
Brealloy flux	120
Oxyd-Stop-PM	128
	128 129
Oxyd-Stop-PM Oxyd-Stop-NPM Oxyd-Stop Silver-Palladium alloy	128 129 129

Brecid pickling agent.....

# bredent - Casting technology

### Manual



The "bredent casting technique" loose-leaf folder (Dental casting, accurate – homogeneous – compatible) is intended to be used as a manual by the user. The folder's purpose is not to "transform" the dental technician into a material scientist or metallurgist using scientific data and chemical formulas. It is rather the objective of the folder to simplify reliable scientific data and experience for the dental technician as a user.

#### bredent casting technique

The dental cast

precisely fitting - homogeneous - compatible

230 pages

REF 992 961G B




# Course program


In the "bredent casting technique" course you will learn the systemati cal procedures.

The know-how conveyed

in this course is to enable you to achieve reproducible high-quality results.




Helpful knowledge of furnace management and the behaviour of the investment material allows many problems to be solved.



Linear preheating

How must the preheating furnace be set for each mould sizes?



How does the galvanic element work and what happens in the mouth?






For easy and secure removal of the modelling, insulate the model die using Isobre wax insulating liquid. Micro-fine application does not alter the accuracy of fit of the model.

Low distortion waxes are required in order to achieve a model that fits exactly. For this, bredent offers a wide range of colours and qualities for every requirement. Various wax aids and the well-known Gnathoflex silicone moulds for occlusals are available for time-saving modelling.

Light-curing materials also enable rapid and secure modelling and facilitate everyday work.



The special dipping waxes facilitate the manufacture of wax caps which fit perfectly. The visual check of the coating thickness provides the highest level of security and reduces expensive processing.





The distortion-free cervical waxes make an exact border seal possible. Ideal processing is ensured due to the uniform heating of the wax with an electric wax knife.



Modelling waxes in various colours and qualities provide a wide range of products. Additional products, such as wax aids or moulds for occlusals facilitate everyday work and reduce the outlay of time.

#### Materials used

Isobre wax insulating liquid	page 86
Elaflex	page 86
Visio-Dip	page 80
Dipping wax	page 87
Cervical disc	page 87
Waxpool duo	page 88
Wax knife bwm 3	page 89
Quick Change	page 90
Wax for outer copings	page 9
Biotec-Wax for outer copings	page 9
K2 exact	page 9
Standard Modelling Wax	page 9
Gecko	page 92
Splendido	page 92
KBI-Wax	page 92
Life-Color-Wax	page 92
Biotec modellling wax	page 93
Cervical wax	page 93
Biotec-Cervical wax	page 93
SERACOLL UV	page 94
compoForm UV	page 9
Biotec metal-ceramic blocks without collar	page 90
In-between pontics bwg	page 97
In-between hollow pontics bwhg	page 97
Betweenblocs bwbl	page 97
Between hollow blocks bwhbl	page 98
Hollow pontic blocks hpbl	page 98
Metal-ceramic blocks mkbl	page 99
Metal-ceramic blocks with shallow collar	
fg-mkbl	page 99
Aesthetic and ergonomic metal-ceramic	
blocks äe-mkbl	page 100
Aesthetic wax veneers ä-wv	page 100
Hollow pontics hpg	page 10
Aesthetic-Gnathoflex	page 102
Gnathoflex Premium	page 10!
Optiguss	page 108
Crystals and beads	page 109



# Isobre wax insulating liquid



Isobre wax insulating liquid 750 ml REF 540 0104 0



Brush pen pk 20 20 ml REF 540 0072 0



2



Isobre wax insulating liquid on organic basis is absolutely reliable, solvent-free and can be washed off easily. Neutral against plastic, ceramic, metal, plaster and painted surfaces. Even when the insulated surface has dried, Isobre wax insulating liquid will produce a highly efficient, micro-fine insulating layer which ensures simple and safe removal of the wax pattern. Highly absorbing surfaces must be insulated 2 to 3 times.

# reliable, exact separation of the wax pattern against all dental materials.

Micro-fine insulating liquid on organic basis for

## Elaflex




Elaflex purple 130 g REF 510 0090 0



Even in inlays, the cavities can be precisely prepared using Elaflex. This way, modelling is simplified.



Elaflex is so flexible that the wax coping is not deformed when it is removed.





Elaflex is self-insulating on all metal parts that are polished to high luster.

# Visio-Dip



Visio-Dip yellow 130 g REF 510 0073 0



At a wax thickness of more than 0.4 mm the die is no longer visible.



The die becomes visible at a layer thickness of 0.3 mm. Less finishing work is required if a precise wax thickness is ensured.

Visual dipping wax with a wax thickness of < 0.3 mm. Ideal coating thickness checks due to the transparency of the wax.



### **Tauchwachs**



# Accurate, precise fitting wax copings with properties similar to resin.

The dipping waxes allow the fabrication of highly precise wax copings with perfect fit.

Thanks to the properties that are similar to those of resins they can be processed on all surfaces without separating. Different colors ensure perfect contrast to the subsurface.

The basis for efficient and precise working!

**green** 130 g **REF 510 0087 0** 

Odentine color 130 g REF 510 0089 0 yellow 130 g REF 510 0085 0

**brown**130 g **REF 510 0088 0** 

Shows the constant wax thickness and excellent marginal fit without having to re-wax the cervical margin. Dipping wax is supplied as beads.

Use the white wax to block out untercuts or to build up the ideal shape of preparation. It does not adhere to the other waxes.

red, 130 g

Accessories:



Cervical disc REF 320 0091 0

REF 510 0086 0



The dipping wax requires no separating medium for removal from smooth preparations or metal work. Hence it is ideal for precision dental technical work.



Tooth-colored dipping wax perfectly suitable for Life-Color wax structures.

### Cervical disc



Cervical disc REF 320 0091 0



Increases precision and reduces working time when making wax or plastic patterns.



Shows a comparison of the cuts made in 0.5 mm thick wax using a scalpel (left) and cervical disc (right).



This precision steel cutting disc is 0.1 mm thick, 3.0 mm in diameter and can be guided exactly when cutting.



Shows an extremely precise pattern, produced using our dipping wax without the need to rewax the cervical margin.



# Waxpool duo



Waxpool duo unit REF 110 0150 0

### Assortment

- 4 parts
- 1 Waxpool duo unit
- 1 Waxpool duo handpiece
- 2 Waxpool duo contouring blades at your choice

REF 110 0152 0

Accessories:

Waxpool duo handpiece REF 110 0151 0



Rest REF 140 0096 5

Wax dipping unit and wax knife all in one – digital control for added comfort

- Stable and easy to clean plastic housing
- Exchangeable plastic lids
- Clear design
- °C or °F can be selected

#### Wax dipping unit

- Precise temperature control of the dipping wax for increased quality
- High-performance heating elements reduce the time for heating the wax
- Uniform wax copings thanks to constant temperature control
- Special, lowered safety dipping wax to avoid burning of fingers
- Melting temperature up to 120° C

#### Wax knife

- A separate wax knife can be connected
- A single unit at the working place
- Non-tiring working thanks to ergonomic design of the handle
- Special insulating elements reduce heating up of the handle
- Simple exchange of blades
- Boost key for quick heating up to the end temperature
- Maximum temperature of 240° C



Contouring blade size 1 REF 320 WP4G 1



Contouring blade size 3 REF 320 WP4G 3



Contouring blade size 5 REF 320 WP4G 5



Contouring blade standard REF 320 WP47 2

### Wax knife bwm 3



Control unit bwm 3 with handpiece and contouring blade size 5

REF 140 0096 3

Control unit bwm 3 REF 140 0096 0

Handpiece bwm 3 REF 140 0096 2



Rest bwm 3 REF 210 0045 1




Footswitch bwm 3 REF 140 0096 1

#### Accessories:

Foam rubber grip lining 4 pieces REF 140 0096 4

Electric wax knife featuring integrated advanced technology and high quality. The ergonomic handpiece allows to take up wax quickly and ensures comfortable working.

- Ergonomically designed handpiece
- Quick heating up with the Rapid-Speed footswitch
- Adjustable temperature control
- Simple and fast exchange of the contouring blades



size 3
Contouring blade bwm 3

size 5

Contouring blade bwm 3
Standard

REF 320 004G 1

REF 320 004G 3

REF 320 004G 5

REF 320 0047 2



Comfortable and quick removal of the contouring blades.



Device for firm, reliable hold of the handpiece at the unit.



Mobile rest for safe depositing of the handpiece.



Blade shapes proven over numerous years allow individual application.



Integration into the grip for quick and simple exchange of the contouring instruments without the risk of injuries.



The special instrument grip avoids twisting of the contouring tip whilst working.



The contouring tips are stored on the control unit in a safe and clearly arranged manner.



If the wax knife is not needed, it can be placed on the rest in the direct reach of the technician.



The footswitch allows to quickly reach a higher temperature than the one that has been set. Activation of the footswitch is indicated by the control lamp.



Logical and clearly arranged control unit for stress-free and safe working.



Handpiece with flexible, stable cable for simple working.



High-tech dental equipment featuring highly useful function and design – for comfortable and simple working.

# Quick Change



The combination of design, function and systematics.

#### Quick change system for instruments for ceramic, model fabrication and prosthetics

- Design carbon handle esthetic and
- Stainless blade holder with magnetic receptacle for perfect fixation of all blades
- All metal components and magnets are corrosion-resistant
- Temperature resistance of instruments inside the handle: 80°C
- Single hand use with quick change system
- Individual indications for ceramic, model fabrication and prosthetic
- Easy control thanks to reduced range of instru-

- Clever system allows to find the suitable instruments quickly (protection of registered design!)
- Troublesome screwing with keys is avoided
- Blades can be adjusted at any position - familiar working position is retained
- High safety thanks to immediate release of the heated blades
- Proper storage of sensitive ceramic blades
- Ceramic brushes can be perfectly stored in a hanging position
- Ceramic blades with high surface quality for outstanding gliding properties



Carbon handle L 101 mm, Ø 8 mm REF 310 0103 1



w 102 x d 100 x h 75 mm Weight approx.. 575 g REF 310 0103 0

Overview of instruments Dimensions in mm



REF 310 0105 6



Fissure tool REF 310 0103 4



Olive REF 310 0105 7



Croco, smooth REF 310 0103 2



Croco, serrated REF 310 0103 3





REF 310 0105 3



MagicContrast size 8 REF 310 0105 4



MagicContrast REF 310 0105 5



size 6 REF 310 0104 4



KoliBrush size 8 REF 310 0104 5



REF 310 0104 6



REF 310 0104 0





REF 310 0104 2



Blade 0308 Fig. 3 REF 310 0103 7



Blade 0408 Fig. 4 REF 310 0103 9



# Wax for outer copings



For secondary metal elements. Special consistency avoids the formation of grooves on the inner side of the coping.

Wax for outer copings 25 g yellow REF 510 0042 0



Easily spread and, whilst cooling, will not shape creases on the surface exposed to the metal. Extremely high precision of fit, thanks to the minimal shrinkage.

# Biotec-Wax for outer copings



The modelling wax for uniform application of coats with minimal shrinkage and unsurpassed burning out properties.

Biotec-Wax for outer copings 28 g violet REF 510 0061 3



Uniform coping thickness due to perfect scraping properties. When the die shines through, a layer thickness of 0.3 to 0.4 mm has been achieved.



Low shrinkage and excellent burning out properties ensure utmost precision of fit and homogeneous castings.

### K2 exact



K2 exact

60 g

green

grey REF 510 0090 2 yellow REF 510 0090 3 beige REF 510 0090 4

REF 510 0090 5

Extraordinary carving qualities for the highest precision in all crown, bridge and inlay work. Low shrinkage and, as a consequence, precise waxing up are distinctive features of this carving wax. The solidification phase is very short; hence K2 exact carving wax is suitable for selective waxing up. The hardness results in good scraping properties and allows perfect smoothening of the wax model.



Marble plinth and dome REF 320 0042 0



The extremely low shrinkage allows high precision even with feather edges.



A short solidification phase makes it easier and quicker for modelling into the correct tooth shape required.



The strength and homogeneity of the wax provide optimal carving qualities and allow to produce smooth polished surfaces

# Standard Modelling wax



Standard Modelling wax beige 70 g Click-Clack jar REF 510 0078 5 Beige modelling wax for crown and bridgework and for the inlay technique. The solidification temperature of 50°C allows to work quickly. The beige color avoids tiring of the eyes during waxing up and supports the three-dimensionality so that deep occlusal surfaces can be perfectly recognized.



### Gecko





The wax pattern can be compared objectively with the adjacent tooth.





The wax reproduces the contours and colour of the model exactly, which is pleasant to work with and will not tire the user.






Opaque wax allows for improved determination of the depth of the fissures in the wax pattern.



The contours are easier to discern thanks to the way in which these pastell shades reflect the light.

Available in different colors for convenient and non-tiring waxing up.

The special wax quality allows highly precise application and perfect sculpturing.




REF 510 0060 2 REF 510 0060 4 REF 510 0060 1 REF 510 0060 3

# Splendido



This wax is suitable for any type of wax-ups: crowns, bridges and inlays. Also suitable for milling. Splendido is also available as summer wax "Splendido hard", which can be milled up to a room temperature of 40 °C.







The light green colour of this wax provides for improved light reflection and facilitates determination of the final contouring. The opacity of this wax allows for improved determination of the depth of the fisures in the wax pattern.

### KBI-wax



Wax for crowns, bridgework, inlays. Minimal shrinkage, high stability, good modelling properties and smooth surfaces after scraping offer ideal possibilities for any waxing-up technique. Suitable for milling

KBI-Wachs 25 g, blue, medium REF 510 0091 0 hard REF 510 0092 0

techniques.





The light blue colour enables the technician to view the contours and surface structure of the pattern in greater detail. "KBI hard" is available for use in summer. Both waxes are, of course suitable for milling.

### Life-Color-Wax



Tooth-colored wax in two consistencies. Particularly low-shrinkage wax especially developed for the waxing-up technique according to M. A. Polz. Life-Color-Wax

25 g

dentine color, medium dentine color, hard REF 510 0081 0

100 g

dentine color, medium dentine color, hard REF 510 0079 0 REF 510 0078 0



Precise application and superior scraping properties are the distinctive features of this wax.



# Biotec modellling wax



Modelling wax that meets highest demands on modelling properties, shrinkage and complete burning out.

burning out.
The excellent scraping properties ensure perfect fit of the wax crown on the die.
Wax residues can be blown away easily.
Low shrinkage leads to high precision of fit.

**Biotec modelling wax,** 60 g green grey

REF 510 0061 1 REF 510 0061 0 Residue-free burning out is the prerequisite for homogeneous casting.
Perfectly suitable for modelling pressed ceramic crowns and inlays.



Easily controllable stability for specific application across small and large areas.

### Cervical wax



In order to achieve a perfect marginal seal on crowns, inlays, onlays etc., the wax must adapt well and be completely shrink-free.

Cervical wax 25 g red REF 510 0060 5



Cervical wax is used for shaping the cervical margin and adheres perfectly to the coping and sculpturing wax.

### Biotec-Cervical wax



Modelling wax for precision-fit crown margins due to low shrinkage and outstanding burning out properties.

Biotec-Cervical wax 28 g red REF 510 0061 2



The special consistency of this wax, minimum shrinkage and extremely low quantity of residues of combustion ensure perfect fit of the crown margin.

# Biotec milling wax



Excellent milling wax with superb modelling properties. Outstanding scraping and milling properties since sticking of wax to the bur is avoided.



Biotec milling wax 28 g blue REF 510 0061 4



Enormous amount of time is saved thanks to good modelling properties since no other wax is required for the shear distributor.



Extremely accurate milling wax to produce smooth and shining surfaces during milling.



The bredent bur system can be found in chapter C!



Residue-free burning allows for use in pressed ceramics.



### SERACOLL UV




### SERACOLL UV connects bridges



The wax-up is prepared separately and thus any stress within the wax structure is avoided. After checking and fine contouring the crown margins, the wax-up is placed onto the model again.



One drop of SERACOLL UV is added into the separating gap using the probe.
Thanks to the good capillary effect of SERACOLL UV the gap is evenly filled. The optimum size of the separating gap is < 0,3 mm. If required, add another drop to restore the original shape.





After each application of material, SERACOLL UV is hardened in standard UV light-curing units for at least 90 seconds.



The wax-up with sprues being attached can be removed from the model without the formation of stress and invested subsequently.



SERACOLL UV light-curing wax adhesive 2 x 3 ml 2 dosing dishes REF 540 0115 1



SERACOLL UV is perfectly suited for stress-free bar structures and other applications in the field of . implant restorations.



SERACOLL UV is suited for all standard light sources (UV and LED). When using hand lamps, each surface of the model must be cured for



By applying a thin layer of SERACOLL UV, all rough model surfaces can be smoothened. More homogeneous casting results are obtained.



Thanks to the enclosed small dosing dish, the required quantity of SERACOLL UV can be taken up with the probe.



Pi-Ku-Plast, compoForm UV and waxes can be connected with each other without any problems.



Attachment elements can be positioned in the parallelometer; one drop of SERACOLL UV is added into the gap and the rounded, clean transition zones are cured using a hand lamp.

# compoForm UV





The stability of compoForm UV renders the material perfectly suitable for the transfer of the jaw situation and, consequently, stress-free working is ensured.

Light-curing composite for modelling, fixation of separated bridges and for quick fabrication of post and core restorations. compoForm UV burns out without leaving any residue and produces homogeneous casting results.

Individual modelling directly from the syringe. Thanks to immediate hardening with a polymerization lamp the model can be built up in a safe and controlled manner.



compoForm UV can be used in conjunction with modelling wax and is perfectly suitable for inter-locking the model prior to investing. This way investing without any deformation is possible.



Modelling and further processing of telescopic and conical crowns can be perfectly controlled by means of a visual check of the layer. The high stability of the hardened composite allows reworking with a bur.

compoForm UV 2 x 3 ml syringes 10 application cannulas REF 540 0115 0



Thanks to low shrinkage and burning without any residue, the composite is ideal for fixation of bridges to be soldered.



Thanks to low shrinkage and burning without any residue, the composite is ideal for fixation of bridges to be soldered.

Accessories:



Application cannulas 25 pieces REF 580 0001 8



Undercuts on dies can be quickly and completely blocked out.



Burning without any residue and reduced swelling behavior provide perfect preconditions for topquality casting results.

### Biotec metal-ceramic blocks without collar b-mkbl



Auxiliary wax elements with properties similar to modelling wax and very limited quantity of residues of burning out. The melting point, hardness and scraping properties are adapted to the modelling waxes to allow simple and specific connecting of the auxiliary wax elements with the crown pattern. The extremely low quantity of residues of burning out of the Biotec auxiliary wax pattern provide perfect preconditions for smooth, homogeneous cast surfaces.

Tooth	17-14 RP	12-22 RP	24-27 RP/
Size C	0000	0000	<b>PPD0</b>
В	0000	0000	0000
А	0000	<b>6000</b>	0000
А	0000		0000
В	0000	HHH	0000
С	0000	HHH	0000
Tooth	47-44	42-32	34-37






Refill packages (RP): Each form and size is available as refill package cont. 10, 25, 50 or 100 pieces each. Please enter the exact number of pieces into the box next to the desired form.



Minikit: 18 forms x 2 parts REF B13 000 MK

Sender (Stamp):	Customer No.
	Date, Signature

### In-between pontics bwg



Design by Jan Langner

Refill packages (RP): containing 50 pieces each

Assortment In-between pontics bwg: containing 540 pcs. (54 different forms with 10 individual pieces each)

REF D00 5401 0


Please enter number of desired packages in the box.

### In-between hollow pontics bwhg

Tooth	17 / 16 RP/	13 RP	12	11 RP	21 RP/	22 RP	23 RP/	26 / 27	UJ premolar RP/
Size						<u> </u>			
С	40				Section 1			-00°	
В	-								
А				0					
А	0	•	-			*	<b>*</b>	4	
В	*	-	<b>*</b>		9	*	<b>*</b>	*	
С	*	<b>*</b>	•	•	9	<b>→</b>	*	*	<b>*</b>
Tooth	47 / 46	45	44	43	33	34	35	36 / 37	LJ-Anteriors

### Betweenblocs bwbl

Refill packages (RP): containing 50 pieces each Anterior assortment In-between hollow pontics bwhg containing 300 pieces (27 different forms)



REF D01 2701 0

Please enter number of desired packages in the box.

Anterior and posterior assortment In-between pontics and In-between hollow pontics bwhg containing 540 pieces. (different forms with 54 individual pieces

each 27 Hollow pontics and 27 Massive pontic) REF D01 5401 0

Please enter number of desired packages in the box.

Refill packages (RP): each 25 pieces Assortment

Between blocks bwbl: containing 180 pieces (18 different forms, each 10 blocks)

REF D00 1801 0

Please enter number of desired packages in the box.


Sender (Stamp):

Illustrations are full size

Customer No.

Further order:

### Between hollow blocks bwhbl





Refill packages (RP): containing 25 blocks each

#### Assortment

Between hollow blocks bwhbl: containing 12 hollow blocks (posterior blocks) and 6 massive blocks (anterior blocks)

#### REF D01 1801 0

Please enter number of desired packages in the box.

### Hollow pontic blocks hpbl

Tooth	17-14 RP	12-22 RP	24-27 RP
Size C	<b>6494</b>	<b>\$866</b>	6666
В	0000	***	6600
А	<b>66.00</b>	***	6660
А	<b>6600</b>		6640
В	G-G-49	7776	6600
С	60M	7799	<b>66.00</b>
Tooth	47-44	42-32	34-37



Refill packages (RP): containing 25 blocks each

#### Assortment

Hollow pontic blocks hpbl: containing 180 blocks (18 different forms with 10 blocks each)

#### REF A11 1801 0

Please enter number of desired packages in the box.

Sender (Stamp):	Customer No.	Furth
	Date, Signature	

her order:

Illustrations are full size

Tooth	17-14 RP	12-22 RP	24-27 RP
Size C			0000
В			0000
А	0000		0000
А	*****	***	0000
В	****	0000	9000
С	****	0-0-0	***************************************
Tooth	47-44	42-32	34-37

Refill packages (RP): 25 each

Assortment

Metal-ceramic blocks mkbl:

containing 180 pieces. (different forms with 18 individual pieces each 10 blocks)

REF A00 1801 0

Please enter number of desired packages in the box.

Metal-ceramic blocks with shallow collar fg-mkbl

Tooth I	17-14	RP	12-22	RP	24-27	RP
Size C	0000		0000		0000	
В	0000		0000		0000	
А	0000		0000		0000	
А	-		4444		1000	
В	9999		9999		1000	
С	9999		0000		1000	
Tooth I	47-44		42-32		34-37	



Refill packages (RP): 25 each

Metal-ceramic blocks with shallow collar fg-mkbl: containing 180 pieces. (different forms with 18 individual pieces each 10 blocks)

REF A01 1801 0

Please enter number of desired packages in the box.

Sender	(Stamp)

Illustrations are full size

Customer No.


Further order:

Date, Signature



### Aesthetic and ergonomic metal-ceramic blocks äe-mkbl

Tooth	17-14 RP	12-22 RP	24-27 RP
Size C	6-6-6-6	6669	<b>6666</b>
В	6666	6666	4444
А	6666		6666
А		HH	
В		HH	
С			
Tooth	47-44	42-32	34-37



Refill packages (RP): containing 25 blocks each


Aesthetic and ceramic metal-ceramic blocks äe-mkbl: containing 180 blocks (18 different forms with 10 blocks each)

REF A02 1801 0

Please enter number of desired packages in the box.

### Aesthetic wax veneers ä-wv

Design by Jan Langner, Master Dental Technician





Aesthetics and ergonomic ceramic pontics as a basis



Aesthetics wax veneers from the palatal side on aesthetic and ergonomic ceramic pontics



Aesthetics wax veneers from the labial side



Assortment Aesthetic wax veneers ä-wv: containing 240 pieces (24 different forms with 10 individual pieces each)

Refill packages (RP): containing 50 pieces each



### REF C13 2401 0

Please enter number of desired packages in the box.




Sender (Stamp):

Customer No.

Further order:

Date, Signature

Illustrations are full size



Sender (Stamp):	Customer No.	Further order:
	Date, Signature	

### Aesthetic-Gnathoflex

16	1 <u>4</u> 5	<b>2</b> <sup>4</sup> / <sub>5</sub>	26	Tooth
16D	1 <del>4</del> D	2 <sup>4</sup> <sub>5</sub> D	26D	Size <b>D</b>
16C	1 <u>4</u> C	24€€	26C	Size C
16B	1 <u>4</u> B	2 <sup>4</sup> / <sub>5</sub> B	26B	Size <b>B</b>

Assortment

36 pieces (12 different moulds in sizes B, C, D)

#### REF 429 Ä003 6

Please enter number of desired packages in the box.

Please enter number of desired parts in the respective box.

47B	46B	45B	44B	348	35B	36B	37B	Size <b>B</b>
47C	46C	45C	44C	34C	35C	36C	37C	Size <b>C</b>
47D	46D	45D	44D	34D	35D	36D	37D	Size <b>D</b>
47	46	45	44	34	35	36	37	Tooth

Sender (Stamp): Further order: Customer No. Date, Signature

### Aesthetic-Gnathoflex



Flexible, re-usable silicone moulds for creating wax, acrylic or ceramic occlusals. Shaping occlusal surfaces within seconds thanks to highly flexible silicone moulds.

- Can be used for many applications with wax, acrylic and porcelain
- Saves time thanks to the teflon coating super smooth surfaces are created immediately
- Only needs to be bought once Aesthetic-Gnathoflex moulds can be reused Achieve increased turnover in less time yet maintain constant high quality!



Gnathoflex is fabricated from high-grade silicone which maintains its stability up to 250 °C. The occlusal path is created by the mould, which is only 0.5 mm thick.



Produce acrylic occlusals

Gnathoflex is extremely flexible and very durable and ensures that its shape is maintained.



Produce porcelain occlusals

Gnathoflex precisely reproduces the anatomical cusps and fissures in wax, acrylic or porcelain. The result is a smooth, glazed surface.

Produce <u>wax</u> occlusals in 40 secs.



The aesthetic, but functional anatomy which Gnathoflex produces in wax patterns, may be modified to suit individual requirements.



in 90 secs.

Gnathoflex provides for high-luster occlusals when using the acrylic or composite of your choice. It is also ideal for temporary bridgework.



in 180 secs.

Porcelain work may be built up using the shade and anatomy of your choice. Any porcelain may be used.



The copings are prepared as usual, using wax or acrylic.



Gnathoflex is filled with the modelling wax of your choice.

You may also use Gnathoflex as the basic mould for fabrcating your laboratorys own individual occlusals.



As soon as the wax begins to gel,



place the Gnathoflex on the die.

Open the articulator 0.5 mm, measured at the surface being waxed-up, in order to allow for the thickness of the Gnathoflex.

### Choose for yourself: Precisely articulated multiple contacts, minimum contact or exactly 0.5 mm out of occlusion.



Close the articulator and attach the Gnathoflex to the wax coping, using a drop of wax.



The second bridge abutment is waxed-up using the same method.



Shows the high-luster precision wax reproduction of the silicone mould.



A pre-formed wax pattern is used when waxing-up the occlusal surface of the pontic.



Shows the completed occlusal aspect of the bridge. Gnathoflex ensures uniform, aesthetic occlusals.



Lingual view of the intercuspidation between the wax pattern and the opposing teeth.

103



Fax (+49) 0 73 09 / 8 72-4 44

The buccal view shows the central cusp-to-fossa relationship.



Terminal occlusion contact areas, made visible by means of articulating film.

### Aesthetic-Gnathoflex

#### Fabricate occlusals using any acrylic, no trimming required.



Prepare the coping as usual and apply the crown and bridge acrylic (dentine) of your choice.



No separating agent is necessary when filling the Gnatoflex with acrylic. First fill the cusps with incisal and then fill the mould completely with dentine. Place the mould on the bridge.



Close the articulator. When using photo-curing acrylic, commence polymerisation now, in order to fix the bite.



Remove the bridge, apply the interdental contact areas and complete the polymerisation. Having carried this out, remove the silicone moulds.

#### Accurately shaded porcelain occlusals, easier than ever before.



Apply the opaque, fire it and build-up the bridge using dentine. Brush a thin coat of Gnathoflex separating agent into the silicone mould.



First fill the Gnathoflex with incisal and then with dentine. The incisal material should be brushed out from the cusps toward the margins. In order to allow for the shrinkage, the articulator should be opened by more than 0.5 mm when building up porcelain occlusals.



Place the Gnathoflex on the bridge and fix it using porcelain. Complete the build-up using incisal, dentine or a mixture of both – depending upon the shade.



Dry the porcelain as usual or use a hair dryer. Carefully remove the Gnathoflex from the bridge. The bridge can now be removed from the model, further porcelain applied to the contact areas and the bridge fired as normal.

#### Accessories:



Isoflex 20 ml REF 540 0101 3

### Aesthetic-Gnathoflex



# Create occlusal surfaces in seconds with these extremely flexible Teflon silicone moulds

- Can be used for many applications with wax, acrylic or porcelain
- Saves time thanks to the Teflon coating super smooth surfaces are created immediately.
- Only has to be bought once.
   Aesthetic Gnathoflex moulds can be re-used.

   Achieve increased turnover in less time yet maintain constantly high quality!
- For wax patterns
- For ceramic
- For acrylic

# **Gnathoflex Premium**



Extremely flexible, reuseable silicone moulds for creating occlusals.

40 to 180 seconds for a gnathological, aesthetically functional occlusal with wax, acrylic and ceramic.





Gnathoflex
Study model
FF1 Set
1 UJ model
1 □ model
REF 992 5027 3

Original size



Gnathoflex Study model FF1 mini Set 1 UJ model mini 1 LJ model mini REF 992 5027 4

Original size






Isoflex - if 20 ml REF 540 0101 3

### **Gnathoflex Premium**

#### Wax occlusals



The copings are prepared as usual, using wax or acrylic.



To compensate the thickness of the Gnathoflex, the bite is raised by 0.5



Fill Gnathoflex with wax and wait until the wax begins to gel.



Once the wax has hardened, place the Gnathoflex occlusal onto the coping.








Close the articulator in the position of maximum intercuspidation and attach the occlusal to the coping using a drop of wax. Depending on the situation, two or more Gnathoflex occlusals can be placed simultaneously or one after the other. The contact can be strongly varied by raising or lowering the antagonist.







High-luster, gnathologically shaped wax occlusals with perfect contact to the antagonist are the perfect basis for smooth and precision-fit casting. Gnathoflex Premium helps to save time during the preparation of the wax model as well as during finishing of crowns and bridges.





Prepare the structure as usual and apply crown and bridge acrylic (dentine).



No separating agent is required when filling the Gnathoflex with acrylic. First fill the cusps with incisal and then fill the mould completely with



Place Gnathoflex on the bridge, close the articulator and polymerize with UV light in order to fix



Then the bridge is removed, interdental contact areas are applied and polymerization is completed.

### Ceramic occlusals




The opaque is fired on the metal structure.



Hold the Gnathoflex with the tweezers and apply Isoflex insulating liquid onto the inner surface.



Tap several times on the Gnathoflex to remove excess Isoflex insulating liquid.







Fill incisal into the cusps and brush out from the cusps toward the margins. Fill the  $\label{thm:continuous} \textbf{Gnathoflex with dentine and place on the bridge \ \ structure. \ Close \ the \ articulator \ and$ turn it. Fix the Gnathoflex occlusals to the bridge using dentine material. Dry the object and carefully remove the Gnathoflex. The other occlusals are prepared accordingly. Complete the bridge.



# **Gnathoflex Premium**

17	16	15	14	24	25	26	27	
17C	16C	15C	14C	24C	25C	26C	27C	С
429 P017 C	429 P016 C	429 P015 C	429 P014 C	429 P024 C	429 P025 C	429 P026 C	429 P027 C	
17B 429 P017 B	429 P016 B	15B 429 P015 B	14B 429 P014 B	24B 429 P024 B	25B 429 P025 B	26B 429 P026 B	27B 429 P027 B	В
17A	16A	15A	14A	24A	25A	26A	27A	Α
429 P017 A	429 P016 A	429 P015 A	429 P014 A	429 P024 A	429 P025 A	429 P026 A	429 P027 A	
47A	46A	45A	44A	34A	35A	36A	37A	A
429 P047 A	429 P046 A	429 P045 A	429 P044 A	429 P034 A	429 P035 A	429 P036 A	429 P037 A	
47B	46B	458	44B	34B	35B	36B	37B	В
429 P047 B	429 P046 B	429 P045 B	429 P044 B	429 P034 B	429 P035 B	429 P036 B	429 P037 B	
47C	46C	45C	44c	34C	35C	36C	37C	С
429 P047 C	429 P046 C	429 P045 C	429 P044 C	429 P034 C	429 P035 C	429 P036 C	429 P037 C	
47	46	45	44	34	35	36	37	

Please select the desired	d parts	from	the	original	il-
lustrations.					

Please	enter	number	of	desired	parts	in	the	respective
box or	highli	ight the	ass	ortment	t.			

Sender (Stamp):	Customer No.
	Date, Signature


Assortments **Gnathoflex Premium:** 16 pcs, containing 16 moulds in size A REF 429 P000 A

16 pcs, containing 16 moulds in size B REF 429 P000 B

16 pcs, containing 16 moulds in size C REF 429 P000 C

48 pcs, containing 16 moulds in 3 sizes ABC each REF 429 P004 8

Please enter number of desired packages in boxes.



# Framework modelling

# **Optiguss**

The solution for increased perfection with less effort.

Optiguss Micro – 5 micron coating – or Optiguss Macro – 10 micron coating – can be applied easily and quickly to the wax pattern to smooth, seal and reinforce it without changing its shape.

The use of Optiguss reduces the finishing time by 50 % compared to a conventional cast surface.



Optiguss-macro 15 ml REF 520 0092 0

Optiguss-micro 15 ml REF 520 0093 0



Optigussschale macromacro 2 pieces REF 390 0035 0



3 Brushes size A + brush holder REF 330 0114 6

3 Brushes size B + brush holder REF 330 0114 7

3 Brushes size C + brush holder REF 330 0114 8



3 Brushes size A

3 Brushes size B 3 Brushes size C

2 Brush cleaning pot

1 Brush cleaner REF 520 0091 0







Brush cleaning pot 2 pieces REF 390 0037 0



Brush cleaner 20 ml REF 520 0094 0



Applying Optiguss creates super smooth surfaces.



Deep fissures, which cannot be reached with a fissure bur, can be smoothed with Optiguss. This simplifies polishing of gnathologically designed occlusals.



The finishing time can be reduced by more than 50 % due to more homogeneous surfaces.



Approximal contact areas are strengthened, yet retain their shape.



Fitting surfaces are built-up properly and smoothed, which reduces the time required for trimming.



# Crystals and beads



Retention crystals 0.2 mm, 20 ml REF 530 0048 0 0.5 mm, 20 ml REF 530 0050 0 0.8 mm, 20 ml REF 530 0051 0



### Retention beads



0.2 mm Ø, 20 ml REF 530 0201 0 0.4 mm Ø, 20 ml REF 530 0220 0 0.6 mm Ø, 20 ml REF 530 0210 0 0.8 mm Ø, 20 ml REF 530 0200 0



Retention adhesive

20 ml

REF 540 0071 1

Thinner

REF 540 0071 2

Retention adhesive, white - no capillary action, solubilizes the crystals and adheres perfectly as well as having a long working time.

Optimum retention leads to the strongest possible acrylic/metal junctures.





Crystals have double the retentive area of beads.

Original size



10x magnification (gold)

















Micro retention beads result in elegant facings.


























# **INNOVATIONS**



For almost 40 years, bredent has offered innovative solutions for use in the dental technology laboratory and this forms a central part of the company's philosophy.

New developments and procedural techniques will significantly influence the future of dental technology and dentistry.

Due to close contact with customers and monitoring of national and international dental markets, ideas and recommendations for modifications find their way into the company's product portfolio. This exchange enables the processes in the laboratory and in the practice to be optimised and costs to be reduced.

bredent desires and strives to achieve the consistent implementation of the ISO 9001 quality standards, thereby offering the user maximum product safety and, as a result, ensuring that patients receive an implant prosthetic restoration of the highest quality.

Efficiency through innovation!







Thanks to their design, the special casting channels improve the casting result. Centrifugal casting and vacuum pressure casting must be differentiated between in order to create the ideal prerequisite for each casting system.


bredent's casting system is completed by the casting

channel system. With the right choice of casting channel and rinsing heads, the casting is more precise and homogeneous. Correct metallurgic work is therefore produced, which can also be better processed and defects in the additional processing can be reduced.



The rinsing heads create melt reservoirs and act as a pressure balance for homogeneous casting. Alloy is saved by reducing the casting channel system and an outstanding casting result is achieved.






The semi-circular casting mould system enables exact placement of the model outside of the heat centre. The casting technology products are matched to one another to ensure that perfect and reproducible casting results can be produced.

Sprues for vacuum pressure casting	page 112
Sprues for centrifugal casting	page 113
Rinsing heads	page 114
Rinsing heads for voluminous castings	page 114
Double rinsing heads	page 115
Double rinsing heads for voluminous	
castinags	page 115
Protek wax patterns cut to size	page 115
Biotec reels of wax pattern	page 116
Reels of wax pattern	page 116
Wax pattern sticks	page 116
Quadro wax profile	page 117
Quadrosticks	page 117
Casting pears	page 117

# Sprues for vacuum pressure casting

Sprues and rinsing heads suitable for all casting techniques to ensure homogeneous, uniform and predictable casting results.








### Assortment

450 pieces

Vacuum pressure casting, containing 30 sprues and 30 rinsing heads each,

REF 430 0146 0

# Sprues for vacuum pressure casting for voluminous castings



Ø mm	Ø mm	Ø mm	pieces	KEF
3.5	6.5	5.0		430 0143 7 430 0146 9
3.5	7.5	5.0		430 0143 8 430 0147 1
3.5	8.5	5.0		430 0143 9 430 0147 2
3.5	9.5	5.0		430 0144 0 430 0147 3




### Assortment

211 pieces

Vacuum pressure casting for voluminous castings, containing 30 sprues and 30 rinsing heads each, 25 g Protek wax patterns (rods) cut to size, Ø 1.0 mm, REF 430 0147 0

# Sprues for centrifugal casting



Sprue channe	l Head	Sprue	Pack.	REF
Ø mm	Ø mm	Ø mm	pieces	
2.5	4.5	3.0	50	430 0144 1
			150	430 0147 7
2.5	5.0	3.5	50	430 0144 2
			150	430 0147 8
3.0	6.0	3.5	50	430 0144 3
			150	430 0147 9
3.5	6.5	4.0	50	430 0144 4
			150	430 0148 1






### Assortment

390 pieces

Centrifugal casting, containing 30 sprues and 30 rinsing heads each, **REF 430 0148 0** 

# Sprues for centrifugal casting for voluminous castings



Sprue cha Ø mm	nnel Head Ø mm	<b>Sprue</b> Ø mm	Pack. pieces	REF
3.5	7.5	4.0	50	430 0144 5
			100	430 0148 2
3.5	8.0	4.0	50	430 0144 6
			100	430 0148 3
3.5	9.5	4.0	50	430 0144 7
			100	430 0148 4



### Assortment

181 pieces

Centrifugal casting for voluminous castings, containing 30 sprues and 30 rinsing heads each, 25 g Protek wax patterns (rods) cut to size, Ø 1.0 mm, REF 430 0148 5

# Rinsing heads


Rinsing heads for vacuum and centrifugal casting. Since the residual air is forced into the rinsing heads, a high density of the structure is obtained to deliver superior casting results.

Sprue channel Ø mm	<b>Head</b> Ø mm	Pack. pieces	REF
2.5	4.0	50 150	430 0144 8 430 0148 6
2.5	5.0	50 150	430 0144 9 430 0148 7
2.5	5.5	50 150	430 0145 0 430 0148 8
3.0	6.0	50 150	430 0145 1 430 0148 9
3.5	6.5	50 150	430 0145 2 430 0149 1



The retainer helps to attach the rinsing heads.


# Rinsing heads for voluminous castings





The retainer is cut off with the wax knife after waxing up the rinsing heads.


# Double rinsing heads






The double rinsing heads are used for two objects with the same volume. Attaching is easier and the amount of work is reduced.

# Double rinsing heads for voluminous castings



# Protek wax patterns cut to size

Protek wax patterns cut to size, rods, for cooling fins, pressure compensation and vent channels





During casting the compressed air is displaced into the air channels and produces castings free from shrinkage cavities, which can be easily fitted.


# Biotec reels of wax pattern



Sprue wax with organic components, highly flexible and burns out perfectly.

# **Biotec reels of wax pattern** 250 g, beige

Cross sec	tion in Ø mm	REF
•	1.2	430 0801 2
•	1.5	430 0801 5
•	2.0	430 0802 0
	2.5	430 0802 5
	3.0	430 0803 0
	3.5	430 0803 5
	4.0	430 0804 0
	5.0	430 0805 0



shaping allow specific, stress-free attaching of the sprues. Residuefree burning out is the basis for perfect casting results. Perfectly suitable for modelling pressed ceramic crowns and inlays.

High flexibility and low

elastic recovery after



# Reels of wax pattern



Various diameters of wax pattern are available in medium and hard consistencies.

### Reels of wax pattern, 250 g

REF

Cross section

in Ø mm		blue (medium hard)	green (hard)
•	1.2	430 0115 0	
•	1.5	430 0115 5	
•	2.0	430 0116 0	430 0111 0
	2.5	430 0116 5	430 0111 5
	3.0	430 0117 0	430 0112 0
	3.5	430 0117 5	430 0112 5
	4.0	430 0118 0	430 0113 0

430 0118 5

REF

430 0113 5



The wax patterns can be bent without recovering elastically or becoming pinched.



# Wax pattern sticks





Wax pattern for sprues made of extremely hard special wax.

Extremely hard special wax to ensure that the wax model will not be deformed when it is removed; result: highly accurate castings, even for large-span work. Dimensional stability at room temperature so that the wax model can be safely removed.

### Wax pattern sticks 250 g, red

5.0

Ø mm x le	REF	
•	2.0 x 115	430 0172 3
•	2.5 x 115	430 0172 1
•	3.0 x 115	430 0168 0
	3.5 x 115	430 0169 0
	4.0 x 115	430 0170 0
	4.5 x 115	430 0172 2
	5.0 x 115	430 0171 0
	6.5 x 115	430 0172 4



Wax model with direct fitting of the sprues. No deformation of the model during removal if wax profile sticks are used.



Wax pattern with attached sprues connected with a bar. The wax pattern sticks can be easily bent by heating them slightly and thus adapted to the bridge shape. Safe removal of the model at room temperature.



# Quadro wax profile



### Square sprues for better casting results.

Studies have shown that all liquids – including liquid metal – flow in drops; that also applies to flowing into a square sprue.

Accordingly, the gas (air) contained in the cavity (casting mould) can escape freely across the unfilled corners. Results:

- no swirling of molten metal due to the back pressure of the residual air
- faster flowing in of the molten metal
- more homogeneous castings
- smoother surfaces
- increased precision of fit



### Quadro wax profile

250 g, green

■ 1.75 x 1.75 mm

2.25 x 2.25 mm

REF 430 0691 0 REF 430 0692 0


3.00 x 3.00 mm

REF 430 0693 0

### Quadrosticks



The Quadrosticks made of extra-hard special wax can not be deformed at room temperature. This way distortion of the wax pattern is avoided when removing it from the model. This is a crucial prerequisite for precision-fit dental work.



# Quadrosticks qust

### Assortment

150 pieces Quadrosticks 65 pieces 1.75 mm 50 pieces 2.25 mm 35 pieces 3.00 mm

REF 430 0164 0

### Quadrosticks, green

- 1.75 X 1.75 X 50 mm, 150 PCSREF 430 0165 0
- 2.25 X 2.25 X 50 mm, 125 PCS REF 430 0166 0
- 3.00 X 3.00 X 50 mm, 85 PCS REF 430 0167 0

# Casting pears



Pointed "Lost head" for fast and specific attaching the casting object with wax; suitable for centrifugal casting.

### Casting pears

100 pieces each

Ø	length	
6 mm	9 mm	REF 430 0740 6
7 mm	10 mm	REF 430 0740 7
8 mm	11 mm	REF 430 0740 8
9 mm	12 mm	REF 430 0740 9



# BETTER QUALITY OF LIFE – BACK TO A HEALTHY SLEEP




Snoring has become the most widespread disease. In 50% of bedrooms, healthy sleep is disturbed by extreme snoring – this can be as loud as 90 decibels, which corresponds to the noise of a passing lorry.

Two forms of treatment are available. The protrusion splint, which is worn in the case of slight or medium obstructive sleep apnoea syndrome, and biofunctional treatment with a vacuum activator, are used in the

case of primary snoring without obstruction.

Become an advisor for your practice and establish your laboratory's position in a new and attractive area of business. Comprehensive information material is available for laboratories, practices and patients.


Added value for your laboratory!



bredent alloys allow for very easy processing due to the low degree of hardness. This reduces expensive processing and reduces the need for a large amount of burs. A high level of quality is very important to bredent. This is also the case when it comes to investment materials, which have been developed for timesaving processing due to perfectly fitting results. The investment materials, which can be inserted exactly, enable stress-free working and significantly reduce adjustment times. Accuracy is aided by the special shape of the mould ring and the exact placement.



Brevest investment materials for results that fit perfectly! The exact controllability of the material enables precise fitting of the expansion and therefore reduces expensive processing or repetitions. This enables stress-free working and saves time.



### Materials used

Wax-Lite surface tension reducing agent page 120 Silicone and wax surface tension reducing page 120 agent Mould release agent page 120 page 121 Metal muffle rings Fleece liners page 121 Wash primer for fleece page 121 Silicone muffle rings page 122 page 123 Vakuum Anrührsystem ecovac Brevest C+B Speed page 124 Brevest Rapid 1 page 124 page 125 Transfuser Investment marker page 126 Casting ring marker page 126 Golden booklet page 126 page 127 Brealloy C+B 270 Brealloy MK page 127 **Brealloy Lot** page 128 page 128 Brealloy flux Oxyd-Stop-PM page 129 Oxyd-Stop-NPM page 129 Oxyd-Stop Silver-Palladium alloy page 130 Heat absorbent paste page 130 Brecid pickling agent page 130



# Wax-Lite surface tension reducing agent



Alcohol-free surface tension reducing agent for bubble-free investing of wax patterns.

Wax-Lite surface tension reducing agent 750 ml REF 520 0100 8





Wax surfaces that are coated with the tension reducing agent allow flowing of the investment material into very small cavities of the model. This results in smooth, homogeneous surfaces and perfect occlusal surfaces. With the spray bottle micro-fine layers of Wax-Lite can be applied on the wax surface.

# Silicone and wax surface tension reducing agent



Improves the flow characteristics of plaster on silicone impressions.

Spraying on silicone and wax tension reducing agent will improve the flow characteristics of plaster for silicone impressions. The impression must be dry before the arch is poured.

Silicone and wax surface tension reducing agent 750 ml REF 540 0070 5




The fine spray head of the plastic spray bottle simplifies spraying uniform quantities of the agent.

Accessories:

Spray bottle, plastic sp 1 piece, 125 ml REF 540 0075 0



The spraying head of the spray bottle simplifies uniform wetting of the surface with silicone and wax surface tension reducing agent.



After the application of the agent onto the surface (left), the flow characteristics of the plaster have been clearly improved.



Silicone and wax surface tension reducing agent produces a homogeneous plaster surface. This will ensure precise dental work

## Mould release agent



Mould release agent 125 ml REF 520 TM12 5



Allows the investment material residue to be removed with from the silicone shape. To preserve and store the silicone moulds, simply spray on release agent and let dry.

Refill package 750 ml REF 520 TM75 0




# Metal muffle rings

Steel ring adapted to the shape of the arch

Semi-round steel rings allow positioning the castings outside the center of heat. The stress-free castings can be easily fitted and allow to continue working quickly.





Steel ring


REF

SX3 SX6 SX9

REF 360 ESR0 3 360 ESR0 6 360 ESR0 9

Base former for steel ring

SX3 SX6 SX9 360 ESS0 3 360 ESS0 6 360 ESS0 9



Steel rings for SX3, SX6 and SX9, compatible with all standard casting machines – familiar procedures can be used.



Metal-reinforced silicone base formers for high stability.

### Assortment

12 pieces

- 1 steel ring SX3, SX6, SX9 each
- 1 base former SX3, SX6, SX9 each
- 1 fleece liner SX3, SX6, SX9

20 ml wash primer for fleece 200 ml wash primer – refill package 125 ml mould release agent

REF 360 ESSE T

#### Accessories:



Mould release agent 125 ml REF 520 TM12 5 750 ml REF 520 TM75 0

### Fleece liners



These fleece liners do not absorb moisture from the investment material; hence uniform expansion of the investment material and high accuracy are ensured.



Precut fleece liners for sizes SX3, SX6 and SX9 for easier processing.

 Sizes
 SX3
 SX6
 SX9

 Quantity
 200 pieces
 100 pieces
 50 pieces

 REF
 360 ESV0 3
 360 ESV0 6
 360 ESV0 9

# Wash primer for fleece



Wash primer 20 ml REF 520 HG02 0



Refill package 200 ml REF 520 HG20 0



The wash primer to glue the fleece liners in the muffle ring avoids the penetration of investment material and simplifies devesting.



# Silicone muffle rings

Silicone rings made from addition-cured silicone are poor heat conductors. During setting of the investment material, the accumulated heat causes quick increase in temperature and leads to a higher final temperature. Stronger expansion movement is obtained. Tolerances are reduced and the precision of investment material is increased. The silicone sleeve ensures uniform expansion pressure, increases precision and produces reliable results. The high-quality silicone can be easily cleaned and features special durability.

A central aspect during cooling and solidification of the liquid molten mass is the position of the prosthetic object towards the outer wall of the muffle. bredent casting technique offers the suitable design of investment muffles to always position the object outside the center of heat.



The bridge is positioned in the center of heat.

Thanks to the mould design, all bridge moulds are positioned outside the center of heat to obtain a homogeneous casting structure.

Silicone sleeve


SX3 SX6 SX9
REF 360 SIM0 3 360 SIM0 6 360 SIM0 9



Base former for silicone sleeve

SX3 SX6 SX9

REF 360 SISO 3 360 SISO 6 360 SISO 9



### Sortiment

7 pieces

1 silicone sleeve SX3, SX6, SX9 each 1 base former SX3, SX6, SX9 each 125 ml mould release agent

REF 360 SISE T

### Accessories:



Mould release agent 125 ml REF 520 TM12 5 750 ml REF 520 TM75 0

### Vakuum Anrührsystem ecovac



### ecovac

# Precision-fit restorations obtained through optimal use of material properties.

The user-friendly and compact design simplifies work and reduces sources of errors. A powerful and maintenance-free vacuum pump, adjustable in two different levels (15 mbars, 200 mbars), ensures bubble-free mixing of materials and results in a perfect casting surface. Stirring time and speed can be adjusted continuously to allow correct processing of different materials.

ecovac (230 V) REF 140 0093 0

(Wall mounting, without mixing cup and base)

- 1 mains cable
- 1 spare filter
- 1 drilling template for wall mounting
- 4 screws and plugs for wall mounting

Accessories

Base ecovac, 1 piece REF 210 0045 0



### ecovac mixing spiral

The mixing spiral takes up the components to be mixed from all areas of the mixing cup and stirrs them horizontally and vertically. No unmixed materials will remain on the bottom of the mixing cup, which may cause different expansion of the material later on.

All features and components listed provide increased reliability, lead to improved fit when preparing dental restorations and avoid time-consuming reworking.

Mixing spiral,	50 ccm	REF 140 0R94 5
Mixing spiral,	250 ccm	REF 140 0R94 0
Mixing spiral,	750 ccm	REF 140 0R94 2
Mixing spiral,	1000 ccm	REF 140 OR94 3



### ecovac mixing cups

The smooth inner surface of the stainless steel mixing cup prevents any material or liquid residues from adhering to or depositing in scratches or undercuts. The conical shape ensures that material which has been taken up will flow back to the center of the mixing cup. Accordingly, the mixing ratio is retained exactly and better results can be achieved with minimal effort.

Mixing cup,	50 ccm	REF 140 0B94 5
Mixing cup,	250 ccm	REF 140 0B94 0
Mixing cup,	750 ccm	REF 140 0B94 2
Mixing cup,	1000 ccm	REF 140 0B94 3



Mixing cup, D (for the use in the Degussa mixing unit), 425 ml

REF 140 0B94 4



# Brevest C+B Speed



Brevest C+B Speed 50 bags, 160 g each REF 570 CBSO 8 125 bags, 160 g each REF 570 CBS2 0

Bresol Speed \*
1000 ml bottle
REF 520 000S 1
5000 ml
REF 520 000S 5

\* frost-resistant

# Assortment 25 bags Browest C. B. Space

Brevest C+B Speed 1000 ml Bresol Speed REF 570 CBS0 4 Very fine-grained, phosphate-bonded investment materials for crowns and bridges made of precious and non-precious metal alloys featuring outstanding reproduction of details.



At a room temperature of 21°C a processing time span of 4 to 6 minutes is obtained for bubble-free pouring out of casting rings.





Dosing bottle REF 520 0101 1



Dosing syringe 6 pieces REF 520 0101 2



Highly accurate and dimensionally precise large-span bridges are produced with Brevest M1 C+B and Brevest C+B Speed.

Exact control of concentrations for precision-fit dentures is possible with the frost-resistant mixing liquids Bresol N and Bresol Speed.

### Brevest Rapid 1



Rapid-heating, universal precision investment material for crowns and bridges as well as the entire field of CoCr work.

**Brevest Rapid 1**40 bags, 200 g each **REF 570 000R 8**100 bags, 200 g each **REF 570 00R2 0** 

Bresol R 1000 ml bottle REF 520 000R 1 5000 ml REF 520 000R 5 Brevest Rapid 1

REF 570 160R 8 125 bags, 160 g each REF 570 16R2 0

50 bags, 160 g each

### Assortment

20 bags Brevest Rapid 1 1000 ml Bresol R REF 570 0002 5



Fine grained, rapid-heating precision investment material for all largespan bridges, can also be used without casting rings.



Brevest Rapid 1 can be placed into the furnace at a temperature of 900 °C already 15 minutes after mixing.

### Accessories:



Dosing bottle REF 520 0101 1



Dosing syringe 6 pieces REF 520 0101 2



### Transfuser



### Gentle and safe transfer of investment compound and gypsum

While investing and casting an impression, quite frequently air gets entrapped – especially in narrow and inverse areas – which may cause undesired bubble formation. So far it was merely possible to utilize instruments and brushes for the preparation. But sharp-edged instruments bear the risk of damaging both modellation and impression. While using brushes, it is possible to transfer excess fluids, and if the brush is rather dry it may withdraw some moisture. Both may affect the

expansion behavior of investment compounds.

The curved and flexible silicone tip of the Transfuser allows
a gentle, damage- and bubble-free application prior to standard filling procedures. The smooth and dense surface provides excellent gliding characteristics for gypsums and investment compounds, and thus a bubble-free material



Highly flexible, soft tip and thus no risk of damaging the wax model.



Narrow areas are thoroughly filled by gentle application without modellation damage.



Optimal material flow and bubble-free condensation with the Transfuser



Smooth and gentle filling of extreme hollows with gypsum in case of silicone impressions.



No risk of damage during the filling process.



Each area is well accessible and allows for optimal and bubble-free material distribution.



### Investment marker



Helps with the positive identification of investment muffles.

Investment marker REF 330 0115 0



The necessary information is noted down quickly and easily.



The marker can be clearly read on all investment materials up to 1100 °C.

# Casting ring marker



Casting ring marker 1 marker with 4 spare cartridges REF 330 0115 1

Refill package with 8 cartridges REF 330 0115 2



Comprehensive information even on the smallest of rings.



The positive identification is assured up to 950° C.


For correct identification of casting rings.

### Golden booklet



Golden booklet DIN A 6 REF 610 0020 0

Thanks to the clear and simple structure of the golden booklet, reliable stock-keeping of precious metal alloys is ensured. The booklet simplifies the control and provides a quick survey on the consumption of alloys.





Golden booklet DIN A 4 REF 610 0010 0

# Brealloy C+B 270



Brealloy C + B 270 cylinder, each 6.3 g 50 g REF 500 CB05 0 200 g REF 500 CB20 0 500 g REF 500 CB50 0 1000 g REF 500 CB00 0 Ceramic bonding alloy with a hardness of 270 HV 10 which can be milled easily. Brealloy C + B 270 is free from nickel, beryllium and gallium. The alloy corresponds to the standard DIN 13912: 1996 for non-precious metals and DIN EN ISO 9693: 1995 for metal-ceramic systems.

# Composition (in mass-%)

Cobalt	66
Chromium	20
Molybdenum	6
Wolfram	6
Silicone	0.9
Carbon	0.02
Manganese	0.7

# Physical values (guide values)

Density (g/cm<sup>3</sup>)

Vickers hardness (HV 10)	270
Solidus point (°C)	1280
Liquidus point (°C)	1350
Casting temperature (°C)	1450
0.2% proof stress (MPa)	600
Modulus of elasticitiy	
(MPa) approx.	200,000
Strain at break (%)	10
Espansion coefficient	
(WAK 20-600 °C)	14.4 μm/mk

8.4



Milling technique: Brealloy C + B 270 can be milled perfectly.



Partial crowns made of Brealloy C + B 270: slender and precise.



Attachment technique with Brealloy C + B 270: precision in the one-piece casting method.

## **Brealloy MK**



brealloy MK 50 g REF 500 MK05 0 200 g REF 500 MK20 0 500 g REF 500 MK50 0 1000 g REF 500 MK00 0 Solderable CoCr based alloy for crown and bridge technology. The low hardness permits easy processing. It is therefore ideal for milling. brealloy MK is free from nickel, beryllium and gallium.

# Composition (in % of mass)

cobalt	65
chrome	20
molybdenum	6.5
tungsten	6.5
silicon	0.8
manganese	0.8
iron	<0.5
carbon	<0.1

# Physical properties (guide values)

density (g/cm³)	8.4
Vickers hardness (HV 10)	265
solidus point (°C)	1280
liquidus point (°C)	1350
casting temperature (°C)	1420
0.2% yield point	480
tensile strength (N/mm2)	790
E-module (mPa)	190,000
Elongation at break (%)	10
Coefficient of expansion	
(WAK 20-600°C) 14	.8 μm/mK



The low hardness facilitates milling, shaping and polishing.



ldeally suited for longer bridges.



Even inlays with fine spring tips can be easily produced.



The low hardness creates ideal conditions for milling and is thus highly suitable for all types of attachments.



# Brealloy Lot



Brealloy Lot 7 g REF 500 0001 0 Solder especially matched with CoCr alloys for chrome cobalt and ceramic bonding techniques to avoid the formation of galvanic elements and undesired reciprocal action with the ceramic material.

# Brealloy flux



Brealloy flux 8 g REF 500 0001 1 **Suitable for all CoCr alloys,** supports the flow characteristics of the solder.

# Oxyd-Stop-PM



Prevents polished, precious and semi-precious alloy

or in a furnace and while firing porcelain.

Oxyd-Stop-PM 20 ml REF 520 0065 0

Thinner 20 ml REF 520 0067 0




Oxyd-Stop-PM maintains the high luster. No refinishing is required.



Refinishing oxidized surfaces reduces the material and may adversely affect the precision of fit.

### Oxyd-Stop-PM - Range of Applications



Oxyd-Stop-PM prevents precious metal crowns and bridges from oxidizing during soldering. Also suitable for Wiron 88.



Safeguards high grade milled attachement restorations against oxidization while porcelain is being fired.



Perfect for post-soldering porcelain bridgework. The high luster metal surface is protected.



Shake the Oxyd-Stop-PM well and brush it onto the areas to be protected. Allow briefly to dry. The surface is now optimally protected against oxida-



Always apply a fresh coat of Oxyd-Stop-PM after the casting has been heated. Only then is optimum protection against oxidation provided.



After firing or soldering, the protective layer is easily removed with hot water or steam cleaner.



Refinishing is reduced to a minimum, which saves time and enhances the quality.

## Oxyd-Stop-NPM



Prevents oxidation while soldering all chrome cobalt and non-precious alloys. This reduces refinishing to a minimum and saves time.

Oxyd-Stop-NPM 2 x 50 ml REF 520 0061 0



Oxyd-Stop-NPM protects non-precious bridges or single crowns against oxidation. The high luster surfaces are maintained.



Oxyd-Stop-NPM is perfect for use on chrome cobalt frameworks. Repairs and extensions can then be completed even faster.



When using Oxyd-Stop-NPM, the metal surface has the same high luster after soldering as it had after polishing.



Apply Oxyd-Stop-NPM direct from the tube onto the areas being protected. The dispensing tip facilitates this procedure.



Use an instrument to spread the Oxyd-Stop-NPM evenly.



Oxyd-Stop-NPM effectively prevents oxidation while soldering with a flame. High luster areas remain free of oxides and retain their precision of fit



After soldering, Oxyd-Stop-NPM can be brushed off easily and quickly under running water. This saves time.



The surfaces only require minimal polishing to restore their high luster. Only the exposed crown has oxidized.



# Oxyd-Stop Silver-Palladium alloy



Oxyd-Stop Silver-Palladium alloy 20 ml REF 520 0033 0



Apply a thin coat of Oxyd-Stop Silver-Palladium onto the areas to be protected and let it dry. A white protective layer is formed.



Oxyd-Stop Silver-Palladium avoids the oxidation of the soldering object while being soldered with a flame or in a furnace.



Remove Oxyd-Stop with the sandblasting unit, in the ultrasonic unit or with a pickling agent. Repolishing of the protected areas is reduced to a minimum degree.

# reduced alloys. No reworking required and thus time is saved.

Avoids oxidizing of silver-palladium and

# Heat absorbent paste



There is no better method of protecting against heat.

Heat absorbent paste 250 g REF 540 0020 0



paste very close to the joint.

Apply the heat absorbent





This paste will not melt when warmed.

# Brecid pickling agent



For the first time ever, precious metal retains its high luster after pickling

Brecid pickling agent 3 x 75 g REF 520 0099 0



being soldered onto a full gold crown, polished to a high-luster.

Shows a contact area



Pickling precious metals in Brecid removes the oxide yet leaves the high-luster intact therefore saving time.



Although conventional pickling agent does remove the oxide, it also dulls the high luster achieved with great effort.



Anchoring elements for secure hold of the prosthesis gives patients a better quality of life. bredent therefore offers a variety of anchoring elements for various indications, from stud and rod attachments, various bar variants, locks and screws for dentist-removable dental prostheses.







R,	ochures for dentists and patients	
וט	Help to make you and your dentist more	
	successful with these brochures!	132
m	portant information1	33
'n	terlock	
	Interlock	
Vá	vks-oc rs 2,2 abutments	
	vks-oc extracoronal use - fixation in the	100
	CoCr structure	140
	vks-oc uni	142
	vks-oc exchangeable stud	148
	vks-oc/sg exchangeable stud with adhesive	1 = 0
/-	sleeve1sleeve	150 1 <b>51</b>
V C	vks-oc extracoronal use	152
	The complete product range	
	vks-sg matrix housing	
	vks-sg bar patrix	
Ro	od attachments 1 Vario-Soft 3	
	Vario-Soft 3 sv	
	Vario-Soft 3 matrix housing	
	Vario-Soft 3 sv bridge-sectioning attachment	
	Vario-Soft 3 mini	166
	Vario-Soft 3 mini sv	
	Vario-Soft 3 zircon sv and zircon sv mini	
Re	Vario-Soft 3 conicalbridge1	
•	Inverto Plus	
	Activatable frictions cylinder	
	Vario Compress 1	
	Vario Compress 2	
	Stud fixator  Cylindrical attachment zg	
R۶	ars	1/0
	Bar connections	183
	Vario-Soft-Bar-Pattern vsp	
	Vario-Soft-Bar vss	
	Wax bars	190
LU	Lock	193
	Swivel-type lock sr	
	Swivel-type lock src	
	Locking Pin bs 1	
	KS-lock	
	Locking Pin Easy-Snap E	201 202
	Locking Pin Easy-Snap A	
	Locking Pin Snap System	
	Locking Pin Snap E	
	Locking Pin Snap	
	Locking Pin Snap A	
50	Locking Pin activatable	
,	Universal screwdriver set	
	Universal screwdriver set for contra-angles	
	Screwdriver	
	Milling and drilling oil	219
	Transverse fixation	220
	Security-Lock	
	Security-Lock-Ceramic	
	Security-Lock adhesive sleeve	224
	Friction Splint FS1	
D	Tool set for individual screw connections	228
or	idge Sectioning Attachment Bridge Sectioning Attachment oc	230
	Custom Bridge Sectioning Attachment	
	Prefabricated screwing set	

## Brochures for dentists and patients


# Help to make you and your dentist more successful with these brochures!



### Vario-Soft product range

The different types of attachments of the Vario-Soft attachment group and their advantages are illustrated in detail in this brochure. This will provide dentists with a better survey on the variety of attachments so that patients will receive more qualified advice and high quality restorations.


REF 000 041G B




### A healthy mouth - for your whole life

The biofilm in which the bacteria are arranged in the mouth, constitutes the main risk to oral health in the long term. In particular if the pathogenic bacteria become out of control and inflammation occurs. HELBO treatment makes a scientifically-proven method of fighting these pathogenic bacteria available to the dentist - the bacteria are stained and killed using a suitably coordinated Theralite laser. This brochure explains the combination and possibilities of the treatment to the patient in a straight forward manner.

REF 000 484G B



bredent medical offers various patient brochures for implant restoration, in order to meet the different requirements of patient groups.



This brochure gives young patients who are missing one or a few teeth an overview of which options implantology offers for them and what they must take into consideration in the case of implant restoration. REF 000 421G B



### tempting ...?

Patients who are on the verge of becoming edentulous are a special customer group. They feel they are too young for a total prosthesis and want a fixed restoration as soon as possible. The options of immediate restorations with the SKY fast & fixed treatment concept are set out in detail in these patient brochures. REF 000 342G B



### ...in the thick of it...

Patients who are already edentulous have different requirements again. These patients can already be helped a great deal by fixing their prosthesis to implants. These patient brochures show the possibilities of securely fixing removable restorations so that the patient is able to eat comfortably and socialise again happily.

REF 000 422G B

The brochures are available free of charge!





### Important information for users of bredent attachments!

To ensure trouble-free and lasting function of the attachments, the stable position of the removeable denture is of utmost importance. A circumferential shoulder with parallel milled interlock at the abutment crown and a corresponding shear distributor at the removeable restorations are essential elements and indispensable. Tilting movements of the denture must be avoided since they result in frequent locking and unlocking of the snap attachments and – in conjunction with crystalline deposits – may cause premature wear and thus affect the proper function of the attachments.

bredent Research Information

# Vario-Stud-Snap attachment vks oc + sg

Latest findings have shown that in a very limited number of cases deposits may be formed on natural teeth, dentures and fixed restorations in the oral environment.

If, due to insufficient oral hygiene, these crystals are not removed, some exceptional cases of inclusion of these crystals in the surface of the plastic matrix might result.

This leads to an abrasive effect on the stud of the patrix resulting in the possible loss of snap. Very rare cases of this unexplained and previously unknown phenomenon have been reported for the Stud-Snap attachments sold (1 of 5000 patients).

Accordingly, we recommend the exclusive use of hard alloys and to clean the teeth, the denture and the fixed restoration two times a day as well as to have them regularly checked by the dentist. To ensure perfect function of the Vario-Stud-Snap attachment it is necessary that the patient acquires the snap point with his finger when inserting the denture and locks it by pressing on it with his finger.



### Interlock



Parallel- and 2°-Interlock made of high-melting special wax.


After determining the direction of insertion, the copings are produced (wax or resin).



- Fast and reliable attaching of the Interlock
- No damage to the die when drilling the Interlock
- Only drill with a groove bur
- Defined wall thickness of just 0.4 mm

The Interlock is integrated into the model using the paralleling mandrel. Then the circular groove is modelled and milled.

### Interlock parallel




8 pieces REF 430 0736 9



Paralleling mandrel Interlock, parallel 1 piece REF 360 0116 6

### Interlock 2°



8 pieces REF 430 0736 8



Paralleling mandrel Interlock 2° 1 piece REF 360 0116 5



Fast and correct attaching of the Interlock with shear distributor ensures quick reworking.



The Interlock is only redrilled (reamed) with the groove bur F538 2H 10. The shear distributor is processed in the usual way.



A reliable and durable removeable denture is obtained, if the correct Interlock is selected.



Product	REF	Ø	Width	Height
Interlock 0°	430 0736 9	0.9 mm	2.2 mm	6.0 mm
Interlock 2°	430 0736 8	1.4 mm	1.0/1.4 mm	6.0 mm



### vks-oc

Stud attachments must be processed very precisely for functionally-secure application. It must be possible to guarantee the positional stability of the prosthesis, in order to prevent overloading caused by micro movements of the matrix and the patrix. The attachment is only designed for secure hold of the prosthesis, it is not, however, for permanent weight-bearing. In the case of extracoronal use, a shear distributor with integrated Interlock must be available. The required stability and force distribution are thus achieved.

Due to its snap function, the stud attachment offers patients more security, as a clear snapping in of the prosthesis is noticeable.

#### Indications

- Root caps
- Extracoronal applications
- Retaining element for bars
- Implant prosthetics (only vks-oc rs)
- Plastic fitting
- Model casting fixation

### **Product variations**

- Vario-Kugel-Snap vks-oc 1.7 uni for root caps and bar solutions
- Vario-Kugel-Snap vks-oc 2.2 uni for root caps and bar solutions
- Vario-Kugel-Snap vks-oc 1.7 for extracoronal application
- Vario-Kugel-Snap vks-oc 2.2 for extracoronal application
- Vario-Kugel-Snap vks-oc rs 2.2 for various implant systems
- Vario-Kugel-Snap vks-oc 1.7 exchangeable stud
- Vario-Kugel-Snap vks-oc 2.2 exchangeable stud

#### **Properties**

- Cost-effective purchase
- Time-saving and secure processing
- Reduction of the alloy variety for high biocompatibility
- Cost-effective solution provided by castable plastic patrices

### Materials

- Patrices
- Combustible plastic
- Castable precious metal alloy
- Matrices
- Biocompatible thermoplast

### Dimensions

The Vario-Kugel-Snap vks-oc stud head attachments are available in sizes 1.7 and 2.2

### The matrices have a retention of









vks-oc 1.7/2.2 uni for root caps and bar solutions



vks-oc 1.7/2.2 for extracoronal application



vks-oc rs 2.2 for various implant systems



vks-oc 1.7 / 2.2 exchangeable stud

### Vario-Kugel-Snap vks-oc



### vks-oc rs 2.2 abutments

### nine abutments suitable for Brånemark, 3i (except 6 mm) and Steri-Oss hex-lock

The patrix of the stud attachment has been produced with utmost precision in an industrial process. It consists of titanium, grade 5. This material guarantees precision of fit, high durability and biocompatibility.

In the area of the head the screw of the vks-oc stud abutment features a recess. This ensures gap-free fit of the abutment on all common implants with hexagon. Only the suitable diameter and the desired distance height have to be observed.





vks-oc rs Ø 2.2 mm abutment Ø 4 mm, distance height 2 mm impression matrix 1 piece each REF 460 0004 2



vks-oc rs Ø 2.2 mm abutment Ø 5 mm, distance height 2 mm impression matrix 1 piece each REF 460 0005 2



vks-oc rs Ø 2.2 mm abutment Ø 6 mm, distance height 2 mm impression matrix 1 piece each REF 460 0006 2



vks-oc rs Ø 2.2 mm abutment Ø 4 mm, distance height 4 mm impression matrix 1 piece each REF 460 0004 4



vks-oc rs Ø 2.2 mm abutment Ø 5 mm, distance height 4 mm impression matrix 1 piece each REF 460 0005 4



vks-oc rs Ø 2.2 mm abutment Ø 6 mm, distance height 4 mm impression matrix 1 piece each REF 460 0006 4



vks-oc rs Ø 2.2 mm abutment Ø 4 mm, distance height 6 mm impression matrix 1 piece each REF 460 0004 6



vks-oc rs Ø 2.2 mm abutment Ø 5 mm, distance height 6 mm impression matrix 1 piece each REF 460 0005 6



vks-oc rs Ø 2.2 mm abutment Ø 6 mm, distance height 6 mm impression matrix 1 piece each REF 460 0006 6



vks-oc rs Ø 2.2 mm implant analogue Ø 4 mm 2 pieces REF 460 0000 4



vks-oc rs Ø 2.2 mm implant analogue Ø 5 mm 2 pieces REF460 0000 5



vks-oc rs Ø 2.2 mm implant analogue Ø6 mm 2 pieces REF 460 0000 6





Screwdriver is elbow 1 piece REF 460 0001 0



Screwdriver is manual, short 1 piece REF 460 0001 1



Angle measuring device 1 piece REF 460 0010 0



vks-oc rs Ø 2.2 mm axle abutment 2 pieces REF 460 0010 2 8 pieces REF 460 0010 8




Impression matrix 2 pieces REF 460 0000 2 8 pieces REF 460 0000 8



The transfer matrices are now in the total impression. Retention grooves ensure the firm position of the transfer matrices



After successful osseointegration of the implant. the mucous membrane is opened and the suitable vks-oc abutment is screwed in. The torque must be limited to 30 Ncm by using a vks-oc screwdriver.





Abutments in the distance heights of 2, 4 and 6 mm are available. Distance sleeves are not required. The stud abutments can also be used as gingiva former. Prior to taking the impression, place the transfer matrices onto the stud abutments.



in the impression mate-



Important information about processing of attachments, see catalogue page 133

### vks-oc rs 2,2 abutments



Prior to the fabrication of the model, laboratory implants are inserted into the tranfer matrices. This way precise transfer of the position of the implants is guaranteed.



The model precisely shows the position of the implants. Normally, the implants have not been inserted in an entirely parallel position to each other so that divergences or convergences between several stud abutments result.



The vks-oc system accepts maximum deviations of 15  $^{\circ}$  to the planned angle of insertion. To determine the deviation from the vertical axle, axle abutments are placed onto the abutment analogue.



A special angle measuring device represents the range of tolerance of the vks-oc system. In this example the deviation of the implants lies within the range of tolerance. Accordingly, supply with the vks-oc rs system is possible.

### Rigid matrices for fixation in an acrylic denture





Rigid matrices vks-oc rs Ø 2.2 mm green - reduced softsnap-in effect 8 pieces REF 440 0070 8

50 pieces REF 440 0075 0

Rigid matrices



Rigid matrices vks-oc rs Ø 2.2 mm red - high soft-snap-in effect 8 pieces REF 440 0090 8

Accessories:

50 pieces REF 440 0095 0

### Assortment

- 13 pieces
- 2 Rigid matrices, green
- 2 Rigid matrices, yellow
- 2 Rigid matrices, red
- 2 Matrix housing for fixation in acrylics
- 2 Screwdriver is
- 1 Angle measuring device
- 1 Matrix inserting instrument
- 1 vks paralleling mandrel

REF 440 0066 4



Fig. 1:1

vks-oc rs Ø 2.2 mm yellow - medium softsnap-in effect 8 pieces REF 440 0080 8 50 pieces

Matrix housing for fixation in acrylics vks-oc rs Ø 2.2 mm

REF 440 0085 0

2 pieces REF 440 0030 2

8 pieces REF 440 0030 8



Matrix inserting instrument vks-oc rs Ø 2.2 mm REF 360 0116 1




Paralleling mandrel universal 2 vks-oc rs Ø 2.2 mm REF 360 0116 0



Matrix pliers vks-oc Ø 2.2 mm + zg REF 310 0000 6



The matrices should be placed onto the divergent abutments so that parallel direction of insertion is ensured. For this purpose the procedure in figures 3 – 9 must be adhered to.



Orientation of the matrices at the axles of the implants would result in a divergent position of the matrices. The function of the matrices would be affected.



Press the yellow plastic matrix (resilient or rigid) into the metal housing using the inserting instrument. The paralleling mandrel is used to achieve parallel placement of all matrices on the abutments.



Fix the position of the matrices with plaster; fill undercuts up to the jaw ridge. The plaster base ensures correct position of all matrices during subsequent working steps.



To set up the teeth, use an acrylic base plate. The matrices are integrated into the base plate using a small amount of acrylic. The teeth are set up on this plate.



For try-in the matrices and the set-up can be removed from the model. Due to the plaster base (cf. fig. 5) the parallel position of the matrices is ensured so that it can be reproduced easily.



Prior to the completion, remove the matrices from the base plate and place them on the abutments. All techniques (e.g. flask pressing technique, casting technique) can be used for completion.



The matrix housing polymerized in the acrylic material: The plastic matrix can be removed with the matrix pliers. The rigid matrices can be exchanged for the resilient ones (and vice versa) at any time.



# Vario-Kugel-Snap vks-oc

# vks-oc rs 2,2 abutments

### Rigid matrices for fixation in a CoCr structure



vks-oc rs Ø 2.2 mm Processing of the matrix on the abutment









Accessories:



Rigid matrices vks-oc rs Ø 2.2 mm green – reduced softsnap-in effect 8 pieces REF 440 0070 8 50 pieces

REF 440 0075 0
Rigid matrices
vks-oc rs Ø 2.2 mm
red – high soft-snap-in
effect
8 pieces
REF 440 0090 8

50 pieces
REF 440 0095 0
Wax matrix housing
vks-oc rs 2.2 mm
8 pieces
REF 440 0100 8
50 pieces
REF 440 0105 0

Matrix pliers vks-oc Ø 2.2 mm + zg 1 piece REF 310 0000 6





Fig. 1:1

yellow – medium softsnap-in effect 8 pieces REF 440 0080 8 50 pieces REF 440 0085 0 Matrix housings vks-oc rs 2.2 mm for glueing or laserwelding 2 pieces

Rigid matrices

vks-oc rs Ø 2.2 mm



Duplicating matrix vks-oc rs 2.2 mm 8 pieces REF 440 0110 8

REF 440 0020 2

### Assortment

13 pieces
2 Rigid matrices, green
2 Rigid matrices, yellow
2 Rigid matrices, red
2 Matrix housings

2 Rigid matrices, red 2 Matrix housings 2 Wax matrix housings 2 Duplicating matrix

2 Screwdriver is1 Angle measuring device

 Matrix inserting instrument
 vks paralleling mandrel


REF 440 0066 5

Important information about processing of attachments, see catalogue page 133

# vks-oc rs 2,2 abutments




The matrices (duplicating matrix housing) should be placed onto the divergent abutments so that parallel alignment is ensured. For this purpose the procedure in figures 3 – 13 must be adhered to.



Orientation of the matrices at the axles of the implants would result in a divergent position of the matrices. The function and the durability of the matrices would be affected.



To ensure parallel position of the duplicating matrices, they are placed onto the patrices using the paralleling mandrel. Then the position of the matrices is fixed using plaster.



Fill undercuts between duplicating matrix and jaw ridge with plaster. The plaster base ensures the correct position of the matrices during all subsequent working steps.



To produce a chrome cobalt framework, the model is blocked out and duplicated in the usual way. The duplicating matrix housing guarantees the correct size of the glueing gap.



In order to wax up the chrome cobalt framework over the matrix housings, special wax matrix housings are available. They ensure correct thickness of the chrome cobalt frame in the area of the matrices.



The chrome cobalt framework is cast, finished and polished according to standard criteria. The cavities to hold the glueing matrices should only be cleaned with a sandblasting unit (do not grind).



Press plastic matrices into the metal matrix housing using the inserting instrument. To glue the metal matrix housings, place them on the abutments.



The plaster bases, which were built up by parallel placement of the duplicating matrices (cf fig. 4) ensure parallel position of the metal matrix housing prior to glueing.



Clean the recesses for the glueing matrices with the sandblasting unit and fill with DTK-adhesive. Caution: Residues in sandblasting units with circulation affect the adhesion of the DTK-adhesive.



Press the chrome cobalt framework onto the glueing matrices which have been fixed in the correct positions. Press excess glue out of the recesses.




Then prepare wax set-up, try-in and complete acrylic denture in the usual way. Resilient and rigid matrices are inter-exchangeable at any time.

### **Dimensions**



	Product	REF	Thread	Distance height	Width	Ø Stud
6	vks-oc rs Ø 2.2 mm	460 0004 2	M2 x 0.4	2 mm	4 mm	2.2 mm
7		460 0004 4	M2 x 0.4	4 mm	4 mm	2.2 mm
		460 0004 6	M2 x 0.4	6 mm	4 mm	2.2 mm
		460 0005 2	M2 x 0.4	2 mm	5 mm	2.2 mm
		460 0005 4	M2 x 0.4	4 mm	5 mm	2.2 mm
		460 0005 6	M2 x 0.4	6 mm	5 mm	2.2 mm
		460 0006 2	M2 x 0.4	2 mm	6 mm	2.2 mm
		460 0006 4	M2 x 0.4	4 mm	6 mm	2.2 mm
		460 0006 6	M2 x 0 4	6 mm	6 mm	2.2 mm



### vks-oc extracoronal use - fixation in the CoCr structure



The extracoronal vks-oc must always be used in conjunction with a milled shear distributor. This way optimal transfer of resulting forces onto the anchor tooth is ensured.

vks-oc is available in two different angles to allow optimal adaptation to the course of the gingiva.



Ø 1.7 mm, 30° 8 pieces REF 430 0734 5 50 pieces REF 430 0734 6

Patrix vks-oc



Patrix vks-oc Ø 2.2 mm 8 pieces REF 430 0539 0 50 pieces REF 430 0556 0



Patrix vks-oc Ø 1.7 mm, 60° 8 pieces REF 430 0734 7 50 pieces REF 430 0734 8



vks Paralleling mandrel oc/sg 1 piece ph-vks 1.7 REF 430 0677 0 ph-vks 2.2

REF 360 0113 0

### Assortment

22 pieces

### Vario-Stud-Snap vks-oc 1.7 30°/60°

- 4 Blocking out discs oc 1.7
- 1 Inserting instrument vks 1.7
- 4 Matrices oc 1.7 each yellow, green, red
- 1 Paralleling mandrel metal ph-vks 1,7
- 2 Patrices oc 1.7 each 30° + 60°

REF 430 0734 9

### Assortment

REF 430 0531 0

Assortment

12 pieces

12 pieces

### Vario-Stud-Snap vks-oc 2.2

Vario-Stud-Snap vks-oc 2.2

- 2 Blocking out discs oc 2.2
- 1 Inserting instrument vks 2.2
- 2 Matrices oc 2.2 each yellow, green, red
- 2 Patrices oc 2.2
- 1 Paralleling mandrel metal ph-vks 2.2

REF 430 0531 0

### Assortment

10 pieces

### Vario-Stud-Snap vks-oc 2.2

- 2 Blocking out discs oc 2.2
- 2 Matrices oc 2.2 each yellow, green, red
- 2 Patrices oc 2.2

REF 430 0534 0



At the beginning a crown is waxed up in the usual way and a milled shear distributor with groove is prepared in wax.



Select the suitable vks-oc according to the course of the papillae and bring it into the correct position using the paralleling mandrel.



Fix vks-oc at the crown using hot wax.



The transition zone of vks-oc/crown must be coated richly with hot wax. vks-oc patrices consist of burnout plastic. They are cast together with the crowns.



The one-piece castin failitates processing. After casting, the vks-oc is only slightly polished to high lustre using a buff.



The extracoronal vks-oc patrices are assembled in the chrome cobalt framework in a very easy



Product	Ø	Angle	Length	Height
Patrix vks-oc	Stud 1.7 mm	30°	5.8 mm	3.9 mm
	Stud 1.7 mm	60°	6.6 mm	6.6 mm
	Stud 2.2 mm	-	6.7 mm	7.5 mm



Important information about processing of attachments, see catalogue page 133

### vks-oc extracoronal use - fixation in the CoCr structure



Accessories:



Matrices vks-oc 1.7 red – high soft-snap-in effect 8 pieces

REF 430 0656 0 50 pieces REF 430 0657 0



Matrices vks-oc 2.2 red – high soft-snap-in effect

8 pieces REF 430 0546 0 50 pieces REF 430 0548 3



Matrices vks-oc 1.7 yellow – medium softsnap-in effect

8 pieces REF 430 0659 0 50 pieces REF 430 0658 0



Matrices vks-oc 2.2 yellow – medium softsnap-in effect

8 pieces REF 430 0545 0 50 pieces REF 430 0549 0



Matrix inserting instrument vks-oc Ø 1.7 mm 1 piece REF 430 0621 0

Metal transfer

vks-oc Ø 1.7 mm

patrices

8 pieces REF 430 0662 0



Matrices vks-oc 1.7 green- reduced snap 8 pieces REF 430 0655 0 50 pieces REF 430 0654 0



Matrices vks-oc 2.2 green- reduced snap 8 pieces REF 430 0544 0 50 pieces REF 430 0548 4



Matrix inserting instrument vks-oc Ø 2.2 mm 1 piece REF 430 0548 0



Blocking out discs vks-oc 1.7 8 pieces REF 430 0652 0 50 pieces

REF 430 0653 0



Blocking out discs vks-oc 2.2 12 pieces REF 430 0540 0 50 pieces REF 430 0548 5



Metal transfer patrices vks-oc Ø 2.2 mm 8 pieces REF 430 0548 2

### vks-oc Ø 1.7 mm and vks-oc Ø 2.2 mm: assembly in chrome cobalt framework



Block out to the basal direction starting from the blocking out disc. This way the perfect recess in the chrome cobalt framework to hold the matrix is obtained.



Then produce duplicate with chrome cobalt investment material.



Wax pattern of the planned chrome cobalt supply: the matrix is coated with a wax layer (thickness approx. 0.4 mm).



The completed chrome cobalt framework is ready for the assembly of the matrix with the inserting instrument.




The matrix is mounted with the special inserting instrument. Retention is ensured due to the conical outer shape. To exchange the matrix use a round bur or the matrix pliers.

	Product	Ø	Height
	Matrix vks-oc 1.7	2.7 mm	2.0 mm
	Matrix vks-oc 2.2	3.3 mm	2.7 mm
C	Blocking out disc vks-oc 1.7	2.8 mm	0.4 mm
	Blocking out disc vks-oc 2.2	3.5 mm	0.4 mm



# Vario-Kugel-Snap vks-oc

# vks-oc uni



### For root caps and bars.

vks oc uni in burnout plastic are cast together with the root cap. They can be processed easily and are particularly biocompatible since there is no electrochemical potential difference caused by a different alloy.

vks-oc uni are also available in a cast-on, highmelting alloy. They are particularly precise since reworking after casting is no longer required.



Patrix vks-oc uni Ø 1.7 mm HL-patrix cast-on 2 pieces REF 430 0701 0



Patrix vks-oc uni Ø 1.7 mm 8 pieces REF 430 0676 0 50 pieces REF 430 0675 0



Patrix vks-oc uni Ø 2.2 mm 8 pieces REF 430 0538 0 50 pieces REF 430 0550 0



Patrix vks-oc uni Ø 2.2 mm HL-patrix cast-on

2 pieces REF 430 0700 0

#### Accessories:



vks Paralleling mandrel oc/sg 1 piece ph-vks 1.7 REF 430 0677 0 ph-vks 2.2 REF 360 0113 0



Wax bars wstg 1.6 1.6 x 8 x 50 mm approx. 65 pieces REF 430 0265 0

Wax bars wstg 1.9 1.9 x 4 x 50 mm approx. 120 pieces REF 430 0266 0

Wax bars wstg 2.2  $2.2 \times 6 \times 50 \text{ mm}$ approx. 65 pieces REF 430 0267 0



Paralleling mandrel 1.9 - 2.2 for wstg 1.9 - 2.2 1 piece REF 430 0270 0



	Product	REF	Ø	Height
	Patrix vks-oc uni/HL	430 0675 0	Stud 1.7 mm	2.2 mm
		430 0550 0	Stud 2.2 mm	3.2 mm

Important information about processing of attachments, see catalogue page 133

### vks-oc uni



Root cap is modelled in the usual way. Bring the vks-oc uni into the correct position using the paralleling mandrel and fix with hot wax.



Remove paralleling mandrel and apply rich coat of hot wax onto the transition zone of vks-oc uni/root cap. The one-piece casting facilitates processing.



Casting is carried out according to standard criteria. After casting, the vks-oc unit is only slightly polished to high lustre using a textile buff.



The blue blocking out disc is placed onto the patrix below the equator. Plug the plastic matrix into the metal matrix housing using the inserting instrument.



Press the metal matrix housing with the plastic matrix onto the patrix. The blocking out disc ensures parallel position of the matrix.



For try-in, fix the metal matrix at the acrylic base plate using a small amount of acrylic.



Matrix housing in the set-up. The picture shows that only little space is required. For try-in, remove the blocking out discs.



For completion, place on the blocking out disc and cover root cap with liquid silicone. Do not cover the occlusal part of the stud.



Press the metal matrix housing with integrated plastic matrix into the silicone that is still soft.



The denture is completed in the usual way after the silicone cover has hardened.



Basal view of completed denture: To change the snap friction, remove the plastic matrix with a round bur and insert a different plastic matrix.

### Bar variation



Place a vks-oc uni in the correct position on a completely waxed-up bar using the paralleling mandrel and fix with hot wax.



Remove paralleling mandrel and apply wax onto the transition zone of vks-oc uni/wax bar.



Casting is carried out in the usual way. vks-oc uni is only slightly polished to high lustre using a textile buff.



# Vario-Kugel-Snap vks-oc

### vks-oc uni



Rigid matrixes for fixation in an acrylic denture.





Fig. 1:1

#### Matrices red

vks-oc Ø 1.7 mm high snap 8 pieces REF 430 0656 0

vks-oc Ø 2.2 mm 8 pieces REF 430 0546 0

Matrices yellow vks-oc Ø 1.7 mm medium snap 8 pieces REF 430 0659 0 vks-oc Ø 2.2 mm 8 pieces

REF 430 0545 0



#### Matrices green

vks-oc Ø 1.7 mm reduced snap 8 pieces REF 430 0655 0

vks-oc Ø 2.2 mm 8 pieces REF 430 0544 0



Metal matrix housing mmg vks-oc Ø 1.7 mm 2 pieces REF 430 0697 0

8 pieces

REF 430 0661 0

mmg vks-oc Ø 2.2 mm 2 pieces

REF 430 0696 0

8 pieces REF 430 0547 0

#### Accessories:



Matrix inserting instrument vks-oc Ø 1.7 mm REF 430 0621 0



Metal transfer patrices vks-oc Ø 2.2 mm 8 pieces REF 430 0548 2



Matrix housings tmg vks-oc Ø 1.7 mm 2 pieces REF 430 0699 0 vks-oc Ø 2.2 mm 2 pieces REF 430 0698 0



Metal transfer patrices vks-oc Ø 1.7 mm 8 pieces REF 430 0662 0



Matrix inserting instrument vks-oc Ø 2.2 mm 1 piece REF 430 0548 0



Blocking out discs vks-oc Ø 1.7 mm REF 430 0652 0 vks-oc Ø 2.2 mm 12 pieces REF 430 0540 0

#### Assortment Vario-Stud-Snap vks universal 1.7

14 pieces

- 2 Matrices each red, yellow, green
- 2 Metal matrix housings
- 2 Blocking out discs
- 2 Patrices
- 1 Matrix inserting instrument
- 1 Paralleling mandrel

REF 430 0674 0

Note: Parallel alignment of vks-oc matrices (only rigid) and vks-oc rs matrices (either rigid or resilient) on root caps is always carried out using a blocking out disc. The processing methods of vks-oc matrices and vks-oc rs matrices are identical. vks-oc and vks-oc rs require different blocking out discs, matrix housings and inserting instruments.

Additional assortments can be found in the price list!

#### Assortment

#### Vario-Stud-Snap vks universal 2.2

14 pieces

- 2 Matrices each red, yellow, green
- 2 Metal matrix housings
- 2 Blocking out discs
- 2 Patrices
- 1 Matrix inserting instrument
- 1 Paralleling mandrel

REF 430 0532 0

	Product	REF	Ø	Height
	Metal matrix housing vks-oc 1.7	430 0697 0	3.5 mm	2.3 mm
	Matrix housings vks-oc 1.7	430 0699 0	3.5 mm	2.3 mm
	Metal matrix housing vks-oc 2.2	430 0696 0	4.3 mm	3.1 mm
	Matrix housings vks-oc 2.2	430 0698 0	4.3 mm	3.1 mm
	Matrices vks-oc 1.7	430 0655 0	2.7 mm	2.0 mm
		430 0659 0	2.7 mm	2.0 mm
		430 0656 0	2.7 mm	2.0 mm
	Matrices vks-oc 2.2	430 0544 0	3.3 mm	2.7 mm
		430 0545 0	3.3 mm	2.7 mm
		430 0546 0	3.3 mm	2.7 mm
C	Blocking out disc vks-oc 1.7	430 0652 0	2.8 mm	0.4 mm
	Blocking out disc vks-oc 2.2	430 0540 0	3.5 mm	0.4 mm



Important information about processing of attachments, see catalogue page 133

### vks-oc uni



Use on root caps for fixation in the CoCr structure.





vks-oc uni Ø 1.7 mm 8 pieces REF 430 0676 0 50 pieces REF 430 0675 0

**Patrix** 

Patrix vks-oc uni Ø 2.2 mm 8 pieces REF 430 0538 0 50 pieces REF 430 0550 0



Patrix vks-oc uni Ø 1.7 mm HL-patrix cast-on 2 pieces REF 430 0701 0



Patrix vks-oc uni Ø 2.2 mm HL-patrix cast-on 2 pieces REF 430 0700 0

#### Accessories:



vks Paralleling mandrel oc/sg 1 piece ph-vks 1.7 REF 430 0677 0 ph-vks 2.2

REF 360 0113 0




Root cap is modelled in the usual way. Bring the vks-oc uni into the correct position using the paralleling mandrel and fix with hot wax.



Remove paralleling mandrel and apply rich coat of hot wax onto the transition zone of vks-oc uni/root cap. The onepiece casting facilitates processing.



Casting is carried out according to standard criteria. After casting, the vks-oc unit is only slightly polished to high lustre using a textile buff.



Fill undercuts between the blocking out disc and the marginal line with blocking out wax and block out and duplicate the chrome cobalt framework in the usual way.



In order to wax up the chrome cobalt framework over the matrix housings, use the special wax matrix housings. They ensure correct thickness of the chrome cobalt frame.

Complete the chrome cobalt framework in the usual way.

Dimensions



Product	REF	Ø	Height
Patrix vks-oc uni/HL	430 0675 0	Stud 1.7 mm	2.2 mm
	430 0550 0	Stud 2.2 mm	3.2 mm

# Vario-Kugel-Snap vks-oc

### vks-oc uni



Rigid matrices vks-oc rs Ø 2.2 mm green – reduced softsnap-in effect 8 pieces REF 440 0070 8



Rigid matrices vks-oc rs Ø 2.2 mm yellow – medium softsnap-in effect 8 pieces REF 440 0080 8



Rigid matrices vks-oc rs Ø 2.2 mm red – high soft-snap-in effect 8 pieces REF 440 0090 8



Matrix housings vks-oc rs 2.2 for glueing or laserwelding 2 pieces REF 440 0020 2



Wax matrix housing vks-oc rs 2.2 8 pieces REF 440 0100 8 50 pieces REF 440 0105 0



Blocking out discs vks-oc rs 2.2 8 pieces REF 440 0010 8



Duplicating matrix vks-oc rs 2.2 8 pieces REF 440 0110 8



Matrix pliers vks-oc Ø 2.2 mm + zg 1 piece REF 310 0000 6

#### Assortment Vario-Stud-Snap vks-oc rs 2.2

18 pieces

- 2 Rigid matrices each, red, yellow, green
- 2 Matrix housings
- 2 Duplicating matrix inkl. 2 Matrices yellow
- 2 Wax matrix housing
- 2 Blocking out discs
- 2 Patrices
- 1 Matrix inserting instrument
- 1 Paralleling mandrel

REF 440 0001 0





Matrix inserting instrument vks-oc rs Ø 2.2 mm 1 piece REF 360 0116 1




universal 2 vks-oc rs Ø 2.2 mm 1 piece REF 360 0116 0

Paralleling mandrel



DTK-adhesive REF 540 0010 6





# 10 YEARS OF THE SKY IMPLANT SYSTEM

# SUMMARISED IN THE 2012 SCIENTIFIC BOOK

# INTERNATIONAL, FUTURE-ORIENTED, SCIENTIFIC AND PRACTICAL:

The three-part, 136-page long Scientific Book contains, in its first part, a summary of the abstracts of the presentations, Excellence courses and workshops that took place during the 2012 SKY Meeting and presents, in its second part, the scientific and clinical posters. The scientific posters show the latest results of the trials involving the SKY implant family (blueSKY, SKY classic, whiteSKY and miniSKY). Furthermore, results achieved using the HELBO treatment (aPDT) and the regeneration materials are presented. The clinical posters show the prosthetic possibilities of the SKY Implant System, miniSKY and whiteSKY, in connection with bredent

materials. The conclusion in the third part provides an overview of specialist articles published nationally and internationally on the SKY Implant System and on the whiteSKY zirconium abutments.

The many authors and the bredent group are convinced that the 2012 Scientific Book will serve as a source of inspiration and invite you to gain an impression of the benefits that the many suggestions could have for you in your daily work.



# Vario-Kugel-Snap vks-oc

# vks-oc exchangeable stud



Safety, precision and biocompatibility due to easily exchangeable titanium



Stud-head screw vks-oc/sg 1.7 titanium 1 piece REF 450 0005 6



Stud-head screw vks-oc/sg 2.2 titanium 1 piece REF 450 0004 7



Thread sleeve vks-oc 1.7 1 piece REF 450 0005 4 platinum-iridium REF 450 0005 5



Thread sleeve vks-oc 2.2 1 piece ΗĹ REF 450 0004 6 platinum-iridium REF 450 0005 3

#### Assortment

vks-oc 1.7 exchangeable stud

- 5 pieces
- 1 Stud-head screw 1 Thread sleeve
- 1 Fixation screw 1 Screwdriver
- 1 Paralleling mandrel

REF 450 0005 8

#### Accessories:



vks Paralleling mandrel oc/sg 1 piece ph-vks 1.7 REF 430 0677 0 ph-vks 2.2 REF 360 0113 0



Screwdriver short, hexagon 1 piece REF 330 0069 0

#### Assortment

vks-oc 2.2 exchangeable stud 5 pieces

- 1 Stud-head screw
- 1 Thread sleeve
- 1 Fixation screw
- 1 Screwdriver 1 Paralleling mandrel

REF 450 0004 5



Fixation screw M 2 1 piece REF 450 0004 8 M 1.6 1 piece

REF 450 0005 7



Screwdriver Stud-head screw vks oc/sg 1.7 1 piece REF 330 0116 4



Tap vks exchangeable stud 1.7 1 piece REF 460 0011 7



Tap vks exchangeable stud 2.2 1 piece REF 460 0012 2

# Vario-Kugel-Snap vks-oc



Important information about processing of attachments, see catalogue page 133

# vks-oc exchangeable stud



The stud-head screw is only slightly screwed into the thread sleeve and held to the root cap waxup using the paralleling mandrel.



The attachment patrix is waxed to the wax model in the determined path of insertion.



Processing is continued using vks-oc 2.2 mm rigid matrices



The stud-head screw is turned out (anticlockwise) of the thread sleeve using the screwdriver.



Prior to investing the model, the stud-head screw must be replaced by the fixation screw.



or vks-oc rs 2.2 mm rigid



Colloid graphite is applied onto the thread area of the fixation screw; then the screw is turned into the thread sleeve exerting minimum



The casting is sandblasted and the fixation screw is turned out. The root cap is finished, the stud-head screw turned in and polished to high lustre using titanium polishing paste.








Product	REF	Ø	Thread	Height
Stud-head screw vks-oc/sg 1.7	450 0005 6	Stud 1.7 mm	M 1.6 x 0.2	2.9 mm
Stud-head screw vks-oc/sg 2.2	450 0004 7	Stud 2.2 mm	M 2 x 0.25	3.5 mm
Thread sleeve vks-oc 1.7	450 0005 4	3.4 mm	_	1.7 mm
Thread sleeve vks-oc 2.2	450 0004 6	3.4 mm	_	1.7 mm

# vks-oc/sg exchangeable stud with adhesive sleeve





**Auxiliary modelling** element 1.7 1 piece REF 450 0007 3



Auxiliary modelling element 2.2 1 piece REF 450 0007 5





Thread sleeve titanium 1.7 2 pieces REF 450 0007 4



Thread sleeve titanium 2.2 2 pieces REF 450 0007 6

#### Accessories:







The shape of the auxiliary modelling element allows to recognize the final alignment of the attachment.



DTK-adhesive REF 540 0010 6



Place the matrix on the stud-head screw and continue processing in the usual way.



Remove the auxiliary modelling element prior to investing.



Processing of vks-oc is carried out using the same auxiliary modelling element.

#### **Dimensions**

The auxiliary modelling

mandrel into the model

according to the path of

After polishing, turn the

stud-head screw into the


thread sleeve and glue in

the sandblasted seating

using DTK adhesive.

element is integrated

with the paralleling





### vks-sg

Stud attachments must be processed very precisely for functionally-secure application. It must be possible to guarantee the positional stability of the prosthesis, in order to prevent overloading caused by micro movements of the matrix and the patrix. The attachment is only designed for secure hold of the prosthesis, it is not, however, for permanent weight-bearing. In the case of extracoronal use, a shear distributor with integrated Interlock must be available. The required stability and force distribution are thus achieved.

Due to its snap function, the stud attachment offers patients more security, as a clear snapping in of the prosthesis is noticeable.

#### Indications

- Extracoronal applications
- Retaining element for bars
- Model casting fixation

#### **Product variations**

- Vario-Kugel-Snap vks-sg sv 1.7
- Vario-Kugel-Snap vks-sg 1.7
- Vario-Kugel-Snap vks-sg 2.2
- Vario-Kugel-Snap vks-sg 1.7 exchangeable stud
- Vario-Kugel-Snap vks-sg 2.2 exchangeable stud
- Vario-Kugel-Snap vks-sg 1.7 bar patrix
- Vario-Kugel-Snap vks-sg 2.2 bar patrix

#### **Properties**

- Cost-effective purchase
- Time-saving and secure processing
- Reduction of the alloy variety for high biocompatibility
- Cost-effective solution provided by castable plastic patrices

#### Materials

- Patrices
  - Combustible plastic
- Castable precious metal alloy
- Matrices
- Biocompatible thermoplast

#### Dimensions


The Vario-Kugel-Snap vks-oc stud head attachments are available in sizes 1.7 and  $^{2}$  2

### The matrices have a retention of

green 4N




yellow 6N




red 8N





vks-sg sv 1.7



vks-sg 1.7/2.2




vks-sg 1.7/2.2 exchangeable stud



vks-sg 1.7/2.2 bar patrix

# Vario-Kugel-Snap vks-sg

# vks-sg exchangeable stud



Cast-on thread sleeve and exchangeable titanium stud for precision, biocompatibility and reliability.



Thread sleeve vks-sg 1.7 1 piece HL REF 450 0005 9 platinum-iridium REF 450 0006 0



Thread sleeve vks-sg 2.2 1 piece HL REF 450 0005 1 platinum-iridium REF 450 0005 2



Stud-head screw vksoc/sg 1.7 1 piece titanium REF 450 0005 6



Stud-head screw vksoc/sg 2.2 1 piece titanium REF 450 0004 7

#### Assortment

vks-sg 1.7 exchangeable stud 5 pieces

- 1 Stud-head screw
- 1 Thread sleeve
- 1 Fixation screw 1 Screwdriver
- 1 Paralleling mandrel

REF 450 0006 1

### Assortment

vks-sg 2.2 exchangeable stud

- 5 pieces
- 1 Stud-head screw 1 Thread sleeve
- 1 Screwdriver 1 Paralleling mandrel

1 Fixation screw

REF 450 0004 9

#### Accessories:



vks Paralleling mandrel oc/sg 1 piece ph-vks 1.7 REF 430 0677 0 ph-vks 2.2 REF 360 0113 0



Screwdriver short, hexagon 1 piece REF 330 0069 0



Fixation screw M 2 1 piece REF 450 0004 8 M 1.6 1 piece REF 450 0005 7



Screwdriver Stud-head screw vks oc/sg 1.7 1 piece REF 330 0116 4



Tap vks exchangeable stud 1.7 1 piece REF 460 0011 7



Tap vks exchangeable stud 2.2 1 piece REF 460 0012 2



Important information about processing of attachments, see catalogue page 133

### vks-sg exchangeable stud



The stud-head screw is only slightly screwed into the thread sleeve and held to the wax model using the paralleling mandrel.



The attachment patrix is waxed to the wax model in the path of insertion of the shear distributor with parallel interlock.



The stud head screw is turned out (anticlock-wise) of the thread sleeve using the screwdriver.



The stud head screw is turned out (anticlock-wise) of the thread sleeve using the screwdriver.



Colloid graphite is applied onto the thread of the fixation screw; then the screw is turned into the thread sleeve exerting minimum force.



The casting is sandblasted and the fixation screw is turned out. The crowns are finished and the stud-head screw is turned in.




The stud-head screw is polished to high lustre using titanium polishing paste.



The yellow matrix is placed on the stud and the model is prepared for duplicating. Further processing with the Vario-Stud-Snap vks-sg.







# The complete product range



With integral shear distributor.



For custom applications.



For use on bars in tooth-bounded gaps when fabricating complex restorations.



Matrices green - reduced snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0668 0 Ø 2.2 mm 8 pieces REF 430 0541 0



Matrices yellow - regular snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0666 0 Ø 2.2 mm 8 pieces REF 430 0542 0



Matrices red - high snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0664 0 Ø 2.2 mm 8 pieces REF 430 0543 0



Vario-Stud-Snap vks-sg/sv 1.7 Refill package: **Patrices** 8 pieces REF 430 0735 3



Patrix sg Refill package: Ø 1.7 mm 8 pieces REF 430 0670 0 Ø 2.2 mm 8 pieces REF 430 0537 0



Patrix sg universal Refill package: Ø 1.7 mm 8 pieces REF 430 0676 0 Ø 2.2 mm 8 pieces REF 430 0538 0



vks-Patrix universal н Refill package: Ø 1.7 mm 2 pieces REF 430 0701 0 Ø 2.2 mm 2 pieces REF 430 0700 0



vks Paralleling mandrel oc/sg 1 piece ph-vks 1.7 REF 430 0677 0 ph-vks 2.2 REF 360 0113 0



Matrix inserting instrument vks-oc Ø 1.7 mm 1 piece REF 430 0621 0 vks-oc Ø 2.2 mm 1 piece REF 430 0548 0



Assortments vks assortment

vks assortment sg/2.2

Metal transfer patrices vks-oc Ø 1.7 mm 8 pieces REF 430 0662 0 vks-oc Ø 2.2 mm 8 pieces REF 430 0548 2

35 pieces **REF 430 0530 0** 

10 pieces REF 430 0533 0



Paralleling mandrel universal for vks-sq/sv REF 360 0115 1




### **Dimensions**



	в	м	١
	я		
×	ы		
	N		
	4		À

Product	REF	Depth	Width	Height
Patrix vks-sg 1.7	430 0670 0	2.7 mm	3.0 mm	4.1 mm
Patrix vks-sg 1.7 sv	430 0735 3	4.3 mm	3.5 mm	4.5 x 5.5 mm
Patrix vks-sg 2.2	430 0537 0	3.5 mm	3.8 mm	5.4 mm
Matrix vks-sg 1.7	430 0668 0	2.3 mm	3.2 mm	3.1 mm
	430 0666 0	2.3 mm	3.2 mm	3.1 mm
	430 0664 0	2.3 mm	3.2 mm	3.1 mm
Matrix vks-sg 2.2	430 0541 0	2.85 mm	4.2 mm	4.1 mm
	430 0542 0	2.85 mm	4.2 mm	4.1 mm
	430 0543 0	2.85 mm	4.2 mm	4.1 mm




vks assortment sg/uni 1.7 25 pieces REF 430 0651 0



Important information about processing of attachments, see catalogue page 133

### The complete product range



Shear distributor vks-sg/sv 1.7

sg/sv 1.7 patrix with a completely new snap-in attachment including a shear distributor with 1.7 mm stud. No additional shear distributor required.



Vario-Stud-Snap vks-sg/sv 1.7 Assortment REF 430 0735 2



Paralleling mandrel for the sg/sv 1.7 patrix. Its slender design simplifies waxing the attachment onto the crown.





After casting, the stud and friction surfaces must not be trimmed.



The high lustre surfaces are perfect for the friction snap-in matrix to fit precisely.



The model must always be duplicated with the yellow matrix in position. This ensures that the housing will be perfect for setting the snap-in friction as required.



The investment model can be cast using standard methods.




The precisely reproduced matrix is integrated into the pattern.



The outer coping, with the matrix pressed into place, guarantees long-term, gentle snap-in-friction.


If the vks attachment is to function perfectly it is essential that the patient finds the "snap-in spot" with the fingers and presses on the restoration to lock it into place.



vks-sq for free-end dentures



sg patrix
The concave waxing surface and mirrorfinish on the resin provide the best possible conditions for producing precise castings.





The concave waxing surface on the sg patrix permits it to be fitted in close proximity to the crown.



The diameter of the stud must not be modified.



No spacer wax should be applied beneath the matrix during blocking out so that the matrix can be fully enveloped in metal.



The framework pattern must cover the matrix completely.



Trim the chrome cobalt framework as usual and fit it down. Coat the matrix housing with wax when polishing the framework.




Use the inserting instrument to insert the matrix with the desired friction.



# Vario-Kugel-Snap vks-sg

# The complete product range



For custom applications



Refill package: Ø 1.7 mm 8 pieces REF 430 0676 0 Ø 2.2 mm 8 pieces REF 430 0538 0

Patrix sg universal



vks-Patrix universal HL Refill package: Ø 1.7 mm 2 pieces REF 430 0701 0 Ø 2.2 mm 2 pieces REF 430 0700 0



Matrices green - reduced snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0668 0 Ø 2.2 mm 8 pieces REF 430 0541 0



yellow - regular snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0666 0 Ø 2.2 mm

Matrices

8 pieces REF 430 0542 0



Matrices red - high snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0664 0 Ø 2.2 mm 8 pieces REF 430 0543 0



Wax bars - save time when waxing up bars/attachments

1.6 x 8 x 50 mm REF 430 0265 0 1.9 x 4 x 50 mm REF 430 0266 0 2.2 x 6 x 50 mm REF 430 0267 0



vks Paralleling mandrel oc/sg 1 piece

ph-vks 1.7 REF 430 0677 0 ph-vks 2.2 REF 360 0113 0



Matrix inserting instrument vks-oc Ø 1.7 mm 1 piece REF 430 0621 0 vks-oc Ø 2.2 mm 1 piece REF 430 0548 0



Paralleling mandrels for wax bars

Paralleling mandrel 1.6 for wstg 1.6 1 piece REF 430 0268 0

Paralleling mandrel 1.9 - 2.2 for wstg 1.9 - 2.2 1 piece REF 430 0270 0




Once the cost-effective, one-piece casting has been completed, it can be blocked out for duplicating.



The bar and matrix are simply coated with wax. The remaining sections of the pattern should be waxed up as required.



The investment model must be fabricated using a yellow matrix.



The Vario-Stud-Snap sg retains every type of denture securely. The dentist can select the friction according to the patient's needs.




Important information about processing of attachments, see catalogue page 133

# vks-sg matrix housing



Matrix housing vks-sg 1.7 8 pieces REF 430 0670 8

The matrix housing ensures reliable hold of the matrix and simultaneously allows to obtain a stress-free CoCr structure thanks to glueing.



Matrix housing vks-sg 2.2 8 pieces REF 430 0680 8



Prepare the crown in the usual way.



Insert the green matrix into the plastic matrix housing and place it on the stud patrix.



Use Pi-Ku-Plast for modelling the shear distributor and connect it with the matrix housing. Remove the matrix prior to casting.



Sandblast the appendix before glueing and attach retentions.



Use DTK adhesive to glue matrix housing and CoCr structure together.



The appendix can also be polymerized directly into the plastic saddle.
A shear distributor is always required.

### Vario-Kugel-Snap vks-sg



Important information about processing of attachments, see catalogue page 133


### vks-sq bar patrix



Vario-Stud-Snap bar patrix. Bar element with three integrated vks-studs in the sizes 1.7 or 2.2 mm. A titanium bar is also available for the exchangeable stud.

vks-sg bar patrix titanium 1.7 1 piece without stud-head screw REF 450 OSA1 7

vks-sq bar patrix titanium 2.2 1 piece without stud-head screw REF 450 OSA2 2



vks-sg bar patrix 1.7 8 pieces REF 430 0800 8

vks-sg bar patrix 2.2 8 pieces REF 430 0810 8

#### Assortment

vks-sq bar patrix 1.7 13 pieces

3 Matrices each green, yellow, red 2 Bar patrixes 1 Matrix inserting instrument 1 Paralleling mandrel

REF 430 0806 0

#### Assortment

vks-sq bar patrix 2.2 13 pieces 3 Matrices each green, yellow, red

2 Bar patrixes 1 Matrix inserting instrument 1 Paralleling mandrel

REF 430 0816 0



Stud-head screw vks-oc/sg 1.7 titanium 1 piece REF 450 0005 6



Screwdriver Stud-head screw vks-oc/sg 1.7 1 piece REF 330 0116 4



Paralleling mandrel 1.6 1 piece REF 430 0623 0



Stud-head screw vks-oc/sg 2.2 titanium 1 piece REF 450 0004 7



Screwdriver short, hexagon 1 piece REF 330 0069 0



Matrix inserting instrument vks 1.7, 1 piece REF 430 0621 0



Matrices green - reduced snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0668 0 Ø 2.2 mm 8 pieces REF 430 0541 0



Matrices yellow - regular snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0666 0 Ø 2.2 mm 8 pieces REF 430 0542 0



vks 2.2, 1 piece REF 430 0548 0 Matrices red - high

snap-in friction Refill package: Ø 1.7 mm 8 pieces REF 430 0664 0 Ø 2.2 mm 8 pieces REF 430 0543 0



Time is saved during waxing up thanks to the integrated vks studs. The bar is cut to the proper length using a separating disc and fitted into the



Use paralleling mandrel for waxing up the bar patrix to the crowns according to the direction of insertion.



The cast bar and any undercuts below the matrix are blocked out (filled) with wax.



The bar and the matrix are simply coated with wax. The remaining sections of the pattern should be waxed up as



The chrome cobalt framework is fitted down and polished to high lustre using Brepol.



Use the inserting instrument to press the corresponding matrix in the housing.



### Rod attachments Vario-Soft 3

The rod attachments from the Vario-Soft group facilitate processing due to their rounded shape. A higher level of precision is thus achieved during manufacture, which offers easier handling for patients. The rod attachments are fricative and can be adjusted individually using the various colour-coded matrices. Soft incorporation and removal is achieved by means of the special plastic of the matrices. Application of a shear distributor with milled Interlock is also necessary for correct distribution of forces. The attachments with the integrated shear distributor are the exceptions.

#### Indications

- Extracoronal applications
- Plastic fitting
- Model casting fixation

#### **Product variations**

- Vario-Soft 3 the classic!
- Vario-Soft 3 sv with integrated shear distributor
- Vario-Soft 3 mini in the case of a small amount of space
- Vario-Soft 3 mini sv
- Vario-Soft 3 matrices housing
- Vario-Soft 3 bridge-sectioning attachment
- Vario-Soft 3 mini sv zircon specially for the application of zirconium frameworks
   Vario-Soft 3 Conicalbridge the small bridge-sectioning attachment

#### **Properties**

- Cost-effective purchase
- Time-saving and secure processing
- Reduction of the alloy variety for a high level of biocompatibility
- Cost-effective solution provided by castable plastic patrices

#### Materials

- Patrices
  - Combustible plastic
- Castable precious metal alloy
- Matrices
  - Biocompatible thermoplast

#### The matrices have a retention of

green 4N




yellow 6N



red 8N





Vario-Soft 3 - the classic!




Vario-Soft 3 sv – with integrated shear distributor



Vario-Soft 3 mini - in the case of a small amount of space



Vario-Soft 3 mini sv



Vario-Soft 3 matrix housing



Vario-Soft 3 bridge-sectioning attachment



Vario-Soft 3 mini sv zircon – specially for the application of zirconium frameworks





Vario-Soft 3 Conicalbridge - the small bridgesectioning attachment



# **Rod attachments**


### Vario-Soft 3





offers numerous application possibilities using the same matrices

vs 3 Patrix 8 pieces REF 430 0520 0






vs 3 Patrix patrix withoutparalleling mandrel 8 pieces REF 430 0737 0


### Soft matrices

Matrixes that have proved their reliability for 15 years provide safety and ensure high comfort of wear for the patient.



Soft soft matrices

Special soft plastic compensates small divergences and minor processing imperfections.



Duplicating matrix 8 pieces REF 430 0737 2



green - reduced friction 8 pieces REF 430 0519 0



green - reduced friction 8 pieces REF 430 0565 0



Wax matrix housing 8 pieces REF 430 0521 0



yellow - regular friction 8 pieces REF 430 0518 0



yellow - regular friction 8 pieces







red - high friction 8 pieces REF 430 0517 0



red - high friction 8 pieces REF 430 0563 0

#### Assortment

13 pieces

#### Vario-Soft 3 2 vs 3 Patrices

- 1 Matrix inserting instrument
- 2 Duplicating matrix
- 2 Wax matrix housing
- 2 Soft Matrices, green reduced friction
- 2 Soft Matrices, yellow regular friction
- 2 Soft Matrices, red high friction

REF 430 0516 0

#### Assortment

13 pieces

#### Vario-Soft 3 Soft

- 2 vs 3 Patrices
- 1 Matrix inserting instrument
- 2 Duplicating matrix
- 2 Wax matrix housing
- 2 Soft Soft Matrices, green reduced friction
- 2 Soft Soft Matrices, yellow regular friction
- 2 Soft Soft Matrices, red high friction

REF 430 0561 0

#### Assortment

13 pieces

# Vario-Soft 3 without integrated paralleling mandrel

- 2 vs 3 Patrices without paralleling mandrel
- 1 Matrix inserting instrument
- 2 Duplicating matrix
- 2 Wax matrix housing
- 2 Soft Matrices, green reduced friction
- 2 Soft Matrices, yellow regular friction
- 2 Soft Matrices, red high friction

REF 430 0738 2

### Vario-Soft 3



After casting, patrices must only be processed using rubber polishers and high-lustre buffs.



The white duplicating matrix that has been adapted from the basal direction provides the perfect precondition for all other types of friction.



Master model prepared for duplicating.



Wax matrix housing on the investment material model guarantees a uniform chrome cobalt housing.



Completed wax pattern of the later chrome cobalt framework.



The use of the inserting instrument ensures precise positioning of the matrices.

#### Accessories:



Paralleling mandrel universal for vks-sg/sv REF 360 0115 1



Product	REF	Ø	Depth	Width	Height	Max. reduction
vs 3 Patrix	430 0737 0	1.8 mm	3.1 mm	3.0 mm	6.0/7.0 mm	3.0 mm
vs 3 Matrix	430 0519 0	-	3.6 mm	3.2 mm	7.0 mm	3.0 mm
	430 0518 0	-	3.6 mm	3.2 mm	7.0 mm	3.0 mm
	430 0517 0	_	3.6 mm	3.2 mm	7.0 mm	3.0 mm

### Vario-Soft 3 sv



with integrated shear distributor

saves time and provides perfect options for esthetic design while ensuring maximum transfer of forces





Patrix with shear distributor 8 pieces REF 430 0737 4



**Duplicating matrix** 8 pieces **REF 430 0737 2** 



Wax matrix housing 8 pieces REF 430 0521 0



Soft Matrix green - reduced friction 8 pieces REF 430 0519 0



Soft Matrix yellow - regular friction 8 pieces REF 430 0518 0



Soft Matrix red - high friction 8 pieces REF 430 0517 0



Matrix inserting instrument 2 pieces REF 430 0736 6



Accessories:

Paralleling mandrel universal for vks-sg/sv REF 360 0115 1

#### Assortment

13 pieces

Vario-Soft 3 sv

2 vs 3 Patrices with integrated shear distributor

- 1 Matrix inserting instrument
- 2 Duplicating matrix
- 2 Wax matrix housing
- 2 Soft Matrices, green reduced friction
- 2 Soft Matrices, yellow regular friction
- 2 Soft Matrices, red high friction

REF 430 0738 3

#### Dimensions



Product	REF	Ø	Depth	Width	Height	Max. reduction
vs 3 sv Patrix	430 0737 4	8 mm	5,3 mm	3.5 mm	6.0/7.0 mm	3.0 mm
vs 3 Matrix	430 0519 0	-	3.6 mm	3.2 mm	7.0 mm	3.0 mm
	430 0518 0	-	3.6 mm	3.2 mm	7.0 mm	3.0 mm
	430 0517 0	_	3.6 mm	3.2 mm	7.0 mm	3.0 mm



The patrix based on computer-aided-design includes all requirements of a modern filigree retaining element.

Precise investment

material model ensures

precision-fit integration of shear distributors.



The white duplicating matrix guarantees precision-fit integration of the various friction matrices.



Due to the integrated shear distributor patient-friendly constructions that protect the periodontium can be



Master model prepared for the production of the investment compound model.



If other friction values are desired, simply exchange the matrices.




tions that pro periodontium achieved.






## Vario-Soft 3 matrix housing



Made of plastic to produce a precision-fit metal matrix housing with any alloy.



Matrix housing 8 pieces REF 430 0737 6 Assortment
6 pieces
Vario-Soft vs 3 Matrix housing
2 Matrix housing
2 Duplicating matrix housing
2 Wax housing
REF 430 0738 4





Duplicating matrix housing 8 pieces REF 430 0737 8

Wax housing 8 pieces REF 430 0738 0




Accessories:

DTK-adhesive REF 540 0010 6





Matrix yellow soft regular friction 8 pieces REF 430 0564 0





The matrix housing is perfectly suitable for all vs 3 patrices. The shear distributor must be integrated.



Fig. 1:1

Slide vs 3 matrix into matrix housing and adapt to the prevailing conditions from the basal direction; then place it onto the patrix.



The shear distributor is coated with Pi-Ku-Plast modelling resin and connected to the matrix housing. Retention crystal must not be spread onto the retention area of the matrix housing.



Prior to investing, remove vs 3 matrix from the matrix housing and cast in the alloy of your choice.



After removing inaccuracies in the cast object, insert the matrix with the inserting instrument.



The duplicating matrix housing is placed onto the retention element. The defined wall thickness of 0.2 mm ensures an optimum gap for glueing.



The outer shape of the duplicating matrix has also been designed in a way to ensure clamping of the adhesive during glueing.



The wax housing is placed onto the retention appendix and connected to the chrom cobalt model.



Prior to glueing, vaseline is applied to the master model and the parts to be glued are sandblasted with 110  $\mu$  aluminium oxide.



A thin coat of DTK adhesive is applied onto the matrix housing and the chrome cobalt framework.



which are then glued exerting uniform pressure.







Product	REF	Width	Height	Max. reduction
vs 3 Matrix housing	430 0737 6	1.8 / 4.7 mm	5.0 / 7.6 mm	individual

# Vario-Soft 3 sv bridge-sectioning attachment

### One attachment ...



Bridge-sectioning attachment with integrated shear distributor suitable for casting-on.



Patrix HL suitable for casting-on 2 pieces REF 430 0730 4

Titanium screw

REF 330 0070 0

REF 330 0071 0

1 piece

10 pieces



Patrix HL cast-on 1 piece REF 450 0000 1

#### Accessories:



Paralleling mandrel universal 1 piece REF 360 0115 1



Screwdriver short 1 piece REF 330 0069 0

#### Assortment

4 pieces, 1 piece each Patrix HL cast-on Titanium screw Patrix HL suitable for casting-on Screwdriver short

REF 450 0000 2



The bridge sectioning attachment that can be cast on is positioned at the wax pattern using the parelleling mandre.



After casting, the crown framework is checked and finished.



After ceramic veneering, the bridge-sectioning attachment is polished with high-lustre buffs.



Fix the cast-on closing ring with titanium screw and cover with Pi-Ku-Plast.



The bridge is waxed up in the usual way.



Completed and fitted bridge framework. Ready for ceramic veneering.

Dimensions



Product	REF	Ø	Depth	Width	Height	Max. reduction
Patrix	450 0000 1	1.8 mm	6.1 mm	3.0 mm	7.6/7.0 mm	2.8 mm
Titanium screw 1.4	330 0070 0	M1.4 x 0.3	_	2.1 mm	4.5 mm	1.2 mm
Patrix HL suitable for casting-on	430 0730 4	2.5 mm	_	_	2.1 mm	1.4 mm

# Vario-Soft 3 sv bridge-sectioning attachment

### ... two indications



In case of loss of the terminal abutment of the bridge the previous patrix becomes the fixation base for the new removable attachment denture.




Fig. 1:1

Matrix red - high friction 8 pieces REF 430 0517 0

Matrix

8 pieces REF 430 0519 0

Matrix

8 pieces REF 430 0518 0



**Duplicating matrix** 8 pieces REF 430 0737 2



Wax matrix housing 8 pieces REF 430 0521 0





Matrix adhesive assortment REF 540 0103 1

If the plastic attachment matrix is not sufficiently retained in the chrome cobalt framework, this tested and approved adhesive system should be used.



After taking the impression and producing the model, the patrix is slid onto the previous bridgesectioning attachment and duplicated.



Wax pattern of the later chrome cobalt framework on the investment material model.



After casting, press in the desired friction matrix.



Completed chrome cobalt framework with attachment to prepare the set-up of teeth.

### Vario-Soft 3 mini



The computerized slender design and three patient-specific soft friction types provide reliable retention even in cases of limited space available.



**Duplicating matrix** 

REF 430 0732 5

**Patrix** 

8 pieces



white 8 pieces REF 430 0732 3



Wax matrix housing REF 430 0732 0



Matrix green - reduced friction 8 pieces REF 430 0731 7



Matrix yellow - regular friction 8 pieces REF 430 0731 5



Matrix red - high friction 8 pieces REF 430 0731 3

#### Assortment

13 pieces

Vario-Soft 3 mini

- 2 Patrices
- 2 Duplicating matrices
- 2 Wax matrix housings
- 2 Matrices, green reduced friction
- 2 Matrices, yellow regular friction
- 2 Matrices, red high friction
- 1 Matrix inserting instrument

REF 430 0731 2



Fig. 1:1

Matrix inserting instrument 2 pieces REF 430 0736 5



#### Accessories:



Paralleling mandrel universal 1 piece REF 360 0115 1



The slender design of the paralleling mandrel ensures safe retention and leaves sufficient space for waxing up.



The duplicating matrix ensures precise fabrication of the metal matrix housing in the chrome cobalt framework.




The working steps are carried out in the usual way. This way quality is assured.

4	Product	REF	Depth	Width	Height	Max. reduction
	vs 3 mini Patrix	430 0732 5	2.3 mm	3.1 mm	6.0 mm	3.0 mm
	vs 3 mini Matrix	430 0731 7	2.0 mm	3.0 mm	6.0 mm	3.0 mm
iii		430 0731 5	2.0 mm	3.0 mm	6.0 mm	3.0 mm
		430 0731 3	2.0 mm	3.0 mm	6.0 mm	3.0 mm
100						




### Vario-Soft 3 mini sv






Patrix 8 pieces REF 430 0734 3



Matrix green - reduced friction 8 pieces REF 430 0733 5



Duplicating matrix white 8 pieces REF 430 0734 1



Matrix yellow - regular friction 8 pieces REF 430 0733 3



Wax matrix housing 8 pieces REF 430 0733 8



Matrix red - high friction 8 pieces REF 430 0733 1



Matrix inserting instrument 2 pieces REF 430 0736 4



Accessories:



Paralleling mandrel universal 1 piece REF 360 0115 1



The pattern is waxed up according to standard criteria; no new techniques have to be learned.



2 Wax matrix housings

REF 430 0733 0

2 Matrices, green - reduced friction


2 Matrices, yellow - regular friction 2 Matrices, red - high friction

1 Matrix inserting instrument

The optimized combustion behaviour of the patrix guarantees the precision in the cast object.



The duplicating matrix can be individually adapted to any situation.



### Vario-Soft 3 zircon sv and zircon sv mini



Regardless whether the model cast is made of Bio HPP or Brealloy, the new idea of creating removable prostheses in connection with zirconium is compelling.

# Special attachments for zirconium constructions in connection with removable prostheses

- two versions: for lateral and anterior area
- · reinforced shear distributor for highest stability
- surface-enlarged connection with retentions for a better hold

An establishedmatrix system-


to use for reliable hold within

formore than 15 years. Convenient

can be individually cut to the desired size

prostheses

- easily tactile due to perfect geometry
- integrable into any system



Vario-Soft 3 zircon sv 13-parts assortment REF 430 0732 8

Refill packages 8 male parts REF 430 0732 2



Vario-Soft 3 zircon sv mini 13-parts assortment REF 430 0732 7

Refill packages 8 male parts REF 430 0732 9



The primary construction is modeled with compoForm UV and prepared for copy milling. The attachment is arranged with the parallel holder.



The attachment is affixed with compoForm UV and immediately cured with bre.Lux N. In addition, the transition attachment/ modellation is wetted with SERACOLL UV and also cured in order to obtain a secure connection.



The modellation is affixed onto a five motions model frame. The perfectly dimensioned shape of the attach-ment simplifies the milling process.



After the milling process, the attachment is cleanly milled in the proportionate magnification ration. The zirconium construction can be prepared for sintering.



After working out and fitting the attachment, the construction is polished to a high gloss with Zi-polish. The model is duplicated as usually.

#### Accessories:



Parallel holder universal 1 piece REF 360 0115 1



compoForm UV 2 x 3 ml syringes 10 application canullae REF 540 0115 0



SERACOLL UV Light curing wax adhesive 2 x 3 ml 2 serving trays REF 540 0115 1

# Vario-Soft 3 conicalbridge



A precisely fitting sectioned bridge is made in a time-saving, economic, tension-free way without individual milling work and independent of alloy situation and.

Bridge sectioning attachment for fixed prosthesis in case of divergent abutment teeth

- · Precisely fitting, full burn-out synthetic mold
- · Conic shape for easy processing
- Integrated parallel holder on male parts and
- matrix save time and expand the application
- Designed for intra- and extra-oral use
- No individual milling work necessary
- Primary and secondary parts are fabricated
- simultaneously to save time and money



Vario-Soft 3 conical bridge 4 females, 4 males REF 430 0734 0



Female and male parts are simply exchanged for intraor extra-coronal use - any application is possible!



The modellation is made according to insertion direction and esthetic requirements.



Female and male parts are assembled and the parallel holder is taken off the not required part.



The existing retention and the height of the attachment are individually adjusted to the situation with a metal



The secondary part is completed by of the modellation of the bridge link. By individualizing the attachment, it adapts to any situation perfectly.



Simply use a finger or an instrument to remove the parallel holder at the "predetermined breaking point".



Primary and secondary part are mo-delled in one step - time saving, material saving, efficient.



The modellation is moun-ted and invested according to the bredent Casting Technique in one step. With Transfuser and Brevest Rapid 1, the slender modellation is invested easily and bubble-free.



Following casting, the attachment is blastpolished with 50 µm pearls. The attachment is now assembled without having to work it over ela-borately.



The attachment distinguishes itself through a special shape and precise fit. The longlasting connection is proof for the success!



The integrated parallel holder allows intracoronal use of the female within the primary part.

#### Accessoires:



Fissure Designer **REF B153 NF 04** 





# IMMEDIATE RESTORATIONS WITH A REDUCED NUMBER OF IMPLANTS

# CONCEPTUAL BACKGROUND AND CLINICAL RESULTS

Over 5,000 happy patients to date have proven that SKY fast & fixed is functioning in a reliable manner and the quality of life of the patients, who are on the verge of becoming edentulous, has been significantly improved by fixed bridges.

The initiators of the ideas, Dr. Georg Bayer, Dr. Frank Kistler, Dr. Steffen Kistler, DT Stephan Adler and PD Dr. Jörg Neugebauer are now able to present the results of their 5 years of experience with SKY fast & fixed. The authors explain how these innovative procedures can also be successfully introduced and used in your practice. In addition, they demonstrate the possibility of temporary and fixed restorations in which the SKY fast & fixed treatment concept can be successfully used in patients, as well as describing the commercial opportunities that this also offers for every practice.



**REF 9929710D (German)** REF 992971GB (English)

# Retaining elements

The special retaining elements provide individual use and enable a holding strength that is exactly adapted to the patient. The retaining elements can be used with all removable prostheses and thus offer a high degree of flexibility.

#### Indications

- Extracoronal applications
- Snap and friction

#### **Product variations**

- Inverto Plus the classic metal attachment
- Activatable friction cylinder individual and biocompatible
- Vario Compress 1 and 2 for individual friction adjustment
- Stud fixator for quick repair after a loss of friction
- Cylinder attachment friction and snap matrices for the individual solution

#### **Properties**

- Cost-effective purchase
- Time-saving and secure processing
- Reduction of the alloy variety for high biocompatibility

#### Materials

Biocompatible precious metal alloys, silicones and thermoplastic plastics are used depending on the type of retaining element.

The retaining elements have individually-adjustable friction.



Inverto Plus - the classic metal attachment




Activatable friction cylinder – individual and biocompatible



Vario Compress 1 and 2 for individual friction adjustment



Stud fixator - for quick repair after a loss of friction



Cylinder attachment – friction and snap matrices for the individual solution

### Inverto Plus



Exchangeable, intracoronal attachment with activating screw.



Patrix 45° incl. activating screw, glue-in sleeve and basal screw 1 piece REF 450 00P4 5

Matrix

1 piece

HL suitable for

REF 450 0004 0

casting-on



Basal screw for 45° and 90° 1 piece REF 450 0004 4







Matrix resin 2 pieces REF 450 0004 1



Patrix 90° incl. activating screw, glue-in sleeve and basal screw 1 piece REF 450 00P9 0



Activating screw for 45° type 1 piece REF 450 00A4 5



Ceramic spacer 1 piece REF 450 0004 3



Activating screw for 90° type 1 piece REF 450 00A9 0



Accessories:



Attachment is attached to the wax model with HL or plastic matrix.

Paralleling mandrel

REF 360 0116 0

universal 2

1 piece



Finish the casting and adjust height of attach-



Replace auxiliary duplicating element with glue-in sleeve and block out undercuts. Duplicate in the usual way.

Glueing in the attachment



Produce CoCr structure and glue the glue-in sleeve onto the model.



Clean glueing areas and remove excess material after the adhesive has hardened.

**Dimensions** 



Product	REF	Ø	Depth	Length	Width	Height
Matrix resin/HL	450 0004 0	_	1.55 mm	_	2,4 mm	5.4 mm
Patrix 45°	450 00P4 5	-	5.1 mm	_	2.5 mm	5.0 x 3.1 mm
Patrix 90°	450 00P90	-	5.1 mm	-	2.5 mm	5.0 x 3.1 mm
Glue-in sleeve	450 0005 0	2.5 mm	_	-	3.1 mm	-
Basal screw	450 0004 4	2.0 mm	_	0.8 mm	-	-
Activating screw 45°	450 00A4 5	1,0 mm	_	1.7 mm	-	-
Activating screw 90°	450 00A9 0	1,0 mm	-	4.0 mm	-	-
Auxiliary duplicating element	450 0004 2	2.9 mm	_	_	_	3.2 mm

# Activatable frictions cylinder



Individually adjustable, biocompatible plastic cylinder with titanium screw.

due to the special shape of the friction cylinder.



#### Assortment

- 4 pieces
- 2 Friction cylinders 2 Titanium screws


REF 440 0068 0



#### Assortment

- 20 pieces 10 Friction cylinders
- 10 Titanium screws

REF 440 0068 1



- Easy integration and safe hold of the denture
- individually adustable friction
- safe hold in the CoCr structure due to the retention stud



Friction cylinders Ceramic spacer REF 440 0068 3



Activatable use of the friction cylinder is possible with the attachment of the VS 3 group or with telescopic crowns.



Always use a shear distributor for attachments.



The plane surface of the friction cylinder is attached to the patrix.



Prior to duplicating, blocking out with wax to the basal direction is carried out and the model is prepared in the usual way.



The precise reproduction of the friction cylinder ensures accurate fit in the CoCr structure.



Prepare the model for investing in the usual



The friction cylinder is pressed into the CoCr structure using a blunt object. The screw seat is facing the basal area.



Due to the adjustment of the titanium screw, the hold of the denture can be individually adjusted to the respective patient.



Perfectly suitable for telescopic crowns



Product	REF	Thread	Depth	Length	Width	Height
Friction cylinders	440 0068 0	_	2.4 mm	-	2.4 mm	3.2 mm
Titanium screw		M 1.4 x 0.3	_	_	_	2.6 mm

### **Retention elements**

### Vario Compress 1

VC 1: Individually adjustable friction with an adjusting screw.

VC 1 adjusting screw and VC 1 friction silicon

individually

guaranteed

- adjusting screw compresses the silicone

due to the individual compression of the friction

silicon the friction of the telescopic anchor is friction silicon features a hollow space inside, under pressure the wall of the friction silicon springs quickly into this hollow space this way soft friction and smooth integration are

#### VC 1 adjusting screw

- precision thread
- titanium grade 5
- can be shortened individually

#### Approximal aperture of the VC 1 friction silicon bearing

- the retention force of the friction silicone is transferred to the primary element through this aperture



#### Thread turn and bearing for compressable VC1 friction silicon

- initial mould is produced with a ceramic mould in the casting procedure
- special tools for finishing the thread turn and bearing of the friction silicon ensure precise guidance of the components

#### Rod attachment

- any type of rod attachment can be used
- VC 1 can also be used at telescopic crowns

Vario Compress VC 1: A cylinder of abrasion-resistant special silicon is individually compressed with an adjusting screw. Due to the compression, this friction silicon exerts gentle pressure to the wall of the telescopic anchor that can be adjusted perfectly. Accordingly, individual adjustment of the static friction of the telescopic anchor is possible. The thread turn for the adjusting screw and the bearing of the friction silicon are shaped with a heat-resistant ceramic mould in the casting procedure. After casting, special tools ensure the precision of the thread turn and the bearing. Perfectly suitable for chrome cobalt supply - safe, economical and precise.



VC 1 Ceramic screw Ø 2 mm length 9,5 mm

1 piece REF 460 0010 3



VC 1 Friction silicon Ø 1.7 mm length 4 mm

1 piece REF 460 0010 4



VC 1 Titanium adjusting screw



Ceramic removing tool 1 piece REF 460 0010 6



Second tap, tungsten carbide

1 piece REF 460 0010 M



Ø 2 mm length 4 mm 1 piece REF 460 0010 5



Last tap, tungsten carbide 1 piece REF 460 0010 F



Tap handwheel REF 330 0115 3



Screwdriver, short 1 piece REF 330 0069 0

Assortment

3 pieces, 1 piece each Vario Compress 1

VC 1 Ceramic screw

VC 1 Friction silicon

VC 1 Titanium adjusting screw

REF 460 0010 7

#### Assortment

8 pieces, 1 piece each Vario Compress 1

VC 1 Ceramic screw

VC 1 Friction silicon

VC 1 Titanium adjusting screw

Ceramic removing tool


Second tap, tungsten carbide

6 Attachments, Locks and Screw connections

Last tap, tungsten carbide

Tap handwheel Screwdriver, short

REF 460 0010 1





# Vario Compress 1



also be used on telescopic crowns. Produce primary elements in the usual way. Secondary elements are waxed up directly or produced in the chrome cobalt technique (see figure 2). After high luster polishing, remove the ceramic screw, recut the thread, position the VC 1 friction silicon and adjust the friction with the adjust-

Vario Compress 1 can



Position ceramic screw so that contact with the primary element is ensured and fix with wax. Then complete the wax pattern of the secondary construction (see figure 3).



Invest and cast as usual; after casting, the secondary construction must be fitted onto the primary element and polished to a high luster.

#### VC 1: Individually adjustable friction for all telescopic attachments.

ing screw.



Prepare wax pattern of the anchor crowns in the usual way. Attach the selected type of extracoronal rod attachment to the wax moulds. Vario Compress 1 can also be used with telescopic crowns.



Cast, polish and finish crowns in the usual way. Parallel surfaces of the primary elements must be milled according to standard techniques. Prepare primary elements to produce the secondary elements.



Block out the model to produce a chrome cobalt duplicate. Prepare duplicating mould and chrome cobalt model. Then model the chrome cobalt structure according to the dental technical rules.



The special VC 1 spacer is fixed with wax in the correct position on the chrome cobalt model. Contact with the friction surface of the primary element must be ensured. The VC 1 ceramic spacer provides precise reproduction of the shape of the screw and the cylindrical VC 1 friction silicon.



Then complete the wax model of the secondary structure. The VC 1 ceramic spacer projects from the wax model. This way safe retention in the investment material of the casting ring is ensured.



After casting, complete the chrome cobalt structure on the primary construction. After polishing the chrome cobalt structure, remove the ceramic spacer by turning the ceramic removing tool gently.



Prethread the thread with the tungsten carbide first tap. The second tap removes small ceramic residues in the thread turn. Then recut the thread with the tungsten carbide last tap. Use sufficient quantities of milling and drilling oil when tapping.



Insert the VC 1 friction silicon into the clean drill hole. The VC 1 features a hollow space inside. Under pressure, the wall of the friction silicon springs into this hollow space so that soft friction is ensured.



Lock the thread with the adjusting screw and fasten the screw slightly. Due to the pressure of the adjusting screw the VC 1 friction silicon is compressed. This way the friction of the attachment is adjusted.


# Vario Compress 2


Friction silicon allows to adjust individual friction for attachments and telescopic crowns.

> Abrasion-resistant special silicon provides extended durability and safe hold of the denture.

The friction silicon is compressed by turning the VC 2 titanium adjusting screw. The friction is adjusted individually. Vario Compress 2 can be integrated from the basal or oral direction.

The primary element can be designed individually. Vario Compress 2 can be used for various indications.





VC 2 Friction silicon Ø 1.9 mm Length 6 mm 1 piece REF 460 0011 5



VC 2 Titanium adjusting screw Ø 2 mm Length 5 mm 1 piece REF 460 0011 4



VC 2 Ceramic spacer Ø 2 mm Length 5 mm 1 piece REF 460 0011 3

#### Assortment

8 pieces, 1 piece each Vario Compress 2 VC 2 Friction silicon VC 2 Adjusting screw titanium VC 2 Ceramic spacer Ceramic removing tool Second tap, tungsten carbide

Last tap, tungsten carbide Tap wheel Screwdriver, short

REF 460 0011 0

Assortment

3 pieces, 1 piece each Vario Compress 2 VC 2 Friction silicon VC 2 Titanium adjusting screw VC 2 Ceramic spacer


REF 460 0011 2



The model is prepared for duplicating.



The position of the ceramic spacer is marked on the investment model using a pen.



A small cavity at the crown and a hole in the basal area are drilled with the Rapidy 2.1 mm until the correct position of the ceramic spacer on the model is ensured.



Complete the model so that only the ceramic spacer needs to be inserted.



Integrate the ceramic spacer completely into the wax model and invest.

The titanium adjusting

screw is turned in after

completion and shortened adequately.



The ceramic removing tool and the taps, see Vario Compress 1.

Dimensions

The ceramic is removed from the thread with the ceramic removing tool. Residual ceramic particles are sandblasted with 50 my glass beads.



Prepare a chamfer at the thread opening using the Rapidy 2.1 mm and recut the thread using the taps.



Product	REF	Ø/Thread	Length	Max. reduction
Titanium adjusting screw	460 0011 4	M 2 x 0.4	5 mm	2.5 mm
Friction silicon	460 0011 5	1.9 mm	6 mm	individual



### Stud fixator



# As a snap element or to increase the friction for new restorations and repairs.

Ceramic stud and cavity-filling silicone as buffers ensure durability and soft integration of the restoration.

- Time is saved thanks to quick and easy integration
- Friction is restored subsequently
- Ceramic stud for prolonged comfort of wearing
- Hygiene-friendly thanks to cavity-filling silicone



Stud fixator 2 pieces REF 440 0265 1

#### Accessories:



DTK adhesive REF 540 0010 6

#### Procedure in the laboratory



To reproduce the oral situation accurately, use Pi-Ku-Plast to fabricate the primary construction



and to produce a working model.



Prepare a matrix before removing the resin saddle.



Drill a hole with a diameter of 2.1 mm into the secondary element and place it back on the model.



Use the drill (Ø 2.1 mm) to carefully prepare a groove with a max. depth of 0.4 mm in the resin saddle.



The stud fixator is fitted in the CoCr structure and fixed with DTK adhesive.



The housing of the stud fixator must be flush with the crown wall. Only the ceramic stud may stand out in the crown. Reattach the resin saddles.



Prepare a coping of the resin die using a thermoforming foil.



Mark the groove on the resin die with a pen.
Drill a hole (Ø 2.1 mm) through the die coping at this point.

#### Procedure in the practice



Place the die coping onto the primary construction in the mouth and transfer the position of the groove accurately.




Integrate the restoration with friction being restored.

Product	REF	Ø	Length
Stud fixator	440 0265 1	2.2 mm	3.7 mm

### **Retention elements**

## Cylindrical attachment zg

Universally suitable attachment featuring snap or friction



#### Friction


Resin matrices in three different

colors, with different pull-off force and easy exchangeability allow fast and specific adjustment of the total pull-off force of the removeable restoration.

Friction matrices or snap matrices are available in 3 different colors each and different retention levels.

The green matrix for reduced, the yellow matrix for normal and the red matrix for strong retention.

Whether the restoration is held by friction or snap can be decided individually and changed by simply exchanging the resin friction matrices or resin snap matrices.



#### Snap

For receiving the resin matrices and for stress-free, low-cost integration into the denture, the titanium matrix housing K for direct fixation in acrylics or the titanium matrix housing M for fixation (glueing) in the chrome cobalt framework are used.

#### Please select:

#### 1. Resin patrix or metal patrix



Resin patrix 8 pieces REF 440 0120 8 50 pieces REF 440 0125 0





Patrix cast-on 2 pieces REF 440 0120 2



# 2. Metal matrix housing for the integration in resin



Titanium matrix housing K
2 pieces
REF 440 0230 2
8 pieces
REF 440 0230 8
50 pieces
REF 440 0235 0



Duplicating matrix 2 pieces REF 440 0250 2





for the integration in metal

Titanium matrix housing M 2 pieces REF 440 0240 2 8 pieces REF 440 0240 8

REF 440 0240 8 50 pieces REF 440 0245 0



Wax matrix housing 8 pieces REF 440 0260 8 50 pieces REF 440 0265 0



# Cylindrical attachment zg



Friction and snap matrices can be exchanged among each other.

# 3. Matrixes with resilience, friction or snap Friction



Matrix green reduced friction 8 pieces REF 440 0150 8 50 pieces REF 440 0155 0



Matrix yellow normal friction 8 pieces REF 440 0140 8 50 pieces REF 440 0145 0



Matrix red high friction 8 pieces REF 440 0130 8 50 pieces REF 440 0135 0

#### Snap



Matrix green reduced friction 8 pieces REF 440 0180 8 50 pieces REF 440 0185 0



Matrix yellow normal friction 8 pieces REF 440 0170 8 50 pieces REF 440 0175 0



Matrix red high friction 8 pieces REF 440 0160 8 50 pieces REF 440 0165 0

#### Accessories:



Paralleling mandrel universal 2 1 piece REF 360 0116 0



Insertion pin 1 piece REF 360 0116 4



Impression transfer set Transfer patrix 2 pieces Transfer matrix 2 pieces REF 440 0116 3



Matrix pliers
1 piece
REF 310 0000 6



DTK adhesive REF 540 0010 6

### Dimensions



Product	REF	Ø	Height	Ø Wax-on area
Resin patrix	440 0120 8	2.5 mm	3.8 mm	4.6 mm
Patrix HL	440 0120 2	2.5 mm	3.7 mm	4.3 mm
Metal matrix housing K	440 0230 2	4.8 mm	4.2 mm	_
Titanium matrix housing M	440 0240 2	4.3 mm	4.2 mm	_
Matrices Friction / Snap	440 0150 8	3.75 mm	3.8 mm	_
	440 0140 8	3.75 mm	3.8 mm	-
	440 0130 8	3.75 mm	3.8 mm	_
	440 0180 8	3.75 mm	3.8 mm	-
	440 0170 8	3.75 mm	3.8 mm	_
	440 0160 8	3.75 mm	3.8 mm	_

## Cylindrical attachment zg

#### Cylindrical attachment and integration into full denture





A base for the wax pattern in which the metal matrix housings are fixed is produced using tray material.



During completion, blocking out with liquid silicone below the metal matrix housing is required to prevent resin from reaching into the matrix during pressing.



The cast-on patrix or the resin patrix is waxed on using the parallel holder.



2.1 Matrix with resilience buffer.
2.2 The matrix is pressed in the metal matrix housing using the insertion pin.
2.3 The functional principle of the resilience matrix.



While processing the resilience matrix the spacer disc must be placed under the metal matrix housing.



The cylindrical attachment provides a simple and economic way of processing that can be used in many application fields.



If required, the resin matrix can be removed with the matrix pliers and replaced by a new matrix with different friction.

#### Basic assortment

12 pieces

#### for the integration in resin \*, Friction

- 2 Resin patrices
- 2 Matrices Friction, green, reduced friction
- 2 Matrices Friction, yellow, normal friction
- 2 Matrices Friction, red, high friction
- 2 Titanium matrix housing K
- 1 Paralleling mandrel universal 2
- 1 Insertion pin

REF 440 0115 5

#### Basic assortment

12 pieces

#### for the integration in resin \*, Snap

- 2 Resin patrices
- 2 Matrices Snap, green, reduced friction
- 2 Matrices Snap, yellow, normal friction
- 2 Matrices Snap, red, high friction
- 2 Titanium matrix housing K
- 1 Paralleling mandrel universal 2
- 1 Insertion pin

REF 440 0115 4

#### Refill package

10 pieces

#### assorted \*, Friction

- 2 Resin patrices
- 2 Matrices Friction, green, reduced friction
- 2 Matrices Friction, yellow, normal friction
- 2 Matrices Friction, red, high friction
- 2 Titanium matrix housing K

REF 440 0115 8

#### Refill package

10 pieces

assorted \*, Snap

- 2 Resin patrices
- 2 Matrices Snap, green, reduced friction 2 Matrices Snap, yellow, normal friction
- 2 Matrices Snap, red, high friction
- 2 Titanium matrix housing K

REF 440 0115 7

\*The cast-on patrices (2 pieces) are available separately REF 440 0120 2.



## Cylindrical attachment zg

# Cylindrical attachment with friction matrices zg-f The ideal combination with other parallel retention elements.





After the try-in of the telescopic crowns the patrix is waxed onto the root cap parallel to the direction of insertion of the other supporting elements



In order not to change the shape of the cylindrical patrix, polishing to high luster should only be carried out using a cotton buff.

#### Basic assortment

16 pieces

for the integration in metal \*, Friction

- 2 Resin patrices
- 2 Matrices Friction, green, reduced friction
- 2 Matrices Friction, yellow, normal friction
- 2 Matrices Friction, red, high friction
- 2 Duplicating matrices
- 2 Wax matrix housings
- 2 Titanium matrix housing M
- 1 Paralleling mandrel universal 2
- 1 Insertion pin

REF 440 0115 2







Basic assortments for the integration in resin, see page 180.



The external form of the metal matrix housing ensures safe retention in the resin.



Individual adjustment of friction even after several years provide maximum comfort of wear.

### Refill package

- 14 pieces, assorted\*, Friction
- 2 Resin patrices
- 2 Matrices Friction, green, reduced friction
- 2 Matrices Friction, yellow, normal friction
- 2 Matrices Friction, red, high friction
- 2 Duplicating matrices
- 2 Wax matrix housings
- 2 Titanium matrix housing M

REF 440 0116 1

#### Cylindrical attachment with snap matrices zg-s Stress-free glueing of the metal matrix housing in the CoCr denture is possible.




Friction and snap matrices can be exchanged

among each other.



The cylindrical attachment provides a wide indication range. After casting, the duplicating matrix (6.1) is placed onto the patrix and




blocked out to the basal direction. Then the wax matrix housing (6.2) is placed on the investment material model and the CoCr structure is modelled.





Due to the small size of the cylindrical attachment an esthetic pattern can be prepared.



The metal matrix housing is glued into the CoCr structure using DTK-adhesive.

#### Basic assortment

16 pieces

for the integration in metal \*, Snap

- 2 Resin patrices
- 2 Matrices Snap, green, reduced friction
- 2 Matrices Snap, yellow, normal friction
- 2 Matrices Snap, red, high friction
- 2 Duplicating matrices
- 2 Wax matrix housings
- 2 Titanium matrix housing M
- 1 Paralleling mandrel universal 2
- 1 Insertion pin

REF 440 0115 1

## Refill package

14 pieces, assorted\*, Snap

- 2 Resin patrices
- 2 Matrices Snap, green, reduced friction
- 2 Matrices Snap, yellow, normal friction
- 2 Matrices Snap, red, high friction
- 2 Duplicating matrices
- 2 Wax matrix housings
- 2 Titanium matrix housing M

REF 440 0116 0



<sup>\*</sup> The cast-on patrices (2 pieces) are available separately REF 440 0120 2.

# www.bredent.com THE WEB ADDRESS FOR DENTAL TECHNICIANS



Current information concerning products, systems and treatment concepts is available on the home page at any time. Extensive news is displayed in an overview on the home page, along with details about the latest innovations in the field of dental technology. Detailed brochures are also available for download. These contain more extensive information on the introductory features about the products.

The event overview guides you through exhibitions, congresses and internal and external professional development events.

Further information from the field of dental technology can be found at www.bredent-medical.com.

The advantage for your knowledge!

## Bar connections

The classic solution for implant works. The stable connection between both abutments provides secure grip for the prosthesis. The various bar profiles and the matrices which accompany them allow the prosthesis to be structured in a way that is suitable for the patient.

#### Indications

- Implant prosthetics
- Root cap solutions
- Definitive bar stub

#### **Product variations**

- Vario-Soft-profile bar for snap, friction and friction-snap
- Vario-Soft-Soft the conical bar with friction matrices
- Wax bars varied choice for every type of use

#### **Properties**

- Cost-effective purchase
- Time-saving and secure processing
- Reduction of the alloy variety for high biocompatibility
- Cost-effective solution provided by castable plastic patrices

#### Materials

- Patrices
  - Combustible plastic or wax
  - Biocompatible titanium
- Matrices
- Biocompatible thermoplast

#### The matrices have a retention of

green 4N yellow 6N red 8N











Vario-Soft-profile bar for snap, friction and friction-snap



Vario-Soft-Soft – the conical bar with friction matrices



Wax bars – varied choice for every use



## Vario-Soft-Bar-Pattern vsp

The undersides of the resin and titanium bars are rounded to facilitate oral hygiene for the patient.



The snap-in retention provides for additional grip in the matrix housing. The well proven matrices are colour coded to enable the dentist to determine the degree of friction currently in use and how it can be changed.




The exterior shape of the matrices are all exactly the same, so that the matrix can be replaced with one providing a different degree of friction.

Titanium bars and high-tech Duroplast matrices, which have been tested for biocompatibility, guarantee the highest possible oral compatibility.

This type of bar has multiple indications if used as an extra-coronal bar stub.








Vario-Soft-Bar-Patterns feature adjustable friction and snap-in effect, matrices with minimal dimensions which can be placed as required, and the possibility to be reduced from beneath.

The reliable and cost-effective bar system with 3 precision matrices each with different degrees of friction for all bar indications.

Matrices for parallel bar restorations



Bar patterns made of non-distorting, fully combustible special high-tech Thermoplast, guarantee optimum castings.



Resin bar vsp-f 4 pieces REF 430 0647 0 25 pieces REF 430 0646 0



The classic parallel bar can be used for a wide range of indications.

Friction matrices vsp-f

8 pieces 50 pieces green 430 0639 0 430 0638 0 yellow 430 0641 0 430 0640 0 red 430 0643 0 430 0642 0 Duplicating matrixes 8 pieces REF 430 0625 1 50 pieces

50 pieces REF 430 0624 1

Accessories:

#### Assortment

20 pieces

Vario-Soft-Bar-Pattern vsp-f, Friction

- 4 Matrices vsp-f each, red, yellow, green
- 2 Bars vsp-f
- 4 Duplicating matrices vsp-f
- 1 Paralleling mandrel metal vsp-f/fs/gs
- 1 Insertion pin vsp-f/fs/gs

REF 430 0650 0



Titanium bar vsp-f REF 560 0001 0



Insertion pin 2 pieces REF 430 0622 0



vsp-f 8 pieces REF 430 0640 8 50 pieces REF 430 0645 0

Matrix housing



Paralleling mandrel 1 piece REF 430 0623 0



# Vario-Soft-Bar-Pattern vsp

#### Matrices for snap-in bar restorations



Friction snap-in matrices vsp-fs

8 pieces 50 pieces 430 0633 0 green 430 0632 0 yellow 430 0635 0 430 0634 0 430 0637 0 430 0636 0

#### Assortment

18 pieces

red

Vario-Soft-Bar-Pattern vsp-fs, Friction-Snap

- 4 Matrices vsp-fs each, red, yellow, green
- 2 Bars vsp-fs
- 1 Paralleling mandrel metal vsp-f/fs/gs
- 1 Insertion pin vsp-f/fs/gs

REF 430 0649 0



Accessories:

REF 430 0694 0 25 pieces REF 430 0695 0

Resin bar

vsp-fs

4 pieces



Implant in the lower jaw with a medium friction snap-in bar.



Titanium bar vsp-fs / gs REF 560 0002 0




Paralleling mandrel 1 piece REF 430 0623 0



Insertion pin 1 piece REF 430 0622 0

#### Matrices for jointed restorations





8 pieces 50 pieces

430 0627 0 430 0626 0 green yellow 430 0629 0 430 0628 0 red 430 0631 0 430 0630 0

**Duplicating matrixes** 

8 pieces REF 430 0625 0

50 pieces

REF 430 0624 0



Resin bar vsp-fs 4 pieces REF 430 0694 0 25 pieces REF 430 0695 0



The special, small, replaceable snap-in jointed matrices result in optimum bar joint restorations.

## Assortment

20 pieces

Vario-Soft-Bar-Pattern vsp-gs, joint snap-in

- 4 Matrices vsp-gs each, red, yellow, green
- 2 Bars vsp-gs
- 4 Duplicating matrices vsp-gs
- 1 Paralleling mandrel metal vsp-f/fs/gs
- 1 Insertion pin vsp-f/fs/gs

REF 430 0648 0



Accessories:

Titanium bar vsp-fs / gs REF 560 0002 0



Insertion pin 2 pieces REF 430 0622 0



Paralleling mandrel 1 piece REF 430 0623 0



## Vario-Soft-Bar-Pattern vsp

#### Dimensions



Product	REF	Length	Width	Height
Bar Friction	430 0646 0	50 mm	1.5 mm	3.5 mm
Bar Friction-Snap/				
joint snap-in	430 0695 0	50 mm	1.5 mm	3.5 mm
Matrix Friction	430 0640 0	6.5 mm	3.0 mm	4.5 mm
Matrix Friction-Snap	430 0634 0	5.6 mm	2.7 mm	2.3 mm
Matrix joint snap-in	430 0628 0	5.7 mm	2.7 mm	4.5 mm

### Fabricating implant-borne restorations using a parallel bar



The bar should be fitted between the implant abutments with a paralleling mandrel. The bar is made of rigid acrylic which can be trimmed easily and quickly.



After casting and trimming, the bar is secured on the abutments with the paralleling mandrel. They should be soldered together to create a non-stressed unit.



Duplicating is always carried out with the yellow matrix intended for the bar system. This provides the optimum conditions for changing the degree of friction later on.



The restoration is blocked out and duplicated using standard methods. No spacer wax should be applied around the matrix.



The matrix is also duplicated and acts as a spacer for the matrix housing in the chrome cobalt framework.



The bar and matrix are simply coated with wax. The remaining sections of the pattern are waxed up as required.



Before pressing the matrix into its housing in the chrome cobalt framework, check the housing for high spots.



The matrix with the desired degree of friction is selected and pressed in with the inserting instrument. The additional snap-in retainers on the matrix provide optimum hold in the housing.



Shows the unterside of the finished restoration with parallel bar and high, firm friction (red matrix). The friction can be increased or reduced as required by replacing the matrix.

## Implant-borne restorations on jointed bars



Once the jointed bar has been soldered and trimmed, the duplicating matrix for the jointed bar snap-in matrix should be placed on it. The underside is blocked out using standard methods. To ensure that the joint matrix fits exactly, the duplicating matrix must not be coated with blocking out wax.



This chrome cobalt framework has been trimmed and checked for high spots and is ready to be fitted with the jointed snap-in matrix with the ideal snapping force for the patient.



The joint snap-in matrix is easily pressed into the chrome cobalt framework with the inserting instrument.

Prior to duplicating, the implant caps and the vertical bar areas are coated with a wax layer with a thickness of 0.3 mm to allow rotational movement of the denture later on. During this process, however, the rounded occlusal end of the bar must not be coated with wax.



# THE GUIDE FOR MODERN IMPLANT PROSTHETICS

A comprehensive guide to the rapid developments in implant prosthetics, including numerous patient cases – some of which from guest authors – and an exciting overview regarding the topic of starting out in dental implantology under the slogan:

"KISS" – keep it simple and successful.

\*\*BURNODER MODER THETTER DER MODER THE TENDER DER MODER THE

REF 9929700D (German) REF 992970GB (English)

## Vario-Soft-Bar vss

Bar system with three interchangeable versions with different degrees of friction. Gentle to the periodontium.

3 precision matrices with different degrees of friction.



Reduced, light friction


Matrices vss green 8 pieces REF 430 0527 0 50 pieces

REF 430 0610 0

Normal, medium friction








High, firm friction

Matrices vss red 8 pieces REF 430 0525 0 50 pieces REF 430 0620 0





The Snap retainers guarantee retention in the matrix housing

The 4 rounded edges of the matrix create guidance grooves to retain it securely in the removeable section of the denture




The double matrix technique provides for reliability

## Assortment

- 2 Patrices vss
- 2 Matrices vss each red, yellow, green
- 1 Insertion pin

REF 430 0523 0

Patrices vss 8 pieces REF 430 0524 0 50 pieces REF 430 0595 0



#### Accessories:



Insertion pin REF 430 0736 3

Press fit pin. Practical, small and cost-effective. Facilitates handling for all users.



The vss attachment can be shortened by up to half its length, making it ideal for many indications even in cases with complex bites.



The patrix can be fitted perfectly, thus ensuring that the gingiva is protected and that the papilla remain free.



Test the excellent, gentle friction for yourself. Both you and your dentist will be enthusiastic about it.



## Dimensions

Product	REF	Length	Width	Height
Bar patrix vss	430 0595 0	48 mm	2.2 / 2°	7.1 mm
Matrix vss	430 0610 0	6.7 mm	3.4 mm	8.0 mm
	430 0594 0	6.7 mm	3.4 mm	8.0 mm
	430 0620 0	6.7 mm	3.4 mm	8.0 mm

The double patrix system requires minimal space, making it perfect for use as an extracoronal attachment.



## Vario-Soft-Bar vss

#### The very gentle friction will impress and enthrall you!



The vss bar patrix can be shortened as required, to suit any particular case. The double patrix can be positioned mesially or distally.



Adapt the underside of the bar to the ridge. The special resin will not distort and is easy and quick to work on.



Shows the bar patrix waxed into place. It can be adjusted with wax at any time.



As the crowns and bar are cast in one piece, no soldering is required and any alloy can be used. This makes the vss ideal for patients with allergies.



The matrix is placed over the double patrices and its underside adapted to fit. The pattern is then blocked out for the chrome cobalt framework, using standard methods.



The model is duplicated with gel or silicone and the investment model is poured. vss can be used with any investment material, thus making it unnecessary to acquire special materials.



The pattern is waxedup over the matrix reproduced in investment material. This guarantees that the chrome cobalt denture base will fit absolutely precisely.



The crome cobalt denture base is fitted down, trimmed and polished. The simple handling and gentle friction will impress you immediately.



Shows the matrix in position. Additional snap retainers guarantee optimum retention in the matrix housing.

#### Individually adjustable friction within reduced time and at low costs even after several years



The patrix is adapted to the situation and then waxed onto the crown. The papilla remain free, as required.



As it is cast in one piece the casting is a homogeneous structure of one alloy, which prevents stresses.



Once the yellow matrix has been fitted, block out the restoration in readiness for fabricating the chrome cobalt denture base. Ensure that no wax is applied around the matrix.



This ensures that the chrome cobalt structure reaches down to the gingiva and the matrix is retained completely in metal.



The press fit pin for the matrix is used to insert it precisely into the chrome cobalt denture base.



Thanks to the guidance grooves, the matrix fits the metall housing perfectly. It can be replaced at any time with a matrix with increased or reduced friction.

# Wax bars

## Wax rod attachments wbgs



vvax rod attachments	Short designation	KEF	Supply form	Order quantity
Head 2 Ø x 50 mm Length	wbgs 2.0	430 0261 0	approx. 170 pieces	
Head 3 Ø x 50 mm Length	wbgs 3.0	430 0262 0	approx. 90 pieces	
Paralleling mandrel 2.0 for wbgs	2.0 ph 2.0	430 0263 0	1 piece	
Paralleling mandrel 3.0 for wbgs	3.0 ph 3.0	430 0264 0	1 piece	

## Wax bar attachments wsgs



Wax bar attachments	Short designation	REF	Supply form	Order quantity
micro 2.2 bar height x 50 mm	wsgs m 2.2	430 0271 0	approx. 250 pieces	
normal 3.0 bar height x 50 mm	wsgs n 3.0	430 0272 0	approx. 125 pieces	
Paralleling mandrel 1.6 for wsgs i	m 2.2 ph 1.6	430 0623 0	1 piece	
Paralleling mandrel 2.2 for wsgs i	n 3.0 ph 2.2	430 0270 0	1 piece	

## Wax T-attachments wtgs



Wax T-attachments	Short designation	REF	Supply form	Order quantity
Wax T-attachments 2.75	wtgs 2.75	430 0275 0	approx. 150 pieces	
Wax T-attachments 3.5	wtgs 3.5	430 0276 0	approx. 90 pieces	
Paralleling mandrel 2.75 for wtgs	2.75 phT 2.75	430 0277 0	1 piece	
Paralleling mandrel 3.5 for wtgs 3	.5 phT 3.5	430 0278 0	1 piece	



Product	REF	Ø	Length	Width	Height	Max. reduction
Wax rod attachments	430 0261 0	Head 2.0 mm	50 mm	2.0 mm	4.5 mm	individual
	430 0262 0	Head 3.0 mm	50 mm	3.0 mm	5.5 mm	individual
Wax bar attachments	430 0271 0	-	50 mm	2.2 mm	1.5 mm	individual
	430 0272 0	-	50 mm	2.2 mm	3.0 mm	individual
Wax T-attachments	430 0275 0	_	50 mm	2.75 mm	3.4 mm	individual
	430 0276 0	-	50 mm	3.5 mm	4.75 mm	individual

Customer No.
Date, signature

Additional order:

# Wax bars

## Wax bars wstg



Wax bars	Short designation	REF	Supply form	Order quantity
1.6 x 8 x 50 mm	wstg 1.6	430 0265 0	approx. 65 pieces	
1.9 x 4 x 50 mm	wstg 1.9	430 0266 0	approx. 120 pieces	
2.2 x 6 x 50 mm	wstg 2.2	430 0267 0	approx. 65 pieces	

Paralleling mandrel 1.6 for wstg 1.6	ph 1.6	430 0623 0	1 piece	
Paralleling mandrel 2.2 for wstg 1.9 and wstg 2.2	ph 2.2	430 0270 0	1 piece	

## Wax bar hinges wsgl



Wax bar hinges	Short designation	REF	Supply form	Order quantity
micro 2.2 bar height x 50 mm	wsgl m 2.2	430 0273 0	approx. 300 pieces	
normal 3.0 bar height x 50 mm	wsgl n 3.0	430 0274 0	approx. 160 pieces	
Paralleling mandrel 1.6 for wsgl i	m 2.2 ph 1.6	430 0623 0	1 piece	
Paralleling mandrel 2.2 for wsgl i	n 3.0 ph 2.2	430 0270 0	1 piece	

#### Round wax bars wstr



Round wax bars	Short designation	REF	Supply form	Order quantity
1.5 Ø x 50 mm	wstr 1.5	430 0279 0	approx. 400 pieces	
1.8 Ø x 50 mm	wstr 1.8	430 0280 0	approx. 300 pieces	
2.0 Ø x 50 mm	wstr 2.0	430 0281 0	approx. 250 pieces	

Dimensions	Product	REF	Ø	Length	Width	Height	Max. reduction
	Wax bars	430 0265 0	_	50 mm	1.6 mm	8.0 mm	individual
		430 0266 0	_	50 mm	1.9 mm	4.0 mm	individual
		430 0267 0	-	50 mm	2.2 mm	6.0 mm	individual
	Wax bar hinges	430 0273 0	_	50 mm	1.4 mm	2.2 mm	individual
		430 0274 0	_	50 mm	2.1 mm	3.0 mm	individual
	Round wax bars	430 0279 0	1.5 mm	50 mm	_	_	individual
		430 0280 0	1.8 mm	50 mm	-	-	individual
		430 0281 0	2.0 mm	50 mm	_	_	individual

Sender (stamp):	Customer No.	Additional order:
	Date, signature	

# **NEW OPPORTUNITIES FOR SUCCESS**

## FURTHER TRAINING AT AND WITH THE BREDENT GROUP ACADEMY



We strive to be amongst the best. In this spirit, we have reviewed our course programme and the results are now being introduced under the bredent group academy banner.

The bredent group academy course programme is aligned to the strategic direction of the bredent group. On the one hand, we have the part that is in line with the market, oriented towards tailored solutions to meet the needs of the market. Such solutions have been developed according to the specific segment. By contrast, there is the market-oriented segment. This is geared towards a changing market, such as economic cycles and trends.

However, it is far more important to recognise users' wishes from the outset. This means being intuitive to what users want and need, and knowing what the future holds - in terms of materials, systems and treatment concepts. Only this way can we as a manufacturer ensure that our company progresses and understands the associated obligations to inform, support and help you, as our client, on the path to commercial success.

Further training is a key component of this duty, since the cornerstone for new technology and a smooth transfer of knowledge are being laid. The instructors

are the life of this course – all of them experts in their field – equipped with the vision and the necessary drive to develop innovations and put them into practice. This programme unites these visions and appeals to those who never cease to develop themselves further.

The individual development and optimisation of your practice, your laboratory and your marketing has a decisive effect on realising success in your practice or laboratory, your income and therefore your abilities, ideas, wishes and dreams.

A new range of events in the areas of patient communication, hygiene, conflict management, practice marketing, data protection and social media, in addition to the exciting beginners' and improvers' courses in implantology and dental technology, can help you become even more successful. Not only is further training at the forefront of your professional qualifications, but so too is positioning your business as a successful brand on the market.

Join us and strive to be amongst the best – for the benefit of you and your patients!

Find out more about our course programme at www.bredent.com.



## Lock

Swivel-type lock or locking pin? bredent offers cost-effective solutions for secure, fixed, removable dental prostheses which creates a secure feeling for the patient. Simple processing and handling of the lock ensures a high level of comfort.

#### Indications

- Extracoronal applications
- Plastic fitting
- Model casting fixation

#### **Product variations**

- Swivel-type lock sr and src easy operation for patients
- Locking Pin bs1 the individual locking pin
- KS-Lock assembled push lock
- Locking Pin Easy-Snap offers secure grip with simple handling
- Locking Pin Snap-System the versatile lock!
- Activatable locking pin individual friction adjustment possible

#### **Properties**

- Cost-effective purchase
- Time-saving and secure processing
- Cost-effective solutions

#### Materials

- Titanium
- Precious metal



KS Lock – assembled push lock



Locking pin Easy-Snap – offers secure grip with simple handling



Swivel-type lock sr and src – easy operation for patients



Locking pin Snap-System – the versatile lock!



Locking pin bs1 – the individual locking pin



Activatable locking pin – individual friction adjustment possible



# Swivel-type lock sr

These particularly slender patterns can be used for a wide range of custom applications for all removable, passive designs – perfect for implants.





Latch retainer with integral shear distributor left, 4 pieces REF 430 0735 8 right, 4 pieces REF 430 0735 9



Latch tongue 4 pieces REF 430 0735 7

# Assortment

REF 430 0736 2

14 pieces Swivel-type lock sr left + right

#### Assortment


14 pieces

Swivel-type lock sr left REF 430 0730 5



14 pieces

Swivel-type lock sr right REF 430 0730 6



Latch tongue and latch box



housing left, 4 pieces REF 430 0730 9 right, 4 pieces REF 430 0731 0

Shear distributor



Latch box 4 pieces REF 430 0735 6



Oxide-steel pins 20 pieces REF 430 0293 0



Latch spring
Guaranteed for 5 years
10 pieces
REF 430 0334 0



Latch tongue and latch box in a shear distributor housing

Cross-section through the Latch System sr



locked



unlocked



Fully assembled Latch System sr

#### Accessories:



Paralleling mandrel universal 1 piece REF 360 0115 1



# Swivel-type lock sr

This cost-effective latch allows you to calculate your prices to optimize your profit.



Classic shear distributor with Interlock and a complete latch system. The chrome cobalt framework and latch system were luted with dtk without creating stresses.



Latch box with individually modelled latch box housing. Ideal when minimal space is available.



Construction of a monoreducer with integrated shear distributor. The usage of the latch box housing avoids modelling of the shear distributor.

Applications for combined fixed/removable appliances with classic shear distribution arms



The latch retainer can be fitted onto the papillae accurately.



The non-soldered, onepiece casting reduces the costs and the number of alloys used in the mouth.



The shear distribution arm pattern is built up with Pi-Ku-Plast brushon resin to guarantee that all details are reproduced.



The latch system provides numerous combinations for fabricating custom restorations.

Applications for combined fixed/removable appliances with a space-saving latch retainer



As the latch retainer is designed to fit around the papillae, it can be waxed close to the crown with a paralleling mandrel.



The shear distributor shoulder on the patrix eliminates the need for labour intensive milling, which saves time and money.



The slender design of the latch system allows the shear distributor housing to be waxed up as required.



The restoration is designed so as not to stress the abutment teeth.

"Monoreducer" with integral shear distributor and custom designed latch box housing



Patrix with integral milled shoulder for the shear distributor - saves time and space.



The latch tongue swivels horizontally into the latch retainer.



The latch system is very easily assembled.



Once the latch has been opened, the partial denture can be released without stressing the abutment tooth.

#### Dimensions



Product	REF	Length	Width	Height	Ø
Latch retainer	430 0735 9	4.0 mm	2.9 mm	4.3 mm	_
Latch tongue	430 0735 7	5.8 mm	3.8 mm	2.9 mm	_
Latch box	430 0735 6	6.2 mm	5.0 mm	2.9 mm	_
Shear distributor housing	430 0731 0	6.4 mm	5.9 mm	4.8 mm	_
Oxide-steel pins	430 0293 0	10.0 mm	_	_	1.0 mm

# Swivel-type lock src

#### Ceramic spacer for the simple fabrication of locks in the one-piece casting technique.

Swivel-type locks for CoCr restorations: low-cost, accurate and time-saving.





#### Wax latch retainer

- is cast together with the anchor crowns, hence reduction of metals in the mouth
- integrated shear distributor provides enhanced esthetics and simplifies the fabrication



#### Ceramic latch blade with lock axle

- Latch box premodelled in wax, hence quick fabrication is possible
- is integrated into the CoCr model
- only sandblasting required after casting
- creates precise fitting surfaces for metal, latch blade and latch axle



#### Latch spring

 determines the position of the lock blade when locked or unlocked and provides additional safety for the patient



#### Titanium latch blade

- shape matched exactly with the ceramic lock blade prefabricated precision lock blade
- ensures efficient processing

#### Latch axle

- stainless-steel, hence orally stable
- matches exactly with the ceramic spacer for the latch axle, simplifies the integration

High-precision ceramic patterns are available which reduce the amount of work tremendously and simplify the fabrication of an individual lock. The ceramic pattern designs are based on the design of the latch blade and the latch axle so that they fit exactly into each other. The latch box is cast in one piece together with the CoCr structure using the one-piece casting technique. Soldering is not required. Accordingly, the amount of alloys used intraorally is reduced and the costs for individual lock restorations are lowered.



Latch blades src ceramic 2 pieces REF 430 0738 5



Latch retainer left 4 pieces REF 430 0735 8



Latch retainer right 4 pieces REF 430 0735 9



Latch spring 10 pieces REF 430 0334 0



Latch tongue titanium 2 pieces REF 430 T735 7



Oxide-steel pins 20 pieces REF 430 0293 0

## Assortment

10 pieces

Latch tongue src left + right

REF 430 0738 8

#### Accessories:



Paralleling mandrel universal 1 piece REF 360 0115 1



# Swivel-type lock src

#### Ceramic spacer for the simple fabrication of a latch retainer.



The latch retainer is waxed with the paralleling mandrel to the primary construction according to the path of insertion. When producing a bar restoration, the integrated shear distributor may be covered with wax.



After casting, prepare the model for duplicating. Block out the lock of the latch retainer so that approx. 0.5 mm of the margin remains visible after duplicating.



The latch retainer can be easily seen on the investment model. The ceramic pattern can be safely positioned in the lock.



Trim the ceramic pattern with a separating disc in a way that it fits precisely into the lock of the latch retainer and ...



... ends exactly at the latch retainer but can still be positioned safely in the lock of the latch retainer.



Fix the ceramic pattern with the axle and attach with wax.



Complete the model in accordance with the situation and integrate the ceramic pattern. The axle should stand out on both sides of the model.



The titanium swivel-type lock is fitted into the sandblasted housing. Put the latch spring behind the titanium swivel-type lock and fix with the axle.

#### **Dimensions**



Product	REF	Length	Width	Height	Ø
Latch retainer	430 0735 9	4.0 mm	2.9 mm	4.3 mm	_
Latch tongue titanium	430 T735 7	5.8 mm	3.8 mm	2.9 mm	_
Oxide-steel pins	430 0293 0	10.0 mm	_	_	1.0 mm

# Locking Pin bs 1

The lock axles can be shortened according to the respective situation and an individual unlocking lens can be added.



The locking pin bs 1 is perfectly suitable for the use in the anterior area. The lock can be opened using a bent wire and then the denture can be removed.



Pin axles 2 pieces REF 450 0006 4

Refill packages



Bolt screws 2 pieces REF 450 0006 5



The locking pin bs 1 can be used individually. Thanks to its size the locking pin bs 1 is perfectly suitable for unilateral removable dentures. An undercut can be integrated into the pin axle to open it.



Wax screws 2 pieces REF 430 0748 2



Wax sleeves 2.0/1 2 pieces REF 450 0007 2



The completely individual solution:

The unlocking lens is prepared individually using denture resin or composite. This way the unlocking device is no longer visible.



Auxiliary modelling elements 2.0 x 3.5 2 pieces REF 450 0007 0



Auxiliary modelling element 2.0 2 pieces REF 450 0006 3



Auxiliary modelling elements 1.3 2 pieces REF 450 0007 1



First tap M 1.6 1 piece REF 330 0116 V



Last tap M 1.6 1 piece REF 330 0116 F



Tap handwheel 1 piece REF 330 0115 3



Tungsten carbide center drill Ø 1.4 1 piece REF 330 0066 0



Diatit-Multidrill spiral drill 2.0 1 piece REF 330 0072 0



Assortment 17 pieces Locking Pin bs 1 REF 450 0006 2



Accessories:

Milling and drilling oil 20 ml REF 550 0000 8




# Locking Pin bs 1

Type 1: Bolt screw in metal framework



Wax up the pattern using standard methods. The prefabricated wax bar (REF 430 0265 0) is perfectly suitable to allow quick fabrication.



Prepare a small depression with the center drill and drill a hole through the bar using the Diatit-Multidrill spiral drill.



Assemble the auxiliary modelling elements 2.0 and 1.3 and position them in the drill hole so that between primary element and auxiliary element 1.3 ...



... a minimum space of 1.5 mm is obtained. Fix the auxiliary modelling elements with Pi-Ku-Plast.



Complete the model in accordance with the situation. Remove the auxiliary modelling elements, invest and then cast.



After casting, cut the thread with the taps whilst adding a rich quantity of milling and drilling oil.



Cut the lock axle according to the respective situation. Add an unlocking device and fix the pin axle with the bolt screw.



The individually fabricated locking pin can also be used if only limited space is available.

Type 2: Bolt screw in resin saddle



Insert the auxiliary modelling element 2.0 x 3.5 in the center of the drill hole as spacer for the investment model. Prepare the model for duplicating.



Place the wax sleeves on the investment model to obtain an accurate, uniform wall thickness of the model.



Complete the model (waxing up). The plugs of the auxiliary modelling element remain visible.



Fix the wax screw in the auxiliary modelling element 2.0 and fit it in the drill hole. There must be a distance of 1 mm between the wax screw and the secondary construction.



After completion the wax screw is removed with steam. Precise has been achieved in the resin which safely accepts the bolt screw. The pin axle is cut individually and an unlocking device is added. If required, the lock axle can be veneered in a suitable shade. A reliable, simple solution for any type of removable dentures.

#### Dimensions




Product	REF	Ø/Thread	Length	Max. reduction
Pin axle	450 0006 4	2.0 mm	15.0 mm	individual
Bolt screw	450 0006 5	M 1.6 x 0.35	4.4 mm	-

## Lock attachments

## KS-lock





Locks 2 pieces REF 450 0007 9



Titanium sleeves 2 pieces REF 450 0007 8

#### Prefabricated sliding lock with individual opener.

- limited number of components
- simple design
- can be shortened individually
- individual design of opener
- for motorically handicapped patients
- snap mechanism when locking and unlocking
- perfectly suitable for implant structures



Use resin to model the primary element with a retainer for the lock. For this purpose attach the lock with titanium sleeve.



Try the lock with titanium sleeve in the cast primary element.



Primary element prepared for the supraconstruction – with or without electroplating.



Supraconstruction fabricated for the veneer.



Glue the titanium sleeve into the supraconstruction.




Fit the cast lock with the customized opener into the construction.



The opened lock with individually designed opener. Perfectly suited for motorically handicapped patients.

#### Dimensions



# Locking Pin Easy-Snap




Perfect locking pin system for secure grip of the prosthesis.

Its small size permits versatile use.

- New mechanism provides secure hold for the prosthesis
- Noticeable snap when opening and closing gives patients more security
- Small size for every situation
- Can be used even in difficult locations
- Simple fitting possible in three variants





The snap ring offers secure grip when opening and closing the locking pin axis. The simple mechanism offers the highest degree of security.



The castable variant allows for use regardless of the alloy.

Castable up to 1800°C.



The polymerisation of the locking pin leaves all options open. Simple use for beautiful teeth!

# Locking Pin Easy-Snap E

The primary construction is always prepared in the same manner for every application. This reduces the need for clarification and accelerates the manufacturing process.



Simple modelling of the primary construction according to the model situation. A set-up for exactly determining the position of the lock is always advantageous.



The hole is measured using the centre drill. Slipping of the Diatit multidrill is thus prevented.



Drilling is simplified using drilling and milling oil.
Drill only once, as this prevents the formation of an oval hole.



The hole is closed with the wax.



The wax is removed again by hand with the Rapidy 2.0. In doing so, a slight hollow is created which positions the investment material model exactly.



The model is completely prepared for doubling. If you are working in the withdrawal procedure, omit this step.



Cast the investment material model using the bredent doubling system. The high level of design accuracy of the investment material facilitates further processing.



The wax guide is fixed lightly in the cast mould. The framework is modelled with a minimum thickness of 0.5 mm.



The wax guide is removed and an image is created for the locking pin.



The exit hole for the locking pin is easily recognisable on the opposite side.



The locking pin is fixed in the correct position for completion with Qu-resin.



The prosthesis is com-

#### Accessories:



Tungsten carbide center drill REF 330 0066 0



Diatit-Multidrill REF 330 0073 0



Rapidy Microbur REF H001 NH 21



Milling and drilling oil 20 ml REF 550 0000 8



DTK adhesive REF 540 0010 6



FGP insulating liquid REF 540 0102 7



Qu-resin dentin 50 ml cartridge REF 540 0116 6

Qu-resin rosa 50 ml cartridge REF 540 0116 5



Pi-Ku-Plast resin REF 540 0017 3 Assortment small REF 540 0019 6



# Locking Pin Easy-Snap A

Sleeves made from platinum-iridium alloy offer time-saving processing for precious metal and non-precious metal constructions, with the exception of titanium.



For patients with less dexterity, the lock can be constructed in such a way that it can even be opened with a thin object from the buccal side. An individual stop must be fitted to this so that the lock is not opened accidentally.



The individual bridge sectioning attachment is suitable for wax modelling. Rapid modelling saves time.



Manufacturing out of brealloy enables a space-saving model to be created for an aesthetic dental prosthesis.



The position of the lock is measured using the centre drill and a small mould is fitted.



The hole is drilled for the locking pin using the Diatit-Multidrill with copious amounts of milling and drilling oil.



The small amount of space required by the Locking Pin Easy-Snap facilitates problem-free positioning of the bore hole in the male part.



The locking pin sleeve which can be cast on, is fitted to modelling pin A and pushed into the locking hole of the patrix until it stops.



The modelling pin A with Locking Pin sleeve which can be cast on is moulded with Pi-Ku-Plast to its largest diameter.



The blue clip forceps hold the modelling pin steady whilst removing this from the model. The sleeve which can be cast on remains within the model.



The modelling is carried out and cast in accordance with the bredent casting technique. The sleeve is fixed in the correct position with the investment material.



In order not to damage the sleeve which can be cast on, the investment material is blasted with glass beads. This means that it maintains its shape and performs the correct function.



The snap ring is set on the impression pin and pressed into the sleeve. It springs into its planned position.



The lock axis is then pushed in. The lock holds firm due to the mechanism and offers the patient the highest level of comfort.

# Locking Pin Easy-Snap

#### Bonding/polymerisation of locking pin Easy-Snap E





Locking Pin Easy-Snap E 1 piece REF 440 0N65 8



#### Moulding of Locking Pin Easy-Snap A





Locking Pin Easy-Snap A 1 piece REF 440 0N65 9



Device for waxing-on 1 piece REF 440 0066 1



Locking Pin Easy-Snap E 1 piece REF 440 0N65 8





Locking Pin Easy-Snap A 1 piece REF 440 0N66 0

Modelling pin A 1 piece REF 440 0N65 5



Assortment
4 pieces
Locking Pin Easy-Snap E
2 Locking Pin Easy-Snap E
2 Modelling pin E
REF 440 0N65 2



Assortment
4 pieces
Locking Pin Easy-Snap A
2 Locking Pin Easy-Snap A
2 Modelling pin E
REF 440 0N65 4



#### Assortment

5 pieces

Locking Pin Easy-Snap E Model casting fitting
2 Locking Pin Easy-Snap E
2 Ceramic spacer E
1 Device for waxing-on
REF 440 0N65 3



## Assortment

A spieces
Locking Pin Easy-Snap E Plastic fitting
Locking Pin Easy-Snap E
Device for waxing-on
REF 440 0N65 1

#### Accessories:

Snap ring 10 pieces REF 440 0N66 3 Matrix inserting instrument 1 pieces REF 440 0N66 2



# Locking Pin Snap System

Suitable for numerous applications in combined work.



The soft, resin-supported guidance results in a soft snap of the locking pin during locking in the closed or open position.



The snap ensures safe locking in closed position.

closed





open The snap informs the patient that the lock is completely open and the denture can be removed.

Locking Pin Snap E The alloy is not relevant








All metal parts are made of titanium. The biocompatible plastic matrix provides long service life and soft snap-friction.





## Locking Pin Snap A Fast and precise integration in precious metal supply






The platinum-iridium-containing alloy of the locking pin sleeve allows to save much time and ensures high precision when casting onto the metal framework.

## Lock attachments

# Locking Pin Snap E

#### Glueing in the chrome cobalt framework





Locking Pin E 1 piece REF 440 0065 8



Ceramic spacer E 2 pieces REF 440 0065 7



Device for waxing-on 1 piece REF 440 0066 1

#### Integration in resin



#### Assortment

5 pieces

Locking pin snap E for integration in chrome cobalt

- 2 Locking Pin E
- 2 Ceramic spacer E
- 1 Device for waxing-on

REF 440 0065 3

## Assortment

3 pieces

Locking pin snap E for integration in resin

- 2 Locking Pin E
- 1 Device for waxing-on

REF 440 0065 1

#### Accessories:



Tungsten carbide center drill REF 330 0066 0



Diatit-Multidrill REF 330 0073 0



Rapidy Microbur REF H001 NH 21



Milling and drilling oil 20 ml REF 550 0000 8



Wax bars wstg 1.6 x 8.0 mm REF 430 0265 0



FGP insulating liquid REF 540 0102 7




DTK adhesive REF 540 0010 6




# Locking Pin Snap E



Waxing up of the pattern and casting is carried out using standard methods.




The position of the locking pin is determined with the centre drill and a small depression is prepared.



The hole for the pin is drilled with the Diatit Multidrill with a diameter of 1.5 mm.



The hole for the pin is filled with wax before duplicating.



A depression is scraped on both sides using a Rapidy Microbur 2.1 mm



The pattern is prepared for duplicating and duplicated in the usual way.



Exact reproduction of the depressions on the bar is required.



The ceramic spacer E is exactly positioned with the device for waxing on.



The ceramic spacer E is integrated into the wax pattern up to its largest diameter.



The spacer is sandblasted with a maximum grain size of 110  $\mu$  at a pressure of 4 bar.



In order to try the function, the pin is inserted into the assembled construction.

## Glueing in of Locking Pin Snap E.

Contact points that must not be glued must be covered with FGP insulating liquid.



The hole in the bar and 2-3 mm in the vicinity.



The contact area of the locking pin lens at the secondary element.



The locking pin stud up to the locking pin sleeve.



The contact area of the locking pin lens at the secondary element.



Primary and secondary element are assembled. A drop of DTK adhesive is evenly spread in the hole in the secondary element.





The locking pin sleeve is covered with a thin layer of DTK adhesive and pressed into the secondary element. Excess adhesive residues are removed after hardening of the DTK adhesive.



# Locking Pin Snap E

#### Locking Pin Snap E for integration in resin.


Easy, fast and secure integration.



The crown is modelled with a shear distributor with interlock and the end of the bar is waxed on.



Casting and polishing are carried out after casting.



A small depression is prepared with the Diatit centre drill and in this way the position of the drillhole is determined.



The Diatit-Multridrill is safely fixed by centering.



The hole for the pin is filled with wax before duplicating.



A depression is scraped on both sides of the Locking Pin patrix using a Rapidy Microbur 2.1 mm.



The pattern is prepared for duplicating and duplicated.



The small depressions are reproduced in the investment compound model.





The plugs of the device for waxing-on lock in position in the depressions.



The cylindrical plugs are integrated in the pattern using modelling wax.



This way two round apertures are obtained on both sides.



The two apertures are parallel to the axis due to the drilled hole.



The holes are prepared using a Diatit-Multidrill with a diameter of 1.5 mm.



The Locking Pin is easily tried in.



The Locking Pin is fixed to the chrome cobalt framework using resin.



The Locking Pin lens is integrated into the wax pattern up to its outer margin and the denture is completed.



# Locking Pin Snap

When the wax pattern is lifted from the model in order to invest it, there are two options:

Glueing in Locking Pin Snap E





Locking Pin Snap E 1 piece REF 440 0065 8



Modelling pin E 1 piece REF 440 0065 6

## Casting in Locking Pin Snap A



Assortment

4 pieces Locking Pin Snap E 2 Locking Pin Snap E 2 Modelling pin E REF 440 0065 2

### Assortment

4 pieces Locking Pin Snap A 2 Locking Pin Snap A 2 Modelling pin A REF 440 0065 4



Locking Pin Snap A 1 piece REF 440 0066 0



Modelling pin A 1 piece REF 440 0065 5



Locking Pin Snap 1 piece REF 440 0065 9

#### Accessories:



Tungsten carbide center drill REF 330 0066 0



Wax bars wstg 1.6 x 8.0 mm REF 430 0265 0



DTK adhesive REF 540 0010 6



Pi-Ku-Plast resin REF 540 0017 3 Assortment small REF 540 0019 6



Milling and drilling oil 20 ml REF 550 0000 8




FGP insulating liquid REF 540 0102 7




## Lock attachments

# Locking Pin Snap E

#### Integration of Locking Pin snap in any alloy.



The wax bars by bredent are fitted in individually.



The position of the locking pin is determined with the centre drill and a small depression is prepared.



The Diatit-Multridrill is safely positioned by centering.



Thanks to little space required by the Locking Pin snap, the hole can be easily positioned in the patrix element.



The modelling pin E is inserted in the patrix hole up to the stop.



The modelling pin E is integrated in the pattern using Pi-Ku-Plast resin and modelling wax.



After completion of the pattern, the modelling pin E is removed by turn-ing it slightly with a pair of pliers.



The investment compound in the Locking Pin housing is sandblasted with a grain size of 110  $\mu$  and a pressure of 4 bar.





Insulating and glueing in are carried out as described on page 199.



# Locking Pin Snap A

Time-saving casting-on to precious metal secondary constructions.





The wax bars by bredent are fitted in individually.



The position of the locking pin is determined with the centre drill and a small depression is prepared.




The bar is perforated with the Diatit-Multi-drill whilst adding rich quantities of milling and drilling oil.



Thanks to the little space required by the locking pin snap, the hole can be easily positioned in the patrix element.



The cast-on locking pin sleeve is put onto the modelling pin A and inserted into the locking pin hole of the patrix up to the stop.



The modelling pin A with the cast-on locking pin sleeve is integrated in the pattern up to its largest diameter using Pi-Ku-Plast resin and modelling



After waxing up, the modelling pin A is removed with a slight turn.



In order not to damage the cast-on locking pin sleeve, the investment compound is removed with glass beads.



The Locking Pin Snap is pressed in the locking pin sleeve that has been cast in.



#### Dimensions



Product	REF	Ø Axle	Ø Ring	Length	Ø	Max. reduction
Locking Pin Snap	440 0065 8	1.5 mm	3.5 mm	6.25 mm	_	-
Locking Pin Snap sleeve	440 0066 0	2.8 mm	-	3.6 mm	2.8 mm	_

# Locking Pin activatable

Universal active. The pin can be located on either the oral or buccal surfaces.





Locking pin: Resistant to the oral environment. Made of special steel, with spark eroded activating slot and wax spacer.

## Assortment

Pack of 2 sets Locking pin system\*

1 blocking out matrix

- 2 locking pin matrices
- 2 locking pin patrices
- 1 steel pin 1.5 mm
- 2 locking pins, activatable

REF 430 0445 0



The metal matrix preformer ensures that the pin guidance is absolutely parallel.



Locking pin matrix: This prefabricated wax matrix reduces the time required when waxing-up chrome cobalt appliances.



Locking pin patrix with concave surface for waxing onto the pattern.



"Mini" locking pin: The smallest in the bredent Locking Pin System.



The "Mini" locking pin matrix, simplifies fitting of the locking pin.



- 1 blocking out matrix
- 2 locking pin matrices
- 2 locking pin patrices
- 1 steel pin 1.5 mm
- 2 locking pins, activatable

REF 430 0460 0



"Mini" locking pin patrix: Saves space, ideal for anterior use.

\* For further information see price list.



The pin viewed from the lingual direction. The pin passes through an extracoronal retaining lug.



This pin can be operated from the buccal aspect.



The spark eroded activating slot is simply activated from the underside.



# Locking Pin activatable

#### Locking Pin System



Wax-up the pattern using standard methods, and then use the paralleling mandrel to wax the patrix into place.



Adapt the underside of the patrix to fit the model and integrate it into the shear distributor.



Trim and polish the framework before applying the porcelain.



Place the matrix pre-former in position and secure it with the oxidized steel pin. Block out the underside.



Remove the matrix preformer and fill the pin aperture with wax, leaving a slight depression.



Position the wax matrix correctly on the investment model.



Wax-up the denture base framework using standard methods.



Shows the chrome cobalt framework after casting and trimming. The apertures for the pin have been aligned accurately.



Shows the try-in, with temporary pin made of clasp wire.



Insert the pin until the wax sleeve touches the chrome cobalt framework.



Shows the completed saddle: The locking pin is pulled to open it.



Viewed from the underside. The slot in the pin enables it to be activated easily.

#### Dimensions



	Product	REF	Ø Axle	Length	Width	Height	Ø Ring
	Locking pin	430 0459 0	1.5	6.2 mm	_	_	2.9 mm
	Locking pin, mini	430 0500 0	1.5	4.5 mm	_	_	2.9 mm
	Locking pin matrix	430 0458 0	_	5.6 mm	2.5 mm	4,1 mm	_
í	Locking pin matrix, mini	430 0490 0	_	4.6 mm	1.9 mm	3.6 mm	_
	Locking pin patrix	430 0458 0	_	5,4 mm	3.7/1.2 mm	3.4 mm	_
١	Locking pin patrix, mini	430 0490 0	_	4.3 mm	3.7/0.9 mm	2.8 mm	_

# **IMPLANT PROSTHETICS**

# ONE-STOP INTERDISCIPLINARY COMPETENCE -WHAT THE BREDENT GROUP STANDS FOR



As long ago as 1974 and with the prime intention of ensuring the long-lasting success of its products when used by customers, the bredent group set out to do more than simply manufacturing a series of stand-alone products.

By providing a comprehensive range of products, system solutions and treatment concepts for dental technology and dentistry that can be combined in a professionally competent manner, the bredent group is amongst the few companies that are held in equally high esteem by customers and professionals alike.

In close cooperation with reference laboratories and practices, well-matched, patient-centred products, system solutions and treatment concepts are developed, which can be implemented efficiently and rapidly in the laboratory and the practice.

The positive effect of ground-breaking development!



#### Screw connections

Using innovative screwing, bredent facilitates the simple manufacture of removable dental prostheses. Secure hold due to screws that cannot be loosened individually achieves a high degree of security.

#### Indications

- Screwing or fixation
- Implant suprastructures
- Individual solutions possible

#### **Product variations**

- Security-Lock-System three different processing options
- Friction Splint FS1 the threadless fixation
- Bridge-sectioning attachment diagonal or transverse screwing
- Toolkit for individual screws endless possibilities
- Transverse fixation and assembled transverse fixation

#### **Properties**

- Time-saving and secure processing
- Cost-effective solution with a high level of individuality

#### Materials

The screws are manufactured from biocompatible titanium.



Bridge-sectioning attachment – diagonal or transverse screwing



Toolkit for individual screws - endless possibilities



Security-Lock-System – three different processing options



Friction Splint FS1 – the threadless fixation



Transverse fixation and assembled transverse fixation

#### Universal screwdriver set



Screwdriver set for 98% of all screws available on the market. To be inserted into the torque ratchet, adjustable from 10 to 40 Ncm. This way screws can be turned in correctly and safely.



Universal screwdriver set to loosen and tighten all types of screwed implant abutments.



Universal screwdriver set with instruments REF 310 0001 2


On the lid you can find important information required for the quick selection of the necessary screwdriver and the torque needed to tighten the screw.







Torque ratchet REF 330 0115 5

Torque adjustable from 10 to 40 Ncm.







## Universal screwdriver set for contra-angles



Screwdrivers with seating for contra-angles. Thanks to the integrated torque they simplify turning in screws with special motors. In conjunction with the adapter, the screwdrivers can also be used with the torque ratchet.



Universal Screwdriver-Set for contra-angles, with instruments REF 310 W001 2 Universal Screwdriver-Set for contra-angles, without instruments REF 310 W001 1



#### Accessories:




Torque ratchet REF 330 0115 5

Torque adjustable from 10 to 40 Ncm.



Ratchet adapter REF 580 0116 8



	_	Screv	vdriver short	
6 Screwdriver	1	short	Torx 6	REF 310 W0K0 1
☐ Screwdriver	2	short	slotted 1.6	REF 310 W0K0 2
Screwdriver	3	short	slotted 2	REF 310 W0K0 3
Screwdriver	4	short	Allen 0.03"	REF 310 W0K0 4
Screwdriver	5	short	Allen 0.05"	REF 310 W0K0 5
0.9 Screwdriver	6	short	Allen 0.9	REF 310 W0K0 6
Screwdriver	7	short	Allen 1.0	REF 310 W0K0 7
1.2 Screwdriver	8	short	Allen 1.2	REF 310 W0K0 8
1.8 Screwdriver	9	short	Allen 1.8	REF 310 W0K0 9
Screwdriver	10		Allen 2.5 onl	y available as long type
1.3 Screwdriver	11	short	Square 1.3	REF 310 W0K1 1
5.5 Screwdriver	12	short	Torx 5.5	REF 310 W0K1 2

### **Screw connections**

## Screwdriver long



Screwdriver long 1 piece REF 330 0081 2

The long screwdriver allows perfect visual control of the horizontal path of screwing in the laboratory. The screw connection can be more easily achieved by the dentist. For screws with 0.9 mm hexagon socket.

## Screwdriver short



Screwdriver short 1 piece REF 330 0069 0 Ideal for practice and laboratory. The grooved handle simplifies turning in of screws since safe hold is ensured. For screws with 0.9 mm hexagon socket

## Screwdriver for contra-angles



Screwdriver for contra-angles 1 piece REF 330 0081 3 For mechanical turning in of screws with 0.9 mm hexagon socket. The use of special motors allows to control the torque.

## Screwdriver-Set



#### Assortment

3 pieces

1 x Screwdriver long

1 x Screwdriver short

1 x Screwdriver for contra-angles

REF 330 0081 0

## Screwdriver is



Screwdriver is for contra-angles 1 piece REF 460 0001 0



Screwdriver is manual short 1 piece REF 460 0001 1 Special screwdrivers for the vks-oc rs abutments. Suitable as manual screwdriver and for contraangles for enhanced control of the torque with special motors.

## Screwdriver for stud-head screw



Screwdriver for stud-head screw 1 piece REF 330 0116 4 Screwdriver for the stud-head screw vks-oc/sg 1.7 exchangeable stud.



## Milling and drilling oil



Milling/drilling oil 20 ml REF 550 0000 8

#### Especially developed for the milling and drilling technique.

This milling and drilling oil does not contain any ethereal additives. Accordingly, the evaporation temperature is increased considerably; gumming of the oil is no longer possible. Due to special components and the particular consistency, the oil film remains between the metal and the milling tool. This results in the fact that metal chips come out of the cutting sections of the burs more quickly and thus easier milling is possible. The cutting performance and the service life of the milling tools is enhanced correspondingly. By using this milling and drilling oil, more material can be removed while exerting less pressure and obtaining a considerably smoother surface. The oil that has been especially developed for dental techniques withdraws the heat during processing of the object more quickly and avoids overheating of the milling and drilling tools.



When tapping, always use a rich quantity of milling and drilling oil. This simplifies turning in of the tap.



The surface of the object becomes clearly smoother if the oil is used.



This milling and drilling oil avoids overheating of the milling and drilling tools; consequently, the service life of the milling tools is increased considerably.

Use:

Always use a rich quantity of milling and drilling oil during centring, drilling, milling and tapping.

## **Screw connections**

## Transverse fixation



In angled implants, diagonal screwing often means that the screw channel is on the surface of the crown, which also cannot always be compensated for using appropriate angled abutments. Aesthetic problems can therefore occur. In the area of the lateral teeth, the opening of the screw channel can lead to problems with diagonal load distribution.





The solution for this is transverse fixation. The bredent group offers various systems for this:

- Assembled transverse screwing to the SKY system
- Individual transverse fixation, suitable for all implant systems



Security-Lock - the screw sits in the secondary part and the prosthetic restoration is bolted in the abutment



Friction Splint – fixation of the prosthetic restoration without a threaded hole



Individual screwing – conical titanium screw with the corresponding toolkit

#### Clinical case

Final restoration with milled NPM framework and veneered with the visio.lign system. 4 implants screwed transversally and 2 implants screwed occlusally. (Stefan Adler, Dental Technician, Landsberg)














## Assembled transverse fixation



For the straight and angled abutments of the SKY fast & fixed system, large and small bridges can be manufactured to a high aesthetic level with prosthetic caps for transverse screwing, as no screw channels impede the aesthetics. All framework materials can be used with this type of screwing – titanium, gold, NPM, ceramic, BioHPP. The "passive fit" of the bridge construction is ensured by oral bonding.



A highly secure and sealed connection is created using three-point fixation and the contact pressure generated by tight screwing of the prosthetic caps into the abutment platform.



The simplicity of use, particularly in the patient's mouth is ensured thanks to the variability of the screw position (360° in the case of straight abutments and 270° in the case of angled abutments), as access to the screw can always be positioned in an optimum manner. Moreover, the screw remains in the bridge so that protracted and difficult threading in the mouth is avoided. Transverse fixation can also be combined simply with diagonal screwing.





SKY fast & fixed Abutment 0° with integrated screw Height 1 mm REF SKYFT001 Height 2 mm REF SKYFT002 Height 4 mm REF SKYFT004



SKY fast & fixed Abutment 17,5° with screw 2,2 Height 3 mm REF SKYFT173 Height 5 mm REF SKYFT175



SKY fast & fixed Prosthetic coping Transversal screwretained REF SKYFTPKS



SKY fast & fixed Abutment 35° with screw 2,2 Height 4 mm REF SKYFT354 Height 5 mm REF SKYFT355

## Security-Lock



Patented threaded units which do not loosen or break. As the non-threaded section is in the inner coping, no micro-movements can be transferred. This guarantees that the threaded rod will not loosen inadvertently.

The matrix sleeve made of a high-melting cast-on alloy can be cast on up to max. 1300 °C. The threaded rods are available in three different sizes (1.0/1.4 and 1.8 mm) and suitable for any situation.



**HM-Centring drill** 1.0 REF 330 0081 5 REF 330 0066 0



Diatit-Multidrill 1.0 REF 330 0061 0 1.4 REF 330 0079 0 1.8 REF 330 0080 0



Threaded rods 2 pieces 1.0 REF 430 0729 3

1.4 REF 430 0729 4 REF 430 0729 5

Matrix sleeve with fixing screw 2 pieces each

REF 430 0729 6 REF 430 0729 7 1.8

REF 430 0729 8

#### Accessories:



Screwdriver short 1 piece REF 330 0069 0

Additional screwdrivers see pages 216-218.



Milling/drilling oil 20 ml, see page 219 REF 550 0000 8

#### Assortment

9 pieces Security-Lock 1.0

2 Threaded rods 2 Matrix sleeves 2 Fixing screws

1 Diatit-Multidrill

1 HM-Centring drill 1 Screwdriver short

REF 430 0729 0

Assortment 9 pieces

Security-Lock 1.4 2 Threaded rods 2 Matrix sleeves

2 Fixing screws

1 Diatit-Multidrill 1 HM-Centring drill

1 Screwdriver short

REF 430 0729 1

Assortment

9 pieces Security-Lock 1.8

2 Threaded rods 2 Matrix sleeves 2 Fixing screws

1 Diatit-Multidrill

1 HM-Centring drill 1 Screwdriver short

REF 430 0729 2



The different sizes can be used for a wide range of applications with implants and bridges etc.

A centring drill is used to

create a purchase point


in the correct position.



In this case, a superstructure is to be retained with a screw. The abutment should be waxed up using standard procedures.



The correct size of Multidrill is used to drill a hole at the correct angle for the threaded rod. It is absolutely essential that bredent milling/drilling oil is used.



been cast, it should be milled and polished.

Once the abutment has



Screw the threaded rod into the threaded sleeve. Both the pin and hexagonal socket (max. reduction: 2.3 mm) can be reduced as required.



Coat the threaded rod and sleeve with Pi-Ku-Plast, REF 540 0017 6.



Pi-Ku-Plast garantees optimum strength for continued processing.



Screw a retention screw coated with colloidal graphite into the threaded sleeve to retain it in the investment material. REF 540 0070 6.

#### **Dimensions**







Product	REF	Ø	Length	Thread	Length/Rod	max. Reduction
Threaded rod titanium 1.0	430 0729 3	Rod 1.0 mm	8.5 mm	M 2 x 0.4	3.5 mm	2.3 mm
Threaded rod titanium 1.4	430 0729 4	Rod 1.4 mm	8.5 mm	M 2 x 0.4	Rod 3.5 mm	2.3 mm
Threaded rod titanium 1.8	430 0729 5	Rod 1.8 mm	8.5 mm	M 2.5 x 0.45	Rod 3.5 mm	2.3 mm
Matrix sleeve HL 1.0	430 0729 6	2.8 mm	5.3 mm	-	-	2.3 mm
Matrix sleeve HL 1.4	430 0729 7	2.8 mm	5.3 mm	-	-	2.3 mm
Matrix sleeve HL 1.8	430 0729 8	3.2 mm	5.3 mm	-	-	2.3 mm



## Security-Lock-Ceramic



Security-Lock-Ceramic 1.4 allows splinting for all alloys without a thread sleeve.

Restorations made of a CoCr alloy which are to be veneered are biocompatible and can be produced without any additional alloy components.



Diatit-Multidrill 1.4 x 6 mm REF 330 0079 0



Threaded rods 1.4 2 pieces REF 430 0729 4







Ceramic removing tool REF 460 0010 6



10 pieces, 1 piece each Auxiliary modelling element Ceramic screw with wax sleeve HM-Centring drill Diatit-Multidrill Threaded rod 1.4

First tap, tungsten carbide Second tap, tungsten carbide Ceramic removing tool Tap handwheel Screwdriver, short

REF 430 0739 1



Auxiliary modelling element 1.4 REF 360 0116 9

Ceramic screws with

wax sleeve1.4

2 pieces REF 360 0117 0



Tap handwheel REF 330 0115 3





Screwdriver short 1 piece REF 330 0069 0



Milling/drilling oil see page 219 REF 550 0000 8



HM-Centring drill 1.4 REF 330 0066 0

The wax model of the

primary constractions is prepared in the usual

way.



Any alloy can be used for casting, even CoCr alloys.



After parallel milling, the secondary element is moulded using Pi-Ku-



To determine the exact position of the screw. the wax-up is modelled according to the situ-



The wax is removed at the specific point to determine the exact drilling position.



A groove is prepared at this point using the tungsten carbide centering drill.



The Diatit-Multidrill 1.4 and milling and drilling oil are used to prepare a hole in the desired direction of screwing.



The auxiliary modelling element is attached to the model using Pi-Ku-Plast and reduced with wax according to the situation.



The wax-up is reduced for ceramic veneering according to the situation.



Using tweezers, the auxiliary modelling element is turned slightly and removed.



After attaching the sprues, the ceramic spacer with wax sleeve is inserted into the opening up to the stop.



The wax sleeve and the model are connected.



The ceramic spacer remains in the metal structure until the ceramic veneer is completed.



The ceramic spacer is removed with the ceramic removing tool - do not remove with the sandblaster.



The thread is recut using the first and the second tap. Use milling and drilling oil when tapping.



The threaded rod is turned in and primary and secondary element are screwed with each other.



RFF

The threaded rod is shortened to the required length (max. reduction: 2.3 mm) using the Tita-Pol prepolishing wheel.



Security-Lock-Ceramic 1.4 can be processed safely and quickly with just a single alloy. There are no temperature-related alloy problems since finished parts are cast in.

**Dimensions** 



Product	
Threaded rod	titanium 1.4

430 0729 4

Ø Rod 1.4 mm

Length 8.5 mm

Thread M 2 x 0.4

Length/Rod Rod 3.5 mm

max. Reduction

2.3 mm

223

## Security-Lock adhesive sleeve



Security-Lock system is perfectly suitable for situations difficult to access such as small jaws or large-span bridges. The titanium threaded sleeve that can be glued in allows processing independent of the alloy.



Auxiliary modelling element 1.4 REF 360 0116 9




**HM-Centring drill** 1.4 mm REF 330 0066 0



Diatit-Multidrill 1.4 x 6 mm REF 330 0079 0



Threaded rods 1.4 2 pieces REF 430 0729 4





Tap handwheel REF 330 0115 3

#### Assortment

5 pieces, 1 piece each Auxiliary modelling element 1.4 HM-Centring drill 1.4

Diatit-Multidrill Threaded rod 1.4 Matrix sleeve titanium REF 430 0739 5

#### Accessories:



Screwdriver short 1 piece REF 330 0069 0

Additional screwdrivers see pages 216-218.



Milling/drilling oil see page 219 REF 550 0000 8



FGP insulating agent REF 540 0102 7



DTK-adhesive REF 540 0010 6

#### **Dimensions**





Product	REF	Ø	Length	Thread	Length/Rod	max. reduction
Threaded rod titanium 1.0	430 0729 3	Rod 1.0 mm	8.5 mm	M 2 x 0.4	3.5 mm	2.3 mm
Threaded rod titanium 1.4	430 0729 4	Rod 1.4 mm	8.5 mm	M 2 x 0.4	Rod 3.5 mm	2.3 mm
Threaded rod titanium 1.8	430 0729 5	Rod 1.8 mm	8.5 mm	M 2.5 x 0.45	Rod 3.5 mm	2.3 mm
Matrix sleeve titanium 1.4	430 0739 7	2.8 mm	5.3 mm	-	-	2.3 mm
Matrix sleeve HL 1.0	430 0729 6	2.8 mm	5.3 mm	-	-	2.3 mm
Matrix sleeve HL 1.4	430 0729 7	2.8 mm	5.3 mm	-	-	2.3 mm
Matrix sleeve HL 1.8	430 0729 8	3.2 mm	5.3 mm	-	_	2.3 mm

## Security-Lock adhesive sleeve



Any alloy is suitable for casting, even CoCr alloys.



After parallel milling and high luster polishing, the secondary element is moulded with Pi-Ku-Plast.



To determine the exact position of the screw, the wax up is modelled according to the situation.



The wax is removed at the specific point to determine the exact drilling position.



A groove is prepared at this point using the tungsten carbide centering drill 1.4.



The Diatit-Multidrill 1.4 and milling and drill oil are used to prepare a hole in the desired direction of screwing.



The auxiliary modelling element is attached to the model using Pi-Ku-Plast and moulded (completed) with wax according to the situation.



Prior to investing, the auxiliary modelling element is removed by turning it slightly with a pair of tweezers.



Any alloy can be used for casting the secondary construction.



All elements that must not be glued – such as the primary construction, the external surfaces of the primary construction and the screw are ...



... separated with FGP insulating liquid (REF 540 0102 7) so that the excessive adhesive can be removed easily.



After separating, the threaded rod is turned into the matrix sleeve.



Primary and secondary element are assembled. A drop of DTK adhesive is filled and spread evenly in the hole in the secondary element.



Matrix sleeve and threaded rod are inserted into the hole and no longer moved until the DTK adhesive has hardened.



The protruding matrix sleeve and the threaded rod are cut off to obtain the required length (max. reduction: 2.3 mm) using the Tita-Pol prepolishing wheel.



Ideal for processing highmelting alloys or titanium in very narrow jaws. The Security-Lock-adhesive sleeve ensures simple and quick processing.


## Friction Splint FS1

Connecting elements for superstructures.

- uncomplicated integration in the mouth
- defective screw connectors can be repaired with FS1
- FS1 is replaceable
- time-saving, no tapping necessary
- variable application for all indications
- can be individually shortened
- no loosening through expansion




Friction Splint FS1 sleeve 1 piece REF 450 0008 0 10 pieces REF 450 0008 4



Friction Splint FS1 pin 1 piece REF 450 0008 1

10 pieces REF 450 0008 5



REF 450 0008 7 Spacer Ø 2,0 mm 1 piece REF 450 0008 2 10 pieces

REF 450 0008 6


Modeling aid

REF 450 0008 3

Ø 2,0 mm

10 pieces

1 piece







The FS1 sleeve is pushed into the congruent splint hole running through the primary and secondary part with the pre-assembled splint screw.



The FS1 sleeve is pushed into the congruent splint hole running through the primary and secondary part with the pre-assembled splint screw.



...the splint screw is screwed in.



No tapping is necessary.



Damaged screw connectors....



....can be re-tooled with the FS1.

#### Accessories:



Diatit-Multidrill Ø 2,0 1 piece REF 330 0072 0



Fixing screw
2 pieces
REF 360 0103 0



Screwdriver short 1 piece REF 330 0069 0



Milling/drilling oil see page 219 REF 550 0000 8

Additional screwdrivers see pages 216-218.

## Friction Splint FS1



Wax-up with matrix.



Remove the wax-up. The pin hole is drilled into the abutment with the Diatit-Multidrill Ø 2,0 mm.



The wax-up is placed back onto the model. The modeling aid is integrated in the wax-up. Holes with a diameter of 2.0 mm are drilled into the full wax-up at the positions for the attachments.



The attachments are milled. The previously prepared matrix serves for orientation. The ceramic spacer can be used to ensure perfect casting of the splint holes.



Using the modeling aids...



... the secondary units are waxed up and prepared for casting.



Sleeve...



...and splint screw are shortened to the same length if required.



In the case of zirconium crowns it must be ensured that ...



...the diameter of the drillhole is 2.0 mm after the sintering process and ...



...the ceramic firings. Stress/tension within the ceramic can only be avoided in this way.



The splint screw which is screwed half way into the sleeve is positioned using tweezers...



...and pressed in. The remaining section of the splint screw is turned in.



The splint can be removed by turning it out with the screwdriver SW 0.9.



The screwed-in fixing screw is removed from the sleeve.



In case of usage of less than 1 year and in undamaged condition, the removed sleeve can be reinserted.

## Tool set for individual screw connections 1.4 and 1.6



For any situations and possibilities of dental technical screw connections.

Fast, inexpensive and tension-free screw connections.



Perfectly suitable for two-section bridges and dentures that are removeable to a limited degree.



Available in two different thread sizes.



Assortment

10 pieces Tool set for individual screw connections M 1.4 REF 330 0060 0



Assortment

10 pieces Tool set for individual screw connections M 1.6 REF 330 0001 6

The screw head is lowered 0.3 mm deep into the primary element. This way maximum tensile strength and protection against acting shear stress are ensured.

The conical screw head provides a self-locking effect. It is not possible for the screw to loosen itself.

Individual screw connections must be prepared for all gold content alloys at the points dictated by the dental-technical conditions. This way new dental-technical indications are obtained.



**HM-Centring drill** Ø 1.4 1 piece for M 1.4 and M 1.6 REF 330 0066 0



First tap 1 piece each M 1.4 REF 330 0067 1 M 1.6 REF 330 0116 V



Titanium screw 1 piece M 1.4 x 0.3 REF 330 0070 0

Head length 2.5 mm

10 pieces M 1.4 x 0.3 REF 330 0071 0

1 piece M 1.6 x 0.35 REF 330 0116 0 Head length 2.5 mm

10 pieces



Diatit-Multidrill 1 piece each M 1.4 REF 330 0063 0 M 1.6 REF 330 0115 7



Second tap 1 piece each M 1.4

REF 330 0067 0 REF 330 0116 F

Auxiliary modelling

element

M 1.4

1 piece each



M 1.6 x 0.35

REF 330 0116 1 Titanium screw

extended head

1 piece M 1.4 x 0.3



Diatit-Multidrill with stop 1 piece each M 1.4 REF 330 0075 0

M 1.6 REF 330 0115 8

Facing cutter

1 piece each



M 1.6 REF 330 0116 3

REF 330 0115 6

Screwdriver short 1 piece REF 330 0069 0

Additional screwdrivers see pages 216-218.




REF 330 0115 9



Accessories:

Milling/drilling oil see page 219 REF 550 0000 8





Tap holder 1 piece REF 330 0068 0



Head length 3.5 mm 10 pieces each M 1.4 x 0.3 REF 330 0K71 0

REF 330 0K70 0

1 piece M 1.6 x 0.35 REF 330 K116 0 Head length 3.5 mm

10 pieces M 1.6 x 0.35 REF 330 K116 1 Head length 3.5 mm



## Tool set for individual screw connections 1.4 and 1.6

## Two possibilities for a successfull screw connection

The quick screw connection without milling machine, only with the handpiece



The patrix of the bridgesectioning attachment features the same direction of insertion as residual abutment teeth.



Wax-up the second bridge element, cast and finish.



Prepare a small groove at the point where the screw is to be placed.



Drill through the secondary element approx. 1.5 mm deep into the primary element using the Diatit-Multidrill.



Remove the secondary element and drill into the primary element up to the stop using the Diatit-Multidrill with stop.



Assemble the primary and secondary element and drill up to the stop using the tungsten carbide facing cutter.



Cut the thread into the primary element. First use the first tap and then the last tap.



Assemble primary and secondary element and turn into the screw.



The screw head with the secondary element is ground flush and polished.

#### Use of the auxiliary modelling element

The safe method once the direction of the screw has been determined



Grind a small groove into the patrix using the center drill.



The Diatit-Multidrill drills down to the exapt depth.



Integrate the auxiliary modelling element into the pattern using the brush resin.



Complete the pattern using modelling wax.



Turn the auxiliary modelling element with apair of pliers and remove it.



After casting, assemble the bridge elements. Drill to the stop using the facing cutter. Further working steps are described in figures 7, 8 and 9.

**Dimensions** 





Product	REF	Ø	Length	Thread	Length/Head	max. Reduction
Titanium screw M 1.4	330 0070 0	2.1 mm	4.5 mm	M 1.4 x 0.3	2.5 mm	1.2 mm
Titanium screw M 1.4 / 3.5	330 0K70 0	2.3 mm	5.5 mm	M 1.4 x 0.3	3.5 mm	1.8 mm
Titanium screw M 1.6	330 0116 0	2.3 mm	5.2 mm	M 1.6 x 0.35	2.5 mm	1.2 mm
Titanium screw M 1.6 / 3.5	330 K116 0	2.6 mm	6.2 mm	M 1.6 x 0.35	3.5 mm	2.0 mm

## Tool set for individual screw connections Additional Set zirconium



#### Assortment

3 pieces, 1 piece each Diatit-Multidrill Facing cutter zirconium Positioning pin

REF 330 2432 4

The preparation of screw connections in zirconium restorations is simplified in conjunction with the tool set for individual screw connections



The tools that are 30 % larger compensate the shrinkage of zirconium and allow precision-fit screw connections



Diatit-Multidrill 1.5 x 8 mm REF 330 0073 0



Facing cutter zirconium REF 330 2432 6



Positioning pin REF 330 2432 7



## Bridge-sectioning Attachment oc



This prefabricated unit facilitates fabrication of a sectioned bridge with occlusal screw.

The titanium screw has a hexagonal socket to facilitate tightening and loosening it.

The circumferential ring marks the maximum level to which it can be shortened.

Made of cast-on alloy.



Fixation screw 2 pieces REF 360 0103 0



Titanium screw 1 piece REF 330 0070 0 10 pieces REF 330 0071 0



Closing ring HL, cast-on 2 pieces REF 430 0730 4



**Bridge-sectioning** studs oc 2 pieces REF 430 0730 3

#### Assortment

6 pieces, 1 piece each Titanium screw Closing ring HL, cast-on Screwdriver short Bridge sectioning studs oc

Fixation screw Paralleling mandrel REF 430 0730 2

Accessories:



Screwdriver short 1 piece REF 330 0069 0

Additional screwdrivers see pages 216-218.



Paralleling mandrel for oc and custom bridge-sectioning attachments 1 piece REF 360 0115 7



The paralleling mandrel positions the sectioning attachment correctly.



The design and minimal dimensions of the threaded sleeve in the sectioning attachment enable it to be adapted to the papillae as required.



The section connecting the attachment to the coping is rounded, has a diameter of 1.0 mm\* and can be trimmed accurately with a 1.0 mm rotary cutter if required.



The exterior design of the locking ring, which consists of a cast-on gold alloy, ensures that it is retained securely in resin.



The threaded sleeve is made of a cast-on alloy and can be used with any gold or semi-precious



on the locking ring marks the level to which the screw and locking ring can be reduced.

outer section.



The fixing screw which is coated with colloidal graphite retains the threaded sleeve precisely in the investment mate-



To ensure that the locking ring is fixed in place securely, the outer section must be moulded with Pi-Ku-Plast brushon resin.



The titanium screw can be around to blend it into the occlusal surface.

Dimensions




Product	REF	Ø	Length	Thread	Length/Head	max. Reduction
Titanium screw 1.4	330 0070 0	2.1 mm	4.5 mm	M 1.4 x 0.3	2.5 mm	1.4 mm
Closing ring HL, cast-on HL	430 0730 4	2.5 mm	2.1 mm	-	-	1.4 mm
Bridge-sectioning studs oc	430 0730 3	3.0 mm	6.9 mm	M 1.4 x 0.3	-	3.3 mm



## Custom Bridge-sectioning Attachment







Custom bridge-sectioning attachments 8 pieces REF 430 0735 0

#### Accessories:



Tool set 10 pieces REF 330 0060 0



Paralleling mandrel for oc and custom bridge-sectioning attachments 1 piece REF 360 0115 7



Milling/drilling oil see page 219 REF 550 0000 8



A paralleling mandrel is used to position the bridge-sectioning attachment as required for the case.



The plastic component can be adapted to the papillae as required.



The section connecting the attachment to the coping is rounded, has a diameter of 1.0 mm\* and can be trimmed with a cylindrical cutter (size 010) if required.



The purchase point is created with a centring drill.



Bredent milling/drilling oil should be used when drilling. All other oils, especially etheric oils, are unsuitable and impede correct drilling.



A Multidrill (1.2 x 5) from the tool set is used to drill an approximately 2 mm deep hole. The use of generous amounts of drilling oil prevents the drill overheating.



Use a stop drill (1.2 x 2) to drill the threaded hole precisely to the required depth. Use Bredent drilling oil to ensure that the hole is drilled neatly and smoothly.



A countersinking drill is used to widen the hole to 1.4 mm for the thread tap and create space for the conical screw head.



The pre-tap taps the first stage of the thread. The final tap taps a high precision thread. Drilling oil prevents the tap



The conical screwhead fits into the inner section by approximately 3/10 mm. It withstands higher shear forces (155 kg) than conventional systems.



The screw should be coated with Pi-Ku-Plast resin and integrated into the pattern. The screw should be reduced after casting.



The minimal dimensions of the screw provide pleasant aesthetics for all screw-retained restorations.


ni	m	01	٠ci	^	n	•
וט	ш	CI	131	U	11.	3



Product	REF	Ø	Length	Thread	Length/Rod Head length	max. Reduction
Custom Bridge-sectioning Attachment	430 0735 0	3.0 mm	7.0 mm	_	_	custom

## Prefabricated screwing set

For occlusal and horizontal screw connections.






Titanium screw
1 piece
REF 330 0070 0
10 pieces
REF 330 0071 0



Closing ring HL, cast-on 2 pieces REF 430 0730 4



Tapped bush HL 2 piece REF 330 0081 1



Closing ring HL, cast-on






Fixation screw 2 pieces REF 360 0103 0



Screwdriver short 1 piece REF 330 0069 0 Additional screwdrivers see pages 216-218.



Tapped bush HL



#### Assortment

5 pieces, 1 piece each Titanium screw Closing ring HL, cast-on Tapped bush HL Fixation screw M 1.4 Screwdriver short

REF 430 0735 1

Dimensions



Product	REF	Ø	Length	Thread	Length/Head	max. Reduction
Titanium screw 1.4	330 0070 0	2.1 mm	4.5 mm	M 1.4 x 0.3	2.5 mm	1.4 mm
Closing ring HL, cast-on HL	430 0730 4	2.5 mm	2.1 mm	_	-	1.4 mm
Tapped bush HL	330 0081 1	2.3 mm	3.5 mm	-	-	_

Harmonized products in the casting technique system offer unsurpassed precision for high-quality dental restorations, which is ensured and can be reproduced at any time through the use of products for the casting technique developed and manufactured by bredent.



The primary structure is fabricated as described in chapter 5. Brevest C+B Speed is used to ensure accurate fit of the restoration.




The block-out and saddle waxes feature exceptional modelling properties for accurate exposure. This way the time required for fitting the secondary structure is reduced.



The duplicating system with the duplicating silicone is an essential component of the system. Accurate duplicating is the key to precision of fit of the CoCr structure.



The investment materials that can be exactly controlled reduce the amount of time for fitting and contribute to fabricating precision-fit castings.



The system of sprues ensures of homogeneous structure of the casting and hence a compatible substructure. Additionally, the dense surface facilitiates polishing.



Surfaces can be finished more quickly through the use of alloys that are easy to process.

nccis of wax pattern	23
Quadro wax profile	23
Wax patterns cut to size wpz	
Protek wax adhesive wk 2	
Pi-Ku-Plast / Pi-Ku-Plast HP 36	
Laser joint	
DTK-adhesive	
Double-T Adhesive Connector	
Double-T Adhesive Mini Connector dtk	
Optiguss	24
Duplicating	
Exaktosil N 15 / N 21	24
Technosil NT duplicating silicone	25
Technolit	
Duplicating system	
Isosil	
Master-Copy	
Bre-Gel 1	
Bre-Gel 2 opaque, Bre-Gel 3 opaque-liquid	25
Investing and casting	
Microkeramik	
Brevest M1	
Brevest Rapid 1	
Brevest exakta M and Brevest exakta Speed	
Brevest ESG	25
Brealloy F 400	25
Brealloy MO	
Brealloy solder	
Brealloy flux	
Duro-Top	
Investment hardener	
Crepe sleeve	
Investment marker	
Sprue	
Golden booklet	26
Implementing the repeatable solution in your	
lab	26
Milling	
Milling unit BF 2	
Milling base	26
Model support BF 1	
Transfer device	
Brenometer surveying system	
5. Cristician surveying system	2 0
Restoring the friction	

Statik-Disc......234

Protek sculpturing wax......235

Clasp patterns......237

Retentions.....

Wax sheets.....

Assortment box .....

Modelling


Activating pliers.....

Friction fit system FGP......265

#### Statik-Disc


Time-consuming design drawings when planning CoCr dentures are no longer required if the Statik-Disc is used.

The statically correct position of the supporting elements is quickly determined by dentists and dental technicians.



- quick determination of correct static
- can be applied individually to any situation
- suitable for all models
- simple handling


The clasp line principle The clasp line principle applies to all denture constructions. The clasp supporting line (a) runs



Statik-Disc REF 360 0126 7 peripheral to the denture body. In the case of saddle dentures it should lie in the center of the jaw ridge. It is always obtained by connecting the clasp supports. The main clasp line (b) separates the jaw halves diagonally. It results from connecting the supports of tooth no. 13 and 27.

## Determination of the tilting axis

Work arm and power arm are vertical on the tilting axis. The lever principle



applies: load (L) x work arm (W) = power (P) x power arm (p). Load and power are given; therefore it must be attempted to keep power x power arm on the same level or higher than load x work arm.

#### Kennedy class I

On both sides the gaps are in the distal area of the residual dentition (bilateral free-end dentures). This type of denture creates the following static situation:

If a saddle is lowered after exposure to masticatory

pressure, diagonal tensile stress is obtained on the opposite side. The rotation axis runs through the support on the same side and the end of the saddle on the opposite Mark planned position of the tooth set up last on the modeled diagonal to the final natural tooth on the shorter arch.
(1). The Statik Disc is placed on the model so that the red line runs through the center of the planned support between the

planned tooth and the final natural tooth on the opposite side. The red line is also the tilting axis. The black line is turned to the planned support of the final natural tooth on the opposite side to determine where the final articifial tooth needs to be set up (2). Simultaneously the blue lines allow to read the power/load relationships. The blue lines have a distance



of 10 mm to allow rapid and simple determina-

In this example it can be recognized that the power/work arm relationship is not perfect when integrating the final tooth. Accordingly, this tooth should not be replaced and thus the arch be shortened (3).



#### Kennedy class II

The gap is on one side in the distal area from the residual dentition (unilateral free end denture) or – in conjunction with a larger gap on the other side. This type of denture cre-



ates the following static situation:


Unless clasps are correctly attached to the denture, incorrect loading and tilting of the denture may be caused. Therefore a supporting element must be attached to avoid tilting .

The Statik-Disc is placed on the model so that the red line is in the center between the final tooth at the shortened arch (mesial support) and the last tooth on the opposite arch



(mesial support)

The black line is turned so that it points to the desired tooth that was set up last. The position of the anti-tilting element is now indicated on the other side. If this element is situated



mainly in the aesthetic area, the black line must be turned further to the mesial direction towards the tooth that is set up next. The position of the the antitilting element will then be shifted to the distal area.



Kennedy class IV defines gaps that are located left or right from the central line and are limited by the residual dentition in the distal area.



In this type of dentures, clasps are attached to the dorsal area. If the Statik-Disc is placed on so that the red line



serves as tilting axis, the blue lines allow to recognized immediately that the power/work arm relationship can be balanced by the clasps attached to the dorsal area. In this case the power arm is missing so that a long work arm is obtained. Open clasps with distal support must be



used. The clasp arms act as retentions in the presence of tensile stress since they are held by the equator if exposed to pull-off movement.

#### Attachment techniques



The Statik-Disc also simplifies correct planning if attachments are used. In this case the red line is placed on the attachment; this line also serves as tilting axis. The power/work arm relationship is read with the help of the blue line and thus the expansion of the teeth to be set up can be determined.



## Protek sculpturing wax



Protek sculpturing wax - emphasizes the contrast for improved viewing and adjusting.

Protek sculpturing wax 25 g, green REF 510 0090 1



The sculpturing wax has the same consistency as all Protek components, which enables junctures to be waxed-up effortlessly and harmoniously. It is no longer necessary to carve from hard into soft wax.

## Biotec blocking out wax



Biotec blocking out wax 28 g, pink REF 510 0061 5



The special components of the blocking out wax ensure perfect blocking out of undercuts.



Easy and quick scraping allows to save time.




No color additives penetrate into the plaster surface after boiling out the model. The master model remains clean.


#### entire field of CoCr work. Blocking out wax with very good scraping properties. No discoloration on the plaster model after boiling out.

Special wax for blocking out undercuts in the

## Protek saddle wax with pre-formed border



Protek saddle wax Size A 0,40 REF 430 \*571 0 90 pieces each, right/left 0.60 REF 430 \*573 0 90 pieces each, right/left



Protek saddle wax Size B 0,40 REF 430 \*572 0 105 pieces each, right/left 0,60 REF 430 \*574 0 105 pieces each, right/left





This self-adhesive saddle wax with pre-formed border, available in 2 sizes and thicknesses, guarantees that the underside of the acrylic is absolutely precise and even. The border matches the lingual bar joiners exactly.

\* Also available as Protek summer wax (ideal wax quality for higher temperatures). Just replace the asterisk in the REF with an "S" for summer wax or with "O" for standard wax quality.

## Protek - Spacer wax



Spacer wax 75x150 mm plaques 15 pieces emp.

0,30 mm REF 430 \*582 0 self-adhesive: 0,40 mm REF 430 \*583 0

0,30 mm REF 430 \*586 0 0,50 mm REF 430 \*584 0 0,40 mm REF 430 \*587 0 0,60 mm REF 430 \*585 0 0,50 mm REF 430 \*588 0 0,60 mm REF 430 \*589 0



The quality of Protek spacer wax is better than ever before - extremely ductile and tear-resistant. After duplicating, it can be off the model without leaving a residue. Simplifies preparation of the model for duplicating and saves a great deal of time.

\* Also available as Protek summer wax (ideal wax quality for higher temperatures). Just replace the asterisk in the REF with an "S" for summer wax or with "O" for standard wax quality.



## Lingual bar patterns



Protek lingual bar wax pattern Ergonomically shaped lingula bar pattern. Adapatation is simplified by the concave

shape matched with the jaw; hence time is saved during finishing.

#### Protek lingual bar wax pattern

REF 430 0743 0 3.6 x 1.85 30 pieces 80 pieces REF 430 0748 0



The conventional bar pattern is more difficult to adapt; a wax knife must be used for coating with wax.




The Protek lingual bar wax pattern adapts to the gingival situation so that reshaping with the wax knife can be omitted.



	1.7 x 4	30 pc. 80 pc.	REF 430 0124 C REF 430 0125 C
	2.0 x 4	30 pc. 80 pc.	REF 430 0124 B REF 430 0125 B
	2.3 x 4	30 pc. 80 pc.	REF 430 0124 A REF 430 0125 A
	2.45 x 4.3	30 pc. 80 pc.	REF 430 012A 0 REF 430 013A 0
Ass	ortment com	REF 430 0124 6	

Protek wax bars are available in 3 sizes to fit every type of jaw. Thanks to the structure of the wax they are easily adapted and exhibit no elastic recovery. The high pressure-resistance of this wax prevents deformation of the patterns and ensures that the framework pattern is shaped aesthetically.

## Protek - Lingual bar joiners



1.7 le A REF 430 517 LA

1.7 ri A REF 430 517 RA

2.0 le A REF 430 520 LA

2.0 ri A REF 430 520 RA

2.3 le A REF 430 523 LA

2.3 ri A REF 430 523 RA

- save work

le B

ri B

le B

le B

1.7

2.0 ri B

2.3

2.3 ri B REF 430 517 LB

REF 430 517 RB


REF 430 520 LB

REF 430 520 RB

REF 430 523 LB

REF 430 523 RB

12 bars of each



50 pc.

50 pc.

50 pc.

50 pc.

50 pc.

50 pc.

15 pc. each size A+B ri/le REF 430 0575 0 2.0 REF 430 0576 0 2.3 REF 430 0577 0

Assortment:



The lower edge of the Protek joiner is shaped to fit the investment finishing line exactly, which was duplicated with Protek saddle wax.



Protek lingual bar joiners are available to fit every size of bar. The joiner is fitted into place and waxed onto the bar thus ruling out the need to wax-up the joint, which is very time-consuming. Protek lingual bar joiners are available in 2 sizes and matched to the three thicknesses of Protek lingual bar, to suit any situation.

## Protek - Clasp / bar joiners



Protek - Clasp / bar joiner size A, 100 pc. each REF 430 0578 0



Protek - Clasp / bar joiner size B, 100 pc. each REF 430 0579 0



Interdental joinments of Protek clasps and bars effortlessly. The Protek joiner is matched to the relevant bar exactly and only requires placing in



lateral view





If Opticast is brushed onto the pattern, the components do not have to be waxed together and even finest of cracks are sealed. The investment material cannot creep under the pattern.



## Clasp patterns



The bent premolar clasp features a shifted ridge to ensure that the chyme is carefully transported away via the papilla. Compared to the previous cross-section of the clasp this results in considerable protection of the periodontium.



When using this pre-bent clasp pattern, no compression or elongation will result whilst bending the wax. This way the casting of the clasp profile will be more homogeneous.

Premolar clasps, bent wlf pmk 10 trays REF 430 0748 1

Premolar clasp, bent, for resin injection moulding

10 sheets of 20 clasps, 10 left + 10 right **REF 430 0748 5** 

Visit the course **"CoCr work is cast information".** Please request the course program!







The equator is marked in the usual way. The position is determined using the surveying plate 2 (undercut depth of 0.35 mm) of the Brenometer surveying system. If a short clasp is used (8 mm), the clasp tip is placed over the determined point (figure 1).







In the case of a premolar clasp (11 mm) the clasp tip is placed on the point (figure 2) and in case of a molar clasp (14 mm) the tip is placed below the point (figure 3). According to the determined position, the bent premolar clasps are placed against the tooth. Do not attach the patterns with wax to avoid changing the pattern.







During finishing, the cast clasps are only smoothed using a rubber polisher and the clasp tip is rounded off. This way the corresponding shape and the length of the pattern will always allow to obtain the same draw-off strength.



The clasp profile is separated from the central rod ("tree") using a knife.



Place the clasp pattern on the template and cut off the desired length.



Accessories:

Wax adapter REF 360 0120 5



The pre-bent clasp pattern is placed against the corresponding tooth and fixed using the wax adapter, REF 360 0120 5. Do not attach with wax to avoid changing the pattern.



Thanks to the ideal clasp pattern design the chyme is directed away from the tooth and the gingiva is protected.



Molar clasp 10 sheets 20 clasps each

REF 430 0157 1





Bonyhard clasp 10 sheets, 12 clasps each REF 430 0157 6





## Clasp patterns



Circumferential clasp, bent

10 sheets of 20 clasps

REF 430 0157 2

Thanks to their resistance to deformation and pressure, all Protek patterns can be fitted easily and quickly. All Protek components stay in place exactly when bent, which obviates the need to wax them into place and, especially in the case of clasp tips, avoids grinding which would otherwise



## Retentions

A wide range of retentions for any situation. Special wax offering high elasticity of bending for precise working.

Perforated retainers 25 pieces 13.5 cm long REF 430 0159 0





Perforated retainers staggered 25 pieces 13.5 cm long REF 430 0159 1





Comb-shaped retainers 25 pieces

13.5 cm long REF 430 0157 5





Finishing bands wral 20 pieces REF 430 0157 7





Perforated mesh 1.5 / 2.0 20 pieces 7x7 cm sheets 1.5 REF 430 0599 0


2.0 REF 430 0158 3





Reinforcing mesh upper, preformed 12 pieces REF 430 0219 0







#### Wax sheets

#### Stippled wax sheets

transparent so that markings and areas to be reduced remain visible.















#### fine stippling plw green

REF 430 \*161 0 0.30 mm 15 pieces 0.35 mm 15 pieces REF 430 \*161 1 0.40 mm 15 pieces REF 430 \*161 2 0.45 mm 15 pieces REF 430 \*161 3 REF 430 \*161 4 0.50 mm 15 pieces 0.60 mm 15 pieces REF 430 \*161 5

### medium stippling plw

0.30 mm 15 pieces REF 430 \*161 6 0.35 mm 15 pieces REF 430 \*161 7 0.40 mm 15 pieces REF 430 \*161 8 REF 430 \*161 9 0.45 mm 15 pieces 0.50 mm 15 pieces REF 430 \*162 0 0.60 mm 15 pieces REF 430 \*162 1

#### coarse stippling plw

REF 430 \*162 2 0.30 mm 15 pieces REF 430 \*162 3 0.35 mm 15 pieces 0.40 mm 15 pieces REF 430 \*162 4 0.45 mm 15 pieces REF 430 \*162 5 0.50 mm 15 pieces REF 430 \*162 6 REF 430 \*162 7 0.60 mm 15 pieces



#### Assortment box



The Protek assortment box provides a clear overview and simplifies working procedures. Can be filled according to your wishes.

Protek assortment box E 12 (empty box without content) REF 640 0084 0

All Protek patterns are available in refill packs which can be fully recycled and are harmless to the environment

## Reels of wax pattern



Various diameters of wax pattern are available in medium and hard consistencies.

Reels of wax	REF	REF		
pattern, 250 g	blue (medium)	green (hard)		
Cross-section in mm				
• 1.2	430 0115 0			
• 1.5	430 0115 5			
2.0	430 0116 0	430 0111 0		
2.5	430 0116 5	430 0111 5		
3.0	430 0117 0	430 0112 0		
3.5	430 0117 5	430 0112 5		
4.0	430 0118 0	430 0113 0		
5.0	430 0118 5	430 0113 5		



The wax patterns can be bent without recovering elastically or becoming pinched.



## 7.0 x 1.5 x 180 mm

REF 430 0156 0 220 g

## Quadro wax profile



#### Square sprues for better casting results.

Studies have shown that all liquids - including liquid metal - flow in drops; that also applies to flowing into a square sprue.

Accordingly, the gas (air) contained in the cavity (casting mould) can escape freely across the unfilled corners. Results:

- no swirling of molten metal due to the back pressure of the residual air
- faster flowing in of the molten metal
- more homogeneous castings
- smoother surfaces
- increased precision of fit



#### Quadro wax profile 250 g, green

- 1.75 x 1.75
- REF 430 0691 0
- 2.25 x 2.25
- REF 430 0692 0
- 3.00 x 3.00 REF 430 0693 0

## Wax patterns cut to size wpz

Cross-section in mm, green					
•	0.8	REF 430 0125 0			
•	1.2	REF 430 0121 0			
•	1.5	REF 430 0121 5			
•	2.0	REF 430 0122 0			
_	1.8 x 0.9	REF 430 0122 5			
•	2.0 x 1.0	REF 430 0123 0			
	3.0 x 1.5	REF 430 0123 5			
	4.0 x 1.5	REF 430 0124 0			
	4.0 x 1.7	REF 430 0124 5			

Wax pattern assortment: 150 g Size 1.2 mm and above wax patterns, cut to size REF 430 0120 0



An assortment of round and semi-round wax patterns in high Protek quality - resistant to deformation and pressure, no elastic recovery which facilitates the attachment of retainers. All patterns are available separately in 55 g packs.

## Protek wax adhesive wk 2 - soaks into the investment material



Wax adhesive wk 2 20 ml REF 540 0099 0 100 ml REF 540 0100 2 Thinner 100 ml REF 540 0100 1



Protek wax adhesive can be applied to the model in a thin film and soaks into the investment material.



The patterns stick securely to the investment model, with no marginal gap whatsoever.

See "Sprues", chapter 5

## Pi-Ku-Plast / Pi-Ku-Plast HP 36



Exceptional material properties such as perfect contouring characteristics and no slumping provide the precondition for top-quality casting results. The brush resin is available in five different colors. Both resins differ only in their contraction. HP 36 has a contraction value of just 0.036 %. Since the resin sets quickly, it is perfectly suitable for the fabrication of resin dies or resin copings in the double crown technique.



The micro-fine grain size allows to reproduce all details and increases the precision.

Advantages of Pi-Ku-Plast HP 36 Five translucent colors simplify control of layer thicknesses so that reworking is minimized.



The flat, pointed shape of the brush which is available in two different sizes allows to take up the desired quantity of material and reduces material consumption.





Wet the brush with Pi-Ku-Plast HP 36 monomer. The amount and firmness of the Pi-Ku-Plast portion can be controlled by the amount of monomer and the time it is immersed in the polymer.

# Assortments big Pi-Ku-Plast

3 vessels 1 brush each, size A + B 1 brush holder

100 ml cleaner 100 ml monomer 85 g polymer blue
 yellow
 orange
 red
 transparent
 REF 540 0017 3
 REF 540 0017 4
 REF 540 0017 5
 REF 540 0017 6
 REF 540 0017 7

# Assortments Pi-Ku-Plast HP 36

3 vessels
1 brush each, size A + B
1 brush holder
100 ml cleaner

1 brush holder 100 ml cleaner 100 ml monomer 85 g polymer blue
 yellow
 orange
 red
 transparent
 REF 540 0021 7
 REF 540 0021 8
 REF 540 0022 0
 REF 540 0021 6

#### Refill package

100 ml	cleaner		REF 540 0016 9
85 g	polymer		REF 540 0016 7
100 ml	monomer	<ul><li>blue</li></ul>	REF 540 0016 8
		<ul><li>yellow</li></ul>	REF 540 0017 8
		orange	REF 540 0017 9
		red	REF 540 0018 0
		<ul><li>transparent</li></ul>	REF 540 0018 1

#### Refill package

100 ml	cleaner		REF 540 0022 4
85 g	polymer		REF 540 0021 5
100 ml	monomer	blue	REF 540 0021 3
		yellow	REF 540 0021 1
		orange	REF 540 0021 2
		red	REF 540 0021 4
		transparent	REF 540 0021 0

#### Refill package

vessel for cleaner, 8 ml	REF 540 0017 2
vessel for monomer, 8 ml	REF 540 0017 1
vessel for polymer, 8 ml	REF 540 0017 0
brush size A + brush holder, 3 pieces	REF 330 0114 6
brush size B + brush holder, 3 pieces	REF 330 0114 7

#### Refill package

vessel for cleaner, 8 ml	REF 540 0020 9
vessel for monomer, 8 ml	REF 540 0020 7
vessel for polymer, 8 ml	REF 540 0020 8
brush size A + brush holder, 3 pieces	REF 330 0114 6
brush size B + brush holder. 3 pieces	REF 330 0114 7

# Assortment small Pi-Ku-Plast

12 g polymer

20 ml cleaner
2 modelling dishes silicone, red
20 ml monomer red
1 brush size B and brush holder

REF 540 0019 6



Pi-Ku-Plast separating agent 10 ml REF 540 0018 2

## Pi-Ku-Plast / Pi-Ku-Plast HP 36



brush size A + holder REF 330 0114 6

#### Optimal control of layer thickness thanks to the transparent colors of Pi-Ku-Plast



brush size B + holder REF 330 0114 7



Dip the flat side of the brush into the polymer to take up large portions.



Dip the narrow side of the brush into the



polymer to take up small portions.

Dip only the brush tip

up very small portions.



The high-lustrous reproduction of the metal surface of the primary element results in a perfect inner surface of the secondary element and thus allows to save precious working time.

Gap-free fit of the outer coping for unsurpassed precision of the cast secondary elements.




Pi-Ku-Plast separating agent, REF 540 00182, allows to produce stable primary elements directly on the plaster die and provides a convincing alternative to wax.



Wax and metal can be connected rigidly using Pi-Ku-Plast HP 36 which renders the material universally suitable.

#### The incineration phase of the resin elements in the casting ring frequently determines whether dental castings could be produced successfully or not.



The competitor's resin and Pi-Ku-Plast HP 36 in the incineration test.



At 275°C the competitor's product foams and expands considerably.



At 300°C the competitor's product reveals distinctive expansion whereas Pi-Ku-Plast HP 36 reduces the volume.



Identical copings produced with brush resin.



The competitor's resin and Pi-Ku-Plast HP 36, prepared for investing.



The considerable expansion of the competitor's resin during the incineration phase resulted in the fracture of the investment material die in the casting ring. After casting, the crown is sealed with a lid and can not be used. A section through the cast crown (figure 8) shows the fractured die.





Owing to friction heat, the competitor's resin may reach the plasticity phase which may result in deformation of the model and require considerable reworking.



Pi-Ku-Plast HP 36 is insensitive to heat, retains dimensional stability and ensures precision of fit which is well above the standard.




Allow the wax element to cool down to obtain a tension-free bridge model, separate using a thin blade and connect using Pi-Ku-Plast HP 36.



The extremely low shrinkage of Pi-Ku-Plast HP 36 allows to obtain a tension-free model and a precision-fit casting.



## Laser joint



#### Assortment

- 30 pieces
- 10 Inner housings 10 Outer duplicating
- housings 10 Outer wax-up housings

REF 440 0000 4

Laser weld joints rationally and precisely. The LV 1 laser joint ensures that the joint is always of the correct size, fits precisely and can be fabricated

quickly.



Wax the inner housing of the laser joint onto the outer housing - It only has to be paralleled if the outer housing is to be welded at two spots. Please note: The approximal "collar" should always face the occlusal aspect.



Before duplicating, place the red outer duplicating housing on the inner housing of the laser joint.



Before casting the investment model, replace the red outer duplicating housing with a blue outer wax-up housing. The outer duplicating housing is red - the outer wax-up housing is blue.



Refill packages:

Inner housing

Measurements in mm:

Outer duplicating housing

Outer wax-up housing

16 pieces	REF 440 0000 5
50 pieces	REF 440 0000 1
16 pieces	REF 440 0000 6
50 pieces	REF 440 0000 2
16 pieces	REF 440 0000 7
50 pieces	REF 440 0000 3
	50 pieces 16 pieces 50 pieces 16 pieces

L 4.6 x W 1.6 x H 2.5

L 5.2 x W 2.6 x H 2.5

L 5.2 x W 2.6 x H 2.5



Paralleling mandrel universal REF 360 0115 1

Custom laser joints are complicated and time-consuming to fabricate. To achieve precision of fit and high strength, the joint must be made to precise dimensions. The LV 1 laser joint ensures that the weld is strong and accurate. The outer housing cannot move due to contraction of the weld seam.



Shows the investment model with a outer wax-up housing: The chrome cobalt framework should be waxed up as usual. The interior dimensions of the outer wax-up housing are slightly larger than those of the outer duplicating housing. Therefore, the cast chrome cobalt framework fits the inner housing without requiring adjusting. The retention grooves can also be used to check the position of the outer housing.



Before welding the outer housing, remove the occlusal bar from the LV 1 laser joint. The outer housing should be fixed in place with two spot-welds placed diagonally above and two beneath the joint. The precision of fit should then be checked. The entire joint should then be welded, placing the welds diagonally.



The outer wax-up housing is minimally oversized. This ensures that the joint fits precisely after welding. If several outer housings are to be welded, proceed consecutively - always weld one joint properly, check the precision of fit and then fix the next housing in place.

## DTK-adhesive



Dual-hardening composite adhesive for the fixation of dental attachment elements.

DTK-adhesive REF 540 0010 6



#### Accessories:

5 g REF 540 0111 K Catalyst paste K, 5 g REF 540 0111 B Base paste B, Mixing block, 10 pieces REF 330 0114 4 Spatula, 100 pieces REF 330 0114 3



## Double-T Adhesive Connector



Tension-free, precise and low-cost metal junctures can be prepared very quickly with only very little space required.



Matrix

**Patrix** 



Patrix with Matrix



**Duplicating matrix** 



Patrix with duplicating matrix

#### Different inclination angles and sizes for all jaw situations

III. 1:1	90° A		90° B		120° A		120° B
Patrix	L 6.0 mm W 4.0 mm H 3.5 mm	id.	L 4.5 mm W 2.5 mm H 3.5 mm	4	L 7.0 mm W 4.0 mm H 5.0 mm	4	L 5.0 mm W 2.5 mm H 3.0 mm
16 pieces 50 pieces	REF 430 0405 A REF 430 0342 A		REF 430 0405 B REF 430 0342 B		REF 430 0402 A REF 430 0422 0		REF 430 0402 B REF 430 0423 0
Matrix	L 5.5 mm W 4.0 mm H 3.0 mm		L 4.0 mm W 2.5 mm H 2.0 mm		L 5.5 mm W 4.0 mm H 3.0 mm		L 4.5 mm W 2.5 mm H 2.5 mm
16 pieces 50 pieces	REF 430 0404 A REF 430 0341 A		REF 430 0404 B REF 430 0341 B		REF 430 0401 A REF 430 0420 0		REF 430 0401 B REF 430 0421 0
Duplicating matrix	L 5.5 mm W 4.0 mm H 3.0 mm		L 4.0 mm W 2.5 mm H 3.0 mm		L 5.5 mm W 4.0 mm H 3.0 mm		L 4.5 mm W 2.5 mm H 2.5 mm
16 pieces 50 pieces	REF 430 0406 A REF 430 0343 A		REF 430 0406 B REF 430 0343 B		REF 430 0403 A REF 430 0424 0		REF 430 0403 B REF 430 0425 0
Paralleling mandrel, 1 piece	REF 430 0345 A		REF 430 0345 B		REF 430 0344 A		REF 430 0344 B

#### Accessories:



DTK-adhesive REF 540 0010 6

#### Assortment

#### Double-T Adhesive Connectors dtk 90°

- 5 Patrices A
- Patrices B
- 10 Matrices A
- 10 Matrices B 5 Duplicating
- matrices A
- Duplicating matrices B
- Paralleling mandrel each, size A + B

REF 430 0340 0

#### Assortment

#### Double-T Adhesive Connectors dtk 90°

- 3 Patrices A
- Patrices B
- 6 Matrices A
- Matrices B
- Duplicating matrices A
- Duplicating matrices B

REF 430 0347 0

#### Assortment

#### Double-T Adhesive Connectors dtk 120°

- 5 Patrices A
- 5 Patrices B
- 10 Matrices A
- 10 Matrices B
- 5 Duplicating matrices A
- Duplicating matrices B
- Paralleling mandrel each, size A + B

REF 430 0408 0

#### Assortment

#### Double-T Adhesive Connectors dtk 120°

#### 3 Patrices A

- 3 Patrices B
- 6 Matrices A
- 6 Matrices B
- 3 Duplicating matrices A
- Duplicating matrices B

REF 430 0407 0



## Double-T Adhesive Connector

#### Double-T Adhesive Connectors at the crowns

Patrixes must be fixed parallel



The patrix with the paralleling mandrel is waxed onto the wax pattern.



Depending on the jaw situation, use the 90° or 120° patrix.



Place the precision-fit duplicating matrix onto the patrix and fix it.



Prepare the model for duplicating in the usual way. The duplicating matrix must not be modified.



Insert the blue matrix in the duplicating mould at the preshaped point



and prepare the investment material model. Now the blue matrix is in the correct position on the matrix.



Prepare the CoCr pattern in the usual way and connect it with the matrix.



After casting, sandblast the matrix and finish and polish the CoCr object.

#### Double-T Adhesive Connectors at the CoCr frame

Patrixes can be adapted to the jaw situations, no parallelism required.



The patrix is waxed onto the CoCr structure and the shape remains unchanged.



After finishing and polishing of the CoCr object, the crowns are modelled and the matrix is waxed on.



Fit the crown, sandblast the joints using 110 µm aluminium oxide and ensure stress-free glueing to the CoCr object.

## Double-T Adhesive Mini Connector dtk



#### dtk mini

Thanks to the 2 different sizes A+B, the 3 different angles 90°, 105° and 120° and the minimal dimensions of the prefabricated wax patterns, the correct type of connector can be fabricated for all cases.

#### Accessories:



DTK-adhesive REF 540 0010 6



#### dtk mini front

dtk-front for tooth-bounded gaps in the anterior region. No problems with space when setting up anterior teeth, even in cases with severe overbites.



#### dtk mini super flat

dtk-super flat: A super flat connector for use in the posterior region. Maximum strength yet requires only a minimal amount of space.

dtk mini	90° A	90° B	105° A	105° B
Patrix	L 11.0 mm W 3.0 mm H 4.0 mm	L 7.5 mm W 2.5 mm H 3.5 mm	L 10.0 mm W 3.0 mm H 4.0 mm	L 7.5 mm W 2.5 mm H 3.5 mm
16 pieces 50 pieces	REF 430 0693 A REF 430 0694 A	REF 430 0693 B REF 430 0694 B	REF 430 0699 A REF 430 0700 A	REF 430 0699 B REF 430 0700 B
Matrix	L 5.5 mm W 3.0 mm H 3.0 mm	L 3.0 mm W 2.5 mm H 3.0 mm	L 5.5 mm W 3.0 mm H 3.0 mm	L 3.0 mm W 2.5 mm H 3.0 mm
16 pieces 50 pieces	REF 430 0691 A REF 430 0692 A	REF 430 0691 B REF 430 0692 B	REF 430 0697 A REF 430 0698 A	REF 430 0697 B REF 430 0698 B
Duplicating matrix	L 5.5 mm W 3.0 mm H 3.0 mm	L 3.0 mm W 2.5 mm H 3.0 mm	L 5.5 mm W 3.0 mm H 3.0 mm	L 3.0 mm W 2.5 mm H 3.0 mm
16 pieces 50 pieces	REF 430 0689 A REF 430 0690 A	REF 430 0689 B REF 430 0690 B	REF 430 0695 A REF 430 0696 A	REF 430 0695 B REF 430 0696 B

50 pieces	REF 430 0690 A	REF 430 0690 B	REF 430 0696 A	REF 430 0696 B
dtk mini	120° A	120° B	dtk mini front	dtk mini super flat
Patrix	L 10.0 mm	L 7.5 mm	L 9.0 mm	L 10.0 mm
	W 3.0 mm	W 2.5 mm	W 2.0 mm	W 5.0 mm
	H 4.0 mm	H 3.5 mm	H 2.0 mm	H 2.0 mm
16 pieces	REF 430 0705 A	REF 430 0705 B	REF 430 0711 0	REF 430 0717 0
50 pieces	REF 430 0706 A	REF 430 0706 B	REF 430 0712 0	REF 430 0718 0
Matrix	L 5.5 mm	L 3.0 mm	L 5.5 mm	L 6.0 mm
	W 3.0 mm	W 2.5 mm	W 2.0 mm	W 5.0 mm
	H 2.5 mm	H 2.5 mm	H 1.5 mm	H 2.0 mm
16 pieces	REF 430 0703 A	REF 430 0703 B	REF 430 0709 0	REF 430 0715 0
50 pieces	REF 430 0704 A	REF 430 0704 B	REF 430 0710 0	REF 430 0716 0
Duplicating matrix	L 5.5 mm	L 3.0 mm	L 5.5 mm	L 6.0 mm
	W 3.0 mm	W 2.5 mm	W 2.0 mm	W 5.0 mm
	H 2.5 mm	H 2.5 mm	H 1.5 mm	H 2.0 mm
16 pieces	REF 430 0701 A	REF 430 0701 B	REF 430 0707 0	REF 430 0713 0
50 pieces	REF 430 0702 A	REF 430 0702 B	REF 430 0708 0	REF 430 0714 0

#### Double-T Adhesive Mini Connector dtk

#### Assortment

dtk mini A + B 90°, 105°, 120°

with 2 connectors each 90°, 105°, 120°

- 1 Paralleling mandrel 90°
- 1 Paralleling mandrel 105°/120°
- 2 anterior connectors2 super flat connector

REF 430 0558 0

Assortment dtk mini A + B

3 Patrices each

909

- 6 Matrices each
- 3 Duplicating matrices each

REF 430 0684 0

Paralleling mandrel REF 430 0623 0

Assortment

- dtk mini A + B 105°
- 3 Patrices each 6 Matrices each
- 3 Duplicating matrices each

REF 430 0685 0

Paralleling mandrel REF 360 0112 0

Assortment

- dtk mini A + B 120°
- 3 Patrices each 6 Matrices each
- 3 Duplicating matrices each

REF 430 0686 0

Paralleling mandrel REF 360 0112 0

Assortment

dtk mini front AI + B

- 3 Patrices each
- 6 Matrices each
  3 Duplicating matrices each

REF 430 0687 0

#### Assortment

dtk mini super flat Al + B

- 3 Patrices each
- 6 Matrices each
  3 Duplicating matrices

each

REF 430 0688 0

dtk mini





Thanks to the 3 different angles of the patrixes, the prefabricated wax patterns can be jawn optimally.



Precisely fitting duplicating matrices blocked out on the patrices - ready for duplicating.



Once the crome cobalt has been polished, adhere - non-stressed - with auto-curing resin cement or composite.

dtk mini front





The dtk-Front in minute, for use in the anterior region. The bar has a notch on its underside to ensure that the papillae remain unimpeded at all times. As this component has been reduced as much as possible, it is only used in tooth-bounded gaps in the anterior region.



Precisely fitting duplicating matrices in the patrices. Block out and duplicate using standard procedures. Can be integrated into the chrome cobalt denture optimally, even if the alveolar ridge is very narrow.



Non-stressed connectors without having to solder – even possible where too little space is available. There is always sufficient space to arrange the anterior teeth as required, even in cases where the teeth are very small.

dtk mini super flat





Super flat adhesive connector - total height only 2 mm for use in the posterior region. Adheres extremely well thanks to the retentive surface being as large as possible. As the waxed surface is relieved to prevent the papillae being impeded, it can be adapted perfectly to the alveolar ridge.



Precisely fitting, super flat, duplicating matrix. This is replaced with the matrix - with an 0.2 mm cement gap - in the duplicating mould.



Super flat, non-stressed adhesive connectors in the posterior region – sufficient space occlusally for setting up denture teeth. Adhere instead of soldering – even if too little space is available.

## **Optiguss**

The solution for increased perfection with less

Optiguss Micro – 5 micron coating – or Optiguss Macro - 10 micron coating - can be applied easily and quickly to the wax pattern to smooth, seal and reinforce it without changing its shape.

The use of Optiguss reduces the finishing time by 50 % compared to a conventional cast surface.



Optiguss-macro 15 ml REF 520 0092 0

Optiguss-micro 15 ml REF 520 0093 0



Optiguss mixing well macro 2 pieces REF 390 0035 0



3 Brushes size A + brush holder REF 330 0114 6

3 Brushes size B + brush holder REF 330 0114 7

3 Brushes size C + brush holder REF 330 0114 8



3 Brushes size A

3 Brushes size B

3 Brushes size C 2 Brush cleaning pot

1 Brush cleaner


REF 520 0091 0



Brush cleaning pot 2 pieces REF 390 0037 0



Brush cleaner 20 ml REF 520 0094 0



Deep fissures, which cannot be reached with a fissure bur, can be smoothed with Optiquss. This simplifies polishing of gnathologically designed occlusals.



waxed-up as carefully as possible, minute scratches and rough areas remain in the wax which have to be trimmed out of the casting.

Even when the pattern is



Approximal contact areas are strengthened, yet retain their shape.

Applying Optiguss creates

super smooth surfaces.



Fitting surfaces are built-up properly and smoothed, which reduces the time required for trimming.



The finishing time can be reduced by more than 50 % due to more homogeneous surfaces.



## Exaktosil N 15 / N 21

**Tests have proven the excellent properties of Exaktosil!** The silicone duplicating materials Exaktosil N15 and N21 with a processing time span of 5 - 6 minutes are highly fluid and hence ensure accurate reproduction of details. Thanks to the exceptional restoring capacity, the high tear resistance and elongation at rupture, Exaktosil silicone duplicating materials protect moulds against damage when removing them and offer technicians a superior level of quality. The suitable silicone duplicating material for all types of indications - Exaktosil!




Linear dimemsional change: 1.8 %0 (according to DIN EN 24 823) A specimen (stylized dental arch) is duplicated with Exaktosil N 21.



Exaktosil N 15 Component A 1000 ml REF 540 0114 A Exaktosil N 15 Component B 1000 ml REF 540 0114 B






Comparison of the dimensions of the specimen and the duplicating mould.
The extraordinary low shrinkage of only 1.8 % ensures precise fit of CoCr objects.



Exaktosil N 15 Component A 5000 ml REF 540 0115 A Exaktosil N 15 Component B 5000 ml REF 540 0115 B

Assortment Exaktosil N 15 5000 ml A 5000 ml B

REF 540 0103 9



The high tear resistance of approx. 1.25 N/mm<sup>2</sup> and an elongation at rupture of approx. 350 % protects duplicating moulds against damage when removing the material from the mould.



Exaktosil N 21 Component A 1000 ml REF 540 0116 A Exaktosil N 21 Component B 1000 ml REF 540 0116 B



REF 540 0114 7







Exaktosil N 21 Component A 5000 ml REF 540 0117 A Exaktosil N 21 Component B 5000 ml REF 540 0117 B

## Assortment Exaktosil N 21

5000 ml A 5000 ml B REF 540 0114 8





## **Duplicating**

## Technosil duplicating silicone



Addition-cured, shrinkage- and filler-free duplicating material for dimensionally accurate duplicates. Technosil NT is mixed in the ratio of 1:1 for simple processing. The shore hardness of 25 makes the material suitable for "ringless" model fabrication with the bredent duplicating system.

Technosil NT duplicating silicone 1000 g of Component A REF 540 TS01 A Component B REF 540 TS01 B



The short setting time allows to continue working quickly. Reduced shrinkage for accurate models.



Technosil NT duplicating silicone 5000 g of Component A REF 540 TS05 A Component B REF 540 TS05 B

# Assortment Technosil NT duplicating silicone component A + B 1000 g each REF 540 TS01 0

Assortment
Technosil NT
duplicating silicone
component A + B
5000 g each
REF 540 TS05 0

## **Technolit**



Surface tension reducing agent avoids the formation of bubbles and improves the flow characteristics of investment material and plaster.

Technolit 125 ml REF 520 ET12 5



After a reaction time of 2 minutes the duplicating mould is blown dry using compressed air. Technolit avoids surface segregation for investment materials and plasters. Consequently, a more homogeneous surface is achieved.



Refill package 750 ml REF 520 ET75 0



# **Duplicating** system

The duplicating method as major element and basis for highly accurate duplicates. The stable plastic components ensure precision during duplicating and reduce errors.



The flask tray serves as basis for the flask sleeve.



The flask sleeve is placed into the flask tray to ensure a stable position.



Flask tray large, REF 520 DBKS G small, REF 520 DBKS K



The spacer - base insert is filled with block-out material to ensure safe hold of the model when duplicating and to exclude shifting.



The block-out kneading material is used to fix the model and to block out undercuts. It will not bond to the silicone and can be reused.



Flask sleeve large, REF 520 DBKM G small REF 520 DBKM K



The model is placed in a central position onto the block-out kneading



The stabilizer is put into the opening of the flask sleeve and the height is adusted according to the model. This protects the silicone mould against undesired deformation when filling the mould.



Spacer - base insert REF 520 DBPE G small. REF 520 DBPE K



The flask sleeve is filled with Technosil.



The duplicating mould is fixed using the aluminium investment aid. Stress-free model fabrication is quaranteed on every type of surface.



Stabilizer large, **REF 520 DBBS** small REF. 520 DBBS K



Small and large sets.

Assortment

large, 5 pieces


1 flask tray

insert

1 stabilizer

1 flask sleeve

1 spacer - base



22 pieces

- flask tray each large and small
- flask sleeve each large and small
- spacers base inserts each large and small
- stabilizers each large and small
- investment aids, aluminium
- block-out kneading materials

125 ml Isosil

je 1000 g Technosil

duplicating silicone A+B

125 ml Technolit REF 520 DBST E



Investment aid, aluminium REF 520 DBAL W



**Block-out kneading** 100 g REF 540 0101 8

### Assortment

small, 5 pieces 1 flask tray

- 1 flask sleeve
- 1 spacer base insert
- 1 stabilizer
- 1 investment aid, aluminium

REF 520 DBST K

### aluminium

REF 520 DBST G

1 investing aid,

### Isosil



Isosil 125 ml REF 520 IS12 5



Refill package 750 ml REF 520 IS75 0



Plastic components that are wetted with Isosil allow easy removal or repositioning of the duplicating mould.



# Master-Copy



The perfect model duplicating system with transfer into the articulator.



Master-Copy base plate 1 piece REF 360 0124 0



Master-Copy stabilizer small 1 piece REF 360 012S K



Master-Copy base plate ring 1 piece REF 360 0124 1



Master-Copy silicone sleeve small 1 piece REF 360 012M K



The investment material model in the articulator features exactly the same occlusion and precision of fit as the master model.



Master-Copy silicone sleeve large 1 piece REF 360 012M G



Master-Copy base former 1 piece REF 360 0124 2



Master-Copy stabilizer large 1 piece REF 360 012S G



Magnetic plates 50 pieces REF 360 0118 1

Assortment large REF 360 0125 6



Master-Copy base plate 1 piece



Master-Copy stabilizer large 1 piece



Master-Copy base plate ring 1 piece



Master-Copy base former 1 piece



Master-Copy silicone sleeve large 1 piece



Magnetic plates 50 pieces


Prerequisite for the function of the Master-Copy system is the fact that the model has been processed with the Master-Split. Please request brochures on the Master model system.

# Master-Copy



### The initial situation ...

A frequently occurring initial situation. The lower jaw model must be duplicated for the CoCr structure



The base plate is the basis for the master model.
The model produced with Master-Split fits exactly on the base plate.



The master model is fixed on the master model with the magnet.



The base ring is placed on the base plate with the master model.



The snap of the silicon sleeve is ensured by catches in the base ring and the sleeve is held safely.




The stabilizer provides absolute stability and a highly uniform silicone layer in the tooth area.



The Master-Copy duplicating mould is filled with silicone up to the openings of the stabilizer.




Once the silicone has hardened, turn the duplicating mould and remove the base plate.



The model can be lifted by blowing in compressed air and removed from the mould.



If the model is difficult to remove, the base ring can be removed temporarily.



The base former features a magnetic plate.



The duplicating mould is filled with investment material up to 5 mm below the rim.



The base former is placed on and the mould is filled up to the base plate. Entrapped air can be easily recognized.



### ... the result

After hardening, the model is removed from the mould and the points for casting-on are ground. The model can be placed into the articulator.



# **Duplicating**

### Bre-Gel 1

Low-viscous agar duplicating gel for precise investment material models, suitable for microwave units



Bre-Gel BG 1 6000 ml REF 540 0103 6



**Low-viscous consistency** Bubble-free casting thanks to excellent flow characteristics.



Transparent color Perfect control during exposure of the model due to transparent color and low viscosity.




High edge stability Stable edges ensure precise reproduction of details of the duplicate models.

Remeltable
Low viscosity to ensure bubble-free casting.



Can be remelted in the duplicating unit or the microwave at least 20 times due to the reversibility.





A low pouring temperature with minimum difference between the gel and model guarantees tension-free, detailed duplicates.

# Bre-Gel 2 opaque, Bre-Gel 3 opaque-liquid

Opaque duplicating gel for the entire duplicating technique, suitable for microwaves.




Bre-Gel BG 2 opaque 6000 ml REF 540 0105 3



Bre-Gel BG 3 opaque-liquid 4 x 400 ml REF 540 0105 4



**Opaque Color.**The bright, opaque color simplifies the evaluation of filigree duplicating areas.



High tensile strength
The high elasticity and
tensile strength allow
easy removal from the
cast even in undercut
areas. Thus precise working is also possible in the
resin casting technique.




The outstanding elasticity ensures recovery of deformed duplicating areas when removing the duplicate model.



### Microkeramik

Perfect cast surfaces thanks to microfine ceramic layers for crowns and bridges and CoCr work.



In the field of crowns and bridges, Microkeramik is especially suitable for NP alloys since very fine cast surfaces are obtained. The Microkeramik is adapted to the expansion of the investment material.



An extended processing time span allows precise application of the Microkeramik. Microscopically fine ceramic particles ensure perfect reproduction of very fine details of wax models.



Microkeramik

with

without

Devesting is simplified since there is no bonding between the investment material and Microkeramik.



Accessories:

The difference after sandblasting with glass beads can be clearly recognized: The entire oxide layer can be asily removed so that less working time is required.



Microkeramik avoids extreme formation of oxide on NPM alloys. Cast objects are only sandblasted with 50  $\mu$  glass beads to obtain almost perfect high luster. Consequently, time for further processing is saved.



Microkeramik 125 g REF 550 0001 2



3 **Brushes** size A + 1 brush holder

3 **Brushes** size B

+ 1 brush holder

3 **Brushes** size C + 1 brush holder

REF 330 0114 6

REF 330 0114 7

REF 330 0114 8

### Brevest M1

Very precise, universal investment material for all CoCr alloys. Precision-fit crowns and bridges, clasps and CoCr attachment work as well as one piece casting work can be produced with two different liquids.



Bresol N \* 1000 ml bottle REF 520 000N 1

5000 ml canister **REF 520 000N 5** 

Brevest M1 40 bags 200 g each REF 570 0000 8

100 bags 200 g each **REF 570 0002 0** 

\* frost-resistant

### Assortment

20 bags 200 g each Brevest M1 1000 ml Bresol C+B \* 1000 ml Bresol M \*

REF 570 0002 2

Please order the course documents for the attachment course VS 3 and one piece casting!

### Accessories:



Dosing bottle REF 520 0101 1




The frost-resistant Bresol C+B liquid which is suitable for expansion control is used for crowns and bridges.



For the precise one piece casting technique different expansion values can be obtained using Brevest M1. The extended reworking time span of 5 to 6 minutes provides the perfect precondition for this purpose.



Dosing syringe 6 pieces REF 520 0101 2



The frost-resistant precision liquids Bresol C+B and Bresol M are perfectly suitable for all types of CoCr work in the entire field of dental techniques.



Due to optimal expansion control perfect fit of attachments and CoCr clasps can be achieved.

### Investing and casting

# Brevest Rapid 1



Rapid-heating, universal precision investment material for crowns and bridges as well as the entire field of CoCr work.

Bresol R 1000 ml bottle REF 520 000R 1 REF 520 000R 5 Brevest Rapid 1 50 bags 160 g each REF 570 160R 8 125 bags 160 g each REF 570 16R2 0

Brevest Rapid 1 40 bags 200 g each REF 570 000R 8 100 bags 200 g each REF 570 00R2 0

Accessories:

Dosing bottle REF 520 0101 1 Dosing syringe 6 pieces REF 520 0101 2



Fine grained, rapid-heating precision investment material for all largespan bridges, can also be used without casting



Perfectly suitable for one piece casting. Precise expansion control with Bresol R



Accurate and precise attachment work and CoCr clasps - even if little time is available.



25 bags 160 g each Brevest Rapid 1 1000 ml Bresol R REF 570 160R 4

20 bags 200 g each Brevest Rapid 1 1000 ml Bresol R REF 570 0002 5



Brevest Rapid 1 can be placed into the furnace . at a temperature of 900 °C already 15 minutes after mixing.

### Brevest exakta M and Brevest exakta Speed

Phosphate-bonded investment materials for gel and silicone duplicating. The expansion for attachment work and CoCr clasps can be precisely controlled with the frost-resistant special mixing liquids.



Brevest exakta M 20 bags 400 g each REF 570 00XM 8 50 bags 400 g each REF 570 0XM2 0

Bresol N \* 1000 ml bottle REF 520 000N 1 5000 ml canister REF 520 000N 5

Assortment

10 bags 400 g each Brevest exakta M 1000 ml Bresol N 3

REF 570 0002 3





Brevest exakta Speed 20 bags 400 g each REF 570 0ES0 8 50 bags 400 g each REF 570 0ES2 0

Bresol Speed \* 1000 ml bottle REF 520 000S 1 5000 ml REF 520 000S 5

\* frost-resistant

### Assortment

10 bags 400 g each Brevest exakta Speed 1000 ml Bresol Speed \*

REF 570 0ES0 4



Dosing bottle REF 520 0101 1





Brevest exakta M and Brevest exakta Speed are particulary suitable for gel duplicating. After devesting, the investment material model is hardened in Duro-Top immersion hardener.

### Silicone duplicating



These investment materials feature good flow characteristics and a processing time span of 2 to 3 minutes. No tension reducing agent is required for silicone duplicating.



Dosing syringe 6 pieces REF 520 0101 2



### **Brevest ESG**

Special investment material in the one piece casting technique for extra-smooth surfaces.



Assortment Brevest ESG 20 Beutel 200 g each 500 ml Bresol ESG 1 Transfuser REF 570 ESGO 4





- Harmonized products to reduce effort and allow for easy to achieve precision castings
- Safe making of dental restorations by simple work
- Reduction of miscasts allows efficient work
- Repeatable results boost your success
- High precision fit minimizes trimming time



A carefully blocked out model shortens wax up time and reduces fitting efforts of Secondary construction.



A correct sized duplicating flask is chosen. This doublicating system reduces waste of silicone and is the basis for precision fit

secondary constructions.



Only the primary parts are poured with special investment ESG using the Transfuser silicone tool. This gives a defined expansion for the precision fit secondary construction.



The entire mould is poured with Brevest Rapid 1 after setting of ESG investment. Both investments join very well.



The correct mixing ratio of each investment is key for a perfect fit. Adapting the investments to your lab equipment is provided by a bredent System Consultant on site. This guarantees smooth workflow and repeatable results.



The special funnel avoids velocities when casting and reduces trapped air inside the cast.



Ready for veneering using visio.lign shells which stand for fit, esthetics and speed



# Investing and casting

# Brealloy F 400





**CoCrMo alloy for clasps and attachments in chrome cobalt restorations.**Brealloy F 400 is nickel-free and complies with the standard
DIN EN ISO 6871 – part 1: 1996.

Brealloy F 400	VPE	100 g	500 g	1000 g
Cylinder, 7.5 g each	REF	500 ML10 0	500 ML50 0	500 ML00 0









The outstanding material properties of Brealloy F 400 allow rapid finishing and polishing.







Brealloy F 400 features a hardness of 400 HV 10. The alloy has been especially developed for non-precious attachment dentures. The chrome cobalt system of bredent offers additional innovative techniques allowing the production of locks and individual screw connections using Brealloy F 400. The combination of the physical values of Brealloy F 400 allows to obtain extremely slender chrome cobalt clasp dentures. Your patients will be enthusiastic about the high comfort of wear of these dentures.

### Accessories:

 Brealloy solder
 Brealloy flux

 7 g
 8 g

 REF 500 0001 0
 REF 500 0001 1

# Physical values (quide values)

Density (g/cm <sup>3</sup> )	8.4
Vickers hardness (HV 10)	400
Solidus point (°C)	1320
Liquidus point (°C)	1380
Casting temperature (°C)	1480
0.2 % proof stress (MPa)	700
Modulus of elasticity (MPa) ap	prox.220.00
Tensile strength (MPa)	900
Elongation at break (%)	4
Expansion coefficient	
(TEC 25 - 600 °C)	15 μm/mk

Composition (in mass %)	
Cobalt	64.7
Chrome	29
Molybdenum	5
Manganese	0.4
Silicone	0.5
Carbon	0.4

# **Brealloy MO**



An alloy designed to meet the requirements of the production of clasp and attachment model castings and for single shot casting technology. Easy shaping reduces milling material use. Brealloy MO is nickel-free.

Brealloy MO

100 g
REF 500 M010 0

500 g
REF 500 M050 0

1000 g
REF 500 M00 0

# Physical properties (guide values)

Density (g/cm³) 8.3	
Vickers hardness (HV 10)	380
Solidus point (°C)	1260
Liquidus point (°C)	1350
Casting temperature (°C)	1420
0.2% proof stress (mPa)	640
Tensile strength (N/mm²)	700
E-modulus (mPa)	210,000
Elongation at break (%)	<6

Composition (in % of mass)	
Cobalt 62.2 Chrome Molybdenum Silicone Manganese Carbon Others	30 5.5 1.0 0.6 0.6



The high e-modulus permits production of delicate brace prostheses.



Easy shaping of brealloy MO facilitates the production of attachments.



Individual bars can be produced to fit perfectly.

### Accessories:

 Brealloy solder
 Brealloy flux

 7 g
 8 g

 REF 500 0001 0
 REF 500 0001 1



# Brealloy solder



Brealloy solder 7 g REF 500 0001 0 Solder especially matched with CoCr alloys for chrome cobalt and ceramic bonding techniques to avoid the formation of galvanic elements and undesired reciprocal action with the ceramic material.

# Brealloy flux



Brealloy flux 8 g REF 500 0001 1 Suitable for all CoCr alloys, supports the flow characteristics of the solder.

# Duro-Top



Duro-Top 1000 ml REF 570 0005 4

Immersion hardener for precise and clean

modelling on duplicate model surfaces.

### For the agar duplicating technique



Immersion hardening liquid especially for agar duplicating – for sealing model surfaces.

### Stabilization of edges



Thin edges and filigree areas withstand increased stress due to the immersion process.

### Surface smoothing



Prefabricated wax elements adhere to the smooth model surface without using any adhesive.

### **Excellent diffusion**



Due to the low-viscous consistency the hardener easily penetrates into the surface.

### High yield



Excellent hardening effect and robust models are obtained even after numerous immersion processes.

### Investment hardener



Improves the hardness and surface texture of all models duplicated in silicone.





The improved strength toughens the edges and prevents damage to the fine wax-coated margins.



The greater scratch resistance allows waxing up without damaging the model surface.



# Investing and casting

# Crepe sleeve



Crepe sleeve 25 m REF 570 0002 1

For individual overbedding of CoCr work.

- Surface enlargement
- Uniform absorption and release of heat
- Investment material is saved



### Investment marker



Helps with the positive identification of investment muffles.

Investment marker REF 330 0115 0



The necessary information is noted down quickly and easily.



The marker can be clearly read on all investment materials up to 1100° C.

# Sprue



made of high-quality plastic 25 pieces REF 360 0002 5



Specially shaped sprue for CoCr work. Made of high-quality plastic for extended durability. Compared to the shape of conventional sprues, this shape improves the filling behavior.

# Golden booklet



breden	Übertrag	Angelo	Verbenuch	Detrers	Delimined

Golden booklet DIN A 6 REF 610 0020 0

Thanks to the clear and simple structure of the golden booklet, reliable stock-keeping of precious metal alloys is ensured. The booklet simplifies the control and provides a quick survey on the consumption of alloys.





Golden booklet DIN A 4 REF 610 0010 0



# Implementing the repeatable solution in your lab!

Competent practical advice and adjustment of equipment by our Special System Consultant – for your immediate success.

- Course content is adapted to daily cases in your lab
- Up to three participating technicians without any fee that's worth it!
- Taking place at your lab saves you money
- Cost effective and prompt on site converting into the new system
- Course content is adapted to individual requirements of each lab

### Content of Workshop:

The Special System Consultant checks and adjusts all corresponding materials and machinery. This means you can achieve reproducible results straight after the course.

On the course a secondary crown is made on a primary one you have prepared prior to the training. After casting the secondary construction is easily fitted on to the primary one.

### Duration of Workshop:

Duration of this workshop is one day. Whilst this beneficial action takes place, lab business can go on as usual. Multiple technicians can participate without extra fees.

#### Workshop order number:

The "Precision Casting" workshop can be ordered under REF 950 0074 0. The specially trained System Consultants then will meet you at your lab by appointment.

Book an individual appointment for this informative and beneficial workshop.





Photo: DL Marco Zelmer, Sondershausen

Highest precision and perfect fit of the denture due to the manufacturing process. By avoiding galvanic caps space is gained for esthetic veneering – minimizing costs.







Photo: DL René Thiere, Gera

Durable dental restorations requiring no maintenance such as fixing the loss of friction. This saves the patient's money.

# Milling unit BF 2



Milling unit BF 2 including 1 Model support BF 2 1 piece

REF 140 0098 0

### Technical Data

230 Volt / 50/60 Hz Power supply 80 Watt Power rating Speed 0 - 30,000 U/min. Ø 2.35 mm Chuck Fuse thermal overload protection Torque 2.6 Ncm 18 kg Weight Width/Depth/Height 250 x 370 x 510 mm

### Accessories:

Chuck 2.35 mm	REF 730 0016 9
Chuck 3 mm	REF 730 0015 3
Tap handwheel	REF 330 0115 4
Model support BF 2	REF 730 0017 0
Milling base	REF 140 0089 3
Adapter airaqua turbine	
16 mm	REF 730 0018 4
18 mm (for BF1)	REF 730 0018 3
28,5 mm	REF 730 0018 5
Transfer device	
3 mm shaft	REF 360 0116 3
2.5 mm shaft	REF 360 0126 5



# Milling base



Milling base with integrated thread for fixation on the milling base of the BF 1 unit. Additionally, plaster is removed completely and without damaging the metal plate by slightly turning the locking bolt.

Milling base 1 piece REF 140 0089 3



### Model support BF 2



The model support can be used for all milling units including units with magnetic circuit. Turning by 90° permits to perform lateral drilling of bars without removing the model.

Model support BF 2 1 piece REF 730 0017 0



### Transfer device



Permits correct transfer of the position of the model to the milling base. Up to 8 units can be transferred at the same time.

Transfer device 3 mm shaft REF 360 0116 3 2,35 mm shaft REF 360 0126 5

# Brenometer surveying system



Brenometer surveying system

Four different surveying plates according to Ney allow accurate positioning of the clasp profiles whilst ensuring correct depth of undercuts. A locating pin and a red marker with a holder ensure correct surveying.



Marking the clasps and surveying with a single unit – this is how time and money can be saved.

### Assortment

- 1 Brenometer marker holder
- 1 Brenometer locating pin
- 1 Brenometer plate 0.25
- 1 Brenometer plate 0.35
- 1 Brenometer plate 0.50
- 1 Brenometer plate 0.75

REF 310 0000 2

### Refill packages:

Brenometer marker holder REF 310 0000 4 REF 310 0000 3 Brenometer plate 0.25 REF 310 0002 5 Brenometer plate 0.35 REF 310 0003 5 Brenometer plate 0.50 REF 310 0005 0 Brenometer plate 0.75 REF 310 0007 5




# Restoring the friction

# Activating pliers



The problem: Conical and telescopic crowns have lost their friction.

The solution:
Activating pliers - Pliers which recreate the retentive forces for telescopes which have lost their function.



bredent Activating pliers REF 320 0043 0

Provide conical and telescopic crowns with "new" friction easily and quickly.



The friction zones in the outer coping create 3 new contact areas between the inner and outer copings. This restores the retentive friction. Should the unit be activated too much, the surface can be trimmed to reduce the friction.




These pliers have a ball and socket for creating one or several new friction zones. The long lever of the pliers enables the forces to be applied as required.



The activating pliers can also be used to reduce a friction zone which is too retentive. If necessary, the facing should be removed for activating the unit and replaced again later.

### Novo-Grip



Sharpenable, sintered diamond coated inserts for conical pliers with "grip".

Novo-Grip pliers
1 pair of pliers
+ 2 standard inserts
+ 1 Allen key
REF 310 0000 8


Novo-Grip pliers 1 pair of pliers + 2 small inserts + 1 Allen key REF 310 0011 3



2 pieces REF 310 0001 A

standard insert

Novo-Grip



Novo-Grip small insert 2 pieces REF 310 0001 B

Accessories:



Diabolo cleaner grindstone for inserts 1 piece REF 340 0100 0

Set screws M3 4 pieces REF 310 0011 2





Exchangeable insert with a diameter of 2.35 mm, also for small primary crowns.

Special shafts



Hardened shafts provide high stability even if strong pressure is exerted.

Rotatable



Worne inserts can be rotated. New diamond grains ensure that inner surfaces of crowns can be held safely again.

Sharpenable



To regain the maximum abrasiveness, the inserts are clamped into the handpiece and new diamond grains are obtained on the surface using the grindstone.



# Friction fit system FGP



### Individual friction for highest demands.

The friction fit system offers the dentist and the dental technician an entirely new perspective during the preparation and the restoration of the friction for all types of telescopic metal restorations. Long service life and simple, time-saving processing render the friction fit system a comfortable solution for your patients.

### Application fields of the FGP system



#### Safety and outstanding quality

The FGP system by bredent offers optimum and individual friction when preparing new conical and telescopic restorations.



# Direct solution instead of extended waiting times

Due to the use of FGP directly in the dental practice. The simple use during the restoration of the friction of telescopic work is the solution for the dentist and the patient.



### Individuality and precision

These requirements can still be fulfilled even in areas difficult to access, whether new dental supply or relining work are concerned.



#### No compromises

During the preparation of new individual attachments. The FGP system allows to obtain results that fulfill highest demands.

### 20 years of experience with the FGP

Discover the personal advantages:

- Saving of time due to fast and simple preparation
- Preparation of individual friction at favourable costs
- No fitting of secondary elements
- Long service life
- Maximum comfort of wear for the patients
- Allows low-cost single-piece casting
- Can be processed in the mouth
- Almost without any wear
- Low susceptibility to plaque thanks to highly compacted resin surface

Up until today these advantages have contributed in more than 50,000 cases to achieving soft integration and removal of the denture.

The principle of the FGP resin is based on the fact that the metal fit that has been common in the telescopic technique so far will now be replaced by a metal-resin fit

The metal-resin fit offers the benefit of a considerably more favourable coefficient of friction than the one of a pure metal fit. Consequently, increased resistance to wear and extended service life are obtained.



# Restoring the friction

# Friction fit system FGP

### New fabrication of telescopic crowns



Thermo-forming or immersion wax copings serve as spacer for the FGP resin.



with a wall thickness of at least 0.2 mm ending 1 mm above the cervical margin.



The investment material model with cervical step is prepared before



the usual outer telescopic and cast pattern.



After casting – made with any alloy –



the cast frame is finished and veneered with resin or ceramic materials.



Due to the preparation during the modellation a gap resulted which is now filled with FGP.



In a preparatory step the pattern is insulated.



FGP bonding agent is applied equally thinly onto the inner surfaces.



The material is hardened at air for 5 minutes; during this time a visible layer is obtained.



The FGP two-component resin is mixed in the ratio of 1:1



and filled into the outer telescopes without any bubbles.



The restoration is placed onto the model exerting uniform pressure.



The hardened FGP resin with a clearly visible border at the cervical margin.



The FGP system offers individual friction with maximum comfort of wear.

### The enhanced friction

Tests and scanning electron microscope studies with FGP reveal clearly better values of friction than those of metal fits.



Conventional metal/ metal fit. Metal fit after completion adjusted to a frictional force of 8 Newton.



FGP resin/metal fit. Resin fit after completion adjusted to a frictional force of 8 Newton.

For this comparison between a classical metal fit and a FGP fit 21,000 integration and removal processes were simulated. This corresponds to a period of wear of approx. 20 years.



Scanning electron microscope picture of the inner side of a telescopic secondary element made of a precious metal alloys with a magnification x 100.



Scanning electron microscope picture of the inner side of a telescopic secondary element made of FGP resin with a magnification x 100

Result:

Residual friction 2 Newton, that is only 25 %. Result:

Residual friction 6 Newton, that is still 75 %.



# Friction fit system FGP

### Restoration in case of loss of friction



Telescopic work after numerous years of wearing.



During the integration there is no sufficient friction.



Primary telescopes in situ prior to friction relining.



The dial caliper is used to measure the thickness of the outer telescopes.



The outer parts are ground to obtain space for the FGP resin.



Any residual grinding particles are removed with compressed air.



Retraction threads are put around the primary elements.



Then the inner telescopes are insulated with a small amount of liquid vaseline.



FGP bonding agent is applied equally thinly onto the inner surface of the outer parts.



The FGP two-component resin is mixed in the ratio of 1:1



and filled into the outer telescopes without any bubbles



After the denture has been integrated, the patient is able to bite evenly exerting normal masticatory pressure.



The resin residues must be removed with the probe. Approx. 120 seconds after beginning of mixing, remove the restoration from the primary elements and place it on again.



The denture is removed after approx. 7 minutes and excess material is removed with a rotating tool.



The result is a functional denture that exhibits excellent comfort of wear within a very short period

### FGP in implantology Absolutely tension-free fit.



The excellent sliding properties of FGP resin ensure gentle, implant-protecting integration and removal of the supraconstructions.



Even very small tensions in the low-cost and biocompatible single-piece casting process are perfectly compensated.



The high resistance to abrasion and non-tilting integrating and removing of the supraconstruction provide patients with a high comfort of wear and simple handling of their dentures.



The friction with FGP resin that will remain stable over many years guarantees the patients' happiness and satisfaction.

# Assortment Friction fit system FGP REF 540 0102 8

- 1 x 2.5 g Friction resin component A
- 1 x 2.5 g Friction resin component B
- 1 x 1.25 ml FGP bonding agent
- 1 x 3.0 ml FGP insulating agent
- 1 Spatula
- 5 Brushes
- 1 Brush holder
- 1 Mixing block



### Refill packages:

Friction resin component A	REF 540 0108 A
Friction resin component B	REF 540 0108 B
FGP bonding agent	REF 540 0102 6
FGP insulating agent	REF 540 0102 7

#### Accessories:

Mixing block		
35 x 50 x 10 mm	10 pieces	REF 330 0114 4
Disposable brushes	100 pieces	REF 330 0114 2
Spatula	100 pieces	REF 330 0114 3
Brush holder, bent	12 pieces	REF 330 0114 1
Application cannulas, black	25 pieces	REF 580 0001 8



# **INNOVATIONS**



For almost 40 years, bredent has offered innovative solutions for use in the dental technology laboratory – this is a central part of the company's philosophy.

New developments and procedural techniques will significantly influence the future of dental technology and dentistry.


Due to close contact with customers and monitoring of national and international dental markets, ideas and recommendations influence modifications in the product portfolio. This exchange enables the processes in the laboratory and in the practice to be optimised and costs to be reduced. bredent strive for and desire consistent implementation of the quality standards in accordance with ISO 9001, which offers the user maximum product safety and, as a result, ensures patients receive implant prosthetic restoration of the highest quality.

### Efficient through innovation!



# Processing of ceramic and zirconium

The proven bonder system for NPM alloys balances the CTE value and hence avoids chipping of the ceramic veneer. More innovative products from bredent, such as liquids or instruments, contribute to simple processing in the field of ceramics.



Ceram-Bond can be easily applied and the components ensure perfect coverage of the framework. After firing, a golden-yellow coat is formed on the framework, which avoids chipping of the ceramic.



Accurate dosing of the ceramic powder enables perfect reproduction of shades. Shrinkage during firing is reduced by mixing the ceramic powder with the ceramic mixing liquid. As a result, time is saved during processing.



The ceramic brushes with black hair provide good contrast to the ceramic material and enable fatigure-free working.



The ceramic processing set that includes several perfectly matched components facilitates surface finishing and reduces the amount of working time needed.

Chrom-Kobalt-Bonding	270
Ceram-Bond	270
Opaquer mixing liquid	271
Ceramic mixing liquid	271
Stain liquid	271
Porcelain liquid set	
Ceramic separating set	
Quicktool	
Ceramix	
MagicContrast	
MagicBrush	
Magic	
KoliBrush	
Unique Brush	
Brush holder	
Ceramico mixing plates	
breCeram	281
anual processing of zirconium	
compoForm UV	
SERACOLL UV	
Vario-Soft 3 zircon sv	
5-motions-zircon	
Zi-polish - the perfect mirror polish!	283

5-motions-glue and 5-motions-active ...... Scanners and grinders .....

Diagen-Turbo-Grinder.....

# **Ceramic veneering**

# Chrom-Kobalt-Bonding



Chrom-Kobalt-Bonding 4,5 g REF 520 0032 1 19 g REF 520 0032 0

The microfine layer of bonding material ensures a perfect bond between the porcelain and chrome cobalt, fully compensating for differences in their thermal coefficients of expansion.

Reduces the problems for alloys with strong tendency to oxide layer formation

The micro-fine intermediate layer that is fired at 980 °C allows to balance the CTE values of the chrome cobalt alloy and the ceramic material.

Chrome Cobalt Bonding protects against spalling and avoids time-consuming remakes. In cases of unfavorable space conditions, CCB allows to fire the ceramic material directly on the CoCr structure.

### Ceram-Bond



Ceram-Bond 30 g REF 520 0032 2 7 q

REF 520 0032 3

### For increased reliability with all alloys.

The premixed, ready-to-use Ceram-Bond allows to omit oxide firing when veneering metal frameworks.

Ceram–Bond is applied immediately after finishing, sandblasting and cleaning the metal framework.

This micro-fine layer improves bonding of the ceramic material to the metal framework, protects against spalling and offers increased reliability.





# Opaquer mixing liquid

Opaquer mixing liquid 18 ml REF 520 0085 0 200 ml REF 520 0012 2 For enhanced wetting and perfect flow characteristics.







Developed and tested by leading ceramists.

### Tip:

Adding a few drops of opaquer liquid to porcelain mixed with mixing liquid prolongs its working time and facilitates building-up of large-size restorations

# Porcelain mixing liquid

Porcelain mixing liquid 30 ml

**REF 520 0086 0** 200 ml

REF 520 0012 3

- Much less shrinkage thanks to improved condensing properties
- Prevents occlusal and interdental contraction cracks in the porcelain
- Easier to condense







#### Tip:

Mix the porcelain slightly thinner; to obtain the ideal consistency leave it for 2 minutes. If building-up takes a long time, spatulate the mixture from time to time; if necessary add a few drops of porcelain mixing liquid because the porcelain already begins to condense on the mixing slab.

### Stain liquid

Stain liquid 7 ml

REF 520 0084 0

30 ml REF 520 0012 1

- $\bullet$  Provides for an absolutely even glaze
- Holds the stains in place perfectly on porcelain
- Thanks to a new formula, this stain liquid can be used for inlay in stains







# Porcelain liquid set

Porcelain liquid set for testing and comparing

30 ml Porcelain mixing liquid 18 ml Opaquer mixing liquid 7 ml Stain liquid REF 520 0087 0



# **Ceramic veneering**

# Ceramic separating set



Plaster sealing liquid gvs 20 ml REF 520 0012 9

For separating ceramic materials against plaster.

- Extremely thin separating film provides outstanding separating effect
- Suitable for all commercial ceramic materials (also low-melting) thanks to optimal composition
- Harmonized components avoid discoloration of the ceramic materials



Apply plaster sealing liquid onto the areas to be separated so that a homogeneous, smooth surface is obtained



Brush plaster sealing liquid also onto approximal areas. Let the plaster sealing liquid dry for 2 minutes.



Ceramic separating liquid kis 20 ml REF 540 0070 3



Apply ceramic separating liquid onto the plaster model so that a wet, shining layer is obtained.



Ceramic separating liquid also applied onto the approximal contacts. The model must not be dried with compressed air.



Ceramic materials are layered directly on the wet ceramic separating liquid.



The special composition of the ceramic separating liquid avoids discoloration of the ceramic materials.



Thinner for ceramic separating liquid REF 550 0000 3



Remove the ceramic structure carefully from the plaster model.



The ceramic separating liquid provides absolute reliability when producing the model. Spalling of ceramic is avoided; accordingly considerable time can be saved.

### Assortment

20 ml Plaster sealing liquid gvs 20 ml Ceramic separating liquid kis

REF 520 0100 0

### Quicktool



Ceramic structures are held safely without any pressure thanks to the three galvano plated diamond tips and the locking mechanism – even galvano formed copings.




The gripping force can be adjusted to the crown size without deforming the crown.



In case of limited space, a diamond tip can be removed - ideal for lower anterior crowns.



The three bud-shaped diamond tips ensure safe and uniform hold of the crown. Accordingly, safe hold is achieved also when condensing.






diamond tips 3 pieces REF 310 0102 1




The integrated condenser condenses the ceramic material in next to no time.

### Ceramix



Reproducing individual shade mixtures quickly and easily. Ceramic material is saved thanks to controlled dosing.

Ceramix REF 360 0119 5



The desired quantity is determined using the scale and the Ceramix is inserted into the ceramic material.



The corresponding filling quantity is determined for additional ceramic materials.



When inserting the Ceramix into the ceramic material, it must be ensured that the material is properly condensed.

To obtain the same shade mixture continuously, note down the ratio of the mixed shades. This way ceramic material is saved.



To obtain individual sample shade tabs, stir the mixed ceramic material thoroughly with a spatula. Take up with the Ceramix, press onto the fibrous pad, wet and fire.



# **Ceramic veneering**

# MagicContrast



Scale 1:1


# MagicBrush



Scale 1:



# Magic...



The MagicContrast brushes = black and MagicBrush = brown feature absolutely identical functional chacteristics!




The pointed shape of dry brush hair is immediately restored by wetting, tapping off or vibrating the brush lightly.



Thanks to the strong contrast, any contamination such as dust or dry ceramic particles can be clearly recognized.







From the spatula shape to the original shape by rinsing the brush in liquid and then tapping if off or vibrating it lightly with a suitable instrument.



High elasticity to pick up specific quantities of ceramic material more easily.



The optimized springiness restores the shape of the brush tip immediately after picking up ceramic material. Additionally, the respective quantity can be easily assessed thanks to the contrast.



The spatula shape can be easily achieved after pressing the tip with two fingers. This way, the brush can be shaped individually.



Large quantities of ceramic material are picked up and time-saving layering is achieved thanks to the elasticity and the springiness of the brush hair.



The spatula shape does not reduce the stability of the brush hair; hence less time is required for applying the ceramic material to the framework.

# **Ceramic veneering**

# KoliBrush



KoliBrush – golden brown natural hair Natural hair brushes made of superior quality Kolinsky hair.



Improved design of the tip of the BigBrush is achieved tanks to the integrated spheres for simpler modelling. Fine, stable tip thanks to carefully selected hair.




The shape and quality of the hair for perfect retention of moisture and improved adhesion and application of ceramic material.



The desired elasticity is obtained by the unique design and combination. This way already applied ceramic layers will not be damaged.

Scale 1:1	Product name	Size	Qty.	REF
	KoliBrush	4, 6, 8 B	1 piece each	390 KSET 1
	KoliBrush	1	2 pieces	390 K001 0
	KoliBrush	2	2 pieces	390 K002 0
	KoliBrush	4	2 pieces	390 K004 0
	KoliBrush	6	1 piece	390 K006 0
	KoliBrush	8	1 piece	390 K008 0
	KoliBigBrush	8 BigBrush	1 piece	390 K008 B
	KoliBrush	1/0	2 pieces	390 KS01 0
	KoliOpakerBrush	5	2 pieces	390 KS03 0

# Unique Brush



Ceramic processing "par excellence".
Unique Brush – the truly esthetic brush.

With newly developed matt black bristles and unique design of each brush – even more benefits when processing ceramic materials.


- Perfectly shaped brush tips enables accurate placement of the precious ceramic materials
- Superior retention of moisture for prolonged modelling
- High elasticity of the brush hair so that large quantities of ceramic material can be picked up
- Spatula function for layering small gaps
- Individual color of each brush renders the tool unique
- Ergonomic shape enables fatigue-free working

The brushes are available in the assortment (in an attractive case) or individually.



Scale 1:1

Product name **VPE** REF Size Unique Brush Set 390 USET 1 1, 4, 6, 8 1 piece each 390 U001 0 Unique Brush 1 1 piece 390 U004 0 Unique Brush 1 piece Unique Brush 1 piece 390 U006 0 Unique Brush 8 1 piece 390 U008 0



The ergonomically designed shape of the brush enables fatigue-free and relaxed working.



The newly developed, matt black brush hair with extremely high elasticity and excellent retention of the mixing liquid ensures prolonged modelling.



Unique Brush

The application-oriented shape of the brush tip enables accurate placement of the precious ceramic material. The contrast between the brush hair and the ceramic material ensures simple determination of the quantity that has been taken up.

BigBrush



1 piece

The optimized, high elasticity of the newly developed brush hair allows also to take up large quantities of ceramic materials.

390 U008 B



The required new shape and function are easily achieved by pressing the brush tip together with the fingers.



The spatula shape allows to build up very small interdental spaces and to separate them. It renders the brush a perfect tool for the preparation of marginal ridges.



The spatula function of the brush allows to apply ceramic material very accurately and to smoothen large areas. As a result, time during finishing is saved after firing.



The perfectly shaped brush tip is restored by rinsing the tip with water and then tapping off or vibrating it gently with a suitable instrument.



# **Ceramic veneering**

### Brush holder



The ingenious storage tool for brushes and instruments.

- High stability
- Stores up to 14 brushes and instruments
- Soft silicone for safe hold
- No restrictions to sizes
- Prolonging the service life of brushes thanks to perfect storage
- Residual substances (for example ceramic materials) will not reach the area of the collars of the brushes
- No damage to sensitive working sections of instruments
- No risk of injuries caused by pointed or sharp instrument inserts

Stainless steel holder (approx. 600 g) with silicone storage rack to hold and store ceramic brushes and instruments in a perfect manner. The size of the 14 clip openings ensures proper fixation of most ceramic brushes. All brushes will be safely held in the ideal section, which may be the section of the handle or the collar as well. As a result, the brushes are not aligned on a single level but arbitrarily and better handling is ensured. On both sides there is a central opening with a larger diameter to hold brushes with large diameter handles.

Brush holder 1 piece REF 310 0102 9



Two diameters for all brush sizes.



Soft, flexible silicone



Improper storage





Holder with brushes / instruments in different sizes - depending on the diameter brushes / instruments are clipped in the area of the handle or the collar. The resulting differences in height (not aligned on one level) allows perfect access. Sensitive instruments can also be fastened in the brush holder!

# Mixing plates



### Perfectly shaped

Mixing plates "Form follows function"

The fanciful mixing plates from bredent. Featuring new shapes and an exquisite design which reflects the uniqueness. Perfectly shaped, individual, functional and unconventional – they add a personal touch to your working place. Created for mixing ceramics, light-curing materials and stains.

### Individuell

### Ceramico GlossOne

GlossOne glazed mixing plates are not self-wetting systems. A special seal has been adhered into the lid

The lid protects the content from drying too quickly and from becoming contaminated. Wet and clean the seal, exert pressure onto the central section of the correctly positioned lid and the air is displaced. The resilience of the lid creates a slight vacuum. The lid adheres to the plate and hence an airtight seal is obtained.



#### Ceramico GlossOne3

White-glazed mixing plate in the shape of an incisor with three indentations and plastic lid with special seal.

3 indentations (approx. 30 x 20 x 4.5 mm)
Dimensions - plate:
approx. 190 x 144 x 15 mm
Dimensions - lid:
approx. 192 x 148 x 18.5 mm
Weight:
approx. 510 g
REF 390 0040 0



#### Ceramico GlossOne14

White-glazed mixing plate in the shape of an incisor with 14 indentations and plastic lid with special seal.

3 indentations
(approx. 43 x 30 x 5 mm)
4 indentations
(approx. 33 x 22 x 4.5 mm)
7 indentations
(approx. 29 x 20 x 4.5 mm)
Dimensions - plate:
approx. 190 x 144 x 15 mm
Dimensions - lid:
approx. 192 x 148 x 18.5 mm
Weight:
approx. 460 g
REF 390 0040 11



The glazed surfaces allow to apply individual markings which can be wiped off with a finger or a cloth if necessary.





When interrupting the working process, the plate can be covered to protect the ceramic from becoming contaminated. A seal has been integrated into the lid.



After exerting pressure onto the lid, excess air will escape and a tight seal is obtained to ensure that the ceramic material will dry less quickly. How long the necessary consistency of the ceramic material can be maintained, depends on the time of coverage, the ambient temperature and exposure to sunlight, which must be avoided.







### **Ceramic veneering**

### Mixing plates

### Ceramico WetOne

The porosity within the plates allows to achieve necessary wetting to maintain the desired consistency of the ceramic materials. Moreover the smooth surface protects bristles against rapid wear. When storing in the closed container and constantly checking the degree of wetness, the ceramic material can be stored almost for an indefinite period and material consumption is reduced considerably. Moreover time is saved since it is not required to keep ceramic materials permanently wet and additional costs for fleece, filters, etc. can be reduced. When closing the lid during interruptions of work, ceramic materials will not dry and become contaminated.



# Ceramico WetOne3 Self-wetting mixing plate in the shape of an incisor made of absorptive special ceramic with 3 indentations in the plastic container.

# 3 indentations (approx. 30 x 20 x 4.5 mm) Dimensions - plate: approx. 190 x 144 x 13 mm Dimensions - container: approx. 192 x 146 x 25 mm Weight: approx. 580 g REF 390 0040 2



Ceramico WetOne14
Self-wetting mixing plate in the shape of an incisor made of absorptive special ceramic with 14 indentations in the plastic container.

3 indentations
(approx. 43 x 30 x 5 mm)
4 indentations
(approx. 33 x 22 x 4.5 mm)
7 indentations
(approx. 29 x 20 x 4.5 mm)
Dimensions - plate:
approx. 190 x 144 x 13 mm
Dimensions - container:
approx. 192 x 146 x 25 mm
Weight: approx. 530 g
REF 390 0040 3

### Funktional thanks to prolonged processing consistency

The photos below show how moisture is withdrawn from the mixing plate by the ceramic material and the material takes on a darker shade. If the brush is pressed lightly into the material, the moisture level is increased and the material can be perfectly taken up for further processing.











Use a soft brush to clean the plates under running water. Alternatively, a steam blasting unit can be used. Better cleaning results can be achieved by treatment of the previously dried plate in the preheating furnace. The plate should be supported on firing trays and slowly heated up to 100 to 150° C; the temperature should be maintained for 1 hour. Then the plate is slowly heated up to 800° C. Finally, the furnace is turned off and the plate is slowly cooled down.

If necessary, finish the surface with waterproof abrasive paper and clean once more under running water.

Attention! If the plate is heated and cooled too quickly, a temperature shock – and consequently fracture – may result.

### Unconventional

### Ceramico

#### BlackMolar and ColorImplant

"Form follows function" - this is part of a famous slogan created by Louis Sullivan, the American architect and main representative of the Chicago School, one of the first leading-edge high-rise architects. The unconventional and unique shapes of Ceramico BlackMolar and ColorImplant are to set new trends for laboratory routines of dental technicians and to appreciate their "works of art".



### Ceramico BlackMolar Black-glazed mixing plate

in the shape of a molar with overglazed nerve canals and a lid impenetrable to UV light for mixing light-curing materials.

5 indentations
(approx. 17 x 12 x 2.5 mm)
1 indentation
(approx. 30 x 22 x 4.5 mm)
Dimensions - plate:
approx. 210 x 115 x 11 mm
Dimensions - lid:
approx. 210 x 117 x 17 mm
Weight: approx. 300 g
REF 390 0040 4




### Ceramico ColorImplant

White-glazed mixing plate in the shape of an implant with 13 indentations and dust lid, for stains.

1 indentation
(approx. 42 x 23 x 4.5 mm)
12 indentations
(approx. 14 x 9 x 3 mm)
Dimensions - plate:
approx. 170 x 67 x 12 mm
Dimensions - lid:
approx. 172 x 69 x 16 mm
Weight:
approx. 190 g
REF 390 0040 5

### breCeram





### for the ceramics specialist

- inverted cone with relief grinding technology for smooth surfaces
- fine but abrasive diamond grinder
- two different abrasion stages of the Abraso-Fix-Roundbrushes permit rapid polishing, as the polish paste is already carried in the bristles
- Ceragum coarse is suitable for use with ceramic and for rubberized metal
- Cerafine adds a very high gloss to ceramic and metal very quickly



The diamond grinder is used for coarse processing of ceramic. But nevertheless a fine surface is achieved thanks to the fine diamonds.



The inverted cone is ideally suited for shaping of the occlusal surfaces. Simultaneously the undercut produces a polished ceramic surface.



Ceragum coarse is a universal product. It removes material fast, leaving an optimal surface structure.



Abraso-Fix green is used to produce the coarse surface structure. It is also highly suitable for smoothing ceramic and metal occlusal surfaces.



Abraso-Fix red already achieves a light polish on the surface. It is used to give the first polish to ceramic and metal.



Cerafine is the high gloss polisher for ceramic and metal. It is particularly suitable for transition areas of metal to ceramic, as after glaze firing the metal polish no longer has a matting effect.



Fast and simple processing of ceramic and metal. breCeram offers the ideal combination.





breCeram processing set for ceramics

- 1 diamond grinder fine 1 Tungsten carbide 1.2
- 1 Abraso-Fix green
- 1 Abraso-Fix red
- 1 Ceragum coarse, wheel
- 1 Cerafine, wheel REF 520 2028 6



Diamond grinder fine 1 piece REF 340 0107 1



Abraso-Fix green 2 pieces REF 350 0059 0 8 pieces REF 350 0075 5



500 104 010006 012 1 piece **REF H010 NH 12** 

Tungsten carbide

ISO-Nr.



Abraso-Fix red 2 pieces REF 350 0060 0 8 pieces REF 350 0075 3



Ceragum coarse, wheel not mounted 12 pieces REF PRK G221 2 50 pieces REF PRK G225 0 100 pieces

REF PRK G220 0



Cerafine, wheel 1 piece REF 520 2028 5

#### Accessories:



Tungsten carbide ISO-No. 500 104 010006 008 1 piece **REF H010 NH 08** 



Tungsten carbide ISO-No. 500 104 010006 010 1 piece **REF H010 NH 10** 



Tungsten carbide ISO-No. 500 104 010006 016 1 piece **REF H010 NH 16** 

# Manual processing of zirconium


# compoForm UV



compoForm UV 2 x 3 ml syringes 10 nozzles REF 540 0115 0



Pin-point application directly from a sysringe for quick and precise modeling.



- stable allows safe modeling to match the demands of profile milling
- low shrinkage for high precision, custom made modelations.
- millable for ease of trimming
- easy syringe application reduces modeling time



Asseccories:

Nozzles 25 pieces REF 580 0001 8

bre.lux Power Unit

REF 140 0097 0



compoForm UV easy trimmable with tungsten

### SERACOLL UV



SERACOLL UV 2 x 3 ml 2 dipping bowels REF 540 0115 1



Asseccories:

Fast and simple – making individual

abutments from

compoForm UV.

bre.lux Power Unit REF 140 0097 0

Light curing wax glue with high grade of capillary attraction means tension free constructions.

- safety by high stability
- speedy work flow by short curing time
- smooth surface of modelation equals smooth zircon surface
- high capillary attraction safe joining of bridge segments
- shrinkage free ideal for custom made construc-



Tension free, stable and safe joining of bridge segments for high precision modelations.




Gaples closure between bridge units through high capillary effect. More safety by applying this procedure.



SERACOLL UV smoothes the surface means smooth probe scanning of the structure.




### Vario-Soft 3 zircon sv



Vario-Soft 3 mini sv females green 8 pieces REF 430 0733 5

females yellow 8 pieces REF 430 0733 3

females red 8 pieces REF 430 0733 1



Vario-Soft 3 zircon sv 8-pieces pack REF 430 0732 2

Vario-Soft 3 zircon sv mini sv 8-pieces pack REF 430 0732 9

The approved attachment system Vario-Soft 3 is now also available for zircon. The joint and integrated stress compensator, especially designed for zircon specifications, clearly indicates Vario-Soft 3 zircon sv for combined denture work. Two sizes available for anterior and posterior areas allow application in low space regions

- approved attachment system with three grades of friction for gaining more individuality
- · solid joints between attachment and crown for more safety
- individually reduceable for low vertical space regions



Joining the components with compoForm UV Completion of the joint with SERACOLL UV to be light cured.



Mirror polishing of the attachment after sintering with Zi-polish. The secondary framework to be processed as usual.

#### Asseccories:

Parallel holder universal REF 360 0115 1

### 5-motions-zircon



For denture applications, the cleanliness of the material ranks first and foremost.

Bi-axial pressed zircon for a steady quality - crosscharge. Different sizes offer full flexibility and high efficiency.



5-motions-zircon zircon blanks small 6/16, 16 mm height 3 pieces REF 360 1061 6

6/20, 20 mm height 3 pieces

REF 360 1062 0



5-motions-zircon zircon blanks medium

10/16, 16 mm height 2 pieces

REF 360 1101 6

10/20, 20 mm height 2 pieces

REF 360 1102 0



5-motions-zircon zircon blanks large

14/16, 16 mm height 1 piece

REF 360 1141 6 14/20, 20 mm height 1 piece

REF 360 1142 0



Asseccories:

5-motions-colors zircon coloring liquid each 100 ml

A1 REF 360 1CA1 0 A3 REF 360 1CA3 0

B1 REF 360 1CB1 0 B3 REF 360 1CB3 0

C3 REF 360 1CC3 0

# Zi-polish - the perfect mirror polish!



Zi-polish

REF 360 1002 5

Zircon polishing paste for high gloss surfaces. Two different diamond grain sizes for pre- and high gloss polishing of marginal edges, basal areas, primary crowns, individual abutments and bars. The special binding agent avoids spraying of the polishing past when dipping the brush into the paste container. This also reduces paste consumption.

- the high content of diamond particles reduces polishing time and leaves a highly glossy surface
- special binding agent reduces paste consumption
- proper wipe off of pre polish residues from dents



Due to a special binding agent the paste sticks to the brush, saving material also.



Perfect and speedy mirror polish is achieved by high diamond particle density.



# Manual processing of zirconium

# 5-motions-glue and 5-motions-active



5-motions-active activator spray for 5-motions-glue 200 ml REF 360 1002 7

Tough-flowing super glue for clean and safe connection of blanks and modelation. Activator-spray for 5-motions-glue provides safe joint between modelation and blank and to other materials.

- non-dripping and clean application
- quick setting with 5-motions-active activator spray saves time
- bonds all materials to one another and is hence universally suitable



5-motions-glue super glue gel 20 g REF 360 1002 6



Pin point fixation with 5-motions-glue. Quick and safe connection by using 5 motions-active spray.



Quick frame mounting of the zircon blanc. By using the activator spray different materials are joint together. This widens the field of indication.


### Scanners and grinders

Multiple fields of applications in dental technique require special geometries for precise zircon processing.

The approved relief grinding technique of bredent created special burs for precise cutting and smooth surfaces.

- applicable for all units using a magnification factor of 1:1,25
- Relief grinding technology stands for smooth surfaces and reduced retouching work.
- 11 different designs for a wide range of applications
- special tool steel stands for long life span and endurance

### Relief grinding technology means smooth surfaces avoiding material ruptures:




# Tripled life-span and endurance compared to conventional bredent-gearing.

In a special process, bredent tungsten carbide burs of the latest development state are provided with relief grinding properties.

This stabilises the cutting edges against ruptural destruction. This also enhances the endurance of those tools by the factor three compared to conventional tools for green state processing.

In addition the relief grinding design optimizes the chip angle for immaculate cutting performance.






# Scanners and grinders

Relief grinding technology for smooth surfaces, avoiding surface ruptures:





Harmonized scanners and burs in various shapes and sizes covering all aspects of processing. The solution for simple and succesful zirconium work. Relief grinding geometry of burs provides an immaculate balance, a clean and smooth surfaced grinding pattern, avoiding ruptures of cutting edges.





All burs and scanners have a spindle diameter of 3 mm and are suitable for units with an enlarge factor of 1:1,25.



Preliminary work by 5-motions-milling cutter 2, speedy and smooth performance.



Deep grindings like with abutments become effortless with long shaped burs.



The right angled, flat headed burs provide a good seating surface on implant platforms.




Drilling and surface work is carried out precisely with the Rapidy 2.



# Manual processing of zirconium

### Diagen-Turbo-Grinder



Highest grinding performance and abrasivity on zircon, ceramic and metal surfaces at low contact pressure.

Improved endurance compared to regular bindings effectively allow a wide field of applications.

The diamond grinder system with extraordinary cutting properties through special Diagen-Diamant-binding. Highly abrasive grinding performance on zircon surfaces. Applying only light pressure for gentle surface processing.

- two grades of abrasion offers more flexibility of use
- special binding (rough) equals 20 % more endurance and saves on tool consumption
- Cooling effect to avoid damages of zircon
- recommended by leading zircon- and ceramic manufacturers

### Zircon grinding patterns



New shapes in 2 grades of abrasion. Boosts trimming through rough diamond grain sizes and special Diagen binding. Surface roughness: 26 µm.



New shapes in wellproven quality. Surface smoothening and fine cut by fine grain diamonds. Surface roughness: 2 µm.

Diagen-Turbo-Grinder, rough New quality of diagen binding with large diamond grain size for speedy trimming.



Cylinder pointy Ø 3.5 x 11 mm 2 pieces REF 340 G015 5




Cylinder rounded Ø 4.5 x 13 mm 2 pieces REF 340 G016 5



**Cylinder** Ø 4.5 x 13 mm, 2 pieces **REF 340 G016 0** 



**Cylinder** Ø 6.0 x 13 mm 2 pieces **REF 340 G017 0** 



Cone Ø 3.5 x 11 mm 2 pieces REF 340 G015 0



Lentil Ø 22 x 2 mm 1 piece REF 340 G021 0

Diagen-Turbo-Grinder, fine Established fine grain quality for perfect surface smoothness.



Cylinder pointy Ø 3.5 x 11 mm 2 pieces REF 340 0015 5



Cylinder rounded Ø 4.5 x 13 mm 2 pieces REF 340 0016 5



Cone Ø 3.5 x 11 mm 2 pieces REF 340 0015 0



**Cylinder**Ø 4.8 x 13 mm
2 pieces **REF 340 0016 0** 



**Cylinder**Ø 6.5 x 13 mm
2 pieces **REF 340 0017 0** 



Top down cone with notch Ø 6 x 8 mm 1 piece REF 340 0025 0



Top down cone with notch Ø 12 x 6 mm 1 piece REF 340 0024 0



Lentil Ø 22 x 2 mm 1 piece REF 340 0021 0



Wheel Ø 22 x 2 mm 1 piece REF 340 0022 0



Wheel
Ø 15 x 3,5 mm
2 pieces
REF 340 0018 0



Wheel Ø 22 x 4,5 mm 1 piece REF 340 0019 0



## visio.lign veneering system / Teeth / Denture resins

Resins and composites are high-tech materials used in modern dentistry. They have a long tradition at bredent and have become established components of complex system solutions for all prosthetic indications. We pay special attention to coordinating material properties, material compound, modelling and handling as well as shading and polishing.



ulala Ilau Manaaulua arakau

Lab. Od. Lazetera Antonio - Savona - Italy, Dott. Vescia Luca -Villa Dossola - Italy

visio.ligii velicerilig systelli	
visio.lign	
Indications - application areas	289
visio.lign system components	290
Processing steps	292
Indications of the visio.lign primers	292
Use of the system components	
novo.lign A veneers	294
Thermo-Pen	
novo.lign P veneers	
Indications of the visio.lign primers	
MKZ Primer and MKZ EM-Aktivator	296
combo.lign Opaquer	296
visio.link	
combo.lign luting composite	297
K-Primer	
Bond strength tests	298
crea.lign veneering material	299
haptosil D	301
visio.sil	301
visio.sil ILT	301
visio.sil fix	
Location matrix drill	302
visio.lign sets	
Red-white esthetics set	
visio.lign tooth cabinet	
visio.lign set 2	
visio.lign set 3	
novo.lign A und P veneers	
neo.lign A ans P (full) teeth	304
Opaquer Set	305
crea.lign Sets	305
visio.lign Toolkit	305
Light-curing units and times	
Suitable light-curing units	306


	bre.Lux polymerization timesbre.Lux Power Unit	
2:	ange of design novo.lign veneers	
10	novo.lign A, upper anterior	300
	novo.lign A, upper/lower anterior	300
	novo.lign P, upper/lower posterior G-design	
	novo.lign P, upper/lower posterior W-design	.311
) (	esign chart neo.lign (full) teeth	
	neo.lign A, upper anterior	312
	neo.lign A, upper/lower anterior	
	neo.lign P, upper/lower posterior, G-design	314
	neo.lign P, upper/lower posterior, L-design	315
P	ease order from your dealer	316
) (	enture resin	
	HIGH.lign	321
	uni.lign	
	uni.lign speed	
	Overview of the uni.lign resins	
	top.lign professional	326
21	eaning	
	Dentaclean impression desinfectant	
	Shipping bags	
	Dentaclean denture cleaning agent	
	Dentaclean ultrasonic cleaning agent	328
	Dentaclean plaster removing agent /	
	Dentaclean plaster removing agent Speed	
	Dentaclean mixing fluid for pumice powder	328
Jı	nits / Instruments	
	Ergonom wax knife	329
	Repositioning tweezers	330
	Thermo-syringe	330
	Posi-boy	
	Articulation paper holder	331
n	sulating agents	
	Wax insulating agent	332
	Plaster insulating agent	332
	Isoplast ip	333
٨	ax	
	Bite blocks	334
	Set-up wax asw	335
	Modelling wax pink Standard mdwst	
	Wax palatal patterns gf	336
26	esins and systems	
	Tray material UV	337
	casting system	
	Multisil-Soft	340
) <sub> </sub>	paker	
	Ropak UV	
	Ropak Kompaktopaker UV	341
	Kompaktopaker tooth-colored UV	341
r	ocessing acrylics	
	Abraso-Gum Acryl	342
	Set-up grinding tool	
	Diacryl grinding tool	
	Tungsten carbide tools	
'n	lishing	
	Acrylic polishing set	345
	High luster polishing with the handpiece	345
	Polishing brushes	
	Daliahina huff	247

## visio.lign veneering system

## visio.lign

The visio.lign system is based on multi-layer veneers (novo.lign) developed in accordance with the example of natural teeth. A bonding system with matched shades and suitable individualization and add-on materials enable the perfect combination of individual and efficient veneering techniques. With abrasion values comparable to natural teeth, visio.lign provides gentle occlusion – an advantage especially for implant-supported restorations. Patients can have a preview of the result thanks to the esthetic try-in with the novo.lign veneers. By using the sili-


cone key based on the try-in, perfect design of the framework can be easily achieved. The high-strength novo.lign veneers can be bonded to all framework materials used in dental techniques and individualized with crea.lign veneering material.














## Indications - Application areas:



Horizontally screw-retained visio.lign bridge on SKY Unicone with reduced number of implants (LI/UJ).

Photo: MDT Oliver Heinzmann



Complex case on SKY fast & fixed, layered using the Inverse Layering Technique with crea.lign, upper and lower immediate definitive restorations.

Photo: Vincenzo Musella



Horizontally screw-retained visio.lign bridge, ZrO2 framework on individual SKY uni.fit ZrO2 abutments.

Photo: DT Andreas Lüdtke



Non-precious metal bar-lock restoration combined with visio.lign and neo.lign teeth.



Fixed bridge made of Bio HPP, partially veneered with visio. lign.



Non-precious veneer bridge, individualized with crea.lign using the cut-back technique.



Onlay- Overlay – Inlay layered with crea.lign using the Inverse Layering Technique.

Photo: Vincenzo Musella



Definitive lab-made crea.lign veneers, adhesively cemented using the non-prep technique.

Photos: Vincenzo Musella



Removable NPM veneer bridge on telescopic crowns, combined with neo.lign denture teeth and novo.lign veneers.



Immediate temporary SKY fast & fixed restoration with novo.lign and top.lign breformance framework resin.

MDT Oliver Heinzmann



## visio.lign veneering system

## visio.lign system components



#### novo.lign veneers

(sets of anterior and posterior teeth)
Multi-layer veneers made of high-impact PMMA
composite, anatomical designs of anterior
and posterior teeth: natural esthetics for all
indications. Available in the classic A-D shades.



#### neo.lign denture teeth

(sets of anterior and posterior teeth)
Anatomical denture teeth with natural shade gradient and designs for implant and fixed/removable restorations made of high-impact PMMA composite. The cross-system occlusal design was developed for all common occlusion concepts.

Available in the classic A-D shades.



#### combo.ligr

Light- and dual-curing luting composite for reliable shade reproductin and superior bonding of the novo. lign veneer to all framework materials. Available in the classic A-D shades.



#### combo.lign Opaquer

Light- and dual-curing opaque material which ensures perfect curing. Bonding system tested at universities, compatible with all metal primers and silanizing techniques; perfectly suitable when using mechanical retentions. Available in three tooth and one gingiva color to ensure color stability in all A-D shades in combination with combo.lign.



#### visio.link

Universal, light-curing PMMA and composite primer. visio.link is used for superior bonding of PMMA denture resins/composites, high-impact PMMA composite materials (novo.lign veneers, neo.lign denture/acrylic teeth) and combo.lign.



#### crea.lign veneering material

- with nano-ceramic fillers only (does not contain dental glass)
- superior polishing properties to produce perfect surfaces, Ra 0.03
- flowable, improved coating transition to the veneer without gaps
- low water absorption, ensure mechanical stability
- no accumulation of plaque
- available in the classic A-D dentin shades, incisal, GLIM and stains
- light-curing crea.lign opaquer available in new tooth shades and a gingiva color for "free-hand layering"



## visio.lign system components



#### crea.lign Modelling Liquid

crea.lign Modelling Liquid can be used to dilute crea.lign and to increase the flowability; it also allows to obtain a "homogeneous" transition from the denture resin to crea.lign during red-white individualization.

Just like crea.lign, crea.lign Modelling Liquid is resistant to discoloration and plaque. When mixing with crea.lign, the maximum proportion is 30 %.



#### K-Primer

K-Primer is used to achieve adhesion of veneering ceramic and press ceramic, e.max press as well as silicate ceramic, such as CAD-Blanks, Mark II, to composites such as crea.lign. K-Primer is not suitable for direct use and can be used for repairing ceramic veneers.



#### MKZ Primer

MKZ Primer is used for conditioning all metal NPM alloys and titanium, zirconium dioxide, aluminium oxide and spinell ceramic. This way chemical bonding to composites such as crea.lign is achieved. Accordingly, MKZ Primer is perfectly suited for luting individual abutments.



#### MKZ EM-Aktivator

MKZ EM-Aktivator can only be used in combination with MKZ Primer - in a ratio of 1:1 for conditioning precious metal frameworks (Au, Ag, Pt, Pd) and to achieve chemical bonding to composites such as crea.lign.



#### visio.sil (transparent silicone for keys)

Transparent, addition-curing silicone for keys; hardness: approx. 60 Shore A; for the fabrication of translucent keys and for using and processing light-curing materials.



#### visio.sil ILT (transparent silicone for keys)

visio.sil ILT was developed especially for the Inverse Layering Technique. This transparent, addition-curing silicone for keys features a high hardness of approx. 75 Shore A. It is flowable to fill the approximal space and to ensure perfect reproduction of the surface texture; however, it is firm and stable and can be applied very precisely. visio.sil ILT must not be applied in the patient's mouth (class I medical device) in the hardened condition.



#### visio.sil fix (high-precision silicone for keys)

visio.sil fix is an addition-curing, high-precision silicone for keys. The use of visio.sil fix enables perfect and detailed reproduction of approximal spaces and the gingival situation of the wax setup. Only high-gloss polishing needs to be carried out after fabricating the restoration so that a lot of time is saved and the amount of work is reduced.




#### haptosil D (kneading silicone)

haptosil D is an addition-curing kneading silicone with a hardness of 90 Shore A for the fabrication of exact and stable keys. haptosil D is mostly used in combination with the other visio.sil silicones to add more stability to these materials



## Versatile polymerization unit for practice and laboratory

The bre.Lux Power Unit was developed for processing veneering and dental materials. This light-curing unit features a wavelength range of 370 nm to 500 nm. The LEDs have a life of 20,000 hours.



# The visio.lign Toolkit has been optimized for processing composite and visio.lign veneers and ensures a perfect finish. The combination of the materials of the visio.lign

veneering system and the visio.lign Toolkit results in surfaces featuring plaque resistance and shade stability and their quality and resistance is identical to the ones of a ceramic material.



## visio.lign veneering system

## **Processing steps**

#### Esthetic try-in



Selecting the suitable color and design.



If required, the novo.lign veneer is milled thinly in the cervical area.



Fixation of the veneer with tooth-colored wax, beauty setup wax.



Completed esthetic setup.




Try-in, control and correction if required



Fixing the esthetic setup with visio.sil (transparent silicone for keys).

#### Modelling the framework



Perfect modelling of the framework is achieved through the use of the key.



The perfectly designed metal framework

#### Conditioning the framework



Conditioning with MKZ Primer.



Apply opaquer and polymerize.



## **Processing steps**

#### **Bonding**



Veneers were sandblasted with 110 µm aluminium oxide.



Apply visio.link thinly and cure with light.



Apply combo.lign to the inner side of the veneer.



Spread excess material and cure with light subsequently.

#### Finish with crea.lign



After bonding, crea.lign is applied in the approximal area.



Contouring the palatal area with crea.lign.



Individual design of the redwhite area.



A hand lamp can be used for intermediate curing of crea.lign.

#### Finishing



Prepolishing with Abraso-Fix brush and Acrypol.



High-gloss polishing with cotton buff and Abraso-Starglanz.



## novo.lign A veneers



## novo.lign A (anterior) veneer

Layering and design of the novo.lign veneers are based on natural teeth. Despite the very thin and slender design of the veneers, the characteristics, such as mamelons and buildup, could be integrated into the veneers.

novo.lign A veneer:

- facet thickness in the cervical and central areas:
- 10 upper designs
- 3 lower designs



All designs of the novo.lign A veneers are available as sets of 6 pieces, sets of incisors and sets of cuspids.

As a result, veneers, for example as a cuspid set, can be combined with the neo.lign incisors (no. 4) depending on the respective case.

Moreover, the cuspid (no. 2) veneers of large designs, for example 147, can be combined with the incisor veneers of other designs, for example 145.

#### Thermo-Pen



#### Thermo-Pen

is a hot-air device with piezo technology without open flame. The veneers can be heated without the risk of burning them. The temperature required to shape the veneers is approx. 250° C. Improved coating of e.g. telescopic crowns is achieved through expansion.




The novo.lign is heated with the Thermo-Pen.






The hot veneer is expanded.

## novo.lign P veneers, G-design







novo.lign P (posterior) multi-functional veneers, G-design

The novo.lign P (posterior) veneers are distinguished by their multi-functional occlusal design which enables the use for all occlusion concepts. The various layers of natural teeth were copied for the novo. lign veneers.

novo.lign P multi-functional veneers:

- facet thickness in the cervical and central areas: 1.2 mm
- 2 sizes G3 and G4 for upper and lower jaws



## novo.lign P veneers, W-design



novo.lign P (posterior) veneer for crown and bridge restorations, W-design The novo.lign P (posterior) veneer, W-design, was developed especially for crown and bridge restorations. These veneers are more voluminous, larger and allow coating large tooth stumps. The layering structure could be maintained despite the layer thickness which is lower than the ones of the G3 and G4 designs.

novo.lign P veneer for crown and bridge restorations:

- facet thickness in the cervical and central areas:
- 3 sizes W3, W4 and W5 for upper and lower jaws

## neo.lign A denture teeth




neo.lign A (anterior) denture tooth The development of visio.lign is based on natural teeth. The lifelike design of the novo.lign veneers was transferred to the neo.lign denture teeth. The shade, design and layering structure of neo.lign anterior and posterior are key aspects within the visio. lign veneering system. Perfect harmony for partial restorations, attachment techniques and implant prosthetics.

neo.lign A denture tooth:

- 9 upper anterior sets
- 3 lower anterior sets

## neo.lign P denture teeth, G-design



neo.lign P (posterior) multi-functional denture tooth, G-design The neo.lign P (posterior) denture tooth is based on the novo.lign P veneer. The multi-functional occlusal design was transferred to the neo.lign P denture teeth and enables the use for all familiar occlusion concepts.

neo.lign P denture tooth, G-design:

- 3 sizes - G2, G3 and G4 - for upper and lower jaws

## neo.lign P denture teeth, L-design



neo.lign P (posterior) denture tooth for lingualized setup, L-design

The neo.lign P (posterior) denture tooth, L-design, was developed especially for lingualized setup techniques. The L-design allows a lingualized setup without grinding adjustments.

neo.lign P denture tooth, L-design:

- 3 sizes - L2, L3, L4 - for upper and lower jaws



## Indications of visio.lign primers



The bond strength values achieved for all primers listed above are clearly high than the minimum value of 5 MPa (DIN EN ISO 10477). See chart, page 298.

#### MKZ Primer and MKZ EM-Aktivator



MKZ Primer 4 ml REF MKZ02004

MKZ EM-Aktivator 4 ml REF MKZEM004

## MKZ Primer and MKZ EM-Aktivator - conditioning of the framework.

The finished metal frameworks (NPM/precious metal-free/CoCr/precious metal) need to be sandblasted at 3-4 bar (zirconium frameworks: max. 2 bar) using 110 µm aluminium oxide. The frameworks must not be steam-cleaned with a steam jet. Use alcohol and a clean brush for cleaning. When blowing dry with compressed air, make sure that the compressed air is free from oil. Use a brush to apply MKZ Primer to the framework and leave to dry completely; then apply the desired opaquer.

When conditioning precious metal/reduced precious metal content frameworks, MKZ Primer and MKZ EM-Aktivator need to be mixed in the ratio of 1:1 and applied subsequently.



Apply MKZ Primer / MKZ EM-Aktivator and leave to dry. Then the opaquer can be applied.

## combo.lign Opaquer



combo.lign Opaquer can be used individually – either as a light-curing or a dual-curing material – by mixing the catalyst with the base paste in a ratio of 1:1. For mechanical retentions, we recommend dual-curing to ensure adequate curing of the opaquer.

REF see Order Form.



Use a suitable brush (opaquer brush) to apply combo.lign Opaquer.



Light-curing of combo.lign Opaquer. Polymerization for 180 sec is recommended after each layer.



#### visio.link



visio.link 10 ml REF VLPMMA10

#### visio.link - conditioning of veneers.

The veneers need to be sandblasted at 2 to 3 bar using 110  $\mu$ m aluminium oxide. Do not use a steam jet for cleaning the veneer. A thin coat of visio.link is applied to the novo.lign veneers; make sure to apply a very thin coat since excessively thick coats may affect the bond strength. After applying visio. link and polymerizing for 90 sec (lightwave range: 370–400 nm), the surface should be dry and exhibit a silky mat gloss.



If required, the novo.lign veneer is milled thinly in the cervical area.



Sandblast novo.lign veneer at 2-3 bar using 110 µm aluminium oxide.



Apply visio.link thinly and polymerize for 90 sec.



A silky mat gloss indicates that visio.link has been applied correctly. Excessive application of visio.link must be avoided.

## combo.lign luting composite



The dual-curing combo.lign hardens chemically and by exposure to light. To achieve the final hardness, combo.lign needs to be polymerized in the light-curing unit for 180 sec. An adhesive gap of 40 µm – 2 mm can be filled with combo.lign. combo.lign should always be covered with crea.lign because it is difficult to process and polish. **REF see Order Form.** 



Applying combo.lign into the conditioned veneer.



Spread or remove excess material.



Harden each tooth for 10 sec., remove the key and use light-curing unit for final polymerization (180 sec.).

#### K-Primer



K-Primer 3 ml REF APK25003

## **K-Primer – conditioning of veneering ceramic.** The veneering ceramic needs to be sandblasted

The veneering ceramic needs to be sandblasted at max. 2 bar using 110 µm aluminium oxide or roughened using a diamond bur (dry). After sandblasting/roughening, the framework must not be cleaned with a steam jet! Any contamination needs to be removed with alcohol and a clean brush. The veneering ceramic is bonded to crea.lign using K-Primer. K-Primer and crea.lign are suitable for direct and indirect use. Damage to ceramic veneers caused by chipping can be quickly and easily repaired.



Sandblast "chipping area" at a pressure of 2 bar using 110 µm aluminum oxide or roughen with diamond abrasive tool (dry). Then K-Primer is generously applied; let it dry for approx. 30 sec.



After applying the suitable crea.lign material and subsequent polishing, no transition is visible at the repaired spot.




## Bond strength tests

Bond strength according to DIN EN ISO 10477 / University of Jena Minimum value: 5 MPa




Bond strength – combo.lign (luting composite) vs veneering composite, MKZ Primer




Bond strength – combo.lign (luting composite) vs veneering composite, MKZ Primer + MKZ II Primer




Bond strength – combo.lign (luting composite) vs veneering composite K-Primer



Bond strength - composite (combo.lign/crea.lign) vs resins visio.link



Bond strength - combo.lign (luting composite) vs competitors



## crea.lign veneering material



The final design/shape of the veneer is prepared with crea.lign. Without intermediate curing, the layer thickness of crea.lign should not exceed 1 mm. crea.lign may be used in situ, e.g. to repair ceramic veneers.

crea.lign can be applied from the syringe or using a brush.

The use of crea.lign facilitates the application

with the brush and optimizes modelling of the approximal spaces. All crea.lign materials can be mixed with one another. Various incisal, neck, dentin and gingiva materials are available for individualization purposes. The purely light-curing crea.lign opaquer is used to cover the frameworks in the suitable tooth shade.

REF see Order Form.



crea.lign Modelling Liquid 10ml REF CLFMOD10

#### Applying crea.lign:



After luting the novo. lign veneers, crea.lign is applied to the approximal and distal areas.



The tooth shape is built up in the palatal area using crea.lign.



crea.lign GUM materials are used for individual contouring in the redwhite area.



Intermediate polymerization is required each time material is applied. A hand lamp may be used. Final polymerization is carried out in the bre.lux Power Unit (light-curing unit).

#### Finishing crea.lign:



Finishing with a tungsten carbide bur.



Prepolishing with Abraso-Fix brush and Acrypol prepolishing paste.



Prepolishing with a goathair brush and Acrypol prepolishing paste.



High-gloss polishing with cotton buff and Abraso Starglanz polishing paste.


#### Accessories:



MagicBrush size 2 2 pieces REF 390 M002 0



bre.Lux Power Unit REF 140 0097 0



Composite finishing set incl. Acrypol polishing paste Abraso-Starglanz high-gloss polishing paste REF VLTOOLKIT



## crea.lign veneering material

#### Inverse Layering Technique Fabrication of veneers (Non-prep technique) with crea.lign composite



The prepared master model with removable dies.



Morphological restoration of the anterior teeth (wax).



Fabrication of the silicone key with visio.sil ILT (separated dies).



A brush is used to apply E2 incisal material into the key.



Light-curing for approx. 5 sec after each layer.



Application of the effect materials to create contrasts. Light-curing: approx. 5 sec.



Opalescent effects with Incisal opal and Incisal blue materials; light-curing: approx. 5 sec.



Application of dentin, then the key is placed on the model for lightcuring.



Final curing: 360 sec. in the bre.Lux Power Unit.




Excellent translucency of the completed crea.lign veneers.



Initial situation.



Veneers in situ after definitive bonding.





## haptosil D





haptosil D Components A and B 1300 g each REF 540 0118 0 7500 g each REF 540 0119 0

# haptosil D addition-curing kneading silicone With a hardness of 90 Shore A, haptosil D is perfectly suited for the fabrication of precise and stable keys. Models for repair and extension can also be produced within a short time. As a result, considerable time can be saved compared to conventional model fabrication techniques.

#### visio.sil



visio.sil 50 ml REF 540 0120 0 visio.sil mixing cannulas 12 pieces REF 320 0045 7

#### visio.sil transparent silicone for keys

Translucent keys are fabricated with visio.sil (hardness: approx. 60 Shore A); these keys are the bases for luting the novo.lign veneers to the respective framework. To add more stability to a visio.sil key, the additional use of haptosil D (hardness: approx. 90 Shore A) silicone for keys is recommended.



The haptosil D key is shaped similar to a picture frame.



The key is built up with visio.sil.



The key is smoothed with a finger dipped in detergent.



Light can penetrate the key, which ensures reliable bonding/luting of the veneers.

#### visio.sil ILT



visio.sil ILT 50 ml REF 540 0140 0 visio.sil mixing cannulas 12 pieces REF 320 0045 7

visio.sil ILT (Inverse Layering Technique) transparent silicone for keys visio.sil ILT (hardness: approx. 75 Shore A) was developed for the Inverse Layering and the flasking technique. This technique is used, for example for the fabrication of veneers. The veneer is layering directly into the key in reverse order. The visio.sil ILT key must not be placed into the patient's mouth for luting the veneers.



Initial situation



Plaster master model



Wax up on the master model ready for fabrication of the key.



The key is placed into the mouth (crea.lign is still soft).



crea.lign is polymerized with the hand lamp.



Completed esthetic preview of the additional veneers in situ.

Technique and photos: Vincenzo Musella.



9 visio.lign veneering system / Teeth / Denture resins

#### visio.sil fix



visio.sil fix 50 ml REF 540 0130 0

#### visio.sil fix high-precision silicone for keys

visio.sil fix can be used in combination with other silicones for keys, such as visio.sil, visio.sil ILT or haptosil D .

visio.sil fix provides perfect reproduction of approximal spaces/surfaces so that only high-gloss polishing is required as a final step. The suction effect of visio.sil fix allows to avoid cyanoacrylate adhesive to fix the veneers in the key; as a result, the surface structure of the veneers is retained.

#### visio.sil fix in combination with haptosil D



visio.sil fix is applied to the setup.



haptosil D is mixed and pressed to visio.sil that is still soft.



Holes are drilled into the key using the matrix drill.



Light penetrates through the drill hole and the veneer.

#### visio.sil fix in combination with visio.sil



visio.sil fix is applied to the setup.



visio.sil is applied to visio. sil fix that is still soft.



The precisely molded approximal spaces and margins can be clearly recognized.




Light penetrates through the key.

#### Location matrix drill



Location matrix drill 1 piece REF 330 0078 0


## No adhesive – no displacement of the teeth in the matrix! conventional method with matrix drill



undesired gap



utmost precision of fit



## The safest and most accurate method for reliable positioning of acrylic teeth in plaster or silicone matrixes.

With the conventional method acrylic teeth are fixed incisally/occlusally in the matrix using adhesive wax. This includes the risk that the tooth is pulled out of the matrix due to the contraction of the wax and an undesired gap results.

The holes that are drilled through the matrix allow the application of adhesive wax onto the acrylic teeth. Due to the contraction of the wax the tooth is pulled into the matrix and held safely.

#### Red-white esthetics set



- nanofilled crea.lign provides shade stability and unsurpassed plaque resistance.
- visio.link ensures permanent bonding of crea.lign to all PMMA denture resins; tested at universities.

#### Assortment

28 pieces

#### Red-white esthetics set

- 1 x aluminium folder
- 8 x crea.lign
- 1 x GUM Opaquer
- 1 x visio.link
- 1 x MKZ Primer
- 1 x crea.lign Modelling Liquid
- 2 x demo denture
- 2 x MagicBrush, size 2
- 1 x layering instructions
- 1 x single-hand grip
- 1 x mixing pad
- 8 x application cannulas

REF CLIGNSETG

- the fully cross-linked uni.lign denture resin enables perfect chemical bonding to crea.lign and avoids visible transitions and adjustments.
- since crealign does not contain dental glass, a surface quality identical to denture resins (Ra value:  $0.03~\mu m$ ) is achieved to guarantee fast and simple polishing.
- high gloss is achieved in only 3 minutes using the same polishing procedure as for denture resin
- bridges veneered with ceramic can also be individualized with the red-white esthetics set.
- individual red-white esthetics in just 6 steps.
- a dental technician needs only approx. 45 minutes to individualize anterior restorations.
- this set contains all necessary materials including step by step layering instructions.
- elegant aluminium folder for presentation to patients and dentists.
- two demo dentures are included to the set; one of them can be individualized to demonstrate patients and dentists the particular advantages.

## visio.lign tooth cabinet



#### Assortment

10 pieces

#### visio.lign tooth cabinet

- 1 x visio.lign tooth cabinet
- 1 x drawer insert, crea.lign
- 8 x drawer insert, novo.lign A

REF VLIGNBOX1

## visio.lign set 2



#### Assortment

43 pieces

#### visio.lign set 2

- 1 x visio.lign tooth cabinet
- 1 x visio.link
- 1 x combo.lign
- 8 x novo.lign P posterior sets
- 10 x novo.lign A anterior sets
- 10 x mixing cannulas, combo.lign
- 8 x drawer insert, novo.lign A
- 4 x drawer insert, novo.lign P

REF VLIGNSET2



## visio.lign sets

## visio.lign set 3



#### Assortment

61 pieces

#### visio.lign set 3

- 1 x visio.lign tooth cabinet
- 1 x visio.link
- 1 x combo.lign
- 1 x Opaquer combo.lign
- 1 x Opaquer catalyst
- 2 x crea.lign
- 8 x novo.lign P posterior sets
- 10 x novo.lign A anterior sets
- 2 x single-hand grip
- 10 x mixing cannulas, combo.lign
- 10 x application cannulas
- 8 x drawer insert, novo.lign A
- 4 x drawer insert, novo.lign P
- 1 x drawer insert, Opaquer combo.lign
- 1 x drawer insert crea.lign/single-hand grip

#### **REF VLIGNSET3**

## novo.lign A and P veneers



#### Assortment

21 pieces

#### novo.lign A and P veneers

- 1 x aluminium folder
- 12 x novo.lign A anterior set
- 8 x novo.lign P posterior set

REF VLIGNPRVN

## neo.lign A and P denture teeth



#### Assortment

25 pieces

#### neo.lign A and P denture teeth

- 1 x aluminium folder
- 12 x neo.lign A anterior set
- 12 x neo.lign P posterior set

#### REF VLIGNPRVTI

#### Assortment

24 pieces

#### neo.lign A and P denture teeth

12 x neo.lign A and P denture to

12 x neo.lign P posterior set

REF NLUPSET2

## Opaquer set



#### Assortment

7 pieces

#### Opaquer set

- 1 x Opaquer combo.lign light
- 1 x Opaquer combo.lign medium
- 1 x Opaquer combo.lign intensive
- 1 x Opaquer combo.lign GUM
- 2 x Opaquer catalyst
- 1 x drawer insert Opaquer combo.lign

**REF OLIGNSET1** 

## crea.lign sets




#### Assortment

46 pieces

#### crea.lign set

- 10 x crea.lign
- 2 x crea.lign Stains
- 1 x crea.lign Modelling Liquid
- 1 x syringe holder
- 12 x single-hand grip
- 20 x application cannulas

**REF CLIGNSET12** 



#### Assortment

13 pieces

#### crea.lign set

- 10 x crea.lign
- 3 x single-hand grip

**REF CLIGNSETN** 

## visio.lign Toolkit



#### Composite processing kit REF VLTOOLKIT

The kit includes: Acrypol polishing paste for veneering composites 170 g Abraso-Starglanz high luster polishing paste

50 ml

The visio.lign Toolkit has been optimized for processing composite and visio.lign veneers and ensures a perfect finish. The combination of the materials of the visio.lign veneering system and the visio.lign Toolkit results in surfaces featuring plaque resistance and shade stability and their quality and resistance is identical to the ones of a ceramic material.

- Burs with relief produce smooth surfaces and reduce the time for polishing
- Prepolishing and high luster polishing paste included in the kit
- Tools arranged at different heights to ensure easy access and grasp
- Removable glass jar to avoid drying out of the high luster polishing paste
- Pictograms and reference numbers printed on the tools for better visual
- Three empty spaces for additional tools



## Light-curing units and times

## Suitable light-curing units

Polymerization times for visio.link, combo.lign and crea.lign

Manufacturer	Product name	Wavelength in nm *	Polymerization time visio.link	Polymerization time combo.lign	Polymerization time crea.lign
bredent	bre.Lux Power Unit	370 - 500	90 s	180 s	6 min
Dentsply / Degudent	Triat, Triat 2000 Eclipse	400 - 500 k.A.	3 min 60 s	6 min 180 s	10 min 6 min
Heraeus Kulzer	Dentacolor XS, Uni XS Heraflash	320 - 520 320 - 520	90 s 90 s	180 s 180 s	6 min 6 min
GC	Laboligth LV-III	380 - 490	2 min	5 min	10 min
Ivoclar Vivadent	Targes Power furnace Lumamat 100	400 - 580 400 - 580	4 min 4 min	180 s 180 s	8 min 6 min
Schütz Dental	Spektra 2000	310 - 500	2 min	180 s	6 min
Shofu Dental	Solitilite EX	400 - 550	90 s	180 s	6 min
Kuraray Dental	CS 110	k.A.	2 min	5 min	8 min
Hager & Werken	Speed Labolight	320 - 550	90 s	180 s	8 min
3M ESPE	Visio BETA (newP1 - P4) Visio BETA (old U0 - U3)**	400 - 500 400 - 500	> 4 min(P2) 7 min (U1, U3)	7 min (P2) 15 min (U0)	15 min (P1) 15 min (U0)

<sup>\*</sup> manufacturer's data \*\* new set of lamps is recommended

## bre.Lux polymerization times

		bre.Lux LED N (	manual lamp)	bre.Lux Power U		
Manufac- turer	Product name	Final polymeri- zation	Fixation/ Prehardening	Intermediate polymerization (layers)	Final polymerization	Prepolymerization function/Reduction
bredent	visio.link	30 s	-	-	90 s	40 s 50 %)
bredent	combo.lign	Х	15 s	120 s	180 s	-
bredent	crea.lign	X	15 s	180 s	360 s	20 s (50 %)
bredent	Opaquer combo.lign	Х	15 s	180 s	180 s	-
bredent	Ropak UV	Х	-	180 s***	360 s	-
bredent	Kompaktopaker	Х	-	180 s***	360 s	-
bredent	Kompaktopaker tooth-colored UV	Х	-	180 s	360 s	-
bredent	compoForm UV	30 s	15 s	-	180 s	-
bredent	Tray material UV*	X	X	90 s	2 x 180 s	40 s (50 %)
bredent	Stumpflack lichthärt.	30 s**	15 s	90 s	180 s	20 s (50 %
bredent	SERACOLL UV	15 s	15 s	-	90 s	-
bredent	Qu-connector	30 s	-	-	90 s	40 s (50 %)
Heraeus	Signum	X	-	180 s	360 s	20 s (50 %)
Heraeus	Palatray XL	Χ	-	90 s	2 x 180 s	40 s (50 %)
Shofu	Solidex	X	-	180 s	360 s	20 s (50 %)
GC	Gradia	Х	15 s	180 s	360 s	20 s (50%)
Wegold	S-Lay	-	-	180 s	360 s	20 s (50 %)
VITA	VITA VM LC Opaque	-	30 s	-	2 x 360 s	-
VITA	VITA VM LC Compos.	-	30 s	180 s	Pontics up to max. 2mm: 360 s	fix up to 1.5 mm, 180s (50%)
Degudent	in:joy	-	-	180 s	360 s	20 s (50%)

<sup>180</sup> s Polymerization time - not intended

<sup>\*</sup> If UV tray material is used, polymerization is carried out from both sides for 180 seconds each. Optionally, prehardening for 90 sec. may be carried out (upper side); during final polymerization, the bottom side is polymerized first \*\* for a single application \*\*\* Apply opaquer in two layers



X Contraindication

#### bre.Lux Power Unit



## Assortment bre.Lux Power Unit

Product package consists of

- 1 polymerization device
- 1 bre.Lux LED N manual lamp with spiral cable
- 1 bre.Lux UpDown 1 flexible hose 1 mains cable

REF 140 0097 0

## Versatile polymerization unit for practices and laboratories.


To date, several devices were necessary for processing facing and dental materials. With bre.Lux, the concept of processing all current materials with one single device becomes reality for the first time.

#### Performance

- Fixation / Hardening / Intermediate polymerization and final thorough hardening directly in the workplace and with one single device
- 370 500 nm covers the required wavelength range, for the manual lamp as well
- Performance delay and reduction minimize the properties and results of dental materials
- Start-up delay and polymerization time can be easily combined
- Large volume for 2 models, optimally and uniformly illuminated

# Keyboard layout of the light polymerization device

- (1) Display
- (2) 20 seconds with 50 % power
- (3) 40 seconds with 50 % power
- (4) Continuous operation with 70 % power
- (5) 90 seconds full power
- (6) 120 seconds full power
- (7) 180 seconds full power
- (8) 360 seconds full power
- (9) Continuous operation at full power



#### Consistency

The membrane keyboard already comes with several parameters. On the left side: Reduced power for the start-up phase with 20 and 40 seconds and for continuous operation at full power as well. The combination (such as start with reduced power) can be combined with the programmed running times of 60 to 360 seconds according to demand. Example: 40 seconds with reduced power combined with a 180-second total running time means that the total running time is lowered from 180 to about 40 seconds, and the running time would still be 140 seconds under full power.

Please take note of the application recommendations for bredent products.

#### Keyboard layout for the manual lamp



(1) Continuous operation, start/stop

(2) 15 seconds

3) 30 seconds



#### Energy

The bre.Lux Power Unit consists of one LED light polymerization device with 21 power LEDs in 3 different capacities, from 370 nm to 500 nm. The LEDs have a useful life of 20,000 hours. The bre.Lux LED N manual lamp (with spiral cable) features a capacity ranging from 370 nm to 500 nm. The flexible hose – with receptacle ring for the manual lamp – serves as third hand and allows two-handed work.



## Range of designs novo.lign veneers


## novo.lign A, upper anterior Facet in the cervical and central area is 1 mm thick

C43 A44 145 **S46** 46.0 147 47.4 D48 48.4 M48 D49


308

48.6

## novo.lign A, upper/lower anterior Facet in the cervical and central area is 1 mm thick



novo.lign A, novo.lign P and combo.lign are available in shades of the classic A-D shade system as well as in the bleached shade BL3.




Design -	- combination	table
anterior		posterior
Upper	Lower	Upper/Lower
C43	T35	G3/W3
A44	T35	G3/W3
145	T35	G3/W3
S46	T35/D38	G3/W3
147	D38	G3/W4
D48	D38	G3/G4/W4
M48	D38	G3/G4/W4
D49	D38/D41	G3/G4/W4
B51	D41	G4/W5
K53	D41	G4/W5


## Range of designs novo.lign veneers


## novo.lign P, upper / lower posterior, G-design Facet in the cervical and central area is 1 mm thick

Multi-functional veneers

G3


















G4


















310


## novo.lign P, upper / lower posterior, W-design Facet in the cervical and central area is 1 mm thick


Veneers for crowns and bridges

W3























W4





















W5



















## neo.lign A, upper anterior



312

49.0


B51







T35



9.6




D38







D41








neo.lign is available in the classic VITA A-D shade system.

ISO 22112:2006

Supply based on quadrant



Design - combination table

anterior	posterior	
Upper	Lower	Upper/Lower
C43	T35	G2
A44	T35	G2/G3
145	T35	G2/G3
S46	T35/D38	G3
147	D38	G3
D48	D38	G3/G4
M48	D38	G3/G4
D49	D38/D41	G3/G4
B51	D41	G4

## neo.lign P, upper/lower posterior, G-design


G2

















G3

















G4

















314

## neo.lign P, upper/lower Posterior, L-design

neo.lign P for lingualized setups (dimension according to pilot series, deviations possible)

L2













L3













L4











novo.lign A		BL3	A1	A2	А3	A3,5				A-D sl B3		systen C1	1) C2	С3	C4	D2	D3	D4
Upper anterior veneers	REF	BL3	A10	A20	A30	A35	A40	B10	B20	B30	B40	C10	C20	C30	C40	D20	D30	D40
C43 Set of 6 (13, 12, 11, 21, 22, 23)	V0C43 S																	
C43 Set of 4 (12, 11, 21, 22)	V0C43 4																	
C43 Set of 2 (13, 23)	V0C43 3																	
A44 Set of 6 (13, 12, 11, 21, 22, 23)	V0A44 S																	
A44 Set of 4 (12, 11, 21, 22)	V0A44 4																	
A44 Set of 2 (13, 23)	V0A44 3																	
l45 Set of 6 (13, 12, 11, 21, 22, 23)	V0I45 S																	
l45 Set of 4 (12, 11, 21, 22)	V0I45 4																	
l45 Set of 2 (13, 23)	V0I45 3																	
S46 Set of 6 (13, 12, 11, 21, 22, 23)	V0S46 S																	
S46 Set of 4 (12, 11, 21, 22)	V0S46 4																	
S46 Set of 2 (13, 23)	V0S46 3																	
l47 Set of 6 (13, 12, 11, 21, 22, 23)	V0I47 S																	
l47 Set of 4 (12, 11, 21, 22)	V0I47 4																	
l47 Set of 2 (13, 23)	V0I47 3																	
D48 Set of 6 (13, 12, 11, 21, 22, 23)	VOD48 S																	
D48 Set of 4 (12, 11, 21, 22)	VOD48 4																	
D48 Set of 2 (13, 23)	VOD48 3																	
M48 Set of 6 (13, 12, 11, 21, 22, 23)	V0M48 S																	
M48 Set of 4 (12, 11, 21, 22)	V0M48 4																	
M48 Set of 2 (13, 23)	V0M48 3																	
D49 Set of 6 (13, 12, 11, 21, 22, 23)	VOD49 S																	
D49 Set of 4 (12, 11, 21, 22)	VOD49 4																	
D49 Set of 2 (13, 23)	VOD49 3																	
B51 Set of 6 (13, 12, 11, 21, 22, 23)	VOB51 S																	
B51 Set of 4 (12, 11, 21, 22)	VOB51 4																	
B51 Set of 2 (13, 23)	VOB51 3																	
K53 Set of 6 (13, 12, 11, 21, 22, 23)	V0K53 S																	
K53 Set of 4 (12, 11, 21, 22)	V0K53 4																	
K53 Set of 2 (13, 23)	V0K53 3																	
Lower anterior veneers																		
T35 Set of 6 (43, 42, 41, 31, 32, 33)	VUT35 S																	
T35 Set of 4 (42, 41, 31, 32)	VUT35 4																	
T35 Set of 2 (43, 33)	VUT35 3																	
D38 Set of 6 (43, 42, 41, 31, 32, 33)	VUD38 S																	
D38 Set of 4 (42, 41, 31, 32)	VUD38 4																	
D38 Set of 2 (43, 33)	VUD38 3																	
D41 Set of 6 (43, 42, 41, 31, 32, 33)	VUD41 S																	
D41 Set of 4 (42, 41, 31, 32)	VUD41 4																	
D41 Set of 2 (43, 33)	VUD41 3																	

ase enter the order quant

Sender (Stamp):

Customer No.

Date, Signature

Shade

REF

VLIGNSET2

							Shade	es (cla	assic A	A-D sl	hade s	systen	1)					
novo.lign P, multi-functional		BL3	A1	A2	А3	A3,5	<b>A4</b>	В1	B2	В3	B4	C1	C2	C3	C4	D2	D3	D4
Veneers, posterior G 3	REF	BL3	A10	A20	A30	A35	A40	B10	B20	B30	B40	C10	C20	C30	C40	D20	D30	D40
1G3 1. quadrant (14, 15, 16, 17)	V01G3																	
2G3 2. quadrant (24, 25, 26, 27)	V02G3																	
3G3 3. quadrant (34, 35, 36, 37)	VU3G3																	
4G3 4. quadrant (44, 45, 46, 47)	VU4G3																	
Veneers, posterior G4																		
1G4 1. quadrant (14, 15, 16, 17)	V01G4																	
2G4 2. quadrant (24, 25, 26, 27)	V02G4																	
3G4 3. quadrant (34, 35, 36, 37)	VU3G4																	
4G4 4. quadrant (44, 45, 46, 47)	VU4G4																	
Veneers posterior W3																		
1W3 1. quadrant (14, 15, 16, 17)	V01W3																	
2W3 2. quadrant (24, 25, 26, 27)	V02W3																	
3W3 3. quadrant (34, 35, 36, 37)	VU3W3																	
4W3 4. quadrant (44, 45, 46, 47)	VU4W3																	
Veneers posterior W4																		
1W4 1. quadrant (14, 15, 16, 17)	V01W4																	
2W4 2. quadrant (24, 25, 26, 27)	V02W4																	
3W4 3. quadrant (34, 35, 36, 37)	VU3W4																	
4W4 4. quadrant (44, 45, 46, 47)	VU4W4																	
Veneers posterior W5																		
1W5 1. quadrant (14, 15, 16, 17)	V01W5																	
2W5 2. quadrant (24, 25, 26, 27)	V02W5																	
3W5 3. quadrant (34, 35, 36, 37)	VU3W5																	
4W5 4. quadrant (44, 45, 46, 47)	VU4W5																	

visio.lign veneer system Set 3	REF	Shade
such as VLIGNSET2, 2 x crea.lign, 2x Opaquer combo.lign, 12 (instead 10) x novo.lign A, 27 pieces, shade at customer's option	VLIGNSET3	
such as VLIGNSET2, 2 x crea.lign, 2x Opaquer combo.lign, 12 (instead 10) x novo.lign A, 27 pieces, shade at customer's option	VLIGNSET3	

10 x novo.lign A, 8 x novo.lign P, 1 x combo.lign, 1 x mixing cannulas, 1 x visio.link, 21 pieces, shade at customer's option

Sender (Stamp):

Customer No.

Date, Signature

visio.lign veneer system Set 2

12 pcs. 320 0045 1

110 0147 0

combo.lign			BL3	<b>A</b> 1	A2	АЗ	A3,5	A4	B1	B2	В3	B4	C1	C2	С3	C4	D2	D3	D4
REF			BL3	A10	A20	АЗ	A35	A40	B10	B20	B30	B40	C10	C20	C30	C40	D20	D30	D40
Fixation composite, 8 g	C02x4 _																		
Mixing cannulas combo.lign, 10 pieces	COMKG2																		
combo.lign tooth-colored	Cont.	REF			Quanti	ty	visio.s	il silic	one f	or key	'S			Con	t.	REF		Qu	antity
Composite, dual-curing	8 g	C02	X4GU	М			visio.s	il trans	sparer	nt				50 ı	nl	540 0	120 0		
	•						visio.s	il ILT t	ranspa	arent				50 ı	ml	540 0	140 0		
combo.lign Opaquer	Cont.	REF			Quanti	tv	Mixing	g cann	ulas (	G4) vis	sio.sil/v	visio.si	LILT	12	ocs.	320 0	045 7		
light for A1 / A2 / A2 / B2	4 a		V40D		Qualiti	СУ	visio.s	il fix, k	nigh-p	recisio	n			50 ı	ml	540 0130 0			

Mixing cannulas (G2) visio.sil fix

Thermo-Pen

combo.lign Opaquer	Cont.	REF	Quantity
light for A1 / A2 / A3 / B2	4 g	CO1X40PL	
medium for A3,5 / B3 / D3	4 g	CO1X40PM	
intensive for A4 / C3	4 g	CO1X40PI	
GUM	4 g	CO1X40PG	
Catalyst for all shades	4 g	CO1X4KAT	
all 4 shades + 2 x cat	6 x 4 g	OLIGNSET1	

beauty setup wax	Cont.	REF	Quantity
tooth-colored, light	25 g	430 0030 0	
tooth-colored, dark	25 g	430 0031 0	

visio.link	Cont.	REF	Quantity
PMMA & Composite Primer	10 ml	VLPMMA10	

visio.lign Toolkit	Cont.	REF	Quantity
Composite processing kit	10 pcs.	VLTOOLKIT	

MKZ Primer	Cont.	REF	Quantity
Metal and zirconium oxide primer	4 ml	MKZ02004	

Accessories	Cont.	REF	Quantity
Dispenser 5 ml 1:1	1 pcs.	320 0044 1	
Grip vor one-hand use crea.lign syringe	1 pcs.	320 0044 2	
crea.lign Syringe holder (Tray)	1 pcs.	320 0044 3	
Applikation cannulas	10 pcs.	320 0094 0	
crealign Modelling Liquid	10 ml	CLFMOD10	

MKZ EM-Aktivator	Cont.	REF	Quantity
Activator for precious metal frameworks	4 ml	MKZEM004	

K-Primer	Cont.	REF	Quantity
Veneering ceramic primer	3 ml	PPK25003	

Sender (Stamp):

Customer No.

Date, Signature

crea.lign	Cont.	REF	Quantity
crea.lign Dentin BL3	5 g	CLFNDBL3	
crea.lign Dentin A1	5 g	CLFNDA10	
crea.lign Dentin A2	5 g	CLFNDA20	
crea.lign Dentin A3	5 g	CLFNDA30	
crea.lign Dentin A3,5	5 g	CLFNDA35	
crea.lign Dentin A4	5 g	CLFNDA40	
crea.lign Dentin B1	5 g	CLFNDB10	
crea.lign Dentin B2	5 g	CLFNDB20	
crea.lign Dentin B3	5 g	CLFNDB30	
crea.lign Dentin B4	5 g	CLFNDB40	
crea.lign Dentin C1	5 g	CLFNDC10	
crea.lign Dentin C2	5 g	CLFNDC20	
crea.lign Dentin C3	5 g	CLFNDC30	
crea.lign Dentin C4	5 g	CLFNDC40	
crea.lign Dentin D2	5 g	CLFNDC10	
crea.lign Dentin D3	5 g	CLFNDC20	
crea.lign Dentin D4	5 g	CLFNDC30	
crea.lign Incisal E1	5 g	CLFN00E1	
crea.lign Incisal E2	5 g	CLFN00E2	
crea.lign Incisal E3	5 g	CLFN00E3	
crea.lign Incisal E4	5 g	CLFN00E4	
crea.lign Incisal opal	5 g	CLFN00I1	
crea.lign Incisal blue	5 g	CLFN00I2	
crea.lign Incisal rose	5 g	CLFN00I3	
crea.lign Incisal universal	5 g	CLFN00I4	
crea.lign GUM light	5 g	CLFN00G1	
crea.lign GUM rosa	5 g	CLFN00G2	
crea.lign GUM pink	5 g	CLFN00G3	
crea.lign Modifier beige	5 g	CLFN00M1	
crea.lign Modifier oliv	5 g	CLFN00M2	
crea.lign Modifier caramel	5 g	CLFN00M3	
crea.lign Modifier lila	5 g	CLFN00M4	
crea.lign Stain orange	1,4 g	CLFN00S1	
crea.lign Stain brown	1,4 g	CLFN00S2	

Accessories	Cont.	REF	Quantity
crea.lign Modelling Liquid	10 ml	CLFMOD10	

crea.lign Opake		Cont.	REF	Quantity
crea.lign Opake 1	A1 / B2	4 g	CLFHOP1	
crea.lign Opake 2	A2	4 g	CLFHOP2	
crea.lign Opake 3	A3 / D3	4 g	CLFHOP3	
crea.lign Opake 4	BL3 / B1 / C1	4 g	CLFH0P4	
crea.lign Opake 5	C2 / C3 / D2 / D4	4 g	CLFHOP5	
crea.lign Opake 6	B3 / B4	4 g	CLFHOP6	
crea.lign Opake 7	A3,5	4 g	CLFH0P7	
crea.lign Opake 8	A4 / C4	4 g	CLFHOP8	
crea.lign Opake GUM		4 g	CLFHOGUM	

crea.l	crea.lign Starter kit														Cont.	R	REF		Quantity		
syring	syringes a 5 g each (at customer's option)												10 pcs	i. C	LIGNSETI	Ν					
DBL3	DA1	DA2	DA3	DA3,5	DA4	DB2	DB3	DB4	DC2	DC3	DD2	DD3	DD4	E2	E3	opal	blue	univ	G1	G2	G3

crea.li	crea.lign Starter kit with syringe holder														Cont.	RI	F	(	Quantity		
syringe	syringes of 5 g each, 2 x Stain, Modelling Liquid, Syringe holder, 12 x single-hand grip 12 pcs.													CLIGNSET12							
DBL3	DA1	DA2	DA3	DA3,5	DA4	DB2	DB3	DB4	E2	E3	opal	blue	univ	rose	M1	M2	M3	M4	G1	G2	G3

syringes a 5 g, 2 x Stain, Modelling Liquid, Syringe holder, 12 x single-hand grip

Sender (	(Stamp)	):

Customer No.

Date, Signature



		_						•			nade s	,	•					
neo.lign A		BL3	A1	A2	А3	A3,5	A4	В1	B2	В3	B4	C1	C2	C3	C4	D2	D3	D4
Anterior teeth upper anterior	REF	BL3	A10	A20	A30	A35	A40	B10	B20	B30	B40	C10	C20	C30	C40	D20	D30	D40
C43 Set of 6 (13, 12, 11, 21, 22, 23)	T0C43 S																	
A44 Set of 6 (13, 12, 11, 21, 22, 23)	T0A44 S																	
l45 Set of 6 (13, 12, 11, 21, 22, 23)	T0I45 S																	
S46 Set of 6 (13, 12, 11, 21, 22, 23)	T0S46 S																	
l47 Set of 6 (13, 12, 11, 21, 22, 23)	T0I47 S																	
D48 Set of 6 (13, 12, 11, 21, 22, 23)	TOD48 S																	
M48 Set of 6 (13, 12, 11, 21, 22, 23)	T0M48 S																	
D49 Set of 6 (13, 12, 11, 21, 22, 23)	TOD49 S																	
B51 Set of 6 (13, 12, 11, 21, 22, 23)	TOB51 S																	
Anterior teeth lower anterior																		
T35 Set of 6 (43, 42, 41, 31, 32, 33)	TUT35 S																	
D38 Set of 6 (43, 42, 41, 31, 32, 33)	TUD38 S																	
D41 Set of 6 (43, 42, 41, 31, 32, 33)	TUD41 S																	

							Shade	es (cla	assic A	A-D sl	hade s	ysten	1)					
neo.lign P multi-functional		BL3	<b>A</b> 1	A2	А3	A3,5	A4	В1	B2	В3	B4	C1	C2	С3	C4	D2	D3	D4
Posterior teeth posterior G 2	REF	BL3	A10	A20	A30	A35	A40	B10	B20	B30	B40	C10	C20	C30	C40	D20	D30	D40
1G2 1. quadrant (14, 15, 16, 17)	T01G2																	
2G2 2. quadrant (24, 25, 26, 27)	T02G2																	
3G2 3. quadrant (34, 35, 36, 37)	TU3G2																	
4G2 4. quadrant (44, 45, 46, 47)	TU4G2																	
Posterior teeth posterior G 3																		
1G3 1. quadrant (14, 15, 16, 17)	T01G3																	
2G3 2. quadrant (24, 25, 26, 27)	T02G3																	
3G3 3. quadrant (34, 35, 36, 37)	TU3G3																	
4G3 4. quadrant (44, 45, 46, 47)	TU4G3																	
Posterior teeth posterior G 4																		
1G4 1. quadrant (14, 15, 16, 17)	T01G4																	
2G4 2. quadrant (24, 25, 26, 27)	T02G4																	
3G4 3. quadrant (34, 35, 36, 37)	TU3G4																	
4G4 4. quadrant (44, 45, 46, 47)	TU4G4																	

								Shade	es (cla	assic A	A-D sl	nade s	systen	1)					
neo.	lign P lingualized		BL3	Α1	A2	А3	A3,5	A4	B1	B2	В3	B4	C1	C2	C3	C4	D2	D3	D4
Post	erior teeth posterior L 2	REF	BL3	A10	A20	A30	A35	A40	B10	B20	B30	B40	C10	C20	C30	C40	D20	D30	D40
L2	UJ (14, 15, 16, 17/24, 25, 26, 27)	TOWL2																	
L2	□ (34, 35, 36, 37/44, 45, 46, 47)	TUWL2																	
Post	erior teeth posterior L 3																		
L3	UJ (14, 15, 16, 17/24, 25, 26, 27)	T0WL3																	
L3	□ (34, 35, 36, 37/44, 45, 46, 47)	TUWL3																	
Post	erior teeth posterior L 4																		
L4	UJ (14, 15, 16, 17/24, 25, 26, 27)	TOWL4																	
L4	□ (34, 35, 36, 37/44, 45, 46, 47)	TUWL4																	

Design chart neo.lign	REF	Shade
24 pcs. (12 x neo.lign A, 12 x neo.lign P), shade at customer' s option	VLIGNPRVTI	

neo.lign Upgrade Kit	REF	Shade	
24 pcs. (12 x neo.lign A, 12 x neo.lign P), shade at customer' s option	NULPSET2		

Sender (	(stamp)	):
----------	---------	----

Customer No.

Date, Signature


## HIGH.lign

#### HIGH.lign - the resin for the future!

The new HIGH.lign series of resins includes uni.lign as a high-quality denture resin and top.lign as a crown and bridge material. The resins of the uni.lign group offer dental technicians a harmonized range of cold-curing denture resins for top-quality dentures, which can be fabricated using teeth from the visio.lign veneering system. The innovative uni.lign speed repair resin enables repairing restorations within a short time. top.lign breformance is used for the

fabrication of high-quality long-term temporaries. top.lign professional is a material to prepare permanent restorations. All resins offer outstanding mechanical properties and performance for the respective indication.





## uni.lign

#### Top quality for sophisticated dentures

Thanks to the innovative formula and a new manufacturing process, the requirements of DIN EN ISO 20795-I are exceeded by far. As a result, denture resins with exceptional mechanical properties are produced which comply with the latest standards. Exclusively high-quality raw materials, which fulfill the stringent requirements of bredent, are processed to achieve utmost and consistent quality and shade stability. Consistent quality is ensured by comprehensive in-house controls for each production lot.

Simple processing in the laboratory and harmonization of the individual

resins ensure ease of use and a high level of comfort for dental technicians. The optimized modulus of elasticity, high flexural strength and resistance to discoloration lead to durable restorations for patients. The low residual monomer content and excellent polishing properties (plaque resistance) guarantee outstanding compatibility.

The use of crea.lign enables esthetic individualization in particular when fabricating implant-supported restorations and underlines the laboratory's competence.



Lab. Od. Lazetera Antonio - Savona - Italy Dott. Vescia Luca -Villa Dossola - Italy



Lab. Od. Lazetera Antonio - Savona - Italy

Three different shades provide more flexibility. Three levels from transparent to opaque and three veined shades enable perfect reproduction of the gingiva for partial dentures.



Information on ordering on page 325

### uni.lign

### Cold-curing denture base material

uni.lign consists of a powder component in the respective shade/color which is either mixed with uni.lign liquid or uni.lign liquid cast in accordance with application-specific characteristics. As a result, the processing time span can

be varied and additional options for the fabrication process of dentures are obtained.




Three different package sizes - from the convenient small package (70 g) to the large refill package (1000 g).



Easily distinguishable and clear labelling of the resins. Systematic and clear designation of the shade prevents mixing up the different colors.

### The various liquids

Compared to uni.lign liquid, the modelling phase can be extended when using uni.lign liquid and a wide range of processing options is obtained. Large dentures or the Opti-Cast flask can be cast without any bubbles and working stress is eliminated.



### Indication - uni.lign with uni.lign liquid



Completion of CoCr dentures.



Relinings and shaping of functional margins.



Any type of repair work, such as cracks, fractures, etc.

### Indication - uni.lign with uni.lign liquid cast



Fabrication of full dentures in the resin casting technique.



Completion of CoCr dentures.




Relinings and shaping of functional margins.




Any type of repair work, such as cracks, fractures, etc.

### uni.lign




#### Cold-curing denture base material

The swelling behavior clearly indicates that the use of uni.lign liquid cast mixing liquid slowly dissolves the resin beads and hence the processing time is extended. The use of uni.lign liquid, however, accelerates swelling of the resin beads so that a shorter processing time for repairs, relinings or small saddles is obtained.



An optimum modulus of the elasticity results in a more flexible denture, minimizes the tendency to fracture and increases the reliability for the patient.



Surface density optimized for the respective indication and very low water absorption contribute to retarding the ageing process of the uni.lign resins and ensuring a high quality standard and extended durability. Moreover, excellent shade stability guarantees that the esthetic appearance is preserved throughout the life of the denture.



# uni.lign speed

### Fast-curing repair resin

Performing repairs, extensions and fabricating small saddles with uni.lign speed – these types of work can be done since the new resin has a short swelling time which does not affect the mechanical values.



Since it has been matched with the uni. lign colors, uni.lign speed enables repairs or extensions with no difference in color.



The short pouring time (1 minute) and the extremely short polymerization time of only 5 minutes allow fast further processing so that considerable time can be saved.



The modelling phase of 2 minutes enables accurate application of the resin and hence reduces the time required for finishing the cured resin.

# Overview of the uni.lign resins





	uni.lign		uni.lign speed		
	70 g	500 g	1000 g	70 g	500 g
PC 10	uniPC101	uniPC102	uniPC103	unispc11	unispc15
PC 20	uniPC201	uniPC202	uniPC203	unispc21	unispc25
PC 30	uniPC301	uniPC302	uniPC303	unispc31	unispc35
PF 10	uniPF101	uniPF102	uniPF103	_	_
PF 20	uniPF201	uniPF202	uniPF203	_	_
PF 30	uniPF301	uniPF302	uniPF303	_	_
TC 10	uniTC101	uniTC102	uniTC103	_	-

uni.lign liquid		uni.lign liquid cast			uni.lign speed liquid		
100 ml	500 ml	1000 ml	100 ml	500 ml	1000 ml	100 ml	500 ml
unil0100	unil0500	unil1000	unilc100	unilc500	unilc000	unisplq1	unisplq5

# top.lign professional

### top.lign professional

 $\label{linear_power_liquid} Discoloration-resistant powder-liquid system for fast and simple fabrication of long-term temporaries or permanent restorations.$ 

Particularly suitable for immediate large-span bridges based on the fast

& fixed system by bredent. The high density enables quick and excellent polishing. The high gloss avoids adhesion of plaque and extends the period of wearing. As a result, the reliability is increased during the healing phase of the implants.

### Application examples



Implant-supported, screw-retained SKY fast & fixed bridge.



Pouring of the implant-supported bar restoration with top. lign professional.



Clasp-retained dentures.



Completed implant-supported bar restoration.



top.lign professional liquid 100 ml REF tlp liq0 1

#### Accessories



Isoplast ip 750 ml REF 540 0101 9




Round brush Rodeo 15 pieces, Ø 18 mm REF 350 0096 0



# Dentaclean impression and denture disinfectant



Disinfecting with Dentaclean impression and denture disinfectant avoids the transmission of viruses, bacteria and fungi – from the patient to the laboratory. The concentrate is mixed to obtain 10 liters of ready-to use solution which is highly effective and has a surprisingly mild odor.



Pathogens can be transmitted to the laboratory with impressions.

Dentaclean impression and denture disinfectant 1000 ml concentrate to obtain 10 liters ready-to-use solution incl. 25 shipping bags REF 520 0100 6

Tested and approved by the Institute for clinical hygiene and infection control, Giessen.



After the use of Dentaclean impression disinfectant, active viruses, bacteria and fungi can no longer be detected.

### Shipping bags



The shipping bags have already been labeled "disinfected". Additionally, a separate bag holds the refte to protect them against moisture.

Shipping bags 200 pieces REF 520 0100 2

### Dentaclean denture cleaning agent

Konzentrat zur mü-

helosen Entfernung

von Plaque, Zahn-

stein und Belägen

auf Prothesen.



Dentaclean denture cleaning agent 1000 ml concentrate to obtain 11 liters ready-to-use solution REF 520 0099 2



Contaminated dentures are unpleasant and require a lot of time for cleaning.

Up to now the removal of

tartar has been difficult

and could often only be

achieved through grind-

ing. This is unpleasant

and takes a lot of time.



in Dentaclean denture cleaning agent remove difficult coatings from dentures safely and quickly within only 15 minutes.

Now high-quality con-

centrate components



Dentures can be quickly and easily cleaned with Dentaclean denture cleaning agent.



### Disinfecting and cleaning

# Dentaclean ultrasonic cleaning agent



Concentrate for removal of polishing paste residues. Mild odor, powerful cleaning capacity.



Cleaning of polishing contaminations takes a lot of time. Therefore aggressive agents that are injurious to health are frequently used.



Matched surfactants and emulsifiers remove contaminations carefully and quickly thus saving time for the technician.

#### Dentaclean ultrasonic cleaning agent 1000 ml concentrate to obtain 11 liters ready-to-use solution REF 520 0099 7

### Dentaclean plaster removing agent / Dentaclean plaster removing agent Speed



# Ready-to-use solution to remove plaster residues

The Dentaclean plaster removing agent is available in two types: normal and Speed. The ready-to-use solution removes plaster residues from all surfaces. If no time is to be wasted, Dentaclean Speed should

Dentaclean plaster removing agent 1000 ml REF 520 0011 9 2500 ml

REF 520 0099 3

Dentaclean plaster removing agent Speed 1000 ml REF 520 0101 0 2500 ml REF 520 0099 4



Hard plaster particles are carefully removed from the mixing bowl without any damage.



Gentle and fast removal of plaster protects the acrylic surface and the color.

### Dentaclean mixing fluid for pumice powder



Dentaclean mixing fluid for pumice powder 5000 ml

Dentaclean mixing fluid for pumice powder 1000 ml REF 520 0099 8 REF 520 0099 9

#### Application:

Simply mix the pumice powder with Dentaclean mixing fluid for pumice powder, do not add water. This way the pumice powder stays moist for two to three weeks.

Protects against germs.

#### Dentaclean mixing fluid for pumice powder

- · Remains moist and free of germs for two to three weeks without having to be remixed.
- Contains skin-care additives to protect employees' hands.
- · Contains natural odours which still smell fresh after several weeks.
- Mixed polish adheres to the brush and restorstion better so that the pumice splatters less. This saves time when polishing as the pumice slurry does not have to be applied repeatedly



In wet punmise powder disease microbs are present. Disinfection action occures within one hour.



The nurturing effect for the skin is acheived by skin nurturing addative.



# Ergonom wax knife



### Modelling knife for dental prosthetics.

Various instruments all in one – hence instruments do not need to be changed any longer so that faster and more efficient processing of the wax model is possible.

Ergonomic design of handle - suitable for right- and left-hand users.

Ergonom Wax knife REF 310 0001 3



Special, ground edge of the knife tip for simple and fast modelling of age-specific papillae shapes.



The spoon which features a ground edge is perfectly suitable for modelling the alveolar area.



The deep spoon perfectly allows to apply large wax quantities within a very short time.



Well-aimed, fast application of wax reduces the time for remodelling in the interdental area.



Transitions towards the functional margin and the functional margin itself can be prepared swiftly and neatly thanks to the curved design of the spoon element.






With the Ergonom wax knife, wax models can be easily and quickly shaped so that a natural appearance is obtained.

### **Units / Instruments**

### Repositioning tweezers








The pointed pincer tips permit secure holding of teeth and other small parts.



special pincer tips for secure holding
 pincer tips of hardened material for a long working life
 no slipping of small parts – no irritating searching
 fine tips for narrow areas

The special and well designed denticulation of the pincers provides optimal hold security.

Repositioning tweezers 1 piece REF 310 0011 5



Secure holding of teeth is not possible with normal tweezers. Time-consuming searching is eliminated.



Small parts such as screws or attachments are gripped easily and securely. A useful instrument particularly for implantology.



On completion there is always a problem – the repositioning of the teeth! The special fine pincers at the tweezers tips permit secure gripping of the teeth.

### Thermo-syringe



Thermo-syringe REF 110 0121 1



After heating, the adhesive acrylic wax is directly applied onto the glueing point using the Thermo-syringe. Firm bonding is ensured.



The adhesive acrylic wax can be applied onto any type of material. Afterwards it can be removed from the objects without leaving any residues.

### Accessories:



**Adhesive acrylic wax** Pack cont. 250 g

Bucket cont. 1000 g

REF 510 0070 1

REF 510 0070 0

Fixing and glueing that can be dissolved without any residues for any type of model situation.

The adhesive acrylic wax can be moulded by heating and easily applied to the models.



# Posi-boy



# The perfect "third hand" to hold any model in the desired position.

Posi-boy simplifies processing of cold-curing acrylics. The solid metal base ensures firm stand and the correct position in the pressure pot. No tilting, no leaking of acrylics, no change of pre-shaped saddles.



The corrosion-resistant V2 A material guarantees a long service life for the Posi-boy and keeps the acrylics "in shape" in any pressure pot. Thanks to its robust design and the individual adjustment options, firm hold of the model is ensured.

#### Posi-boy REF 360 0101 0

# Articulation paper holder



Repeated taking up and placing down the handpiece and articulation paper are no longer required. Grinding in a flick of the wrist!



Articulation paper holder size 1 1 piece REF 360 0121 7



Articulation paper holder size 2 1 piece REF 360 0122 0

### Acessories:



Set-up grinding tool 1 piece REF 340 0101 0

# **Insulating agents**

# Wax insulating agent



Microfeine Isolierflüssigkeit für alle Wachsmodellationen. Isoliert Gips, Kunststoff, Metall und sogar Wachs gegen Wachs.

Wax insulating liquid wis with brush pen pk 20 750 ml REF 540 0070 4


#### Accessories:




Plastic spray bottle sp 125 ml REF 540 0075 0

20 ml

REF 540 0072 0



The brush pen allows to apply defined quantities of the wax insulating liquid to the desired areas of the model.



# Plaster insulating agent



For reliable insulation of plaster against plaster. Alginate-based plaster insulating liquid which ensures gap-free fit. For utmost precision and separating of sawcut models without any damage.

Plaster insulating agent 750 ml REF 540 0013 5

### Accessories:



Plastic spray bottle sp 125 ml REF 540 0075 0



Brush pen pk 20 20 ml REF 540 0072 0



The plaster insulating liquid soaks into the plaster and seals the surface without layering. The brush pen allows quick application.



The plaster insulating liquid allows separating the two flask halves without any damage.



# Isoplast ip



Plaster-acrylic insulating liquid for all cold- and hot-curing acrylics.



Isoplast ip 750 ml

REF 540 0101 9

Accessories:



The desired quantity of Isoplast can be applied to ensure economic use.



The plaster-acrylic insulating liquid seals the surface. This way precise impressions are ensured.



Isoplast allows to obtain extremely smooth, shining acrylic surfaces. The finishing time is reduced.



Brush pen pk 125 125 ml REF 390 0033 0



### Bite blocks



Prefabricated wax bite blocks – available in the shape of jaws or rods featuring different degrees of hardness.



The basal profile of the bite blocks allows time-saving adaptation on the base plate.



Bite blocks bw rods medium, red 104 pieces 14 x 8 x140 mm REF 430 0023 0



medium, red, UJ/⊔
74 pieces
REF 430 0022 0
medium, red, UJ
74 pieces
REF 430 0020 0
medium, red, ⊔
74 pieces
REF 430 0021 0



hard, yellow, UJ/LJ 74 pieces REF 430 0017 0 hard, yellow, UJ 74 pieces REF 430 0015 0 hard, yellow, LJ 74 pieces REF 430 0016 0



soft, pink 104 pieces 14 x 8 x140 mm REF 430 0028 0



soft, pink, UJ/LJ 74 pieces REF 430 0027 0 soft, pink, UJ 74 pieces REF 430 0025 0 soft, pink, LJ 74 pieces REF 430 0026 0



super-hard, white, UJ/LJ 74 pieces REF 430 0012 0 super-hard, white, UJ 74 pieces REF 430 0010 0 super-hard, white, LJ 74 pieces REF 430 0011 0



hard, yellow 104 pieces 14 x 8 x140 mm REF 430 0018 0



super-hard, white 104 pieces 14 x 8 x140 mm REF 430 0013 0

The height and the width of prefabricated bite blocks are suitable for the use on partial dentures.



Prepare situation model in the usual way.



The basal profile of the bite block simplifies adapting on the base plate.



The consistency of the bite blocks allows simple reduction of the height and width using the wax knife.



The prefabricated wax bite block can be easily integrated.



Since additional application of wax to the buccal and lingual area is no longer required, a considerable amount of time and material can be saved.



The high stability and functional processing of the bite blocks ensures precise bite-taking.



# Set-up wax asw



For setting up and changing the position of acrylic teeth without heating.



Set-up wax asw 4 pink 220 g REF 430 0157 4



Set-up wax asw 5 pink 220 g REF 430 0152 0



Set-up wax asw 3 pink 220 g REF 430 0151 0



Three different sizes of the pink set-up wax allow the individual use.



Thanks to its consistency the set-up wax can be perfectly processed without being heated.



The set-up wax allows quick adapting on the base plate.



No additional wax is required for flushing of the set-up wax.



Due to the adhesive capacity of the set-up wax, acrylic teeth are fixed prior to waxing on.



Even after waxing on, acrylic teeth can be brought into any individual position.

# Modelling wax pink Standard mdwst



Modelling waxes in sheets are used for a large number of applications in denture work.

Modelling wax pink Standard mdwst sheets.

Two thicknesses and three different qualities provide the technician with individual processing options.

Sheet thickness 1.25 mm quantity 1000 g 75 x 150 x 1.25 mm soft, pink medium, pink hard, pink

REF 430 0164 3 REF 430 0164 2 REF 430 0164 1

Sheet thickness 1.50 mm quantity 1000 g 75 x 150 x 1.5 mm soft, pink medium, pink hard, pink

REF 430 0164 6 REF 430 0164 5 REF 430 0164 4



Due to the particular stability of the pink modelling wax sheets, sufficient stability for the base plates is provided.



By rolling up the wax sheets and waxing them to the base plate, acrylic teeth can be set up immediately.



Bite blocks can be easily produced by rolling up and kneading this modelling wax.



The structure of the pink modelling wax sheets allows easy blocking-out for individual trays.




Even during extended try-in, the original stability of this modelling wax is maintained.

# Wax palatal patterns gf



More quality, function and esthetics within a short time.

The recesses for the acrylic teeth simplify adapting of the pre-shaped wax palatal patterns to the situation



Assortment 0.5 mm A, B Wax palatal patterns gf 1.5 mm A, B 120 pieces, REF 430 02 30 pieces each

0.5 mm A, B 1.5 mm A, B REF 430 0218 0

The use of wax palatal patterns for the wax setup simplifies modelling and saves time.



Cut out the pink modelling wax and replace it by wax palatal patterns.




The pre-shaped wax palatal patterns can be easily adapted.



The transition of the wax palatal patterns to the approximal area can be perfectly designed with the fine modelling tip of the wax knife.



The natural function and esthetics of the palate is restored.



size A

0.5 mm 100 pieces REF 430 0214 A

1.5 mm 60 pieces REF 430 0211 A



size E

\_\_\_\_ 0.5 mm 110 pieces REF 430 0215 B

1.5 mm 70 pieces **REF 430 0212 B** 

### Tray material UV



Highly stable light-curing resin for trays and base plates.

The flexibility of the material allows easy placement onto the model without tearing. The required shape can be cut with an instrument. The pink color provides the perfect basis for the set-up.



Tray material UV 50 pieces UJ REF 540 0011 0



25 x Tray material UM - UJ 25 x Tray material UV - LJ

Assortment

REF 540 0011 2

50 pieces

Tray material UV 50 pieces ⊔ REF 540 0011 1



Tray material UV band 2.5 mm x 90 mm 1350 g REF 540 0016 6

Tray material UV block 1000 g REF 540 0011 3

### Accessories:



Polylux pl 20 Polylux polymerization unit with material container (see page 251)

REF 140 0088 0



The high flexibility of the material simplifies the placement onto the model. The material will not be damaged.



The tray material can be precisely cut with any instrument. Accordingly, the amount of work is reduced.



Perfect adaptation to any situation guarantees uniform wall thicknesses.



Due to the high stability the position of the handle which has been determined will not be changed during the polymerization process.



The tray material has hardened after only 10 minutes in the Polylux unit.



The high stability of the tray material avoids deformation during impression taking. Precise models will be obtained.



The pink color offers the perfect basis for any type of set-up.



As a basic material for bite patterns or functional trays with bite rims, the resin ensures that the work will not be deformed.

### Resins and systems

### casting system

### Casting resin set



Flask 1 piece REF 360 0125 7

Sealing plugs 20 pieces REF 360 0125 8 Producing, finishing and polishing dentures with the least effort. The perfect flow characteristics and the well-matched processing times render Opti-Cast casting resin the first choice product. Due to the carefully selected raw materials, inaccuracy of fit is a matter of the past. The optimized material combination minimizes residual monomer.



Silicone plugs 3 pieces REF 360 0125 9

Small punching tube REF 360 0126 0



Mixing cup maxi 1 80 ml 1 piece REF 320 004M 1



Isoplast ip 750 ml REF 540 0101 9



Measuring cup liquid 25 ml, 1 piece REF 360 0126 2



Measuring cup powder 50 ml, 1 piece REF 360 0126 3



Bre-Gel BG 3 4 x 400 ml REF 540 0105 4

### Assortment uni.lign

#### 1 flask

- 20 Sealing plugs Silicone plugs
- Small punching tube
- Mixing cup maxi 1
- Measuring cup liquid

Measuring cup powder 500 g uni.lign

casting system powder 500 ml uni.lign

casting system liquid

30 ml Isoplast ip

REF 360 0126 4

### Assortment

### uni.lign

- 1 flask
- 20 Sealing plugs
- Silicone plugs Small punching tube
- Mixing cup maxi 1
- Measuring cup liquid

Measuring cup powder

500 g uni.lign

casting system powder

500 ml uni.lign

casting system liquid 30 ml Isoplast ip

REF 360 0126 8

### Bre-Gel BG 3 opaque liquid

Special duplicating gel for the economical production of dentures.



The model with the waxup is soaked.



To avoid low pressure when deflasking, the plug is inserted in the outside of the lower flask element.



The silicone plugs are pressed into the charging holes.



To ensure optimal positioning of the sprues, the model and the upper flask element are assembled as shown. A magnet in the base plate helps to hold the model.



Shake duplicating gel to obtain homogeneous consistency



and then melt in the microwave for 3 min at 600 to 800 watt.



# casting system



Stir Bre-Gel to achieve uniform heating. Melt two more minutes.



Excess pressure caused by boiling is avoided by the opened lid.



The flask is cooled in the cold water bath to 40-45 °C while stirring.



Pour Bre-Gel into the flask until the vents are slightly overfilled.



Final strength is reached after 45 minutes in the cold water bath.



The circumferential groove simplifies the removal of the upper flask element.



The model is carefully deflasked using compressed air.



The complex gingival model is reproduced in precise details.



Soak the model for 10 minutes before it is completed.



The charging hole and the vent are neatly punched with the small punching tube.



Before the teeth are placed back into the gel mould, they require circumferential ...



... and basal roughening with the setup grinding tool (REF 340 0101 0).



Thinly applied Isoplast (REF 540 0101 9) with short drying time ensures a perfect insulating film.



The model is placed back into the gel mould.



Die Silikonpfropfen verbleiben zur Stabilisierung bis zum Schließen der Küvette in den Einfüllöffnungen.



The flask is closed in the correct position using a centering snap.



The flask is placed onto the flattened lower flask element. Opti-Cast casting resin can now be poured in from above.





Durch ein Hin- und Herschwenken der Küvette treten Lufteinschlüsse aus.



A delay in polymerization can be achieved with cold water. This way resin can flow during the polymerization phase.



The resin is polymerized in the pressure pot for 30 min at 40 to 50 °C and a pressure of 2-6 bars.

# Resins and systems

### Multisil-Soft



The permanently soft relining system.

Multisil-Primer 5 ml REF 520 0100 4



reliable bonding

The bonding agent is matched to the bond of denture resin and silicone.



ready-to-use

The flow characteristics allow rapid processing and applying through the dosing device.



Silicone burs REF S187 QG 23 REF S263 QG 60 REF S237 QG 65



grindable

Silicone burs with a special cutting geometry simplify grinding of functional margins and transition zones.



Multisil sealing liquid 10 ml REF 520 0100 5



plaque-resistant

Multisil sealing agent avoids the accumulation of plaque to the surface and penetration of bacteria into the surface.



permanently elastic

Prolonged comfort of wear is a distinctive feature of the



tear-resistant

Highly cured materials create exceptional wear resistance and special tear resistance.



Dosing device REF 320 0044 0

### Assortment

Multisil-Soft

2 x 50 ml Multisil-Soft

in cartridges

5 ml Multisil-Primer 10 ml Multisil sealing liquid 12 pieces Mixing cannulas 1 piece Silicone burs

S237 QG 65

REF 540 0104 5



#### Refill packages:

50 ml Multisil-Soft cartridges REF 540 0104 6

5 ml Multisil-Primer REF 520 0100 4

10 ml Multisil sealing liquid REF 520 0100 5

12 St. Mixing cannulas yellow REF 320 0045 1

0104 5

# Ropak UV



Light curing acryliccolored opaque to coat CoCr objects.

Ropak UV F - liquid 10 ml REF 520 0016 4



Ropak UV P - powder 10 g REF 520 0016 5




The viscosity of Ropak UV can be adjusted to the individual requirements.



Mix powder and liquid on a mixing tray to obtain a homogeneous consistency.





Esthetics beyond compare – pink opaque shows perfection



Use disposable brush to apply the material. Ropak UV will coat the object even if it exhibits a thin consistency.



Apply thinly using the disposable brush; even dark metal elements will be coated in an aesthetic way.

### Ropak Kompaktopaker UV



The ready-to-use alternative for convenient coating of CoCr objects.

Ropak Kompaktopaker UV 20 ml REF 540 0013 3



Apply Ropak Kompaktopaker with the integrated brush directly onto the clean metal surface.



All metal elements applied with Ropak are perfectly covered.



The use of Ropak provides the future acrylic area with a more pleasant look.

# Kompaktopaker tooth-colored UV



To enhance esthe-tics in the area of acrylic teeth.

Kompaktopaker tooth-colored UV 10 ml REF 540 0010 5



Kompaktopaker toothcolored is particularly suitable for the anterior area.



The tooth-colored opaque that features a fine coating capacity is applied to the desired area.



This way perfect coating of the metal and thus esthetic restorations are obtained.



# **Processing acrylics**

#### Abraso-Gum Acryl processing set for acrylics

#### Ready to hand for minor Diatit bur adjustments of dentures.

bredent tungsten carbide burs and acrylic polishers are helpful tools for repairs, remove all tender spots, smoothen surfaces and produce high luster.



1 piece REF D200 KF 23





A smooth surface is achieved by exerting slight pressure.

Acrylic polisher coarse, green



6 pieces **REF P243 HG 10** 



The coarse acrylic polisher removes traces of the bur and shapes the surface.



5 pieces

- 1 Diatit bur D263 KG 60
- 1 Diatit bur D200 KF 23
- 1 Abraso-Gum Acryl coarse, green 1 Abraso-Gum Acryl
- medium, grey 1 Abraso-Gum Acryl fine, red

REF 350 0099 2

Acrylic polisher medium, grey



6 pieces REF P243 HM 10



The grey acrylic polisher features a slight abrasive capacity and smoothens the surface in a single working step.



6 pieces **REF P243 HF 10** 



The fine acrylic polisher produces a perfect highluster on all acrylic materials in next to no

### Set-up grinding tool



Set-up grinding tool REF 340 0101 0

Two grinding tools in one. Grinding without exchanging tools in a single working step

- quick adaptation of the underside of the tooth to be set up
- grinding in of occlusal stops

### Two grinding tools in one.



### occlusal



The small, precisely shaped grinding tip with fine, perfectly cutting diamond grains provides the ideal precondition for well-aimed and rapid grinding in of occlusal contacts.



The large grinding area with its optimized shape and selected natural abrasive diamonds ensures maximum removal of material and thus accurate and quick

#### Accessories:



Articulation paper holder size 1 1 piece REF 360 0121 7



Articulation paper holder size 2 1 piece REF 360 0122 0

# Diacryl grinding tool





Coarse grinding tool
1 piece
REF 340 0103 0




The special diamond grain size and the hollow shape of the grinding tool result in excellent grinding properties and ensure maximum cooling.



Margin grinding tool, round 1 piece REF 340 0106 0



Due to tapering in the middle of the grinding tool uniform margins of functional trays can be produced.



Due to the uniform coarse-grain diamonds with sharp cutting edges and the special shape, Diacryl grinding tools are perfectly suitable for finishing acrylic dentures and tray material in a quick and purposeful manner.



Universal grinding tool 1 piece REF 340 0104 0



Can be universally used for coarse and large papillae as well as for root bases.



Margin grinding tool, pointed 1 piece REF 340 0102 0



Recesses of labial and buccal frenula can be perfectly finished with this Diacryl grinding tool.



Papilla grinding tool 1 piece REF 340 0105 0



The fine, pointed flame design allows filigree finishing of alveolar and papillary bases.



Rubber grinding tool 1 piece REF 340 0090 0



Thanks to the fine grinding performance the object is prepared for polishing within a short time The rubber grinding tool is used instead of sandpaper.



Assortment

REF 340 0107 0

1 piece each



# **Processing acrylics**

### Tungsten carbide tools

### Tungsten carbide burs

For processing of acrylics.

Special types with cross cut for smooth surfaces on all acrylics.



Rapidy Microbur with relief 1 piece REF H001 NH 10



The microbur with relief ensures quick removal of material even at inaccessible spots.

### Diatit burs

With longer service life and increased grinding performance.



1 piece REF D194 KS 70



The coarse cross cut allows the quick removal of material across large areas.


A bredent bur that features Diatit wear protection reaches a degree of hardness of up to 3700 according to Vickers (HV).



1 piece REF D468 GG 16



Acrylic can be cut precisely and almost without any chips with the Diatit bur.





1 piece REF D237 KG 65



The coarse cross cut produces smooth surfaces in a quick and pressure-free manner.



1 piece REF D263 KG 60



Due to the universal design of the bur time-consuming exchange of tools is no longer required.

A hardness that is 100 % higher thanks to the Diatit wear protection results in a service life that is three times longer than the one of uncoated bredent burs.



1 piece REF D194 KG 23



The acute wedge angle of the individual cutting edges ensure precise milling with a high cutting performance.

Further information on burs in chapter 9!



1 piece REF D274 KG 60



The diagonal cutting edges of this bur allow to produce extremely smooth surfaces in no time

### Acrylic polishing set

The complete range of polishing products for acrylics – systematic polishing.





Abraso-Buff Acryl Ø 80 mm 1 piece REF 350 0078 0

High luster buff Acryl Ø 100 mm 1 piece REF 350 0082 0

for polishing of acrylics and metal 3 x 500 g REF 520 0016 0

Pumice polishing paste



Abraso Star K50 slightly abrasive 320 g REF 520 0016 1



The pumice polishing paste diffuses into the brushes and allows particularly long prepolishing.



The fine abrasive components of the pumice polishing paste simplify careful polishing of the acrylic denture material.

### Acrylic polishing set

REF 350 0084 0

1 x 150 g Abraso-Star K50 slightly abrasive 1 x 500 g Pumice polishing rpaste 1 piece Abraso-Soft Acryl

Abraso-Buff Acryl 1 piece 1 piece High luster buff Acryl

### High luster polishing with the handpiece



Polishing of acrylics with the handpiece. Brushes and buffs for handpieces provide brilliant high luster on all dental acrylics.



Star brushes goat hair white Ø 19 mm 15 pieces REF 520 0015 1



Cotton buff Ø 22 mm 15 pieces REF 350 0091 0

Leather buff

REF 350 0066 0

Ø 22 mm

15 pieces

Star brushes

Ø 13 mm

15 pieces

goat hair white

REF 520 0014 1



Velvet-soft small cotton threads polish palatal patterns excellently so that smooth surfaces are obtained to which coatings can not adhere.



The star design allows to increase the polishing performance up to 50 % and reduces the working time considerably.



All filigree areas are pre-polished with the smaller star brush. Perfectly suitable in the approximal area; protects acrylic teeth.



The dimensionally stable linen buff produces a mirror-like finish even on the hardest veneering materials.



Polishing with the leather buff avoids damage to thin transitions towards metal.



Accessories:

Abraso-Starglanz asg Universal high luster polishing paste 2 x 50 ml REF 520 0016 3

Acrypol polishing paste

for veneering materials

REF 520 0017 0

170 g



Polishing buff felt, three layers Ø 22 mm 15 pieces REF 350 0064 0



The three felt layers are perfectly suitable for any type of structure. Extremely fine polishing results are achieved.



### Polishing brushes



### Abraso-Soft Acryl

Due to the polishing heat the open-pore special fleece and the bleached Chungking bristles absorb more polishing paste and therefore up to 50 % of working time can be saved.

Unlike conventional brushes, the open-pore structure of the fibre fleece allows to take up considerably larger quantities of pumice or polishing paste. Accordingly, less polishing paste needs to be applied.

The fleece is able to absorb more air so that the polishing temperature is reduced and gentle polishing is ensured. Overheating of the surface is avoided.

Abraso-Soft Acryl Ø 80 mm 1 piece REF 350 0080 0



The combination of fibre fleece and bleached Chungking bristles let the pumice polishing paste diffuse deeply into the brush





Mixed pumice diffuses into the brush and the fibre fleece. The polishing agent remains longer on the brush and is gradually applied onto the surface in uniform quantities.

The brush hair are made of bleached Chungking bristles. Bleaching roughens the bristles, makes them softer and increases the absorbing capacity. On the one hand the rough surface holds the pumice paste more easily and on the other hand acrylic is polished more actively without overheating the surface.



Round polishing brushes with plastic core. Round brushes – Chungking white for abrasive polishing.

Due to its small width, the large brush is suitable for polishing areas that are difficult to access.



Chungking white Ø 80 mm 4 rows 12 pieces REF 350 0034 0



Chungking white Ø 65 mm 4 rows 12 pieces REF 350 0074 0



Chungking white Ø 70 mm
3 rows
12 pieces
REF 350 0030 0



Chungking white Ø 60 mm 3 rows 12 pieces REF 350 0075 0



Chungking white Ø 50 mm 2 rows 12 pieces REF 350 0027 0



Narrow brush

White goat hair with metal core for polishing that protects the structure.

The soft goat hair brush avoids the abrasion of the surface structure of acrylic teeth and thus simplifies polishing of approximal areas.

Narrow brush - white goat hair with metal core Ø 48 mm 10 pieces REF 350 0061 0



# Polishing buff



### Abraso-Buff Acryl

Three rows of high-quality Chungking bristles and special fabric liners guarantee prepolishing with an excellent result.

#### Abraso-Buff Acryl Ø 80 mm 1 piece REF 350 0078 0



The small width of the Abraso-Buff Acryl allows polishing in filigree approximal areas. Exchanging the brush is no longer required.





After finishing, the optimum prehigh-luster is easily achieved in no time.

The 2 x 2 special textile layers of the Abraso-Buff acrylic retain polishing pastes or pumice considerably longer than conventional brushes. They gradually spread abrasive materials and thus simplify polishing.



### Prepolishing buff Acryl

Aggressive polishing behaviour – working time reduced.

The stable layers of the prepolishing buff Acryl consist of silicone-treated linen. Accordingly, particularly aggressive polishing is possible.

# Prepolishing buff Acryl Ø 80 mm

24 layers 1 piece REF 350 0099 1

Ø 60 mm, 24 layers 1 piece REF 350 0098 0



The buff allows time-saving polishing at reduced temperatures.



Due to the different sizes palatal areas can also be polished easily.



# Polishing buff



High luster buff Acryl
No formation of fuzz and only reduced evolution
of heat.

# High luster buff Acryl 1 piece each,

Ø 60 mm, 40 layers REF 350 0094 0 Ø 100 mm, 35 layers REF 350 0082 0



The outer layers with reinforced fibres provide the buff with a stability never achieved before.



The 35 resp. 40 layers of linen have been welded with an ultrasonic unit to protect them against twisting and produce a unique high luster due to the high stability.





The particularly loosely woven linen ensures circulation of air during high luster polishing so that overheating of the acrylic is avoided. This results in a very gentle polishing process.

Special linen consequently avoids excessive evolution of heat on the acrylic surface.



**Leather buff**Produces high luster in a quick and gentle manner.

### Leather buff for acrylics

1 piece each

Ø 80 mm, 5 layers REF 350 0036 0 Ø 100 mm, 5 layers REF 350 0035 0

### Leather buff for metal

1 piece each

Ø 60 mm, 9 layers REF 350 0099 0



Acylics can be polished at lower temperatures in a very gentle manner using the leather buff for speeds up to 1500 rpm.



Polishing at lower temperatures produces a high luster even in the approximal area so that coatings will not remain there.



The leather buff produces perfect high luster without any retentions for bacteria. This way cleaning of dentures is simplified.




The strong rise in allergies and intolerances to dental alloys, amalgams, residual monomers, and the resulting allergic reactions have led to a rethink in the choice of dental materials for restorative dentistry.

The use of highly-compatible thermoplastic resins enables the manufacture of high-quality, metal-free and biocompatible dentures. An appropriate thermoplastic is available for every area of application of dental technology.

Using the *fol* 2 press vacuum press system, the highperformance polymer BioHPP is produced for the manufacture of highly-stable framework structures for permanent fixed and removable dentures. The elasticity of this high-performance polymer has been specially adapted to human bones.

There are five different thermoplastics available in a variety of colours for the thermopress 400 injection moulding system. Thanks to these different types of plastic, there are many areas of use.



### for 2 press vacuum press system

System steps	35 I
for 2 press vacuum press device	352
for 2 press - The system components	353
for 2 press and BioHPP	354
fixed dentures	355
removable dentures	356
Implantology	356

### thermopress 400 plastic injection moulding system

System steps	357
thermopress 400	358
bre.crystal HP	360
bre.dentan HP	
Bio Dentaplast	361
bre.flex	
bre.flex 2 <sup>nd</sup> Edition	362

# BREDENT GROUP ACADEMY INTERNATIONAL



Flexibility is the concept that sets the bredent group academy apart from other training providers on an international level! The modular structure of the individual courses is based on indications – whether related to dental technology or dental medicine. This offers you an ideal combination of logically coordinated and interlinked modules relating to fixed prostheses, removable prostheses, combined prostheses, prostheses on implants or natural teeth, and also the issue of prevention.

All modules can be booked separately and highlight the individual advanced modules as well as transverse solutions for each module. This diverse and self-contained concept is an option for anyone interested in further and personal development.

This flexibility allows each bredent partner to respond to the wishes and knowledge of their clients and further develop these according to their level of knowledge. What is more, the concept is highly targeted and seeks to

achieve the greatest impact possible by integrating the full range of bredent and bredent medical products, systems and treatment concepts. Within the individual modules, the aim of the course, the course description and the materials list can be viewed in detail by potential participants in the form of text and images.

The benefit lies in the fast and easy navigation through the concept. The bredent medical partners will quickly identify the indication-related path and bredent partners will identify a product and system-related route.

Courses provided on site are delivered by qualified and trained bredent group partners. Those partners will be familiar with all products, systems and treatment concepts within the bredent group academy, and will be regularly trained and kept up to date.

Realise your own unique potential with the bredent group academy!



### for 2 press vacuum press system



A wax model is created using a plaster or plastic master model under typical laboratory conditions. The dentist should carry out a tangential or chamfer preparation, just as with ceramic crowns or bridges. Fully anatomical occlusion surfaces are possible because of the white colour of the BioHPP. The framework structures can be veneered with traditional composite veneers.

In the for 2 press vacuum press system, metal-free and biocompatible dentures are manufactured to the highest standard. The range of indications stretches to include fixed and removable dentures with all kinds of friction elements. The material to be processed is a high-performance polymer "BioHPP", which has been specially developed as a framework material for the field of dentistry. BioHPP can be veneered with all composite veneers through the use of a special bonding agent, visio.link. The purchase prices for the whole systems are very low, meaning the investment will pay for itself, even if only a small number of jobs are carried out.



Once the specific wax weight of the wax model has been determined, and the minimum amount of Bi-OHPP necessary for the pressing procedure has been established, the wax model is invested in a special silicone mould. After a setting time of 20 minutes, the mould can be placed directly together with the disposable press plunger in the pre-heating oven.



The investment material mould can be heated conventionally or using high-speed heating. The BioHPP is also melted in the pre-heating oven at 400°C within 20 minutes. After that, the mould is used with the attached disposable press plunger in the for 2 press vacuum press device. By manually raising the pressing table, the pressing procedure is triggered and automatically completes within 38 minutes. After that, the mould is cooled at room temperature and devested.



The BioHPP can very easily be processed with little contact pressure using cross-toothed carbide mills. The speed of rotation may not exceed 8,000 rpm. Additional mechanical retention for a necessary adhesion can be applied in the form of retention beads or splitters. However, this is not a stringent requirement.



The use of the bonding agent visio.link is crucial for a high level of adhesion. This must also be used for composite veneers foreign to the system. After applications, visio.link is polymerised in a light-hardening device (e.g. bre.lux) for 90 seconds. Veneering then follows afterwards.



The BioHPP can be very well polished by using the right rotating instruments.

1st step: Carbide mill (cross-tooth)
2nd step: Diagen-Turbo-Grinder

- 3. Ceragum rubber-polishing cylinder
- Goat-hair brush with pumice stone powder (0-90μm particle size)
- Goat-hair brush with high-gloss polishing paste (Abraso-Starglanz)

#### The materials used

for 2 press vacuum press device	page 352
for 2 press - the system components	page 353
for 2 press and BioHPP	page 354
fixed dentures	page 355
removable dentures	page 356
implantology	page 356

bredent

### for 2 press vacuum press device



### for 2 press vacuum press device for the manufacture of metal-free and biocompatible dentures.

Using the for2 press vacuum press device, the thermoplastic high-performance polymer BioHPP is processed into a mould made of phosphate-bonded investment material. The end result is a metal-free, white framework structure, which bonds extremely well with traditional composite veneers individually and in a highly-aesthetic manner. Dentures manufactured from BioHPP are certified for permanent

The melting procedure for the BioHPP high-performance polymer is carried out in your pre-heating oven that is already available. The subsequent pressing procedure is completed fully automatically in a vacuum.

for 2 press device REF 140 0060 0

### for 2 press Basic Set


1 x for 2 press mould (consisting of mould plates) 3/16 mm and silicone ring, 35 x 210 g Brevest for 2 press EBM incl. 2 litres Bresol for 2 press Liquid, 25 x for 2 press filler 16 mm (disposable press plunger for pressing the materials into the mould) 20 g BioHPP, 1x processing instructions REF 140 0060 1



Once the investment material mould has been pre-heated in the pre-heating oven, and the wax or plastic of the model has melted away, the BioHPP is also melted at 400°C in the pre-heating oven. Immediately after that, the mould is placed in the for 2 press vacuum press device with the melted BioHPP and the mounted disposable press plunger. By raising the lift, the automatic pressing procedure is triggered in a vacuum atmosphere.



After the 3-minute pressing procedure including vacuum admission, the cooling process begins while maintaining the pressing pressure. This is the only way to guarantee that the material properties of Bio HPP are fully exploited.



After a total of 35 minutes, the entire pressing procedure is ended and is signalled to the user by means of an acoustic signal and an optical LED display. The Bio HPP framework structure can then be immediately devested and processed further.



Example of application on the basis of a bridge with five splinted individual crowns from BioHPP veneered with viso.lign veneers.

Image: Zahntechnischer Meisterbetrieb Harald

#### Technical data for 2 press

Power supply Power consumption Venturi nozzle vacuum performance Weight Dimensions (W x H x D)

approx. 760 mbar 13 kg 250 x 600 x 290 mm IP 34

Protection class Noise level Input air pressure

< 70 dB min. 4.5 to max. 6 bar

90 - 250 Volt, 50 - 60 Hz

Fuse

T 2.5 A

15 Watt

### for 2 press - The system components



### The material BioHPP

Bio HPP is sealed against humidity and comes in a transparent tube. The required quantity can be dosed by the gram, looked up from the wax to Bio HPP conversion chart. Only that amount of material

### BioHPP

REF 540F2PB2 20 g 100 g REF 540F2PB3 10 x 15 g Pellet for mould XXII REF 540F2PB4



### Mould system for 2 press mold

Consisting of sprueing base and silicone ring. Available in 3 different: 3, 9 and 9 XXL.

Sprueing base and silicone ring

Size 3 REF 360F2P16 REF 360F2P20 Size 9 Size 9 XXL REF 360F2P30



Silicone ring suitable for the forl 2 press mould. Easy and quick demolding of set investment. Checkered inside, enlarged surface for effective evaporation during pre-heat process.

### Silicone ring

Size 3, REF 360F2PR3 Size 9, REF 360F2PR9\*

\* Also suitable for the 9 XXL mould plate



#### Sprueing base single

REF 360F2PT1 Size 3 REF 360F2PT2 Size 9 Size 9 XXL REF 360F2PT3\*

\* Only for processing the BioHPP pellet



#### Disposable plunger for 2 press filler for safe pressing results

Disposable plungers for pressing the high performance polymer into the mould.

Rounded at one end for smooth forward movement. Syncronized cooling process between plunger and investment material, no risk of chipping or cracks of polymer. Pressure resistant.

### Disposable plunger

25 pieces 16 mm REF 570F2P16 REF 570F2P20 25 pieces 20 mm 14 pieces 30 mm conical REF 570F2P20\*

\* Only for processing the BioHPP pellet



#### Investment brevest for 2 press

Micro grain special investment for the for 2 press System.

Appropriate for speed or conventional burn out processes.

### Brevest for 2 Press

Box of (approx. 7,35 kg) 35 x 210 g bags incl. 1000 ml Bresol for 2 press

REF 570F2PV1



### Bresol for 2 press

Investment liquid Brevest for 2 press



REF 520F2PL1



Accessories:



"Generation M" relief burs REF H274 M5 16



"Generation M" relief burs REF H263 M5 40



Diagen-Turbo-Grinder REF 340 0020 0



Ceragum

Rubber polishing cylinder **REF PWKG0600** 



Goat-hair brush REF 350 0061 0



Abraso-Starglanz REF 520 0016 3



Acrypol REF 520 0017 0



Wax casting bulbs REF 430 0144 7



### for 2 press and BioHPP

### Advantages and benefits of BioHPP

### Reproducible manufacturing process

Advantage

Consistent quality thanks to automatic and electronically-monitored pressing procedure

Benefit

Consistent material properties and avoidance of complaints

### Shock-absorbing effect (peak-off)

Advantage

Protection of the implant against high masticatory forces

Durability and increased comfort for the patient

### Abrasion-resistant tooth-like material



Occlusal surfaces keep their shape over long periods of use

Benefit

Increases quality of life

### White framework material that can be veneered



Can be veneered individually with composite

veneers

Benefit

Enables individual adjustment to the remaining teeth and prevents chipping

### Low density (1.3 to 1.5 g/cm<sup>3</sup>)

Advantage



Very simple dentures

Benefit



Increases comfort for patients

### Steady friction in connecting elements

Advantage



Prevents loss of friction

Benefit


Increases comfort and saves on replacing the

### Homogeneity



Uniformly distributed fillers in the semicrystalline polymer matrix

Benefit



### **Biocompatibility**

Advantage

No harmful substances such as metals or residual monomers are released

Benefit



Offers body-friendly and healthy dentures

### Material properties particular to BioHPP

#### Mechanical properties according to ...

**DIN EN ISO 10477** 

E-modulus Flexural strength 4,000 MPa >150 MPa

(no material failure)

Water absorption 6.5 μg/mm<sup>3</sup> Water solubility – < 0.3 μg/mm³</p> Thermocycling 10,000 cycles 5°C / 55°C in accordance with DIN EN ISO 10477

E-modulus - 4,000 MPa

Flexural strength - >150 MPa (no material failure)

Breaking load tests on 3-part bridges

Maximum load without failure

Maximum load without failure

>1,200 N

(after 24 hr immersion in water, 37 °C)

>1,200 N

(after mechanical and

thermal alternating load

1.2 million x 50 N, 10,000 x 5 °C / 55 °C)

Other

properties

Melting range (DSC) Bond strength

- approx. 340 °C > 25 MPa

Density - 1.3 to 1.5 cm<sup>3</sup>

Hardness (HV) - 110 HV 5/20



3-part anterior bridge with a very high-quality aesthetic veneer.



### Indications of BioHPP - fixed dentures

### Posterior area



To enable the manufacture of an aesthetically-pleasing and durable prosthetic with BioH-PP, the bond between BioHPP and the composite veneer is of the utmost importance.



The visio.link adhesive from the visio.lign veneer system incorporates all outstanding adhesive properties to materials such as PMMA and composites. In this way, visio.link is a primer and bonder in one.

### Anterior area



The aesthetic zone in the anterior area offers the best possible range of indications for this biocompatible and metal-free framework material.



As far as aesthetics are concerned, this bridge is comparable to a zirconium oxide bridge.



The use of BioHPP dentures means that a flexible material adapted to bone is used in the mandible.

### Molars



Thanks to its colouring, BioHPP can be used as a fully-anatomical restoration in the molar area and offers an ideal environment thanks to the abrasion resistance particular to the material.



The BioHPP bridge is securely fitted to the uni. fit adhesive sleeves using DTK adhesive.

#### Red-white aesthetic



Using visio.link adhesive, it is possible to manufacture gum fillings.

Photos: Dental-Labor Schwindt, Landau/Pfalz



### Indications of BioHPP - removable dentures

### Bridgework



The BioHPP offers outstanding friction properties in combination with metallic and ceramic primary structures. As far as the patient is concerned, the high comfort and ease of insertion and removal stand to the fore. The inert properties preclude any adverse effects.

### Telescopic bridges/crowns



There is no loss of friction thanks to the flexibility of the material and the ease with which the secondary structure can be veneered. No abrasion to the primary part.

Photo: Dentallabor Fiedler, Neulußheim

# Indications of BioHPP - implantology

### Individual abutments



The individual abutment "SKYuni.HPP" has been developed for bredent medical's SKY and blueSKY implants. In the future, additional abutments will also become available for other companies' implant systems



The SKYuni.HPP is blasted with aluminium dioxide with a grain of 110µm and max. 3 bar pressure before wax modelling.



Final wax modelling on the SKYuni.HPP



The wax model is sprued with the SKYuni.HPP on the base mould of the *for* 2 press system. A casting bulb is used for this.



The individual abutment is over-pressed with BioHPP.



The final abutment can now be veneered or fitted with a ceramic crown directly or finished with a BioHPP crown and a visio.lign veneer.

### thermopress 400 plastic injection moulding system



The master model prepared for the silicone duplication is secured in the middle of the duplicating flask. The mixed duplicating silicone Exaktosil N21 is then poured into the flask without any bubbles. The edges of the lining wax should be cut off in undercuts like in a watch seam. This provides an outstanding mechanical bond for the resin saddle material.

Various thermoplastic materials can be processed using the thermopress 400 injection moulding system. The advanced injection moulding device has been adjusted to the ever-higher material-specific requirements. Practical implementation in the manufacture of metal-free dentures can be achieved quickly for the dental technician.



The silicon duplicate can be cast immediately after the hardening period using the expansion plaster Expando-Rock. Here, it is important to pay attention to the necessary mixing ratio of the expansion plaster Expando Rock for the thermoplastic resin to be used. These mixing ratios can be found in the user manual. In the thermoplastic bre.crystal HP, a class 4 plaster, Exakto-Rock S, is used instead of the Expando-Rock.



The wax modelling should be done using modelling waxes with low melting temperatures. The edges of the resin saddles should be moulded for mechanical retention and an undercut end strip in the form of a watch seam.



The duplicate model is invested into the lower half of the flask with a type 3 dental stone. In order to easily separate both halves of the flask for the wax boil out without any damage, the undercut area should be blocked off with plaster. During the injection process, the fused plastic is guided into the wax model via a 10 mm injection sprue and then a 1.5 mm thick film sprue.




All thermoplastic materials can be processed very easily using cross-toothed carbide mulls, such as mills for silicone, for example. The surface is then smoothed using sandpaper cloth and pre-polished on the polishing lathe with pumice stone. Finally, the surface is polished to a high gloss with a cloth buffing wheel using a little high-gloss polishing paste.



After finishing, the thermoplastic framework structure, here in the form of a telescopic bridge, can be finished with gum-coloured resins. Biocompatible and residual monomer-reduced bre.crystal HP can likewise be used by means of a second investing and injection.

#### The materials used

thermopress 400	page 358
bre.crystal HP	page 360
bre.dentan HP	page 360
Bio Dentaplast	page 361
bre.flex	page 362
bre.flex / bre.flex 2 <sup>nd</sup> Edition	page 362





### thermopress 400 plastic injection moulding system

### thermopress 400



#### Injection moulding device for processing thermoplastic resins with a melting temperature up to 400°C

- No additional equipment such as CO₂bottles or a compressed air supply are required. This saves on additional costs and time. Consistent quality due to avoidance of a drop in pressure.
- High operational comfort when handling the device.
- Time-saving function thanks to simultaneous operation of both heating chambers.
- The injection process can only be performed with the lid closed, providing additional safety.

Convenient removal of the flask thanks to automatic cartridge ejection if the bracket is unlocked.

#### thermopress 400

- 1 unit with power cord
- 2 allen keys
- 1 cleaning brush 1 special tool
- REF 110 0040 0

The melting process is specially adapted to the demands of plastic thanks to high-performance heating elements. In this way the best possible material properties are



The forces of up to 2 kN exerted on the special flask are accommodated in the device by means of a special input and fastening system. In this way, it is ensured that the fused thermoplastic resin is free of air pockets and is injected into the hollow part of the flask in the correct shape.



All 6 different processing parameters are stored in the device software. In total, up to 30 melting programmes can be used in the device in a userfriendly manner.

#### Technical data thermopress 400

Width 650 mm Height 250 mm Depth 300 mm Weight 40 kg 220 - 230 V Voltage Power 0.5 - 1.6 kW max. 2.2 kW Accessories thermopress 400:

REF 140 0090 4 1 press-out device and punch\* 1 pair of cartridge pliers\* REF 140 0090 6 REF 140 0091 2 1 flask hook with hex\* 1 flask N small\* (I 122 mm, w 102 mm, h 72 mm) REF 140 0N90 3 1 flask N large (I 140 mm, w 102 mm, h 72 mm) REF 140 0N90 5 REF 110 0040 2 1 cleaning brush 1 thermopaste 400 special paste, 50 g\* REF 540 0105 1 Expando-Rock-Set 5 kg expansion plaster, 500 ml Expandosol REF 570 OERS 5

### thermopress 400 accessories assortment (Unit)

7 pieces

REF 110 0040 1

View above, accessories thermopress 400: articles marked with \* are included in thermopress 400 accessories assortment!

**Brochure Patient Information** thermopress sensitive REF 000 415G B

thermopress 400 Patienten passport REF 000 628G B

thermopress 400 Information Brochure REF 992 945G B

thermopress 400 Patient Information REF 000 135G B

Training program 2-day courses (Senden) REF 950 0020 0

#### thermopress introductory assortment (20 cartridges)

5 x 30 q bre.crystal HP pink 1 2 x 20 g bre.flex pink bre.crystal HP clear 2 x 16 g 3 x 24 g bre.dentan HP A bre.flex 2nd Edition Bio Dentaplast A2 5 x 24 q 3 x 16 g pink veined

REF 140 0090 2

### thermopress 400 assortment S1 (20 cartridges)

Bio-Dentaplast A3 bre.crystal HP pink 1 1 x 20 g 2 x 30 g bre.crystal HP pink 1 x 20 g Bio-Dentaplast A3 2 x 24 g bre.flex 2nd Edition 1 x 250 ml Acryl Sep 5 x 24 g pink veined 1 x 20 ml Die varnish, light-curing bre.crvstal HP clear opaque 1 x 24 a bre.dentan HP A 1 x 50 g thermopaste 400 REF 540 S000 1

### thermopress 400 assortment S2 (20 cartridges)

5 x 16 g Bio Dentaplast A2 5 x 16 g Bio Dentaplast B3 5 x 16 g Bio Dentaplast A3 1 x 250 ml Acrylic Sep acrylic 5 x 16 g Bio Dentaplast B2 separating liquid REF 540 S000 2



### thermopress 400

Thanks to the large choice of thermoplastic resins, the thermopress 400 device can be used for almost all areas of application.



#### bre.crystal HP

features long-term stability, provides a dense and hence smooth surface.

This results in enhanced comfort of wearing of full dentures.

- Almost no residual monomer high biocompatibility
- limited water absorption constant suction effect, lasting precision of fit
- available in the shades: crystal clear, pink 1, pink 2, pink 3, pink stippledh
- Melting temperature 280° C



#### bre.dentan HP

is an industrially polymerised thermoplastic; this increases the fracture resistance and biocompatibility in temporary crowns and bridges.

- can be veneered with conventional C+B resins
- available in the three common dentine shades A, B, C
- Melting temperature 280° C



#### bre.flex

Prosthetic base material for partial dentures.

The indication range also includes splints and sports mouthguards.

- available in the shades: translucent, pink 1, pink 2, pink 3 and tooth shade B
- Melting temperature 222° C

#### bre.flex 2nd Edition

- In the colours transparent, pink 2 and pink veined
- Melting temperature 280° C



#### Bio Dentaplast

Clasps and attachments which are normally made of metal can be produced using tooth-colored Bio Dentaplast. The esthetic appearance of teeth at which retaining clasps have been attached is improved.

Additional indications are

- crowns and bridges (temporary)
- telescopic and attachment work
- tooth-colored clasps
- shades A1, A2, A3, B2, B3 based on the VITA shade guide
- to be processed at 220° C



### bre.crystal HP



#### bre.crystal HP (Thermoplastic resins - partial and full dentures)

color	Quantity	REF
crystal-clear	20 x 24 g	540 OP32 4
crystal-clear	20 x 30 g	540 OP33 O
crystal-clear	1 x 500 g	540 OP30 5
pink 1	20 x 24 g	540 OP12 4
pink 1	20 x 30 g	540 OP13 O
pink 1	1 x 500 g	5400P10 5
pink veined	20 x 24 g	540 OPO2 4
pink veined	20 x 30 g	540 OPO3 0
pink veined	1 x 500 g	540 OPOO 5
pink 2	20 x 24 g	540 OP22 4
pink 2	20 x 30 g	540 OP23 O
pink 2	1 x 500 g	540 OP20 5
pink 3	20 x 24 g	540 OP42 4
pink 3	20 x 30 g	540 OP43 O
pink 3	1 x 500 g	540 OP40 5

#### Esthetically appealing thermoplastic resin.

	bre.crystal HP	Heat-/cold-curing resin
Shrinkage	0.6 %	5-7 %
Water absorption	0.35 %	2 %
Residual monomer	0.2 %	2-7 %

#### Wax model



The flexible plastic casting channel wax with the required diameter of 10 mm is fed via a 1.5 mm thick wax plate (film sprue) into the wax model.



Dentasil tooth protection silicone (REF 520 0029 6) allows to achieve a high final hardness (Shore 65) so that the teeth are perfectly protected against the high press-in pressure.

#### Investing



The wax model is invested into the flask using class III stone and then boiled out.

#### Retentions



The roughened acrylic teeth with retentions are wetted with the Haftconnector (bonding agent) for 5 min. The circular retention produced using Vb2 ensures safe hold.

#### Polishing



The highly compressed bre.crystal HP denture can be processed and polished just like conventional acrylics.

# bre.dentan HP



Thermoplastic resin featuring high resistance to fracture for tooth-colored crowns, bridges and long-term temporary restorations - simple, fast and favorably-priced.



Quantity bre.dentan HP A 20 x 16 g bre.dentan HP B 20 x 16 g 3 tooth shades (crowns and bridges)

bre.dentan HP in

REF 540 ODA1 6 540 ODB1 6 bre.dentan HP C 20 x 16 g 540 ODC1 6



Outstanding material properties for accurate positioning and perfect fit. Ideal to prepare exact, high-quality temporary restorations.

### Bio Dentaplast



.... offers a wide range of applications in the area of attachment and chrome cobalt work. Reinforced pre-bent clasp pattern for resin injection moulding. No application of additional wax required, hence time is saved and correct cross-section design is ensured.

Expando-Rock-Set 5 kg expansion plaster 500 ml Expandosol REF 570 OERS 5



Premolar clasps, bent for resin injection moulding 10 sheets with 10 clasps each left + right REF 430 0748 5

#### **Bio Dentaplast**

(tooth shade based on the Vita shade guide)

	Quantity	REF
A1	20 x 16 g	540 BA11 6
A1	20 x 20 g	540 BA12 0
A2	20 x 16 g	540 BA21 6
A2	20 x 20 g	540 BA22 0
A2	1 x 500 g	540 BA20 5
А3	20 x 16 g	540 BA31 6
А3	20 x 20 g	540 BA32 0
А3	1 x 500 g	540 BA30 5
B2	20 x 16 g	540 BB21 6
B2	20 x 20 g	540 BB22 0
B2	1 x 500 q	540 BB20 5
В3	20 x 16 g	540 BB31 6
В3	20 x 20 g	540 BB32 0
В3	1 x 500 g	540 BB30 5

#### Preparatory work (for a telescopic prosthesis, for example)



The alveolar ridges for the resin saddle are covered with preparation wax (REF 430 584 0). The finished edges should be cut off like a watch seam so that the saddle resin can be seated mechanically. The model is then cast using Exaktosil N21.

#### Duplicate



So that the cast primary parts are not damaged by the high injection forces, metal reinforcements should be placed in the primary parts straight after pouring in the Expando-Rock.

Model




The wax modelling should be carried out using softer modelling waxes with low melting temperatures. Points and sharp edges must be avoided without fail. Due to the extremely simple processing of the thermoplastic resin, the wax model can turn out stronger.

#### Investing



In order to obtain the best reproducible material properties, a film sprueing should be carried out. This will ensure that the liquid resin is guided into the framework geometry in the shortest and quickest way.

#### Separating



Both halves of the flask must be separated from each other using a special separating liquid, Acrylic Sep (REF 520 0029 1). It is important to ensure that both halves of the flask are first screwed together shortly before the injection process and placed in the thermopress plastic injection moulding machine. Otherwise there is the risk of blistering in the moulded objects.

#### Finishing



All thermoplastic materials can be processed very easily using cross-toothed carbide mulls, such as mills for silicone, for example. The surface is then smoothed using sandpaper cloth and pre-polished on the polishing lathe with pumice stone. Finally, the surface is polished to a high gloss with a cloth buffing wheel using a little high-gloss polishing paste.

#### Finishing



After finishing, the thermoplastic framework structure, here in the form of a telescopic bridge, can be finished with gum-coloured resins. Biocompatible and residual monomer-reduced bre.crystal HP can likewise be used by means of a second investing and injection.

### bre.flex and bre.flex 2nd Edition



bre.flex is a flexible, highly compatible polyamide and has proved its suitability for dentures for allergic persons.

Silicone burs are perfectly suitable for processing bre.flex.

bre.flex and bre.flex 2nd Edition (Telescopic prostheses, splints and sports mouthguards, chrome cobalt and attachment work)

bre.flex		
color	Quantity	REF
translucent	20 x 24 g	540 0F12 4
translucent	20 x 20 g	540 0F12 0
translucent	20 x 16 g	540 0F11 6
translucent	1 x 500 g	540 0F10 5
tooth shade B	20 x 24 g	540 0F22 4
tooth shade B	20 x 20 g	540 0F22 0
tooth shade B	20 x 16 g	540 0F21 6
tooth shade B	1 x 500 g	540 0F20 5
pink 1	20 x 24 g	540 0F02 4
pink 1	20 x 20 g	540 0F02 0
pink 1	20 x 16 g	540 0F01 6
pink 1	1 x 500 g	540 0F00 5
pink 2	20 x 24 g	540 0F42 4
pink 2	20 x 20 g	540 0F42 0
pink 2	20 x 16 g	540 0F41 6
pink 2	1 x 500 g	540 0F40 5
pink 3	20 x 24 g	540 0F32 4
pink 3	20 x 20 g	540 0F32 0
pink 3	20 x 16 g	540 0F31 6
pink 3	1 x 500 g	540 0F30 5

Colour Quantity pink 2 1 x 500 g 540 0F50 5 540 0F51 6 pink 2 20 x 16 g pink 2 20 x 24 g

bre.flex 2nd Edition

540 0F52 4 540 0F60 5 1 x 500 g pink veined pink veined 20 x 16 g 540 0F61 6 540 0F62 4 pink veined 20 x 24 g clear 1 x 500 g 540 0F80 5 20 x 16 g 540 0F81 6 clear clear 20 x 24 g 540 0F82 4


Accessories:



Aluminium Catridges empty 18 pieces REF 540 KL01 8



Diatit-Multidrill 1.5 Ø x 8 mm REF 330 0073 0



Silicone burs REF S187 QG 23 REF S263 QG 60



Diamond point Veneering technique Vb 2 1 piece REF 340 0083 0



Plaster insulating agent 750 ml REF 540 0013 5



Austragegerät REF 320 0044 0



Exaktosil N 21 Componente A 1000 ml Componente B 1000 ml REF 540 0114 7



Abraso-Gum Acryl 6 pieces REF P243 HG 10



Abraso-Gum Acryl 6 pieces REF P243 HM 10



Rodeo round brushes 15 pieces, Ø 21 mm REF 350 0097 0



Catridges empty 18 pieces REF 540 KL01 8



Qu-resin pink 50 ml cartridge REF 540 0116 5



Acrylic Sep Acrylic/plaster separating liquid 250 ml REF 520 0029 1 750 ml REF 520 0029 4



Dentasil tooth protection silicone 2 x 50 ml with 24 mixing cannulas size 2, yellow REF 520 0029 6

Mixing cannulas size 2, yellow 12 pieces REF 320 0045 1



Flexible acrylic sprue wax 275 g, Ø 10 mm REF 430 0741 0

**Bonding agent** REF 520 0029 2



Snoring has become the number 1 endemic illness. Extreme snoring disrupts restorative sleeping in 50 % of all bedrooms – it can reach a volume of 90 decibels, which is the same as the noise produced by a lorry driving by.


Two forms of treatment are available. The protrusion splint, which is worn in the case of slight or medium obstructive sleep apnoea syndrome, and biofunctional treatment with a vacuum activator, are used in the case of primary snoring without obstruction.



Snoring therapy  <i>  Sleep</i> Plus mandibular advancement device	
Snoring therapy	364
SleepPlus mandibular advancement device	364
Snoring therapy / Biofunctional therapy	
Biofunctional therapy	366
Vacuum activator silencos	366
Vacuum activator silencos kids	367
Silencos night device	367

### Snoring therapy / SleepPlus mandibular advancement device

### Snoring therapy




#### Snoring - a common disease

Extreme snoring (up to 90 decibel) disrupts restorative sleeping in 50 % of all bedrooms! With increasing age, more than half of the population suffers from snoring problems:

- 40 60 % experience simple snoring without disrupted breathing (primary snoring)
- 5 10 % suffer from pathologic snoring with health risks caused by blockage of the airway (obstructive sleep apnea)

But snoring does not simply generate annoying sounds or sleep disruption – it may cause serious diseases, such as high blood pressure, stroke or heart attacks.




The efficiency of intraoral mandibular advancement devices to combat minor or medium sleep apnea syndromes is already established through several studies. Their use is recommended by the Deutsche Gesellschaft Zahnärztliche Schlafmedizin – DGZS (German Society of Dental Sleep Medicine) for this area of indication. Specific dental, oral and functional diagnostics are required prior to use.

The individually designed intraoral mandibular ad-vancement device *SleepPlus* keeps the lower jaw toward the median line at night, and thus prevents from snoring and reduces a disrupted breathing significantly.

### SleepPlus mandibular advancement device



A system for the treatment of minor or medium obstructive sleep apnea syndrome (OSAS).





#### Mandibular Advancement Device

The mandibular advancement device SleepPlus is a novel and comfortable treatment device with an especially high wearing comfort to combat obstructive snoring efficiently.

This advanced development of the already known mandibular advancement device offers a very good lateral and vertical mobility, due to telescopic ball and socket joints. Since the telescopic rails of the mandibular advancement device *Sleep*Plus are occlusally integrated, the buccal mucosa will not become irritated.

A temporomandibular joint support, which is modeled into the synthetic material, provides an additional relief of the musculature during sleep. The design of the mandibular advancement device *SleepPlus* was developed in the dental practice under clinical conditions and warrants high patient acceptance. The treatment concept has proven itself for several years. The mandible advancement device *SleepPlus* allows you to expand your range of products and to supply your patients with an established and proven treatment device that was developed in the dental practice.

The unique comfort of *Sleep*Plus is going to make your laboratory even more successful! Offer your dentists a product that makes for satisfied patients.

#### We are going to support you with:

- Workshops for practical production (certification)
- Advanced training courses about snoring therapy (certification)
- Marketing support through
- Patient flyers
- Waiting room posters
- Patient information on the internet, including reference to certified den tists and dental laboratories, www.sleepplus.de





# Snoring therapy / SleepPlus mandibular advancement device

### Sleep Plus mandibular advancement device



### Sleep Plus

#### Mandibular Advancement Device

Ball and socket joints on telescopic rails provide high lateral and vertical mobility.

#### The system components and their advantages:

- The telescopic rails are delivered with individu ally adjustable springs, which are variable to suit any bite situation.
- Ball and socket joints allow for high three-dimensional mobility.
- Occlusally incorporated telescopes offer optimal wearing comfort. The buccal mucosa will not
- The bite elevation relieves the temporomandibu lar joints at night.
- The optional frontal immobilization keeps the mandible in a comfortable position.



Visit our workshop and get acquainted with the professional production of the intraoral mandibular advancement device SleepPlus.

The clinically required construction criteria and system components match each other, and the production methods are conveyed by qualified instructors.

Ask for current course offers.

#### Your benefit:

- Acquisition of new customers
- Expansion of offers for the private insurant

Our office and field consultants would be pleased to inform you about further details.

**Sleep**Plus Waiting room poster 2 pieces REF OPOOO8GB

manufacturing manual

**Sleep**Plus

REF 000322GB

1 piece

Patient flyer 20 pieces REF 000 299G B

SleepPlus

#### Assortment


- 2 Telescopic rails
- 2 Pressure springs
- 1 Key
- 1 Spherical-head screw
- 1 Matrix, green
- 1 Matrix case
- 1 Threaded sleeve

REF 580 0119 0

#### Refill package:

10 Telescopic rails	
Incl. 10 pressure springs	REF 580 0119 1
10 Pressure springs	REF 580 0118 1
1 Key	REF 580 0119 2
1 Spherical-head screw	
vks-oc/sg 2,2 titan	REF 450 0004 7
8 Matrices green	
vks-oc Ø 2,2 mm	REF 430 0544 0
2 Matrix case	
vks-oc Ø 2,2 mm	REF 430 0696 0
2 Threaded sleeve	REF 580 0119 3
1 Screwdriver	
short, Hexagon	REF 330 0069 0
1 Inserting Instrument	
vks-oc Ø 2,2 mm	REF 430 0548 0
•	






365

### Snoring therapy / Biofunctional therapy

### Biofunctional therapy



Daily training with the vacuum activator helps to keep the lips closed and to create a self-contained state of rest within the mouth through swallowing.

This procedure is visualized by the pressure gauge of the device. The bio-functional therapy uses the vacuum activator's pressure gauge as bio-feedback signal for training the proper position of a self-contained state of rest.

The treatment steps are systematically controllable. Training procedure and training time are displayable via pressure monitoring. Therapeutic goal is a preferably permanent system stabilization of the self-contained state of rest.

This method was developed at the University of Göttingen/Germany under Prof. Dr. Dr. W. Engelke and is applied since 2003.

Literature: Engelke, W.: Systematische Ronchopathiebehandlung in der zahnärztlichen Praxis, Cuvillier Verlag, Göttingen.

The application is simple and may be integrated into the dental or orthodontic practice at any time

An "oral-friendly" design and high class materials grant a high wearing comfort and encourage your patients' active cooperation.

We would, of course, be more than happy to send you a detailed information leaflet.

#### Basic exercise with silencos:

Low pressure is created in the mouth when collecting saliva and swallowing. After swallowing, breathing takes place steadily through the nose and the low pressure created is maintained. Active muscle tension is not required. The membrane responds to the pressure in the mouth. During basic exercise, the membrane should be drawn into the funnel; this indicates the correct adoption of the closed position at rest. The exercises should be carried out several times a day for 15–30 minutes. The aim of the basic exercise is to learn the closed position at rest and maintain this in the long-term. The duration and frequency of the exercises is agreed between the patient and his/her therapist.

# Vacuum activator silencos for adult therapy



silencos provides an effective primary snoring therapy without obstruction.

Nasal respiration, the self-contained oral state of rest, and the velum are being trained through regular exercise and daily practice.

The mouth guard is made from soft, flexible silicon and can be worn with comfort.

#### Additional therapeutic possibilities are

- Practice of the self-contained state of rest
- Development of a nasal respiration habit
- Stabilization of tongue and velum
- Immobilization support for X-ray taking (OPG or CT)
- · Closed mouth training
- Mobilization of the mouth base after tumor surgical measures
- Exercises to restore the oral functions within neurologic rehabilitation

silencos incl. stora ge box 1 piece REF 580 0600 0 silencos Waiting room poster 2 pieces REF 0P0004GB silencos Patient flyer 20 pieces REF 000 274G B silencos Membrane 5 pieces REF 580 M600 0

Accessories:

silencos Funnel 2 pieces REF 580 T600 0 silencos Stopper 5 pieces REF 580 S600 0

### Vacuum activator silencos kids for early child therapy



Malocclusions and jaw malformations may be caused by incorrect swallowing patterns, pathologic mouth breathing or habits, such as thumbsucking. These should be recognized and treated at an early stage.

The vacuum activator offers all functions of an oral vestibule plate (OVP), yet reaches far beyond the possibilities of the OVP, because the exercises are controllable via membrane.

silencos kids was especially developed for the preschool child.

The mouth guard is made from soft, flexible silicon and can be worn with comfort.

#### Additional therapeutic possibilities are

- Controlled habit manipulation, such as thumbsucking, lip-biting
- Balance of forces within the orofacial system
- Adaption of a natural rest tongue position
- Myofunctional dysfunctions
- Treatment of habitually conditioned incorrect tongue positions
- Frontal open bite therapy
- · Adenotonsillectomy aftercare
- Habitual dysfunction therapy
- Play therapy

silencos kids incl. storage box 1 piece REF 580 0600 K silencos kids Waiting room poster 2 pieces REF OPO 005G B silencos kids Patient flyer 20 pieces REF 000 278G B silencos Membrane 5 pieces REF 580 M600 0

Accessories:

silencos Funnel 2 pieces REF 580 T600 0 silencos Stopper 5 pieces REF 580 S600 0

### silencos night device



silencos night device incl. storage box 1 piece REF 580 N600 0

# For use in the case of velar snoring or in the case of habitual breathing through the mouth.

The silencos night device is a silicon mouth guard. When combined with an individually manufactured maxilla splint, this closes the mouth from the outside and promotes breathing through the nose at night. This prevents flapping of the velum as a result of the air current produced when breathing through the mouth.

- Snoring sounds are clearly reduced or are prevented altogether. The device can only be used if the person using this can breathe through his/her nose adequately. The silencos night device is a favourable solution for snoring therapy and can be manufactured comparatively quickly. It works like a curtain in the vestibule of the mouth and promotes breathing through the nose as a result of this. The elasticity of this allows for a sufficient level of movement of the jaw and offers the patient a high level of comfort.

Due to the connection to the maxilla splint, the guard is always correctly positioned in the vestibule of the mouth, even if the mouth is open. Before using for the first time, it should be checked on the model or by the dentist on the patient whether there are any irregularities in the area between the lips, cheeks and dental arch, such as e.g. a deepset labial frenulum or an unfavourably shaped oral vestibule. In case the labial frenulum or oral mucous membranes are compressed by the device, the edges can be reduced appropriately with a scissors. We recommend using the silencos Night device in combination with the silencos Vacuum activator and biofunctional treatment. Regular therapeutic training stabilises the mouth and throat muscles and leads to habitual correctness of the tongue position.



Retaining elements made from wire are inserted in the molar region...



...on an individually manufactured maxilla splint (e.g. in deep-drawn procedures).

The retaining elements are individually bent from e.g. 0.9 or 1.0 mm hard wire and fixed to the maxilla splint using plastic.



Thea mouth guard is then attached to the openings of the silencos night device in the molar region. The soft, smooth material can be shortened using scissors or a silicon shaper if required.

bredent

# BETTER QUALITY OF LIFE – BACK TO A HEALTHY SLEEP



Snoring has become the most widespread disease. In 50% of bedrooms, healthy sleep is disturbed by extreme snoring – this can be as loud as 90 decibels, which corresponds to the noise of a passing lorry.

Two forms of treatment are available. The protrusion splint, which is worn in the case of slight or medium obstructive sleep apnoea syndrome, and biofunctional treatment with a vacuum activator, are used in the

case of primary snoring without obstruction.

Become an advisor for your practice and establish your laboratory's position in a new and attractive area of business. Comprehensive information material is available for laboratories, practices and patients.

Added value for your laboratory!

All individual components have been exclusively developed by experienced epitheticians. The technique that can be learned in courses opens new business fields to your laboratory.

Future developments which allow lasting bonding of metal/ silicone and silicone/silicone will become milestones in the field of epithetics.

The advantages of soft silicone and hard epithelial resin combined in a soft resin provide the basis for additional new developments.







Starter set for silicone epithetics	
Multisil-Epithetics Set	370
mpression material	
Multisil-Epithetics soft-form and hard-form	370
Modelling wax	
Modelling wax for epithetics	.371
Epithelial material	
Multisil-Epithetics transparent	.371
Multisil-Epithetics city / country / beach	.371
Coloring and characterization of epithetics	
Multisil stains	372
Multisil-Epithetics thickener	372
Multisil sealing agent	372
Multisil-Primer	372
Multisil intensive colors	373
Multisil fibers	373
Surface sealing agent	
Matt sealing agent for epithetics	374

### Starter set for silicone epithetics / Impression material

### Multisil-Epithetics Set



#### Starter set to produce a silicone epithesis.

- Content of the case prepared in cooperation with experienced epitheticians
- Contains all materials required to produce a silicone epithesis
- Robust aluminium case and clearly arranged comparments to find the necessary materials quickly

#### Multisil-Epithetics Set

- 1 aluminium case with foam lining
- 20 mixing cannulas, pink
- 1 dosing device
- 1 brush
- 1 mixing spatula for epithetics
- 1 sliding caliper
- 30 ml Isoplast ip
- 5 ml Multisil-Epithetics thickener
- 5 ml Multisil-Primer
- 10 different intensive colors, 5 g each
- 10 different fibers, 2.5 g each
- 10 vdifferent stains, 5 g each
- 2 x 50 ml Multisil-Epithetics soft-form
- 2 x 50 ml Multisil-Epithetics hard-form
- 1 x 50 ml Multisil-Epithetics city
- 1 x 50 ml Multisil-Epithetics country
- 1 x 50 ml Multisil-Epithetics beach
- 3 x 50 ml Multisil-Epithetics transparent

#### REF 540 0106 0





#### Accessories:

12 Mixing cannulas, pink

1 Dosing device

Brush, size A + holder

1 Brush, size C + holder1 Mixing spatula, epithetics

1 Sliding caliper

750 ml Isoplast ip

1 Mixing block 80 PE-foil cut-outs REF 330 0114 8 REF 320 0045 3 REF 320 0045 4

REF 540 0101 9

REF 320 0045 2

REF 320 0044 0

REF 330 0114 6

REF 320 0045 5

REF 320 0045 6

### Multisil-Epithetics soft-form and hard-form

#### Soft-form



Impression material for epithetics on 1:1 silicone basis in time- and material-saving double mixing cartridges.



- Due to the low hardness of 25 Shore A it is particularly suitable for undercut areas
- Extreme firmness ensures reliable impression-taking



Multisil soft-form to reproduce undercut areas.





Multisil-Epithetics hard-form 2 x 50 ml REF 540 0106 2

- The hardness of 45 Shore A provides stability for larger surfaces and for covering Multisil-Epithetics soft-form
- Extreme firmness simplifies reliable impressiontaking




Multisil hard-form – for covering and stabilizing Multisil soft-form.



### Modelling wax for epithetics mdwe



Skin-colored plate wax for epithetics.





- thickness of 2.8 mm adjusted softening temperature optimal hardness, improved plasticity well-balanced stickiness, which is perfectly matched with the epithetics, ensure quick and reliable modelling
- After heating, the modelling wax can be shaped for an extended time and so modelling of the epithesis is simplified.

### Multisil-Epithetics transparent



Multisil-Epithetics transparent 2 x 50 ml REF 540 0106 3

#### transparent


#### Transparent, soft epithelial material on 1:1 silicone basis.

- Convenient double mixing cartridge for consistent mixing quality
- Crystal-clear silicone for optimal, individual coloring with Multisil intensive colors
- Extended processing time of two hours at room temperature provides sufficient time for individualizations
- Simple polymerizing at 60° C does not require special equipment
- Final hardness of 35 Shore A and high tear strength ensure comfort of wear for the patient
- Fine flow behavior of the silicone allow most accurate reproduction of details of the model

### Multisil-Epithetics city / country / beach



Multisil-Epithetics city 50 ml REF 540 0106 4



Colored, soft epithelial material on 1:1 silicone basis.


• 3 different skin tones simplify coloring of the epithesis:

city = bright skin type country = normal skin type beach = dark skin type

- Mixing in Multisil-Epithetik transparent saves time during invididualizing
- Individual coloring with Multisil stains adds new possibilities of shade adaptation
- Final hardness of 35 Shore A and high tear strength ensure comfort of wear for the patient
- Fine flow behavior of the silicone allow most accurate reproduction of details of the model




Multisil-Epithetics country 50 ml REF 540 0106 5



Multisil-Epithetics

beach 50 ml **REF 540 0106 6** 





371

Fax (+49) 0 73 09 / 8 72-4 44

# Coloring and characterization of epithetics

### Multisil stains



### Stains for surface characterization. 10 different stains provide all options of optimal adaptation to the patient situation.



#### Multisil stains Glass jar, cont. 5 g



REF

# Multisil-Epithetics thickener



To change the viscosity of addition-linked silicones.

Multisil-Epithetics thickener transparent 5 ml REF 540 0106 8



Silicone without Multisil-Epithetics thickener.



Silicone with Multisil-Epithetics thickener renders the silicone firm and simplifies layering of the epithesis.

# Multisil-sealing agent



Sealing varnish for silicone surfaces.

Multisil-sealing agent transparent 10 ml REF 520 0100 5



Sealing the base of the epithesis avoids the accumulation of dirt and secretion and thus simplifies cleaning.

### Multisil-Primer



Multisil-Primer **Bonding agent** 

Multisil-Primer REF 520 0100 4



Optimal bonding of silicone and acrylic resin.



# Coloring and characterization of epithetics

### Multisil intensive colors



#### Silicone colors for coloring addition-linked silicones.

- 10 different intensive colors provide comprehensive possibilities of color characterization
- High color stability avoids discoloration of the epithesis



### Multisil fibers



#### Viscose fibers for characterizing epithetics.

- 10 different intensive colors provide comprehensive possibilities of color adaptation and characterization
- Special, thin fibers allow perfect reproduction of blood vessels, downy hair, etc.




Viscose fibers with different colors especially matched with the requirements of epithetics.

### Multisil fibers Plastic jar, cont. 2.5 g

brown

, , ,		
	color	REF
	silver	530 0060 0
	white	530 0060 1
	beige	530 0060 2
	signal red	530 0060 3
	ruby-colored	530 0060 4
	purple	530 0060 5
	bordeaux	530 0060 6
	blue	530 0060 7
	mocha	530 0060 8
	ochre	530 0060 9

540 0107 9



373

# Surface sealing agent

### Matt sealing agent for epithetics



Matt sealing agent for epithetics Primer 10 ml REF 540 0109 1

Matt sealing agent for epithetics Coating

20 ml REF 540 0109 2

Matt sealing agent for epithetics Matting powder REF 540 0109 3

Assortment

Matting powder

REF 540 0109 4

1 Matt sealing agent for epithetics, Primer

1 Matt sealing agent for epithetics

1 Matt sealing agent for epithetics, Coating

3 pieces





A thin coat of primer is applied on the surface of the epithesis.



Multisil stains allow further characterization of the epithesis.



Uniform, thin layers of



coating are dabbed on using a brush.



Allow to dry for approximately two minutes at 65° C (e.g. hot-air blower).



The matting powder is spread on after 15 minutes.



Place the epithesis in hot water for two minutes.

### Matt surface sealing for addition-curing

- Creates a matt surface on silicone epithetics and thus ensures a natural appearance of the epithesis
- Prevents the stain coat from coming off the sur face and thus offers extended comfort of wearing
- Contains a UV protective varnish for the silicone colors and stains and provides lasting protection against fading
- Simple handling ensures reliable and permanent sealing of the surface of the epithesis

# Devices – the best quality for demanding jobs



The devices produced in Germany offer the highest quality and reduce the amount of work required, which saves time. Specially developed by dental technicians for the fields of work relevant to dental technicians.



ecovac vacuum mixing system	376
Master pin drill unit mpb 1	377
KoEx Measuring Device	378
ore.Lux Power Unit	379
Polylux 2	380
Protective chamber	
forl 2 press and BioHPP	381
hermopress 400	382
CPS Cordless Prosthodontic Screwdriver	384
Naxpool duo	385
Wax knife bwm 3	386
Milling unit BF 2	387
Milling base	388
Model support BF 2	
Fransfer device	388
Brenometer surveying system	388
airaqua turbine	389
Polierjet	390

### ecovac vacuum mixing system



#### ecovac

#### Precision-fit restorations obtained through optimal use of material properties.

The user-friendly and compact design simplifies work and reduces sources of errors. A powerful and maintenance-free vacuum pump, adjustable in two different levels (15 mbars, 200 mbars), ensures bubble-free mixing of materials and results in a perfect casting surface. Stirring time and speed can be adjusted continuously to allow correct processing of different materials.

ecovac (230 V)

REF 140 0093 0

(Wall mounting, without mixing cup and base)

- 1 mains cable
- 1 spare filter
- 1 drilling template for wall mounting
- 4 screws and plugs for wall mounting

Accessories:

base ecovac, 1 piece REF 210 0045 0



### ecovac mixing spiral

The mixing spiral takes up the components to be mixed from all areas of the mixing cup and stirrs them horizontally and vertically. No unmixed materials will remain on the bottom of the mixing cup, which may cause different expansion of the material later on.

All features and components listed provide increased reliability, lead to improved fit when preparing dental restorations and avoid time-consuming

Mixing spiral,	50 ccm	REF 140 OR94 5
Mixing spiral,	250 ccm	REF 140 0R94 0
Mixing spiral,	750 ccm	REF 140 0R94 2
Mixing spiral,	1000 ccm	REF 140 0R94 3



### ecovac mixing cups

The smooth inner surface of the stainless steel mixing cup prevents any material or liquid residues from adhering to or depositing in scratches or undercuts. The conical shape ensures that material which has been taken up will flow back to the center of the mixing cup. Accordingly, the mixing ratio is retained exactly and better results can be achieved with minimal effort.

Mixing cup,	50 ccm	REF 140 0B94 5
Mixing cup,	250 ccm	REF 140 0B94 0
Mixing cup,	750 ccm	REF 140 0B94 2
Mixing cup,	1000 ccm	REF 140 0B94 3



Mixing cup, D (for the use in the Degussa mixing unit),

REF 140 0B94 4

### Master pin drill unit mpb 1

The powerful, high quality and maintenance-free motor features high true running accuracy.

Accordingly, the precision of the drilled hole and the accuracy of the models are increased. Working is simplified thanks to the easy-to-operate lifting mechanism

Master pin drill unit mpb 1 REF 140 0092 0


(without 15° base)
1 spare fuse

1 flat wrench

Master-Pin Diatit
 tungsten carbide bur

1 plaster collecting tray

r 1 plug axle 1 power cord




#### Accessories:



Adapter base 15° inclination / precious wood REF 210 0044 0

15° inclination adapter / precious wood








Tungsten carbide drill Special drill for Master-Pin Radix-K Ø 2,0 mm 3 mm shaft REF 360 0123 3



# Master-Pin Diatit tungsten carbide step drill special/yellow RFF 360 0119 3





# Master-Pin Diatit tungsten carbide step drill special/red REF 360 0119 4

If the drilled hole is too large to receive the Master-Pin, the special/red Master-Pin Diatit tungsten carbide drill can be used to prepare a smaller drillhole. The diameter of this drill is 0.01 mm smaller than the one of the standard/green Master-Pin Diatit tungsten carbide drill.

### KoEx Measuring Device

Why are there fitting differences between the model and the actual situation in the mouth? Thanks to the KoEx measuring device, material measurements of contraction and expansion are now possible for the first time.



Facilitating contraction and expanding measurements for the first time

Why do discrepancies in fit exist between the cast and the intraoral situation?

**KoEx Measuring Device** 1 piece including 2 contraction inserts REF 110 0148 0

#### Impression Materials, Contraction

Studies have indicated that impression materials differ greatly in their contraction (shrinkage) behavior, reproducing the oral situation inaccurately. The brecision impression material provides




stable values after two hours, permitting further processing to be performed rapidly.

#### Dental stone, Expansion

Check the expansion values for your dental stone materials and compare them to those of bredent's Exakto-Rock S. Exakto-Rock S expands by a maximum of 0.06 % after two hours; after 48 hours, the expansion is still less than 0.08%.



#### Comparison of contraction values for A silicones and polyether materials 0.140 0.120 0.100 0.080 0.000 Weau + 0.020 0.000 0,5 h p=1.00 p<0.001 p<0.001




### Silicone duplicating materials, Contraction Contraction measurements of different silicone

duplicating materials have shown substantial differences between these. Exaktosil N15 was stable at 30 minutes, at 0.02 %. The values for other



duplicating silicones changed after 24 hours, adversely affecting the fit of the restoration.



#### Investment compound, Expansion

Investment compounds that can be controlled exactly and individually are a prerequisite for nonprecious alloy precision one-piece attachment



casting as well as for K+B plastic injection molding using thermo press 400.

#### bre.Lux



#### Polymerization Multi Talent for Practice and Laboratory

To date, several devices were necessary for processing facing and dental materials. With bre.Lux, the concept of processing all current materials with one single device becomes reality for the first time.

#### Assortment

bre.Lux Power Unit Product package consists of

1 bre.Lux UpDown 1 flexible hose

1 polymerization device

1 bre.Lux LED N manual lamp with spiral cable

1 mains cable REF 140 0097 0

- Fixation / Hardening / Intermediate polymerization and final thorough hardening directly in the workplace and with one single device
- 370 500 nm covers the needed wavelength range, for the manual lamp as well
- Performance delay and reduction minimize the properties and results of dental materials
- Start-up delay and polymerization time can be easily combined
- Large volume for 2 models, optimally and uniformly



#### Keyboard layout of the light polymerization device

- (1) Display
- (2) 20 seconds with 50 % power
- (3) 40 seconds with 50 % power
- (4) Continuous operation with 70 % power
- (5) 90 seconds full power
- (6) 120 seconds full power
- (7) 180 seconds full power
- (8) 360 seconds full power
- (9) Continuous operation at full power

#### Consistency

The membrane keyboard already comes with several parameters. On the left side: Reduced power for the start-up phase with 20 and 40 seconds and for continuous operation at full power as well. The combination (such as start with reduced power) can be combined with the programmed running times of 60 to 360 seconds according to demand. Example: 40 seconds with reduced power combined with a 180-second total running time means that the total running time is lowered from 180 to about 40 seconds, and the running time would still be 140 seconds under full power.

Please take note of the application recommendations for bredent products.

#### Keyboard layout for the manual lamp



bredent







- (1) Continuous operation, start/stop
- (2) 15 seconds
- 3) 30 seconds



The bre.Lux Power Unit consists of one LED light polymerization device with 21 power LEDs in 3 different capacities, from 370 nm to 500 nm. The LEDs have a useful life of 20,000 hours. The bre.Lux LED N manual lamp (with spiral cable) features a capacity ranging from 370 nm to 500 nm. The flexible hose – with receptacle ring for the manual lamp - serves as third hand and allows two-handed work.

### **Devices**

### Polylux 2

The powerful and universally applicable light-hardening device for materials with a wave length range of 350 to 500 Nm.



Polylux 2 with drawer, 230 V Polylux 2 with drawer, 115 V REF 140 0099 0 REF 140 0099 1 Two different and energy-saving special lamps guarantee excellent polymerization. Moreover, perfect distribution of light ensures careful and appropriate curing of the materials.

- Convenient and simple operation thanks to clearly arranged buttons
- Fully mirrored polymerization area for selective illumination with compact drawer
- Device can be accessed from three sides to enable polymerization of large objects
- lower energy consumption for a cost-effective service life

#### Technical data

Number of light sources Wavelength range Triple timer function Mains voltage Dimensions of unit Dimensions drawer Weight 2 fluorescent lamps, 9 watts each 350 – 500 Nm 180 sec/360 sec/continuous operation 115 V, 50 Hz / 230 V, 50 Hz approx. 250 x 120 x 90 mm approx. 140 x 110 x 55 mm approx. 1500 g

### Protective chamber



The protective chamber avoids inhaling of dust, protects your eyes and, consequently, protects your health. Available with or without extraction nozzle. The extraction nozzle can be directly connected with the extraction system.

Protective chamber with extraction nozzle Dimensions: approx. w 410 x d 350 x h 260 mm  $\emptyset$  35 mm

REF 220 0010 0

**Protective chamber without extraction nozzle** Dimensions: approx. w 410 x d 350 x h 260 mm

REF 220 0011 0

### for 2 press and BioHPP



# for 2 press vacuum-press system for the manufacture of metal-free, biocompatible dental prosthesis.

With the for 2 press vacuum-press system, the thermoplastic high-performance polymer BioHPP is processed in a mould made from a phosphate-bonded investment material. The end result is a metalfree, white framework construction, which can be veneered perfectly with conventional veneering

composites in a customised, aesthetically pleasing manner. The dental prosthesis made from BioHPP is certified for permanent applications.

The melting process of the BioHPP high-performance polymer is carried out in the pre-heating furnace already available to you. The subsequent pressing process is fully automated and takes place in a vacuum.

for 2 press device 1 piece REF 140 0060 0

#### for 2 press Basic Set


1 x foil 2 press mould (consisting of a mould plate 3/16 mm and silicone ring, 35 x 210 g
Brevest foil 2 press EBM incl.
2 litre Bresol foil 2 press liquid,
25 x foil 2 press filler 16 mm (one-way extrusion die for moulding the material in the mould
20 g BioHPP, 1x processing instructions
REF 140 0060 1



After the investment material mould has been pre-heated in the pre-heating furnace and the modelling wax or plastic has melted, the BioHPP is melted at 400° C, also in your pre-heating furnace. Immediately afterwards, the mould with the melted BioHPP and with the fitted one-way extrusion die is moved into the *for* 2 press vacuum-press device. The raising of the lift triggers the automatic pressing process in a vacuum atmosphere.



After the 3-minute pressing process, which includes the application of a vacuum, the cooling process begins, while maintaining the pressing pressure. This is the only way to guarantee that the material properties of the BioHPP are exploited fully.



After 35 minutes in total, the entire pressing process is completed and the user is informed of this by an acoustic signal and an optical LED display. Now the BioHPP framework structure can immediately be divested and processed further.



Sample application using a 3-part bridge made of BioHPP veneered with visio.lign veneers.

Image: Dental technical certified company Harald Schwindt.

#### Technical data for 2 press

Power supply
Power consumption
Vacuum performance
Venturi nozzle
Weight
Size (W x H x D)
Protection class
Sound level
Input air pressure

Hardware protection

90 - 250 volts, 50 - 60 Hz 15 watts

approx. 760 mbar 13 kg 250 x 600 x 290 mm IP 34 < 70 dB min. 4.5 to 6 bar max. T 2.5 A

#### Accessories:

High Performance Polymer
BioHPP, 20 g REF 540F2PB2
BioHPP, 100 g REF 540F2PB3

# Mould system for 2 press mould set,

16 mm, third mould plate 16 mm and third silicon ring REF 360F2P16 for 2 press mould set,

20 mm, ninth mould plate 20 mm and ninth silicon ring REF 360F2P20

for 2 press silicone ring, Gr. 3REF 360F2PR3for 2 press mould plate, Gr. 3REF 360F2PR9for 2 press mould plate, Gr. 3REF 360F2PR9for 2 press mould plate, Gr. 3REF 360F2PT1for 2 press mould plate, Gr. 9REF 360F2PT2

#### One-way extrusion die

for 2 press filler, 16 mm, 25 pieceREF 570F2P16for 2 press filler, 20 mm, 25 pieceREF 570F2P20

#### Investment material

 Brevest for 2 press, 35 x 210g, 1000 ml Bresol
 REF 570F2PV1

 Bresol for 2 press liquid 1000 ml
 REF 520F2PL1

### thermopress 400



Injection moulding device for the processing of thermoplastic plastics with a melting temperature of up to  $400^\circ$  C

- No additional equipment such as CO<sub>2</sub> cylinder or a compressed-air connection required. This way, additional costs and time are saved. The same level of quality is maintained by avoiding a drop in pressure.
- Higher user comfort during handling of the device.
- Time-saving function thanks to simultaneous operation of both heating chambers.
- The injection process can only be performed with the lid closed; additional safety is provided.
   Convenient removal of the flask thanks to automatic cartridge ejection if the bracket is unlocked.

thermopress 400

1 unit with power cord
2 allen keys

1 cleaning brush 1 special tool REF 110 0040 0



The melting process is adapted to the requirements of the plastic thanks to special high-performance heating elements. This ensures that the best possible mechanical material properties are achieved.





The compression forces of up to 2 kN occurring on the special cuvette are absorbed by a special conduction and fastening system into the device. This guarantees that the melted thermoplastic plastic is injected into the cuvette hollow mould without any air pockets and without deformation.



### thermopress 400

#### Technical data thermopress 400

Width 650 mm

Height 250 mm

Depth 300 mm

Weight 40 kg

Voltage 220 - 230 V

Power 0.5 - 1.6 kW

max. 2.2 kW

#### thermopress 400 accessories assortment (device)

7 pieces

#### REF 110 0040 1

View above, accessories thermopress 400: articles marked with \* are included in thermopress 400 accessories assortment!

#### Patient information brochure thermopress sensitive REF 000 415G B

thermopress 400 patient pass REF 000 628G B

thermopress 400 information brochure REF 992 945G B

thermopress 400 patient information REF 000 135G B

#### Training and further education programme

2-day course (Senden) REF 950 0020 0

#### Accessories thermopress 400:

REF 140 0090 4 1 press-out device and punch\* 1 pair of cartridge pliers\* REF 140 0090 6 1 flask hook with hex\* REF 140 0091 2 1 flask N, small\* (I 122 cm, w 102 cm, h 72 cm) REF 140 0090 3 1 flask N, large (I 140 cm, w 102 cm, h 72 cm) REF 140 0090 5 1 cleaning brush REF 110 0040 2 1 thermopaste 400 special paste, 50 g\* REF 540 0105 1 Expando-Rock-Set 5 kg expansion plaster, 500 ml Expandosol REF 570 OERS 5

#### 3 1 1 1

5 x 30 g bre.crystal HP pink 1 2 x 20 g bre.flex pink 3 x 24 g bre.crystal HP clear 2 x 16 g bre.dentan HP A 5 x 24 g bre.flex 2nd Edition 3 x 16 g Bio Dentaplast A2 pink veined PEE 140 0000 3

pink veined REF 140 0090 2

thermopress introductory assortment (20 cartridges)

#### thermopress 400 range S1 (20 cartridges)

2 x 30 g bre.crystal HP pink 1 1 x 20 q Bio-Dentaplast A3 1 x 20 g 2 x 24 g bre.crystal HP pink Bio-Dentaplast A3 5 x 24 g bre.flex 2nd Edition 1 x 250 ml Acryl Sep 1 x 20 ml light-hardening Stumplack pink veined bre.crystal HP clear 1 x 24 g transparent bre.dentan HP A 1 x 50 g thermopaste 400 5 x 16 g REF 540 S000 1

### thermopress 400 range S2 (20 cartridges)

5 x 16 g Bio Dentaplast A2 5 x 16 g Bio Dentaplast B3
5 x 16 g Bio Dentaplast A3 1 x 250 ml Acrylic Sep
6 Acrylic/plaster separating liquid
7 REF 540 S000 2

### CPS Cordless Prosthodontic Screwdriver



#### How much time do you spent loosening and fixing implant screws?

You can save app. 50 % of the time using the CPS.

#### Are you sure that your screws are always tight correctly?

With CPS there will be no screw loosenings due to wrong torque.

### Can you reach your palatinal screws easily?

With CPS you have a perfect cordless access to all sites, because it is cordless.

#### **CPS Cordless Prosthodontic Screwdriver**

- Cordless motor
- Battery runtime for 40 screws
- Charging point
- Handpiece 80:1
- Screw drivers:

REF 580 CPS4 0

#### Hexagon
























# Waxpool duo



#### Waxpool duo unit REF 110 0150 0

#### Assortment

- 4 parts
- 1 Waxpool duo unit
- 1 Waxpool duo handpiece
- 2 Waxpool duo contouring blades at your choice REF 110 0152 0

#### Accessories:

Waxpool duo handpiece REF 110 0151 0 Rest REF 140 0096 5

# Wax dipping unit and wax knife all in one – digital control for added comfort • Stable and easy to clean plastic housing

- Exchangeable plastic lids
- Clean design
- Clear design
- °C or °F can be selected

#### Wax dipping unit

- Precise temperature control of the dipping wax for increased quality
- High-performance heating elements reduce the time for heating the wax
- Uniform wax copings thanks to constant temperature control
- Special, lowered safety dipping wax to avoid burning of fingers
- Melting temperature up to 120° C

#### Wax knife

- A separate wax knife can be connected
- A single unit at the working place
- Non-tiring working thanks to ergonomic design of the handle
- Special insulating elements reduce heating up of the handle
- Simple exchange of blades
- Boost key for quick heating up to the end temperature
- Maximum temperature of 240° C



Contouring blade size 1 REF 320 WP4G 1



Contouring blade size 3 REF 320 WP4G 3



Contouring blade size 5 REF 320 WP4G 5



Contouring blade standard REF 320 WP47 2

### Wax knife bwm 3



Control unit bwm 3 with handpiece and contouring blade size 5 REF 140 0096 3

Control unit bwm 3 REF 140 0096 0

Handpiece bwm 3 REF 140 0096 2



Rest bwm 3 REF 210 0045 1




Footswitch bwm 3 REF 140 0096 1

#### Accessories:

Foam rubber grip lining 4 nieces REF 140 0096 4

Electric wax knife featuring integrated advanced technology and high quality. The ergonomic handpiece allows to take up wax quickly and ensures comfortable working.

- Ergonomically designed handpiece
- Quick heating up with the Rapid-Speed footswitch
- Adjustable temperature control
- Simple and fast exchange of the contouring blades



Contouring blade bwm 3 size 1

Contouring blade bwm 3 size 3

Contouring blade bwm 3 size 5

Contouring blade bwm 3 Standard

REF 320 004G 1

REF 320 004G 3

REF 320 004G 5

REF 320 0047 2



Comfortable and quick removal of the contouring blades.



Device for firm, reliable hold of the handpiece at the unit.



Mobile rest for safe depositing of the handpiece.



Blade shapes proven over numerous years allow individual application.



Integration into the grip for quick and simple exchange of the contouring instruments without the risk of injuries.



The special instrument grip avoids twisting of the contouring tip whilst working.



The contouring tips are stored on the control unit in a safe and clearly arranged manner.



If the wax knife is not needed, it can be placed on the rest in the direct reach of the technician.



The footswitch allows to quickly reach a higher temperature than the one that has been set. Activation of the footswitch is indicated by the control lamp.




Logical and clearly arranged control unit for stress-free and safe working.



Handpiece with flexible, stable cable for simple working.



High-tech dental equipment featuring highly useful function and design - for comfortable and simple working.



Milling unit BF 2 including 1 Model support BF 2 1 piece

REF 140 0098 0

#### Technical Data

230 Volt / 50/60 Hz Power supply 80 Watt Power rating Speed 0 - 30,000 U/min. Ø 2.35 mm Chuck Fuse thermal overload protection Torque 2.6 Ncm Weight 18 kg Width/Depth/Height 250 x 370 x 510 mm

#### Accessories

Chuck 2.35 mm	REF 730 0016 9
Chuck 3 mm	REF 730 0015 3
Tap handwheel	REF 330 0115 4
Model support BF 2	REF 730 0017 0
Milling base	REF 140 0089 3
Adapter airaqua turbine	
16 mm	REF 730 0018 4
18 mm (for BF1)	REF 730 0018 3
28,5 mm	REF 730 0018 5
Transfer device	
3 mm shaft	REF 360 0116 3
2.5 mm shaft	REF 360 0126 5

### Milling base



Milling base with integrated thread for fixation on the milling base of the BF 1 unit. Additionally, plaster is removed completely and without damaging the metal plate by slightly turning the locking bolt.

Milling base 1 piece REF 140 0089 3



### Model support BF 2



The model support can be used for all milling units including units with magnetic circuit. Turning by 90° permits do perform lateral drilling of bars without removing the model.

Model support BF 2 1 piece REF 730 0017 0



### Transfer device



Permits correct transfer of the position of the model to the milling base. Up to 8 units can be transferred at the same time.

Transfer device 3 mm Schaft REF 360 0116 3 2,35 mm Schaft REF 360 0126 5

# Brenometer surveying system



Brenometer surveying system

Four different surveying plates according to Ney allow accurate positioning of the clasp profiles whilst ensuring correct depth of undercuts. A locating pin and a red marker with a holder ensure correct surveying.



Marking the clasps and surveying with a single unit – this is how time and money can be saved.

#### Assortment

- 1 Brenometer marker holder
- 1 Brenometer locating pin
- 1 Brenometer plate 0.25
- 1 Brenometer plate 0.35
- 1 Brenometer plate 0.50
- 1 Brenometer plate 0.75

REF 310 0000 2

#### Refill packages:

Brenometer marker holder REF 310 0000 4 REF 310 0000 3 Brenometer plate 0.25 REF 310 0002 5 Brenometer plate 0.35 REF 310 0003 5 Brenometer plate 0.50 REF 310 0005 0 Brenometer plate 0.75 REF 310 0007 5

### airaqua turbine



#### Technical data:

Speed 300,000 rpm
Energy supply compressed air
Operating pressure 2.8 - 3.2 bars
Air consumption 40 l/min.
Water reservoir 350 ml
Collet 1.6 mm
Lubrication manual

Width approx. 190 mm Height approx. 190 mm Depth approx. 125 mm airaqua turbine is a handy, compact unit with a light-weight handpiece for precise processing of hard materials such as high-performance ceramics (sintered zirconium oxide), press and metal ceramics. The airaqua turbine features a spraying device to spray an air/water mixture onto the processing area. Water cooling avoids overheating of the material.

The formation of microcracks is reduced considerably so that safe processing of materials is ensured. The water spray traps the grinding particles, protects the grinding tools and thus extends their service life. As an option, an adapter is available or using the turbine handpiece in milling units.

The spray can be switched on and off quickly with the switch on the handpiece. Fine adjustment is achieved with the two regulators in the table unit. A very small rotor allows extremely comfortable working and perfect view on the workpiece. The lubricant is directly fed into the bearings. The handpiece features a special adapter (midwest) and thus can be used with angle handpieces, turbines and air motors with the same standard.

airaqua turbine REF 110 0146 0 Scope of delivery:

Table unit with filter, controller, manometer, water reservoir and regulators, footswitch, handpiece with rotor, special oil 30 ml and

#### Accessories:



16 mm REF 730 0018 4 18 mm (für BF1) REF 730 0018 3 28,5 mm REF 730 0018 5

Adapter for airaqua turbine





Refill package airaqua turbine oil 30 ml REF 520 0033 5

### Polierjet

Further development of an industrial polishing technique rationalizes dental surface processing.



Quadro-Finish

Quadro-Finish polishing unit incl. starter kit polishing material 4 polishing drums REF 130 0046 0 Technical data
Height 860 mm
Width 830 mm
Depth 600 mm
Motor power 0.75 KW
Current consumption 2.7 A /230 V
Weight 152 kg





Duo-Finish polishing unit incl. starter kit polishing material 2 polishing drums REF 130 0045 0 Technical data
Height 670 mm
Depth 755 mm
Width 600 mm
Motor power 0.75 KW
Current consumption 2.7 A /230 V
Weight 120 kg

Accessories:

Foot rack REF 730 0016 8 Polishing drum REF 730 0016 7 Leasing at request



Polishing porcelain 3/3, 8000 g REF 730 0015 7 Polishing porcelain 6/12, 7200 g REF 730 0015 8 Polishing porcelain 9/9, 7800 g REF 730 0015 9 Polishing porcelain coarse, 5500 g REF 730 0016 2 2. Main polishing



Polishing porcelain spheres, 9800 g REF 730 0016 0 Polishing porcelain pins, 8800 g REF 730 0016 1 Polishing powder, fine, 3500 g REF 730 0016 3 Compared to conventional polishing, the biocompatibility of chrome cobalt castings is increased due to a compacted surface:

• no mechanical irritation of the mucosa by

- no mechanical irritation of the mucosa by sharp edges
- no mechanical irritation of the mucosa due to pores
- constant polishing performance
- · no deformed chrome cobalt objects due to manual polishing
- improves and standardizes the precision of fit of chrome cobalt work
- time-consuming, unpleasant polishing with rubber polishers which is also injurious to health is no longer required
- less dirt and more pleasant working conditions in the laboratory
- surface hardening of non-precious alloys clasps become clearly more elastic
- perfect polishing of the inside of the clasp

Due to the rolling effect of the polishing materials the structure of chrome cobalt claps is enhanced. A hard shell is formed on the outside and a soft core in the inside so that the clasps become more elastic and flexible similar to a cornstalk.




Chrome cobalt supply, crowns and bridges made of non-precious alloys and titanium are polished to pre-high luster in three working steps. Compared to conventional polishing, approx. 20 minutes can be saved for each chrome cobalt object.

3. High luster polishing



Polishing granulate, 4000 g REF 730 0016 4 Polishing cream, 290 g REF 730 0016 5 Polishing stainless steel pins, 2500 g REF 730 0016 6 Innovative instruments simplify procedures during day-to-day work. In accordance with this principle, bredent develops sophisticated products that generate time-saving and cost-reducing workflows.



Abdruck-Cut	392
Cervical disc	392
Thermo-syringe	392
Piezo-Blitz pb1	393
Plaster knife	393
Ergonom wax knife	394
Quick-Mandrell-System	394
Quick Change	395
Spot Clip	396
Mamelon cutter	396
Quicktool	397
Ceramix	397
MagicContrast	
MagicBrush	
KoliBrush	400
Unique Brush	
Blue-Clip	402
Ball-Clip	
Repositioning tweezers	404
Activating pliers	405
Novo-Grip	405
Pollygrip	406
Transfuser	406
Wax adapter	
Universal screwdriver set	
Universal screwdriver set for contra-angles	
Screwdriver	410

### Abdruck-Cut



Undercuts can be removed from the impression material with ease and in a targeted manner, using the loop blade that is as sharp as a scalpel.

Abdruck-Cut 1 piece REF 360 0114 0



The scalpel-sharp loop blade allows cutting even in areas difficult to access.

#### Accessories:



Loop knife 1 piece REF 360 0115 0

### Cervical disc



Increases precision and reduces working time when making wax or plastic patterns.

Cervical disc REF 320 0091 0



This precision steel cutting disc is 0.1 mm thick, 3.0 mm in diameter and can be guided exactly when cutting.



Shows a comparison of the cuts made in 0.5 mm thick wax using a scalpel (left) and cervical disc (right).



Shows an extremely precise pattern, produced using our dipping wax without the need to rewax the cervical margin.

### Thermo-syringe



Fixing and glueing, that can be dissolved quickly without any residues, for any type of model situation.

The adhesive resin wax can be moulded by heating and easily placed on the models.

Thermo-syringe REF 110 0121 1



After heating, the adhesive resin wax is directly applied onto the glueing point using the thermosyringe. Firm bonding is ensured.



The adhesive resin wax can be applied onto any type of material. Afterwards it can be removed from the objects without leaving any residues.

#### Accessories:



Adhesive resin wax Pack cont. 250 g Bucket cont. 1000 g

REF 510 0070 1 REF 510 0070 0

# Piezo-Blitz pb 1




Piezo-electric ignitor for all gas burner types. Suitable for all burner types (even old ones)!



Piezo-Blitz pb1 REF 360 0126 6









The main and the economy flame can be ignited by turning the ignition electrode.

### Plaster knife




Multi-purpose knife with ergonomically shaped plastic handle for optimum transfer of force, simplifies your daily work.

- Long blade made of stainless hardened steel.
- Dimensionally stable, easy-to-clean hard plastic handle. Ergonomic shape for right and left hand
- Multi-purpose element for easy removal of im pression tray. Features impact surface with opposing chisel.



Plaster knife REF 310 0011 4







When opening flasks, the lateral chisel ensures improved transmission of force thanks to the high leverage effect of the knife handle.

Plaster edges can be

perfectly trimmed with

the permanently sharp

and stable blade.



A separate impact surface has been added opposite the chisel to protect the back and the blade of the knife.



### Ergonom wax knife



#### Modelling knife for dental prosthetics.

Various instruments all in one - hence instruments do not need to be changed any longer so that faster and more efficient processing of the wax model is

Ergonomic design of handle - suitable for right- and left-hand users.

Ergonom Wax knife REF 310 0001 3



Special, ground edge of the knife tip for simple and fast modelling of age-specific papillae shapes.



The spoon which features a ground edge is perfectly suitable for modelling the alveolar area.



The deep spoon perfectly allows to apply large wax quantities within a very



Well-aimed, fast application of wax reduces the time for remodelling in the interdental area.



Transitions towards the functional margin and the functional margin itself can be prepared swiftly and neatly thanks to the curved design of the spoon element.





With the Ergonom wax knife, wax models can be easily and quickly shaped so that a natural appearance is obtained

### Quick-Mandrel-System



Save 40 seconds every time you change a separating disc, wheel or polishing disc. Every second counts!

The magnetic screwdriver holds the hexagonal nut.

Hexagonal nut with large support, made of magnetizable, hardened steel.



Stainless steel Quick-Gradle 1 piece REF 360 0115 5

up to 1 mm thick discs 10 pieces REF 360 0115 4

size 1



1-3 mm thick discs 10 pieces REF 360 0115 3



size 3 3-5 mm thick discs 10 pieces REF 360 0115 2

#### Save 40 seconds

Assortment

1 Quick-Gradle 2 Quick-Mandrels size 1

2 Quick-Mandrels size 2 2 Quick-Mandrels

size 3

REF 360 0115 6





Changing the disc with pliers and an instrument wastes a lot of time!



in future

The Quick-Mandrel-System - Unbeatably fast and pratical.

### Be faster and get ahead



Place the magnetic screwdriver on the Quick-Mandrel and turn it to loosen the hexagonal nut. The hexagonal nut is retained in the magnetic screwdriver.



Tighten the magnezitable hexagonal nut to position the new separating disc on the Quick-Mandrel, securely and in the centre.



## **Quick Change**



The combination of design, function and systematics.

#### Quick change system for instruments for ceramic, model fabrication and prosthetics

- Design carbon handle esthetic and
- Stainless blade holder with magnetic receptacle for perfect fixation of all blades
- All metal components and magnets are corrosion-resistant
- Temperature resistance of instruments inside the handle: 80°C
- Single hand use with quick change system
- Individual indications for ceramic, model fabrication and prosthetic
- Easy control thanks to reduced range of instru-

- Clever system allows to find the suitable instruments quickly (protection of registered design!)
- Troublesome screwing with keys is avoided
- Blades can be adjusted at any position - familiar working position is retained
- High safety thanks to immediate release of the heated blades
- Proper storage of sensitive ceramic blades
- Ceramic brushes can be perfectly stored in a hanging position
- Ceramic blades with high surface quality for outstanding gliding properties



Carbon handle L 101 mm, Ø 8 mm REF 310 0103 1



w 102 x d 100 x h 75 mm Weight approx.. 575 g REF 310 0103 0

Overview of instruments Dimensions in mm



REF 310 0105 6



Fissure tool REF 310 0103 4



Olive REF 310 0105 7



Croco, smooth REF 310 0103 2



Croco, serrated REF 310 0103 3





REF 310 0105 3



REF 310 0105 4



MagicContrast REF 310 0105 5



size 6 REF 310 0104 4



KoliBrush size 8 REF 310 0104 5



size 8B REF 310 0104 6



to Zahle REF 310 0104 0



Probe 0.8 REF 310 0104 1



REF 310 0104 2



Blade 0308 Fig. 3 REF 310 0103 7



Blade 0408 Fig. 4 REF 310 0103 9

### **Instruments**

# Spot Clip



Artery clip with spot-shaped holding area simplifies veneering - in two variations



Every ceramic specialist is familiar with the problems of metal-ceramic crowns without a metal margin: the holding spot of the artery clip is not sufficiently covered by base material (opaque). The problem can be solved with Spot Clip.



This clip covers only a tiny spot of the surface to be veneered. The base material can be easily applied around the holding spot of the clip.



**Spot Clip** 1 piece REF 310 0000 5





After removing the Spot Clip, the aqueous base material fills the holding spot of the clip. This way a uniform smooth layer of base material is obtained.

Large mamelon cutter

for upper incisors.



Spot Clip simplifies the application of stains and glaze material. No smearing of stains, no subsequent application of stains in the area of the holding spot of the clip ist required.

### Mamelon cutter



Simplifies the incisal design of ceramic crowns.

Mamelon cutter 1 piece REF 310 0000 1







The incisal edge is cut back using the mamelon cutter.

Small mamelon cutter

for lower incisors.



The contoured dentine core after firing: A base for incisal design options is obtained in a fast, safe and easy way.



The individual shades can be applied onto the dentine core – regardless whether firing has been carried out or not.



The incisal edges of the finished crowns exhibit a vivid play of colours.



## Quicktool



Ceramic structures are held safely without any pressure thanks to the three galvano plated diamond tips and the locking mechanism even galvano formed copings.



The gripping force can be adjusted to the crown size without deforming the crown.



In case of limited space, a diamond tip can be removed - ideal for lower anterior crowns.



achieved also when condensing.

uniform hold of the crown. Accordingly, safe hold is

Quicktool REF 310 0102 0



Accessories:

3 diamond tips REF 310 0102 1




The integrated condenser condenses the ceramic material in next to no time.

### Ceramix



Reproducing individual shade mixtures quickly and easily. Ceramic material is saved thanks to controlled dosing.

Ceramix REF 360 0119 5



The desired quantity is determined using the scale and the Ceramix is inserted into the ceramic material.



The corresponding filling quantity is determined for additional ceramic materials.



When inserting the Ceramix into the ceramic material, it must be ensured that the material is properly condensed.

To obtain the same shade mixture continuously, note down the ratio of the mixed shades. This way ceramic material is saved.



To obtain individual sample shade tabs, stir the mixed ceramic material thoroughly with a spatula. Take up with the Ceramix, press onto the fibrous pad, wet and fire.

### **Instruments**

## MagicContrast




Scale 1:1

# MagicBrush



Scale 1:1

# Magic...



The MagicContrast brushes = black and MagicBrush = brown feature absolutely identical functional chacteristics!




The pointed shape of dry brush hair is immediately restored by wetting, tapping off or vibrating the brush lightly.



Thanks to the strong contrast, any contamination such as dust or dry ceramic particles can be clearly recognized.







From the spatula shape to the original shape by rinsing the brush in liquid and then tapping if off or vibrating it lightly with a suitable instrument.



High elasticity to pick up specific quantities of ceramic material more easily.



The optimized springiness restores the shape of the brush tip immediately after picking up ceramic material. Additionally, the respective quantity can be easily assessed thanks to the contrast.



The spatula shape can be easily achieved after pressing the tip with two fingers. This way, the brush can be shaped individually.



Large quantities of ceramic material are picked up and time-saving layering is achieved thanks to the elasticity and the springiness of the brush hair.



The spatula shape does not reduce the stability of the brush hair; hence less time is required for applying the ceramic material to the framework.

## **Instruments**

## KoliBrush



KoliBrush – golden brown natural hair Natural hair brushes made of superior quality Kolinsky hair.



Improved design of the tip of the BigBrush is achieved tanks to the integrated spheres for simpler modelling. Fine, stable tip thanks to carefully selected hair.




The shape and quality of the hair for perfect retention of moisture and improved adhesion and application of ceramic material.



The desired elasticity is obtained by the unique design and combination. This way already applied ceramic layers will not be damaged.

Scale 1:1	Product name	Size	Qty.	REF
	KoliBrush	4, 6, 8 B	1 piece each	390 KSET 1
	KoliBrush	1	2 pieces	390 K001 0
	KoliBrush	2	2 pieces	390 K002 0
	KoliBrush	4	2 pieces	390 K004 0
	KoliBrush	6	1 piece	390 K006 0
No.	KoliBrush	8	1 piece	390 K008 0
	KoliBigBrush	8 BigBrush	1 piece	390 K008 B
	KoliBrush	1/0	2 pieces	390 KS01 0
	KoliOpakerBrush	5	2 pieces	390 KS03 0

## Unique Brush



Ceramic processing "par excellence".

Unique Brush – the truly esthetic brush.

With newly developed matt black bristles and unique design of each brush – even more benefits when processing ceramic materials.


- Perfectly shaped brush tips enables accurate placement of the precious ceramic materials
- Superior retention of moisture for prolonged modelling
- High elasticity of the brush hair so that large quantities of ceramic material can be picked up
- Spatula function for layering small gaps
- Individual color of each brush renders the tool unique
- Ergonomic shape enables fatigue-free working

The brushes are available in the assortment (in an attractive case) or individually.



Scale 1:1

**VPE** REF Product name Size Unique Brush Set 390 USET 1 1, 4, 6, 8 1 piece each 390 U001 0 Unique Brush 1 piece 390 U004 0 Unique Brush 1 piece 390 U006 0 Unique Brush 6 1 piece Unique Brush 8 1 piece 390 U008 0 Unique Brush BigBrush 1 piece 390 U008 B



The ergonomically designed shape of the brush enables fatigue-free and relaxed working.



The newly developed, matt black brush hair with extremely high elasticity and excellent retention of the mixing liquid ensures prolonged modelling.



The application-oriented shape of the brush tip enables accurate placement of the precious ceramic material. The contrast between the brush hair and the ceramic material ensures simple determination of the quantity that has been taken up.



The optimized, high elasticity of the newly developed brush hair allows also to take up large quantities of ceramic materials.



The required new shape and function are easily achieved by pressing the brush tip together with the fingers.



The spatula shape allows to build up very small interdental spaces and to separate them. It renders the brush a perfect tool for the preparation of marginal ridges.



The spatula function of the brush allows to apply ceramic material very accurately and to smoothen large areas. As a result, time during finishing is saved after firing.



The perfectly shaped brush tip is restored by rinsing the tip with water and then tapping off or vibrating it gently with a suitable instrument.



### Blue-Clip

The three clips with diamond-coated working sections - suitable for numerous small items.

- avoid time-consuming searching
- hold numerous small items safely
- functional design and appealing color
- three-point contact to hold spherical and cylindrical small items safely
- diamond-coated working section
- color contrast to items that have been taken up
- self-clamping or locking functions
- double function (spreading and clamping)

Screws, nuts, attachment components, items in the latch technique and numerous others are becoming smaller and smaller and more and more delicate. To remove these items from the packaging, clips and tweezers without diamond-coating are being used today, which include the problem that these tools do not provide sufficient grip and the items to be gripped may slip off and fall down.

Blue-Clip – the highly special clips with two functional principles. Three different 2-point and 3-point tip shapes. All working sections are fully coated with a fine diamond grain, which ensures firm grip of all small items that can be taken up with Blue-Clip using the clamping or spreading function. To achieve enhanced contrast, all three clips were coated in blue using plasma coating. As a result, the position of shiny metal items can be perfectly recognized.



Blue-Clip "Shape 1" Angled and diamondcoated tips (3-point) with cross-clamping function

REF 310 0008 1






This type of tweezers allows to put down the clamped items and the clip after taking up the small items. The force of pretension has been matched with this type of handling.



Blue-Clip "Shape 2"

Double bent and diamond-coated tips (2-point) with spreading and clamping function REF 310 0008 2





Very short and double force is sufficient for the easily hold anterior teeth clamping function of the





bent tips. The spring spreading behavior to in the inner side. The clip also allows to take up small items.



Blue-Clip "Shape 3" Single bent and diamond-coated tips (3-point) provide a spreading and clamping function REF 310 0008 3









The spring force for the spreading behavior is sufficient to easily take up small bridges in the inner side. The clamping function of the clip also allows to take up small items. Spherical or cylindrical small items can be safely grasped with the threepoint contact tips.





Self-clamping "Shape 1" thanks to cross-clamping



Locking function "Shapes 2 + 3"

Raster function with quickrelease lock. To reduce the stress on your fingers, the clamping lock of both clamps "Shape 2 + 3" snaps into place in the perfect position. When the working process is completed, the securing function can be deactivated by pressing the clip with the thumb; the object that has been clamped in before, can be removed again.

## Ball Clip



The holding system for crowns and bridges..

Simple and reliable fixation of crowns and bridges to be veneered with composite or ceramic. Easy-touse locking function. No reworking or remakes of dropped crowns and bridges.

#### Additional applications

- · Application of oaque
- Adding of ceramic
- · Staining and glazing
- Sandblasting and steam cleaning
- Reliable hold while condensing
- Removal of the ceramic structure from the model
- Safe placement on the firing try
- Removal tool for frameworks and raw fired structures

#### Assortment

1 Ball Clip holding device (ca. 155 x 15 mm) 1 Ball Clip holding pins, ring with 25 pieces

REF 310 0008 4

100 Ball Clip holding pins 4 rings with 25 pieces each REF 310 00H8 4





25 ball heads on each moulding... Removing a ball head from the moulding...)









Firm grip of single crowns and bridges in all positions. Simple release of the clip function and avoidance of the risk of damage to the restoration caused by dropping.

















Since the diamond-coated boring in the functional area has been matched with the diameters of the head shapes, perfectly safe hold of all single crowns and bridges is ensured. The sophisticated working area offers an almost unlimited number of different positions.






# Repositioning tweezers



Secure holding of plastic and ceramic teeth during repositioning in the cuvette and boiling-out. The special pincer tips of the tweezers take secure hold of teeth and other small parts and permit fast working.







The pointed pincer tips permit secure holding of teeth and other small parts.



 special pincer tips for secure holding
 pincer tips of hardened material for a long working life
 no slipping of small parts – no irritating searching
 fine tips for narrow areas

The special and well designed denticulation of the pincers provides optimal hold security.



Secure holding of teeth is not possible with normal tweezers. Time-consuming searching is eliminated.



Small parts such as screws or attachments are gripped easily and securely. A useful instrument particularly for implantology.




On completion there is always a problem – the repositioning of the teeth! The special fine pincers at the tweezers tips permit secure gripping of the teeth.

## Activating pliers



The problem: Conical and telescopic crowns have lost their friction.

The solution:
Activating pliers - Pliers which recreate the retentive forces for telescopes which have lost their function.



bredent Activating pliers REF 320 0043 0

Provide conical and telescopic crowns with "new" friction easily and quickly.



The friction zones in the outer coping create 3 new contact areas between the inner and outer copings. This restores the retentive friction. Should the unit be activated too much, the surface can be trimmed to reduce the friction.




These pliers have a ball and socket for creating one or several new friction zones. The long lever of the pliers enables the forces to be applied as required.



The activating pliers can also be used to reduce a friction zone which is too retentive. If necessary, the facing should be removed for activating the unit and replaced again later.

## Novo-Grip



Sharpenable, sintered diamond coated inserts for conical pliers with "grip".



Novo-Grip pliers 1 pair of pliers + 2 small inserts + 1 Allen key REF 310 0011 3



Novo-Grip small insert 2 pieces REF 310 0001 B

Novo-Grip

2 pieces

standard insert

REF 310 0001 A





1:1

Diabolo cleaner grindstone for inserts 1 piece REF 340 0100 0

Set screws M3 4 pieces REF 310 0011 2





Exchangeable insert with a diameter of 2.35 mm, also for small primary crowns.





Hardened shafts provide high stability even if strong pressure is exerted.

Rotatable



Worne inserts can be rotated. New diamond grains ensure that inner surfaces of crowns can be held safely again.

Sharpenable



To regain the maximum abrasiveness, the inserts are clamped into the handpiece and new diamond grains are obtained on the surface using the grindstone.

### **Instruments**

## Pollygrip



Grips all crowns, bridges and inlays firmly for finishing and polishing precisely.

Crown holder, wide 1 piece

REF 360 0100 0

Crown holder, narrow 1 piece

REF 360 0099 0

Replacement parts:

Special rubber sleeves 100 pieces

REF 360 0096 0



The rubber sleeves can be exchanged and grip firmly during all procedures.



Even minute inlays can be held firmly and without causing harm.

### Assortment

23 pieces

- 1 Pollygrip
- 1 Crown holder, wide
- 1 Crown holder, narrow
- 20 Special rubber sleeves

REF 360 0095 0

#### Accessories:



Crown holder, wide 1 piece REF 360 0098 0



Crown holder, narrow 1 piece REF 360 0097 0

## Transfuser



Transfuser 1 piece. REF 390 S000 1 4 pieces, REF 390 S000 4

#### Gentle and safe transfer of investment compound and gypsum.

While investing and casting an impression, quite frequently air gets entrapped - especially in narrow and inverse areas - which may cause undesired bubble formation. So far it was merely possible to utilize instruments and brushes for the preparation. But sharp-edged instruments bear the risk of damaging both modellation and impression. While using brushes, it is possible to transfer excess fluids, and if the brush is rather dry it may withdraw some moisture. Both may affect the expansion behavior of investment compounds. The curved and flexible silicone tip of the Transfuser allows

a gentle, damage- and bubble-free application prior to standard filling procedures. The smooth and dense surface provides excellent gliding characteristics for gypsums and investment compounds, and thus a bubble-free material



Highly flexible, soft tip and thus no risk of damaging the wax model.



Narrow areas are thoroughly filled by gentle application without modellation damage.




Optimal material flow and bubble-free condensation with the Transfuser.



Smooth and gentle filling of extreme hollows with gypsum in case of silicone impressions.



No risk of damage during the filling process.



Each area is well accessible and allows for optimal and bubble-free material distribution.



# Wax adapter



Quick and safe adapting of prefabricated wax patterns in the CoCr technique.

Wax adapter REF 360 0120 5



With the flat side, clasps and sublingual clasps can be perfectly and correctly placed on the investment model and pressed against it.





When using the wax adapter, prefabricated wax patterns will no longer be damaged or deformed by the special silicone but safely attached to the model.





The rounded side is perfectly suited for stippled maxillary plates or retentions. The fine wax patterns will not be deformed either.



### Universal screwdriver set



Screwdriver set for 98% of all screws available on the market. To be inserted into the torque ratchet, adjustable from 10 to 40 Ncm. This way screws can be turned in correctly and safely.



Universal screwdriver set to loosen and tighten all types of screwed implant abutments.



Universal screwdriver set with instruments REF 310 0001 2

On the lid you can find important information required for the quick selection of the necessary screwdriver and the torque needed to tighten the screw.

Universal screwdriver set, without instruments REF 310 0001 1



Torque ratchet REF 330 0115 5


Torque adjustable from 10 to 40 Ncm.



Screwdriver long Screwdriver Torx 6 REF 310 0010 1 [6] Screwdriver slotted 1.6 REF 310 0010 2 REF 310 0010 3 3 slotted 2 Screwdriver Screwdriver 0.03" only available as short type Screwdriver Allen 0.05" REF 310 0010 5 0.9 Screwdriver Allen 0.9 REF 310 0010 6 Screwdriver Allen 1.0 REF 310 0010 7 Screwdriver Allen 1.2 REF 310 0010 8 REF 310 0010 9 Screwdriver 9 Allen 1.8 Screwdriver 10 Hexagon 2.5 REF 310 0011 0 REF 310 0101 1 Screwdriver 11 Square 1.3 REF 310 0101 2 Screwdriver 12 Torx 5.5

Screwdriver short Screwdriver REF 310 00K0 1 1 short Torx 6 Screwdriver REF 310 00K0 2 2 short slotted 1.6 Screwdriver REF 310 00K0 3 3 short slotted 2 Screwdriver Allen 0.03" REF 310 00K0 4 4 short Screwdriver Allen 0.05" REF 310 00K0 5 5 short Screwdriver 6 Allen 0.9 REF 310 00K0 6 short Screwdriver 7 Allen 1.0 REF 310 00K0 7 short Screwdriver 8 short Allen 1.2 REF 310 00K0 8 Screwdriver 9 short Allen 1.8 REF 310 00K0 9 Screwdriver 10 Allen 2.5 only available as long type REF 310 00K1 1 Screwdriver 11 short Square 1.3 REF 310 00K1 2 Screwdriver 12 short Torx 5.5

# Universal screwdriver set for contra-angles



Screwdrivers with seating for contra-angles. Thanks to the integrated torque they simplify turning in screws with special motors. In conjunction with the adapter, the screwdrivers can also be used with the torque ratchet.

Universal Screwdriver-Set for contra-angles, with instruments REF 310 W001 2 Universal Screwdriver-Set for contra-angles, without instruments REF 310 W001 1



#### Accessories:



Torque ratchet REF 330 0115 5

Torque adjustable from 10 to 40 Ncm.

Ratchet adapter REF 580 0116 8



_			Screwdriver long	
6	Screwdriver	1	Torx 6	REF 310 W010 1
	Screwdriver	2	slotted 1.6	REF 310 W010 2
	Screwdriver	3	slotted 2	REF 310 W010 3
	Screwdriver	4	0.03" only availab	ole as short type
0.05"	Screwdriver	5	Allen 0.05"	REF 310 W010 5
0.9	Screwdriver	6	Allen 0.9	REF 310 W010 6
1	Screwdriver	7	Allen 1.0	REF 310 W010 7
1.2	Screwdriver	8	Allen 1.2	REF 310 W010 8
1.8	Screwdriver	9	Allen 1.8	REF 310 W010 9
2.5	Screwdriver	10	Hexagon 2.5	REF 310 W011 0
1.3	Screwdriver	11	Square 1.3	REF 310 W101 1
5.5	Screwdriver	12	Torx 5.5	REF 310 W101 2

		-	Screwo	lriver short	
6 9	Screwdriver	1	short	Torx 6	REF 310 W0K0 1
<u> </u>	Screwdriver	2	short	slotted 1.6	REF 310 W0K0 2
	Screwdriver	3	short	slotted 2	REF 310 W0K0 3
0.03"	Screwdriver	4	short	Allen 0.03"	REF 310 W0K0 4
0.05"	Screwdriver	5	short	Allen 0.05"	REF 310 W0K0 5
0.9	Screwdriver	6	short	Allen 0.9	REF 310 W0K0 6
1 9	Screwdriver	7	short	Allen 1.0	REF 310 W0K0 7
1.2	Screwdriver	8	short	Allen 1.2	REF 310 W0K0 8
1.8	Screwdriver	9	short	Allen 1.8	REF 310 W0K0 9
9	Screwdriver	10		Allen 2.5 only av	ailable as long type
1.3	Screwdriver	11	short	Square 1.3	REF 310 W0K1 1
5.5	Screwdriver	12	short	Torx 5.5	REF 310 W0K1 2

# Screwdriver long



Screwdriver long 1 piece REF 330 0081 2

The long screwdriver allows perfect visual control of the horizontal path of screwing in the laboratory. The screw connection can be more easily achieved by the dentist. For screws with 0.9 mm hexagon socket.

## Screwdriver short



Screwdriver short 1 piece REF 330 0069 0

Ideal for practice and laboratory. The grooved handle simplifies turning in of screws since safe hold is ensured. For screws with 0.9 mm hexagon socket.

## Screwdriver for contra-angles



Screwdriver for contra-angles 1 piece REF 330 0081 3

For mechanical turning in of screws with 0.9 mm hexagon socket. The use of special motors allows to control the torque.

### Screwdriver-Set



#### Assortment

3 pieces

1 x Screwdriver long

1 x Screwdriver short

1 x Screwdriver for contra-angles

REF 330 0081 0

## Screwdriver is



Screwdriver is for contra-angles 1 piece REF 460 0001 0



Screwdriver is manual short 1 piece REF 460 0001 1

Special screwdrivers for the vks-oc rs abutments. Suitable as manual screwdriver and for contraangles for enhanced control of the torque with special motors.

## Screwdriver for stud-head screw



Screwdriver for stud-head screw 1 piece REF 330 0116 4

Screwdriver for the stud-head screw vks-oc/sg 1.7 exchangeable stud.



Rotating tools

Only selected, high-quality raw materials are used for the manufacture of burs. This selection makes it possible to create first-class burs for surface processing of the most varied materials.







The bredent Order Number System for Diatit	
and tungsten carbide tools	
Finding the desired tool quickly	
The Cut Types of the bredent Diatit and tungste	
carbide tools	
Cut overview	. 414
Generation M	
"Generation M" relief burs	. 416
Microburs	
Fissure tool	. 418
Tungsten carbide tools with relief	
Comparison:	
bredent tools with and without relief	
Diatit wear protection	. 419
Microbur with relief	
Rapidy Microbur	
Tungsten carbide tools with relief	
Diatit and tungsten carbide tools	
Diatit power burs	
Processing of titanium	. 438
Processing of titanium surfaces	
Titanium Finishing Set	. 439
Tools with relief for the milling technique	
Wax bur, Profile bur, Polishing bur	. 440
Milling wax / Bur application	
Biotec milling wax	
Bur application	. 441
Burs	
Parallel burs	
Parallel burs for titanium, non-precious	
Conical burs	. 446
Conical burs for titanium, non-precious	
and precious metal	
Groove-shoulder attachment	
Wax burs	
Cross-cut burs	
Wax burs and diamond grinding tools	. 451
Milling and drilling oil	
Milling and drilling oil	. 449
Diamond grinding wheels	
Survey	
Giflex-TR	
Giflex-TR Master x-tray	
Diamond grinding wheels	
Ceraflex	
Microflex	
Transflex-T	
Transflex	
Ultraflex, Superflex, Flexibel, Elastisch	. 456
Diamond grinding	
Survey	
Diacryl Grinding Instruments dcs	
Diagen-Turbo-Grinder, fine	
Diagen-Turbo-Grinder, coarse	
Set-up grinding tool	
Special Diamonds for the Veneering Technique	
Diamond grinding tool dsl	
Diabolo	
FG-Diabolo for uirconium oxide	. 470
Those was stone ton diamend animaline to all	

Diabolo Cleaner.....

## The bredent Order Number System for Diatit and tungsten carbide tools

#### Color coding

SH

orange

Finding the cut type quickly with the help of the color code on the shaft of the bur.

NF KF red KS black none NH KM blue GG none orange MH QM light-blue KC purple orange N/MH orange/blue **QG** white KT silver grey orange/green KG green GH M gold

Tool shape, ISO-Number

Three numbers indicate the tool shape according to ISO.

# bredent D263 KG 60

#### Letter in initial position

N = special tool for non-precious metal alloys

**H** = Tungsten carbide

**D** = Diatit wear protection<sup>1</sup>

B = special tools (drills), e.g. fissure tool

**F** = special tools for the milling technique

**S** = silicone bur

<sup>1</sup> For details on the Diatit wear protection see page 329

#### Cut

bredent offers the tools mentioned above in 14 different cut designs. The cut types are marked by the combination of two capital letters. Size Diam

Diameter at the largest point of the working element in tenth of millimeters.

For details on the cut designs see page 325

## Finding the desired tool quickly

This catalogue offers the possibility to always find the desired tools in the fastest possible way. The method of determination is either based on the shape or on the cut of the tool.

#### Determination based on the shape

The outer two columns of the double page 326/327 show all bredent tool shapes.

The desired shape can be selected there. Then a bredent cut type is selected in the row of the desired tool shape.

A page number is indicated in the box of the selected cut. Further information on the selected tool is provided on this page.

Picture	REF								Cut						
Scale 1:1	•••	• NF	NH	МН	GH	SH	KF	KM	QM	QG	KG	KS	GG	КС	KT
	D13723 H13723						426	429	432		434				
•	• • •	••••	• • • •	• • • •	• • •	•	•	•	•	• • • •	•		•	• • • •	• • • •
Picture of the tool in the original size.	Here are th letters of tl		, -	oe. ind ava	ler numb ication o ilable wi mm.	f the cut	. This is	cut t For r	tool is a types KF, more det see page	KM, QM ailed info	, KG. orma-	burs	low findi quickly t been ind	he color	codes

#### Orientation based on the cut

From page 328 all tools are arranged according to the cut.

The arrangement includes fine and coarse cuts as well as special cuts for chrome-cobalt alloys titanium.

#### ISO numbers

are indicated for all tools to ensure enhanced comparability. These internationally standardized numbers feature 15 digits. The numbers include the following information: - 3. digit:
 Materials of the working element

7. - 9. digit: Shape of the working element 13. - 15. digit: Diameter of the working element

509 104 001215 023

4. - 6. digit: Shaft type 10. - 12. digit: Cut



## The Cut Types of the bredent Diatit and tungsten carbide tools

• for processing of any dental material

• single cut instead of "double" cross cut



Normal cut






NH: Normal cut with relief



MH: Central cut with relief



GH: Coarse cut with relief



SH: Super coarse cut with relief



KF: Cross cut Fine



KM: Cross cut Medium



OM: Horizontal cut Medium



OG: Cross cut Coarse



Cross cut Coarse KS:

Cross cut

KG:



Super coarse GG:

Straight cut

Coarse



KC: Cross cut Chrome-Cobalt



KT: Cross cut Titanium





• for processing of precious metals, non-precious metals, resin and ceramics • fine removal of material, very smooth surface of object, low vibration running protects the wrist of the technician and the drive

• relief: wide, stable cutting edge, extended service life

• easy removal of material with perfect control, smooth surface of object

• for processing of precious metals, non-precious metals, resin, plaster

• excellent removal of material and very smooth rotation, smooth surface of object

- relief: wider, stabler cutting edge for extended service life, enhanced cutting performance
- for coarse treatment of precious metals, non-precious metals, resins; in individual cases also for treatment of plaster
- excellent removal of material, low-vibration running and extended service life due to relief
- for processing of plaster and carrying out particulary coarse work on resin surfaces
- excellent removal of material and particularly smooth material surface due to relief
- no loading with shavings due to larger cut spaces
- mainly for more delicate types of work on precious and non-precious metals, resins and ceramics
- moderate and accurate removal of material, smooth surface of object
- for finishing of larger surfaces on precious metals, non-precious metals and resins, in individual cases also on plaster
- efficient removal of material, smooth surface of object, smooth running of tool
- universal application possibilities, therefore reduced frequency of tool exchange
- suitable for finishing of larger surfaces as well as for more delicate work on precious and non-precious metals and resin, therefore reduced frequency of tool exchange
- very fine, economic removal of material, smooth surface
- high smoothness of running protects drive and wrist
- especially for processing of silicones
- very efficient and accurate removal of soft materials
- for coarse and efficient pretreatment of large surfaces on precious metals, non-precious metals and resins, in individual cases also on plaster
- extensive removal of material, larger surface roughness than the finer bredent cut types
- especially for processing of plaster, also suitable for very coarse types of work on resin
- extensive removal of material
- the size of the individual cut space avoids loading with shavings
- to perform cuts in resin or shellack plates
- very economic cutting of plates
- · single, straight cutting edges
- especially for processing of chrome-cobalt alloys
- excellent removal of material, smooth surface
- the characteristic feature of this tool: the resulting metal swarf cause fewer irritations to the skin since they are larger and exhibit a coarse structure
- especially for processing of titanium
- the special dent of this cut increases the cutting volume which reduces the friction. Overheating of titanium is avoided.
- economic, careful removal of material, smooth surface
- offers multiple applications
- for plaster, plastic and high-end thermoplastic
- · quick removal for efficient working



- creates smooth surfaces and reduces the amount of work required
- for precious, non-precious metals, plastics



Medium serration

M3:


- for time-saving surface processing
- · good material removal with very smooth surfaces

fine serration creates very smooth surfaces and makes polishing easier

for all materials



- Long service life for cost-effective working





# **Cut overview**

Picture								C	ut						
Scale 1:1	REF	NF	NH	МН	GH	SH	KF	KM	QM	QG	KG	KS	GG	KC	KT
Size 06	B15302-06 only available in tungsten carbide	418													
Size 23	H001 NH 04-31 only available in tungsten carbide		423												
•	D001 14 only available in Diatit														438
•	D001 23 H001 23 H010 08-16			420				429			434				438
Size 16	H010 08-16		423												
	D137 23 H137 23						426	429	432		432				
	D141 23 H141 23 N141 23			420+424				429							
	H161 60									423					
	D184 16 H184 16			420			426	429							
	D187 23 H187 23 S187 23						426	430		433	434				
	D194 23 H194 23						426	430			434				438
	D194 40 H194 40 N194 40				421+425	422	426	430			434			437	438
	D194 50 H194 50				421		426	430			434				438
	D194 60 H194 60					422						436			
	D194 70 H194 70					422						436			
	D198 23 H198 23 N198 23			424			427	430							438



# **Cut overview**

Picture								Cı	ut						
Scale 1:1	REF	NF	NH	МН	GH	SH	KF	KM	QM	QG	KG	KS	GG	KC	KT
	D200 23 H200 23						427	430			435				
	D225 23 H225 23						427	430			100				
	D237 23 H237 23			420			427	431							
	D237 65 H237 65 S237 65									433	435				
	H244 23				421										
	D251 60 only available in Diatit													437	
Size 16	D257 16/23 H257 16/23							431							
	H263 30 D263 40 H263 40 N263 40				421		428	431							438
	D263 60 H263 60 S263 60 N263 60				421+425	422				433	435				
	D274 60 H274 40/60 N274 40				421+425	422					435				
	D277 14 H277 14 N277 14			420+424				421							
•	D277 23 H277 23			420+424				421							
-	D289 23 H289 23			421			428	432							
	D292 23 H292 23						428	432			435			437	
Size 23	D468 16/23 H468 16/23												436		


## "Generation M" relief burs



#### Relief burs with new multifunctional cut.

Thanks to a significant clearance angle, the "Generation M" relief burs provide extended service life and reduce the number of repurchases. Adequate support of the cutting edges avoids chipping and leads to improved cutting performance. A smooth surface of the object to be processed is obtained and subsequent timeconsuming polishing work is reduced.

Thanks to the optimized geometry, chips can not penetrate into the skin when processing metals. This way the dental technician's health is protected.



#### High-quality materials offer reliability.

Only specially selected, high-quality raw materials are used for the production of burs to obtain topclass tools and enable surface processing of various materials.

Quality-oriented production of the elaborately ground tools results in excellent ease of use. In combination with the relief cut technology, the exact concentricity avoids chipping of the cutting edges and protects the drive of the micromotor to enable perfect and efficient processing of high-quality dental restorations.

#### Generation M - the progressive bur concept

- suitable for soft and hard materials;
   efficiency is increased
- reduces the number of burs and facilitates organisation and control of the workplace
- considerably extended service life thanks to multifunctional cut
- reduction of the amount of work thanks to smooth surfaces

A new level of surface treatment is achieved by the change of the clearance angle and the proven support of the cutting edges of the new "Generation M" bur concept. The variety of materials to be processed in the laboratory requires a number of different processing tools, which is minimized by the "Generation M" burs and their range of applications. As a result, the efficiency during processing is increased and the number of tools to be used is reduced!



#### M7

The super coarse cut is suitable for various applications on plaster, denture acrylics and thermoplastic resins. Also suitable for fast removal of material on NPM alloys.



#### M5

The coarse cut produces a smooth surface and enables extensive removal of material on precious metals, NPM alloys and acrylics. The special clearance angle extends the service life of the bur.



#### М3

The medium cut produces a very smooth surface and hence reduces the amount of work during subsequent surface treatment. Processing in areas difficult to access is facilitated.



#### М1

The fine cut produces a very smooth surface and facilitates polishing. The relief cut extends the service life of the tool, which can be used over a longer period.



## "Generation M" relief burs



**REF H263 M7 40** ISO-No. 500 104 263220 040



REF H263 M5 40 ISO-No. 500 104 263220 040



**REF H274 M7 40** ISO-No. 500 104 274220 040



REF H274 M5 40 ISO-No. 500 104 274220 040



**H274 M7 16** ISO-No. 500 104 274220 016



**H274 M5 16** ISO-No. 500 104 274220 016



**H277 M3 60** ISO-No. 500 104 277190 060



H263 **M7** 40 H263 **M5** 40 The proven, universal design with a wide indication range for metals and resins/acrylics.



H274 M7 40 H274 M5 40 Fast removal of material enables efficient working. The tapering tip allows perfect finishing even in areas difficult to access.



H274 M7 16 H274 M5 16 Surface finishing is completed with the small flame to obtain a very smooth surface.



H277 M3 60 Time-saving processing is guaranteed by the combination of smooth surface and fast finishing with the medium cut bur.

#### Accessories:



#### For perfect high gloss

The proven polishing paste filled with natural diamonds produces perfect high gloss in no time. The ideal complement to the "Generation M".

#### Zi-polish

Pre- and high gloss polishing paste 5 g
REF 360 1002 5



Rodeo round brushes 15 pieces Ø 18 mm REF 350 0096 0

## **Microburs**

## Fissure tool





Due to the shape of the tool smoothening of cusp "slopes" at inaccessible spots is possible. The extremely small diameter allows excellent smoothening in the deep area of the fissure so that polishing of occlusal surfaces is simplified. Well polished occlusal surfaces reduce the accumulation of plaque. This tool offers the dental technician excellent design possibilities.

#### Assortment

REF 330 0082 6

6 pieces, 2 pieces each Fissure tool ISO-No. 500 104 153006 002 ISO-No. 500 104 153006 004 ISO-No. 500 104 153006 006



Perfect fissures with the smallest fissure tool in the world Diameter 0.2 mm



Additionally, the special cutting edge geometry allows recontouring of ceramic occlusal surfaces prior to glaze firing. Therefore it offers new possibilities of designing occlusal surfaces to the ceramic specialist.



The bredent fissure tool in a magnification x 100


Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm			10-20,000	10-20,000	15-20,000	15-20,000


## Comparison: bredent tools with and without relief



# Triple service life comparated to conventional bredent cut.

The latest tungsten carbide tools by bredent are provided with a relief during an additional manufacturing process. The relief supports the sharp cutting edge to avoid breakage of the edges. This way the service life of the relief tools is three times longer than the one of conventional tools. Additionally, the relief allows to optimize the machining angle so that an excellent cutting performance is achieved.





## Diatit wear protection

enhanced cutting performance



#### Smooth running from the very beginning

Diatit tools feature particular wear protection: Diatit. This is a special material which is added into the surface of the bur after it has been produced. It hardens the tool surface and reduces the surface friction. This comprehensive hardening process results in a tool which features very smooth rotation and precise cutting performance from the very beginning – and this is provided over a considerably extended period. Accordingly, accurate removal of material is ensured. Additionally, the service life of the tool (compared to uncoated tungsten carbide burs) is increased considerably by the hardening process.


Tungsten carbide structure bredent tungsten carbide tools consist of a metal sintering material with a very fine grain size. Additionally, Diatit tools



are subject to a hardening process after the cut has been completed. This hardening process reaches into the gaps between the crytals in a depth of up to 100 mm.

# bredent tool with Diatit wear protection.

The surface of the tool is smoothened so that the friction is reduced. The swarf comes



off the tools more easily. This results in smoother running of the tool.

#### bredent tool without Diatit wear protection. Additionally, - compared to

uncoated bredent tungsten carbide tools - jagging of the





cutting edges of Diatit tools is avoided due to the wear protection.
Compared to uncoated bredent tungsten carbide burs the hardness rises up to 3700 HV (compared to 1850 HV) and results in an increased service life of the tool.

C Rotating tools

Fax (+49) 0 73 09 / 8 72-4 44

### Cut: MH







This bur offers numerous possibilities of application; in this picture it is used in the CoCr technique.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	10-20,000	10-20,000	10-20,000	10-20,000	15-20,000




Material Tungsten carbide REF H141 MH 23 ISO-No. 500 104 141190 023



The tool H141 MH 23 used for grinding a ditch. Smooth, chatter-free running of the relief cut increases the reliability during the application.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	10-20,000	10-20,000	10-20,000	10-20,000	15-20,000



Material Tungsten carbide REF H184 MH 16 ISO-No. 500 104 184190 016



The high cutting performance of the relief cut allows more efficient treatment of ceramics; the picture shows smoothening of the transitions of metal/ ceramic.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Tungsten carbide Material REF H237 MH 23 ISO-No. 500 104 237190 023



Tools with relief produce a particularly smooth micrograph. During milling of ceramic materials a silky-mat surface is obtained which is suitable for glaze-firing without any additional treatment. Therefore relief tools ensure high efficiency when processing ceramic

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide REF H277 MH 14 ISO-No. 500 104 277190 014 REF H277 MH 23 500 104 277190 023 ISO-No.



Slender designs allow the use of relief tools even for highly precise work.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000

### Cut: MH and GH




Material	lungsten carbide
REF	H289 MH 23
ISO-No.	500 104 289190 023



The tool H289 MH 23 is particularly suitable for grinding VMK veneer surfaces.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000







Due to the extended service life the relief cut saves material costs Accordingly, costs can be reduced considerably particularly when processing VMK frameworks.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide

REF H244 GH 23

ISO-No. 500 104 244220 023

The shape allows highly accurate finishing of metal structures.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide

REF H263 GH 30

ISO-No. 500 104 263220 030

REF H263 GH 60

ISO-No. 500 104 263220 060



Larger amounts of plaster are quickly removed using the large surface of the bur. A smooth surface is created simultaneously.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-15,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide

REF H274 GH 40

ISO-No. 500 104 274220 040

REF H274 GH 60

ISO-No. 500 104 274220 060




A smooth object surface can be produced with the relief cut. This provides considerable advantages when processing resins.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-15,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000

### Cut: SH and NH








The relief cut produces a very smooth object surface. Smooth and chatter-free cutting of the bur leads to safer use by the technician and protects the joints.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-15.000	8-12.000				








Tools with SH cut have been especially developed for processing of plaster. Smooth running of the relief avoids coarse edges in the plaster.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	8-12.000	8-12.000				








Bur with relief cut for quick removal of denture resin. Also perfectly suitable for plaster.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20.000	10-20.000	10-20.000	10-20.000	10-20.000	15-20.000







The inverted cone is perfectly suitable for shaping occlusal surfaces. Simultaneously, a brilliant ceramic surface is achieved thanks to the relief cut.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	10-20,000	10-20,000	10-20,000	15-20,000	15-20,000

### Cut: NH

# Rapidy Microbur



#### Optimum cutting performance and long service life due to relief.

Material


Tungsten carbide

The Rapidy Microbur also features a relief. Due to this modern edge geometry the Rapidy exhibits a particularly high cutting performance as well as extraordinarily smooth running.

The dental technician is able to benefit from these properties especially when extremely hard materials have to be processed in a fast and precise manner, e.g. when shaping ceramic or non-precious metal alloys. Even on these materials the Rapidy ensures extensive removal of material and creates a particularly smooth object surface. Additionally, the triple service life that is ensured by the relief allows to save costs.



H001 NH 04: The fine cutting performance of the Rapidy Microbur offers excellent possibilities of design to the ceramic specialist.



		. 5					
	QTY	1 piece		5 pieces	10	pieces	
	REF	H001 NF	l 04	330 00	50 4 33	30 0100 4	
	ISO-No.	. 500 104	001006 004				
	REF	H001 NF	l 05	330 00	50 5 33	30 0100 5	
	ISO-No.	. 500 104	001006 005				
	REF	H001 NF	1 06	330 00	50 6 33	30 0100 6	
	- ISO-No.	. 500 104	001006 006				
	REF	H001 NF	l 07	330 00	50 7 33	30 0100 7	
	ISO-No.	. 500 104	001006 007				
	REF	H001 NF	l 08	330 00	50 8 33	30 0100 8	
	ISO-No.	. 500 104	001006 008				
	REF	H001 NF	1 09	330 00	50 9 33	30 0100 9	
	ISO-No.	. 500 104	001006 009				
	REF	H001 NF	ł 10	330 00	51 0 33	30 0101 0	
	ISO-No.	. 500 104	001006 010				
9	REF	H001 NF	l 12	330 00	51 2 33	30 0101 2	
	ISO-No.	. 500 104	001006 012				
•	REF	H001 NF	114	330 00	51 4 33	30 0101 4	
	ISO-No.	. 500 104	001006 014				
g	REF	H001 NF	ł 16	330 00	51 6 33	30 0101 6	
	ISO-No.	. 500 104	001006 016				
9	REF	H001 NF	ł 18	330 00	51 8 33	30 0101 8	
	ISO-No.	. 500 104	001006 018				
-	REF	H001 NF	ł 21	330 00	52 1 33	30 0102 1	
	ISO-No.	. 500 104	001006 021				
g====	REF	H001 NF	1 23	330 00	52 3 33	30 0102 3	
	ISO-No.	. 500 104	001006 023				
Q.	REF	H001 NF	ł 31	330 00	53 1 33	30 0103 1	
	ISO-No.	. 500 104	001006 031				
Application field		Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr den- ture/NPM	Ceramic
Working speed rpm		10-20,000	10-20,000	10-20,000	10-20,000	15-20,000	15-20,000
Troiking specu ipin		10-20,000	10 20,000	10-20,000	10-20,000	13-20,000	13-20,000

### Cut: MH/NPM special burs



Material	lungsten carbide		
REF	N141 MH 23		
ISO-No.	500 104 141190 023		



The tool N141 MH 23 during grinding of a ditch. Smooth, chatter-free running of the relief cut increases the reliability during the application.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	



Material Tungsten carbide

REF N198 MH 23

ISO-No. 500 104 198190 023



The tool N198 MH 23 during grinding of a ditch. The smooth, chatter-free running of the relief cut increases the reliability during the application.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	



Material Tungsten carbide

REF N277 MH 14

ISO-No. 500 104 277190 014



Slender designs allow the use of the relief tools even for very precise work and create an extremely smooth surface.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	




Assortment
Package contains 7 burs
REF 330 0117 0

Thanks to the special relief cut, extended service life of these burs for non-precious metal alloys is achieved. The modified angle of the relief increases the abrasiveness and leads to a better surface quality which allows to save a considerable amount of time.

## Cut: GH/NPM special burs








Thanks to the increased service life the relief cut saves tool costs. Accordingly, costs can be considerably reduced when processing VMK frameworks made of non-precious metal alloys.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	





Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	



Material Tungsten carbide

REF N263 GH 60

ISO-No. 500 104 263220 060

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	



Material Tungsten carbide

REF N274 GH 40

ISO-No. 500 104 274220 040

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					15-20,000	



### Cut: KF



	Material	Tungsten carbide	Diatit
arantmee -	REF	H137 KF 23	D137 KF 23
	ISO-No.	500 104 137140 023	509 104 137140 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H184 KF 16
 D184 KF 16

 ISO-No.
 500 104 184140 016
 509 104 184140 016



Due to the fine cut a smooth surface of the object is obtained. This slender tool is particularly suitable for finishing of veneers.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H187 KF 23
 D187 KF 23

 ISO-No.
 500 104 187140 023
 509 104 187140 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide Diatit REF H194 KF 23 D194 KF 23 ISO-No. 500 104 194140 023 509 104 194140 023 D194 KF 40 REF H194 KF 40 ISO-No. 500 104 194140 040 509 104 194140 040 REF H194 KF 50 D194 KF 50 ISO-No. 500 104 194140 050 509 104 194140 050



The D194 KF 23 is particularly suitable for finishing of metal-ceramic frames.



The fine micrograph of the KF cut simplifies polishing of the metal surface.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000

## Cut: KF




Material	Tungsten carbide	Diatit
REF	H198 KF 23	D198 KF 23
ISO-No.	500 104 198140 023	509 104 198140 023



The slender design and the smooth micrograph of the D198 KF 23 ensure excellent suitability for processing of partial frameworks.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000





Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	H225 KF 23	D225 KF 23
ISO-No.	500 104 225140 023	509 104 225140 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	H237 KF 23	D237 KF 23
ISO-No.	500 104 237140 023	509 104 237140 023



D237 KF 23: Due to the fine cut a smooth surface can be achieved on hard alloys as well.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000

## Cut: KF



iviatciiai	rungsten caroluc	Diatit
REF	H263 KF 40	D263 KF 40
ISO-No.	500 104 263140 040	509 104 263140 040



The D263 KF 40 can be used for a wide range of applications in the CoCr technique.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H289 KF 23
 D289 KF 23

 ISO-No.
 500 104 289140 023
 509 104 289140 023



The KF cut is perfectly suitable for finishing of metal ceramic veneers.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H292 KF 23
 D292 KF 23

 ISO-No.
 500 104 292140 023
 509 104 292140 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		12-18,000	15-20,000	15-20,000	15-20,000	15-20,000


Cut: KM



Material	Tungsten carbide	Diatit
REF	H001 KM 23	D001 KM 23
ISO-No.	500 104 001190 023	509 104 001190 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000







Due to the smooth running and the excellent cutting performance, this tool is particularly suitable for precise and efficient finishing.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



	Material	lungsten carbide	Diatit
(22)22	REF	H141 KM 23	D141 KM 23
	ISO-No.	500 104 141190 023	509 104 141190 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide Diatit

REF H184 KM 16 D184 KM 16

ISO-No. 500 104 184190 016 509 104 184190 016



Application example.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000

### Cut: KM





Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000





Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H198 KM 23
 D198 KM 23

 ISO-No.
 500 104 198190 023
 509 104 198190 023



Efficient and controlled finishing due to the slender tool with KM cut.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H200 KM 23
 D200 KM 23

 ISO-No.
 500 104 200190 023
 509 104 200190 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H225 KM 23
 D225 KM 23

 ISO-No.
 500 104 225190 023
 509 104 225190 023



D225 KM 23 for regrinding exact transitions of metal/resin.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Cut: KM



Material	lungsten carbide	Diatit
REF	H237 KM 23	D237 KM 23
ISO-No.	500 104 237190 023	509 104 237190 023



This tool is particularly suitable for designing filigree CoCr frameworks.


Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



	Material	Tungsten carbide	Diatit
•	REF	H257 KM 16	D257 KM 16
	ISO-No.	500 104 257190 016	509 104 257190 016
	REF	H257 KM 23	D257 KM 23
	ISO-No.	500 104 257190 023	509 104 257190 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000







Controlled removal of material with the bredent KM tool.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000




	Material	Tungsten carbide	Diatit
•	REF	H277 KM 14	D277 KM 14
	ISO-No.	500 104 277190 014	509 104 277190 014
	REF	H277 KM 23	D277 KM 23
	ISO-No.	500 104 277190 023	509 104 277190 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000

#### Cut: KM and QM







Finishing of precious metal alloys: D289 KM 23

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



 Material
 Tungsten carbide
 Diatit

 REF
 H292 KM 23
 D292 KM 23

 ISO-No.
 500 104 292190 023
 509 104 292190 023



Fine cutting performance and smooth surface of the object: D292 KM 23 - application example in the precious metal technique

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20,000	12-18,000	15-20,000	15-20,000	15-20,000	15-20,000



Material Tungsten carbide Diatit

REF H137 QM 23 D137 QM 23

ISO-No. 500 104 137134 023 509 104 137134 023



Precious metal technique: Smooth surface of the object, tool is running steadily.




The fine cutting performance and the smooth running of this tool ensure efficient working in the field of resins as well.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	15-20.000	15-18.000	15-20.000	15-20.000	15-20.000	15-20.000

Cut: QG







The tapering tip is particularly suited to the delicate and precise processing of plastics. It is also perfect for use in the splint technique.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm		10-20,000				



Material	lungsten carbide
REF	S187 QG 23
ISO-No.	500 104 187 023



Due to the different shapes, these burs can be used to shape silicone and also in hardly accessible areas, such as in the alveoli of the gingival mask.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic	Silicone
Working speed rpm							30.000



Material	lungsten carbide
REF	S237 QG 65
ISO-No.	500 104 237 065



Due to controlled removal of material, a smooth transition zone from the silicone to the denture resin is achieved. The ground surface is smooth and free from grooves.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic	Silicone
Working speed rpm				vascu			20.000



Material	lungsten carbide
REF	S263 QG 60
ISO-No.	500 104 263 060




Thermoformed, soft thermoplastic plates with different hardness degrees are rapidly and safely ground with these silicone burs, e.g. for sports mouth guards.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic	Silicone
Working speed rpm							18.000

### Cut: KG







D001 KG 23: Removal of bubbles in the precious metal technique.


Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000





Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000







The slender design and the fine cutting performance of the D187 KG 23 render this tool indispensable for the use in the CoCr technique.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



	Material	Tungsten carbide	Diatit
	REF	H194 KG 23	D194 KG 23
	ISO-No.	500 104 194220 023	509 104 194220 023
	REF	H194 KG 40	D194 KG 40
	ISO-No.	500 104 194220 040	509 104 194220 040
	REF	H194 KG 50	D194 KG 50
The state of the s	ISO-No.	500 104 194220 050	509 104 194220 050



D194 KG 23 for efficient CoCr processing



The coarse cuts of the D194 KG 40 (picture 3) and the D194 KG 50 (picture 2) guarantee fast and perfect finishing of resin.



Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



### Cut: KG



	Material	Tungsten carbide	Diatit
and the same of th	REF	H200 KG 23	D200 KG 23
	ISO-No.	500 104 200220 023	509 104 200220 023



D200 KG 23 for finishing of CoCr frameworks.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



	Material	Tungsten carbide	Diatit
2000	REF	H237 KG 65	D237 KG 65
ISO-N	ISO-No.	500 104 237220 065	509 104 237220 065



D237 KG 65 excellent removal of material and smooth running for efficient processing of resin

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	H263 KG 60	D263 KG 60
ISO-No.	500 104 263220 060	509 104 263220 060




Due to its shape and the coarse cut the D263 KG 60 is well suited for grinding of plaster dies..

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	H274 KG 60	D274 KG 60
ISO-No.	500 104 274220 060	509 104 274200 060



The D274 KG 60 is suitable for all materials due to this coarse cut.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	lungsten carbide	Diatit
REF	H292 KG 23	D292 KG 23
ISO-No.	500 104 292220 023	509 104 292220 023

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000

### Cut: KS and GG



	Material	lungsten carbide	Diatit
	REF	H194 KS 60	D194 KS 60
	ISO-No.	500 104 194223 060	509 104 194223 060
	REF	H194 KS 70	D194 KS 70
-	ISO-No.	500 104 194223 070	509 104 194223 070



D194 KS 60: The supercoarse cut guarantees a particularly powerful cutting performance.



D194 KS 70: The supercoarse cut is particularly effective on plaster and

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	8-12,000	8-12,000				



Material	Tungsten carbide	Diatit
REF	H468 GG 16	D468 GG 16
ISO-No.	500 104 468211 016	509 104 468211 016
REF	H468 GG 23	D468 GG 23
ISO-No.	500 104 468211 023	509 104 468211 023



Simple and precise cutting of plate material is possible due to the straight cut.



Smearing of the plate material due to overheating is avoided so that fast and reliable working is ensured.


Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm	8-12.000	8-12.000				

Cut: KC

#### KC cut: Special cut for CoCr and non-precious metal alloys.

The KC cut provides high cutting performance on hard alloys. Coarse metal swarf results which can not penetrate into the technician's skin. The KC cut rationalizes surface working and simultaneously avoids injuries to the skin.







The KC cut ensures gentle and pressure-free milling with high cutting performance.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					10-20,000	



Material Diatit

REF D251 KC 60

ISO-No. 509 104 251190 060



The KC cut by bredent guarantees efficient finishing of CoCr frames.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					10-20,000	



Material Diatit

REF D292 KC 23

ISO-No. 509 104 292190 023



The high cutting performance ensures efficient finishing of crowns and bridges made of non-precious metal alloys.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd- based	CoCr den- ture/NPM	Ceramic
Working speed rpm					10-20,000	


### Processing of titanium

#### Cut: KT

#### KT cut: Especially for processing of titanium.

At 850 °C titanium reacts with the oxygen in the ambient air and forms a surface with undesired material characteristics (e.g. discolourations, insufficient polishing capacity, embrittlement, etc.). Due to the special diagonal cut the tools with the KT cut offer a larger swarf-cutting volume so that the swarf comes off the tool more easily and the friction is reduced. This special cut avoids overheating of the titanium caused by friction heat. Consequently, this cut ensures efficient and careful removal of material and produces a smooth surface.







The different shapes of the titanium burs by bredent guarantee efficient and reliable processing of titanium frames.

Application	field

Plaster Denture resin

Veneer

Titanium

Ceramic


Working speed rpm

resin

Precious metal/pd-based

10-15,000







D194 KT 23: For controlled processing of areas difficult to access.

ISO-No. 509 104 194190 050

Application field
Working speed rpm

Plaster Denture

Precious metal/pd Veneer

Titanium

Ceramic

size 23+40: 20-25,000 size 50: 20,000



Material D198 KT 23 REF ISO-No. 509 104 198190 023



Precious metal/pd-based Application field Plaster Denture Veneer Titanium Ceramic resin resin Working speed rpm 25-30,000



Material Diatit D263 KT 40 REF ISO-No. 509 104 263190 040



Precious metal/pd-based Application field Plaster Denture Veneer Titanium Ceramic resin resin Working speed rpm 20-25,000

# Titanium Finishing Set



Finish titanium rationally using cutters, polishers, brushes and pastes developed specifically for use on titanium.

Extra-sharp blades, special blade geometry and cutter blades combined with Diatit wear-resistance, which has been proven for many years, guarantee that titanium can be finished quickly, without harming the material yet reducing heat development.





Recommended speed 20,000 r.p.m. REF D 194 KT 50





Recommended speed 20-25,000 r.p.m. REF D 194 KT 40

The titanium-Diatit-cutter grinds exceptionally abrasively yet runs extremely smoothly on the titanium surface. Thus, one can work accurately and quickly to achieve a uniformly smooth ground surface.





Recommended speed 25-30.000 r.p.m. REF D 198 KT 23





Recommended speed 10-15.000 r.p.m. REF D 001 KT 14

439

Thanks to the various shapes and sizes, even areas which are narrow and difficult to access can be finished precisely. When used at the correct speed (refer to range of speeds) and only minimal pressure is exerted, the titanium-Diatit-cutter grinds exceptionally well and lasts a very long time.



**Titapol** 150 g REF 520 0015 3 350 g REF 520 0015 4



Abraso-Star Glaze REF 520 0016 3



Titapol pre high-lustre polishing paste and Abraso-Star universal high-lustre polishing paste - the perfect combination for excellent polishing.



Beech wood stand 8Bo/HP REF 210 0043 0

#### Finishing set for titanium

REF 350 0089 0

1 Diatit tungsten carbide cutter, D194 KT 50	REF D194 KT 50
1 Diatit tungsten carbide cutter, D194 KT 40	REF D194 KT 40
1 Diatit tungsten carbide cutter, D198 KT 23	REF D198 KT 23
1 Diatit tungsten carbide cutter, D001 KT 14	REF D001 KT 14
1 Titapol pre-polishing wheel	REF 350 0087 0

1 Titapol pre-polishing cylinder 1 Round brush zwm db 19 Ø goat hair, white, mounted, double rows

1 Cotton buff, for handpiece

1 Abraso-Soft Metal CSF 2/80 chunking, black, white textile insert

1 High-luster buff, metal, 50 L/100

1 Titapol pre-polishing paste, 150 g

1 Abraso-Star asg universal high-luster polishing

1 Beech wood stand 8Bo/HP


Refill packs:

KEF D194 K1 40
REF D198 KT 23
REF D001 KT 14
REF 350 0087 0
REF 350 0088 0
REF 350 0054 0
REF 350 0065 0
REF 350 0081 0
REF 350 0083 0
REF 520 0015 3
REF 520 0016 3
REF 210 0043 0


Fax (+49) 0 73 09 / 8 72-4 44 C Rotating tools

# Tools with relief for the milling technique


### Wax bur



### Profile bur



# Polishing bur



### Biotec milling wax



Excellent milling wax with superb modelling properties. Outstanding scraping and milling properties since sticking of wax to the bur is avoided.



Biotec milling wax 28 g REF 510 0061 4



Enormous amount of time is saved due to good modelling properties since no other wax is required for the shear distributor.




Extremely accurate milling wax to produce smooth and shining surfaces during milling.




Can be used for press ceramics since the wax burns out almost entirely.

# Systematic preparation of a groove-shoulder attachment with the milling technique tool set by bredent



It is recommended to model the entire, planned crown in wax prior to starting the milling work.




In the first step a semiround shoulder with a marginal step is prepared with the wax bur F137 3W 23.




Then the approximal grooves are prepared with the groove bur F538 3H 10.



Finally, the occlusal shoulder is completed with the shoulder bur F205 3H 27.



After casting and finishing of the crown, the parallel surface is remilled with the profile bur F137 3H 23.



Remilling of the groove is carried out with the tool F 538 3H 10. During this process the groove bur should only be moved up and down in the vertical

axis.



The occlusal shoulder is reworked with the tool F205 3H 27.



Finally, a high lustre is achieved on the parallel surface using the polishing bur F137 3P 23.



The milling tools with relief by bredent produce a perfect high lustre so that additional polishing is not required.



A secondary element is modelled with the pattern resin Pi-Ku-Plast and – if required – shaped with rotating tools.




The secondary element is cast and placed onto the primary element.



The low shrinkage of the pattern resin ensures excellent precision of fit of the secondary element.

441

### Wax bur, parallel, round face






The wax bur F137 3W 23 features rounded face cut. Therefore it is perfectly suitable to prepare precise ditches in the marginal area

- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm

### Profile bur, parallel, round face






The round face cut of the profile bur F137 3H 23 corresponds to the face cut of the wax bur shown above. Precisely designed ditches in wax can be milled additionally and easily with the suitable profile bur.

- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm

### Polishing bur, parallel, round face





A mirror-finish surface is prepared with the polishing bur F137 3P 23. The identical face cuts of all wax, profile and polishing burs of the same size simplify designing of a perfect marginal ditch.

- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm






Milling and drilling oil see page 449 REF 550 0000 8




### Wax bur, parallel, straight face





Wax bur F116 3W 23: Smooth and precise wax surfaces due to modern cutting edge geometry. Milling tools with straight face are particularly suitable for tangential borders in the marginal area.

### Profile bur, parallel, straight face






Rapid and precise shaping with the profile bur F116 3H 23: The relief technology allows to obtain excellent cutting performance.

- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm

### Polishing bur, parallel, straight face






The high-lustre polished land allows to prepare high-lustrous milling surfaces. The secondary elements can be attached to these surfaces without subsequent polishing. This way precision is increased and working time is saved.

- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft  $\emptyset$  3.00 short: total length of tool 30 mm



### Parallel burs for titanium, precious and non-precious metal alloys


### Profile bur abrasive, parallel round face



working speed of non-precious metal 20,000 - 25,000 rpm

- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm

# Polishing bur abrasive, parallel round face



- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm

#### Abrasive burs for titanium, precious and non-precious metal and milling work.

The cutting edge geometry has been especially designed for rapid removal of material. The relief ensures smooth running and allows to obtain very smooth surfaces on the object to be milled. The wide relief extends the edge life and avoids breaking of the sharp cutting edges.




Milling and drilling oil see page 449 REF 550 0000 8




### Parallel burs for titanium, precious and non-precious metal alloys

### Profile bur abrasive, parallel straight face



- shaft Ø 2.35 short: total length of tool 34 mm
- \*\* shaft Ø 3.00 short: total length of tool 30 mm

### Polishing bur abrasive, parallel straight face



- \* shaft Ø 2.35 short: total length of tool 34 mm
  - \*\* shaft Ø 3.00 short: total length of tool 30 mm

445

A rich quantity of milling and drilling oil is applied onto the milling surface and a speed of 20,000 – 25,000 rpm is used for milling. Whilst adding a copious quantity of milling and drilling oil, the milling surface is prepolished with the polishing bur at 20,000 rpm and then polished to high lustre at 10,000 rpm.



Milling and drilling oil see page 449 REF 550 0000 8

Fax (+49) 0 73 09 / 8 72-4 44 C Rotating tools

# Wax bur, conical, round face



	Material	Tungsten carbide	Tungsten carbide
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short
0-6-6-6	REF	F200 2W 23	F200 3W 23
conical 2°	ISO No.	500 103 200362 023	500 123 200362 023
	REF	F200 2W 31	F200 3W 31
conical 4°	ISO No.	500 103 200362 031	500 123 200362 031
	REF	F200 2W 40	F200 3W 40
conical 6°	ISO No.	500 103 200362 040	500 123 200362 040
	working spee	d on wax	2,500 - 5,000 rpm



total length: 32 mm

Preparation of a tapered crown with marginal ditch: The wax bur F200 3W 40 features a rounded face. First a wax pattern is prepared with this wax bur. The rounded face features a sharp cutting edge and produces a precise ditch.

### Profile bur, conical, round face



	Material	Tungsten carbide	Tungsten carbide
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short
	REF	F200 2H 23	F200 3H 23
conical 2°	ISO No.	500 103 200103 023	500 123 200103 023
anananan S	REF	F200 2H 31	F200 3H 31
conical 4°	ISO No.	500 103 200103 031	500 123 200103 031
	REF	F200 2H 40	F200 3H 40
conical 6°	ISO No.	500 103 200103 040	500 123 200103 040
	working spe	ed on precious metal	15,000 - 20,000 rpm



The object is shaped with the profile bur after casting. The radius of the face cut has been precisely matched with the corresponding wax bur. Accordingly, the ditch can be shaped quickly.

total length: 32 mm

# Polishing bur, conical, round face



	Material	Tungsten carbide	Tungsten carbide
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short
	REF	F200 2P 23	F200 3P 23
conical 2°	ISO No.	500 103 200133 023	500 123 200133 023
	REF	F200 2P 31	F200 3P 31
conical 4°	ISO No.	500 103 200133 031	500 123 200133 031
	REF	F200 2P 40	F200 3P 40
conical 6°	ISO No.	500 103 200133 040	500 123 200133 040
	working spe	eed on precious metal	15.000 - 20.000 rpm



surface is finished with the polishing bur. Slight reworking with the polishing bur allows to obtain high lustre polishing even in the area of the ditch, since the radiuses of the polishing, profile and wax burs are absolutely identical.

total length: 32 mm



Milling and drilling oil see page 449 REF 550 0000 8



# Wax bur, conical, straight face



	Material	Tungsten carbide	Tungsten carbide
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short
	REF	F186 2W 23	F186 3W 23
conical 2°	ISO No.	500 103 186362 023	500 123 186362 023
	REF	F186 2W 31	F186 3W 31
conical 4°	ISO No.	500 103 186362 031	500 123 186362 031
	REF	F186 2W 40	F186 3W 40
conical 6°	ISO No.	500 103 186362 040	500 123 186362 040
	working sp	eed on wax	2,500 - 5,000 rpm



marginal design: The wax bur F186 3W 40 is perfectly suitable for this task. First the wax pattern is prepared with the wax bur. The sharp cutting edges ensure particularly smooth wax surfaces.

Preparation of a tapered crown with tangential

total length: 32 mm

# Profile bur, conical, straight face



	Material	Tungsten carbide	Tungsten carbide
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short
	REF	F186 2H 23	F186 3H 23
conical 2°	ISO No.	500 103 186103 023	500 123 186103 023
***************************************	REF	F186 2H 31	F186 3H 31
conical 4°	ISO No.	500 103 186103 031	500 123 186103 031
	REF	F186 2H 40	F186 3H 40
conical 6°	ISO No.	500 103 186103 040	500 123 186103 040
	working spee	d on precious metal	15,000 - 20,000 rpm



The friction surfaces of the primary crown are precisely shaped with the profile bur. The optimized cutting angle of the profile bur ensures efficient profile milling.

total length: 32 mm

# Polishing bur, conical, straight face



	Material	Tungsten carbide	Tungsten carbide
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short
	REF	F186 2P 23	F186 3P 23
conical 2°	ISO No.	500 103 186133 023	500 123 186133 023
	REF	F186 2P 31	F186 3P 31
conical 4°	ISO No.	500 103 186133 031	500 123 186133 031
	REF	F186 2P 40	F186 3P 40
conical 6°	ISO No.	500 103 186133 040	500 123 186133 040
	working spee	d on precious metal	15,000 - 20,000 rpm



After profile milling, a mirror-like high lustre is produced on the friction surface using the polishing bur. Additional polishing of the surface is no longer required.

C Rotating tools

total length: 32 mm



Fax (+49) 0 73 09 / 8 72-4 44

Milling and drilling oil see page 449 REF 550 0000 8

### Conical burs for titanium, precious and non-precious metal alloys

### Profile bur, abrasive, conical, round face



	Material	Tungsten carbide	Tungsten carbide	
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short	
	REF	N200 2H 23	N200 3H 23	
conical 2°	ISO No.	500 103 200 103 023	500 123 200 103 023	
	REF	N200 2H 31	N200 3H 31	
conical 4°	ISO No.	500 103 200 103 031	500 a123 200 103 031	
	REF	N200 2H 40	N200 3H 40	
conical 6°	ISO No.	500 103 200 103 040	500 123 200 103 040	
	working spee	d of non-precious metal	20,000 - 25,000 rpm	

working speed of non-precious metal 20,000 - 25,000 rpm total length 32 mm

Abrasive burs for titanium, precious and non-precious metal and milling work.

The cutting edge geometry has been especially designed for rapid removal of material. The relief ensures smooth running and allows to obtain very smooth surfaces on the object to be milled. The wide relief extends the edge life and avoids breaking of the sharp cutting edges.

### Profile bur, abrasive, conical, straight face



	Material Tungsten carbide		Tungsten carbide	
	Shaft	Ø 2.35 mm short	Ø 3.00 mm short	
	REF	N186 2H 23	N186 3H 23	
conical 2°	ISO No.	500 103 186 103 023	500 123 186 103 023	
	REF	N186 2H 31	N186 3H 31	
conical 4°	ISO No.	500 103 186 103 031	500 123 186 103 031	
	REF	N186 2H 40	N186 3H 40	
conical 6°	ISO No.	500 103 186 103 040	500 123 186 103 040	
	working spee	d of non-precious metal 1	20 000 = 25 000 rpm	


working speed of non-precious metal 20,000 - 25,000 rpm total length 32 mm

A rich quantity of milling and drilling oil is applied onto the milling surface and a speed of 20,000 - 25,000 rpm is used for milling. Whilst adding a copious quantity of milling and drilling oil, the milling surface is prepolished with the polishing bur at 20,000 rpm and then polished to high lustre at 10,000 rpm.



Milling and drilling oil see page 449 REF 550 0000 8

### Groove bur





The groove bur F538 3H 10 exclusively serves to prepare grooves. For this purpose the groove bur is only moved up and down in the vertical axis. The cut on the face simplifies extending the grooves to the cervical direction. At speeds of 15.000 – 20,000 rpm a fine cutting performance and a mirror-like high lustre on the milling surface are obtained.

### Shoulder bur



Material	Tungsten carbide	Tungsten carbide
Shaft	Ø 2.35 mm short	Ø 3.00 mm short
REF	F205 2H 27	F205 3H 27
ISO No.	500 103 205175 027	500 123 205175 027
REF	F205 2H 29	F205 3H 29
ISO No.	500 103 205175 029	500 123 205175 029
working spec	ed on precious metal	15,000 - 20,000 rpm

total length: 34 mm



The occlusal shoulder is prepared with a special shoulder bur. The shoulder bur features a face cut which smoothens the bottom of the shoulder. Optimum use of this tool is ensured at speeds of approx. 15.000 – 20,000 rpm. A mirror-like lustre on the milling surface is achieved with the relief cut. Additional polishing is not required.

### Milling and drilling oil



Milling and drilling oil REF 550 0000 8

The milling and drilling oil was especially developed to be used with the milling and drilling tools by bredent. The special consistency produces a reliable oil film between the metal and the drill so that the metal swarf slides out of the cut sections of the tool. This way the cutting performance and the service life of the milling tools are increased. Gumming of the milling and drilling oil is excluded thanks to the high evaporation temperature.

#### Use:

Always add sufficient quantities of milling and drilling oil during centring, drilling resp. milling.

# Wax burs 0°, 1°, 2°, 4°, 6°





Material	Tungsten ca	Tungsten carbide		
Shaft	Ø 2.35 mm short			
	REF	F137 2W 10		
parallel 0°	ISO No.	500 103 137366 010		
	REF	F137 2W 15		
parallel 0°	ISO No.	500 103 137366 015		
	REF	F137 2W 23		
parallel 0°	ISO No.	500 103 137366 023		
200	REF	F200 2W 29		
konisch 1°	ISO No.	500 103 200362 029		
8-6-6-	REF	F200 2W 23		
conical 2°	ISO No.	500 103 200362 023		
8-6-6-6	REF	F200 2W 31		
conical 4°	ISO No.	500 103 200362 031		
	REF	F200 2W 40		
conical 6°	ISO No.	500 103 200362 040		
working speed on wax		2,500 - 5,000 rpm		

The new 1° wax burs are perfectly suitable for primary conical crowns with a slightly conical angle which are to provide maximum friction below a 0° telescopic crown.

Perfectly suited for electroplating.

### Cross-cut burs





Thanks to the cross cut, high abrasive capacity is achieved to allow quick processing of precious and non-precious metals and titanium.



Material	Tungste	n carbide
Shaft	Ø 2.35 r	nm short
	REF	F 200 2K 29
conical 1°	ISO No.	500 103 200190 029
en e	REF	F 200 2K 23
conical	ISO No.	500 103 200190 023
	REF	F 200 2K 31
conical 4°	ISO No.	500 103 200190 031
	REF	F 200 2K 40
conical 6°	ISO No.	500 103 200190 040
working speed on preciou	ıs metal	20,000 rpm

When using milling and drilling oil, the milled surfaces can be shaped with the same bur and using lower pressure just like with a relief bur. Smooth surfaces can be achieved with the same bur without changing the bur.




Milling and drilling oil see page 449 REF 550 0000 8



### Wax burs and diamond grinding tools, parallel


### Wax bur, straight cutting edge, parallel, round face



Assortment 4 pieces, 1 piece each REF 320 0087 0 All tools feature a total tool length of 36 mm and a shaft diameter of 2.35 mm.

Working speed on wax 5,000 rpm

### Diamond grinding tool, coarse grain, parallel, round face



All tools feature a total tool length of 36 mm and a shaft diameter of 2.35 mm.

The radius of the face of the parallel diamond grinding tool has been matched with the face of the parallel wax burs shown above. The use of diamond grinding tools ensures efficient surface working especially on hard alloys.

Working speed on CrCo/NPM 10,000 - 20,000 rpm

### Diamond grinding tool, fine grain, parallel, round face




If correct usage is ensured, the diamond grinding tools offer high dimensional stability, functionality and edge life thanks to galvanic diamond coating All tools feature a total tool length of 36 mm and a shaft diameter of 2.35 mm.

Working speed on CrCo/NPM 10,000 - 20,000 rpm

### Wax burs and diamond grinding tools, conical

### Wax bur, straight cutting edge, conical, round face



	REF	320 0080 2
conical 2°	ISO No.	330 103 200382 023
	REF	320 0081 4
conical 4°	ISO No.	330 103 200382 031
	REF	320 0082 6
conical 6°	ISO No.	330 103 200382 040

working speed on wax 5,000 rpm

Assortment
3 pieces, 1 piece each
REF 320 0086 0

All tools feature a total tool length of 36 mm and a shaft diameter of 2.35 mm.

### Diamond grinding tool, coarse grain, conical, round face




	REF	340 0088 G
conical 2°	ISO No.	806 103 200534 023
	REF	340 0089 G
conical 4°	ISO No.	806 103 200534 031
and the control of th	REF	340 0090 G
conical 6°	ISO No.	806 103 200534 040

Working speed on CrCo/NPM 10,000 - 20,000 rpm All tools feature a total tool length of 36 mm and a shaft diameter of 2.35 mm.

The radius of the face of the parallel diamond grinding tool has been matched with the face of the parallel wax burs shown above. The use of diamond grinding tools ensures efficient surface preparation especially on hard alloys.

### Diamond grinding tool, fine grain, conical, round face





Working speed on CrCo/NPM 10,000 - 20,000 rpm If correct usage is ensured, the diamond grinding tools offer high dimensional stability, functionality and edge life thanks to galvanic diamond coating.

All tools feature a total tool length of 36 mm and a shaft diameter of 2.35 mm.

453

#### Giflex-TR



Time-saving through more rapid and more precise separation of the saw models than is possible with a hand saw.

### Giflex-TR Master x-tray



Special diamond disc for processing acrylics.

### Diamond grinding wheels



Small diameters, for reliable and specific grinding.

### Ceraflex



Diagonal toothing and abrasive diamond for a high cutting efficiency.

### Microflex



The thinnest diamond grinding wheel coated on both sides – only 0.08 mm thick.

### Transflex-T



The highly flexible grinding wheel with transparency for safe, concerted grinding.

### Transflex



Diagonally arranged cutouts for running transparency with high breaking strength and optimal grinding capacity.

### Ultraflex, Superflex, Flexibel, Elastisch



Special diamond graining in various thicknesses and coatings – the right diamond grinding wheel for every area of application.

Fax (+49) 0 73 09 / 8 72-4 44 C Rotating tools

# Diamond grinding wheels

### Giflex-TR



Time-saving through more rapid and more precise separation of the saw models than is possible with a hand saw.





Due to the wide application range of the Giflex-TR it is possible to the separate the arch from the basal direction if preparation margins are close to the each other. The segmented design of the Giflex-TR allows rapid removal of grinding dust so that jamming of the disc is avoided.

### Giflex-TR Master x-tray



Thickness:

Coating:
on both sides
Version:

Ø 25 mm

REF

340 00M2 5

Special diamond disc for processing acrylics. Giflex-TR Master x-tray features a coarse diamond grit; accordingly, a cooling effect is achieved already in the diamond-coated area when separating acrylics.

### Diamond mini



Small diameters, for reliable and specific grinding.

Thickness: Coating: Version:	0.23 mm on both sides mounted
Ø 8 mm	<del></del> -
REF	340 0014 3
Ø 10 mm	E-Q-3
REF	340 0014 4
Ø 12 mm	
REF	340 0014 5
Ø 14 mm	
REF	340 0014 6



The small diameter is particularly suitable for finishing interdental spaces of ceramically veneered bridges.

#### Ceraflex



Diagonal toothing and abrasive diamond for a high cutting efficiency. The cooling effect of the saw toothing with the abrasive diamond graining creates ideal conditions for rapid, concerted finishing of synthetic and ceramic veneers.

Thickness: Coating: Version:	<b>0.25 mm</b> on both sides mounted
Ø 16 mm	
REF	340 0013 0
Ø 22 mm	
REF	340 0003 0



As a result of its cooling effect, Ceraflex is particularly well suited for finishing composites and other synthetic materials.

#### Microflex



The thinnest diamond grinding wheel coated on both sides - only 0.08 mm thick. The flexibility and thin structure of Microflex allow very fine separation in the anterior and posterior

Thickness:	0.08 mm
Coating:	on both sides
Version:	mounted
Ø 12 mm	
REF	340 0014 2
Ø 16 mm	
Ø 16 mm	340 0014 1
2	П

Very thin incisions create natural vividness in the veneering technique.

### Transflex-T



The longish cutouts in the grinding wheel ensure cooling, high flexibility and maximum running transpar-Thickness: 0.20 mm 0.20 mm 0.25 mm ency. This permits on both sides Coating: on one side on one side working in a way Version: outer coating inner coating mounted which is gentle on materials, clear and Ø 16 mm targeted and allows REF 340 0010 0 Ø 22 mm

340 0008 0



The special arrangement of the longish cutouts makes the rotating grinding wheel transparent. Being able to see the area of application during grinding is a considerable advantage.

### **Transflex**



Diagonally arranged cutouts for running transparency with high breaking strength and optimal grinding capacity. Transflex is particularly well suited for diving and separating in the areas of the front and side teeth and in approximate shaping.

to save time.

REF

Thickness: Coating: Version:	<b>0.20 mm</b> on one side outer coating	<b>0.20 mm</b> on one side inner coating	<b>0.25 mm</b> on both sides mounted
Ø 22 mm			
REF	340 0005 0	340 0006 0	340 0004 0

340 0009 0

340 0007 0



The specially diagonally arraned cutouts ensure running transparency with high stability and abrasiveness of the grinding wheel.

### Diamond grinding wheels

### Ultraflex, Superflex, Flexibel, Elastisch

Special diamond graining in various thicknesses and coatings - the right diamond grinding wheel for every area of application.













#### **Ultraflex**




#### Superflex



#### **Flexibel**



#### Elastisch



### Diacryl Grinding Instruments dcs



Save time and improve quality by grinding acrylic with diamond coated Diacryl rotating instruments.

# Diagen-Turbo-Grinder dtg



The diamond grinder system with the extraordinary grinding properties due to special Diagen diamond binding material.

### Set-up grinding tool



Two grinding tools in one. Grinding without exchanging tools in a single working step.

### Special Diamonds for the Veneering Technique



Perfect finishing of acrylic and ceramic veneers.

# Diamond grinding tool dsl



The all-rounder among diamond grinding tools, available in the most common shapes.

### FG-Diabolo



Economic system of grinding tools with razor sharp diamonds, self-regenerating grit and extended durability.

Fax (+49) 0 73 09 / 8 72-4 44

457

### **Diamond grinding**

### Diacryl Grinding Instruments dcs



Save time and improve quality by grinding acrylic with diamond coated Diacryl rotary instruments. Thanks to their uniform, coarse grit diamond particles with sharp edges and their specific shapes, Diacryl diamond instruments are excellent for trimming acrylic dentures quickly and accurately.



Universal diamond instrument REF 340 0104 0



Papillae diamond instrument REF 340 0105 0

Coarse diamond instrument

REF 340 0103 0



Round diamond instrument for peripheries REF 340 0106 0



Round diamond instrument for peripheries REF 340 0102 0



Rubber grinder REF 340 0090 0





























Thanks to their extra coarse diamond grit and large diameter, these instruments grind aggressively and create a perfectly ground surface. They are perfect for reducing large areas of acrylic. Their hollow shape enables them to be used at high speeds, with maximum cooling effect.

Can be used in lingual and palatal areas, as required. For grinding large papillae and root attachment or lingual bars.


The pointed flame shape permits the papillae and alveolar attachments to be ground interdentally.

Thanks to the tapered central section of the instrument, uniformly thick peripheries can be created easily and quickly on functional impression trays and partial dentures. Acrylic beads and rough areas on the fitting surface of the denture are easily removed with the round head of the instrument.

Narrow frenae can be rounded and perfected with this Diacryl instrument.

This abrasive rubber instrument creates smooth surfaces on acrylic dentures instead of using sandpaper. Thanks to the smooth surfaces, the denture can be prepared for optimum polishing in the shortest possible time.

### Diagen-Turbo-Grinder, fine



The diamond grinder system with the extraordinary grinding properties due to special Diagen diamond binding material.

Maximum grinding power and abrasive capacity on metal and ceramic surfaces at reduced pressure. Increased service life compared to conventional binding material allows a wide range of applications and thus high efficiency.



Cylinder, pointed Ø 3.5 x 11 mm, 2 pcs REF 340 0015 5



Cone Ø 3.5 x 11 mm 2 pieces REF 340 0015 0

Cylinder

2 pieces



Cylinder, round Ø 4.5 x 13 mm, 2 pcs REF 340 0016 5



Cylinder Ø 4.8 x 13 mm 2 pieces REF 340 0016 0



Inverted cone with recess Ø 6 x 8 mm 1 piece REF 340 0025 0



Inverted cone with recess Ø 12 x 6 mm 1 piece REF 340 0024 0



Lens Ø 22 x 2 mm 1 piece REF 340 0021 0



Disc Ø 22 x 2 mm 1 piece REF 340 0022 0



Disc Ø 15 x 3,5 mm, 2 pieces REF 340 0018 0



Disc Ø 22 x 4,5 mm, 1 piece REF 340 0019 0



Assortment 5 pieces, 1 piece each Diagen-Turbo-Grinder REF 340 0020 0



Non-precious and precious metal alloys The dtg grinders are perfectly suitable for processing of metal



Assortment 5 pieces, 1 piece each Diagen-Turbo-Grinder ceramic REF 340 0020 5



Non-precious and precious metal alloys Gentle and pressure-free grinding allows perfect finishing and efficient removal of material.



surfaces of ceramic frameworks thanks to the diamond coating.



Ceramic and resins The grinder system is perfectly suitable for processing of ceramic due to grinding at reduced heat and extensive removal of material.

#### Ideal for processing zirconium oxide - recommended by leading implant manufacturers



The cylindrical shape is perfectly suitable for parallel processing.



The inverted cone with recess provides a cooling effect and thus the perfect precondition for processing zirconium oxide.



Thanks to their diamond coating, Diagen-Turbo-Grinders represent the ideal tools for processing materials such as glass ceramic or zirconium oxide



Rapid and convenient shaping of implant suprastructures. The fine diamond particles produce a smooth surface.




The variety of different shapes allow to obtain a wide indication range.



Extremely hard alloys and be processed as easily

### **Diamond grinding**

### Diagen-Turbo-Grinder, coarse



\*The companies Ivoclar Vivadent, Vita und Teamziereis recommend the Diagen-Turbo-Grinders for processing zirconium and ceramic. Diagen-Turbo-Grinder, coarse – new quality of the Diagen binding with coarse diamond grains for quick preparation.

- two abrasion levels for more flexibility in the usage
- special binding material (coarse) for longer tool life (20 %) and reduced consumption of grinders
- cooling effect avoids damage to the zirconium
- 11 different shapes for all applications ensure perfect processing results
- can be used for zirconium, ceramic and metal
- continuously exposed diamonds guarantee fast removal of material

- recommended by leading zirconium and ceramic manufacturers; \*
- to make sure that the right decision has been made \*

#### New shapes:

- Cylinder, pointed, for easy accesss to narrow areas
- Cylinder, round, for finishing chamfers

#### Comparison of microsection Zirconium



New shapes in additional abrasion level. Fast removal thanks to coarse diamond grain and special diamond binding material. Surface roughness: 26 µm



New shapes in proven quality. Smoothening of surfaces and finishing with fine diamond grain: Surface roughness 2 µm



**Cylinder, pointed** Ø 3.5 x 11 mm 2 pcs **REF 340 G015 5** 



**Cylinder, round** Ø 4.5 x 13 mm 2 pcs **REF 340 G016 5** 



Cylinder Ø 4.5 x 13 mm 2 pcs REF 340 G016 0



**Cylinder** Ø 6.0 x 13 mm 2 pcs **REF 340 G017 0** 



Cone Ø 3.5 x 11 mm 2 pcs REF 340 G015 0



Cone Ø 3.5 x 11 mm 2 pcs REF 340 G015 0



Quick finishing of closely positioned crowns with the pointed cylinder saves time.



Processing without pressure and simultaneous cooling due to the diagen binding material avoid damage to the zirconium.



The round shape facilitates finishing of chamfers for individual abutments. This way the time for finishing is re-

duced and the efficiency is increased.



The unique diamond binding material ensures longer tool life and the coarser diamond grain facilitates efficient processing.



The coarse diamond grain leads to fast removal of material and reworking with the fine Diagen-Turbo-Grinder produces perfect surfaces. Ideal for processing IPS e.max.



Accurate finishing in the shoulder area with the fine Diagen-Turbo-Grinder for precise restorations made of press ceramic, layering ceramic or zirconium. A perfect surface for subsequent polishing with Zi-polish is obtained.

#### Accessories



**Zi-polish** 5 g **REF 360 1002 5** 



The abrasive polishing paste Zi-polish in combination with the Rodeo round brush produces a perfect high luster on zirconium surfaces in a shorter time.



### Set-up grinding tool



Set-up grinding tool 1 piece REF 340 0101 0

Two grinding tools in one. Grinding without exchanging tools in a single working step.

- quick adaptation of the underside of the tooth to be set up
- grinding in of occlusal stops

#### Two grinding tools in one



#### occlusal



The small, precisely shaped grinding tip with fine, perfectly cutting diamond grains provides the ideal precondition for well-aimed and rapid grinding in of occlusal contacts.

basal



The large grinding area with its optimized shape and selected natural abrasive diamonds ensures maximum removal of material and thus accurate and quick grinding.


### Special Diamonds for the Veneering Technique



Diamond grinding tool for veneering techniques Vb 1
1 piece

REF 340 0084 0 ISO No. 806 104 033524 029

Diamond grinding tool for veneering techniques Vb 2
1 piece
REF 340 0083 0
ISO No. 806 104 000524 032



Diamond grinding tool for veneering techniques Vb 3
1 piece
REF 340 0085 0
ISO No. 806 104 171524 033

Perfect finishing of acrylic and ceramic veneers.



The concave part of this tool ensures perfect contouring of the approximal surfaces in the area of the neck of the tooth. Enhanced aesthetics in reduced



The very thin tip of this tool allows to design tooth necks of single crowns and bridges in an efficient way.



In a very simple way this tool creates uniform enamel bulges at the tooth neck. All tools can be used on ceramic and acrylic veneering materials.



The conical grinding tool with coated face provides a universal application range. The entire process of contouring can be carried out with this tool.

# Diamond grinding tool dsl

	Designation	REF	Supply form	Order quantity
	KA 4,0	340 0070 0	1 piece	
	KA 2,5	340 0071 0	1 piece	
	KI 2,5	340 0072 0	1 piece	
CALCULATE TO A STATE OF THE STA	SR 1,0	340 0073 0	1 piece	
	KF 0,5	340 0074 0	1 piece	
	KS 2,0	340 0075 0	1 piece	
90	KS 1,5	340 0076 0	1 piece	
8	RU 2,0	340 0077 0	1 piece	
	RU 1,5	340 0078 0	1 piece	
8	RU 1,0	340 0079 0	1 piece	
<b></b>	LZ 2,0	340 0080 0	1 piece	
1	LZ 1,5	340 0081 0	1 piece	

Sender (stamp):	Customer No.	Additional order
	Date, signature	

### Diabolo



Economic system of grinding tools with razor sharp diamonds, self-regenerating grit and extended durability.



Thanks to embedding the diamond crystals in a special mixture of binding material ceramic surfaces can be processed without leaving any residues.



The variety of shapes of the Diabolo tools offers an individual range for each surface type and all dental materials.



The outstanding hardness of the sintered diamonds allows particularly efficient use and material removal for resistant materials through self-sharpening diamond crystals.




The Diabolo sintered diamond tools ensure reliable and efficient processing of sensitive areas down to the last diamond grain.



Various materials can be processed with a Diabolo diamond grinding tool without time-consuming, repeated exchanging of tools.

### Sintered diamond grinding tools



# Schematic comparison of electroplated and sintered diamond grinding tools:

In the case of electroplated grinding tools, the diamond crystals have been attached to the bur blank in a metal bond.



Diamond crystals in an electrodeposited metal bond.

In the bredent sintering method the razor sharp diamond crystals are embedded in an adapted mixture of binding material.



Diamond crystals embedded in a special mixture of binding material.


# The superior class of diamond grinding tools. For fast and efficient processing of extremely hard dental materials.

Carefully selected natural diamonds are entirely integrated into a mixture of metal and binding material designed for the individual application.

Due to the special manufacturing process worne out diamond grit is automatically removed and replaced by razor sharp diamond crystals so that automatic sharpening is achieved.

The range of different shapes ensures individual selection and permits the use for numerous applications.

The self-sharpening effect allows to perform highly difficult processing of dental material in a simple, fast and efficient way.



Summary: When comparing sintered/electroplated diamond grinding tools in the long-term test, the bredent Diabolo tools excelled by their efficiency and extended service life. Due to the self-sharpening effect of the diamond crystals, the sintered diamond tools feature high cutting performance and ensure perfect removal of material until the end of service life.

bredent §

### Diamond grinding

#### Diabolo

#### Color coding

Diabolo grinding tools feature color coding.


This system indicates the different grit size of the grinding tool and thus simplifies the selection of the suitable tool.

Color coding	Grit size	Grit area	Marking
	200 μm		extra coarse / black
	130 μm		coarse/ green
	100 μm	=	normal
	80 μm		fine / red

### Work recommendations

Always adapt speed and pressure (approx. 0.3 - 2N) to the material to be processed.

#### Recommended speed





#### Order Number


To simplify reordering, the order number of the respective Diabolo grinding tool is engraved on the shaft.

#### The color code: Order number: Diabolo diamond crystals constantly form new cutting The order number is included on the shaft of From fine to extra coarse grit - at a single glance! edges during grinding. This way extremely high A separate color for each of the four grit sizes ensures every tool to exclude errors when ordering. resistance and extended service life are ensured. that you select the correct Diabolo. Simple selection of the desired grit size with the help of the color rings. SF = Sintered diamond, fine 199 = shape of the working tip (C) 050 = largest diameter bredent SF 199 050 of the working tip (E) Precise: Guaranteed bredent quality: Very easily changed: Every Diabolo sintered diamond is absolutely Every sintered diamond undergoes the Diabolo sintered diamonds feature rounded concentric and therefore wears down evenly. bredent quality assurance test. We guarantee ends on their shafts so that they can be Accordingly, restorations fit precisely. This even optimum, uniform cutting performance right inserted quickly into any handpiece. applies to complex milling of non-precious metal down to the last particle of diamond. objects. ISO number C Ε ISO numbers are indicated for all tools 7. - 9. digit: 13. - 15. digit: 1. - 3. digit: to ensure enhanced comparability. Material of the Shape of the Diameter of the These internationally standardized numbers feature working tip working tip working tip 15 digits. The numbers include the following information:

#### Diagram Grinding tool







# Diabolo


Conical, pointed		extra coarse 200 μm	coarse 130 μm	normal 100 μm	fine 80 μm	
	REF	SS 165 023	SG 165 023	SN 165 023	SF 165 023	
A	ISO No.	807 104 165543 023	807 104 165533 023	807 104 165523 023	807 104 165513 023	10
	REF	SS 167 050	SG 167 050	SN 167 050	SF 167 050	
	ISO No.	807 104 167543 050	807 104 167533 050	807 104 167523 050	807 104 167513 050	
	REF	SS 213 016	SG 213 016	SN 213 016	SF 213 016	
7	ISO No.	807 104 213543 016	807 104 213533 016	807 104 213523 016	807 104 213513 016	Highly accurate shaping
	REF	SS 213 031	SG 213 031	SN 213 031	SF 213 031	of occulusal surfaces of ceramic veneers.
	ISO No.	807 104 213543 031	807 104 213533 031	807 104 213523 031	807 104 213513 031	ccianne veneers.
Conical, round		extra coarse 200 μm	n coarse 130 μm	normal 100 μm	fine 80 μm	
	REF	SS 400 005	SG 198 025	CN 100 025	CF 100 025	3
	ISO No.	SS 198 025 807 104 198543 025		SN 198 025	SF 198 025	0
			807 104 198533 025	807 104 198523 025	807 104 198513 025	1 200
5 3	REF	SS 198 037	SG 198 037	SN 198 037	SF 198 037	
	ISO No.	807 104 198543 037	807 104 198533 037	807 104 198523 037	807 104 198513 037	
	REF	SS 199 031	SG 199 031	SN 199 031	SF 199 031	Grinding the incide of
	ISO No.	807 104 199543 031	807 104 199533 031	807 104 199523 031	807 104 199513 031	Grinding the inside of chrome cobalt clasps.
	REF	SS 199 040	SG 199 040	SN 199 040	SF 199 040	
	ISO No.	807 104 199543 040	807 104 199533 040	807 104 199523 040	807 104 199513 040	
	REF	SS 199 050	SG 199 050	SN 199 050	SF 199 050	
	ISO No.	807 104 199543 050	807 104 199533 050	807 104 199523 050	807 104 199513 050	
Conical		extra coarse 200 μm	coarse 130 μm	normal 100 μm	fine 80 μm	
	REF	SS 171 016	SG 171 016	SN 171 016	SF 171 016	
	ISO No.	807 104 171543 016	807 104 171533 016	807 104 171523 016	807 104 171513 016	
<b>4</b> 1	REF ISO No.	SS 171 031	SG 171 031	SN 171 031	SF 171 031	and the same of th
		807 104 171543 031	807 104 171533 031	807 104 171523 031	807 104 171513 031	
	REF ISO No.	SS 172 031	SG 172 031	SN 172 031	SF 172 031	Removing sharp edges
		807 104 172543 031	807 104 172533 031	807 104 172523 031	807 104 172513 031	and flattening of rete
	REF	SS 173 031	SG 173 031	SN 173 031	SF 173 031	tion beads.
	ISO No.	807 104 173543 031	807 104 173533 031	807 104 173523 031	807 104 173513 031	
	REF	SS 173 040	SG 173 040	SN 173 040	SF 173 040	
	ISO No.	807 104 173543 040	807 104 173533 040	807 104 173523 040	807 104 173513 040	
	REF	SS 173 050	SG 173 050	SN 173 050	SF 173 050	
	─ ISO No.	807 104 173543 050	807 104 173533 050	807 104 173523 050	807 104 173513 050	
ud, round		extra coarse 200 μm	coarse 130 μm	normal 100 μm	fine 80 μm	Y
	DEE	SS 201 050	SG 261 050	SN 261 050	SE 261 050	\$
and the same of th	REF ISO No.	SS 261 050	SG 261 050	SN 261 050	SF 261 050	
2 %		807 104 261543 050	807 104 261533 050	807 104 261523 050	807 104 261513 050	
	REF	SS 263 050	SG 263 050	SN 263 050	SF 263 050	
	ISO No.	807 104 263543 050	807 104 263533 050	807 104 263523 050	807 104 263513 050	Smoothing of the sublingual bar.
ud		extra coarse 200 μm	n coarse 130 μm	normal 100 μm	fine 80 μm	
						3
<b>A</b>	REF	SS 254 060	SG 254 060	SN 254 060	SF 254 060	
	ISO No.	807 104 254543 060	807 104 254533 060	807 104 254523 060	807 104 254513 060	



Grinding the retention area of chrome cobalt objects.

### **Diamond grinding**

#### Diabolo





Grinding the surface of all-ceramic frameworks and veneering ceramic.






Opening the interdental spaces.

	REF	SS 141 025	SG 141 025	SN 141 025	SF 141 025
	ISO No.	807 104 141543 025	807 104 141533 025	807 104 141523 025	807 104 141513 025
	REF	SS 141 031	SG 141 031	SN 141 031	SF 141 031
	ISO No.	807 104 141543 031	807 104 141533 031	807 104 141523 031	807 104 141513 031
	REF	SS 143 050	SG 143 050	SN 143 050	SF 143 050
	ISO No.	807 104 143543 050	807 104 143533 050	807 104 143523 050	807 104 143513 050
	REF	SS 143 080	SG 143 080	SN 143 080	SF 143 080
	ISO No.	807 104 143543 080	807 104 143533 080	807 104 143523 080	807 104 143513 080
	REF	SS 153 016	SG 153 016	SN 153 016	SF 153 016
	ISO No.	807 104 153543 016	807 104 153533 016	807 104 153523 016	807 104 153513 016
	REF	SS 153 031	SG 153 031	SN 153 031	SF 153 031
	ISO No.	807 104 153543 031	807 104 153533 031	807 104 153523 031	807 104 153513 031



Cervical finishing of veneering ceramic.





Smoothing the surface structure in the incisal area of ceramic veneers.

# Diamond grinding / dressing stone for diamond grinding tools

### Diabolo



REF ISO No. REF ISO No.




Smoothing of surfaces treated with a solder



extra coarse 200 μm coarse 130 μm		normal 100 μm	fine 80 μm
SS 227 016	SG 227 016	SN 227 016	SF 227 016
807 104 227543 016	807 104 227533 016	807 104 227523 016	807 104 227513 016
SS 227 031	SG 227 031	SN 227 031	SF 227 031
807 104 227543 031	807 104 227533 031	807 104 227523 031	807 104 227513 031



Finishing of deeper surfaces of metal frameworks



	extra coarse 200 μm coarse 130 μm		normal 100 μm	fine 80 μm
REF	SS 233 016	SG 233 016	SN 233 016	SF 233 016
ISO No.	807 104 233543 016	807 104 233533 016	807 104 233523 016	807 104 233513 016
REF	SS 233 031	SG 233 031	SN 233 031	SF 233 031
ISO No.	807 104 233543 031	807 104 233533 031	807 104 233523 031	807 104 233513 031



Smoothing of chrome cobalt areas difficult to access.



	extra coarse 200 μm	coarse 130 µm	normal 100 μm	fine 80 μm
REF	SS 014 018	SG 014 018	SN 014 018	SF 014 018
ISO No.	807 104 014543 018	807 104 014533 018	807 104 014523 018	807 104 014513 018
REF	SS 014 021	SG 014 021	SN 014 021	SF 014 021
ISO No.	807 104 014543 021	807 104 014533 021	807 104 014523 021	807 104 014513 021
REF	SS 014 050	SG 014 050	SN 014 050	SF 014 050
ISO No.	807 104 014543 050	807 104 014533 050	807 104 014523 050	807 104 014513 050
REF	SS 014 080	SG 014 080	SN 014 080	SF 014 080
ISO No.	807 104 014543 080	807 104 014533 080	807 104 014523 080	807 104 014513 080



Finishing of a precise marginal edge.

### Diabolo Cleaner





Diabolo Cleaner 1 piece REF 340 0100 0


The indispensable tool for removing contaminations so that constant cutting performance of the Diabolo tools is ensured.

Contaminated material is removed easily and quickly and new diamond cutting edges are exposed from the bronze binding material.

bredent §

## Diamond grinding

### Diabolo



#### Note

The photos illustrate only  $\underline{\text{one}}$  common application field of the grinding tool.



### Diabolo

#### Discs





Separating and grinding back the sprues.

	extra coarse 200 μm coarse 130 μm		normal 100 μm	fine 80 μm
REF	SS 327 004	SG 327 004	SN 327 004	SF 327 004
ISO No.	807 104 327543 004	807 104 327533 004	807 104 327523 004	807 104 327513 004
Ø 25 x 0,4 mm				
REF	SS 327 010	SG 327 010	SN 327 010	SF 327 010
ISO No.	807 104 327543 010	807 104 327533 010	807 104 327523 010	807 104 327513 010
Ø 25 x 1,0 mm				
REF	SS 327 020	SG 327 020	SN 327 020	SF 327 020
ISO No.	807 104 327543 020	807 104 327533 020	807 104 327523 020	807 104 327513 020
Ø 25 x 2,8 mm				
REF	SS 327 030	SG 327 030	SN 327 030	SF 327 030
ISO No.	807 104 327543 030	807 104 327533 030	807 104 327523 030	807 104 327513 030
Ø 25 x 3,0 mm				
REF	SS 327 080	SG 327 080	SN 327 080	SF 327 080
ISO No. Ø 7 x 0,4 mm	807 104 327543 080	807 104 327533 080	807 104 327523 080	807 104 327513 080

### Disc



Supra Disc



Duo Disc




Supra Disc





Highly accurate separation and preparation of interdental areas.




## **Diamond grinding**

### Diabolo

# Fast, efficient processing of any type of zirconium oxide and hard dental materials.

Friction grip – first-class diamond grinders. Economic system of grinding tools with razor sharp diamonds, self-regenerating grit and extended durability. FG-Diabolo are sintered diamond grinding tools and are used in the turbine or in the T-hand-piece with FG adapter for grinding zirconium oxide and extremely hard materials. FG-Diabolo reduces processing times and renders working processes highly efficient.




FG adapter 1.6 to 2.35 for the use in the handpiece Pack cont. 10 pieces REF 340 0100 2

not included in the assortment

#### Assortment

6 pieces, 1 piece each FG-Diabolo, fine grit Bud, large Flame Cylinder, round head Inverted cone Torpedo Bud, small

REF 330 0116 6



Bud, fine, large REF FF 263 023



Flame, fine REF FF 250 016







Cylinder, round head fine REF FF 141 023






Inverted cone, fine REF FF 227 023





Torpedo, fine REF FF 289 023





Bud, fine, small REF FF 263 014



#### All the tools you need even when processing zirconium oxide

Indispensable tool for removing contaminations. Diabolo Cleaner guarantees constant cutting performance. Contaminated material is removed easily and quickly and new diamond crystals are exposed from the bronze binding material and integrated into the surface of the tool to enhance the cutting performance and reduce the grinding time.

Leading implant manufacturers recommend Diagen-Turbo-Grinders dtg. These abrasive tools have proved their suitability for reworking sintered zirconium frameworks and minimized the amount of work thanks to the special Diagen diamond binding material.



see also page 284

Sortiment 5 pieces, 1 piece each Diagen-Turbo-Grinder REF 340 0020 0





Diabolo Cleaner 1 piece REF 340 0100 0



airaqua turbine REF 110 0146 0



# Polishing brushes / Polishing pastes

D

Specially selected raw materials guarantee perfect prepolishing and a brilliant high luster on all dental materials.

The wide range of shapes and materials offer new possibilities for polishing with dental motors and dental handpieces.





Polishing metal	
Survey	472
Standard polishing brushes	473
Abraso-Soft Metal	474
Abraso-Buff Metal	474
Abraso-Buff Metal Mini	474
Abraso-Buff Polipast Metal	475
High Luster Buff Metal	
The metal polishing set	476
Prepolishing with a handpiece	477
High luster polishing with a handpiece	478
Cotton mandrel	478
Polishing acrylic	
Survey	
Standard polishing brushes	480
Abraso-Sil Acrylic	480
Abraso-Buff Acrylic	481
Prepolishing Buff Silicone	481
Abraso-Soft Acrylic	481
Leather buff	481
High luster Buff Acrylic	482
Acrylic finishing set	482
The acrylic polishing set	483
Prepolishing with a handpiece	484
High luster polishing with a handpiece	484
Universal polishing / Polishing porcelain	
Abraso-Fix	485
Polishing porcelain	485
Polishing pastes	
Metal polishing pastes	
Metal and acrylic polishing pastes	487
Acrylic polishing paste	
Porcelain polishing paste	487
Cleaning	
Dentaclean mixing fluid for pumice powder	488
Instruments	
Pollvarip	488

### Polishing metal - Survey

### Standard polishing brushes



The standard line of polishing brushes includes a wide range of different shapes and sizes for various applications when polishing all types of metal.

### Abraso-Soft Metal



The combination of selected natural, Chungking bristles and an open-pore special fabric absorb up to five times more polishing paste than conventional bristles.

### Abraso-Buff Metal



As polishing requires less effort, the user can work in a more relaxed manner and save up to 50% of his time.

### Abraso-Buff Metal Mini



The fabric layers absorb considerably more polishing paste than conventional polishing brushes.

### Abraso-Buff Polipast Metal



Fabric layers impregnated with polishing paste and selected, especially stable, natural Chungking bristles reduce the time required for polishing by up to 60 %.

## High Luster Buff Metal



Fifty particularly closely woven layers create a previously unattainable high luster on all alloys.

## Brushes for handpieces



A large range of brushes for handpieces allows specific polishing of surface details on all alloys.

### Cotton mandrel



No more spinning and punching of the cotton.



## Standard polishing brushes



Chungking black converging Ø 80 mm, 4 rows 12 pieces REF 350 0033 0



Chungking black straight Ø 80 mm, 4 rows 12 pieces REF 350 0031 0

The standard line of polishing brushes includes a wide range of different shapes and sizes for various applications when polishing all types of metal.



The selected Chungking bristles are very rigid and guarantee the long service life of the brushes.



Chungking black converging Ø 70 mm, 3 rows 12 pieces REF 350 0029 0



Chungking black straight Ø 65 mm, 4 rows 12 pieces REF 350 0072 0



The single tufts of bristles penetrate the fissures to create an optimum pre-luster in the shortest po-ssible time.



Chungking black converging Ø 65 mm, 2 rows 12 pieces REF 350 0028 0



Chungking black converging Ø 60 mm, 3 rows 12 pieces REF 350 0073 0



As the star-shape beats the surface gently, this brush polishes all stippled areas quickly and precisely. Abraso-Star K80 (page 486 with integrated polishing particles adds to the effect.



Chungking black converging Ø 50 mm, 2 rows 12 pieces REF 350 0026 0



Chungking black Ø 48 mm 10 pieces REF 350 0047 0



Different versions of these brushes are available for various applications. This makes polishing much less labour intensive.



Chungking black converging Ø 44 mm, 1 row 12 pieces REF 350 0025 0



Chungking black Ø 42 mm 10 pieces REF 350 0048 0



The small brush with a metal hub and short bristles abrades the metal surface and eliminates any traces left after trimming.



Hexagonal brush Chungking black Ø 48 mm 10 pieces REF 520 0004 8



Chungking black tapering Ø 36 mm 10 pieces REF 350 0063 0



This brush (REF 350 0063 0) has a metal hub and tapering bristles. It is particularly hard due to the short bristles. This enables it to polish very slender metal components and junctures aggressively.

### Polishing metal

### Abraso-Soft Metal



The combination of selected natural, Chungking bristles and an open-pore special fabric absorbs up to five times more polishing paste than conventional bristles. Therefore, polishing paste does not have to be applied to the brush as often. This reduces the time required considerably. The fabric is impregnated with abrasive polishing grit and does not require polishing paste for polishing soft alloys. Scratches are eliminated from surfaces faster than when using conventional brushes. This rules out the need for prepolishing with a rubber polisher.

#### Abraso-Soft Metal

Ø 50 mm 1 piece REF 350 0102 1 Ø 80 mm 1 piece REF 350 0081 0



The special nonwoven fabric with integrated abrasives saves time when reducing/polishing all alloys.



The special fabric adapts to every surface optimally to polish large palatal surfaces of chrome cobalt frameworks and create a perfect prehigh-luster finish.

### Abraso-Buff Metal



The 2 x 3 fabric layers absorb very much more polishing paste (e.g. Abraso-Star K80, page 486) and polish effectively. As polishing requires less effort, the user can work in a more relaxed manner and save up to 50 % of his time.

The three rows of high grade Chungking bristles enhance the polish and create an optimum prehigh-luster. Any traces of trimming can be identified immediately and eliminated quickly. This brush is welded together using a special ultrasonic technique which guarantees that the bristles and fabric are gripped firmly.

#### Abraso-Buff Metal

Ø 50 mm 1 piece REF 350 0102 5 Ø 80 mm 1 piece REF 350 0079 0



The 2 x 3 fabric layers and selected, natural Chungking bristles are for polishing all dental alloys.



The slender shape is ideal for reaching areas which are difficult to access and polishing them to a pre-high-luster.

### Abraso-Buff Metal Mini



When combined with selected, natural Chungking bristles, they produce a perfect pre-high-luster on all slender components such as clasps, crowns, inlays etc.





This small metal hub brush, with 4 layers of special fabric, facilictates abrasive polishing of all alloys.



By using different polishing pastes (e.g. Abraso Star K50 REF 520 0016 1, page 487), a perfect prehighluster can be achieved quickly on all alloys.



### Abraso-Buff Polipast Metal



Fabric layers impregnated with polishing paste and selected, especially stable, natural Chungking bristles reduce the time required for polishing by up to 60 %.

The combination of abrasive fabric and high grade Chung-king bristles creates surfaces with no scratches whatsoever, in one single session.

The eight parts of the brush are welded together permanently using ultrasonics and high pressure. This prevents the impregnated fabric layers from redating loose and thus guarantees high strength and a long service life.



Abraso-Buff Polipast Metal Ø 50 mm 1 piece REF 350 0102 6

Ø 80 mm 1 piece REF 350 0086 0 Selected, especially stable, natural Chungking bristles polish more effectively and facilitate prepolishing to create non-streaky surfaces. This saves time and permits the user to work in a relaxed, non-stressed manner. High grade fabric layers store polishing paste and smooth the metal surface. Strong metal retainers grip the bristles securely.




Accurate, abrasive polishing saves time and reduces costs. Scratches etc. caused by trimming are eliminated.



The Abraso-Star K80 polishing paste (REF 520 0016 2) enhances and reactivates the polishing effect as required.



The outcome is impressive: the surface on the left has been sandpapered and that on the pre-polished. All traces left after trimming are eleminated from all chrome cobalt alloys without using rubber polishers.





### High Luster Buff Metal



Fifty particularly closely woven layers create a previously unattainable high luster on all alloys.

High Luster Buff Metal Ø 60 mm, 50 layers 1 piece REF 350 0093 0

Ø 100 mm, 50 layers 1 piece **REF 350 0083 0** 



The hub is welded ultrasonically and grips the 50 fabric layers to prevent them from redating or becoming loose. This guarantees that the buff remains stable during polishing.



This detailed view of the fabric indicates how closely it is woven. The ends of the single threads are free and the millions of minute threads produce a silky soft surface for creating a radiant high luster.



Abraso Star Glaze high luster polishing paste enhances the polishing effect with selected material components. The high luster buff is shaped in the factory so that labour intensive shaping is no longer required before polishing.



50 layers of specially selected fabric form an extremely dense surface for polishing to a high luster. The high density surface generates a high polishing temperature and guarantees an extremely high luster on all alloys.



# **Polishing metal**

# Metal polishing set

#### A complete range of polishing products for all alloys.



Abraso-Soft Metal

Abraso-Buff Metal

High Luster Buff Metal

Pumice polishing paste for polishing acrylic and metal.

### Metal polishing set

Contents:

1 x 150 g Abraso-Star K50, low abrasion 1 x 150 g Abraso-Star K80, high abrasion

1 piece Abraso-Soft metal 1 piece Abraso-Buff Metal

1 piece High Luster Buff Metal Metal 1 x 500 g Pumice Polishing Paste

Abraso Star Glaze 50 ml

REF 350 0085 0



Abraso-Star K80 high abrasion

Abraso-Star K50 low abrasion



Abraso Star Glaze Universal high luster polishing paste for precious metals, non-precious alloys and acrylics.

## Prepolishing with a handpiece

A large range of brushes for handpieces allows specific polishing of surface details on all alloys.



#### Hexagonal brushes Chungking, black 15 pieces each

Ø 13 mm REF 520 0013 0 Ø 19 mm REF 520 0019 0



The light beating effect caused by the star-shape polishes right into the deepest fissures and eliminates scratches in the shortest possible time



#### Pen-shaped brushes Chungking, black, 7mm long

15 pieces each

 Ø 2 mm
 REF 350 0043 0

 Ø 4 mm
 REF 350 0041 0



When used with Abraso Star K80 (page 486), these pen-shaped brushes facilitate polishing areas which are difficult to get at, e.g. the inner surfaces of telescopic crowns.



#### Round brush Chungking, black, double the bristles

15 pieces

 Ø 19 mm
 REF 350 0049 0

 Ø 22 mm
 REF 350 0056 0

 Ø 25 mm
 REF 350 0050 0



The double row of bristles is very stable when polishing wide areas.



#### Round brush Chungking, black

15 pieces Ø 19 mm REF 350 0051 0 Ø 22 mm REF 350 0052 0 Ø 25 mm REF 350 0053 0



The small diameter of the brush is perfect for polishing fragile components and saves time when polishing all slender areas.



# Linen buff coated

15 pieces Ø 22 mm REF 350 0091 0



The high polishing performance and minimal dimensions provide for brilliantly polished surfaces.



#### Round brush Rodeo

15 pieces each

 Ø 15 mm
 REF 350 0095 0

 Ø 18 mm
 REF 350 0096 0

 Ø 21 mm
 REF 350 0097 0



Fabric discs impregnated with polishing paste polish smoother. This saves time because no polishing paste has to be applied.



#### Hexagonal brushes Rodeo

15 pieces each

Ø 13 mm REF 520 0R13 0 Ø 19 mm REF 520 0R19 0



Specially selected tail hairs from wild horses, which vary in hardness between Chungking bristles and goat-hair, are especially suitable for prepolishing soft alloys.

# High luster polishing with a handpiece

#### Produce a radiant high luster, even in the tiniest areas.



Cotton buff 15 pieces Ø 22 mm

REF 350 0065 0



Fluffy, soft cotton fibres create a mirror-like finish on soft alloys.



Linen buff 15 pieces Ø 22 mm

REF 350 0067 0



When used with Abraso Star Glaze, these stable fabric buffs create a brilliant high luster in areas which are difficult to access.



Felt polishing buff

15 pieces Ø 22 mm

REF 350 0064 0



This three layer felt buff polishes crown/facing junctures gently.



Pen-shaped brushes

Goat-hair, white, 7 mm long

15 pieces

 Ø 2 mm
 REF 350 0044 0

 Ø 4 mm
 REF 350 0042 0



The soft goat hairs create a high luster on outer crowns which is gentle to the surface and produces optimum friction.

### Cotton mandrel



#### No more spinning and punching of the cotton.

- special shape of the retaining eyelet ensures that the cotton is safely held
- easy attaching of the cotton saves time during polishing



Cotton mandrel
Ø 2.35 mm
2 pieces
REF 360 0126 9



Attach a piece of cotton in the area of the retaining eyelet and press it on slightly using the finger.



Wind the cotton around the mandrel at a low speed (< 1000 rpm).



Apply polishing paste onto the cotton.



Simple and fast high luster polishing of bars and milled supports.



Occlusal areas and crowns are polished to high luster in a time-saving manner.



High luster finish of friction surfaces in the double crown technique.



## Standard polishing brushes



The soft Chungking brushes simplify polishing of acrylics and produce surfaces without striae.

## Abraso-Sil Acrylic



This brush absorbs particularly much polishing paste or pumice and only releases it very slowly – for efficient prepolishing.

### Abraso-Buff Acrylic



The special textile layers retain the pumice polishing paste longer so that less polishing paste needs to be added.

### Prepolishing Buff Silicone



The silicone coating increases the stability of the buff. This results in increased abrasion capacity – particularly effective during prepolishing.

## Abraso-Soft Acrylic



These materials absorb more pumice and retain it for a longer time. The fibre fabric reduces the friction heat.

### Leather Buff



This leather buff creates a perfect high luster which prevents bacteria and deposits being trapped. Dentures are then easier to clean.

## High Luster Buff Acrylic



The air is circulated continually to polish acrylic coolly and gently.

## Acrylic finishing set



Abrasive finishing and accurate polishing right up to a brilliant high luster quickly and easily – especially in areas which are difficult to get at, such as on orthodontic appliances.

## Brushes for handpieces



These four handpiece buffs produce a brilliant high luster on all dental acrylics.

# Standard polishing brushes



Chungking white Ø 80 mm 4 rows 12 pieces REF 350 0034 0



Chungking white Ø 70 mm 3 rows 12 pieces REF 350 0030 0



The wet, mixed pumice diffuses into the brush and nonwoven fibre fabric (Abraso-Soft Acrylic).



Chungking white Ø 65 mm 4 rows 12 pieces REF 350 0074 0



Chungking white Ø 60 mm
3 rows
12 pieces
REF 350 0075 0



The reduced width of the Abraso-Soft Acrylic allows perfect polishing of interdental spaces.



Chungking white Ø 50 mm 2 rows 12 pieces REF 350 0027 0



Goat-hair metal core Ø 48 mm 10 pieces REF 350 0061 0



Mandrel for polishing brush 1 piece REF 360 0116 8



Chungking white Ø 44 mm 1 row 12 pieces REF 350 0024 0



Chungking white Ø 24 mm 1 row 12 pieces REF 350 0102 3

## Abraso-Sil Acryl



This buff consists of a nonwoven fibre fabric in the centre between two layers of silicone-coated cotton fabric. On the outside there are two rows of bleached Chungkink bristles. This brush absorbs particularly much polishing paste or pumice and only releases it very slowly - for efficient prepolishing.

Abraso-Sil Acryl Ø 80 mm 1 piece REF 350 0099 3 Ø 50 mm REF 350 0102 2



## Abraso-Buff Acrylic



This buff consists of two special textile layers and three rows of bleached Chunking bristles. The special textile layers retain the pumice polishing paste longer so that less polishing paste needs to be added.

Abraso-Buff Acrylic Ø 50 mm 1 piece REF 350 0102 4 Ø 80 mm 1 piece REF 350 0078 0



The reduced width of the Abraso-Soft Acrylic allows perfect polishing of interdental spaces.



## Prepolishing Buff Silicone



The buff consists of 24 layers of a silicone-coated cotton fabric. The silicone coating increases the stability of the buff. This results in increased abrasion capacity – particularly effective during prepolishing.

Additionally, the silicone coating results in considerably extended service life of the buff.

Prepolishing Buff Silicone Ø 80 mm 1 piece REF 350 0099 1 Ø 60 mm 1 piece REF 350 0098 0



### Abraso-Soft Acrylic



This brush consists of a central, nonwoven fibre fabric and bleached Chungking bristles on the outside. These materials absorb more pumice and retain it for a longer time. The fibre fabric reduces the friction heat.

Abraso-Soft Acrylic Ø 50 mm 1 piece REF 350 0102 0 Ø 80 mm 1 piece REF 350 0080 0



The wet, mixed pumice diffuses into the brush and nonwovern fibre fabric.



### Leather Buff




This leather buff creates a perfect high luster which prevents bacteria and deposits from being trapped. Dentures are then easier to clean.

REF 350 0099 0

Leather Buff 1 piece each Ø 60 mm

Ø 80 mm
 Ø 100 mm
 REF 350 0036 0
 REF 350 0035 0





This leather buff can be used at speeds of up to 1,500 r.p.m. to polish acrylics gently and keep them especially cool.



Cool polishing even creates a high luster interdentally, which prevents deposits being trapped.

## High Luster Buff Acrylic



The air is circulated continually to polish acrylic coolly and gently.



This high luster buff is ready for immediate use on a polishing motor and can be used easily and without fraying. Specially selected fabric prevents the acrylic overheating.



The fibre reinforced outer layers provide this buff with a previously unattainable stability.



The 35 layers of textile have been welded into place ultrasonically to prevent them redating and, due to their high strength, create a previously unheard of high luster



The loose woven textile circulates the air during high luster polishing and prevents the acrylic from overheating. Therefore, it polishes very gently.

### High Luster Buff Acrylic 1 piece each

REF 350 0094 0 Ø 60 mm 40 layers Ø 100 mm 35 layers REF 350 0082 0

## Acrylic finishing set

Abrasive finishing and accurate polishing right up to a brilliant high luster quickly and easily - especially in areas which are difficult to get at, such as on orthodontic appliances.

The Diatit coated cutter is especially long lasting and as cost-effective as never before.

Three different abrasive grits provide for accurate polishing right up to a brilliant high luster.



Tungsten carbide cutter 1 piece REF D 200 KF 23



Pressure can be exerted as necessary to reduce the material as required.




Tungsten carbide cutter 1 piece **REF D 263 KG 60** 



Acrylic polisher coarse green 1 piece REF P 243 HG 10



The green, coarse polisher removes all traces of finishing effortlessly.



REF 350 0099 2

Assortment

Abraso-Gum Acrylic 2 Tungsten carbide cutter

5 pieces



Acrylic polisher medium grey 1 niece REF P 243 HM 10



The grey polisher polishes slightly abrasively and prepolishes in one stage.



Shows the polished surface. A brilliant high luster with no scratches



Acrylic polisher fine red 1 piece REF P 243 HF 10



The high luster polisher creates an excellent high luster on all acrylics in the shortest possible



# Polishing acrylic set

A complete range of polishing products for all acylic techniques.



Abraso-Soft Acrylic

Abraso-Buff Acrylic

High Luster Pumice polishing Buff Acrylic paste for acrylic

#### Polishing acrylic set

Content:

1 x 150 g Abraso-Star K50, low abrasive
1 x 500 g Pumice polishing paste
1 piece Abraso-Soft Acrylic
1 piece Abraso-Buff Acrylic
1 piece High luster Buff Acrylic

REF 350 0084 0



Abraso-Star K50 low abrasion

## Polishing acrylic

## Prepolishing with a handpiece

#### Satin soft goat-hair brushes prepolish gently.



Hexagonal brushes Goat-hair white 15 pieces Ø 19 mm

REF 520 0015 1



The star-shape increases the polishing performance by up to 50 % and saves a great deal of time.



Hexagonal brushes Goat-hair white 15 pieces Ø 13 mm

REF 520 0014 1



The smaller brush prepolishes all narrow areas gently.



Round brush Goat-hair white double the bristles 15 pieces Ø 22 mm

REF 350 0055 0



The large diameter is ideal for polishing large surfaces such as facings.



Round brush Goat-hair white double the bristles 15 pieces Ø 19 mm

REF 350 0054 0



Doubling the number of bristles has provided the necessary stability, even when polishing interdental spaces.

## High luster polishing with a handpiece

#### These four handpiece buffs produce a brilliant high luster on all dental acrylics.



Linen buff 15 pieces Ø 22 mm

REF 350 0067 0



This stable fabric buff creates a brilliant high luster, even on the hardest facing acrylics.



Cotton buff 15 pieces Ø 22 mm

REF 350 0065 0



Super soft cotton threads polish palatal rugae optimally and do not leave rough areas which would trap deposits.



Leather Buff 15 pieces Ø 22 mm

REF 350 0066 0



Polishing cool with a leather buff prevents harm to thin metal margins.



Polishing buff felt, 3 layers 15 pieces Ø 22 mm

REF 350 0064 0



The three rows of felt adapt to every structure optimally, which simplifies polishing intricate details.



# Universal polishing / Polishing porcelain

### Abraso-Fix

























yellow - extrafine

Pen-shaped brushes

REF

4 pieces Pen-shaped brushes 1 piece each: extra fine, fine, regular, coarse REF 350 0075 8

2 pieces

2 pieces

2 pieces

350 0046 0

350 0060 0

350 0045 0

350 0058 0

350 0077 0

2 pieces

2 pieces

350 0075 6

350 0075 7

green - coarse

Ø 4 mm REF

Round brush Ø 22 mm REF

blue - regular

Ø 4 mm REF

Round brush

Ø 22 mm REF

Ø4 mm REF

yellow - extra fine

Pen-shaped brushes

REF

Round brush Ø 22 mm REF

Ø 4 mm

Ø 2 mm

Round brush Ø 22 mm REF

red - fine

Pen-shaped brushes

Pen-shaped brushes

Pen-shaped brushes

8 pieces

8 pieces

8 pieces

350 0076 0

350 0075 3


8 pieces

8 pieces

350 0070 0

350 0076 1

350 0076 2



Perfect for quickly prepolishing stippled chrome cobalt denture bases.



The gentle abrasion enables all outer telescopes to be polished to a high luster accurately.



The soft abrasiveness makes it possible for an exact high luster finish even on secondary attachments.



The fine polishing particles create an optimum pre-high luster on all facing acrylics in the shortest possible time.



The particularly slender shape polishes occlusal surfaces right into the smallest fissures.



## Polishing porcelain

1 piece each: extra fine, fine, regular, coarse

Assortment

REF 350 0075 1

4 pieces

Round brush



The yellow Abraso-Fix brush creates a perfect prehigh luster on all porcelains.



Felt wheels unmounted Ø 12 mm 100 pieces REF 350 0071 0



Extremely long lasting due to pressure impregnation and the hardness.



## Metal polishing pastes

For optimum, abrasive prepolishing right up to high luster polishing – specially developed polishing pastes enhance the properties of all polishing brushes. This saves time, allows the user to work in a relaxed, non-stressed manner and improves the quality of the work.



Abraso-Star K80 high abrasion 320 g REF 520 0016 2



Abraso-Star K80 is highly abrasive, which simplifies polishing all non-precious alloys.



As K50 and K80 stick to all polishing brushes well, abrasive polishing can be carried out longer than when using conventional polishing pastes.



Titapol Polishing Paste
150 g
REF 520 0015 3
350 g
REF 520 0015 4



This titanium polishing paste prepolishes abrasively, almost up to a perfect high luster.



A handpiece brush and Titapol prepolish accurately, even in areas which are difficult to access and in the shortest possible time.



Abraso Star Glaze High Luster Polishing Paste 2 x 50 ml REF 520 0016 3



Abraso Star Glaze creates an optimum high luster quickly and easily.



The excellent polishing properties reduce the effort required when polishing with a handpiece.




Brepol 50 g REF 540 0103 7



The round goat-hair brush and Brepol provide perfect high luster on all non-precious metal alloys.



Safe polishing of clasp dentures with the handpiece



Accessories:



Crowns and bridges made of non-precious metal alloys are polished as easily as gold.



After milling, polish telescopic and conical crowns to high luster without prepolishing. Perfectly suitable for the inner surfaces of secondary crowns.



Round brush Goat-hair white double the bristles Ø 19 mm, 15 pieces REF 350 0054 0



Pen-shaped brushes Chungking, black 7 mm long 15 pieces REF 350 0041 0





Attachment, shear distributor and CoCr structure are quickly and neatly polished to high luster.



## Metal and acrylic polishing pastes



Pumice Polishing Paste for metal and acrylic 3 x 500 g REF 520 0016 0



The gentle polishing properties allow all traces left from finishing to be removed from all soft alloys in seconds.



The low abrasion constituents of this pumice polishing paste simplify polishing metal/acrylic junctures.



Abraso-Star K50 slightly abrasive 320 g REF 520 0016 1

### Acrylic polishing paste



Acrypol High Luster Paste for facing acrylics 170 g REF 520 0017 0



Slightly abrasive materials create a virtually perfect high luster. Simply polishing over the surface with a cotton buff is all that's needed to produce a perfect high luster.

## Porcelain polishing paste



Diamond Polishing Paste 5 g REF 540 0014 0

Polishing paste with a high diamond content and long lasting, impregnated, hard felt wheels provide for the best possible finish on all porcelains.

The high percentage of diamond particles provides for maximum abrasion and the best possible high luster.



The special consistency of the polishing paste enables it to diffuse into the felt and polish for up to five times longer.



The paste liquefies during polishing and can be pushed backwards and forwards on the facing without splashing.

## Zi-polish - the perfect mirror polish!



**Zi-polish** 5 g **REF 360 1002 5** 

Zircon polishing paste for high gloss surfaces. Two different diamond grain sizes for pre- and high gloss polishing of marginal edges, basal areas, primary crowns, individual abutments and bars. The special binding agent avoids spraying of the polishing past when dipping the brush into the paste container. This also reduces paste consumption.

- the high content of diamond particles reduces polishing time and leaves a highly glossy surface
- special binding agent reduces paste consumption
- proper wipe off of pre polish residues from dents



Due to a special binding agent the paste sticks to the brush, saving material also.



Perfect and speedy mirror polish is achieved by high diamond particle density.



# Disinfecting and cleaning / Instruments

## Dentaclean mixing fluid for pumice powder



Dentaclean mixing fluid for pumice powder 5000 ml REF 520 0099 8 Dentaclean mixing fluid for pumice powder 1000 ml REF 520 0099 9 Protects against germs.

#### Dentaclean mixing fluid for pumice powder

- Remains moist and free of germs for two to three weeks without having to be remixed.
- Contains skin-care additives to protect employees' hands.
- Contains natural odours which still smell fresh after several weeks.
- Mixed polish adheres to the brush and restorstion better so that the pumice splatters less. This saves time when polishing as the pumice slurry does not have to be applied repeatedly



In wet punmise powder disease microbs are present. Disinfection action occures within one hour.



The nurturing effect for the skin is acheived by skin nurturing addative.

#### Application:

Simply mix the pumice powder with Dentaclean mixing fluid for pumice powder, do not add water. This way the pumice powder stays moist for two to three weeks.

### Pollygrip



Grips all crowns, bridges and inlays firmly for finishing and polishing precisely.

Crown holder, wide

1 piece **REF 360 0100 0** 

Crown holder, narrow

1 piece REF 360 0099 0

Replacement parts:

Special rubber sleeves

100 pieces **REF 360 0096 0** 



The rubber sleeves can be exchanged and grip firmly during all proce-



Even minute inlays can be held firmly and without causing harm.

### Assortment

23 pieces

- 1 Pollygrip
- 1 Crown holder, wide
- 1 Crown holder, narrow
- 20 Special rubber sleeves

REF 360 0095 0





Crown holder, wide 1 piece REF 360 0098 0



Crown holder, narrow 1 piece REF 360 0097 0



	Accessory
3	Chuck262, 387
	Cleaner241
3D planning	Cleaner, Pi-Ku-Plast HP 3637
mini <sup>1</sup> SKY54, 56	Diabolo-Cleaner467
mini¹SKY laboratory analog56	Disposable brush32, 267
mini <sup>1</sup> SKY planning matrix56	Foot rack
3D-resin58-59	Friction silicone
	Inserting instrument141, 144
5	Inserts
	Latch spring194
5-motion-active284	Lock spring196
5-motion-colors	Master-Split metal magnetic plate26
5-motion-glue284	Matrix inserting instrument137, 146, 154,
5-motion-scanner	156, 162, 166-167, 179, 184-185, 188
5-motions-milling cutter285	Matrix pliers
5-motions-zircon283	Measuring cup
	Measuring syringes
A	Metal transfer patrices vks43, 141, 144, 154
	Mixing block
Abdruck-Cut	Mixing cannulas42
Abraso-Buff Acryl345, 347, 479, 481	
Abraso-Buff Metal472, 474	Mixing cup
Abraso-Buff Metal Mini472, 474	Model support BF1262-263, 387-388
Abraso-Buff Polipast Metal472, 475	Optiguss mixing well macro/micro
Abraso-Fix281, 485	Paralleling mandrel for oc and custom bridge sectioning
Abraso-Gum Acryl342, 362	attachments
Abraso-Sil Acrylic479-480	Paralleling mandrel Interlock
Polishing acrylic	Paralleling mandrel metal vsp184-185
Abraso-Soft	Paralleling mandrel metal wstg
Metal polishing set476	142, 145-146, 148, 154-156
Abraso-Soft Acryl345-346, 479, 481	Paralleling mandrel vks
Abraso-Soft Metal472, 474	Paralleling mandrel wstg190-191
Abraso-Star487	PML surveyor
Abraso-Star K50345	PMP surveyor
Abraso-Starglanz345, 353	Spatulas
Abutment	Spray bottle, plastic
SKY fast & fixed221	Steel pins194, 196
SKY fast & fixed Laboratory analog69	Stirring sticks39
Abutments	Thinner for spacer varnish33
Angle measuring device136	Titanium adjusting screw174, 176
Axle abutment	Tungsten carbide drill20, 202, 206, 209, 344, 377
Implant analogue136	Universal paralleling mandrel
Vario-Stud-Snap vks-oc rs 2.2 Abutments	162, 164, 166-167, 194, 196
	Universal paralleling mandrel 2172, 179
Accessories Milling unit	Vessels ppk241
PML surveyor	Acid
PMP surveyor	Brecid pickling agent130
Accessory	Acrylic
•	Abraso-Gum Acryl342, 362
adapter base	Abraso-Sil Acrylic479-480
Angle measuring device	Abraso-Soft Acryl345-346, 479, 481
Auxiliary modelling198, 223-224, 228	Abraso-Buff Acryl345, 347, 479, 481
Base	Adhesive acrylic wax330
Brealloy flux	bre.crystal358-359
Brealloy solder259	bre.crystal HP360
Bresol Speed	bre.dentan
Brush + brush holder108, 248, 255	bre.dentan HP
Brush cleaning pot, Optiguss108, 248	bre.flex
Brush holder, bent	bre.flex 2 <sup>nd</sup> Edition
Brush holder, straight32	High lustre buff Acryl345, 479, 482, 484
Brush pen12, 86, 332-333	Location matrix drill
Ceramic spacer176	Polishing acrylic
	1 Shoring doryno

Special diamonds for the veneering technique	
Vario-Soft 3 zircon sv	
Acrylic polishing set	
Acrylic polishing set	
Acrylic Sep	zetioning to cut
Acrylic veneering Felt polishing buff	
Felt polishing buff	
Retention Crystals and Beads	
Activating pliers	
Activating of telescopes Activating pliers	
Activating pliers	
Activating pliers	
Activating screw	
Activating telescopical crowns Activating pliers	
Activating pliers	25-27, 40 12628221172172230 ectioning attachment230
Activation spray 5-motion-active	
5-motion-active	
Adapter Assembled transverse fixation	
Wax adapter.237, 407AttachmentAdapter base.20, 377Activating screwAdhesiveBasal screw5-motion-active.284Bridge sectioning attachment oc5-motion-glue.284Bridge sectioning stud oc, bridge sectioning stud o	
Adapter base	
Adhesive Basal screw	
5-motion-active	230 sectioning attachment230
5-motion-glue	ectioning attachment230
	230
DTK-adhesive	
138, 146, 150, 179, 202, 206, 209, 224, 243 Ceramic spacer	172
Retention adhesive	204
Adhesive (glueable) pin Custom bridge sectioning attachm	
Master-Pin Radix-K19 Cylindrical attachment	
Adhesive acrylic wax	
Adhesive connections Inverto Plus	
Double t-adhesive connector244, 245 Vario Compress 1	
Double t-adhesive mini connector246, 247 Vario Compress 2	
Adhesive resin wax  Vario-Soft 3	
Thermo-syringe29, 330, 392 Vario-Soft 3 matrix housing	
Adhesive sleeve Vario-Soft 3 mini	
Security-Lock-adhesive sleeve224-225 Vario-Soft 3 mini sv	
äe-mkbl Vario-Soft 3 sv	
Aesthetic and ergonomic pontics for ceramic, Auxiliary  Vario-Soft 3 sv bridge-sectioning a	
wax elements wht	
Aesthetic and ergonomic pontics for ceramic, Auxiliary wax  Vario-Soft-Bar vss	
elements wht	
Aesthetic wax veneers, auxiliary wax parts wht50, 100 Vario-Stud-Snap vks-oc	
Aesthetic-Gnathoflex102-104 Vario-Stud-Snap vks-oc exchange	able stud with
Isoflex	
Airaqua TurbineVario-Stud-Snap vks-sg	
Air-drying spacer varnishvario-Stud-Snap vks-sg exchange	able stud with
Alloy adhesive sleeve	150, 153
Brealloy C + B 270127 Wax attachment elements	191
Brealloy F 400	
Brealloy flux	255
Brealloy MK	
Brealloy MO258 Biotec blocking out wax	65, 235
Brealloy solder	93 441
Aluminium Catridges empty	160-167
Anchoring Attachments	
	100
9	
Locking Pin Easy-Snap A	190
Locking Pin Easy-Snap E	470
Locking Pin Snap A	
Locking Pin Snap E	
Vario-Stud-Snap vks-sg 151-152, 154-158 Swivel-type lock system src	196-197

Auxiliary	Battery-powered screwdriver
Auxiliary modelling198, 223-224, 228	Cordless Prosthodontic Screwdriver 38-
Auxiliary element (aid)	Beads
Biotec metal-ceramic blocks without collar96	Retention Crystals and Beads109
Auxiliary elements (aids)	Beauty setup5
Gnathoflex-Premium105-107	Between
Auxiliary modelling198, 223-224, 228	Between hollow blocks9
Auxiliary modelling element	Between hollow blocks9
Friction Splint FS1 Splint226	Betweenblocs9
Auxiliary modelling element	BF 2
vks-sg matrix housing157	Milling unit BF 2262, 387-38
Wax matrix housing zg178	BG 1
Auxiliary wax elements	Bre-Gel BG 125
Betweenblocs97	BG 2
Protek clasp/bar joiners236	Bre-Gel BG 225
Auxiliary wax elements wht96-101	BG 3
Aesthetic and ergonomic pontics for ceramic, Auxiliary	Bre-Gel BG 3254, 33
wax elements wht100	BigBrush
Between hollow blocks98	Unique Brush277, 40
Hollow pontic blocks, Auxiliary wax elements wht98	Bio Dentaplast358-359, 361, 38
Hollow pontics, Auxiliary wax elements wht101	Biofunctional therapy
In-between hollow pontics, auxiliary wax elements wht97	BioHPP353-356, 38
In-between pontics, auxiliary wax elements wht97	Biotec blocking out wax65, 23
Metal-ceramic blocks, Auxiliary wax elements wht99	Biotec cervical wax9
Auxiliary wax parts	Biotec metal-ceramic blocks without collar9
Pontics with shallow collars for ceramic, Auxiliary wax	Biotec milling wax
elements wht99	Biotec modelling wax grey/green9
Auxiliary wax parts wht	Biotec reels of wax pattern110
Aesthetic wax veneers, auxiliary wax parts wht	Biotec-Wax for outer copings9
50, 100	Bite blocks
äwv	Bite registration
Aesthetic wax veneers, auxiliary wax parts wht	Security-bite blue
	Bite-taking A0. 22
Axle abutment136	Bite blocks
_	Black hair brush MagicBrush274-275, 299, 398-39
В	Blades for wax knife
Ball Clip403	Blanks
Bar	5-motions-zircon28
Protek wax pattern, lingual bar235-236	breCAM.BioHPP
Vario-Soft-Bar vss188-189	breCAM.cutter
Wax attachment elements191	breCAM.resin
Wax bar attachments190	breCAM.resin milling blanks7
Wax bars142, 156, 190-191, 206, 209	Blanks
Bar attachement	breCAM.wax7
Vario-Soft 3 zircon sv168, 283	Blocking out wax
Bar pattern .	Biotec blocking out wax65, 23
Vario-Soft-Bar-Pattern vsp184-186	Block-out kneading material25
Bar, Vario-Soft-Bar-Pattern	Blocks
Vario-Soft-Bar-Pattern vsp184-186	Between hollow blocks9
Basal screw172	Blooking out at dies
Base16, 123, 376	Undercut wax3
Adapter base20, 377	Blooking out discs141, 14
Base former(s)	Blooking out of dies
Master-Split base former(s)26	Undercut wax3
Base material	Blooking out pattern
HIGH.lign321	Protek236-237, 239-24
Base plate252	Blooking-out material
Base plate ring252	Transblock1
Base stone	Blue-Clip40
Fluid-Rock28, 74	Bonder
Base stone Arti-Rock28	Multisil-Primer34
Basic (base) material	Bonding
Opaquer mixing liquid271	Ceram-Bond27

Bonding		Brevest ESG	
Chrom-Kobalt-Bonding	270	Brevest exakta M	256
Bonding agent		Brevest exakta speed	
FGP Bonding agent		Bresol Speed124	
FGP insulating liquid202, 206, 209, 22		Brevest exakta Speed	256
Multisil-Primer34	10, 372	Brevest for 2 press	353
visio.link	297	Brevest M1	255
X-connector bonding agent	.55. 58	Brevest Rapid 1124	1, 256
Bonyhard clasp		Bridge sectioning attachement	•
Protek clasp patterns2		Vario-Soft 3 conicalbridge	169
Book	0. 200	Bridge sectioning attachment	100
Manual bredent casting technique	8/	Bridge sectioning attachment  Bridge sectioning stud oc, bridge sectioning attached	hment
Booklet		0C	
	000		
Golden booklet12	26, 260	Vario-Soft 3 sv bridge-sectioning attachment 16	
Bow, wax pattern	00.040	Bridge sectioning attachment oc	230
Protek236-237, 2		Bridge sectioning stud oc, bridge sectioning attachn	
bre.crystal3			
bre.crystal HP		Bridge, Crown-bridge-inlay wax, KBI	92
bre.dentan358-35		Bridges	
bre.dentan HP		Quadrosticks	117
bre.flex358-359, 36	52, 383	Wax pattern sticks	116
bre.flex 2 <sup>nd</sup> Edition	362	Bridges and crowns	
bre.lux Power Unit65, 291, 299, 306-30	07. 379	Brealloy C + B 270	127
Brealloy C + B 270		Brevest Rapid 1124	
Brealloy F 400		Brochures	., 200
Brealloy flux12		Patient brochures	132
Brealloy Lot		Brush	102
		Goat-hair brush	252
Brealloy MK			
Brealloy MO		KoliBrush	
Brealloy solder		MagicBrush274-275, 299, 39	
breCAM.BioHPP		MagicContrast274-275, 39	
breCAM.cutter	,	Transfuser125	5, 169
breCAM.resin		Unique Brush277	
breCAM.resin milling blanks		Brush + brush holder108, 248	
breCAM.wax	78	Brush cleaner, Optiguss108	3, 248
breCeram	281	Brush cleaning pot, Optiguss108	3, 248
Brecid pickling agent	130	Brush holder	278
breciform D impression tray		Brush holder, bent32	
brecision implant heavy		Brush holder, straight	
brecision implant light		Brush pen12, 86, 33	
brecision impression materials		Brush resin	_ 000
·		Pi-Ku-Plast HP 36 ppk	27
brecision putty soft	50	Brushes	31
bredent activating pliers	24 405		470
Activating pliers26	04, 405	Pen-shaped brushes	4/8
preformance		Brushes for polishing with the handpiece	
top-lign breformance		346-348, 472, 477, 479	), 484
Bre-Gel BG 1	254	Buff	
Bre-Gel BG 2		Abraso-Buff Acryl345, 347, 479	
Bre-Gel BG 325		Abraso-Buff Metal472	2, 474
Brenometer lead holder26	52, 387	Abraso-Buff Metal Mini472	2, 474
Brenometer locating pin		Abraso-Buff Polipast Metal472	2, 475
Brenometer plate		Cotton buff for handpiece345	
Brenometer surveying system26	63. 388	Fabric buff	
Brepol		Felt polishing buff	
Bresol for 2 press		High lustre buff Acryl	
Bresol M		High lustre buff metal472	
		<del>-</del>	
Bresol N		Prepolishing Buff Silicone	40 I
Bresol R		Buff, leather buff	. 404
Bresol Speed12	24, 256	Leather buff	ı, 481
Brevest C+B speed		Bunsen burner	
Bresol Speed12		Piezo-Blitz	393
Brevest C+B Speed	124	Bur	
Brevest ceram speed		Conical bur44	6-448
Bresol Speed12	24, 256	Groove burs	449

Bur	Casting ring
Parallel bur442-445	Casting ring marker126
Processing of titanium438-439	Investment marker126, 260
Rapidy Microbur with relief202, 206, 423	Casting ring former
Shoulder bur449	Crêpe sleeve260
Silicone bur340, 362	Casting ring marker126
Tungsten carbide tools with relief 419-422, 424-425	Casting ring system
Burnisher	Wash primer for fleece121
Gloss and hardening agent for plasters34	Casting ring system
Burn-out resin	Fleece liners121
Pi-Ku-Plast HP 36 ppk37, 241-242	Casting rings
Pi-Ku-Plast ppk241-242	Fleece liners121
Bush	Wash primer for fleece
Threaded bush, refractory alloy, cast-on232	Casting system338-339
Betweenblocs97	Casting technique
	Block-out kneading material
bwg	Double rinsing heads
In-between pontics, auxiliary wax elements wht97	Double rinsing heads for voluminous castings115 Duplicating system251
bwhbl Between hollow blocks98	Flask sleeve
bwhq	Flask tray251
In-between hollow pontics, auxiliary wax elements wht97	Isosil
bwm 3	Manual bredent casting technique84
foam rubber grip lining89	Rinsing heads114
rest for handpiece89	Rinsing heads for voluminous castings114
Wax knife bwm 389	Spacer – base insert251
bwm 3	Sprues for centrifugal casting113
footswitch89	Sprues for centrifugal casting for voluminous castings
	113
С	Stabilizer251
	Technolit250
C + B 270	Technosil duplicating silicone250
Brealloy C + B 270127	Sprues for vacuum pressure casting112
CAD/CAM	Sprues for vacuum pressure casting for voluminous
breCAM.BioHPP79	castings112
breCAM.cutter73, 81	Casting technique
breCAM.resin80	Brealloy MO258
breCAM.resin milling blanks72	Optiguss108, 248
breCAM.wax78	Wax knife, electrically heated89, 386
CAELO76-77	Casting wax
e.cad	Cervical wax93
e.order	Dipping wax87
e.scan	Gecko
CAELO76-77 Cannulas	K2 exact carving wax91
Application cannulas95	KBI, Crown-bridge-inlay wax92
• •	Life-Color wax50, 92
Cast supply  Rlock-out kneeding material 251	Modelling waxes
Block-out kneading material251  Duplicating system251	Splendido92
Flask sleeve251	Cast-on
Flask tray251	Threaded bush, refractory alloy, cast-on
Isosil251	Cast-on matrix sleeves, Security-Lock-System
Spacer – base insert251	Matrix sleeves, Security-Lock-System
Stabilizer251	Centrifugal casting Sprues for centrifugal casting113
Technolit	Sprues for centrifugal casting for voluminous castings
Technosil duplicating silicone250	Sprides for centificing casting for voluminous castings
Casting	Centring drill
Fleece liners121	Tungsten carbide centring drill224, 281
Manual bredent casting technique84	Cerafine
Metal muffle rings121	Ceraflex diamond disc453, 455
Reels of wax patterns116, 239	Ceragum281, 353
Silicone muffle rings122, 353	Ceram-Bond270
Wash primer for fleece121	Ceramic
Casting pears117	Auxiliary modelling198, 223-224, 228

Ceramic		Clasp patterns	
Ceram-Bond	270	Molar clasp	237
Ceramic separating liquid	272	Protek	236-237, 239-240
Ceramic separating set	272	Clasp/bar joiners, Protek	
Ceramix	273, 397	Protek clasp/bar joiners	236
Chrom-Kobalt Bonding	270	Clean	
Mamelon cutter		Dentaclean impression and denture	e disinfectant
Microkeramik	255		
Opaquer mixing liquid		Cleaner	24
Pontics with shallow collars for ceramic		Cleaner, Pi-Ku-Plast HP 36	
elements wht		Diabolo Cleaner	
Quicktool		Cleaner, Pi-Ku-Plast HP 36	
Special diamonds for the veneering ted		Cleaning	
		Brush cleaner, Optiguss	108, 248
Spot Clip		Dentaclean plaster removing agent	
Spot Clip with supporting ring		Dentaclean plaster removing agent	
Stain liquid		Cleaning of dentures	оросси
Threaded rods		Dentaclean denture cleaning agent	327
Ceramic (porcelain)		Domacican domaio cicaring agone	
Brevest ESG	257	Cleaning of impressions	
Ceramic bonding alloy	201	Dentaclean impression and denture	a disinfectant
Brealloy C + B 270	127		
Ceramic brush	121	Cleaning solution	
Brush holder	278	Dentaclean ultrasonic cleaning age	ant 329
KoliBrush		Clip	320
MagicBrush 274-27		Spot Clip	206
MagicContrast			
•		CNC milling techniques breCAM.wax	70
Unique Brush	277, 401		
Ceramic elements	106 107	CAELO	
Swivel-type lock system src	196-197	e.cad	
Ceramic instruments	00 205	e.order	
Quick Change	90, 395	s breCAM.cutter	
Ceramic mixing plates	270 200	CNC milling techniques Round blanks	
Mixing plates	279-280	breCAM.BioHPP	
Ceramic polishing	204 252	breCAM.resin	80
Ceragum	281, 353	Co Cr work	05/
Ceramic primer	004 000 007	Magnetic plates	252
K-Primer		Coarse grinding tool	457.45
MKZ Primer		Diacryl grinder dcs	457-458
Ceramic removing tool	174, 223	Coated line buff for the handpiece	47/
Ceramic screw	000	Fabric buff	
Security-Lock-Ceramic		Coating	3/4
Ceramic screws		Coating of CoCr elements	0.4
Ceramic separating liquid		Kompaktopaker tooth-coloured UV	
Ceramic separating set		Ropak compact opaquer	34′
Ceramic spacer		Ropak UV powder + liquid	34′
Ceramic spacer E		Cobalt	
Ceramix		Chrom-Kobalt-Bonding	270
Cervical disc		CoCr clasp supply	
Cervical wax		Brenometer lead holder	
Biotec cervical wax	93	Brenometer locating pin	
Change system		Brenometer plate	263
Quick Change	90, 395	CoCr clasp work	
Check bites		Statik-Disc	234
security-bite blue		CoCr supply	
Chrom-Kobalt-Bonding	270	Bio Dentaplast	358-359, 361, 383
Chuck	262, 387	Brealloy F 400	258
Circumferential clasp, bent		Brealloy flux	
Clasp		Brealloy solder	
Thermo-Pen	66, 294	Brevest M1	
Thermoplastic clasp		Brevest Rapid 1	
Clasp patterns		Crêpe sleeve	
Bonyhard clasp	237	Diatit power burs	
Circumferential clasp, bent		Double t-adhesive connector	

CoCr supply	0.40.047	Conical
Double t-adhesive mini connector		Interlock 2°
Exaktosil N 15		Conical and telescopic crowns
Exaktosil N 21		Brenometer locating pin263
Exaktosil silicone duplicating materials		Brevest Rapid 1124, 256
Interlock 2°		Transfer device263, 388
Investment hardener		Conical bur446-448
Kompaktopaker tooth-coloured UV		Conical crowns
LV 1 laser joint	243	Biotec blocking out wax65, 235
Metal polishing set	476	Biotec milling wax93, 44
Microkeramik	255	Conical bur446-448
Optiguss	108, 248	Conical crowns, restoring the friction
Polierjet Duo-Finish	390	Activating pliers264, 405
Polierjet Quadro-Finish	390	Conical drilling
Protek wax adhesive		Tungsten carbide facing cutter228
Quick mandrel system	394	Conical screw
Ropak compact opaquer		Titanium screws164, 228, 230, 232
Ropak UV powder + liquid		Titanium screws, extended head228
Stippled wax sheets		Conicalbridge
Wax patterns cut to size, Protek		Vario-Soft 3 conicalbridge169
CoCr techniques	240	Connections
Casting ring marker	126	Double t-adhesive connector244-245
CoCr work	120	Double t-adhesive connector
Base former	252	
		LV 1 laser joint
Base plate		Contarction measuring device
Base plate ring		KoEx Measuring Device
Bre-Gel BG 3		Contrast  Magic Contract  274 275 200 200
Brevest exakta M		MagicContrast274-275, 398-399
Brevest exakta Speed		Сору
Duro-Top		Base former
Inverto Plus		Base plate252
Master-Copy		Base plate ring252
Silicone sleeve	252	Magnetic plates252
Stabilizer		Master-Copy252-253
Wax adapter	237, 407	Silicone sleeve252
Cold-curing resin		Stabilizer252
top.lign breformance liquid cold	68	Cordless Prosthodontic Screwdriver384
Coloring liquid		Cotton buff345, 478
5-motion-colors	283	Brushes for polishing with the handpiece
Colors		346-348, 472, 477, 479, 484
Multisil intensive colors	373	High lustre polishing with the handpiece345
combo.lign		Cotton mandrel472, 478
Opaquer set	305	CPS
combo.lign luting composite	297	Cordless Prosthodontic Screwdriver384
combo.lign Opaquer		Cradle
Combustible resin		Quick cradle (handle), Quick mandrel system 394
Pi-Ku-Plast HP 36 ppk	37 241-242	Cr-Co alloy
Pi-Ku-Plast ppk		Brealloy MO258
Communication software		•
CAELO	76 77	Cr-Co supply
		Block-out kneading material
e.cad		Duplicating system
e.order		Flask sleeve25
e.scan		Flask tray25
compoForm UV	•	Investment aid25
Application cannulas	95	Isosil25
Composite		Spacer – base insert25
compoForm UV	95, 282	Stabilizer25
DTK-adhesive		Technolit250
138, 146, 150, 179, 202, 206,	209, 224, 243	Technosil duplicating silicone250
Composite finishing set		crea.lign291, 299-300, 305
Compress		Crêpe sleeve260
Vario Compress 1	174-175	Cross cut, rotating tools with cross cut
Vario Compress 2		Tungsten carbide tools with relief419-422, 424-425
1		Cross-cut burs

Crown an bridge technique	Crowns and bridges
bre.dentan358-359, 383	Special diamonds for the veneering technique
bre.dentan HP360	457, 461
Crown and bridge material	Splendido92
top.lign breformance liquid cold68	Sprues for centrifugal casting113
top.lign breformance liquid hot	Sprues for centrifugal casting for voluminous castings
top.lign professional326	113
Crown margin wax	Sprues for vacuum pressure casting
Cervical wax	Sprues for vacuum pressure casting for voluminous
Crown wax K2 exact carving wax91	castings
Crown-bridge-inlay wax, KBI92	Stain liquid271 Wax attachment elements191
Crowns and bridges	Wax bur
Between hollow blocks98	Wax knife, electrically heated
Betweenblocs97	Wax-Lite surface tension reducing agent 10, 120
Biotec blocking out wax65, 235	Wax pattern sticks116
Biotec cervical wax93	Wax patterns cut to size, Protek240
Biotec metal-ceramic blocks without collar96	Crystal
Biotec milling wax93, 441	bre.crystal358-359
Biotec modelling wax grey/green93	bre.crystal HP360
Biotec reels of wax pattern116	bre.flex358-359, 362, 383
Brealloy C + B 270127	bre.flex 2 <sup>nd</sup> Edition362
Brevest C+B Speed124	Crystals and beads
Brevest M1255	Retention adhesive109
Brevest Rapid 1124, 256	CT
Bridge sectioning stud oc, bridge sectioning attachment	X-connector bonding agent55, 58
oc230	X-resin55-56
Cross-cut burs450	X-resin flow57
Custom bridge sectioning attachment231	Cup
Diamond grinding tools451-452, 457, 462	Measuring cup338
Dipping wax87	Mixing cup
Dipping wax Elaflex	Curing material
Dipping wax Visio-Dip86	combo.lign Opaquer296
Double rinsing heads	Curing unit, Light-curing unit Polylux 2
Double rinsing heads for voluminous castings115	Polylux 213-14 Curing unit, Light-curing unit Polylux PL 20
Double t-adhesive mini connector246-247 Felt polishing buff478	Polylux380
Gnathoflex-Premium105-107	Custom bridge sectioning attachment231
Inverto Plus172	Cut
In-between hollow pontics, auxiliary wax elements wht97	Abdruck-Cut
In-between pontics, auxiliary wax elements wht97	Cutters
Isobre wax insulating (separating) liquid86	5-motion-scanner284-285
Life-Color wax50, 92	5-motions-milling cutter285
Mamelon cutter396	breCAM.cutter81
Master-Pin23	Generation M353, 416-417
Master-Pin sleeve23	Mamelon cutter396
Master-Pin system22-24	Cyanoacrylate adhesive
Modelling tips386	5-motion-active284
Modelling waxes69, 91-93	5-motion-glue284
Multisil-Mask hard43	Cylindrical attachment178-181
Multisil-Mask, soft42	Impression transfer set179
Opaquer mixing liquid271	
Optiguss	D
Pollygrip406, 488	Dehubblizing of cilicons and way surfaces
Porcelain mixing liquid271	Debubblizing of silicone and wax surfaces Silicone and wax surface tension reducing agent
Processing of titanium	10, 120
Quadro wax profile	Dentaclean
Quadrosticks	Shipping bags9, 327
Reels of wax patterns	Dentaclean denture cleaning agent327
Retention adhesive	Dentaclean impression and denture disinfectant 9, 327
	Dentaclean plaster removing agent
Rinsing heads114 Rinsing heads for voluminous castings114	Dentaclean plaster removing agent Speed328
Rodeo326, 362	Dentaclean mixing fluid for pumice powder 328, 488
	- · · · · · · · · · · · · · · · · · · ·

Dentaclean ultrasonic cleaning agent	328	Diamond tools	
Dentalos - Accessories		Airaqua Turbine	389
Acrylic Sep	362	Diamond tools	
Dentan		Diabolo sintered diamonds	463-470
bre.dentan358-3	59, 383	Diacryl grinder dcs	457-458
bre.dentan HP	360	Diacryl grinding tool	
Dentasil tooth protection silicone		Diagen-Turbo-Grinder286, 3	
Denture acrylic		Diamond grinding tools	
bre.crystal	358-359	Diatit	
bre.crystal HP		Processing of titanium	438-439
bre.flex358-359, 3	62. 383	Diatit and tungsten carbide tools	
bre.flex 2 <sup>nd</sup> Edition	362	Diatit bur	
Denture cleaning solution		Diatit and tungsten carbide tools	423, 426-439
Dentaclean denture cleaning agent	327	Diatit burs, tungsten carbide burs	
Denture disinfectant	027	Tungsten carbide burs for processing	of plaster 31
Dentaclean impression and denture disinfectan	t	Diatit power burs	
		Diatit tungsten carbide	
Denture repair resin	5, 521	Master-Pin Diatit tungsten carbide ste	an drill
Qu-resin65, 70, 2	02 362		
	.02, 302	Diatit-Multidrill with stop	
Denture resin HIGH.lign	224	Diatit-Multidrill, twist drill	220
•			
uni.lign		198, 202, 206, 209, 2	<u>1</u> 23-224, 228, 362
uni.lign speed	325	Die material	07 044 040
Denture teeth		Pi-Ku-Plast HP 36 ppk	
neo.lign56, 65,		Pi-Ku-Plast ppk	241-242
neo.lign A295, 304, 3		Die varnish	
neo.lign P295, 304, 3		diephos dentine	
Diabolo Cleaner264, 4		Die varnish, light-curing, opaque	32
Diabolo sintered diamonds4		Die varnishes, various	
Diacryl grinder dcs		Light-curing die varnish	
Diacryl grinding tool		Diephos dentine	33
Diacryl rubber grinding tool		Dies	
Diacryl grinder dcs		Retention pins	36
Diagen-Turbo-Grinder286, 353, 457, 4	159-460	Digital	
Diamond		CAELO	76-77
FG-Diabolo	457	e.cad	77
Diamond disc, electroplated		e.order	76
Elastic4	53, 456	e.scan	77
Flexible4	53, 456	Dipping wax	87
Superflex4	53, 456	Dipping wax Elaflex	86
Transflex4		Dipping wax Visio-Dip	
Transflex-T4	53, 455	Dipping wax Elaflex	
Ultraflex4		Dipping wax unit	
Diamond disc, electroplated, Ceraflex4		Waxpool duo	385
Diamond disc, electroplated, Giflex-TR30,		Dipping wax Visio-Dip	
Diamond for turbine		Disc	
FG-Diabolo	457	Blooking out discs	141 144
Diamond grinding disc Giflex TR Master x-tray		Cervical disc	
30, 4		Statik-Disc	
Diamond grinding tool	100 101	Disinfection	207
Set-up grinding tool331, 342, 4	57 /61	Dentaclean impression and denture of	disinfectant
Diamond grinding tools451-452, 4		Dentacioni impression and dentare c	
Diamond grinding wheels		Disinfection bag	
Diamond point Veneering technique			0 227
	302	Shipping bags Disinfection bath 3I	
Diamond polishing paste	00 407		
Zi-polish	.os, 401	Disposable brush	32, 267
Diamond polishing paster, felt wheels	OF 407	DKZ Mixing consults	40
Polishing porcelain	oo, 40/	Mixing cannulas	
Diamond tips	04 405	Dosing device	
Inserts	04, 405	Double rinsing heads	
Diamond tool, sintered	457	Double rinsing heads for voluminous ca	
FG-Diabolo	45/	Double t-adhesive connector	
		Double t-adhesive mini connector	246-247

Dowel pins	E
Retention pins36	e.cad77
Drill	e.order
Location matrix drill	e.scan
Master-Pin Diatit tungsten carbide step drill	Easy-Snap
	Locking Pin Easy-Snap201, 204
Rapidy Microbur with relief202, 206, 423	Easy-Snap A
Tungsten carbide centering174, 223, 228	Locking Pin Easy-Snap A203
Tungsten carbide centring drill224, 281	Easy-Snap E
Orill sleeves	Locking Pin Easy-Snap E202
SKYplanX drilling sleeves59	ecovac
Drilling	base16, 123, 376
Diatit-Multidrill with stop	ecovac Vacuum Mixing unit
Diatit-Multidrill, twist drill	ecovac Vacuum Mixing unit
198, 202, 206, 209, 223-224, 228, 362	Elaflex
Milling unit BF 2262, 387-388	Dipping wax87
Tungsten carbide facing cutter228	Dipping wax Elaflex86
Orilling oil	Elastic
Milling and drilling oil	Electric wax knife
202, 206, 209, 219, 222-224, 228, 231, 449	Foam rubber grip lining89
Drilling templates	Footswitch 89
3D-resin58-59	Rest for handpiece
Orilling unit	Wax knife bwm 389
adapter base20, 377	Wax knife, electrically heated89, 386
Master pin drill unit20-21, 377	·
Drum	Waxpool duo
Polishing drum390	Electroplated diamond
dtk	Special diamonds for the veneering technique
Double t-adhesive connector244-245	Electroplated diamond disc
Double t-adhesive mini connector246-247	Electroplated diamond disc
DTK-adhesive	Ceraflex diamond disc
138, 146, 150, 179, 202, 206, 209, 224, 243	Elastic
Duo-Finish	Microflex
Polierjet Duo-Finish390	Superflex
Duplicate model	Transflex
Base former252	Transflex-T
Base plate252	Electroplated diamond disc, Giflex-TR30, 453-454
Base plate ring252	Electroplated diamond discs Flexible453, 456
Magnetic plates252	Ultraflex
Silicone sleeve252	- · · · · · · · · · · · · · · · · · · ·
Stabilizer252	Electroplated diamond tools
Duplicating gel	Diacryl grinder dcs
Bre-Gel BG 1254	Diacryl grinding tool
Bre-Gel BG 2254	
Bre-Gel BG 3254, 338	Aesthetic-Gnathoflex
Duplicating matrixes zg178	
Duplicating silicone	Epithetics
Technosil duplicating silicone250	coating
Duplicating system	Isoplast
Isosil251	Matt sealing agent for epitheses
Duplicating system251	Matting powder
Block-out kneading material251	Modelling wax for Epithetiks
Flask sleeve251	Multisil Epithetik
Flask tray251	Multisil Epithetik set
Investment aid251	Multisil-Epithetik transparent
Spacer – base insert251	Multisil-Epithetik hard-form
Stabilizer251	Multisil-Epithetik soft-form
Duro-Top259	Multisil fibers
DVT	Multisil intensive colors
X-connector bonding agent55, 58	Multisil sealing agent
X-resin55-56	Multisil sealing liquid340
X-resin flow57	Multisil stains
	Multisil-Epithetik beach
	Multisil-Epithetik city
	Multisil-Epithetik country371
	Marian Epitricux country

Epithetics		FGP insulating liquid	202, 206, 209, 224, 267
Multisil-Epithetik thickener	372	Fibers	
Multisil-Primer		Multisil fibers	373
primer		Fine investment material	
Equipment		Microkeramik	255
Waxpool duo	88	Finish	
Ergonom wax knife		Polierjet Duo-Finish	390
Ergonomic		Finishing of ceramic and acryl	
Aesthetic and ergonomic pontics for cera	amic, Auxiliary	Special diamonds for the ve	
wax elements wht			457, 461
ESG		First tap, tungsten carbide	
Brevest ESG	257	Fissure tool	
Esthetic SetRed-white esthetics set	303	Fit	
Exactoform models		FGP friction fit system	265, 266-267
Gnathoflex Study model FF1 mini	105	Fix, Abraso-Fix	
Gnathoflex Study model FF1 Set		Fixation of a matrix in acrylics	
Exakta M		Metal matrix housing vks-od	
Brevest exakta M	256	Titanium matrix housing vks	
Exakta Speed		Fixing screw, bridge sectioning	
Brevest exakta Speed	256	Fixing screw, Screw connection	
Exakto-Form insulating liquid		Flame	
Exakto-Form, model resin		Piezo-Blitz	393
Exakto-Rock S		Flask	
Exaktosil N 15		Flask sleeve	
Exaktosil N 21		Flask tray	
Exaktosil silicone duplicating materials		Fleece liners	
Expansion measuring device		Flexible	
KoEx Measuring Device	17 45 378	Flexible Acrylic sprue wax	
Extraction system	17, 10, 070	Flexible elements for occlusal	
Protective chamber with/without extraction	n nozzle 380	Aesthetic-Gnathoflex	
1 Totootivo onambor with without oxtraotio	71 1102210000	Fluid-Rock	
_		Flux	20, 7
F		Brealloy flux	128 259
F 400		Foam rubber grip lining	89
Brealloy F 400	258	Foot rack	
Fabric buff		Footswitch	
Fabric buff, coated		for 2 press	
Brushes for polishing with the handpiece		Forceps	
346-348, 472,		Ball Clip	403
Facial prosthetics (Epithetik)	• •	Blue-Clip	
Multisil Epithetik	371-374	Spot Clip	
Multisil Epithetik set		Former(s)	
Facing cutter		Master-Split base former(s)	26
Tungsten carbide facing cutter	228		m/small26
fast & fixed	_	Friction fit system	1// 311141120
SKY fast & fixed	221	FGP Bonding agent	267
SKY fast & fixed Laboratory analog	69	FGP insulating liquid	
Felt			202, 200, 203, 224, 207
Felt polishing buff	478	Friction silicone	
Felt polishing buff		Friction Splint FS1	
Felt wheels	_	Friction Splint FS1 sleeve	226
Polishing porcelain	390, 485, 487	Friction Splint FS1 Spacer	
FG diamond, sintered		Friction Splint FS1 Splint	
FG-Diabolo	457	Front teeth	220
FG diamonds		novo.lign A	204 304 308-300
Airaqua Turbine	389	<u> </u>	
FG-Diabolo		visio.lign FS1	200-303, 300-320
fg-mkbl			226 227
Pontics with shallow collars for ceramic,	Auxiliary wax	Friction Splint FS1Friction Splint FS1 sleeve	
elements wht		Friction Splint FS1 Spacer .	
FGP			
FGP friction fit system	265-267	Friction Splint FS1 Splint FSFriction Splint FS1 Splint	
FGP Bonding agent			
FGP friction fit system		Functional margin wax	50

G	Handpiece, polishing with the handpiece, leather buff
Gecko92	478 High lustre polishing with the handpiece345
Gel	Handpiece, polishing with the handpiece, polishing buff478
Bre-Gel BG 1254	High lustre polishing with the handpiece345
Bre-Gel BG 2254	Handwheel
Bre-Gel BG 3254, 338	Tap handwheel174, 198, 223-224, 262, 387
Gel duplicating	
Brevest exakta M256	haptosil D
Duro-Top259	Hard gingival mask
Generation M353, 416-417	Multisil-Mask, hard43
Giflex:	Hardener
Giflex-TR diamond disc30, 453-454	Gloss and hardening agent for plasters34
	Investment hardener259
Gingiva former	Hardware
Vario-Stud-Snap vks-oc rs 2.2 Abutments136-139	CAELO76-77
Gingival mask	e.cad77
Multisil-Mask hard43	e.order
Multisil-Mask, soft42	e.scan77
Silicone bur340, 362	Heat absorbent paste130
Gloss and hardening agent for plasters34	Hexagon head screwdriver
Gnathoflex	Screwdriver218, 410
Aesthetic-Gnathoflex102-104	Short screwdriver
Gnathoflex Study model FF1 mini105	
Gnathoflex Study model FF1 Set105	148, 152, 164, 174, 222-224, 228, 230, 232
Isoflex	Hexagon insert bit
Gnathoflex Study model FF1 mini105	Screwdriver
Gnathoflex Study model FF1 Set105	Hexagon socket
	Titanium screws164, 228, 230, 232
Gnathoflex-Premium105-107	Titanium screws, extended head228
Goat-hair brush353	Hexagon socket screws, screw connections
Golden booklet126, 260	Screw connection set, partly prefabricated 232
Granulate	High luster polisher
Polymer, Pi-Ku-Plast HP 3637, 241	Cerafine281
Grinder	High luster polishing
Diagen-Turbo-Grinder286, 353, 457, 459-460	Brepol486
Grinding	
Diabolo sintered diamonds463-470	Cotton mandrel
Diacryl grinder dcs457-458	High lustre
Diagen-Turbo-Grinder286, 353, 457, 459-460	Cotton buff for handpiece345, 478
Grinding tool	High lustre buff Acryl345, 479, 482, 484
	High lustre buff metal472, 475
Diamond grinding tools	High lustre paste Acrypol353
Set-up grinding tool331, 342, 457, 461	Acrypol polishing paste, high lustre345, 487
Grip	High lustre polishing
Pollygrip406, 488	Polishing cream390
Grip lining	Polishing granulate390
foam rubber grip lining89	Polishing stainless steel pins
Groove burs449	High lustre polishing with the handpiece345, 478
Guiding sleeves	
SKYplanX guiding sleeves59	HIGH.lign321
	top.lign professional
	uni.lign322-324
Н	uni.lign speed325
Handpiece for milling unit BF 2	High-gloss polishing paste
Chuck262, 387	Zi-polish283, 487
	High-performance polymers
Handpiece, polishing with the handpiece	BioHPP353-356, 381
346-348, 472, 477, 479, 484	for 2 press351-356, 381
Brushes for polishing with the handpiece	High-precision silicone for keys
Cotton buff for handpiece345, 478	visio.sil fix291, 302
Fabric buff478	•
Felt polishing buff478	Hinges
Pen-shaped brushes478	Wax bar hinges191
Rodeo326, 362	Holder
Handpiece, polishing with the handpiece, cotton buff	Articulation paper holder331
478	Posi-boy331
	Tap holder, tool set228
High lustre polishing with the handpiece345	

Holder for crowns and bridges	Implant planning system
Pollygrip406, 488	SKYplanX guiding sleeves59
Holder for discs	SKYplanX sleeve drill59
Quick mandrel394	SKYplanX sleeve seating tool59
Holder for polishers	Implant technique
Quick mandrel394	Multisil-Mask hard43
Holder for tap	Multisil-Mask, soft
Tap handwheel174, 198, 223-224, 262, 387	Impression
Hollow blade	breciform D impression tray8, 54
Blades for wax knife89, 386	
Hollow blocks	brecision implant heavy
	brecision implant light
Between hollow blocks	brecision impression materials8
Hollow pontic blocks, Auxiliary wax elements wht98	brecision putty soft56
Hollow pontics, Auxiliary wax elements wht101	Impression disinfectant
horizontal screw retentions	Dentaclean impression and denture disinfectant
Assembled transverse fixation221	9, 327
SKY fast & fixed221	Impression materials
SKY fast & fixed Laboratory analog69	brecision implant heavy54
Transverse fixation	brecision implant light54
Hot air device	brecision impression materials8
Thermo-Pen66, 294	brecision putty soft56
,	Impression taking
Hot-curing resin	
top.lign breformance liquid hot68	breciform D impression tray
Housing	Impression transfer set179
Metal matrix housing vks-oc144	Impression trays
Titanium matrix housing for acrylic178	breciform D impression tray8, 54
Titanium matrix housing for metal178	brecision implant heavy54
Titanium matrix housing vks-oc144	brecision implant light54
vks-sg matrix housing157	brecision impression materials8
Wax matrix housing zg178	brecision putty soft56
hpbl	Tray material UV13, 337
Hollow pontic blocks, Auxiliary wax elements wht98	In the articulator
hpg	Master-Split metal magnetic plate
Hollow pontics, Auxiliary wax elements wht101	In-between hollow pontics, auxiliary wax elements wht 97
	In-between pontics, auxiliary wax elements wht97
1	Individual impression trays
	Tray material UV13, 337
Ignitor	Individual locking pin
Piezo-Blitz393	Locking pin bs1198-199
Immediate implant	Individual overbedding
mini <sup>1</sup> SKY54, 56	Crêpe sleeve260
mini¹SKY laboratory analog56	Individual screw connection
mini <sup>1</sup> SKY planning matrix56	Tungsten carbide facing cutter228
Immersion (dipping) hardener	Individual screw connections
Duro-Top259	
	Tool set228-229, 231
Implant	Individually adjustable friction
Axle abutment	Vario Compress 1174-175
FGP friction fit system265-267	Information material
Implant analogue136	Patient brochures132
mini SKY OP-Tray54	Injection device362
mini <sup>1</sup> SKY54, 56	Injection moulding ceramic
mini <sup>1</sup> SKY laboratory analog56	Brevest ESG257
mini¹SKY planning matrix56	Injection moulding resin
O-Ring54	Thermopress357-362, 382-383
Implant	
Universal screwdriver set216, 408	Injection resin
Universal screwdriver set for contra-angles217, 409	Bio Dentaplast358-359, 361, 383
	Inlay
Vario-Stud-Snap vks-oc rs 2.2 Abutments136-139	Brealloy C + B 270
Implant analogue136	Inlay, Crown-bridge-inlay wax KBI92
Implant impression	Inserting instrument141, 144
brecision implant heavy54	Inserts
brecision implant light54	Quick Change90, 395
Implant planning system	Inserts
SKYplanX drilling sleeves	

Instruments		Investment model
Quick Change	90, 395	Stabilizer
Repositioning tweezers		Isobre wax insulating (separating) liquid
Insulating	,	Isoflex104-
Isoflex	104-105	Isoplast12, 59, 65, 69, 326, 333,
Multisil-Sep		Isoplast, separating plaster against resin
Plaster insulating liquid		12, 59, 65, 69, 326,
Wax insulating liquid		Isosil
Insulating (separating) liquid		100011
	۰ و	
Isobre wax insulating (separating) liquid		J
Isoplast	330	Joiners, Protek
Insulating liquid	00	Protek clasp/bar joiners
Exakto-Form insulating liquid	39	Protek clasp/bar joiners
Integrating/fixing		
Master-Split metal magnetic plate	26	K
Integrating/fixing in the articulator		K2 exact carring way
Master-Split base former(s)		K2 exact carving wax
Master-Split model system		K50, Abraso-Star K50, polishing paste
Model formers large/medium/small	26	K80, Abraso-Star K80, polishing paste
Interim dentures		KBI, Crown-bridge-inlay wax
Qu-base UV	64-65	Kneading silicone
Thermo-Pen		haptosil D44,
Thermoplastic clasp	,	Knife
top-lign breformance		Plaster knife29,
Interlock		KoEx Measuring Device17, 45,
Paralleling mandrel Interlock	13/	KoliBrush276,
Interlock 2°		Kolinski
Interlock 2		KoliBrush276,
	134	Kompaktopaker tooth-coloured UV
Intracoronal attachment	470	K-Primer291, 296-2
Inverto Plus		KS-lock
Inverto Plus		10 1001
Activating screw		
Basal screw		L
Ceramic spacer		Laboratory analog
Ceramic spacer E	204	mini <sup>1</sup> SKY laboratory analog
Investing		Last tap, tungsten carbide
Manual bredent casting technique	84	
Metal muffle rings	121	Latch spring
Silicone muffle rings	122, 353	Lead holder
Transfuser	125. 169	Brenometer lead holder
Investment aid		Leather buff
Investment hardener		High lustre polishing with the handpiece345,
Investment marker		Life-Color wax50
Investment material		Light-curing adhesive
Bresol for 2 press		SERACOLL UV59, 94, 2
		Light-curing die varnish32
Bresol M		diephos dentine
Bresol N		Light-curing resin
Bresol R		compoForm UV95,
Bresol Speed		Qu-base UV64
Brevest C+B Speed	124	Light-curing resin for blooking out
Brevest ESG	257	Litebloc UV
Brevest exakta M	256	
Brevest exakta Speed	256	Light-curing tray material
Brevest for 2 press		Tray material UV13,
Brevest M1		Light-curing unit
Brevest Rapid 1		bre.lux Power Unit65, 291, 299, 306-307,
	124, 200	Light-curing unit Polylux 2
Investment material for rapid heating up	050	Polylux 213
Brevest exakta Speed	256	Light-curing unit Polylux PL 20
Investment model	2=2	Polylux
Base former		Limitedly removeable constructions
Base plate		Security-Lock-System222-:
Base plate ring		Lineal
Magnetic plates		OPG-ruler
Silicone sleeve		OI O-10161

Liners		
Fleece liners		121
Wash primer for fleece		121
Lingual bar		
Protek wax pattern, lingual bar	225	226
	.233	-230
Link		
visio.link		297
Liquid		
Bresol for 2 press		353
Bresol M		
Bresol N		
Bresol R		
Bresol Speed		
Cleaner		241
Cleaner, Pi-Ku-Plast HP 36		37
Exakto-Form insulating liquid		39
Measuring cup		
Monomer, Pi-Ku-Plast HP 36		
Opaquer mixing liquid		2/1
Porcelain mixing liquid		271
Stain liquid		271
Thinner for ceramic separating liquid		
Thinner for Oxyd-Stop PM		129
Wax insulating liquid		222
		332
Lite		
Wax-Lite surface tension reducing agent		
Litebloc UV		31
Locating pin		
Brenometer locating pin		263
Location matrix drill		303
		302
Lock	040	~4~
Locking Pin activatable		
Locking Pin Easy-Snap		
Locking Pin Easy-Snap A		203
Locking Pin Easy-Snap E		
Locking Pin Snap A205,		
Locking Pin Snap E		
Swivel-type lock system sr		
Swivel-type lock system src	.196-	-197
Lock attachements		
KS-lock		200
Lock spring		
Lock techniques		
Locking pin bs1	100	100
Locking pin bs i	. 190-	-199
Locking Pin activatable	.212-	-213
Locking pin bs1	.198-	-199
Locking Pin Easy-Snap	201,	204
Locking Pin Easy-Snap A		
Locking Pin Easy-Snap E		
Locking Pin Snap A205,		
Locking Pin Snap E	.205	-210
Looking		
Matrix sleeves, Security-Lock-System		224
Looking ring, refractory alloy, cast-on		
Loop knife	,	
Abdruck-Cut	Ω	300
Lustre, high lustre paste Acrypol		
Acrypol polishing paste, high lustre	345,	487
LV 1		
LV 1 laser joint		243
LV 1 laser joint		

## М

M : D	074 075 00	
MagicBrush		
MagicContrast	274-27	5, 398-399
Magnet		
Magnetic plates		252
Quick Change		90, 395
Magnetic plate		•
Master-Split metal magnetic pla	ate	26
Magnetic plates		
Magnetic plates		252
Magnetic screwdriver		
Quick mandrel system		394
Main polishing		
Polishing porcelain spheres		390
Polishing powder, fine		
Mamelon cutter		
Mandibular advancing device		
		204 205
SleepPlus		364-365
Mandrel		
Cotton mandrel		
Paralleling mandrel Interlock		134
Paralleling mandrel metal vsp		184-185
Paralleling mandrel metal wstg		
142	145-146 14	8 15 <i>1</i> -156
Paralleling mandrel wstg		
Quick mandrel system		
Universal paralleling mandrel		
162,		
Universal paralleling mandrel 2		172, 179
Manipulating implant		
Implant analogue		136
Manual		
		126
Screwdriver		
Manual bredent casting technique	<del>)</del>	84
Manual processing of zirconium		
5-motion-active		
5-motion-colors		283
5-motion-glue		
5-motion-scanner		
5-motions-milling cutter		285
5-motions-zircon		200
		203
Margin grinding tool		4== 4=0
Diacryl grinder dcs		457-458
Marker		
Casting ring marker		126
Investment marker		126, 260
Mask – gingival mask		-,
Multisil-Mask hard		43
Multisil-Mask, soft		42
Master model system		
Master-Sep		19, 23
Master pin drill unit		20-21, 377
Master x-tray		
Diamond grinding disc Giflex T	R Master x-tra	av
	3	0 453-454
Master-Copy		252-253
Base former		
Base plate		
Base plate ring		
Magnetic plates		252
Silicone sleeve		252
Stabilizer		
Master-Pin		
Master-Pin Diatit tungsten carbide	etan drill '	23 20 22 277
master i in Platit turigateri carbiut	, otop utili z	_0, _0, 011

Master-Pin Radix-K		Metal polishing pastes	486-487
Tungsten carbide drill20, 202, 206	, 209, 344, 377	Metal polishing set	476
Master-Pin Radix-S	18	Metal polishing, brushes for metal po	lishing
Master-Pin sleeve	23	Metal polishing	472-473
Master-Pin system	22-24	Metal primer	
Master-Sep	19, 23	MKZ Primer	291, 296
Master-Split base former(s)	26	Metal transfer patrices vks	43, 141, 144, 154
Master-Split metal magnetic plate	26	Metal-ceramic blocks	
Master-Split model system		Biotec metal-ceramic blocks witho	ut collar96
Matrix material		Metal-ceramic blocks, Auxiliary wax e	elements wht99
visio.sil	291, 301	Microbur	
visio.sil ILT	291	Rapidy Microbur with relief	202, 206, 423
Matrix housing		Microflex	
Metal matrix housing vks-oc	144	Microkeramik	
Titanium matrix housing for acrylic		Milling	
Titanium matrix housing for metal		Brealloy C + B 270	127
Titanium matrix housing vks-oc		Milling and drilling oil	
vks-sg matrix housing		202, 206, 209, 219, 22	
Wax matrix housing zg		Milling base	
Matrix inserting instrument137,		Milling technique	200, 300
		5-motion-scanner	284-285
Matrix material	7, 10 <del>4</del> -105, 100	5-motions-milling cutter	
	44 204		
haptosil D		breCAM.BioHPP	
Matrix pliers		breCAM.resin	
Matrix sleeves, Security-Lock-System	224	breCAM.resin milling blanks	
Matrixes		breCAM.wax	
Duplicating matrixes zg		Conical bur	
mini¹SKY planning matrix		Cross-cut burs	
Matt sealing agent for epitheses	374	Generation M	,
coating		Groove burs	
matting powder	374	Milling and drilling oil	
primer	374	202, 206, 209, 219, 22	2-224, 228, 231, 449
Matting powder	374	Milling base	263, 388
Measuring abutment		Milling unit BF 2	262, 387-388
Axle abutment	136	Parallel bur	442-445
Measuring cup	338	Shoulder bur	449
Measuring device		Tools for Milling technique	440
Angle measuring device	136	Transfer device	
KoEx Measuring Device		Wax bur 440, 442-44	
Measuring instrument		Milling tools	o,,
Angle measuring device	136	Generation M	353, 416-417
Measuring syringes		Milling unit BF 2	
Metal		Chuck	
Abraso-Buff Metal	472 474	Model support BF1	
Abraso-Buff Metal Mini		Milling wax	202 200, 007 000
Abraso-Buff Polipast Metal		Biotec milling wax	93 441
		Mini	
Abraso-Soft Metal		Abraso-Buff Metal Mini	470 474
High lustre buff metal		Double t-adhesive mini connector.	246 247
Metal muffle rings	121		
Metal attachment	470	Vario-Soft 3 mini	
Activating screw		Vario-Soft 3 mini sv	
Basal screw		mini SKY OP-Tray	
Ceramic spacer		mini¹SKY	
Ceramic spacer E		mini¹SKY laboratory analog	
Inverto Plus	172	mini <sup>1</sup> SKY planning matrix	56
Metal casting ring		Mixing	
Casting ring marker	126	ecovac Vacuum Mixing unit	
Metal connection (joint)		Mixing block	32, 267
LV 1 laser joint	243	Mixing cannulas	
Metal magnetic plate		mixing cup1	
Master-Split metal magnetic plate	26	Mixing liquid	
Metal matrix housing vks-oc		Bresol for 2 press	353
Metal muffle rings		Bresol M	
Metal polishing		Bresol N	

Mixing liquid	Modelling
Bresol R124, 2	
Bresol Speed124, 2	
Monomer, Pi-Ku-Plast HP 3637, 2	
Porcelain mixing liquid2	
Mixing plates279-2	
Mixing system	Wax bur
base16, 123, 3	
Mixing unit	Modelling crowns
ecovac Vacuum Mixing unit16, 38, 123, 3	
Mixing-stirring spiral	
MK	Modelling instruments
Brealloy MK1	
mkbl	Modelling of the crown margin
Metal-ceramic blocks, Auxiliary wax elements wht	
MKZ EM-Aktivator	
MKZ Primer	
mmg Matal matrix hausing uka as	Modelling secondary elements
Metal matrix housing vks-oc1	
MO Breelley MO	Wax for outer copings9
Brealloy MO2	
Model fabrication	Cervical disc
Adapter base20, 3	• .
Base former2	
Base plate2	
Base plate ring2	52 Gecko
Exakto-Form, model resin	
Exakto-Rock S17, 38, 74,	
Functional margin wax	36 Modelling wax for crown margin
Magnetic plates2	
Master-Pin	
Master pin drill unit20-21, 3	
Master-Pin Diatit tungsten carbide step drill	
20, 23, 3	
Master-Pin Radix-K	
Master-Pin Radix-S	
Master-Pin sleeve	
Master-Sep19,	
Master-Pin system22-	
Plaster insulating liquid3	
Retention pins	
Silicone sleeve2	52 Mould release agent120-122
Spacer wax2	
Stabilizer2	52 Diatit-Multidrill with stop228
Model formers large/medium/small	
Model hardener	198, 202, 206, 209, 223-224, 228, 362
Investment hardener2	59 Multisil
Model holder	Modelling wax for Epithetiks37
Model support BF 2 262-263, 387-3	88 Multisil fibers373
Posi-boy3	
Model resin	Multisil stains372
Exakto-Form, model resin	
Model support BF 2262-263, 387-3	
Model system	Multisil-Epithetik country37
Master-Split base former(s)	
Master-Split metal magnetic plate	
Master-Split model system25-27,	
Model formers large/medium/small	
Model table	Multisil Epithetik
Model support BF 2 262-263, 387-3	88 Multisil Epithetik set
Model transfer system	Multisil fibers
Master-Copy252-2	
Modelling	Multisil sealing agent373
Blades for wax knife89, 3	
Diagos for wax will	oo iyaation oodiing liquid

Multisil soft	Opaquer set305
Dosing device56-57, 59, 340	
Multisil stains372	OPG
Multisil-Epithetik beach371	OPG-ruler52
Multisil-Epithetik city371	X-resin OPG52
Multisil-Epithetik country371	OPG-ruler52
Multisil-Epithetik hard-form370	Opti-Cast
Multisil-Epithetik soft-form370	Measuring cup338
Multisil-Epithetik thickener372	Mixing cup
Multisil-Epithetik transparent371	Optiguss
Multisil-Mask hard43	Optiguss mixing well macro/micro108, 248
Multisil-Mask, soft42	OP-Tray
Multisil-Primer340, 372	mini SKY OP-Tray54
Multisil-Sep42	O-Ring54
Multisil-Soft340	Orthodontics
	Silicone bur
N	Overbedding
IN .	Crêpe sleeve260
N 15	Oxyd-Stop PM
Exaktosil N 15249	Oxyd-Stop Silver-palladium130
N 21	Oxyd-Stop-NPM
Exaktosil N 21249, 362	Oxyd Glop 141 W
Natural hair brush	D.
KoliBrush276, 400	Р
neo.lign56, 65, 69, 320	Palatal ridges
neo.lign A295, 304, 312-313	Wax palatal patterns
neo.lign P295, 304, 314-315	Palladium
Neutralizing agent	Oxyd-Stop Silver-palladium130
Wax-Lite surface tension reducing agent10, 120	Paper holder
Non-precious metal alloys	Articulation paper holder
Brealloy MK127	Papilla grinding tool
Diatit power burs437	Diacryl grinder dcs457-458
Non-precious metal polishing paste	Diacryl grinding tool
Brepol486	Parallel
novo.lign A294, 304, 308-309	Interlock parallel134
novo.lign P294, 296, 304, 310-311	Parallel attachment
Novo-Grip264, 405	Cylindrical attachment178-181
Inserts	Parallel bur442-445
NPM alloy	Paralleling mandrel
Brealloy C + B 270127	Milling unit BF 2262, 387-388
Brealloy F 400	
Dieality F 400230	Paralleling mandrel Interlock
	Paralleling mandrel for oc and custom bridge sectioning
0	attachments
oc	Paralleling mandrel Interlock
Bridge sectioning attachment oc230	Paralleling mandrel metal vsp
Vario-Stud-Snap vks-oc 135-141, 144-146, 148-149	Paralleling mandrel metal wstg
Vario-Stud-Snap vks-oc exchangeable stud with	142, 145-146, 148, 154-156
adhesive sleeve150	Paralleling mandrel vks
Occlusal bridge sectioning attachment	Paralleling mandrel wstg190-191
	Partly prefabricated
Bridge sectioning attachment oc	Screw connection set, partly prefabricated
Occlusal surface tool Fissure tool418	Threaded bush, refractory alloy, cast-on
	Paste
Occlusal surfaces	Abraso-Star K50345
Aesthetic-Gnathoflex102-104	Heat absorbent paste
Occlusals (occlusal surfaces)	Paste, high lustre paste Acrypol353
Gnathoflex-Premium105-107	Acrypol polishing paste, high lustre345, 487
Offset screwdriver	Patient brochures
Cordless Prosthodontic Screwdriver384	Patrices
Oil	Metal transfer patrices vks43, 141, 144, 154
Milling and drilling oil	Pattern
202, 206, 209, 219, 222-224, 228, 231, 449	Wax pattern sticks116
Opaque die varnish, light-curing32	PEEK
Opaquer mixing liquid271	breCAM.BioHPP79
Porcelain liquid set 271	

Pen-shaped brushes478	Pliers for conical crowns	
Perlform	Inserts	
Dipping wax87	Novo-Grip	264, 405
Permanent restoration	Plug-type pin	
top.lign professional326	Master-Pin Radix-S	18
Pickling	PM	
Brecid pickling agent130	Oxyd-Stop PM	129
Piezo-Blitz393	PML surveyor	
Pi-Ku-Plast HP 36 ppk37, 241-242	PMMA-Blanks	,
Pi-Ku-Plast ppk241-242	breCAM.resin	80
Pi-Ku-Plast separating varnish37	PMP surveyor	
Pin	Polierjet Duo-Finish	
	Polierjet Quadro-Finish	
Inserting instrument		390
Master-Pin	Polipast	470 475
Master-Pin Radix-K	Abraso-Buff Polipast Metal	472, 475
Master-Pin Radix-S18	Polisher for resin	
Master-Pin sleeve23	Diacryl rubber grinding tool	
Master-Pin system22-24	Polisher, high lustre paste Acypol	
Matrix inserting instrument 137, 146, 154, 156,	Acrypol polishing paste, high lustre	345, 487
162, 166-167, 179, 184-185, 188	Polishing	
Pin Diatit tungsten carbide Pin	Abraso-Sil Acrylic	479-480
Master-Pin Diatit tungsten carbide step drill	Abraso-Starglanz	
20, 23, 377	Acrylic polishing set	
Pin drill unit	Ceragum	
adapter base20, 377	Cotton buff for handpiece	
Master pin drill unit20-21, 377	Cotton mandrel	
Pins	Fabric buff	
Retention pins	Leather buff3	
PL 20 light-curing unit, Polylux PL 20	Metal polishing set	
Polylux380	Pen-shaped brushes	
Planning matrixes	Polierjet Duo-Finish	
mini <sup>1</sup> SKY planning matrix56	Polierjet Quadro-Finish	390
Planning system	Polishing acrylic	479-480, 482-483
SKYplanX drilling sleeves59	Polishing bur	440, 442-447
SKYplanX guiding sleeves59	Polishing porcelain	390, 485, 487
SKYplanX sleeve drill59	Pollygrip	
SKYplanX sleeve seating tool59	Rodeo	
Plaster	Polishing acrylic	
Arti-Rock28	Polishing brush	170 100, 102 100
Fluid-Rock	Abraso-Buff Acryl3	2/5 3/7 /70 /81
Plaster against wax	Abraso-Buff Metal	
Wax insulating liquid332	Abraso-Buff Metal Mini	
Plaster bur	Abraso-Buff Polipast Metal	
Tungsten carbide burs for processing of plaster 31	Abraso-Soft Acryl	
Tungsten carbide tools with relief 419-422, 424-425	Abraso-Soft Metal	
Plaster insulating liquid332, 362	High lustre buff Acryl3	345, 479, 482, 484
Plaster knife29, 393	Prepolishing Buff Silicone	481
Plaster sealing liquid272	Polishing brushes for metal polishing	
Plaster separating agent (liquid)	Metal polishing	472-473
Master-Sep19, 23	Polishing brushes for polishing with the	
Plaster solvent	Brushes for polishing with the handpi	
Dentaclean plaster removing agent29, 328	346-348, 4	
	Doliobing bur	112, 411, 419, 404
Dentaclean plaster removing agent Speed328	Polishing bur	
Plaster-resin separating agent/liquid	Parallel bur	
Acrylic Sep362	Polishing cream	390
Plastic pin	Polishing drum	
Master-Pin Radix-K19	Polishing granulate	390
Master-Pin Radix-S18	Polishing of ceramic	
Plate	Abraso-Fix	281. 485
Master-Split metal magnetic plate26	Polishing of metal	5., .50
Pliers	Abraso-Fix	281 ⊿85
Inserts264, 405	Polishing paste	201, 400
		100
Matrix pliers	Brepol	
Novo-Grip264, 405	Pumice polishing paste	345

Polishing paste	Prefabricated wax elements
Zi-polish283, 487	Protek clasp patterns237-238
Polishing paste, diamond polishing paste	Prehigh lustre
Polishing porcelain390, 485, 487	Rodeo326, 362
Polishing pastes, Abraso-Star K50, Titapol, Abraso-	Premium
Starglanz, pumice polishing paste486-487	Gnathoflex-Premium105-107
Polishing porcelain390, 485, 487	Premolar clasp
Polishing porcelain spheres390	Protek clasp patterns237-238
Polishing porcelain, felt wheels, diamond polishing paste	Preparation margin
Polishing porcelain390, 485, 487	Gloss and hardening agent for plasters34
Polishing powder	Preparation of a thread
Polishing powder, coarse390	First tap, tungsten carbide198, 223, 228
Polishing powder, coarse390	Last tap, tungsten carbide228
Polishing powder, fine390	Second tap, tungsten carbide223
Polishing stainless steel pins390	Preparatory work
Polishing to high lustre	adapter base
Abraso-Starglanz345, 353	Master pin drill unit
Polishing with the handpiece	Plaster sealing liquid
Brepol486	Pre-polish
Pollygrip406, 488	Polierjet Quadro-Finish390
	•
Pollygrip set	Prepolishing
Pollygrip	Polishing powder, coarse
Polyan - Accessories	Prepolishing Buff Silicone
Acrylic Sep362	Press-in pin
Polylux380	Matrix inserting instrument
Polylux 213-14	162, 166-167, 179, 184-185, 188
Polylux, light-curing unit, Polylux 2	Pressing matrixes into the housing
Polylux 213-14	Inserting instrument141, 144
Polylux, light-curing unit, Polylux PL 20	Press-System
Polylux380	BioHPP353-356, 381
Polymer, Pi-Ku-Plast HP 3637, 241	for 2 press351-356, 381
Polymerization unit	Primer374
bre.lux Power Unit65, 291, 299, 306-307, 379	K-Primer291, 296-297
Polymerization unit, light-curing unit	MKZ Primer291, 296
Polylux380	Primer
Polylux 213-14	Multisil-Primer372
Pontic	Processing of acrylics
Aesthetic and ergonomic pontics for ceramic, Auxiliary	Abraso-Gum Acryl
wax elements wht100	Abraso-Sil Acrylic479-480
Pontics with shallow collars for ceramic, Auxiliary wax	Diacryl grinder dcs457-458
elements wht99	Tungsten carbide tools with relief419-422, 424-425
Porcelain liquid set271	Processing of ceramic
Porcelain mixing liquid271	breCeram281
Porcelain liquid set271	Cerafine281
Posi-boy331	Diagen-Turbo-Grinder286, 353, 457, 459-460
Positioning tweezers	Diabolo sintered diamonds463-470
Repositioning tweezers330, 404	Tungsten carbide tools with relief419-422, 424-425
Posterior teeth	Processing of metal
novo.lign P294, 296, 304, 310-311	Tungsten carbide tools with relief419-422, 424-425
visio.lign	Processing of resin
Powder	visio.lign Toolkit291, 305
Polymer, Pi-Ku-Plast HP 3637, 241	Processing of resin (acrylics)
opk Classer	Diacryl grinding tool
Cleaner Di Ku Pleat HP 26	Processing of titanium
Cleaner, Pi-Ku-Plast HP 3637	Processing of zirconium
Pi-Ku-Plast HP 36 ppk37, 241-242	Airaqua Turbine389
Pi-Ku-Plast ppk241-242	Processing set
Vessels ppk241	breCeram
Precious metal	Cerafine
Oxyd-Stop PM129	Processing VMK structure
Precision investment material	Tungsten carbide tools with relief419-422, 424-425
Brevest M1255	Profile
Prefabricated	Quadrosticks117
Screw connection set, partly prefabricated232	

Profile bur	440, 442-448	Quick mandrel, Quick mandrel system.	394
Conical bur	446-448	Quick mandrel system	
Parallel bur	442-445	Quicktool	273, 397
Profile bur, Milling technique		Qu-resin	. 65, 70, 202, 362
Profile bur	440, 442-448		
Prosthetic cap		R	
SKY fast & fixed Prosthetic cap	69, 221		
Prosthetic instruments		Rack	
Quick Change	90, 395	Foot rack	390
Prosthetics		Radiopaque resin	
Acrylic Sep	362	OPG-ruler	
Articulation paper holder		X-connector bonding agent	55, 58
Bio Dentaplast	358-359, 361, 383	X-resin	55-56
Dentasil tooth protection silicone		X-resin flow	57
Ergonom wax knife		X-resin OPG	52
Isoplast		Radix-K	
Kompaktopaker tooth-coloured UV		Master-Pin Radix-K	19
Measuring cup		Radix-S	
Mixing cup		Master-Pin Radix-S	18
Multisil sealing liquid		Rapid	
Multisil-Soft		Brevest Rapid 1	124, 256
Ropak compact opaquer		Rapid-heating investment material	
Ropak UV powder + liquid		Brevest Rapid 1	124, 256
Set-up grinding tool		Rapidy Microbur with relief	
Set-up wax		Red-white esthetics set	
Wax palatal patterns		Reels of wax patterns	
Protective chamber with/without extra		Refractory alloy, cast-on	-,
Protective chamber with without extract Protective paste	CHOT HOZZIE 300	Looking ring, refractory alloy, cast-on	164, 232
Heat absorbent paste	130	Relief	,
Protective wax	130	Conical bur	446-448
Functional margin wax	36	Groove burs	
Protek		Parallel bur	
		Profile bur	
Bonyhard clasp		Rapidy Microbur with relief	
Circumferential clasp, bent		Shoulder bur	
Molar clasp		Tools for Milling technique	
Protek clasp patterns		Relining	
Protek clasp/bar joiners		Multisil sealing liquid	340
Protek sculpturing wax Protek wax adhesive		Multisil-Soft	
		Removable restorations	
Protek wax pattern		Locking Pin activatable	212-213
Protek wax pattern, lingual bar Protek, wax patterns cut to size	235-236	Locking Pin Easy-Snap	
	240	Locking Pin Easy-Snap A	
Wax patterns cut to size, Protek Pumice disinfectant	240	Locking Pin Easy-Snap E	
	220 400	Locking Pin Snap A	
Dentaclean pumice disinfectant		Locking Pin Snap E	
Pumice polishing paste	345	Removal of a matrix	205-210
Pumice polishing paste, polishing pas		Matrix pliers	137-138 170
Metal polishing pastes	486-487	Removing metal	137-130, 173
Putty	50	Profile bur	440 442 449
brecision putty soft	56	FIGHE but	440, 442-440
		Demoving of plactor	
Q		Removing of plaster	20
Quadro		Dentaclean plaster removing agent	
	300	Dentaclean plaster removing agent S	peeu320
Polierjet Quadro-Finish		Removing oxide	400
Quadro wax profile		Brecid pickling agent	130
Quadrosticks		Removing tartar	207
Quadro wax profile		Dentaclean denture cleaning agent	
Quadrosticks		Repositioning tweezers	330, 404
Qu-base UV		Resin	50.50
Quick Change	90, 395	3D-resin	
Quick change system	00 005	Casting system	
Quick Change		compoForm UV	
Quick cradle (handle), Quick mandrel	system394	Exakto-Form, model resin	39

Resin		Rotating tools
HIGH.lign	321	Cross-cut burs450
Litebloc UV		Diabolo sintered diamonds463-470
Pi-Ku-Plast HP 36 ppk		Diacryl grinder dcs457-458
Pi-Ku-Plast ppk		Diacryl grinding tool343
Qu-resin6	5, 70, 202, 362	Diagen-Turbo-Grinder286, 353, 457, 459-460
Thermopress		Diatit and tungsten carbide tools423, 426-439
top.lign breformance liquid cold		Diatit power burs437
top.lign breformance liquid hot		Elastic453, 456
top-lign breformance		Fissure tool418
top.lign professional		Flexible453, 456
uni.lign		Groove burs449
uni.lign speed		Location matrix drill
Resin blanks		Parallel bur442-445
breCAM.BioHPP	79	Polishing bur440, 442-447
breCAM.cutter		Processing of titanium438-439
breCAM.resin		Profile bur440, 442-448
Resin insulating (separating) liquid		Rapidy Microbur with relief202, 206, 423
Isoplast	338	Shoulder bur449
Resin press system		Superflex
BioHPP	353-356 381	Tools for Milling technique430
for 2 press		Transflex453, 455
Resin to be applied with a brush	331-330, 301	Tungsten carbide tools with relief419-422, 424-425
	244 242	•
Pi-Ku-Plast HP 36 ppk		Ultraflex
Pi-Ku-Plast ppk		Round blanks
Rest for handpiece	09	
Restoring the friction	204 405	breCAM.wax78
Activating pliers	204, 405	Round brush
Retaining element for acrylic dentures	470	Abraso-Fix
Titanium matrix housing for acrylic	178	Round face, conical bur, Milling technique
Retaining element for plastic matrixes	4.4.4	Conical bur
Metal matrix housing vks-oc		Round face, parallel bur, Milling technique
Titanium matrix housing vks-oc		Parallel bur442-445
Retention adhesive		Rubber grinding tool
Retention Crystals and Beads		Diacryl grinder dcs457-458
Retention pins	36	Diacryl rubber grinding tool
Retentions	0.007.000.040	Rubber polisher
Protek	6-237, 239-240	Cerafine
Retention Crystals and Beads	109	Rubber polisher for resin
Ridges	000	Diacryl rubber grinding tool343
Wax palatal patterns	336	
Right-angle, screwdriver for right-angle	400 040 440	<b>S</b>
Screwdriver	136, 218, 410	Cayyout madel
Ring		Sawcut model
Looking ring, refractory alloy, cast-on	164, 232	adapter base
Ring clasp		Master pin drill unit20-21, 377
Protek clasp patterns		Master-Pin23
Rinsing heads		Master-Pin sleeve
Rinsing heads for voluminous castings	114	Master-Pin system22-24
Rod attachments		Sawcut models, diamond disc
Wax rod attachments	190	Giflex-TR diamond disc30, 453-454
Rodeo	326, 362	Scanner
Root cap anchor		5-motion-scanner
Vario-Stud-Snap vks-oc	142-143	5-motions-milling cutter285
Root caps		e.scan77
Vario-Stud-Snap vks-oc	142-143	Screw
Ropak compact opaquer	341	Activating screw172
Ropak UV liquid		Basal screw172
Ropak UV powder		Fixing screw, bridge sectioning attachment oc 230
		Fixing screw, Screw connection set232
Rotary instruments		Security-Lock-System222-225
Generation M	353, 416-417	Threaded rods223-224
Rotating tools		Titanium adjusting screw174, 176
Conical bur	446-448	

Screw connection	See above
Ceramic screws174, 223	Security-Lock-adhesive sleeve224-225
Friction Splint FS1226-227	Self-adhesive patterns
Friction Splint FS1 sleeve226	Protek236-237, 239-240
Friction Splint FS1 Spacer226	Separarting
Friction Splint FS1 Splint226	Plaster sealing liquid272
Security-Lock-Ceramic223	Separate plaster against resin 12, 59, 65, 69, 326, 333
Security-Lock-System222-225	Separating
Threaded rods223-224	Acrylic Sep362
Titanium screws164, 228, 230, 232	Ceramic separating liquid272
Titanium screws, extended head228	Ceramic separating set272
Universal screwdriver set216, 408	Separating agent (liquid)
Universal screwdriver set for contra-angles 217, 409	Master-Sep19, 23
Screw connection set, partly prefabricated232	Separating discs
Screw retentions	Quick mandrel394
Assembled transverse fixation221	Separating liquid
SKY fast & fixed221	Exakto-Form insulating liquid39
SKY fast & fixed Laboratory analog69	Separating liquid for plaster
Transverse fixation220	Plaster insulating liquid332, 362
Screwdriver136	Separating liquid for resin, Isoplast, separating plaster
Cordless Prosthodontic Screwdriver384	against resin
Short screwdriver	Separating plaster against resin 12, 59, 65, 69, 326, 333
148, 152, 164, 174, 222-224, 228, 230, 232	SERACOLL UV59, 94, 282
Screwdriver (long, short, right-angle)218, 410	Set
Screwdriver set	Acrylic polishing set345
Universal screwdriver set216, 408	Composite finishing set299
Universal screwdriver set for contra-angles 217, 409	Impression transfer set179
Screws	neo.lign56, 65, 69, 320
Universal screwdriver set216, 408	neo.lign A295, 304, 312-313
Universal screwdriver set for contra-angles217, 409	neo.lign P295, 304, 314-315
Sculpturing	Opaquer set305
Protek sculpturing wax235	Red-white esthetics set303
Sealing (liquid)	Screw connection set, partly prefabricated 232
Multisil sealing liquid340	Universal screwdriver set216, 408
Sealing agent	Universal screwdriver set for contra-angles 217, 409
Matt sealing agent for epitheses374	Metal polishing set476
Multisil sealing agent372	Set, Pollygrip set
Sealing margin	Pollygrip406, 488
Protek	Set-up
Sealing of surfaces	Articulation paper holder331
Optiguss108, 248	beauty setup50
Sealing plugs	Set-up grinding tool
Opti-Cast338	Set-up grinding tool
Second tap, tungsten carbide223	Set-up wax51, 335
Sectioning attachment	sg
Bridge sectioning attachment oc230	Vario-Stud-Snap vks-sg151-152, 154-158
Custom bridge sectioning attachment231	Vario-Stud-Snap vks-sg exchangeable stud with
Vario-Soft 3 sv bridge-sectioning attachment164-165	adhesive sleeve150, 153
Sectioning of bridges	Shipping bags
Security-Lock-System222-225	Short screwdriver
Sectioning stud	148, 152, 164, 174, 222-224, 228, 230, 232
Bridge sectioning stud oc, bridge sectioning attachment	Shoulder bur449
oc230	Sil
Security-bite blue48	Abraso-Sil Acrylic479-480
Security-Lock	Silencos night device367
Auxiliary modelling198, 223-224, 228	Silicon sealing agent
Ceramic removing tool174, 223	Matt sealing agent for epitheses374
Ceramic screws	Silicone
Threaded rods	Block-out kneading material251
Security-Lock-adhesive sleeve224-225	Duplicating system251
Security-Lock-Ceramic223	Exaktosil silicone duplicating materials249
Security-Lock-System222-225	Exaktosil N 15
Matrix sleeves, Security-Lock-System224	Exaktosil N 21
3.33.33, 3334.Nj 200N Ojokolili	Flask sleeve

Silicone		Sleep apnea	
Flask tray	251	Biofunctional therapy	366-36
Friction silicone	174, 176	silencos night device	
haptosil D		SleepPlus	
Investment aid		Vacuum activator silencos	36
Isosil		Vacuum activator silencos kids.	
Location matrix drill		SleepPlus	
Material		Sleeve	
visio.sil ILT	301	Crêpe sleeve	26
		Friction Splint FS1 sleeve	
Multisil Epithetik			
Multisil-Epithetik beach		Master-Pin sleeve	
Multisil-Epithetik city		Silicone sleeve	25
Multisil-Epithetik country		Sleeve drill	_
Multisil Epithetik set		SKYplanX sleeve drill	5
Multisil-Epithetik transparent		Sleeve seating tool	
Multisil-Mask hard		SKYplanX sleeve seating tool	5
Multisil-Mask, soft	42	Small punching tube	
Multisil-Primer	340	Opti-Cast	33
Multisil-Soft	340	Smallest fissure tool	
Prepolishing Buff Silicone		Fissure tool	41
Silicone muffle rings		Smoothing	
Spacer – base insert		Investment hardener	25
Stabilizer		Smoothing of wax	20
Technolit		Optiguss	100 24
		. •	106, 24
Technosil duplicating silicone		Snap	105 110 110 11
visio.sil		Vario-Stud-Snap vks-oc	
visio.sil fix		Vario-Stud-Snap vks-oc exchange	
visio.sil ILT		adhesive sleeve	
Silicone and wax surface tension reducin	g agent .10, 120	Vario-Stud-Snap vks-sg	151-152, 154-15
Silicone brush		Snap	
Transfuser	125, 169	Vario-Stud-Snap vks-sg exchan	geable stud with
Silicone bur	340, 362	adhesive sleeve	
Silicone colors	,	Snap A	, -
Multisil intensive colors	373	Locking Pin Snap A	205 209 21
Silicone duplicating material		Snap E	
Exaktosil N 15	2/10	Locking Pin Snap E	205-21
Exaktosii N 21			203-21
		Snoring therapy	200 20
Exaktosil silicone duplicating materials	249	Biofunctional therapy	
Silicone elements for occlusal surfaces	100.101	silencos night device	
Aesthetic-Gnathoflex	102-104	SleepPlus	
Silicone for keys		Vacuum activator silencos	
visio.sil fix		Vacuum activator silencos kids.	36
visio.sil ILT	301	Soft	
Silicone muffle rings	122, 353	Abraso-Soft Acryl	345-346, 479, 48
Silicone plugs	,	Vario-Soft 3	160-16
Opti-Cast	338	Vario-Soft 3 matrix housing	
Silicone polisher		Vario-Soft 3 sv	
Ceragum	204 252	Vario-Soft 3 sv bridge-sectioning	
Silicone sleeve	252	Vario-Soft attachment group	
Silver		Vario-Soft-Bar vss	
Oxyd-Stop Silver-palladium	130	Vario-Soft-Bar-Pattern vsp	184-18
Sintered diamond tips		Soft gingival mask	
Inserts	264, 405	Multisil-Mask, soft	4
Sintered diamonds		Soft relining	
Diabolo sintered diamonds	463-470	Silicone bur	340. 36
SKY fast & fixed		Software	
SKY fast & fixed Laboratory analog		CAELO	7 <b>6</b> ₋7
SKY fast & fixed Prosthetic cap		e.cad	
SKY X-ray template		e.order	
SKYplanX drilling sleeves		e.scan	7
SKYplanX guiding sleeves		Solder	
SKYplanX sleeve drill		Brealloy solder	25
SKYplanX sleeve seating tool		Soldering	
. 3		Brealloy flux	128. 25

Soldering	Stains
Brealloy solder259	5-motion-colors
Heat absorbent paste130	Multisil intensive colors
Oxyd-Stop PM129	Multisil stains372
Oxyd-Stop Silver-palladium130	Stain liquid271
Solving plaster	Standard blade
Dentaclean plaster removing agent328	Blades for wax knife
Spacer	Standard Modelling wax91
Ceramic spacer172, 176	Starglanz
Ceramic spacer E204	Abraso-Starglanz345, 353
Friction Splint FS1 Spacer226	Statik-Disc234
Spacer – base insert251	Steel pins
Spacer (space retaining element)	Sticks
Ceramic spacer172	Quadrosticks117
Ceramic spacer E204	Wax pattern sticks116
Spacer varnish, air-drying33	Stippled
Spacer wax235	Stippled wax sheets
Spatulas267	Stippled wax sheets
Special diamonds for the veneering technique 457, 461	Protek236-237, 239-240
Speed	Stirrer
Bresol Speed124, 256	ecovac Vacuum Mixing unit16, 38, 123, 376
Brevest exakta Speed256	Mixing cup
Speed investment material	Mixing-stirring spiral
Brevest C+B Speed124	Stirring sticks
Speed investment material(s)	Stirring system
Brevest ESG257	base16, 123, 376
Spiral	Stone
mixing-stirring spiral16, 123, 376	Exakto-Rock S
Splendido92	Fluid-Rock
Splint	Stone for arches
Friction Splint FS1226-227	Exakto-Rock S
Friction Splint FS1 sleeve226	Stop
Friction Splint FS1 Spacer226	Oxyd-Stop PM129
Friction Splint FS1 Splint226	Straight face, conical bur, Milling technique
Security-Lock-Ceramic223	Conical bur446-448
Split-Cast	Straight face, parallel bur, Milling technique
Master-Split base former(s)26	Parallel bur442-445
Master-Split metal magnetic plate26	Stress-free bridges
Master-Split model system25-27, 40	SERACOLL UV59, 94, 282
Model formers large/medium/small26	Stud
Spot Clip396	Vario-Stud-Snap vks-oc135-146, 148-149
Spot Clip with supporting ring396	Vario-Stud-Snap vks-oc exchangeable stud with
Spray	adhesive sleeve150
Investment hardener259	Vario-Stud-Snap vks-sg151-152, 154-158
Spray bottle, plastic10, 120, 332	Vario-Stud-Snap vks-sg exchangeable stud with
Sprue wax	adhesive sleeve
Flexible Acrylic sprue wax362	Stud attachment
Sprueing	Cylindrical attachment178-181
Reels of wax patterns116, 239	Stud fixator177
Sprueing of wax-ups	Stud snap
Wax patterns cut to size, Protek240	Vario-Stud-Snap vks-oc 135-141, 144-146, 148-149
Sprues for centrifugal casting113	Vario-Stud-Snap vks-oc exchangeable stud with
Sprues for centrifugal casting for voluminous castings113	adhesive sleeve150
Sprues for vacuum pressure casting112	Vario-Stud-Snap vks-oc rs 2.2 Abutments 136-139
Sprues for vacuum pressure casting for voluminous	Super hard stone
castings112	Exakto-Rock S
sr	Fluid-Rock
Swivel-type lock system sr194-195	Super-flat
src	Double t-adhesive mini connector246-247
Swivel-type lock system src196-197	Superflex453, 456
Stabilizer251-252	Supporting Ring
Stain liquid271	Spot Clip with supporting ring396
Porcelain liquid set271	Surface conditioning agent
	Multisil-Primer340

Surface detensionizer		Teeth	
Technolit	250	visio.sil ILT	291
Surface treatment		Telescopes	
Abraso-Gum Acryl	342, 362	FGP friction fit system	265-267
Ceragum		Telescopic and conical crowns	
Diatit power burs		Biotec-Wax for outer copings	91
Fissure tool		Brenometer locating pin	263
Gloss and hardening agent for plas		Brevest Rapid 1	
Investment hardener		Cross-cut burs	
Polierjet Duo-Finish		Double t-adhesive connector	
Polierjet Quadro-Finish			
	390	LV 1 laser joint	
Surface treatment (processing)	0.40	Milling unit BF 2	
Multisil sealing liquid		Pen-shaped brushes	
Plaster sealing liquid	272	Pi-Ku-Plast HP 36 ppk	
Surveying system		Pi-Ku-Plast ppk	
Brenometer lead holder		Polishing bur	
Brenometer locating pin		Profile bur	
Brenometer plate	263	Transfer device	263, 388
Surveying system for CoCr clasp supp	oly	Wax bur 440, 442-443	, 446-447, 450-452
Brenometer surveying system		Wax for outer copings	
Surveyor	,	Telescopic crowns	
PML surveyor	262.387	Biotec blocking out wax	65. 235
PMP surveyor		Biotec milling wax	
SV	202, 307	Parallel bur	
Vario-Soft 3 mini	166		
		Vario Compress 2	170
Vario-Soft 3 mini sv		Telescopic pliers	004 405
Vario-Soft 3 sv		Inserts	
Swivel-type lock system sr		Novo-Grip	264, 405
Swivel-type lock system src	196-197	Template	
Sword blade		SKY X-ray template	52
Blades for wax knife	89, 386	Temporaries	
Synthetic hair		top-lign breformance	68-69
MagicContrast	. 274-275, 398-399	top.lign professional	
Syringe	,	Temporary crown and bridge material	
Measuring syringes	39	top.lign breformance liquid cold	
Thermo-syringe		top.lign breformance liquid hot	
Systhetic hair	20, 000, 002	Temporary implant	
MagicBrush 274	275 200 209 200	mini¹SKY	54 56
Magicbrusii 214	-273, 299, 390-399	mini¹SKY laboratory analog	
_			
T		mini¹SKY planning matrix	
T adhasive somestions		The bredent Order Number System fo	
T-adhesive connections	240 247	carbide tools	
Double t-adhesive mini connector	246-247	The Cut Types of the bredent Diatit ar	d tungsten carbide
Tap		tools	• • • • • • • • • • • • • • • • • • • •
Tool set	228-229, 231	Thermo-Pen	66, 294
Tap (drill) holder		Thermoplastic clasp	66
Tap holder, tool set	228	Thermoplastic resin	
Tap handwheel174, 198	, 223-224, 262, 387	Bio Dentaplast	358-359, 361, 383
Tap holder, tool set	228	Thermopress	
T-attachments		Acrylic Sep	
Wax T-attachments	190	Aluminium Catridges empty	
Technolit			
Technosil NT	200	bre.crystal	
Technosil NT duplicating silicone	250	bre.crystal HP	
		bre.dentan	
Technosil NT duplicating silicone	250	bre.dentan HP	
Teeth	000	bre.flex	.358-359, 362, 383
combo.lign Opaquer		bre.flex 2 <sup>nd</sup> Edition	362
crea.lign		Thermopress - Accessories	
crea.lign sets		Bio Dentaplast	.358-359, 361, 383
novo.lign A		Dentasil tooth protection silicone	
novo.lign P294		Thermopress Injection Moulding Syste	
visio.lign	288-305, 308-326	Flexible Acrylic sprue wax	
visio.link		Thermo-syringe	
visio.sil		Adhesive acrylic wax	
	,	rianiconto acigno waxiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiiii	

Thickener		Transverse force distributor
Multisil-Epithetik thickener		Vario-Soft 3 mini sv167
Thinner for ceramic separating liquid	272	Vario-Soft 3 sv162
Thinner for Oxyd-Stop PM		Vario-Soft 3 sv bridge-sectioning attachment 164-165
Thinner for spacer varnish		Vario-Soft attachment group160-167
Threaded bush, refractory alloy, cast-on		Tray material UV13, 337
Threaded rods	223-224	Trays
Titanium		breciform D impression tray
Threaded rods		brecision implant heavy54
Titanium adjusting screw		brecision implant light54
Titanium adjusting screw	174, 176	brecision impression materials8
Titanium coping		brecision putty soft56
SKY fast & fixed Prosthetic cap		Tungsten carbid burs
Titanium matrix housing for acrylic		Diatit and tungsten carbide tools423, 426-427439
Titanium matrix housing for metal		Tungsten carbide burs for processing of acrylics14
Titanium matrix housing vks-oc		Tungsten carbide burs for processing of plaster31
Titanium screws164, 228		Tungsten carbide burs, Diatit burs
Titanium screws, extended head	228	Tungsten carbide burs for processing of plaster31
mg		Tungsten carbide centering174, 223, 228
Titanium matrix housing vks-oc	144	Tungsten carbide centring drill224, 281
Tool set228	3-229, 231	Tungsten carbide drill20, 202, 206, 209, 344, 377
Tools for Milling technique	440	Master-Pin Diatit tungsten carbide step drill
Tooth completion material		20, 23, 377
crea.lign291, 299	9-300, 305	Tungsten carbide facing cutter228
Tooth protection silicone	·	Tungsten carbide tools for Milling technique
Dentasil tooth protection silicone	362	Tools for Milling technique440
Tooth-colored crown and bridge material		Tungsten carbide tools with relief419-422, 424-425
top.lign breformance liquid cold	68	Turbine
top.lign breformance liquid hot		Airaqua Turbine389
Tooth-colored die varnish		Turbo grinder
diephos dentine	33	Diagen-Turbo-Grinder286, 353, 457, 459-460
Tooth-colored resin		Tweezers
top-lign breformance	68-69	Ball Clip
Тор		Blue-Clip
Duro-Top	259	Repositioning tweezers
op.lign breformance liquid cold		Twist drill
op.lign breformance liquid hot		Diatit-Multidrill with stop228
op.lign professional		Diatit-Multidrill, twist drill
op-ligh breformance		198, 202, 206, 209, 223-224, 228, 362
Torque ratchet	00 00	twp, Perlform
Universal screwdriver set	216 408	Dipping wax87
Fransblock		Dipping wax
Transfer device		11
Transfer of investment model	200, 300	U
Master-Copy	252-253	UJ (upper jaw) wax sheets
Fransfer patrices	202-200	Stippled wax sheets
Metal transfer patrices vks43, 141	1 1//-15/	Ultraflex
Transfer set	1, 177-107	Ultrasonic cleaning solution
Impression transfer set	170	Dentaclean ultrasonic cleaning agent
	179	Undercut wax31
Transflex, electroplated diamond disc	4EO 4EE	uni.lign
Transflex		uni.lign speed
Transflex-T	453, 455	Unique Brush
Transflex-T, electropated diamond disc	450 455	Unit
Transflex-T		Milling unit BF 2262, 387-388
Transfuser	125, 169	Waxpool duo
Transparent blooking-out material	40	
Transblock	12	Unit for polymerizing, light-curing unit 2
Transparent silicone for keys	664	Polylux 2
visio.sil ILT		Unit for polymerizing, light-curing unit PL 20
Transverse fixation	220	Polylux380
Transverse force distributor		Universal 2
Interlock 2°		Universal paralleling mandrel 2172, 179
Vario-Soft 3 matrix housing		Universal grinding tool
Vario-Soft 3 mini	166	Diacryl grinder dcs457-458

Universal investment material	VC 2
Brevest M1255	Vario Compress 2
Brevest Rapid 1124, 256	Veneer shells
Universal paralleling mandrel	combo.lign Opaquer296
162, 164, 166-167, 194, 196	crea.lign291, 299-300
Universal paralleling mandrel 2172, 179	crea.lign sets305
Universal polishing paste	novo.lign A294, 304, 308-309
Pumice polishing paste345	novo.lign P294, 296, 304, 310-311
Universal screwdriver set216, 408	visio.lign
Universal screwdriver set for contra-angles217, 409	visio.link
UV, Litebloc UV31	visio.sil
UV-unit, light-curing units	visio.sil ILT
Polylux380	Veneering
Polylux 213-14	Diamond grinding tools
	Opaquer mixing liquid
V	Retention Crystals and Beads109
Vocassar activator allegace	Veneering Ceramic
Vacuum activator silencos366	K-Primer291, 296-297
Vacuum activator silencos kids367	Veneering CoCr alloys
Vacuum press system	Ceram-Bond270
BioHPP353-356, 381	Chrom-Kobalt-Bonding270
for 2 press351-356, 381	Veneers
Vacuum pressure casting	Retention adhesive109
Sprues for vacuum pressure casting112	visio.lign Toolkit291, 305
Sprues for vacuum pressure casting for voluminous	Vessels ppk241
castings112	Visio Dip
Vacuum unit	Dipping wax87
mixing cup16, 123, 376	visio.lign288-305, 308-326
mixing-stirring spiral16, 123, 376	visio.lign Toolkit
Vario Compress 1174-175	visio.link
Vario Compress 2176	visio.sil
Ceramic spacer176	
Friction silicone	visio.sil fix
Titanium adjusting screw174, 176	visio.sil ILT
	Visio-Dip
Vario-Soft 3160-161	Dipping wax Visio-Dip86
Vario-Soft 3 conicalbridge	vks
Vario-Soft 3 matrix housing163	Vario-Stud-Snap vks-oc 135-141, 144-146, 148-149
Vario-Soft 3 mini	Vario-Stud-Snap vks-oc exchangeable stud with
Vario-Soft 3 mini sv167	adhesive sleeve150
Vario-Soft 3 sv162	Vario-Stud-Snap vks-oc rs 2.2 Abutments 136-139
Vario-Soft 3 sv bridge-sectioning attachment164-165	Vario-Stud-Snap vks-sg151-152, 154-158
Vario-Soft 3 zircon sv168, 283	Vario-Stud-Snap vks-sg exchangeable stud with
Vario-Soft attachment group160-167	adhesive sleeve 150, 153
Vario-Soft-Bar vss188-189	vks-oc
Vario-Soft-Bar-Pattern vsp184-186	Blooking out discs141, 144
Vario-Stud-Snap vks-oc 135-146, 148-149	Inserting instrument141, 144
Vario-Stud-Snap vks-oc exchangeable stud with adhesive	Matrix inserting instrument 137, 146, 154, 156
sleeve150	Metal matrix housing vks-oc144
Vario-Stud-Snap vks-oc rs 2.2 Abutments136-139	Titanium matrix housing vks-oc144
Vario-Stud-Snap vks-sg 151-152, 154-158	vks-sg matrix housing
Vario-Stud-Snap vks-sg exchangeable stud with adhesive	VS
sleeve	
Vario-Stud-Snap vks-sg free end dentures	Vario-Soft attachment group160-167
Vario-Stud-Snap vks-sg 151-152, 154-158	VS 3
	Matrix inserting instrument
Vario-Stud-Snap vks-sg exchangeable stud with	Vario-Soft 3
adhesive sleeve150, 153	Vario-Soft 3 matrix housing163
Varnish	Vario-Soft 3 sv bridge-sectioning attachment 164-165
Die varnish, light-curing, opaque32	vs 3 mini
Light-curing die varnis h32, 43	Matrix inserting instrument166-167
Multisil sealing liquid340	Universal paralleling mandrel166-167, 194, 196
Spacer varnish, air-drying33	Vario-Soft 3 mini166
VC 1	vs 3 mini sv
Vario Compress 1174-175	Universal paralleling mandrel166-167, 194, 196
	Vario-Soft 3 mini sv

vs 3 sv		Wax blanks	
Universal paralleling mandrel	162	breCAM.wax	
Vario-Soft 3 sv	162	Wax bur 440, 442-443, 446-447, 450-4	52
vs 3 sv btg		Conical bur446-4	48
Universal paralleling mandrel	164	Parallel bur442-4	45
vsp		Wax casting bulbs3	53
Matrix inserting instrument	.184-185	Wax debubblizer	
Paralleling mandrel metal vsp	.184-185	Silicone and wax surface tension reducing agent	
Vario-Soft-Bar-Pattern vsp	.184-186	10, 1	20
VSS		Wax dipping unit	
Matrix inserting instrument	188	Waxpool duo	
Vario-Soft-Bar vss	.188-189	Wax for outer copings	91
		Wax hinges	
W		Wax bar hinges1	91
Mach muiman for flagge	404	Wax instruments	
Wash primer for fleece	121	Quick Change90, 3	95
Wax	000	Wax insulating (separating) liquid	
Adhesive acrylic wax		Isobre wax insulating (separating) liquid	
Auxiliary wax elements wht		Wax insulating liquid3	32
beauty setup	50	Wax knife	
Biotec blocking out wax		Ergonom wax knife329, 3	
Biotec cervical wax		Foam rubber grip lining	
Biotec metal-ceramic blocks without collar		Footswitch	
Biotec milling wax		Modelling tips3	
Biotec modelling wax grey/green		Rest for handpiece	
Biotec-Wax for outer copings		Waxpool duo88, 3	
Bite blocks		Wax knife bwm 3	
Cervical wax		Wax knife, electrically heated89, 3	86
Dipping wax		Wax matrix housing zg1	78
Dipping wax Elaflex	86	Wax neutralizing agent	
Dipping wax Visio-Dip		Wax-Lite surface tension reducing agent 10, 1	20
Functional margin wax		Wax palatal patterns3	36
K2 exact carving wax		Wax pattern	
Life-Color wax	50, 92	Biotec reels of wax pattern1	16
Modelling wax for Epithetiks		Protek wax pattern1	15
Modelling wax, wax sheets	51, 335	Protek wax pattern, lingual bar235-2	36
Modelling waxes6		Wax pattern sticks1	
Quadrosticks		Wax patterns cut to size, Protek2	
Set-up wax		Wax patterns cut to size, Protek2	40
Spacer wax		Wax rod attachments1	90
Splendido		Wax sheets	
Undercut wax	-	Modelling wax, wax sheets51, 3	35
Wax for outer copings		Wax sheets, stippled, Protek	
Wax palatal patterns		Stippled wax sheets2	38
Wax pattern sticks	116	Wax sprues	
Wax adapter	237, 407	Double rinsing heads1	15
Wax adhesive		Double rinsing heads for voluminous castings1	
Protek wax adhesive	240	Rinsing heads1	
SERACOLL UV59	, 94, 282	Rinsing heads for voluminous castings1	
Wax attachment elements	191	Sprues for centrifugal casting1	
Wax attachments		Sprues for centrifugal casting for voluminous casting	
Wax bar attachments	190	1	
Wax rod attachments	190	Sprues for vacuum pressure casting1	
Wax T-attachments		Sprues for vacuum pressure casting for voluminous	
Wax bands		castings1	12
Protek	239-240	Wax T-attachments	
Wax bands, Protek		Wax wire, Protek	50
Reels of wax patterns	116, 239	Reels of wax patterns116, 2	30
Wax bar attachments		Wax, Gecko modelling wax	.03
Wax bar hinges		GeckoGecko modelling wax	a၁
Wax bars142, 156, 190-191,		Wax, Adhesive resin wax	32
Wax blanks		Thermo-syringe29, 330, 3	ເດວ
breCAM.cutter	73. 81	Wax, Crown-bridge-inlay wax, KBI29, 330, 3	
breCAM.resin milling blanks		vvan, Glowii-biluge-ililay wax, NDI	52
2.22. IIII Con Timing Diamoni			

Wax, joiners, Protek
Protek clasp/bar joiners236
Wax, reels of patterns, Protek
Reels of wax patterns116, 239
Wax-Lite surface tension reducing agent10, 120
Waxpool duo88, 385
Wax-up
Waxpool duo385
Welding
LV 1 laser joint243
Wheel
Tap handwheel174, 198, 223-224, 262, 387
WHT
Between hollow blocks98
wk wax adhesive
Protek wax adhesive240
Work preparation
Master-Pin Diatit tungsten carbide step drill
20, 23, 377
Master-Pin Radix-K19
Master-Pin Radix-S18
Master-Sep19, 23
Master-Split base former(s)26
Master-Split model system25-27, 40
Model formers large/medium/small26
wstg
Wax bars142, 156, 190-191, 206, 209

X	
X-connector bonding agent	55, 58
X-connector bonding agent	55, 58
X-resin	
X-resin flow	
X-ray template	
SKY X-ray template	52
X-resin	
X-resin flow	
X-resin OPG	
X-tray template	
X-resin OPG	52
z	
zg	
Cylindrical attachment	
Cylindrical attachment Matrix inserting instrument	
	179
Matrix inserting instrument	179
Matrix inserting instrumentZi-polish	179 283, 487
Matrix inserting instrument	179 283, 487 284
Matrix inserting instrument	179 283, 487 284 283
Matrix inserting instrument	179 283, 487 284 283 284
Matrix inserting instrument Zi-polish Zirconium 5-motion-active 5-motion-colors 5-motion-glue	179 283, 487 284 283 284 284-285
Matrix inserting instrument Zi-polish Zirconium 5-motion-active 5-motion-colors 5-motion-glue 5-motion-scanner	179 283, 487 284 283 284 284-285
Matrix inserting instrument Zi-polish Zirconium 5-motion-active 5-motion-colors 5-motion-glue 5-motion-scanner 5-motions-milling cutter	
Matrix inserting instrument Zi-polish	