

breedent buyer



Issue 7 • Aug/Sept 2006



Diacryl Grinding Tools Kit

Diamond electroplated trimming burs for quick processing of acrylic dentures. There are 5 differently shaped burs all with coarse-grain sharp edged diamonds for effective denture adjustment from any angle. Burs need to be used at 10,000-15,000 rpm to remain effective.

(340 01070) RRP £154.95

with
**FREE Set-Up
Grinding
Tool**
RRP£35.96

Two diamond grinding tools in one:

- Coarse large area for quick acrylic tooth adjustment
- Fine small precisely shaped grinding tip for delicate occlusal adjustment



(340 01010)



Transblock 250g

Transblock is a transparent, re-usable soft material suitable for various applications:

- Quick, malleable blocking out material that can be easily removed
- Transparent, re-useable spacer for special trays so depth can be easily judged
- An ideal adapter for thermoplastic clasps while they are being heated
- An excellent base for quick tooth set up that will trim perfectly with a wax knife around the gingival margin and can be easily adapted to the model for putty duplication
- And many more...

(540 01149) RRP £29.95

with
**FREE Brush
Pen pk20**
RRP£6.77

Screw top reservoir with a brush tip suitable for accurate application of liquids.



(540 00720)



Kit of 7 Burs for non-precious alloy

This kit was designed entirely with non-precious alloy in mind. Burs have a special relief cut ensuring a long service life. Modified angles of the relief cut increase the abrasiveness and lead to a better surface quality resulting in a considerable time saving. The 7 different shapes cover all possible angles of non-precious alloy trimming.

(330 01170) RRP £129.95

**with
FREE Black
Chungking
Star Brushes
RRP£21.37**

Star shaped brushes cause a light beating effect ensuring that the bristles polish into the deepest fissures and eliminate roughness in the shortest possible time (15 pieces).



(520 00130)



Diagen-Turbo-Grinder Kit

Diamond grinding system with extraordinary grinding properties due to special 'diagen' diamond binding material. Perfect for porcelain and zirconium oxide trimming with no contamination!

(340 00200) RRP £119.95

**with
FREE
Microflex
Ø 12mm disc
RRP£17.06**

Microflex Ø 12mm disc with only 0.08mm thickness.



(340 00142)



Acrylic Polishing Set

This is the complete range for polishing acrylics! It includes a slightly abrasive paste, 3 polishing mops for a lathe and a pumice polishing paste. All you need to get your dentures sparkly.

(350 00840) RRP £34.95

**with
FREE
Leather Buff
Ø 60mm
RRP£10.71**

This Leather Buff Ø 60mm creates a perfect lustre which prevents bacteria and deposits being trapped.



(350 00990)