



Protective Chamber

A portable chamber that can be used for protection and/or extraction of dust and fumes. It can be supplied with or without the extraction nozzle. The extraction nozzle outlet fits an ordinary vacuum hose pipe. Chamber replacement glass is also available for purchase separately.

Without extraction RRP £165.85 with extraction RRP £176.35

receive **FREE** Plaster Tungsten Carbide bur

OFFER
when you buy
the Protective
Chamber



H194 KS70 RRP £20.42

OR



H263 SH60 RRP £20.00

Spend only 1 minute... & save time!!!



Dentaclean impression and denture disinfectant

This is the strongest impression and denture disinfectant on the market. It will kill all viruses, bacteria and fungi on an impression or denture in 1 minute. It is effective for disinfection of 30 impressions or dentures after which effectiveness is reduced. Disinfectant is available in 1 litre concentrate (1litre of concentrate will make 10 litres of solution ~ mixing ratio 1:9 liquid to water). Once concentrate is mixed with water it is effective for 3 days.

1 litre Impression and denture disinfectant 520 01006 RRP £62.79



Dentaclean plaster removing agent

Ready-to-use solution that removes plaster from all surfaces (dentures, orthodontic appliances, mixing bowls, etc.) Speed plaster removing agent will remove plaster in an ultrasonic cleaner within 1 minute. Plaster removing agents (normal and Speed) are available in 1 and 2.5 litre bottles and the mixing ratios are 1:1 liquid to water (creating 2 or 5 litre solutions). Both solutions will not affect plastic, resin or colour of the restoration.

1 litre Speed plaster removing agent 520 01010 RRP £22.00


2.5 litre Speed plaster removing agent 520 00994 RRP £49.56

OFFER
buy Impression
solution
(520 010 06
RRP £62.79)


receive **FREE** litre of giploser speed
(520 01010 worth £22.00)

breedent buyer

Issue 10 • May-July 2007



Come and see us at the Denturist Show at Coventry on
May 11-12th on stand C22



Breciform D

Breciform D impression trays reduce the most frequent sources of errors during impression taking and thus produce better results. Breciform D is suitable for all impression techniques (mono and two phase) eg alginate, hydro-alginate combination, polyether and silicones. Breciform D can be customized by the dentist using stops and kneading silicone and will frequently replace special trays.



Brecision Base Putty Shore Hardness 65

Brecision Putty Soft is a non sticky and kneadable base impression material. It excels by convenient and simple mixing and can be cut and adjusted easily.

OFFER
when you buy
Breciform D
(580 000 35
RRP £258.70)

receive **Brecision Putty FREE!**
(worth £71.12)



Thermopress: Brand new machine

Thermopress is the Injection Moulding System for:

- Metal-free dentures
- Long-term temporary crowns & bridges
- Unbreakable/flexible denture bases for partial dentures, splints and mouthguards
- Tooth-coloured clasps, attachments and telescopic work

Come & see the new Thermopress 400 at stand C22

Advantages of the system

- All 5 materials can be processed in 1 machine
- No need for compressed air
- 2 flasks can be processed at one time
- Injection process takes 1 minute
- There is no excess heat or steam produced

Course

Course 4/5 June 2007 ~ £650.00*

- 2 Day practical course at Bredent, Germany
- Absolutely essential and unmissable
- Tailor made for individual needs
- All 5 system components covered
- Course fee, flights, transfers, accommodation and meals included

Our competitors do Flexible dentures...




Flexiplast

- Nylon based unbreakable denture material
- Suitable for.. **flexible dentures, splints, mouthguards and attachment work**
- Pink 1, 2, 3, tooth shade B and translucent
- Flexible material

...so do we!

...however, Bredent dont just stop there!



Polyan

- Thermoplastic resin based on conventional acrylics but with NO residual monomer
- Suitable for **full & partial dentures, splints, stents, bite-raisers, etc**
- Shades: Crystal clear, Pink 1,2,3 and Pink stippled
- Not a flexible material



Dentalos

- Quick temporary C&B material
- Can be veneered with conventional composites
- Available in A, B & C shades
- Not a flexible material



Bio Dentaplast

- Tooth-coloured material for production of **clasps, tooth-coloured denture frameworks, telescopes and attachments**
- 5 shades available: A1, A2, A3, B2 and B3
- Flexible material



Bio XS

- High-melting, very stable thermoplastic resin in cream colour
- Suitable for C & B work, attachments and precision-fit denture frameworks
- Can be veneered with conventional C&B materials

thermopress 400

* Course fee will be refunded upon purchase of the system or material discounts may apply if used in conjunction with other systems