

bredent

Prosthetics



July/August 2005

Ergonom wax knife



This is the tool that combines various instruments into one hence saving time as instruments do not need to be changed so that faster and more efficient processing of the wax model is possible.

- **Ergonomic design of handle ~ suitable for right- and left-handed**
- **Wax knife is burn-free**

(310 00013) £24.50

OFFER
£19.50

Wax palatal patterns gf



These patterns allow the technician to restore natural function and aesthetics to the palate. They are simple and easy to adapt and allow for individuality within your prosthetic work. Wax palatal patterns are available in 3 sizes (A,B,C) and 2 thicknesses (0.5mm & 1.5mm).

Assortment 430 02180 (120 pieces) £84.75

OFFER
£74.95


Ropak UV



This is a light curing pink coloured acrylic opaquer designed for coating Chrome Cobalt denture frameworks. It is a mixture of powder and liquid and can be used with a heat curing method of acrylic processing. The aesthetics created are beyond compare.

(520 00165 – powder) £29.95

(520 00164 – liquid) £17.95



Ropak Kompaktopaker UV

A convenient pre-mixed pink coloured opaquer applied to Chrome Cobalt frameworks to mask out the darkness of the alloy. The pre-supplied brush allows for ease of application. The opaquer should be used with a bonder and is suitable for cold-cure acrylics.

(540 00133) £46.95

OFFER
£39.95



Kompaktopaker tooth-coloured UV

This is a tooth-coloured opaquer for Chrome Cobalt denture frameworks designed to be used in the area of acrylic teeth. It will mask out the darkness of the alloy and enhance the aesthetic appearance of the denture. There is an applicator brush in the bottle making it easier to apply. The opaquer should be used with a bonder and is suitable for cold-cure acrylics.

(540 00105) £24.95

OFFER
£19.95



Posi-boy

It provides the perfect 'third hand'. The instrument is designed to simplify processing of cold-curing acrylics. The solid metal base ensures firm stand and the correct position in the pressure pot. There is no tilting, no flowing out of acrylic or change of pre-shaped saddles. Posi-boy is also corrosion-resistant in any pressure pot so long service life is guaranteed.

(360 01010) £69.75

OFFER
£59.95



Plaster removing agent Speed

This is a ready-to-use solution that removes plaster from all surfaces. It will help with de-flasking dentures as the plaster is removed without harming the acrylic surface or affecting the shade. Hard plaster is also removed quickly from the mixing bowls without damaging the surface. This is an extremely effective solution. Will remove all plaster off a denture within 1 minute if used in the ultrasonic bath.

(520 01010) 1litre **£19.95** (520 00994) 2.5 litre **£44.95**



Diacryl grinding tools dsc

This is the only kit you will ever need for finishing dentures and special trays. All the burs are diamond-coated and will last for a very long time. 6 different grinding instruments will make any job simple and quick to finish.

Assortment (340 01070) £144.85

OFFER
£129.95