


SKY temp



The ideal basis for temporary restorations

The ideal basis for temporary restorations

Demands of patients on implant restorations are increasing continuously. Owing to professional and social obligations numerous patients are not willing to accept gaps in their dentition over longer periods. Temporary restorations provide an attractive alternative for such patients.

The SKY temp plastic abutments are also ideal for modern implantological concepts such as progressive loading or immediate loading.



SKY temp S
Order No. SKYTEMPS

SKY temp M
Order No. SKYTEMPM

SKY temp L
Order No. SKYTEMPL

Advantages:

- ▶ Perfect gingival reproduction with three prefabricated designs
- ▶ The extra-smooth surface avoids irritation of the mucosa
- ▶ Temporary restorations can be easily and quickly fabricated

SKY temp – as individual gingiva former

SKY temp is also the perfect basis for individual gingiva formers to fulfill high demands on red-white aesthetics.

A bonding agent is sufficient for direct veneering of the SKY temp abutments. It is not required to incorporate mechanical retentions.




Clinical case



After customizing the height of the abutment, the suitable facing can be selected.



The suitable crown and bridge material is used to veneer the palatal area of the facing.



Adequate fit of temporary crowns and SKY temp ensures the success of the restoration. Frequently, luting is not required thanks to the high level of friction that is achieved.



The patient will receive a highly esthetic restoration within just a few minutes.

