

compoForm UV

Light-curing composite for modelling, fixation of separated bridges and for quick fabrication of post and core restorations. compoForm UV burns out without leaving any residue and produces homogeneous casting results.

Individual modelling directly from the syringe.
Thanks to immediate hardening with a polymerization lamp the model can be built up in a safe and controlled manner.



The stability of compoForm UV renders the material perfectly suitable for the transfer of the jaw situation and, consequently, stress-free working is ensured.



1
Modelling and further processing of telescopic and conical crowns can be perfectly controlled by means of a visual check of the layer. The high stability of the hardened composite allows reworking with a bur.



2
Thanks to low shrinkage and burning without any residue, the composite is ideal for fixation of bridges to be soldered.



3
Model cast extensions can be fabricated in a simple and time-saving manner. compoForm UV can be easily removed from the plaster.



4
Undercuts on dies can be quickly and completely blocked out.



5
Burning without any residue and reduced swelling behavior provide perfect preconditions for top-quality casting results.

compoForm UV
2 x 3 ml syringes
10 application cannulas
Order No. 540 0115 0



Accessories

Application cannulas
25 pieces
Order No. 580 0001 8

bredent

GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany · Tel. (+49) 0 73 09 / 8 72-4 40 · Fax (+49) 0 73 09 / 8 72-4 44
www.bredent.com · e-mail info@bredent.com

03/07 197 GB