


Magic...



.... the magic brushes among dental brushes



bredent

GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany · Tel. (+49) 0 73 09 / 8 72-4 40 · Fax (+49) 0 73 09 / 8 72-4 44
www.bredent.com · e-mail info-lab@bredent.com

11/07 239 GB 3

MagicContrast - black hair

Fatigue-free working thanks to perfect contrast between ceramic material and brush hair.

The synthetic hair provides lasting springiness.


The „transformers“ among the brushes which obtain their original pointed shape after they have been washed out or by vibrating them lightly.



1 The pointed shape of dry brush hair is immediately restored by wetting, tapping off or vibrating the brush lightly.



2 Thanks to the strong contrast, any contamination such as dust or dry ceramic particles can be clearly recognized.



3 From the spatula shape to the original shape by rinsing the brush in liquid and then tapping it off or vibrating it lightly with a suitable instrument.



6 High elasticity to pick up specific quantities of ceramic material more easily.



7 The optimized springiness restores the shape of the brush tip immediately after picking up ceramic material. Additionally, the respective quantity can be easily assessed thanks to the contrast.



8 The spatula shape can be easily achieved after pressing the tip with two fingers. This way, the brush can be shaped individually.




9 Large quantities of ceramic material are picked up and time-saving layering is achieved thanks to the elasticity and the springiness of the brush hair.



10 The spatula shape does not reduce the stability of the brush hair; hence less time is required for applying the ceramic material to the framework.












MagicContrast

Scale 1:1	Product name	Size	Qty.	Order No.
	MagicContrast	4, 6, 8	1 piece each	390 CSET 1
	MagicContrast	1	2 pieces	390 C001 0
	MagicContrast	2	2 pieces	390 C002 0
	MagicContrast	4	2 pieces	390 C004 0
	MagicContrast	6	1 piece	390 C006 0
	MagicContrast	8	1 piece	390 C008 0
	MagicContrastBigBrush	8 BigBrush	1 piece	390 C008 B
	MagicContrast	1/0	2 pieces	390 CS01 0
	MagicContrast-Opaker	5	2 pieces	390 CS03 0

MagicBrush - golden brown hair

High springiness of durable synthetic hair for simpler layering of the ceramic material. MagicBrush and MagicContrast only differ by the color of the hair.

Restore the original pointed shape by tapping off or vibrating the brush lightly, just like with the MagicContrast brushes.

Scale 1:1	Product name	Size	Qty.	Order No.
	MagicBrush	4,6,8	1 piece each	390 MSET 1
	MagicPaintBrush	00 000	1 piece each	390 MS23 0
	MagicBrush	1	2 pieces	390 M001 0
	MagicBrush	2	2 pieces	390 M002 0
	MagicBrush	4	2 pieces	390 M004 0
	MagicBrush	6	1 piece	390 M006 0
	MagicBrush	8	1 piece	390 M008 0
	MagicBigBrush	8 BigBrush	1 piece	390 M008 B
	MagicBrush	1/0	2 piece	390 MS01 0
	MagicBrush	2/0	2 pieces	390 MS02 0
	MagicBrush-Opaker	5	2 pieces	390 MS03 0

KoliBrush - golden brown natural hair

Natural hair brushes made of superior quality Kolinsky hair.



Improved design of the tip of the BigBrush is achieved thanks to the integrated spheres - for simpler modelling. Fine, stable tip thanks to carefully selected hair.












1 The shape and quality of the hair for perfect retention of moisture and improved adhesion and application of ceramic material.



2 The desired elasticity is obtained by the unique design and combination. This way already applied ceramic layers will not be damaged.

KoliBrush

Scale 1:1	Product name	Size	Qty.	Order No.
	KoliBrush	4, 6, 8 B	1 piece each	390 KSET 1
	KoliBrush	1	2 pieces	390 K001 0
	KoliBrush	2	2 pieces	390 K002 0
	KoliBrush	4	2 pieces	390 K004 0
	KoliBrush	6	1 piece	390 K006 0
	KoliBrush	8	1 piece	390 K008 0
	KoliBigBrush	8 BigBrush	1 piece	390 K008 B
	KoliBrush	1/0	2 pieces	390 KS01 0
	KoliOpakerBrush	5	2 pieces	390 KS03 0

