

Waxpool duo

Two units in one – advanced technology separately controls two dipping units and a wax knife

Two separate lids ensure reduced contamination

Digitally controlled temperature regulator from 0 – 120° C for each dipping unit

Speed-switch for quick heating of the wax knife to the highest temperature (250° C)

Connection for the wax knife

Lowered safety dipping basin


Waxpool duo
Order No. 110 0151 0

Quick melting of the wax is possible through new electronic controller

Stable, easy to clean plastic housing

Dipping wax with different characteristics for different demands in the bredent DT-Catalogue or at www.bredent.com

Accessories



Foam rubber grip lining
4 pieces
Order No. 140 0096 4



Rest bwm 3
Order No. 140 0096 5

bredent

GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany · Tel. (+49) 0 73 09 / 8 72-4 40 · Fax (+49) 0 73 09 / 8 72-4 44
www.bredent.com · e-mail info-lab@bredent.com

03/07 223 GB